

Contents of Appendix J

TDCS Excavation Spoils – PCB Testing

Text Summary

Table TDCS Spoils

TDCS EXCAVATION SPOILS

The excavated shot rock from the TDCS Shaft and Tunnels was removed from the TDCS after each round of blasting and placed in piles in designated bins on the spoils structure for PCB sampling. The shot rock spoil from every blast was sampled and analyzed for PCB. Five samples were taken at locations from each spoil pile. The results of spoils analyses determined the classification of the spoil materials and ultimate disposition location. Refer to TDCS Spoils Table for a summary of blasts and PCB spoil test results. Spoils with PCB concentrations less than 10 ppm were allowed to remain onsite in the Onsite Spoil Placement Area in the back lot. Spoil materials with PCB concentrations greater than 10 ppm were required to be properly disposed of offsite. Based on the results of the spoils analyses, PCBs were first encountered in the shaft after blast #12 (elevation 120 to 110 ft NGVD). The spoil PCB concentrations were greater than 10 ppm from shaft blast #12 to tunnel blast #82 (tunnel 2 sta 2:1+37 and tunnel 3 sta 3:1+54). The Spoils PCB concentrations were below 10 ppm for the remainder of the excavation from tunnel blast #83 (tunnel 2 sta 2: 1+37 and tunnel 3 sta 3:1+54) to the end of both tunnels 2 and 3 and the final tunnel blast #117. The two blasts for the TDCS Sump at the bottom of the shaft were both well above 10 ppm.

PCB Concentrations

		Shaft Blast 1	Shaft Blast 6.3	Shaft Blast 8.2	Shaft Blast 9.2	Shaft Blast 11.1	Shaft Blast 12.1	Shaft Blast 12.2	Shaft Blast 13.1	Shaft Blast 14.1	Shaft Blast 14.2	Shaft Blast 15.1	Shaft Blast 15.2	Shaft Blast 16.1
Parameter	Units	9/20/07	11/2/07	11/14/07	11/27/07	12/12/07	12/18/07	12/27/07	1/7/08	1/9/08	1/10/08	1/16/08	1/17/08	1/22/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	0.336 AB	0.0564 AB	ND							
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	0.1260 AD	ND	ND	0.165 AD	ND	42.9 AD	214.0	65.9 AD	26.4 AD	4.37 AD	55.7 AD	42.2 AD	2.52 AD
Aroclor-1248	μg/g	ND	ND	0.116 AE	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	0.0407 JAF	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	0.0117 AG,J	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB														
Amount	μg/g	0.1667	0.0117	0.116	0.501	0.0564	42.9	214.0	65.9	26.4	4.37	55.7	42.2	2.52

											Bellout	Bellout	Bellout		
		Shaft Blast	Shaft Blast	Shaft Blast	Shaft Blast	Shaft Blast	Shaft	Shaft Blast	Shaft Blast	Shaft Blast	Blast 1.1 &	Blast 1.2 &	Blast 1.3 &	Bellout	Bellout
		16.2	16.3	17.1	18.1	19.1	Blast 19.2	19.3	20.1	20.2	2.1	2.2	2.3	Blast 3.1	Blast 3.2
Parameter	Units	1/25/08	1/25/08	1/29/08	2/5/08	2/11/08	2/18/08	2/20/08	2/22/08	2/26/08	3/5/08	3/5/08	3/7/08	3/17/08	3/19/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	72.1 AD	42.7 AD	20.5 AD	58.8 AD	3.82 AD	13.9 AD	14.4 AD	19.6 AD	23.2 AD	4.16 AD	8.25 AD	70.8 AD	74.8 AD	24.2 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB															
Amount	μg/g	72.1	42.7	20.5	58.8	3.82	13.9	14.4	19.6	23.2	4.16	8.25	70.8	74.8	24.2

Notes

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.

1

PCB Concentrations

		Bellout Blast 3.3	Tunnel Test Blast 1	Tunnel Test Blast 2	Tunnel Test Blast 3	Tunnel Test Blast 4.1	Tunnel Test Blast 5	Tunnel Test Blast 6.1	Tunnel Test Blast 6.2	Tunnel Test Blast 7.1	Tunnel Test Blast 8.1	Tunnel Test Blast 9.1	Tunnel Test Blast 10.1	Tunnel Test Blast 11.1
Parameter	Units	3/20/08	3/25/08	3/31/08	4/2/08	4/4/08	4/11/08	4/15/08	4/15/08	4/17/08	4/18/08	4/21/08	4/22/08	4/23/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	17.7 AD	59.6 AD	32.7 AD	4.4 AD	16.3 AD	17.0 AD	35.4 AD	3.43 AD	21.0 AD	25.4 AD	16.3 AD	49.5 AD	59.1 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB														
Amount	μg/g	17.7	59.6	32.7	4.4	16.3	17.0	35.4	3.4	21.0	25.4	16.3	49.5	59.1

		Tunnel Test Blast 12.1	Tunnel Test Blast 12.1 Duplicate	Tunnel Test Blast 13.1	Tunnel Test Blast 14.1	Tunnel Test Blast 15.1	Tunnel Test Blast 16.1	Tunnel Test Blast 16.1 Duplicate	Tunnel Test Blast 17.1	Tunnel Test Blast 18.1	Tunnel Test Blast 19.1	Tunnel Test Blast 19.1 Duplicate	Tunnel Test Blast 20.1	Tunnel Test Blast 21.1
Parameter	Units	4/24/08	4/24/08	4/25/08	4/29/08	4/29/08	5/1/08	5/2/08	5/2/08	5/5/08	5/6/08	5/6/08	5/7/08	5/7/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	26.3 AD	30.7 AD	24.0 AD	36.8 AD	46.9 AD	24.5 AD	25.6 AD	23.5 AD	27.9 AD	24.1 AD	37.6 AD	33.9 AD	42.0 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB														
Amount	μg/g	26.3	30.7	24.0	36.8	46.9	24.5	25.6	23.5	27.9	24.1	37.6	33.9	42.0

<u>Notes</u>

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.

		Tunnel Test Blast 22.1		Tunnel Test Blast 23.1 Duplicate	Tunnel Test Blast 24.1	Tunnel Blast 25.1	Tunnel Blast 26.1	Tunnel Blast 26.1 Duplicate	Tunnel Blast 27.1 & 28.1	Tunnel Blast 29.1	Tunnel Blast 30.1	Tunnel Blast 30.1 Duplicate	Tunnel Blast 31.1	Tunnel Blast 32.1
Parameter	Units	5/9/08	5/12/08	5/12/08	5/14/08	5/16/08	5/19/08	5/19/08	5/22/08	5/27/08	5/28/08	5/28/08	5/29/08	5/30/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	70.1 AD	33.6 AD	32.8 AD	38.5 AD	33.5 AD	37.9 AD	28.0 AD	27.2 AD	10.0 AD	29.8 AD	26.9 AD	20.7 AD	7.81 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB Amount	μg/g	70.1	33.6	32.8	38.5	33.5	37.9	28.0	27.2	10.0	29.8	26.9	20.7	7.8

					Tunnel						Tunnel			
		Tunnel	Tunnel	Tunnel	Blast 35.1	Tunnel	Tunnel	Tunnel	Tunnel	Tunnel	Blast 40.1	Tunnel	Tunnel	Tunnel
		Blast 33.1	Blast 34.1	Blast 35.1	Duplicate	Blast 36.1	Blast 37.1	Blast 38.1	Blast 39.1	Blast 40.1	Duplicate	Blast 41.1	Blast 42.1	Blast 43.1
Parameter	Units	5/30/08	6/3/08	6/3/08	6/3/08	6/4/08	6/5/08	6/6/08	6/9/08	6/9/08	6/9/08	6/10/08	6/11/08	6/11/08
Aroclor-1016	μg/g	ND												
Aroclor-1221	μg/g	ND												
Aroclor-1232	μg/g	ND												
Aroclor-1242	μg/g	14.0 AD	25.1 AD	12.0 AD	14.6 AD	81.0 AD	76.6 AD	20.7 AD	16.6 AD	18.1 AD	24.4 AD	103 AD	6.80 AD	19.9 AD
Aroclor-1248	μg/g	ND												
Aroclor-1254	μg/g	ND												
Aroclor-1260	μg/g	ND												
Total PCB														
Amount	μg/g	14.0	25.1	12.0	14.6	81.0	76.6	20.7	16.6	18.1	24.4	103.0	6.8	19.9

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.

		Tunnel Blast 44.1	Tunnel Blast 45.1	Tunnel Blast 46.1	Tunnel Blast 46.1 Duplicate	Tunnel Blast 47.1	Tunnel Blast 48.1	Tunnel Blast 49.1	Tunnel Blast 50.1	Tunnel Blast 51.1	Tunnel Blast 52.1	Tunnel Blast 53.1	Tunnel Blast 54.1	Tunnel Blast 54.1 Duplicate
Parameter	Units	6/12/08	6/13/08	6/16/08	6/16/08	6/16/08	6/17/08	6/18/08	6/19/08	6/20/08	6/23/08	6/24/08	6/25/08	6/25/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	21.3 AD	31.7 AD	39.9 AD	25.2 AD	80.5 AD	10.7 AD	12.3 AD	32.6 AD	45.3 AD	14.4 AD	23.1 AD	40.9 AD	33.3 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB Amount	μg/g	21.3	31.7	39.9	25.2	80.5	10.7	12.3	32.6	45.3	14.4	23.1	40.9	33.3

				Tunnel				Tunnel						Tunnel
		Tunnel	Tunnel	Blast 56.1	Tunnel	Tunnel	Tunnel	Blast 59.1	Tunnel	Tunnel	Tunnel	Tunnel	Tunnel	Blast 64.1
		Blast 55.1	Blast 56.1	Duplicate	Blast 57.1	Blast 58.1	Blast 59.1	Duplicate	Blast 60.1	Blast 61.1	Blast 62.1	Blast 63.1	Blast 64.1	Duplicate
Parameter	Units	6/27/08	7/1/08	7/1/08	7/2/08	7/7/08	7/8/08	7/8/08	7/10/08	7/11/08	7/14/08	7/16/08	7/18/08	7/18/08
Aroclor-1016	μg/g	ND												
Aroclor-1221	μg/g	ND												
Aroclor-1232	μg/g	ND												
Aroclor-1242	μg/g	27.5 AD	23.3 AD	21.8 AD	55.0 AD	12.2 AD	39.2 AD	30.6 AD	22.3 AD	11.7 AD	10.6 AD	19.8 AD	16.0 AD	14.4 AD
Aroclor-1248	μg/g	ND												
Aroclor-1254	μg/g	ND												
Aroclor-1260	μg/g	ND												
Total PCB														
Amount	μg/g	27.5	23.3	21.8	55.0	12.2	30.6	30.6	22.3	11.7	10.6	19.8	16.0	14.4

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.

PCB Concentrations

		Tunnel Blast 65.1	Tunnel Blast 66.1	Tunnel Blast 66.1 Duplicate	Tunnel Blast 67.1 (see note 9)	Tunnel Blast 68.1 (see note 10)	Tunnel Blast 69.1 (see note 11)	Tunnel Blast 70.1 (see note 12)	Tunnel Blast 70.1 Duplicate (see note 12)	Tunnel Blast 72.1	Tunnel Blast 73	Tunnel Blast 74.1	Tunnel Blast 75.1
Parameter	Units	7/21/08	7/24/08	7/24/08	7/25/08	7/25/08	7/28/08	7/29/08	7/29/08	7/30/08	7/31/08	8/1/08	8/4/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	23.2 AD	13.7 AD	14.6 AD	22.5 AD	6.58 AD	7.61 AD	17.4 AD	20.6 AD	0.86 AD	3.79 AD	14.8 AD	2.71 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB Amount	μg/g	23.20	13.70	14.60	22.50	6.58	7.61	17.40	20.60	0.86	3.79	14.80	2.71

		Tunnel										Tunnel	
		Blast 75.1	Tunnel	Blast 84.1	Tunnel								
		DUP	Blast 76.1	Blast 77.1	Blast 78.1	Blast 79.1	Blast 80.1	Blast 81.1	Blast 82.1	Blast 83.1	Blast 84.1	Duplicate	Blast 85.1
Parameter	Units	8/4/08	8/5/08	8/6/08	8/6/08	8/7/08	8/8/08	8/11/08	8/12/08	8/13/08	8/14/08	8/14/08	8/15/08
Aroclor-1016	μg/g	ND											
Aroclor-1221	μg/g	ND											
Aroclor-1232	μg/g	ND											
Aroclor-1242	μg/g	3.32 AD	2.28 AD	19.8 AD	10.8 AD	18.5 AD	14.7 AD	14.3 AD	21.1 AD	8.92 AD	6.07 AD	2.04 AD	1.93 AD
Aroclor-1248	μg/g	ND											
Aroclor-1254	μg/g	ND											
Aroclor-1260	μg/g	ND											
Total PCB	μg/g	3.32	2.28	19.8	10.8	18.5	14.7	14.3	21.1	8.92	6.07	2.04	1.93

Notes

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 9. Sample ID Tunnel Blast 67.1 contained rock spoils from Tunnel Blasts 66, 67 and 68.
- 10. Sample ID Tunnel Blast 68.1 contained rock spoils from Tunnel Blast 69 only.
- 11. Sample ID Tunnel Blast 69.1 contained rock spoils from Tunnel Blasts 66 and 70 only.
- 12. Sample ID Tunnel Blast 70.1 contained rock spoils from Tunnel Blast 71 only. There is no sample named Tunnel Blast 71.1.

		Tunnel Blast 86.1	Tunnel Blast 87.1	Tunnel Blast 88.1	Tunnel Blast 89.1	Tunnel Blast 89.1 Duplicate	Tunnel Blast 90.1	Tunnel Blast 91.1	Tunnel Blast 93.1	Tunnel Blast 94.1	Tunnel Blast 95.1	Tunnel Blast 95.1 Duplicate	Tunnel Blast 96.1
Parameter	Units	8/18/08	8/18/08	8/20/08	8/21/08	8/21/08	8/21/08	8/22/08	8/26/08	8/26/08	8/27/08	8/27/08	8/28/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	0.28 AD	3.08 AD	0.511 AD	0.372 AD	0.706 AD	0.755 AD	0.101 AD	0.382 AD	0.329 AD	0.344 AD	0.343 AD	1.01 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB													
Amount	μg/g	0.28	3.08	0.511	0.372	0.706	0.755	0.101	0.382	0.329	0.344	0.343	1.010

								Tunnel			Tunnel		
		Tunnel	Tunnel	Tunnel 2	Tunnel	Tunnel	Tunnel	Blast 101.1	Tunnel	Tunnel	Blast 103.1	Tunnel	Tunnel
		Blast 97.1	Blast 98.1	Blast 99	Blast 99.1	Blast 100.1	Blast 101.1	Duplicate	Blast 102.1	Blast 103.1	Duplicate	Blast 104.1	Blast 105.1
Parameter	Units	8/29/08	8/29/08	9/2/08	9/3/08	9/4/08	9/4/08	9/4/08	9/8/08	9/9/08	9/9/08	9/10/08	9/11/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	1.18 AD	0 .519 AD	ND	0.594 AD	2.75 AD	0.663 AD	0.553 AD	1.91 AD	0.768 AD	0.520 AD	0.777 AD	0.434 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB	μg/g	1.18	0.519	ND	0.594	2.75	0.663	0.553	1.91	0.768	0.520	0.777	0.434

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 9. Sample ID Tunnel Blast 67.1 contained rock spoils from Tunnel Blasts 66, 67 and 68.
- 10. Sample ID Tunnel Blast 68.1 contained rock spoils from Tunnel Blast 69 only.
- 11. Sample ID Tunnel Blast 69.1 contained rock spoils from Tunnel Blasts 66 and 70 only.
- 12. Sample ID Tunnel Blast 70.1 contained rock spoils from Tunnel Blast 71 only. There is no sample named Tunnel Blast 71.1.

		Tunnel Blast 106.1	Tunnel Blast 107.1	Tunnel Blast 107.1 Duplicate	Tunnel Blast 108.1	Tunnel Blast 108.2	Tunnel Blast 109.1 & 110.1	Tunnel Blast 110.2	Tunnel Blast 111.1	Tunnel Blast 111.2	Tunnel Blast 112.1	Tunnel Blast 113.1	Tunnel Blast 114.1
Parameter	Units	9/12/08	9/15/08	9/15/08	9/16/08	9/16/08	9/18/2008	9/18/2008	9/22/2008	9/22/2008	9/22/2008	9/23/2008	9/24/2008
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	0.285 AD	0.687 AD	0.802 AD	0.224 AD	1.59 AD	1.73 AD	0.459 AD	0.398 AD	0.526 AD	1.67 AD	0.273 AD	0.907 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Total PCB													
Amount	μg/g	0.285	0.687	0.802	0.224	1.59	1.73	0.459	0.398	0.526	1.67	0.273	0.907

		Blast 114.1	Tunnel	Tunnel	Tunnel	Sump	Sump Blast	Sump Blast
		Duplicate	Blast 115.1	Blast 116.1	Blast 117.1	Rock/Concrete	1	2
Parameter	Units	9/24/08	9/24/08	9/26/08	9/29/08	12/3/08	12/5/08	12/11/08
Aroclor-1016	μg/g	ND	ND	ND	ND	ND	ND	ND
Aroclor-1221	μg/g	ND	ND	ND	ND	ND	ND	ND
Aroclor-1232	μg/g	ND	ND	ND	ND	ND	ND	ND
Aroclor-1242	μg/g	1.35 AD	0.263 AD	0.466 AD	2.12 AD	300 AD	108	43.7 AD
Aroclor-1248	μg/g	ND	ND	ND	ND	ND	ND	ND
Aroclor-1254	μg/g	ND	ND	ND	ND	ND	ND	ND
Aroclor-1260	μg/g	ND	ND	ND	ND	ND	ND	ND
Total PCB								
Amount	μg/g	1.35	0.263	0.466	2.12	300	108	43.7

- 1. PCB analyses by EPA method 8082
- 2. ND = Non Detect
- 3. J- Indicates an estimated value. The analyte was detected in the sample at a concentration greater than the Method Detection Limit (MDL) but less than the Practical Quantitation Limit (PQL).
- 4. AD-Aroclor 1242 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 5. AF-Aroclor 1254 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 6. AG-Aroclor 1260 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 7. AE-Aroclor 1248 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 8. AB-Aroclor 1221 is being reported as the best Aroclor match. The sample exhibits an altered PCB pattern.
- 9. Sample ID Tunnel Blast 67.1 contained rock spoils from Tunnel Blasts 66, 67 and 68.
- 10. Sample ID Tunnel Blast 68.1 contained rock spoils from Tunnel Blast 69 only.
- 11. Sample ID Tunnel Blast 69.1 contained rock spoils from Tunnel Blasts 66 and 70 only.
- 12. Sample ID Tunnel Blast 70.1 contained rock spoils from Tunnel Blast 71 only, There is no sample named Tunnel Blast 71.1. \Geoschs-sherman\geoschvol1\common\DOCS\GE\Hudson Falls\2011\Phase 3 Complétion and Hydraulic MonitorngReport\Appendix\Appendix G TDCS Spoils Sample Results Summary\TableSpoils Results