

SCOPE OF WORK - SEPTEMBER 1993 SUPPLEMENTAL WORK PLAN

**Niagara Cold Drawn Corporation
Buffalo, New York**

NOVEMBER 1993

REF. NO. 3829 (3)

This report is printed on recycled paper.

CONESTOGA-ROVERS & ASSOCIATES

TABLE OF CONTENTS

	<u>Page</u>
1.0 INTRODUCTION	1
2.0 PROPOSED ACTIVITIES.....	3
2.1 PCB SAMPLING.....	3
2.2 ABANDONED SULFURIC ACID STORAGE TANKS.....	4
2.3 ANALYTICAL SERVICES PROTOCOL	4
3.0 SUMMARY.....	6

1.0 INTRODUCTION

On May 28, 1993, Mr. John T. Kolaga of Whiteman, Osterman & Hanna, representing the Niagara Cold Drawn Corporation, sent a petition to delist to Mr. Thomas C. Jorling of the New York State Department of Environmental Conservation (NYSDEC) requesting the removal of the Niagara Cold Drawn Site (No. 915046A) (the Site) from the NYSDEC's Registry of Inactive Hazardous Waste Sites. The petition to delist was accompanied by a report entitled "January 1993 Data Collection Program" which summarized historical investigative data and presented data collected in January 1993 to assess those areas of the Site which had previously been identified as areas of possible concern. The January 1993 Data Collection Program report concluded that no further investigatory or remedial activities were required at the Site prior to delisting.

On September 13, 1993, Mr. Michael J. O'Toole, Jr. of the NYSDEC sent a response to Mr. Kolaga indicating that the NYSDEC had denied Niagara Cold Drawn's petition to delist. The letter indicated that there was insufficient information to properly assess the petition to delist and that several specific areas of concern needed to be addressed in order to complete the assessment of the petition. The specific areas of concern identified by the NYSDEC were as follows:

1. The investigation did not sufficiently investigate the potential for PCB exposure. No subsurface soil samples were collected and analyzed for PCBs in the oil stained area or near the north fence where PCBs were detected in the past.
2. Discrete surface and subsurface soil samples need to be collected near the abandoned sulfuric acid storage tank. Composite sampling covering a large area is not the preferred method when specific areas of contamination or potentially contaminated areas are known.
3. It appears that the proper extraction procedure was not followed for the Target Compound List (TCL) analytes. Therefore, the analytical data is

questionable. The TCL parameters should have been analyzed according to the NYS Analytical Services Protocol (ASP).

In response to the above concerns identified by the NYSDEC, Niagara Cold Drawn Corporation has prepared this draft work plan (the September 1993 Supplemental Work Plan) to identify the efforts Niagara Cold Drawn Corporation will make to supply additional information to the NYSDEC to address the specific areas of concern set forth in the September 13, 1993 NYSDEC letter. This work plan presents the proposed scope of work for a supplemental investigation to be conducted immediately upon approval by the NYSDEC in order to attain the delisting of the Site.


2.0 PROPOSED ACTIVITIES

2.1 PCB SAMPLING

In 1984, NUS Corporation, on behalf of the United States Environmental Protection Agency (USEPA), conducted some testing at the Site. No hazardous waste was identified at the Site, with the exception of one soil sample collected near the north fence which was reported to contain PCB (Arochlor 1254 at 62 ppm and Arochlor 1260 at 123 ppm) along with several metals and organics at low levels. There was also an oil-stained area observed at the Site which was not sampled by NUS.

No evidence of PCB presence was identified by Conestoga-Rovers & Associates (CRA) at the Site during prior sampling events. Specifically, in January 1991, CRA collected three discrete surface soil samples (upper six inches) at three locations along the north fence line and one composite surface soil sample (upper three inches at three locations) within the oil-stained area. In January 1993, CRA collected composite surface soil samples (upper six inches within 10-foot by 5-foot area) at five locations along the north fence line with no additional sampling within the oil-stained area. PCBs were not detected in eight of these nine samples, with the remaining result being only 0.03 µg/g above the detection limit.

Because no evidence of PCB presence was identified in the surface soils during previous investigations, CRA did not perform subsequent subsurface soil sampling. Nonetheless, in response to the NYSDEC concerns regarding PCB presence in the subsurface soils, Niagara Cold Drawn Corporation is prepared to collect two discrete subsurface soil samples within each of the areas of concern at NYSDEC-approved locations adjacent to the north fence line and within the oil-stained area, as shown on Figure 2.1. Each sample will be collected by driving a split-spoon sampler twelve inches into the ground and subsequently collecting a sample for analysis from the bottom six inches of each retrieved soil sample. All sampling techniques and cleaning procedures will be identical to those used during the January 1993 Data Collection Program. Each of the four subsurface samples collected will be submitted for analysis of PCBs. In addition, one of



LEGEND

- SURFACE SOIL SAMPLE COLLECTED IN MARCH 1991 FOR PCB ANALYSIS
- COMPOSITE SOIL SAMPLE COLLECTED IN JANUARY 1993 FOR PCB ANALYSIS
- ▲ PROPOSED SUBSURFACE SOIL SAMPLE FOR PCB ANALYSIS (ONE SAMPLE IN OIL-STAINED AREA FOR TCL/TAL)
- ◼ PROPOSED SUBSURFACE SOIL SAMPLE FOR METALS/PCBs/pH/SULPHATE ANALYSIS

figure 2.1
SAMPLE LOCATIONS
NIAGARA COLD DRAWN CORP.
Buffalo, New York

CRA

the two samples collected within the oil-stained area will also be analyzed for the full Target Compound List (TCL)/Target Analyte List (TAL) parameters. The analytical results will be compiled, validated and submitted to the NYSDEC.

2.2 ABANDONED SULFURIC ACID STORAGE TANKS

Another identified area of concern is an area along the west side of the main building where sulfuric acid storage tanks were temporarily staged following their removal from the building. These tanks are no longer present in this area and were never used as sulfuric acid storage tanks in this area outside the building.

Niagara Cold Drawn Corporation will conduct an investigation of the soils in this area and collect one soil sample at a location approved by the NYSDEC. The sample will be collected by driving a split-spoon sampler twelve inches into the ground at the location shown on Figure 2.1 and subsequently collecting a sample for analysis from the bottom six inches of the retrieved soil sample. If the material encountered in the 6 to 12-inch layer is crushed stone fill, the sampler will be driven deeper to the uppermost layer of soil material which may be analyzed. In addition, a subsurface soil sample will be collected in a general use area of the Site to determine background conditions, the location of which to be field-selected with NYSDEC approval. All sampling techniques and cleaning procedures will be identical to those previously used by CRA and each of the two subsurface samples will be submitted for analysis of TAL Metals, EP Tox Metals, PCBs, pH and sulfate. The analytical results will be compiled, validated and submitted to the NYSDEC.

2.3 ANALYTICAL SERVICES PROTOCOL

The NYSDEC has questioned the validity of the previous analytical data for TAL parameters, as the analysis was not performed in accordance with the methods specified by NYS in its Analytical Services

Protocol (ASP). However, the extraction and analysis of the previously collected samples was performed using Method 3050 contained in SW846, which is the specific method referenced in ASP. Therefore, the analytical data collected are valid, notwithstanding the absence of certain deliverables required by ASP. In order to satisfy NYSDEC's concerns, however, CRA has contacted the analytical facility and acquired all available deliverables (attached hereto as Appendix A).

The analysis of sampling events proposed in this work plan will be performed in accordance with ASP and all deliverables will be made available to the NYSDEC in a timely fashion

During this supplemental sampling event, NYSDEC may request split samples collected for its own analysis. Niagara Cold Drawn Corporation will supply the NYSDEC with the required glassware for any samples, upon notification of its intent to analyse such split samples.

3.0 SUMMARY

Niagara Cold Drawn Corporation is prepared to conduct additional sampling in the areas of potential concern identified by the NYSDEC in its September 13, 1993 letter. These samples will be analyzed in accordance with ASP. Following completion of these activities, and assuming all the results are below action levels, Niagara Cold Drawn Corporation will resubmit its petition to delist for approval by the NYSDEC.

APPENDICES

APPENDIX A

PROJECT

DIGESTIONS

NOTEBOOK NO. 12-223

Continued From Page

SAMPLE NUMBER	DATE TIME	PREP BY	TYPE DIGEST	SPIKE SIK USED	INIT SAMPLE	FINAL VOLUME
SAND BLANK	1-26-93	mm	3050F		1.00	200
LCS-5	10A				1.00	
0029-01					1.02	
02					1.02	
03					1.03	
04					1.00	
05					1.04	
06					1.01	
07					1.00	
08					1.00	
09					1.04	
10					1.05	
11					1.01	
12					1.03	
13					1.00	
05-1/2 spk					1.00	
05-1/2 spk	✓	✓	✓	2 ml Furnace SIK	1.00	✓
0029-05	1-26-93	mm	3050P		1.05	200
05	10A				1.00	
05	✓			2 ml SPIKE-1	1.02	
SAND BLANK	1-25-93	mm	3050P		1.00	200
LCS-5	1-26-93 mm				1.00	
0029-05					1.00	
01-1/2 spk					1.00	
02					1.01	
03					1.01	
04					1.00	
05					1.00	
01-1/2 spk					1.03	
01-1/2 spk				2 ml SPIKE-1	1.01	

Continued on Page

Read and Understood By

Signed

Date

Signed

Date


PROJECT MERCURY 1-27-93

Continued From Page _____

SAMPLE NUMBER	Read	Conc	LS	UNITS	AMT SAMPLE	FINAL CONC DL	FINAL CONC	
							CORR	.99937
							SLOPE	42.000
							INTCPT	1.315
BLANK	2							
0.2	10							
0.4	18							
0.8	36							
1.0	40							
2.0	87							
3.0	127							
ICV	49	1.135		ug/L	100ml	0.0002mg/L	0.00114	113.5%
ICB	1	-0.008				0.0002mg/L	0.0002	
LCS-W 10TV	30	0.683				0.0002mg/L	0.00048	68.3%
LCS-S 13.7TV	56	1.322		ug/g	0.0002g	0.100mg/kg	10.85mg/kg	79.2%
CCV	46	1.064		ug/L	100ml	0.0002mg/L	0.00106	106.4%
CCB	2	0.016					0.0002	
ICV-2 15TV	57	1.326					0.00133	88.7%
S th 0103-01	2	0.016		ug/g	0.21g	0.100mg/kg	0.0100	
0103-01 R _{sp}	4	0.064			0.23g		0.0120	
S _{pk}	34	0.778			0.21g		0.37mg/kg	77.8%
0104-01	10	0.207			0.21g		0.0100	
R _{sp}	9	0.159			0.20g		0.0100	
S _{pk}	49	1.135			0.21g		0.54mg/kg	113.5%
02	4	0.064			0.24g		0.0100	
03	3	0.040			0.23g		0.0100	
04	4	0.064			0.23g		0.0100	
05	4	0.064			0.20g		0.0100	
SAND BLANK	1	-0.008			0.20g		0.0100	
CCV	35	0.802		ug/L	100ml	0.0002mg/L	0.00080	80.2%
CCB	2	0.010					0.0002	

Continued on Page _____

Read and Understood By

$$\text{conc} = \frac{\text{Read} - 1.315}{42.000}$$


Signed

1-28-93

Date

Signed

Date

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-A
 HAS Sample #93-0104-01
 Date Sampled: 1/19/93
 Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL VOL.(L)	(ml) (gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC	SPIKE %REC	RPD
ALUMINUM	6010	1/28/93	2910	2	0.698	30.0	86000	8340000	*95		
ANTIMONY	6010	1/28/93	9.91	2	0.698	50.0	143000	<DL	*95		<1.0
ARSENIC	7060	2/02/93	9.5	2	0.698	10.0	28700	<DL	*95	96.8	
BARIUM	6010	1/28/93	42.8	2	0.698	10.0	28700	123000	*95		<1.0
BERYLLIUM	6010	1/28/93	1.76	2	0.698	5.00	14300	<DL	*95		<1.0
CADMIUM	6010	1/28/93	2.6	2	0.698	5.00	14300	<DL	*95		<1.0
CALCIUM	6010	1/28/93	26100	2	0.698	20.0	57300	74800000	*95		
CHROMIUM	6010	1/28/93	13.1	2	0.698	10.0	28700	37500	*95		
COBALT	6010	1/28/93	7.15	2	0.698	10.0	28700	<DL	*95		<1.0
COPPER	6010	1/28/93	18.3	2	0.698	10.0	28700	52400	*95		
IRON	6010	1/28/93	22400	2	0.698	20.0	57300	64200000	*95		
LEAD	7421	2/01/93	7.5 ✓	20 ✓	0.698	5.00	143000	215000	*95		
MAGNESIUM	6010	1/28/93	2920	2	0.698	40.0	115000	8370000	*95		3.5
MANGANESE	6010	1/28/93	475	2	0.698	10.0	28700	1360000	*95		
MERCURY	7471	1/27/93	0.207	0.1	0.147	0.20	136	141	*95	92.8	<1.0
NICKEL	6010	1/28/93	-6.02	2	0.698	40.0	115000	<DL	*95		<1.0
POTASSIUM	6010	1/28/93	1330	2	0.698	3000	8600000	<DL	*95		<1.0
SELENIUM	7740	2/01/93	1	0.2	0.698	5.00	1430	<DL	*95	90.2	<1.0
SILVER	6010	1/28/93	-1.14	2	0.698	10.0	28700	<DL	*95		<1.0
SODIUM	6010	1/28/93	-39.7	2	0.698	50.0	143000	<DL	*95		<1.0
THALLIUM	7841	2/02/93	0.7	0.2	0.698	10.0	2870	<DL	*95		<1.0
VANADIUM	6010	1/28/93	15.5	2	0.698	20.0	57300	<DL	*95		<1.0
ZINC	6010	1/28/93	159	2	0.698	20.0	57300	456000	*95		

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN
 EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR
 SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON
 DRY WEIGHT.

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-B
 HAS Sample #93-0104-02
 Date Sampled: 1/19/93
 Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL VOL.(L)	(ml) (gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC
ALUMINUM	6010	1/28/93	15000	2	0.834	30.0	71900	36000000	*95
ANTIMONY	6010	1/28/93	-15.7	2	0.834	50.0	120000	<DL	*95
ARSENIC	7060	2/02/93	0.5	2	0.834	10.0	24000	<DL	*95
BARIUM	6010	1/28/93	113	2	0.834	10.0	24000	271000	*95
BERYLLIUM	6010	1/28/93	5.42	2	0.834	5.00	12000	13000	*95
CADMIUM	6010	1/28/93	0.0001	2	0.834	5.00	12000	<DL	*95
CALCIUM	6010	1/28/93	91200	2	0.834	20.0	48000	219000000	*95
CHROMIUM	6010	1/28/93	19.5	2	0.834	10.0	24000	46800	*95
COBALT	6010	1/28/93	-2.44	2	0.834	10.0	24000	<DL	*95
COPPER	6010	1/28/93	6.03	2	0.834	10.0	24000	<DL	*95
IRON	6010	1/28/93	11200	2	0.834	20.0	48000	26900000	*95
LEAD	7421	2/01/93	1.1 ✓	20 ✓	0.834	5.00	120000	<DL	*95
MAGNESIUM	6010	1/28/93	22700	2	0.834	40.0	95900	54400000	*95
MANGANESE	6010	1/28/93	930	2	0.834	10.0	24000	2230000	*95
MERCURY	7471	1/27/93	0.064	0.1	0.192	0.20	104	<DL	*95
NICKEL	6010	1/28/93	-2.93	2	0.834	40.0	95900	<DL	*95
POTASSIUM	6010	1/28/93	1770	2	0.834	3000	7190000	<DL	*95
SELENIUM	7740	2/01/93	1.7	2	0.834	5.00	12000	<DL	*95
SILVER	6010	1/28/93	-1.14	2	0.834	10.0	24000	<DL	*95
SODIUM	6010	1/28/93	626	2	0.834	50.0	120000	1500000	*95
THALLIUM	7841	2/02/93	-0.1	0.2	0.834	10.0	2400	<DL	*95
VANADIUM	6010	1/28/93	3.12	2	0.834	20.0	48000	<DL	*95
ZINC	6010	1/28/93	13.1	2	0.834	20.0	48000	<DL	*95

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN
 EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR
 SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON
 DRY WEIGHT.

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-C
HAS Sample #93-0104-03
Date Sampled: 1/19/93
Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL (ml) VOL.(L)	(gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC
ALUMINUM	6010	1/28/93	6860	2	0.797	30.0	75300	17200000	*95
ANTIMONY	6010	1/28/93	-3.3	2	0.797	50.0	125000	<DL	*95
ARSENIC	7060	2/02/93	2.6	2	0.797	10.0	25100	<DL	*95
BARIUM	6010	1/28/93	56.2	2	0.797	10.0	25100	141000	*95
BERYLLIUM	6010	1/28/93	2.85	2	0.797	5.00	12500	<DL	*95
CADMIUM	6010	1/28/93	0.434	2	0.797	5.00	12500	<DL	*95
CALCIUM	6010	1/28/93	91500	2	0.797	20.0	50200	230000000	*95
CHROMIUM	6010	1/28/93	241	2	0.797	10.0	25100	605000	*95
COBALT	6010	1/28/93	8.43	2	0.797	10.0	25100	<DL	*95
COPPER	6010	1/28/93	32.5	2	0.797	10.0	25100	81600	*95
IRON	6010	1/28/93	56300	2	0.797	20.0	50200	141000000	*95
oh LEAD	7421	2/01/93	2.6	20	0.797	5.00	125000	<DL	*95
MAGNESIUM	6010	1/28/93	9850	2	0.797	40.0	100000	24700000	*95
MANGANESE	6010	1/28/93	3890	2	0.797	10.0	25100	9760000	*95
MERCURY	7471	1/27/93	0.04	0.1	0.189	0.20	106	<DL	*95
NICKEL	6010	1/28/93	-27	2	0.797	40.0	100000	<DL	*95
POTASSIUM	6010	1/28/93	1640	2	0.797	3000	7530000	<DL	*95
SELENIUM	7740	2/01/93	0.4	2	0.797	5.00	12500	<DL	*95
SILVER	6010	1/28/93	-10.3	2	0.797	10.0	25100	<DL	*95
SODIUM	6010	1/28/93	191	2	0.797	50.0	125000	479000	*95
THALLIUM	7841	2/02/93	-1	0.2	0.797	10.0	2510	<DL	*95
VANADIUM	6010	1/28/93	81	2	0.797	20.0	50200	203000	*95
ZINC	6010	1/28/93	22.7	2	0.797	20.0	50200	57000	*95

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN
EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR
SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON
DRY WEIGHT.

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-D

HAS Sample #93-0104-04

Date Sampled: 1/19/93

Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL (ml) VOL.(L)	(gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC
ALUMINUM	6010	1/28/93	11600	2	0.827	30.0	72600	28100000	*95
ANTIMONY	6010	1/28/93	7.43	2	0.827	50.0	121000	<DL	*95
ARSENIC	7060	2/02/93	3.1	2	0.827	10.0	24200	<DL	*95
BARIUM	6010	1/28/93	90.8	2	0.827	10.0	24200	220000	*95
BERYLLIUM	6010	1/28/93	4.47	2	0.827	5.00	12100	<DL	*95
CADMIUM	6010	1/28/93	-0.434	2	0.827	5.00	12100	<DL	*95
CALCIUM	6010	1/28/93	76700	2	0.827	20.0	48400	185000000	*95
CHROMIUM	6010	1/28/93	54	2	0.827	10.0	24200	131000	*95
COBALT	6010	1/28/93	2.98	2	0.827	10.0	24200	<DL	*95
COPPER	6010	1/28/93	15.5	2	0.827	10.0	24200	37500	*95
IRON	6010	1/28/93	37200	2	0.827	20.0	48400	90000000	*95
LEAD	7421	2/01/93	1.2 ✓	20 ✓	0.827	5.00	121000	<DL	*95
MAGNESIUM	6010	1/28/93	17500	2	0.827	40.0	96700	42300000	*95
MANGANESE	6010	1/28/93	1190	2	0.827	10.0	24200	2880000	*95
MERCURY	7471	1/27/93	0.064	0.1	0.198	0.20	101	<DL	*95
NICKEL	6010	1/28/93	-16.6	2	0.827	40.0	96700	<DL	*95
POTASSIUM	6010	1/28/93	1670	2	0.827	3000	7260000	<DL	*95
SELENIUM	7740	2/01/93	-0.1	2	0.827	5.00	12100	<DL	*95
SILVER	6010	1/28/93	-4.56	2	0.827	10.0	24200	<DL	*95
SODIUM	6010	1/28/93	329	2	0.827	50.0	121000	796000	*95
THALLIUM	7841	2/02/93	0.3	0.2	0.827	10.0	2420	<DL	*95
VANADIUM	6010	1/28/93	16.3	2	0.827	20.0	48400	<DL	*95
ZINC	6010	1/28/93	8.54	2	0.827	20.0	48400	<DL	*95

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON DRY WEIGHT.

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-OS

HAS Sample #93-0104-05

Date Sampled: 1/19/93

Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL VOL.(L)	(ml) (gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC
ALUMINUM	6010	1/28/93	2720	2	0.853	30.0	70300	6380000	*95
ANTIMONY	6010	1/28/93	7.43	2	0.853	50.0	117000	<DL	*95
ARSENIC	7060	2/02/93	3.4	2	0.853	10.0	23400	<DL	*95
BARIUM	6010	1/28/93	42.8	2	0.853	10.0	23400	100000	*95
BERYLLIUM	6010	1/28/93	1.63	2	0.853	5.00	11700	<DL	*95
CADMIUM	6010	1/28/93	0.0001	2	0.853	5.00	11700	<DL	*95
CALCIUM	6010	1/28/93	76900	2	0.853	20.0	46900	180000000	*95
CHROMIUM	6010	1/28/93	476	2	0.853	10.0	23400	1120000	*95
COBALT	6010	1/28/93	1.92	2	0.853	10.0	23400	<DL	*95
COPPER	6010	1/28/93	22.3	2	0.853	10.0	23400	52300	*95
IRON	6010	1/28/93	48500	2	0.853	20.0	46900	114000000	*95
LEAD	7421	2/01/93	6.3 ✓	20 ✓	0.853	5.00	117000	148000	*95
MAGNESIUM	6010	1/28/93	5530	2	0.853	40.0	93800	13000000	*95
MANGANESE	6010	1/28/93	6480	2	0.853	10.0	23400	15200000	*95
MERCURY	7471	1/27/93	0.064	0.1	0.171	0.20	117	<DL	*95
NICKEL	6010	1/28/93	-28	2	0.853	40.0	93800	<DL	*95
POTASSIUM	6010	1/28/93	624	2	0.853	3000	7030000	<DL	*95
SELENIUM	7740	2/01/93	4	0.2	0.853	5.00	1170	<DL	*95
SILVER	6010	1/28/93	-10.3	2	0.853	10.0	23400	<DL	*95
SODIUM	6010	1/28/93	52.2	2	0.853	50.0	117000	122000	*95
THALLIUM	7841	2/02/93	-0.9	0.2	0.853	10.0	2340	<DL	*95
VANADIUM	6010	1/28/93	95.3	2	0.853	20.0	46900	223000	*95
ZINC	6010	1/28/93	43.6	2	0.853	20.0	46900	102000	*95

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN
EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR
SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON
DRY WEIGHT.

HUNTINGDON ANALYTICAL SERVICES

Sample ID: SOIL BLANK
 HAS Sample #93-0104-SB
 Date Sampled: NA
 Date Prepared: 1/26/93

ANALYTE	EPA METHOD	DATE ANALYZE	RESULT	FINAL VOL.(L)	(ml) (gm) DRY WT.	METHO DL	DET.LIMIT (ug/kg)	RESULT ug/kg	QC
ALUMINUM	6010	1/28/93	22.4	0.2	1.000	30.0	6000	<DL	*95
ANTIMONY	6010	1/28/93	0.0001	0.2	1.000	50.0	10000	<DL	*95
ARSENIC	7060	2/02/93	-1.8	0.2	1.000	10.0	2000	<DL	*95
BARIUM	6010	1/28/93	0.222	0.2	1.000	10.0	2000	<DL	*95
BERYLLIUM	6010	1/28/93	1.49	0.2	1.000	5.00	1000	<DL	*95
CADMIUM	6010	1/28/93	0.434	0.2	1.000	5.00	1000	<DL	*95
CALCIUM	6010	1/28/93	-5.61	0.2	1.000	20.0	4000	<DL	*95
CHROMIUM	6010	1/28/93	-2.84	0.2	1.000	10.0	2000	<DL	*95
COBALT	6010	1/28/93	0.344	0.2	1.000	10.0	2000	<DL	*95
COPPER	6010	1/28/93	0.928	0.2	1.000	10.0	2000	<DL	*95
IRON	6010	1/28/93	-13.3	0.2	1.000	20.0	4000	<DL	*95
LEAD	7421	2/01/93	1.1	0.2	1.000	5.00	1000	<DL	*95
MAGNESIUM	6010	1/28/93	10.6	0.2	1.000	40.0	8000	<DL	*95
MANGANESE	6010	1/28/93	-0.345	0.2	1.000	10.0	2000	<DL	*95
MERCURY	7471	1/27/93	-8E-06	0.1	0.200	0.20	100	<DL	*95
NICKEL	6010	1/28/93	-6.72	0.2	1.000	40.0	8000	<DL	*95
POTASSIUM	6010	1/28/93	67.2	0.2	1.000	3000	600000	<DL	*95
SELENIUM	7740	2/01/93	-1	0.2	1.000	5.00	1000	<DL	*95
SILVER	6010	1/28/93	1.14	0.2	1.000	10.0	2000	<DL	*95
SODIUM	6010	1/28/93	-138	0.2	1.000	50.0	10000	<DL	*95
THALLIUM	7841	2/02/93	-0.7	0.2	1.000	10.0	2000	<DL	*95
VANADIUM	6010	1/28/93	0.321	0.2	1.000	20.0	4000	<DL	*95
ZINC	6010	1/28/93	-0.296	0.2	1.000	20.0	4000	<DL	*95

*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN
 EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR
 SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON
 DRY WEIGHT.

PRQ

PERCENT SOLID - DRY WEIGHT - 160.3

Book 654
Page 34

DATE/TIME/INIT. IN: 1/28/93 1:30pm PB
DATE/TIME/INIT. OUT: 1/29/93 9:00am PB

HAS ID	QC	Dish Wt. (grams)	ID	Wet Weight Sample only (g)	Dry Dish+ Residue (g) 103	Residue (g) 103	% Solid dry weight	
0104-01		0.9900	8	6.5600	5.5700	4.5800	69.8	70.0
0104-01	dup	0.9900	9	8.0700	6.6500	5.6600	70.1	0.5
0104-02		0.9900	10	5.6100	5.6700	4.6800	83.4	
0104-03		0.9900	11	7.8700	7.2000	6.2100	78.9	
0104-04		0.9900	12	6.8800	6.6800	5.6900	82.7	
0104-05		0.9900	13	6.2500	6.3200	5.3300	85.3	

~~Titles, Window, Status, Page~~ PERCENT SOLID - DRY WEIGHT - 160.3
~~Book 654~~

~~Page 34~~

1-29-93
HS

ICP RAW DATA

#	Sample Name	File	Method	Date	Time	OpID	Type	Mode
1	CAL BLANK	P112831	LOW-MIP	01/28/93	10:21	mip	S	CONC
2	ICV-1	P112831	LOW-MIP	01/28/93	10:27	mip	S	CONC
3	ICV-2	P112831	LOW-MIP	01/28/93	10:56	mip	S	CONC
4	ICV-7	P112831	LOW-MIP	01/28/93	11:04	mip	S	CONC
5	0028.01 NIOSH-P.E.	P112831	LOW-MIP	01/28/93	11:15	mip	S	CONC
6	0028.02	P112831	LOW-MIP	01/28/93	11:30	mip	S	CONC
7	0028.03	P112831	LOW-MIP	01/28/93	11:50	mip	S	CONC
8	0028.04	P112831	LOW-MIP	01/28/93	12:08	mip	S	CONC
9	0028.05	P112831	LOW-MIP	01/28/93	12:28	mip	S	CONC
10	CAL BLANK	P112831	LOW-MIP	01/28/93	12:47	mip	S	CONC
11	CCV	P112831	LOW-MIP	01/28/93	12:58	mip	S	CONC
12	CAL BLANK	1 P112831	LOW-MIP	01/28/93	14:49	mip	S	CONC
13	CCV	1 P112831	LOW-MIP	01/28/93	15:00	mip	S	CONC
14	PRP.BLNK 3050P 01/26	1 P112831	LOW-MIP	01/28/93	15:27	mip	S	CONC
15	LCS-S 01/26	1 P112831	LOW-MIP	01/28/93	15:38	mip	S	CONC
16	0104.01 (1/10)	1 P112831	LOW-MIP	01/28/93	16:00	mip	S	CONC
17	0104.01 (1/10) REP.1	1 P112831	LOW-MIP	01/28/93	16:16	mip	S	CONC
18	0104.02 (1/10)	1 P112831	LOW-MIP	01/28/93	16:32	mip	S	CONC
19	0104.03 (1/10)	1 P112831	LOW-MIP	01/28/93	16:49	mip	S	CONC
20	0104.04 (1/10)	1 P112831	LOW-MIP	01/28/93	17:05	mip	S	CONC
21	0104.05 (1/10)	1 P112831	LOW-MIP	01/28/93	17:21	mip	S	CONC
22	CAL BLANK	1 P112831	LOW-MIP	01/28/93	17:49	mip	S	CONC
23	CCV	1 P112831	LOW-MIP	01/28/93	17:59	mip	S	CONC
24	0102.01 (1/10)	1 P112831	LOW-MIP	01/28/93	18:21	mip	S	CONC
25	0102.02 (1/10)	1 P112831	LOW-MIP	01/28/93	18:37	mip	S	CONC
26	0102.02 (1/10) REP.1	1 P112831	LOW-MIP	01/28/93	18:54	mip	S	CONC
27	0102.03 (1/10)	1 P112831	LOW-MIP	01/28/93	19:10	mip	S	CONC
28	0102.04 (1/10)	1 P112831	LOW-MIP	01/28/93	19:26	mip	S	CONC
29	0102.05 (1/10)	1 P112831	LOW-MIP	01/28/93	19:43	mip	S	CONC
30	0102.06 (1/10)	1 P112831	LOW-MIP	01/28/93	19:59	mip	S	CONC
31	CAL BLANK	1 P112831	LOW-MIP	01/28/93	20:26	mip	S	CONC
32	CCV	1 P112831	LOW-MIP	01/28/93	20:37	mip	S	CONC
33	0102.07 (1/10)	1 P112831	LOW-MIP	01/28/93	20:59	mip	S	CONC
34	0102.08 (1/10)	1 P112831	LOW-MIP	01/28/93	21:15	mip	S	CONC
35	0102.09 (1/10)	1 P112831	LOW-MIP	01/28/93	21:31	mip	S	CONC
36	0102.10 (1/10)	1 P112831	LOW-MIP	01/28/93	21:48	mip	S	CONC
37	0102.11 (1/10)	1 P112831	LOW-MIP	01/28/93	22:04	mip	S	CONC
38	0102.12 (1/10)	1 P112831	LOW-MIP	01/28/93	22:20	mip	S	CONC
39	0102.13 (1/10)	1 P112831	LOW-MIP	01/28/93	22:37	mip	S	CONC
40	CAL BLANK	1 P112831	LOW-MIP	01/28/93	23:04	mip	S	CONC
41	CAL BLANK	1 P112831	LOW-MIP	01/28/93	23:15	mip	S	CONC
42	CAL BLANK	1 P112831	LOW-MIP	01/28/93	23:25	mip	S	CONC
43	CCV	1 P112831	LOW-MIP	01/28/93	23:36	mip	S	CONC
44	CCV	1 P112831	LOW-MIP	01/28/93	23:47	mip	S	CONC
45	CCV	1 P112831	LOW-MIP	01/28/93	23:58	mip	S	CONC

#	Sample Name	Al3082	Sb2068	As1936	Ba4934	Be3130	B_2496
1	CAL BLANK	L-10.8	L4.95	L14.3	L.222	L.813	L1.91
2	ICV-1	1830.	550.	97.6	2010.	51.2	30.6
3	ICV-2	L24.8	L.000	512.	L.222	L.813	L10.5
4	ICV-7	991.	L-5.78	L18.5	1050.	L.407	1120.
5	0028.01 NIOSH-P.E.	L6.56	L5.78	L11.3	L.222	L.497	L13.0
6	0028.02	L7.97	L9.91	L8.90	L.222	L.723	L8.91
7	0028.03	L8.78	L10.5	L10.1	L.222	L.768	L8.58
8	0028.04	L1.70	L.826	L4.69	L.222	L.587	L11.1
9	0028.05	L5.61	L.275	L14.2	L.222	L.632	L10.5
10	CAL BLANK	L-2.07	L-1.65	L6.95	L.222	L.632	L8.90
11	CCV	954.	988.	980.	998.	1010.	1060.
12	CAL BLANK	1 L17.0	L-9.08	L16.0	L.222	L.813	L4.78
13	CCV	1 926.	943.	936.	950.	967.	998.
14	PRP.BLNK 3050P 01/26	L22.4	L.000	L8.83	L.222	L1.49	L8.59
15	LCS-S 01/26	1 19700.	228.	L8.92	462.	290.	L14.4
16	0104.01 (1/10)	1 2910.	L9.91	L1.57	42.8	L1.76	L8.44
17	0104.01 (1/10) REP.	1 3420.	L4.95	L1.85	42.8	L1.22	L12.0
18	0104.02 (1/10)	1 15000.	L-15.7	L-85.8	113.	L5.42	L25.7
19	0104.03 (1/10)	1 6860.	L-3.30	L-49.9	56.2	L2.85	L4.18
20	0104.04 (1/10)	1 11600.	L7.43	L-69.0	90.8	L4.47	L12.9
21	0104.05 (1/10)	1 2720.	L7.43	L-21.9	42.8	L1.63	L9.17
22	CAL BLANK	1 L24.6	L-4.13	L5.87	L.222	L.678	L-9.56
23	CCV	1 925.	956.	942.	948.	973.	989.
24	0102.01 (1/10)	1 657.	L-.826	L-.099	387.	L3.25	L15.7
25	0102.02 (1/10)	1 10500.	L-30.6	L-65.5	211.	L1.22	L15.0
26	0102.02 (1/10) REP.	1 10600.	L-29.7	L-64.1	174.	L.813	L16.2
27	0102.03 (1/10)	1 127.	L-9.08	L12.6	L19.5	L.813	L-11.2
28	0102.04 (1/10)	1 2170.	L1.65	L-7.80	728.	L.813	L15.2
29	0102.05 (1/10)	1 2460.	L.826	L-5.99	146.	L.813	L12.4
30	0102.06 (1/10)	1 240.	L5.78	L7.78	30.4	L.813	L-6.59
31	CAL BLANK	1 L27.4	L-.826	L21.6	L.222	L.813	L-9.55
32	CCV	1 899.	961.	921.	930.	962.	976.
33	0102.07 (1/10)	1 639.	L1.65	L2.18	85.5	L2.44	L13.9
34	0102.08 (1/10)	1 151.	L4.13	L21.8	L18.2	L.813	L-2.74
35	0102.09 (1/10)	1 175.	L-9.08	L7.91	L12.4	L.813	L-11.1
36	0102.10 (1/10)	1 70.8	L-14.0	L11.5	52.8	L.813	L-11.1
37	0102.11 (1/10)	1 124.	L-5.78	L14.2	L10.2	L.813	L-3.81
38	0102.12 (1/10)	1 282.	L10.7	L36.8	201.	L.813	L15.1
39	0102.13 (1/10)	1 186.	L-4.13	L5.30	L14.0	L.813	L-11.5
40	CAL BLANK	1 L25.0	L1.65	L15.4	L.222	L.813	L-9.53
41	CAL BLANK	1 L29.4	L5.78	L16.2	L.222	L.813	L-9.53
42	CAL BLANK	1 L29.9	L4.95	L14.1	L.222	L.678	L-9.53
43	CCV	1 955.	973.	929.	967.	979.	1020.
44	CCV	1 925.	965.	941.	951.	977.	993.
45	CCV	1 952.	979.	955.	973.	995.	1010.

#	Sample Name	Cd2288	Ca3179	Cr2677	Co2286	Cu3247	Fe2599
1	CAL BLANK	L-1.30	L2.75	L-.875	L-1.03	L.464	L6.21
2	ICV-1	51.6	50000.	93.4	501.	243.	1030.
3	ICV-2	L4.34	27.8	L-3.50	L.688	L.464	L-17.5
4	ICV-7	L3.47	L10.7	L-3.94	L-.687	L-.696	L-17.5

#	Sample Name		Cd2288	Ca3179	Cr2677	Co2286	Cu3247	Fe2599
5	0028.01 NIOSH-P.E.		29.4	L15.6	454.	L-.344	L.155	L-12.2
6	0028.02		44.0	26.1	637.	L.917	L.851	L-7.21
7	0028.03		93.6	33.4	1130.	L.344	L.155	L14.7
8	0028.04		62.4	L17.0	779.	L-.917	L-.155	L-7.61
9	0028.05		L.000	L19.2	L-2.33	L.802	L.232	L-3.46
10	CAL BLANK		L-.868	L-.165	L-.729	L-1.49	L-1.08	L10.7
11	CCV		1020.	1020.	1000.	1030.	979.	1010.
12	CAL BLANK	1	L.000	L-7.87	L-3.94	L.344	L-1.16	L-3.56
13	CCV	1	973.	973.	952.	968.	941.	973.
14	PRP.BLNK 3050P 01/26	1	L.434	L-5.61	L-2.84	L.344	L.928	L-13.3
15	LCS-S 01/26	1	389.	11000.	816.	629.	454.	45200.
16	0104.01 (1/10)	1	L2.60	26100.	13.1	L7.15	18.3	22400.
17	0104.01 (1/10)	REP.1	L.000	13900.	21.2	L6.44	26.5	30000.
18	0104.02 (1/10)	1	L.000	91200.	19.5	L-2.44	L6.03	11200.
19	0104.03 (1/10)	1	L.434	91500.	241.	L8.43	32.5	56300.
20	0104.04 (1/10)	1	L-.434	76700.	54.0	L2.98	15.5	37200.
21	0104.05 (1/10)	1	L.000	76900.	476.	L1.92	22.3	48500.
22	CAL BLANK	1	L.000	97.8	L-4.81	L-3.09	L-1.62	42.4
23	CCV	1	978.	995.	953.	978.	940.	982.
24	0102.01 (1/10)	1	L4.77	39100.	35.9	L7.49	38.3	23600.
25	0102.02 (1/10)	1	L.000	5320.	215.	65.6	10000.	36000.
26	0102.02 (1/10)	REP.1	L-.434	2800.	219.	72.4	12700.	37400.
27	0102.03 (1/10)	1	L.000	212.	L-1.75	L.334	94.5	3170.
28	0102.04 (1/10)	1	L4.34	13700.	87.3	14.3	78.9	35100.
29	0102.05 (1/10)	1	L.868	12100.	22.3	L6.78	72.6	32900.
30	0102.06 (1/10)	1	L.000	780.	L4.59	L1.69	57.3	10500.
31	CAL BLANK	1	L.000	L1.95	L-5.25	L-1.38	L-.928	L5.05
32	CCV	1	970.	975.	952.	974.	918.	974.
33	0102.07 (1/10)	1	L4.34	1790.	15.3	L5.12	44.6	11400.
34	0102.08 (1/10)	1	L.434	364.	L-1.09	L-.697	31.6	3330.
35	0102.09 (1/10)	1	L.868	336.	L-2.84	L-.696	34.1	2860.
36	0102.10 (1/10)	1	L.434	831.	L-1.31	L1.37	37.4	2860.
37	0102.11 (1/10)	1	L.868	684.	L-3.06	L-2.08	68.2	5320.
38	0102.12 (1/10)	1	L.868	1370.	L7.87	L8.23	28.3	5510.
39	0102.13 (1/10)	1	L.000	388.	L-.875	L-1.04	19.7	3730.
40	CAL BLANK	1	L.434	L-5.18	L-5.47	L-2.06	L-1.62	L-8.90
41	CAL BLANK	1	L.000	L-3.69	L-4.81	L-.687	L-.928	L-14.8
42	CAL BLANK	1	L.868	L.353	L-2.41	L2.06	L1.39	L-17.8
43	CCV	1	989.	980.	964.	987.	968.	959.
44	CCV	1	980.	980.	962.	983.	942.	973.
45	CCV	1	997.	993.	976.	993.	969.	970.

#	Sample Name		Pb2203	Mg2790	Mn2576	Mo2020	Ni2316	K 7664
1	CAL BLANK		L-6.79	L-28.5	L-.345	L.647	L.827	L-62.0
2	ICV-1		48.0	49100.	154.	L-2.16	402.	55700.
3	ICV-2		529.	L6.36	L-.345	L1.29	L-8.40	L244.
4	ICV-7		L13.1	L-5.27	L-.345	L.647	L-5.04	11000.
5	0028.01 NIOSH-P.E.		425.	L-.342	L-.345	L-.647	L-8.97	L34.9
6	0028.02		175.	L3.53	L-.345	L1.58	L-5.88	L122.
7	0028.03		299.	L-.012	L.805	L.719	L-5.34	L-23.0
8	0028.04		128.	L-6.69	L-.345	L.287	L-6.45	L239.
9	0028.05		L2.40	L.707	L-.345	L.431	L-2.52	L-32.5

Sample Name	Pb2203	Mg2790	Mn2576	Mo2020	Ni2316	K_7664
10 CAL BLANK	L-17.7	L-10.2	L-.345	L-.359	L1.11	L-343.
11 CCV	1040.	981.	1010.	1020.	1030.	9330.
12 CAL BLANK	1 L-19.1	L-18.0	L-.345	L-.862	L.837	L-80.9
13 CCV	1 965.	956.	965.	960.	981.	9040.
14 PRP.BLNK 3050P 01/26	(L-8.77)	L10.6	L-.345	L3.67	L-6.72	L67.2
15 LCS-S 01/26	1 244.	7380.	715.	268.	682.	12200.
16 0104.01 (1/10)	1 64.1	2920.	475.	(L2.16)	L-6.02	L1330.
17 0104.01 (1/10) REP.	1 111.	2820.	726.	(L5.17)	L-9.10	L1500.
18 0104.02 (1/10)	1 L23.2	22700.	930.	(L-.647)	L-2.93	L1770.
19 0104.03 (1/10)	1 L16.0	9850.	3890.	11.4	L-27.0	L1640.
20 0104.04 (1/10)	1 L5.20	17500.	1190.	L1.51	L-16.6	L1670.
21 0104.05 (1/10)	1 68.8	5530.	6480.	L7.76	L-28.0	L624.
22 CAL BLANK	1 L-8.78	L-3.21	L7.59	L-.862	L-15.2	L131.
23 CCV	1 976.	928.	974.	963.	981.	9130.
24 0102.01 (1/10)	1 243.	1630.	115.	L3.88	L7.92	L178.
25 0102.02 (1/10)	1 177.	611.	425.	22.6	2920.	L100.0
26 0102.02 (1/10) REP.	1 153.	486.	369.	37.5	2820.	L294.
27 0102.03 (1/10)	1 L.165	L37.6	16.2	L.647	L16.6	L-97.4
28 0102.04 (1/10)	1 593.	1920.	207.	L5.17	L17.7	L1160.
29 0102.05 (1/10)	1 237.	2680.	163.	L3.02	L16.0	L737.
30 0102.06 (1/10)	1 43.7	48.5	46.6	L.647	L.255	L274.
31 CAL BLANK	1 L-30.7	L24.3	L-.345	L-1.72	L-5.88	L511.
32 CCV	1 991.	970.	962.	960.	981.	8820.
33 0102.07 (1/10)	1 317.	277.	71.1	L3.67	L3.99	L-18.3
34 0102.08 (1/10)	1 L2.93	L22.6	16.2	L-.431	L-5.41	L386.
35 0102.09 (1/10)	1 L19.0	L28.3	16.2	L.431	L3.29	L627.
36 0102.10 (1/10)	1 96.8	L36.8	15.2	L.647	L4.14	L573.
37 0102.11 (1/10)	1 L26.0	53.8	24.5	L-.862	L-.737	L98.0
38 0102.12 (1/10)	1 62.0	149.	32.8	L2.80	L6.75	L997.
39 0102.13 (1/10)	1 L32.8	54.0	22.4	L-.647	L3.61	L-184.
40 CAL BLANK	1 L.410	L-15.8	L-.345	L-2.59	L-3.37	L-172.
41 CAL BLANK	1 L-27.1	L6.35	L-.345	L1.08	L-1.67	L264.
42 CAL BLANK	1 L-18.3	L22.2	L-.345	L1.08	L-1.66	L550.
43 CCV	1 971.	991.	978.	969.	981.	10200.
44 CCV	1 983.	964.	978.	972.	972.	9200.
45 CCV	1 986.	957.	992.	986.	1010.	9690.

Sample Name	Se1960	Si2516	Ag3280	Na5889	Tl1908	Sn1899
1 CAL BLANK	L39.7	L.967	L5.70	L-17.9	L28.6	L3.54
2 ICV-1	81.1	82.4	96.9	45900.	177.	L2.90
3 ICV-2	548.	35.4	L7.98	L7.61	483.	L2.58
4 ICV-7	L58.6	848.	1020.	1020.	L53.2	L13.2
5 0028.01 NIOSH-P.E.	L48.8	49.6	L7.22	L14.8	L39.6	L-2.90
6 0028.02	L54.8	35.8	L4.94	L6.15	L33.1	L36.1
7 0028.03	L53.3	69.3	L2.66	L-5.06	L37.6	L-8.70
8 0028.04	L53.4	109.	L3.04	L-17.0	L42.4	L-6.98
9 0028.05	L56.8	L25.2	L3.80	L-18.2	L39.8	L27.0
10 CAL BLANK	76.3	L1.63	L4.56	L-109.	L36.6	L7.09
11 CCV	1060.	5020.	988.	901.	985.	974.
12 CAL BLANK	1 L53.9	L-2.92	L3.42	L-124.	L19.4	L6.12
13 CCV	1 1010.	4840.	954.	940.	940.	924.
14 PRP.BLNK 3050P 01/26	76.2	L4.62	L1.14	L-138.	L27.8	L5.48

Sample Name	Se1960	Si2516	Ag3280	Na5889	Tl1908	Sn1899
15 LCS-S 01/26	1 204.	1490.	198.	1720.	406.	L18.7
16 0104.01 (1/10)	1 L57.7	149.	L-1.14	L-39.7	L-60.0	L11.6
17 0104.01 (1/10) REP.	1 L47.1	121.	L-5.70	L-51.6	L-93.7	L3.87
18 0104.02 (1/10)	1 L53.1	139.	L-1.14	626.	L6.62	L7.73
19 0104.03 (1/10)	1 L25.5	146.	L-10.3	191.	L-163.	L7.73
20 0104.04 (1/10)	1 L37.5	106.	L-4.56	329.	L-92.4	L11.9
21 0104.05 (1/10)	1 L31.8	112.	L-10.3	52.2	L-124.	L5.15
22 CAL BLANK	1 L59.4	L-3.24	L3.42	L-188.	L37.4	L5.15
23 CCV	1 1020.	4840.	954.	917.	926.	920.
24 0102.01 (1/10)	1 L51.0	84.1	L-1.14	230.	L-61.0	51.2
25 0102.02 (1/10)	1 L37.6	96.8	L-5.70	L-61.6	L-98.3	301.
26 0102.02 (1/10) REP.	1 L50.5	105.	L-1.14	L-85.3	L-105.	104.
27 0102.03 (1/10)	1 64.5	31.9	L1.14	145.	L14.0	L1.93
28 0102.04 (1/10)	1 L32.4	79.3	L-3.42	1070.	L-109.	84.7
29 0102.05 (1/10)	1 L35.5	68.6	L-5.70	1180.	L-106.	110.
30 0102.06 (1/10)	1 L43.1	68.8	L.000	11300.	L-13.2	L28.4
31 CAL BLANK	1 84.1	L-2.25	L6.84	L-169.	L45.2	L9.34
32 CCV	1 1030.	4780.	948.	860.	913.	917.
33 0102.07 (1/10)	1 64.4	86.3	L2.28	5290.	L-4.02	65.1
34 0102.08 (1/10)	1 69.2	55.2	L3.42	2860.	L19.3	L23.8
35 0102.09 (1/10)	1 78.9	78.8	L-3.42	6140.	L13.5	L16.1
36 0102.10 (1/10)	1 L55.3	47.9	L1.14	3480.	L25.5	L5.48
37 0102.11 (1/10)	1 L58.9	55.5	L5.70	400.	L8.71	L11.0
38 0102.12 (1/10)	1 78.7	98.3	L7.98	1300.	L12.9	L26.7
39 0102.13 (1/10)	1 L53.6	73.3	L.000	4100.	L9.45	L14.5
40 CAL BLANK	1 71.1	L-.948	L3.42	L-193.	L40.4	L.322
41 CAL BLANK	1 79.0	L.696	L4.56	L-163.	L45.4	L12.9
42 CAL BLANK	1 60.7	L1.35	L2.28	L-189.	L25.2	L5.15
43 CCV	1 1030.	4970.	971.	1070.	953.	931.
44 CCV	1 1010.	4880.	968.	875.	938.	940.
45 CCV	1 1040.	4970.	984.	996.	965.	954.

#	Sample Name	Ti3349	V_2924	Zn2138	N+1999
1	CAL BLANK	L1.13	L-.319	L-.302	.000
2	ICV-1	L1.41	508.	197.	.000
3	ICV-2	L1.69	L1.91	L.305	.000
4	ICV-7	L.846	L1.27	L.107	.000
5	0028.01 NIOSH-P.E.	L1.50	L.320	L.302	.000
6	0028.02	L2.73	L1.48	L2.10	.000
7	0028.03	L1.32	L.315	L2.90	.000
8	0028.04	L1.88	L1.38	L1.90	.000
9	0028.05	L2.73	L1.17	L2.90	.000
10	CAL BLANK	L.376	L.316	L-.101	.000
11	CCV	992.	1010.	1000.	.000
12	CAL BLANK	1 L-.846	L1.27	L-.897	.000
13	CCV	1 952.	946.	964.	.000
14	PRP BLNK 3050P 01/26	L2.26	L.321	L-.296	.000
15	LCS-S 01/26	1 1060.	179.	427.	.000
16	0104.01 (1/10)	1 23.1	L15.5	159.	.000
17	0104.01 (1/10) REP.	1 33.3	L17.1	234.	.000
18	0104.02 (1/10)	1 246.	L3.12	L13.1	.000
19	0104.03 (1/10)	1 501.	81.0	22.7	.000

#	Sample Name		Ti3349	V_2924	Zn2138	N+1999
20	0104.04 (1/10)	1	232.	L16.3	L8.54	.000
21	0104.05 (1/10)	1	541.	95.3	43.6	.000
22	CAL BLANK	1	L.282	L3.17	L-.902	.000
23	CCV	1	956.	951.	965.	.000
24	0102.01 (1/10)	1	31.9	L8.86	1860.	.000
25	0102.02 (1/10)	1	68.2	L4.10	1510.	.000
26	0102.02 (1/10)	REP.1	83.8	L6.68	1640.	.000
27	0102.03 (1/10)	1	L5.36	L1.26	123.	.000
28	0102.04 (1/10)	1	93.3	L6.50	2110.	.000
29	0102.05 (1/10)	1	40.9	L5.99	965.	.000
30	0102.06 (1/10)	1	L7.33	L1.35	459.	.000
31	CAL BLANK	1	L-.846	L3.81	L1.20	.000
32	CCV	1	947.	943.	962.	.000
33	0102.07 (1/10)	1	13.0	L1.80	400.	.000
34	0102.08 (1/10)	1	L5.36	L4.72	64.3	.000
35	0102.09 (1/10)	1	L4.79	L2.28	168.	.000
36	0102.10 (1/10)	1	L.564	L5.77	85.1	.000
37	0102.11 (1/10)	1	L1.69	L4.32	403.	.000
38	0102.12 (1/10)	1	L7.33	L3.64	196.	.000
39	0102.13 (1/10)	1	L4.23	L.829	57.1	.000
40	CAL BLANK	1	L-.282	L-.316	L-.896	.000
41	CAL BLANK	1	L.000	L2.86	L.305	.000
42	CAL BLANK	1	L.846	L2.23	L-.894	.000
43	CCV	1	970.	960.	973.	.000
44	CCV	1	963.	955.	975.	.000
45	CCV	1	980.	970.	986.	.000

ID Weight File: 0059.IDW
Sample Volume: 100 mL

Analyst: LEAD
Nominal Weight: 1.0 g

pc.	Sample ID	Weight	Dilution
1	PB 3020 1/28/93		
2	LCSW 1/28/93		
3	0059-06		
4	0059-06REP		
5	0059-06 SPK		
6	0059-07		
7	0059-08		
8	0059-09		
9	0059-10		
10	PBS 1/26/93		
11	LCSS 1/26 1-10		10
12	0104-01 1-100		100
13	0104-01 REF 100		100
14	0104-01 SPK 100		100
15	0104-02 100		100
16	0104-03 1-100		100
17	0104-04 1-100		100
18	0104-05 1-100		100

*All
analysis.
took place
in 2/1/93.*

Element File: PBBRIEF.GEL Element: Pb Wavelength: 283.3
Date: 01/30/93 Time: 13:14 Slit: 0.70 L
Data File: F320131.DAT ID/Wt File: 0059.IDW Lamp Current: 10
Technique: HGA Calib. Type: Nonlinear Energy: 64
Remark 1: Matrix modifier: 4% NH4PO4 .2% MgNO3 Integration time: 3 s
Remark 2: Lamp: Temp. program:
Remark 3: Energy: Current:
Remark 4: Char. Mass: 20 ppb:
Remark 5: 60 ppb: 100 ppb:

~~~~~  
Pb      ID: Calib. Blank      Seq. No.: 00001      A/S Pos.: 37      Date: 01/30/93  
uL dispensed: 10 from 0. 5 from 38. 20 from 37

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.001  
Background Pk Area (A-s): 0.005  
Blank Corrected Pk Area (A-s): 0.001

Time: 13:16  
Peak Height (A): 0.003  
Background Pk Height (A): 0.006

Auto-zero performed.

~~~~~  
Pb ID: 20 ppb Seq. No.: 00002 A/S Pos.: 36 Date: 01/30/93

uL dispensed: 26 from 0.5 from 38.4 from 36

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.053
Background Pk Area (A-s): 0.014
Blank Corrected Pk Area (A-s): 0.052

Time: 13:19
Peak Height (A): 0.182
Background Pk Height (A): 0.034

Standard number 1 applied. [20.0]
Correlation coefficient: 1.00000

Slope: 0.0026

~~~~~  
Pb ID: 60 ppb Seq. No.: 00003 A/S Pos.: 36 Date: 01/30/93

uL dispensed: 18 from 0.5 from 38.12 from 36

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.130  
Background Pk Area (A-s): 0.033  
Blank Corrected Pk Area (A-s): 0.130  
Concentration (ug/L): 49.8

Time: 13:23  
Peak Height (A): 0.242  
Background Pk Height (A): 0.048

Standard number 2 applied. [60.0]  
Correlation coefficient: 1.00000

Slope: 0.0029

~~~~~  
Pb ID: 100ppb Seq. No.: 00004 A/S Pos.: 36 Date: 01/30/93

uL dispensed: 10 from 0.5 from 38.20 from 36

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.203
Background Pk Area (A-s): 0.047
Blank Corrected Pk Area (A-s): 0.202
Concentration (ug/L): 115.7

Time: 13:26
Peak Height (A): 0.379
Background Pk Height (A): 0.078

Standard number 3 applied. [100.0]
Correlation coefficient: 1.00000

Slope: 0.0045

~~~~~  
Pb ID: ICV Seq. No.: 00005 A/S Pos.: 35 Date: 01/30/93

uL dispensed: 10 from 0.5 from 38.20 from 35

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.106  
Background Pk Area (A-s): 0.033  
Blank Corrected Pk Area (A-s): 0.105  
Concentration (ug/L): 46.7

Time: 13:29  
Peak Height (A): 0.205  
Background Pk Height (A): 0.047

QC sample is within range 45 - 55

~~~~~  
Pb ID: ICB Seq. No.: 00006 A/S Pos.: 37 Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 37
Replicate 1 (Peak Stored) Time: 13:32
Peak Area (A-s): 0.001 Peak Height (A): 0.003
Background Pk Area (A-s): 0.006 Background Pk Height (A): 0.006
Blank Corrected Pk Area (A-s): -0.000
Concentration (ug/L): -0.1

QC sample is within range -10 - 10

~~~~~  
Pb ID: PB 3020 1/28/93 Seq. No.: 00007 A/S Pos.: 1 Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 1  
Replicate 1 (Peak Stored) Time: 13:36  
Peak Area (A-s): 0.002 Peak Height (A): 0.005  
Background Pk Area (A-s): 0.010 Background Pk Height (A): 0.006  
Blank Corrected Pk Area (A-s): 0.001  
Concentration (ug/L ): 0.2

~~~~~  
Pb ID: PB 3020 1/28/93 Seq. No.: 00008 A/S Pos.: 1 Date: 01/30/93

µL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 1
Replicate 1 (Peak Stored) Time: 13:39
Peak Area (A-s): 0.023 Peak Height (A): 0.039
Background Pk Area (A-s): 0.013 Background Pk Height (A): 0.011
Blank Corrected Pk Area (A-s): 0.022
Concentration (ug/L): 6.8

Recovery is 66.0% (outside of specified limits)

~~~~~  
Pb ID: LCSW 1/28/93 Seq. No.: 00009 A/S Pos.: 2 Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 2  
Replicate 1 (Peak Stored) Time: 13:42  
Peak Area (A-s): 0.046 Peak Height (A): 0.075  
Background Pk Area (A-s): 0.018 Background Pk Height (A): 0.017  
Blank Corrected Pk Area (A-s): 0.045  
Concentration (ug/L ): 16.9

~~~~~  
Pb ID: LCSW 1/28/93 Seq. No.: 00010 A/S Pos.: 2 Date: 01/30/93

µL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 2
Replicate 1 (Peak Stored) Time: 13:46
Peak Area (A-s): 0.070 Peak Height (A): 0.121
Background Pk Area (A-s): 0.022 Background Pk Height (A): 0.027
Blank Corrected Pk Area (A-s): 0.069
Concentration (ug/L): 28.2

Recovery is 113.4%

~~~~~  
Pb ID: 0059-06 Seq. No.: 00011 A/S Pos.: 3 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 3  
Replicate 1 (Peak Stored) Time: 13:49  
Peak Area (A-s): 0.006 Peak Height (A): 0.014  
Background Pk Area (A-s): 0.040 Background Pk Height (A): 0.031  
Blank Corrected Pk Area (A-s): 0.005  
Concentration (ug/L ): 1.3

~~~~~  
Pb ID: 0059-06 Seq. No.: 00012 A/S Pos.: 3 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 3
Replicate 1 (Peak Stored) Time: 13:52
Peak Area (A-s): 0.029 Peak Height (A): 0.064
Background Pk Area (A-s): 0.038 Background Pk Height (A): 0.032
Blank Corrected Pk Area (A-s): 0.028
Concentration (ug/L): 9.2

Recovery is 78.6% (outside of specified limits)

~~~~~  
P ID: 0059-06REP Seq. No.: 00013 A/S Pos.: 4 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 4  
Replicate 1 (Peak Stored) Time: 13:55  
Peak Area (A-s): 0.006 Peak Height (A): 0.015  
Background Pk Area (A-s): 0.034 Background Pk Height (A): 0.025  
Blank Corrected Pk Area (A-s): 0.005  
Concentration (ug/L ): 1.4

~~~~~  
P ID: 0059-06REP Seq. No.: 00014 A/S Pos.: 4 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 4
Replicate 1 (Peak Stored) Time: 13:58
Peak Area (A-s): 0.072 Peak Height (A): 0.167
Background Pk Area (A-s): 0.045 Background Pk Height (A): 0.053
Blank Corrected Pk Area (A-s): 0.071
Concentration (ug/L): 29.5

Recovery is 281.5% (outside of specified limits)

~~~~~  
Pb ID: 0059-06 SPK Seq. No.: 00015 A/S Pos.: 5 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 5  
Replicate 1 (Peak Stored) Time: 14:01  
Peak Area (A-s): 0.053 Peak Height (A): 0.120  
Background Pk Area (A-s): 0.042 Background Pk Height (A): 0.049  
Blank Corrected Pk Area (A-s): 0.052  
Concentration (ug/L ): 20.0

~~~~~  
Pb ID: 0059-06 SPK Seq. No.: 00016 A/S Pos.: 5 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 5
Replicate 1 (Peak Stored) Time: 14:05
Peak Area (A-s): 0.072 Peak Height (A): 0.182
Background Pk Area (A-s): 0.047 Background Pk Height (A): 0.058
Blank Corrected Pk Area (A-s): 0.071

Concentration (ug/L): 29.3

Recovery is 92.4%

~~~~~  
Pb ID: CCV Seq. No.: 00017 A/S Pos.: 35 Date: 01/30/93

uL dispensed: 10 from 0. 5 from 38. 20 from 35

Replicate 1 (Peak Stored)

Time: 14:08

Peak Area (A-s): 0.108

Peak Height (A): 0.265

Background Pk Area (A-s): 0.036

Background Pk Height (A): 0.056

Blank Corrected Pk Area (A-s): 0.107

Concentration (ug/L ): 48.1

QC sample is within range 45 - 55

~~~~~  
Pb ID: CCB Seq. No.: 00018 A/S Pos.: 37 Date: 01/30/93

uL dispensed: 10 from 0. 5 from 38. 20 from 37

Replicate 1 (Peak Stored)

Time: 14:11

Peak Area (A-s): 0.002

Peak Height (A): 0.005

Background Pk Area (A-s): 0.013

Background Pk Height (A): 0.016

Blank Corrected Pk Area (A-s): 0.001

Concentration (ug/L): 0.3

QC sample is within range -10 - 10

~~~~~  
Pb ID: 0059-07 Seq. No.: 00019 A/S Pos.: 6 Date: 01/30/93

uL dispensed: 10 from 0. 5 from 38. 20 from 6

Replicate 1 (Peak Stored)

Time: 14:14

Peak Area (A-s): 0.005

Peak Height (A): 0.012

Background Pk Area (A-s): 0.032

Background Pk Height (A): 0.020

Blank Corrected Pk Area (A-s): 0.004

Concentration (ug/L ): 0.9

~~~~~  
Pb ID: 0059-07 Seq. No.: 00020 A/S Pos.: 6 Date: 01/30/93

uL dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 6

Replicate 1 (Peak Stored)

Time: 14:17

Peak Area (A-s): 0.027

Peak Height (A): 0.072

Background Pk Area (A-s): 0.039

Background Pk Height (A): 0.034

Blank Corrected Pk Area (A-s): 0.026

Concentration (ug/L): 8.4

Recovery is 74.2% (outside of specified limits)

~~~~~  
Pb ID: 0059-08 Seq. No.: 00021 A/S Pos.: 7 Date: 01/30/93

uL dispensed: 10 from 0. 5 from 38. 20 from 7

Replicate 1 (Peak Stored)

Time: 14:20

Peak Area (A-s): 0.005

Peak Height (A): 0.012

Background Pk Area (A-s): 0.038

Background Pk Height (A): 0.026

Blank Corrected Pk Area (A-s): 0.004

Concentration (ug/L ): 1.0

~~~~~  
Pb ID: 0059-08 Seq. No.: 00022 A/S Pos.: 7 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 7
Replicate 1 (Peak Stored) Time: 14:23
Peak Area (A-s): 0.025 Peak Height (A): 0.069
Background Pk Area (A-s): 0.042 Background Pk Height (A): 0.036
Blank Corrected Pk Area (A-s): 0.024
Concentration (ug/L): 7.8

Recovery is 68.7% (outside of specified limits)

~~~~~  
Pb ID: 0059-09 Seq. No.: 00023 A/S Pos.: 8 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 8  
Replicate 1 (Peak Stored) Time: 14:26  
Peak Area (A-s): 0.006 Peak Height (A): 0.019  
Background Pk Area (A-s): 0.040 Background Pk Height (A): 0.029  
Blank Corrected Pk Area (A-s): 0.005  
Concentration (ug/L ): 1.4

~~~~~  
Pb ID: 0059-09 Seq. No.: 00024 A/S Pos.: 8 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 8
Replicate 1 (Peak Stored) Time: 14:30
Peak Area (A-s): 0.029 Peak Height (A): 0.069
Background Pk Area (A-s): 0.044 Background Pk Height (A): 0.035
Blank Corrected Pk Area (A-s): 0.028
Concentration (ug/L): 9.4

Recovery is 80.7% (outside of specified limits)

~~~~~  
Pb ID: 0059-10 Seq. No.: 00025 A/S Pos.: 9 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 9  
Replicate 1 (Peak Stored) Time: 14:33  
Peak Area (A-s): 0.009 Peak Height (A): 0.020  
Background Pk Area (A-s): 0.038 Background Pk Height (A): 0.032  
Blank Corrected Pk Area (A-s): 0.008  
Concentration (ug/L ): 2.0

~~~~~  
Pb ID: 0059-10 Seq. No.: 00026 A/S Pos.: 9 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 9
Replicate 1 (Peak Stored) Time: 14:36
Peak Area (A-s): 0.033 Peak Height (A): 0.077
Background Pk Area (A-s): 0.039 Background Pk Height (A): 0.033
Blank Corrected Pk Area (A-s): 0.032
Concentration (ug/L): 11.2

Recovery is 91.6%

~~~~~  
Pb ID: ~~PBS 1726793~~ Seq. No.: 00027 A/S Pos.: 10 Date: 01/30/93

2/1/93  
CA.

uL dispensed: 10 from 0, 5 from 38, 20 from 10  
Replicate 1 (Peak Stored) Time: 14:39  
Peak Area (A-s): 0.005 Peak Height (A): 0.011  
Background Pk Area (A-s): 0.014 Background Pk Height (A): 0.016  
Blank Corrected Pk Area (A-s): 0.005  
Concentration (ug/L ): 1.1

Pb ID: PBS 1/26/93 Seq. No.: 00028 A/S Pos.: 10 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 10  
Replicate 1 (Peak Stored) Time: 14:42  
Peak Area (A-s): 0.028 Peak Height (A): 0.052  
Background Pk Area (A-s): 0.020 Background Pk Height (A): 0.019  
Blank Corrected Pk Area (A-s): 0.027  
Concentration (ug/L ): 8.8

Recovery is 76.9% (outside of specified limits)

Pb ID: CCV Seq. No.: 00029 A/S Pos.: 35 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 35  
Replicate 1 (Peak Stored) Time: 14:45  
Peak Area (A-s): 0.111 Peak Height (A): 0.252  
Background Pk Area (A-s): 0.034 Background Pk Height (A): 0.050  
Blank Corrected Pk Area (A-s): 0.110  
Concentration (ug/L ): 49.4

QC sample is within range 45 - 55

Pb ID: CCB Seq. No.: 00030 A/S Pos.: 37 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 37  
Replicate 1 (Peak Stored) Time: 14:48  
Peak Area (A-s): 0.003 Peak Height (A): 0.006  
Background Pk Area (A-s): 0.016 Background Pk Height (A): 0.018  
Blank Corrected Pk Area (A-s): 0.002  
Concentration (ug/L ): 0.4

QC sample is within range -10 - 10

Pb ID: LCSS 1/26 1-10 Seq. No.: 00031 A/S Pos.: 11 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 11  
Replicate 1 (Peak Stored) Time: 14:51  
Peak Area (A-s): 0.054 Peak Height (A): 0.102  
Background Pk Area (A-s): 0.029 Background Pk Height (A): 0.027  
Blank Corrected Pk Area (A-s): 0.053  
Concentration (ug/L ): 20.4 Corrected Conc (ug/L ): 204.

40.8 mg/kg 91.7% Acc.

Pb ID: LCSS 1/26 1-10 Seq. No.: 00032 A/S Pos.: 11 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 11  
Replicate 1 (Peak Stored) Time: 14:55  
Peak Area (A-s): 0.075 Peak Height (A): 0.137

Background Pk Area (A-s): 0.034  
Blank Corrected Pk Area (A-s): 0.074  
Concentration (ug/L ): 30.9

Background Pk Height (A): 0.036  
Corrected Conc (ug/L ): 309.

Recovery is 104.6%

~~~~~  
P1 ID: 0104-01 1-100 Seq. No.: 00033 A/S Pos.: 12 Date: 01/30/93

UL dispensed: 10 from 0, 5 from 38, 20 from 12
Replicate 1 (Peak Stored) Time: 15:00
Peak Area (A-s): 0.026 Peak Height (A): 0.043
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.019
Blank Corrected Pk Area (A-s): 0.025
Concentration (ug/L): 8.3 Corrected Conc (ug/L): 830.

~~~~~  
P1 ID: 0104-01 1-100 Seq. No.: 00034 A/S Pos.: 12 Date: 01/30/93

UL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 12  
Replicate 1 (Peak Stored) Time: 15:03  
Peak Area (A-s): 0.042 Peak Height (A): 0.091  
Background Pk Area (A-s): 0.028 Background Pk Height (A): 0.024  
Blank Corrected Pk Area (A-s): 0.041  
Concentration (ug/L ): 14.9 Corrected Conc (ug/L ): 1489.

Recovery is 65.9% (outside of specified limits)

~~~~~  
Pb ID: 0104-01 REP 100 Seq. No.: 00035 A/S Pos.: 13 Date: 01/30/93

UL dispensed: 10 from 0, 5 from 38, 20 from 13
Replicate 1 (Peak Stored) Time: 15:06
Peak Area (A-s): 0.029 Peak Height (A): 0.051
Background Pk Area (A-s): 0.022 Background Pk Height (A): 0.019
Blank Corrected Pk Area (A-s): 0.028
Concentration (ug/L): 9.2 Corrected Conc (ug/L): 925.

2/1/93
E.H.

~~~~~  
Pb ID: 0104-01 REP 100 Seq. No.: 00036 A/S Pos.: 13 Date: 01/30/93

UL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 13  
Replicate 1 (Peak Stored) Time: 15:10  
Peak Area (A-s): 0.048 Peak Height (A): 0.080  
Background Pk Area (A-s): 0.028 Background Pk Height (A): 0.022  
Blank Corrected Pk Area (A-s): 0.048  
Concentration (ug/L ): 17.9 Corrected Conc (ug/L ): 1791.

Recovery is 86.6%

~~~~~  
P ID: 0104-01 SPK-100 Seq. No.: 00037 A/S Pos.: 14 Date: 01/30/93

UL dispensed: 10 from 0, 5 from 38, 20 from 14
Replicate 1 (Peak Stored) Time: 15:13
Peak Area (A-s): 0.040 Peak Height (A): 0.076
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.024
Blank Corrected Pk Area (A-s): 0.039
Concentration (ug/L): 14.1 Corrected Conc (ug/L): 1413.

2/1/93
E.H.

1/6.9% Rec.

~~~~~  
Pb ID: 0104-01 SPK 100 Seq. No.: 00038 A/S Pos.: 14 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 14  
Replicate 1 (Peak Stored) Time: 15:16  
Peak Area (A-s): 0.056 Peak Height (A): 0.100  
Background Pk Area (A-s): 0.028 Background Pk Height (A): 0.024  
Blank Corrected Pk Area (A-s): 0.055  
Concentration (ug/L ): 21.4 Corrected Conc (ug/L ): 2142.

Recovery is 72.9% (outside of specified limits)

~~~~~  
Pb ID: 0104-02 100 Seq. No.: 00039 A/S Pos.: 15 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 15
Replicate 1 (Peak Stored) Time: 15:19
Peak Area (A-s): 0.005 Peak Height (A): 0.011
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.023
Blank Corrected Pk Area (A-s): 0.004
Concentration (ug/L): 1.1 Corrected Conc (ug/L): 106.

~~~~~  
Pb ID: 0104-02 100 Seq. No.: 00040 A/S Pos.: 15 Date: 01/30/93

uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 15  
Replicate 1 (Peak Stored) Time: 15:22  
Peak Area (A-s): 0.026 Peak Height (A): 0.048  
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.026  
Blank Corrected Pk Area (A-s): 0.025  
Concentration (ug/L ): 8.2 Corrected Conc (ug/L ): 822.

Recovery is 71.6% (outside of specified limits)

~~~~~  
Pb ID: CCV Seq. No.: 00041 A/S Pos.: 35 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 35
Replicate 1 (Peak Stored) Time: 15:25
Peak Area (A-s): 0.119 Peak Height (A): 0.255
Background Pk Area (A-s): 0.034 Background Pk Height (A): 0.051
Blank Corrected Pk Area (A-s): 0.118
Concentration (ug/L): 53.7

Q sample is within range 45 - 55

~~~~~  
Pb ID: CCB Seq. No.: 00042 A/S Pos.: 37 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 37  
Replicate 1 (Peak Stored) Time: 15:28  
Peak Area (A-s): 0.002 Peak Height (A): 0.005  
Background Pk Area (A-s): 0.017 Background Pk Height (A): 0.017  
Blank Corrected Pk Area (A-s): 0.001  
Concentration (ug/L ): 0.3

Q sample is within range -10 - 10

~~~~~  
Pl ID: 0104-03 1-100 Seq. No.: 00043 A/S Pos.: 16 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 16
Replicate 1 (Peak Stored) Time: 15:32
Peak Area (A-s): 0.011 Peak Height (A): 0.020
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.021
Blank Corrected Pk Area (A-s): 0.010
Concentration (ug/L): 2.6 Corrected Conc (ug/L): 259.

~~~~~  
Pl ID: 0104-03 1-100 Seq. No.: 00044 A/S Pos.: 16 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 16  
Replicate 1 (Peak Stored) Time: 15:35  
Peak Area (A-s): 0.031 Peak Height (A): 0.051  
Background Pk Area (A-s): 0.030 Background Pk Height (A): 0.022  
Blank Corrected Pk Area (A-s): 0.030  
Concentration (ug/L ): 10.3 Corrected Conc (ug/L ): 1032.

Recovery is 77.3% (outside of specified limits)

~~~~~  
Pl ID: 0104-04 1-100 Seq. No.: 00045 A/S Pos.: 17 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 17
Replicate 1 (Peak Stored) Time: 15:38
Peak Area (A-s): 0.006 Peak Height (A): 0.011
Background Pk Area (A-s): 0.024 Background Pk Height (A): 0.023
Blank Corrected Pk Area (A-s): 0.005
Concentration (ug/L): 1.2 Corrected Conc (ug/L): 122.

~~~~~  
Pl ID: 0104-04 1-100 Seq. No.: 00046 A/S Pos.: 17 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 17  
Replicate 1 (Peak Stored) Time: 15:41  
Peak Area (A-s): 0.028 Peak Height (A): 0.047  
Background Pk Area (A-s): 0.029 Background Pk Height (A): 0.021  
Blank Corrected Pk Area (A-s): 0.027  
Concentration (ug/L ): 8.9 Corrected Conc (ug/L ): 889.

Recovery is 76.7% (outside of specified limits)

~~~~~  
Pl ID: 0104-05 1-100 Seq. No.: 00047 A/S Pos.: 18 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 18
Replicate 1 (Peak Stored) Time: 15:44
Peak Area (A-s): 0.021 Peak Height (A): 0.043
Background Pk Area (A-s): 0.027 Background Pk Height (A): 0.021
Blank Corrected Pk Area (A-s): 0.020
Concentration (ug/L): 6.3 Corrected Conc (ug/L): 632.

~~~~~  
Pl ID: 0104-05 1-100 Seq. No.: 00048 A/S Pos.: 18 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 18  
Replicate 1 (Peak Stored) Time: 15:48

Peak Area (A-s): 0.043  
Background Pk Area (A-s): 0.032  
Blank Corrected Pk Area (A-s): 0.042  
Concentration (ug/L ): 15.2

Peak Height (A): 0.072  
Background Pk Height (A): 0.024  
Corrected Conc (ug/L ): 1520.

Recovery is 88.8%

~~~~~  
Pb ID: CCV Seq. No.: 00049 A/S Pos.: 35 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 35

~~~~~  
Pb ID: CCV Seq. No.: 00050 A/S Pos.: 35 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 35

Replicate 1 (Peak Stored) Time: 15:52  
Peak Area (A-s): 0.206 Peak-Height (A): 0.428  
Background Pk Area (A-s): 0.080 Background Pk Height (A): 0.104  
Blank Corrected Pk Area (A-s): 0.205  
Concentration (ug/L ): 101.6

*double injected w/o run*  
*TV = 100 ug/L - in range.*

~~QC sample is out of range 45 - 55~~

~~~~~  
Pb ID: CCB Seq. No.: 00051 A/S Pos.: 37 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 37

~~~~~  
Pb ID: PB 3020 1/28/93 Seq. No.: 00052 A/S Pos.: 1 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 1

Replicate 1 (Peak Stored) Time: 15:56  
Peak Area (A-s): 0.003 Peak Height (A): 0.008  
Background Pk Area (A-s): 0.015 Background Pk Height (A): 0.014  
Blank Corrected Pk Area (A-s): 0.002  
Concentration (ug/L ): 0.6

~~~~~  
Pb ID: PB 3020 1/28/93 Seq. No.: 00053 A/S Pos.: 1 Date: 01/30/93

ul dispensed: 5 from 0. 5 from 38. 5 from 40. 20 from 1

Replicate 1 (Peak Stored) Time: 15:59
Peak Area (A-s): 0.025 Peak Height (A): 0.054
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.017
Blank Corrected Pk Area (A-s): 0.024
Concentration (ug/L): 7.9

Recovery is 73.4% (outside of specified limits)

~~~~~  
Pb ID: LCSW 1/28/93 Seq. No.: 00054 A/S Pos.: 2 Date: 01/30/93

ul dispensed: 10 from 0. 5 from 38. 20 from 2

Replicate 1 (Peak Stored) Time: 16:02  
Peak Area (A-s): 0.047 Peak Height (A): 0.125  
Background Pk Area (A-s): 0.029 Background Pk Height (A): 0.036  
Blank Corrected Pk Area (A-s): 0.046  
Concentration (ug/L ): 17.4

ID: LCSW 1/28/93

Seq. No.: 00055

A/S Pos.: 2

Date: 01/30/93

µL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 2

Replicate 1 (Peak Stored)

Time: 16:05

Peak Area (A-s): 0.067

Peak Height (A): 0.151

Background Pk Area (A-s): 0.036

Background Pk Height (A): 0.043

Blank Corrected Pk Area (A-s): 0.066

Concentration (ug/L ): 26.6

Recovery is 91.8%

ID: 0104-01 1-100

Seq. No.: 00056

A/S Pos.: 12

Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 12

Replicate 1 (Peak Stored)

Time: 16:08

Peak Area (A-s): 0.021

Peak Height (A): 0.042

Background Pk Area (A-s): 0.020

Background Pk Height (A): 0.015

Blank Corrected Pk Area (A-s): 0.020

Concentration (ug/L ): 6.3

Corrected Conc (ug/L ): 632.

ID: 0104-01 1-100

Seq. No.: 00057

A/S Pos.: 12

Date: 01/30/93

µL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 12

Replicate 1 (Peak Stored)

Time: 16:12

Peak Area (A-s): 0.041

Peak Height (A): 0.079

Background Pk Area (A-s): 0.027

Background Pk Height (A): 0.023

Blank Corrected Pk Area (A-s): 0.040

Concentration (ug/L ): 14.5

Corrected Conc (ug/L ): 1453.

Recovery is 82.1% (outside of specified limits)

ID: CCV

Seq. No.: 00058

A/S Pos.: 35

Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 35

Replicate 1 (Peak Stored)

Time: 16:15

Peak Area (A-s): 0.109

Peak Height (A): 0.237

Background Pk Area (A-s): 0.035

Background Pk Height (A): 0.050

Blank Corrected Pk Area (A-s): 0.108

Concentration (ug/L ): 48.4

sample is within range 45 - 55

ID: 0104-01 1-100

Seq. No.: 00059

A/S Pos.: 12

Date: 01/30/93

µL dispensed: 10 from 0, 5 from 38, 20 from 12

Replicate 1 (Peak Stored)

Time: 16:18

Peak Area (A-s): 0.024

Peak Height (A): 0.045

Background Pk Area (A-s): 0.022

Background Pk Height (A): 0.016

Blank Corrected Pk Area (A-s): 0.023

Concentration (ug/L ): 7.5

Corrected Conc (ug/L ): 745.

ID: 0104-01 1-100

Seq. No.: 00060

A/S Pos.: 12

Date: 01/30/93



uL dispensed: 5 from 0, 5 from 38, 5 from 40, 20 from 12  
Replicate 1 (Peak Stored) Time: 16:21  
Peak Area (A-s): 0.046 Peak Height (A): 0.085  
Background Pk Area (A-s): 0.028 Background Pk Height (A): 0.023  
Blank Corrected Pk Area (A-s): 0.045  
Concentration (ug/L ): 16.6 Corrected Conc (ug/L ): 1663.

Recovery is 91.8%

~~~~~  
Pb ID: CCV Seq. No.: 00061 A/S Pos.: 35 Date: 01/30/93

uL dispensed: 10 from 0, 5 from 38, 20 from 35
Replicate 1 (Peak Stored) Time: 16:24
Peak Area (A-s): 0.110 Peak Height (A): 0.229
Background Pk Area (A-s): 0.034 Background Pk Height (A): 0.046
Blank Corrected Pk Area (A-s): 0.110
Concentration (ug/L): 49.2

QC sample is within range 45 - 55

ID/Weight File: 0059.IDW
Sample Volume: 100 mL

Analyst: ARSENIC
Nominal Weight: 1.0 g

Loc.	Sample ID	Weight	Dilution
1	PBS 1/26/93 ✓		
2	LCSS 1/26 1-10 ✓		10
3	0104-01 1-10 ✓		10
4	0104-01 REP 10 ✓		10
5	0104-01 SPK 50 ✓		50
6	0104-02 10 ✓		10
7	0104-03 1-10 ✓		10
8	0104-04 1-10		10
9	0104-05 1-10		10

Element File: ASBRIEF.GEL

Element: As

Analyst:

Print Data: Main+Suppl.

Peak Storage: 1 Repl./Sample

Remarks:

Matrix Modifier: Palladium/ Magnesium Int time = 3 sec

Lamp:

Energy:

Current

Char. Mass: 20 ppb

0 ppb: 100 ppb:

INSTRUMENT: 4100 ZL

Technique: HGA

Wavelength: 193.7 Peak

Slit: 0.7 Low

Signal Type: Zeeman AA

Signal Measurement: Peak Area

Read Time: 3.0

Read Delay: 0.0

BOC Time: 2

Sample Replicates: 1

Standard Replicates: 1

Spike Replicates: Same as Sample

CALIBRATION:

Solutions	ID	Conc	Location	Volume	Diluent	Modifier	
						#1	#2
Calib. Blank	CAL BLANK		37	20	10		5
Standard 1	20 ppb	20.0	36	4	26		5
Standard 2	60 ppb	60.0	36	12	18		5
Standard 3	100 ppb	100.0	36	20	10		5
Samples				20	10		5

Diluent Location: 0

Modifier #1 Location:

Modifier #2 Location: 39

Calibration Units: ug/L

Sample Units: ug/L

Calibration Type: Nonlinear

Furnace Time/Temperature Program:

Step|Temp Ramp Hold|Gas Flow|Read|Gas Type|

1 | 110 1 20 | 250 | | Norm |

2 | 130 5 30 | 250 | | Norm |

3 | 1000 10 20 | 250 | | Alt |

4 | 2000 0 5 | 0 | * | Norm |

5 | 2400 1 2 | 250 | | Norm |

Injection Temp: 20

Pipette Speed: 100%

Extraction System: On

SEQUENCE:

Step Action and Parameters

1 Pipet diluent + modifier 2 + spike + sample/std

2 Run HGA steps 1 to End

CHECKS:

Recalibration Type: Completely Recalibrate

Locations: None

Conc. Above Calibration Action: Dilute & Reanalyze After All Reps

Alternate Sample Volumes (uL): 10,5,4

QC Alternate Volume Blanks: No

If %RSD > 20.0 and Concentration > .05 then Retry 1 times

Check %RSD on: Samples + Standards

Recovery Measurements:

10 uL of 40 ug/L Standard at Location 40 Gives 20.0 ug/L

Measure Recovery on Samples: All Samples

Add to QC Samples: No % Recovery Limits: 85 to 115

QC

#	A/S	QC Sample	Conc. Limits	After	Periodic	At	Count As
	Loc.	ID	Lower	Upper	Calib	Check	End Sample
1	35	ICV	45.0	55.0	X		
2	37	ICB	-10.0	10.0			
3	38	CRA	1.0	20.0			
4	35	CCV	45.0	55.0		X	X
5	37	CCB	-10.0	10.0			

Run Periodic QC Samples: Every 10

Out of Limit Action: Stop Element and Continue With Next

Matrix Check Calculations:

% Difference for Dupls: No

Locations:

% Recovery for Spike: No

Locations:

Conc:

Element File: ASBRIEF.GEL Element: As Wavelength: 193.7
Date: 02/02/93 Time: 09:26 Slit: 0.7 L
Data File: F320231.DAT ID/Wt File: 0059.IDW Lamp Current: 0
Technique: HGA Calib. Type: Nonlinear Energy: 50
Remark 1: Matrix Modifier: Palladium/ Magnesium Int time = 3 sec
Remark 2: Lamp:
Remark 3: Energy: Current
Remark 4: Char. Mass: 20 ppb
Remark 5: 60 ppb: 100 ppb:

~~~~~  
A      ID: CAL BLANK      Seq. No.: 00001      A/S Pos.: 37      Date: 02/02/93

Replicate 1 (Peak Stored)      Time: 09:26  
Peak Area (A-s): 0.005      Peak Height (A): 0.007  
Background Pk Area (A-s): 0.027      Background Pk Height (A): 0.032  
Blank Corrected Pk Area (A-s): 0.005

Auto-zero performed.

~~~~~  
A ID: 20 ppb Seq. No.: 00002 A/S Pos.: 36 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 09:28
Peak Area (A-s): 0.041 Peak Height (A): 0.127
Background Pk Area (A-s): 0.023 Background Pk Height (A): 0.026
Blank Corrected Pk Area (A-s): 0.036

Standard number 1 applied. [20.0]
Correlation coefficient: 1.00000 Slope: 0.0018

~~~~~  
As      ID: 60 ppb      Seq. No.: 00003      A/S Pos.: 36      Date: 02/02/93

Replicate 1 (Peak Stored)      Time: 09:31  
Peak Area (A-s): 0.114      Peak Height (A): 0.324  
Background Pk Area (A-s): 0.030      Background Pk Height (A): 0.029  
Blank Corrected Pk Area (A-s): 0.110  
Concentration (ug/L ): 60.2

Standard number 2 applied. [60.0]  
Correlation coefficient: 1.00000      Slope: 0.0018

~~~~~  
A ID: 100 ppb Seq. No.: 00004 A/S Pos.: 36 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 09:34
Peak Area (A-s): 0.191 Peak Height (A): 0.544
Background Pk Area (A-s): 0.031 Background Pk Height (A): 0.037
Blank Corrected Pk Area (A-s): 0.186
Concentration (ug/L): 101.5

Standard number 3 applied. [100.0]
Correlation coefficient: 1.00000 Slope: 0.0018

~~~~~  
As      ID: ICV      Seq. No.: 00005      A/S Pos.: 35      Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.099  
Background Pk Area (A-s): 0.025  
Blank Corrected Pk Area (A-s): 0.094  
Concentration (ug/L ): 51.5

Time: 09:37  
Peak Height (A): 0.315  
Background Pk Height (A): 0.029

QC sample is within range 45.0 - 55.0

As ID: PBS 1/26/93 Seq. No.: 00006 A/S Pos.: 1 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.001  
Background Pk Area (A-s): 0.017  
Blank Corrected Pk Area (A-s): -0.003  
Concentration (ug/L ): -1.8

Time: 09:39  
Peak Height (A): 0.009  
Background Pk Height (A): 0.024

As ID: PBS 1/26/93 Seq. No.: 00007 A/S Pos.: 1 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.041  
Background Pk Area (A-s): 0.013  
Blank Corrected Pk Area (A-s): 0.036  
Concentration (ug/L ): 19.9

Time: 09:42  
Peak Height (A): 0.128  
Background Pk Height (A): 0.022

Recovery is 108.6%

As ID: LCSS 1/26 1-10 Seq. No.: 00008 A/S Pos.: 2 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.031  
Background Pk Area (A-s): 0.032  
Blank Corrected Pk Area (A-s): 0.026  
Concentration (ug/L ): 14.4

Time: 09:45  
Peak Height (A): 0.074  
Background Pk Height (A): 0.032  
Corrected Conc (ug/L ): 144.

*28.8 mg/Kg 84.29 ug.*

As ID: LCSS 1/26 1-10 Seq. No.: 00009 A/S Pos.: 2 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.069  
Background Pk Area (A-s): 0.024  
Blank Corrected Pk Area (A-s): 0.064  
Concentration (ug/L ): 35.4

Time: 09:47  
Peak Height (A): 0.191  
Background Pk Height (A): 0.022  
Corrected Conc (ug/L ): 354.

Recovery is 104.7%

As ID: 0104-01 1-10 Seq. No.: 00010 A/S Pos.: 3 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.022  
Background Pk Area (A-s): 0.043  
Blank Corrected Pk Area (A-s): 0.017  
Concentration (ug/L ): 9.5

Time: 09:50  
Peak Height (A): 0.036  
Background Pk Height (A): 0.034  
Corrected Conc (ug/L ): 95.

As ID: 0104-01 1-10 Seq. No.: 00011 A/S Pos.: 3 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 09:53  
Peak Area (A-s): 0.058 Peak Height (A): 0.092  
Background Pk Area (A-s): 0.039 Background Pk Height (A): 0.033  
Blank Corrected Pk Area (A-s): 0.053  
Concentration (ug/L ): 29.0 Corrected Conc (ug/L ): 290.

Recovery is 97.8%

As ID: 0104-01 REP 10 Seq. No.: 00012 A/S Pos.: 4 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 09:55  
Peak Area (A-s): 0.029 Peak Height (A): 0.044  
Background Pk Area (A-s): 0.040 Background Pk Height (A): 0.036  
Blank Corrected Pk Area (A-s): 0.024  
Concentration (ug/L ): 13.2 Corrected Conc (ug/L ): 132.

As ID: 0104-01 REP 10 Seq. No.: 00013 A/S Pos.: 4 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 09:58  
Peak Area (A-s): 0.069 Peak Height (A): 0.096  
Background Pk Area (A-s): 0.045 Background Pk Height (A): 0.038  
Blank Corrected Pk Area (A-s): 0.065  
Concentration (ug/L ): 35.6 Corrected Conc (ug/L ): 356.

Recovery is 112.0%

As ID: 0104-01 SPK 50 Seq. No.: 00014 A/S Pos.: 5 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:01  
Peak Area (A-s): 0.075 Peak Height (A): 0.194  
Background Pk Area (A-s): 0.027 Background Pk Height (A): 0.029  
Blank Corrected Pk Area (A-s): 0.071  
Concentration (ug/L ): ~~38.7~~ Corrected Conc (ug/L ): 1936.

96.8% Rec.

As ID: 0104-01 SPK 50 Seq. No.: 00015 A/S Pos.: 5 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:04  
Peak Area (A-s): 0.121 Peak Height (A): 0.224  
Background Pk Area (A-s): 0.033 Background Pk Height (A): 0.031  
Blank Corrected Pk Area (A-s): 0.116  
Concentration (ug/L ): 63.5 Corrected Conc (ug/L ): 3176.

Recovery is 124.0% (outside of specified limits)

As ID: CCV Seq. No.: 00016 A/S Pos.: 35 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:06  
Peak Area (A-s): 0.085 Peak Height (A): 0.123  
Background Pk Area (A-s): 0.101 Background Pk Height (A): 0.215  
Blank Corrected Pk Area (A-s): 0.081  
Concentration (ug/L ): 44.1

QC sample is out of range 45.0 - 55.0

~~~~~  
AL ID: ICV Seq. No.: 00017 A/S Pos.: 35 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:35
Peak Area (A-s): 0.094 Peak Height (A): 0.256
Background Pk Area (A-s): 0.034 Background Pk Height (A): 0.032
Blank Corrected Pk Area (A-s): 0.090
Concentration (ug/L): 49.0

QC sample is within range 45.0 - 55.0

~~~~~  
As ID: 0104-02 10 Seq. No.: 00018 A/S Pos.: 6 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:37  
Peak Area (A-s): 0.006 Peak Height (A): 0.012  
Background Pk Area (A-s): 0.076 Background Pk Height (A): 0.062  
Blank Corrected Pk Area (A-s): 0.001  
Concentration (ug/L ): 0.5 Corrected Conc (ug/L ): 5.

~~~~~  
AL ID: 0104-02 10 Seq. No.: 00019 A/S Pos.: 6 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:40
Peak Area (A-s): 0.047 Peak Height (A): 0.070
Background Pk Area (A-s): 0.068 Background Pk Height (A): 0.065
Blank Corrected Pk Area (A-s): 0.042
Concentration (ug/L): 23.1 Corrected Conc (ug/L): 231.

Recovery is 113.0%

~~~~~  
As ID: 0104-03 1-10 Seq. No.: 00020 A/S Pos.: 7 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:43  
Peak Area (A-s): 0.009 Peak Height (A): 0.014  
Background Pk Area (A-s): 0.074 Background Pk Height (A): 0.051  
Blank Corrected Pk Area (A-s): 0.005  
Concentration (ug/L ): 2.6 Corrected Conc (ug/L ): 26.

~~~~~  
AL ID: 0104-03 1-10 Seq. No.: 00021 A/S Pos.: 7 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:46
Peak Area (A-s): 0.054 Peak Height (A): 0.062
Background Pk Area (A-s): 0.079 Background Pk Height (A): 0.055
Blank Corrected Pk Area (A-s): 0.049
Concentration (ug/L): 27.1 Corrected Conc (ug/L): 271.

Recovery is 122.6% (outside of specified limits)

~~~~~  
As ID: 0104-04 1-10 Seq. No.: 00022 A/S Pos.: 8 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 10:48  
Peak Area (A-s): 0.010 Peak Height (A): 0.017  
Background Pk Area (A-s): 0.076 Background Pk Height (A): 0.062



Blank Corrected Pk Area (A-s): 0.006

Concentration (ug/L ): 3.1

Corrected Conc (ug/L ): 31.

As ID: 0104-04 1-10 Seq. No.: 00023 A/S Pos.: 8 Date: 02/02/93

Replicate 1 (Peak Stored)

Time: 10:51

Peak Area (A-s): 0.055

Peak Height (A): 0.067

Background Pk Area (A-s): 0.082

Background Pk Height (A): 0.069

Blank Corrected Pk Area (A-s): 0.051

Concentration (ug/L ): 27.8

Corrected Conc (ug/L ): 278.

Recovery is 123.6% (outside of specified limits)

As ID: 0104-05 1-10 Seq. No.: 00024 A/S Pos.: 9 Date: 02/02/93

Replicate 1 (Peak Stored)

Time: 10:53

Peak Area (A-s): 0.011

Peak Height (A): 0.016

Background Pk Area (A-s): 0.063

Background Pk Height (A): 0.047

Blank Corrected Pk Area (A-s): 0.006

Concentration (ug/L ): 3.4

Corrected Conc (ug/L ): 34.

As ID: 0104-05 1-10 Seq. No.: 00025 A/S Pos.: 9 Date: 02/02/93

Replicate 1 (Peak Stored)

Time: 10:56

Peak Area (A-s): 0.051

Peak Height (A): 0.060

Background Pk Area (A-s): 0.063

Background Pk Height (A): 0.046

Blank Corrected Pk Area (A-s): 0.046

Concentration (ug/L ): 25.3

Corrected Conc (ug/L ): 253.

Recovery is 109.4%

As ID: CCV Seq. No.: 00026 A/S Pos.: 35 Date: 02/02/93

Replicate 1 (Peak Stored)

Time: 10:59

Peak Area (A-s): 0.101

Peak Height (A): 0.299

Background Pk Area (A-s): 0.023

Background Pk Height (A): 0.026

Blank Corrected Pk Area (A-s): 0.096

Concentration (ug/L ): 52.7

QC sample is within range 45.0 - 55.0

As ID: ICV Seq. No.: 00027 A/S Pos.: 35 Date: 02/02/93

Replicate 1 (Peak Stored)

Time: 11:08

Peak Area (A-s): 0.099

Peak Height (A): 0.272

Background Pk Area (A-s): 0.023

Background Pk Height (A): 0.026

Blank Corrected Pk Area (A-s): 0.095

Concentration (ug/L ): 52.0

QC sample is within range 45.0 - 55.0

As ID: IDL Seq. No.: 00028 A/S Pos.: 10 Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.020  
Background Pk Area (A-s): 0.019  
Blank Corrected Pk Area (A-s): 0.015  
Concentration (ug/L ): 8.4

Time: 11:11  
Peak Height (A): 0.061  
Background Pk Height (A): 0.025

As ID: IDL

Seq. No.: 00029

A/S Pos.: 11

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.021  
Background Pk Area (A-s): 0.019  
Blank Corrected Pk Area (A-s): 0.016  
Concentration (ug/L ): 8.8

Time: 11:13  
Peak Height (A): 0.071  
Background Pk Height (A): 0.023

As ID: IDL

Seq. No.: 00030

A/S Pos.: 12

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.021  
Background Pk Area (A-s): 0.020  
Blank Corrected Pk Area (A-s): 0.016  
Concentration (ug/L ): 8.9

Time: 11:16  
Peak Height (A): 0.063  
Background Pk Height (A): 0.023

As ID: IDL

Seq. No.: 00031

A/S Pos.: 13

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.020  
Background Pk Area (A-s): 0.023  
Blank Corrected Pk Area (A-s): 0.016  
Concentration (ug/L ): 8.6

Time: 11:18  
Peak Height (A): 0.069  
Background Pk Height (A): 0.025

As ID: IDL

Seq. No.: 00032

A/S Pos.: 14

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.021  
Background Pk Area (A-s): 0.037  
Blank Corrected Pk Area (A-s): 0.016  
Concentration (ug/L ): 9.0

Time: 11:21  
Peak Height (A): 0.053  
Background Pk Height (A): 0.046

As ID: IDL

Seq. No.: 00033

A/S Pos.: 15

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.023  
Background Pk Area (A-s): 0.028  
Blank Corrected Pk Area (A-s): 0.018  
Concentration (ug/L ): 9.9

Time: 11:23  
Peak Height (A): 0.062  
Background Pk Height (A): 0.029

As ID: IDL

Seq. No.: 00034

A/S Pos.: 16

Date: 02/02/93

Replicate 1 (Peak Stored)  
Peak Area (A-s): 0.020  
Background Pk Area (A-s): 0.036  
Blank Corrected Pk Area (A-s): 0.015  
Concentration (ug/L ): 8.5

Time: 11:26  
Peak Height (A): 0.056  
Background Pk Height (A): 0.036

~~~~~  
As ID: CCV Seq. No.: 00035 A/S Pos.: 35 Date: 02/02/93

Replicate 1 (Peak Stored) Time: 11:29
Peak Area (A-s): 0.091 Peak Height (A): 0.263
Background Pk Area (A-s): 0.032 Background Pk Height (A): 0.039
Blank Corrected Pk Area (A-s): 0.086
Concentration (ug/L): 47.1

QC sample is within range 45.0 - 55.0

$\bar{X} = 8.81$
 $SD = .5023$
 $IDL = 1.501$

P.E. AA GRAPHITE FURNACE 5100Z RUN LOG

P.E. AA GRAPHITE FURNACE 5100Z RUN LOG

Date: 2/1/93

Element: Se

Date Storage File:
F220131

Lamp: P.E EDL

Energy:

Current:

Char. Mass @ 20 ppb:
60 ppb:
100 ppb:

Temp. Program:

POS

SAMPLE ID

1	PBS 1/26/93	21
2	LCSS 1/26/93 1-10	22
3	0104-01	23
4	0104-01 REP	24
5	0104-01 SPK	25
6	0104-02	26
7	0104-03	27
8	0104-04	28
9	0104-05	29
10		30
11		31
12		32
13		33
14		34
15		35
16		36 100 ppb standard
17		37 Blank
18		38 NH4PO4 Mg(NO3)2 modifier
19		39 Pd/Mg
20		40 40 ppb spk sol

Element File: SEBRIEF.GEL Element: Se
Date: 01/29/93 Time: 13:57
Data Storage File: F220131.DAT ID/Weight File: 0104.A60
Technique: HGA Calibration Type: Nonlinear
Remark 1: Analysis was run with Palladium/Magnesium modifier.
Remark 2: Integration time = 5 seconds.

ID: Blank Seq. No.: 00001 A/S Pos.: 37 Date: 01/29/93

Replicate 1 Time: 13:57
Peak Area (A-s): 0.002 Peak Height (A): 0.005
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.009
Blank Corrected Pk Area (A-s): 0.000
Concentration (ug/L): 0.1

*all analysis performed.
2/1/93. GL*

auto-zero performed.

ID: Standard 1 Seq. No.: 00002 A/S Pos.: 36 Date: 01/29/93

Replicate 1 Time: 14:00
Peak Area (A-s): 0.053 Peak Height (A): 0.112
Background Pk Area (A-s): 0.039 Background Pk Height (A): 0.053
Blank Corrected Pk Area (A-s): 0.050
Concentration (ug/L): 17.5

Standard number 1 applied. [20.0]
Correlation coefficient: 1.00000 Slope: 0.0025

ID: Standard 2 Seq. No.: 00003 A/S Pos.: 36 Date: 01/29/93

Replicate 1 Time: 14:03
Peak Area (A-s): 0.150 Peak Height (A): 0.291
Background Pk Area (A-s): 0.084 Background Pk Height (A): 0.130
Blank Corrected Pk Area (A-s): 0.147
Concentration (ug/L): 58.6

Standard number 2 applied. [60.0]
Correlation coefficient: 1.00000 Slope: 0.0025

ID: Standard 3 Seq. No.: 00004 A/S Pos.: 36 Date: 01/29/93

Replicate 1 Time: 14:06
Peak Area (A-s): 0.223 Peak Height (A): 0.445
Background Pk Area (A-s): 0.118 Background Pk Height (A): 0.217
Blank Corrected Pk Area (A-s): 0.221
Concentration (ug/L): 91.6

Standard number 3 applied. [100.0]
Correlation coefficient: 1.00000 Slope: 0.0025

ID: EP-8 Seq. No.: 00005 A/S Pos.: 35 Date: 01/29/93

Replicate 1 Time: 14:09
Peak Area (A-s): 0.120 Peak Height (A): 0.272
Background Pk Area (A-s): 0.074 Background Pk Height (A): 0.124
Blank Corrected Pk Area (A-s): 0.118
Concentration (ug/L): 47.4 Corrected Conc (mg/L): 0.0474

Check sample is within range 45.0 - 55.0

~~~~~  
ID: PBS 1/26/93 Seq. No.: 00006 A/S Pos.: 1 Date: 01/29/93

Replicate 1 Time: 14:12  
Peak Area (A-s): -0.000 Peak Height (A): 0.004  
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.010  
Blank Corrected Pk Area (A-s): -0.003  
Concentration (ug/L ): -1.0 Corrected Conc (mg/L ): -0.0010

~~~~~  
ID: PBS 1/26/93 Seq. No.: 00007 A/S Pos.: 1 Date: 01/29/93

Replicate 1 Time: 14:15
Peak Area (A-s): 0.020 Peak Height (A): 0.045
Background Pk Area (A-s): 0.029 Background Pk Height (A): 0.027
Blank Corrected Pk Area (A-s): 0.017
Concentration (ug/L): 6.9 Corrected Conc (mg/L): 0.0069

Recovery is 79.2%

~~~~~  
ID: LCSS 1/26 1-10 Seq. No.: 00008 A/S Pos.: 2 Date: 01/29/93

Replicate 1 Time: 14:18  
Peak Area (A-s): 0.047 Peak Height (A): 0.092  
Background Pk Area (A-s): 0.047 Background Pk Height (A): 0.043  
Blank Corrected Pk Area (A-s): 0.045  
Concentration (ug/L ): 17.7 Corrected Conc (mg/L ): 0.0177

~~~~~  
ID: LCSS 1/26 1-10 Seq. No.: 00009 A/S Pos.: 2 Date: 01/29/93

Replicate 1 Time: 14:21
Peak Area (A-s): 0.072 Peak Height (A): 0.135
Background Pk Area (A-s): 0.054 Background Pk Height (A): 0.060
Blank Corrected Pk Area (A-s): 0.070
Concentration (ug/L): 27.9 Corrected Conc (mg/L): 0.0279

Recovery is 102.2%

~~~~~  
ID: 0104-01 1-10 Seq. No.: 00010 A/S Pos.: 3 Date: 01/29/93

Replicate 1 Time: 14:24  
Peak Area (A-s): 0.001 Peak Height (A): 0.005  
Background Pk Area (A-s): 0.033 Background Pk Height (A): 0.017  
Blank Corrected Pk Area (A-s): -0.001  
Concentration (ug/L ): -0.4 Corrected Conc (mg/L ): -0.0004

~~~~~  
ID: 0104-01 1-10 Seq. No.: 00011 A/S Pos.: 3 Date: 01/29/93

Replicate 1
Peak Area (A-s): 0.027
Background Pk Area (A-s): 0.042
Blank Corrected Pk Area (A-s): 0.025
Concentration (ug/L): 9.8
Time: 14:27
Peak Height (A): 0.051
Background Pk Height (A): 0.026
Corrected Conc (mg/L): 0.0098

Recovery is 102.3%

~~~~~  
ID: 0104-01 REP 10      Seq. No.: 00012      A/S Pos.: 4      Date: 01/29/93

Replicate 1  
Peak Area (A-s): 0.000  
Background Pk Area (A-s): 0.034  
Blank Corrected Pk Area (A-s): -0.002  
Concentration (ug/L ): -0.9  
Time: 14:30  
Peak Height (A): 0.004  
Background Pk Height (A): 0.019  
Corrected Conc (mg/L ): -0.0009

~~~~~  
ID: 0104-01 REP 10 Seq. No.: 00013 A/S Pos.: 4 Date: 01/29/93

Replicate 1
Peak Area (A-s): 0.025
Background Pk Area (A-s): 0.046
Blank Corrected Pk Area (A-s): 0.022
Concentration (ug/L): 8.9
Time: 14:34
Peak Height (A): 0.053
Background Pk Height (A): 0.027
Corrected Conc (mg/L): 0.0089

Recovery is 97.6%

~~~~~  
ID: 0104-01 SPK 10      Seq. No.: 00014      A/S Pos.: 5      Date: 01/29/93

Sample abs. exceeds the range of the calibration function.

Replicate 1  
Peak Area (A-s): 0.332  
Background Pk Area (A-s): 0.199  
Blank Corrected Pk Area (A-s): 0.330  
Concentration (ug/L ): -----  
Time: 14:37  
Peak Height (A): 0.586  
Background Pk Height (A): 0.343  
Corrected Conc (mg/L ): -----

~~~~~  
ID: 0104-01 SPK 10 Seq. No.: 00015 A/S Pos.: 5 Date: 01/29/93

Replicate 1
Peak Area (A-s): 0.114
Background Pk Area (A-s): 0.073
Blank Corrected Pk Area (A-s): 0.112
Concentration (ug/L): 45.1
Time: 14:40
Peak Height (A): 0.238
Background Pk Height (A): 0.109
Corrected Conc (mg/L): 0.1802

96.27% Rec.

~~~~~  
ID: 0104-01 SPK 10      Seq. No.: 00016      A/S Pos.: 5      Date: 01/29/93

Replicate 1  
Peak Area (A-s): 0.135  
Background Pk Area (A-s): 0.085  
Blank Corrected Pk Area (A-s): 0.133  
Concentration (ug/L ): 53.7  
Time: 14:43  
Peak Height (A): 0.279  
Background Pk Height (A): 0.129  
Corrected Conc (mg/L ): 0.2148

Recovery is 86.4%

ID: 0104-02 1-10 / Seq. No.: 00017 A/S Pos.: 6 Date: 01/29/93

Replicate 1 Time: 14:46  
Peak Area (A-s): 0.007 Peak Height (A): 0.012  
Background Pk Area (A-s): 0.067 Background Pk Height (A): 0.034  
Blank Corrected Pk Area (A-s): 0.004  
Concentration (ug/L ): 1.7 Corrected Conc (mg/L ): 0.0017

ID: 0104-02 1-10 Seq. No.: 00018 A/S Pos.: 6 Date: 01/29/93

Replicate 1 Time: 14:49  
Peak Area (A-s): 0.027 Peak Height (A): 0.065  
Background Pk Area (A-s): 0.083 Background Pk Height (A): 0.042  
Blank Corrected Pk Area (A-s): 0.024  
Concentration (ug/L ): 9.6 Corrected Conc (mg/L ): 0.0096

Recovery is 78.7%

ID: 0104-03 1-10 / Seq. No.: 00019 A/S Pos.: 7 Date: 01/29/93

Replicate 1 Time: 14:52  
Peak Area (A-s): 0.003 Peak Height (A): 0.010  
Background Pk Area (A-s): 0.060 Background Pk Height (A): 0.029  
Blank Corrected Pk Area (A-s): 0.001  
Concentration (ug/L ): ~~0.4~~ #41 Corrected Conc (mg/L ): 0.0004

0.4

*repeal this result 2/1/93*

ID: 0104-03 1-10 Seq. No.: 00020 A/S Pos.: 7 Date: 01/29/93

Replicate 1 Time: 14:55  
Peak Area (A-s): 0.023 Peak Height (A): 0.059  
Background Pk Area (A-s): 0.066 Background Pk Height (A): 0.034  
Blank Corrected Pk Area (A-s): 0.021  
Concentration (ug/L ): 8.3 Corrected Conc (mg/L ): 0.0083

Recovery is 79.4%

ID: 0104-04 1-10 / Seq. No.: 00021 A/S Pos.: 8 Date: 01/29/93

Replicate 1 Time: 14:58  
Peak Area (A-s): 0.002 #43 Peak Height (A): 0.009  
Background Pk Area (A-s): 0.056 Background Pk Height (A): 0.033  
Blank Corrected Pk Area (A-s): -0.000  
Concentration (ug/L ): -0.1 Corrected Conc (mg/L ): -0.0001

*repeal this result 2/1/93*

ID: 0104-04 1-10 Seq. No.: 00022 A/S Pos.: 8 Date: 01/29/93

Replicate 1 Time: 15:01  
Peak Area (A-s): 0.024 Peak Height (A): 0.060  
Background Pk Area (A-s): 0.074 Background Pk Height (A): 0.037  
Blank Corrected Pk Area (A-s): 0.022  
Concentration (ug/L ): 8.7 Corrected Conc (mg/L ): 0.0087



Recovery is 88.9%

~~~~~  
ID: 0104-05 1-10 Seq. No.: 00023 A/S Pos.: 9 Date: 01/29/93

Replicate 1 Time: 15:04
Peak Area (A-s): 0.001 Peak Height (A): 0.007
Background Pk Area (A-s): 0.044 Background Pk Height (A): 0.020
Blank Corrected Pk Area (A-s): -0.001
Concentration (ug/L): -0.4 Corrected Conc (mg/L): -0.0004

~~~~~  
ID: 0104-05 1-10                      Seq. No.: 00024                      A/S Pos.: 9                      Date: 01/29/93

Replicate 1                                      Time: 15:07  
Peak Area (A-s): 0.019                      Peak Height (A): 0.044  
Background Pk Area (A-s): 0.052                      Background Pk Height (A): 0.029  
Blank Corrected Pk Area (A-s): 0.017  
Concentration (ug/L ):                      6.8                      Corrected Conc (mg/L ):                      0.0068

Recovery is 71.7%

~~~~~  
ID: EP-8 Seq. No.: 00025 A/S Pos.: 35 Date: 01/29/93

Replicate 1 Time: 15:11
Peak Area (A-s): 0.112 Peak Height (A): 0.237
Background Pk Area (A-s): 0.071 Background Pk Height (A): 0.107
Blank Corrected Pk Area (A-s): 0.110
Concentration (ug/L): 44.3 Corrected Conc (mg/L): 0.0443

Check sample is out of range 45.0 - 55.0

~~~~~  
ID: EP-8                                      Seq. No.: 00026                      A/S Pos.: 35                      Date: 01/29/93

Replicate 1                                      Time: 15:14  
Peak Area (A-s): 0.112                      Peak Height (A): 0.234  
Background Pk Area (A-s): 0.065                      Background Pk Height (A): 0.107  
Blank Corrected Pk Area (A-s): 0.110  
Concentration (ug/L ):                      44.3                      Corrected Conc (mg/L ):                      0.0443

Check sample is out of range 45.0 - 55.0

~~~~~  
ID: EP-8 Seq. No.: 00027 A/S Pos.: 35 Date: 01/29/93

Replicate 1 Time: 15:19
Peak Area (A-s): 0.120 Peak Height (A): 0.253
Background Pk Area (A-s): 0.066 Background Pk Height (A): 0.113
Blank Corrected Pk Area (A-s): 0.117
Concentration (ug/L): 47.3 Corrected Conc (mg/L): 0.0473

Check sample is within range 45.0 - 55.0

~~~~~  
ID: PBS 1/26/93                      Seq. No.: 00029                      A/S Pos.: 1                      Date: 01/29/93

Replicate 1                                      Time: 15:23

Peak Area (A-s): 0.002 Peak Height (A): 0.004  
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.011  
Blank Corrected Pk Area (A-s): -0.001  
Concentration (ug/L ): -0.2 Corrected Conc (mg/L ): -0.0002

~~~~~  
ID: PBS 1/26/93 Seq. No.: 00030 A/S Pos.: 1 Date: 01/29/93

Replicate 1 Time: 15:27
Peak Area (A-s): 0.021 Peak Height (A): 0.050
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.024
Blank Corrected Pk Area (A-s): 0.019
Concentration (ug/L): 7.4 Corrected Conc (mg/L): 0.0074

Recovery is 76.4%

~~~~~  
ID: 0104-01 1-10 *2/1/93* Seq. No.: 00031 A/S Pos.: 3 Date: 01/29/93

Replicate 1 Time: 15:30  
Peak Area (A-s): 0.005 Peak Height (A): 0.011  
Background Pk Area (A-s): 0.162 Background Pk Height (A): 0.078  
Blank Corrected Pk Area (A-s): 0.003  
Concentration (ug/L ): 1.0 Corrected Conc (mg/L ): 0.0010

~~~~~  
ID: 0104-01 1-10 Seq. No.: 00032 A/S Pos.: 3 Date: 01/29/93

Replicate 1 Time: 15:33
Peak Area (A-s): 0.020 Peak Height (A): 0.041
Background Pk Area (A-s): 0.190 Background Pk Height (A): 0.098
Blank Corrected Pk Area (A-s): 0.018
Concentration (ug/L): 7.1 Corrected Conc (mg/L): 0.0071

Recovery is 61.2%

~~~~~  
ID: 0104-01 REP Seq. No.: 00033 A/S Pos.: 4 Date: 01/29/93

Replicate 1 Time: 15:36  
Peak Area (A-s): 0.006 Peak Height (A): 0.014  
Background Pk Area (A-s): 0.198 Background Pk Height (A): 0.085  
Blank Corrected Pk Area (A-s): 0.004  
Concentration (ug/L ): 1.5 Corrected Conc (mg/L ): 0.0015

~~~~~  
ID: 0104-01 REP Seq. No.: 00034 A/S Pos.: 4 Date: 01/29/93

Replicate 1 Time: 15:39
Peak Area (A-s): 0.020 Peak Height (A): 0.045
Background Pk Area (A-s): 0.173 Background Pk Height (A): 0.079
Blank Corrected Pk Area (A-s): 0.018
Concentration (ug/L): 7.1 Corrected Conc (mg/L): 0.0071

Recovery is 55.7%

~~~~~  
ID: 4 ppb / 40 ppb Seq. No.: 00035 A/S Pos.: 40 Date: 01/29/93

Replicate 1 Time: 15:43  
Peak Area (A-s): 0.089 Peak Height (A): 0.190  
Background Pk Area (A-s): 0.070 Background Pk Height (A): 0.087  
Blank Corrected Pk Area (A-s): 0.087  
Concentration (ug/L ): 34.6 Corrected Conc (mg/L ): 0.0346

~~~~~  
ID: 4 ppb / 40 ppb Seq. No.: 00036 A/S Pos.: 40 Date: 01/29/93

Replicate 1 Time: 15:49
Peak Area (A-s): 0.091 Peak Height (A): 0.186
Background Pk Area (A-s): 0.055 Background Pk Height (A): 0.082
Blank Corrected Pk Area (A-s): 0.089
Concentration (ug/L): 35.5 Corrected Conc (mg/L): 0.0355

~~~~~  
ID: FBS 1/26/93 Seq. No.: 00037 A/S Pos.: 1 Date: 01/29/93

Replicate 1 Time: 15:53  
Peak Area (A-s): -0.001 Peak Height (A): 0.004  
Background Pk Area (A-s): 0.016 Background Pk Height (A): 0.009  
Blank Corrected Pk Area (A-s): -0.003  
Concentration (ug/L ): -1.3 Corrected Conc (mg/L ): -0.0013

~~~~~  
ID: FBS 1/26/93 Seq. No.: 00038 A/S Pos.: 1 Date: 01/29/93

Replicate 1 Time: 15:57
Peak Area (A-s): 0.020 Peak Height (A): 0.048
Background Pk Area (A-s): 0.025 Background Pk Height (A): 0.024
Blank Corrected Pk Area (A-s): 0.018
Concentration (ug/L): 7.1 Corrected Conc (mg/L): 0.0071

Recovery is 83.8%

~~~~~  
ID: 0104-02 Seq. No.: 00039 A/S Pos.: 6 Date: 01/29/93

Replicate 1 Time: 16:00  
Peak Area (A-s): 0.030 Peak Height (A): 0.048  
Background Pk Area (A-s): 0.624 Background Pk Height (A): 0.711  
Blank Corrected Pk Area (A-s): 0.027  
Concentration (ug/L ): 10.8 Corrected Conc (mg/L ): 0.0108

~~~~~  
ID: 0104-02 Seq. No.: 00040 A/S Pos.: 6 Date: 01/29/93

Replicate 1 Time: 16:03
Peak Area (A-s): 0.039 Peak Height (A): 0.064
Background Pk Area (A-s): 0.634 Background Pk Height (A): 0.621
Blank Corrected Pk Area (A-s): 0.037
Concentration (ug/L): 14.5 Corrected Conc (mg/L): 0.0145

Recovery is 37.1%

~~~~~  
ID: 0104-03 Seq. No.: 00041 A/S Pos.: 7 Date: 01/29/93

Replicate 1 Time: 16:06

Peak Area (A-s): 0.030                      Peak Height (A): 0.052  
Background Pk Area (A-s): 0.409           Background Pk Height (A): 0.239  
Blank Corrected Pk Area (A-s): 0.028  
Concentration (ug/L ):                      11.2           Corrected Conc (mg/L ):                      0.0112

~~~~~  
ID: 0104-03 Seq. No.: 00042 A/S Pos.: 7 Date: 01/29/93

Replicate 1 Time: 16:09
Peak Area (A-s): 0.048 Peak Height (A): 0.085
Background Pk Area (A-s): 0.413 Background Pk Height (A): 0.236
Blank Corrected Pk Area (A-s): 0.046
Concentration (ug/L): 18.3 Corrected Conc (mg/L): 0.0183

Recovery is 71.7%

~~~~~  
ID: 0104-04                      Seq. No.: 00043                      A/S Pos.: 8                      Date: 01/29/93

Replicate 1                                      Time: 16:12  
Peak Area (A-s): 0.034                      Peak Height (A): 0.050  
Background Pk Area (A-s): 0.555           Background Pk Height (A): 0.441  
Blank Corrected Pk Area (A-s): 0.031  
Concentration (ug/L ):                      12.5           Corrected Conc (mg/L ):                      0.0125

~~~~~  
ID: 0104-04 Seq. No.: 00044 A/S Pos.: 8 Date: 01/29/93

Replicate 1 Time: 16:15
Peak Area (A-s): 0.048 Peak Height (A): 0.068
Background Pk Area (A-s): 0.566 Background Pk Height (A): 0.457
Blank Corrected Pk Area (A-s): 0.046
Concentration (ug/L): 18.2 Corrected Conc (mg/L): 0.0182

Recovery is 57.4%

~~~~~  
ID: 0104-05                      Seq. No.: 00045                      A/S Pos.: 9                      Date: 01/29/93

Replicate 1                                      Time: 16:18  
Peak Area (A-s): 0.012                      Peak Height (A): 0.032  
Background Pk Area (A-s): 0.273           Background Pk Height (A): 0.104  
Blank Corrected Pk Area (A-s): 0.010  
Concentration (ug/L ):                      4.0           Corrected Conc (mg/L ):                      0.0040

~~~~~  
ID: 0104-05 Seq. No.: 00046 A/S Pos.: 9 Date: 01/29/93

Replicate 1 Time: 16:21
Peak Area (A-s): 0.028 Peak Height (A): 0.068
Background Pk Area (A-s): 0.296 Background Pk Height (A): 0.110
Blank Corrected Pk Area (A-s): 0.025
Concentration (ug/L): 10.0 Corrected Conc (mg/L): 0.0100

Recovery is 60.4%

~~~~~  
ID: EP-8                                      Seq. No.: 00047                      A/S Pos.: 35                      Date: 01/29/93

Replicate 1 Time: 16:24  
Peak Area (A-s): 0.120 Peak Height (A): 0.255  
Background Pk Area (A-s): 0.080 Background Pk Height (A): 0.106  
Blank Corrected Pk Area (A-s): 0.117  
Concentration (ug/L ): 47.3 Corrected Conc (mg/L ): 0.0473

Check sample is within range 45.0 - 55.0

~~~~~  
ID: Seq. No.: 00048 A/S Pos.: 10 Date: 01/29/93

Replicate 1 Time: 16:27
Peak Area (A-s): -0.000 Peak Height (A): 0.005
Background Pk Area (A-s): 0.021 Background Pk Height (A): 0.013
Blank Corrected Pk Area (A-s): -0.002
Concentration (ug/L): -1.0 Corrected Conc (mg/L): -0.0010

~~~~~  
ID: EP-8 Seq. No.: 00050 A/S Pos.: 35 Date: 01/29/93

Replicate 1 Time: 16:32  
Peak Area (A-s): 0.119 Peak Height (A): 0.252  
Background Pk Area (A-s): 0.065 Background Pk Height (A): 0.102  
Blank Corrected Pk Area (A-s): 0.117  
Concentration (ug/L ): 47.0 Corrected Conc (mg/L ): 0.0470

Check sample is within range 45.0 - 55.0

~~~~~  
ID: IDL Seq. No.: 00051 A/S Pos.: 11 Date: 01/29/93

Replicate 1 Time: 16:35
Peak Area (A-s): 0.012 Peak Height (A): 0.030
Background Pk Area (A-s): 0.028 Background Pk Height (A): 0.020
Blank Corrected Pk Area (A-s): 0.010
Concentration (ug/L): 4.0 Corrected Conc (mg/L): 0.0040

~~~~~  
ID: IDL Seq. No.: 00052 A/S Pos.: 12 Date: 01/29/93

Replicate 1 Time: 16:38  
Peak Area (A-s): 0.015 Peak Height (A): 0.030  
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.019  
Blank Corrected Pk Area (A-s): 0.012  
Concentration (ug/L ): 4.9 Corrected Conc (mg/L ): 0.0049

~~~~~  
ID: IDL Seq. No.: 00053 A/S Pos.: 13 Date: 01/29/93

Replicate 1 Time: 16:41
Peak Area (A-s): 0.011 Peak Height (A): 0.030
Background Pk Area (A-s): 0.025 Background Pk Height (A): 0.015
Blank Corrected Pk Area (A-s): 0.008
Concentration (ug/L): 3.3 Corrected Conc (mg/L): 0.0033

~~~~~  
ID: IDL Seq. No.: 00054 A/S Pos.: 14 Date: 01/29/93

Replicate 1 Time: 16:44

Peak Area (A-s): 0.014                      Peak Height (A): 0.032  
Background Pk Area (A-s): 0.025           Background Pk Height (A): 0.019  
Blank Corrected Pk Area (A-s): 0.011  
Concentration (ug/L ):                      4.4           Corrected Conc (mg/L ):                      0.0044

~~~~~  
ID: IDL Seq. No.: 00055 A/S Pos.: 15 Date: 01/29/93

Replicate 1 Time: 16:47
Peak Area (A-s): 0.012 Peak Height (A): 0.030
Background Pk Area (A-s): 0.026 Background Pk Height (A): 0.019
Blank Corrected Pk Area (A-s): 0.010
Concentration (ug/L): 3.9 Corrected Conc (mg/L): 0.0039

~~~~~  
ID: EP-8                      Seq. No.: 00056                      A/S Pos.: 35                      Date: 01/29/93

Replicate 1                      Time: 16:50  
Peak Area (A-s): 0.133                      Peak Height (A): 0.270  
Background Pk Area (A-s): 0.069           Background Pk Height (A): 0.107  
Blank Corrected Pk Area (A-s): 0.131  
Concentration (ug/L ):                      52.8           Corrected Conc (mg/L ):                      0.0528

Check sample is within range 45.0 - 55.0

~~~~~  
ID: IDL Seq. No.: 00057 A/S Pos.: 16 Date: 01/29/93

Replicate 1 Time: 16:53
Peak Area (A-s): 0.014 Peak Height (A): 0.033
Background Pk Area (A-s): 0.022 Background Pk Height (A): 0.019
Blank Corrected Pk Area (A-s): 0.012
Concentration (ug/L): 4.8 Corrected Conc (mg/L): 0.0048

~~~~~  
ID: IDL                      Seq. No.: 00058                      A/S Pos.: 17                      Date: 01/29/93

Replicate 1                      Time: 16:56  
Peak Area (A-s): 0.013                      Peak Height (A): 0.030  
Background Pk Area (A-s): 0.023           Background Pk Height (A): 0.018  
Blank Corrected Pk Area (A-s): 0.011  
Concentration (ug/L ):                      4.4           Corrected Conc (mg/L ):                      0.0044

$\bar{X} = 4.24$        $SD = .556$        $IDL = 1.669$

~~~~~  
ID: EP-8 Seq. No.: 00059 A/S Pos.: 35 Date: 01/29/93

Replicate 1 Time: 17:26
Peak Area (A-s): 0.133 Peak Height (A): 0.294
Background Pk Area (A-s): 0.066 Background Pk Height (A): 0.114
Blank Corrected Pk Area (A-s): 0.131
Concentration (ug/L): 52.8 Corrected Conc (mg/L): 0.0528

Thallium

P.E. AA GRAPHITE FURNACE 5100Z RUN LOG

P.E. AA GRAPHITE FURNACE 5100Z RUN LOG

Date: 2/2/93

Element: Tl

Date Storage File:

Lamp: Hamamatsu HCL

F220231

Energy: 51

Current: 5

Char. Mass @ 20 ppb: 54

Temp. Program:

60 ppb: 38

100 ppb: 32

POS

SAMPLE ID

1	PBS 1/26/93	21
2	LCSS 1/26/93 1-20	22
3	0104-01	23
4	0104-01 REP	24
5	0104-01 SPK 1-50	25
6	0104-02	26
7	0104-03	27
8	0104-04	28
9	0104-05	29
10		30
11		31
12		32
13		33
14		34
15		35
16		36 100 ppb standard
17		37 Blank
18		38 NH4PO4 Mg(NO3)2 modifier
19		39 Pd/Mg
20		40 40 ppb spk sol

Element File: TLBRIEF.GEL Element: Tl
Date: 02/02/93 Time: 09:27
Data Storage File: F220231.DAT ID/Weight File: 0104.A60
Technique: HGA Calibration Type: Nonlinear
Remark 1: Analysis was run with 4% ammonium phosphate/.2% magnesium
Remark 2: nitrate matrix modifier.
Remark 3: Integration time = 5 seconds.

~~~~~  
ID: Blank                               Seq. No.: 00001                    A/S Pos.: 37                    Date: 02/02/93

Replicate 1                               Time: 09:27  
Peak Area (A-s): 0.004                    Peak Height (A): 0.009  
Background Pk Area (A-s): 0.009           Background Pk Height (A): 0.014  
Blank Corrected Pk Area (A-s): 0.004

Auto-zero performed.

~~~~~  
ID: Standard 1 Seq. No.: 00002 A/S Pos.: 36 Date: 02/02/93

Replicate 1 Time: 09:30
Peak Area (A-s): 0.078 Peak Height (A): 0.147
Background Pk Area (A-s): 0.077 Background Pk Height (A): 0.120
Blank Corrected Pk Area (A-s): 0.074

Standard number 1 applied. [20.0]
Correlation coefficient: 1.00000 Slope: 0.0037

~~~~~  
ID: Standard 2                            Seq. No.: 00003                    A/S Pos.: 36                    Date: 02/02/93

Replicate 1                               Time: 09:33  
Peak Area (A-s): 0.194                    Peak Height (A): 0.340  
Background Pk Area (A-s): 0.176           Background Pk Height (A): 0.322  
Blank Corrected Pk Area (A-s): 0.190  
Concentration (ug/L ):                    51.5

Standard number 2 applied. [60.0]  
Correlation coefficient: 1.00000           Slope: 0.0040

~~~~~  
ID: Standard 3 Seq. No.: 00004 A/S Pos.: 36 Date: 02/02/93

Replicate 1 Time: 09:35
Peak Area (A-s): 0.277 Peak Height (A): 0.437
Background Pk Area (A-s): 0.269 Background Pk Height (A): 0.452
Blank Corrected Pk Area (A-s): 0.273
Concentration (ug/L): 97.9

Standard number 3 applied. [100.0]
Correlation coefficient: 1.00000 Slope: 0.0040

~~~~~  
ID: EP - 8                                Seq. No.: 00005                    A/S Pos.: 34                    Date: 02/02/93

Replicate 1                               Time: 09:38

Peak Area (A-s): -0.003      Peak Height (A): 0.007  
Background Pk Area (A-s): 0.012      Background Pk Height (A): 0.010  
Blank Corrected Pk Area (A-s): -0.007  
Concentration (ug/L ): -1.9

Check sample is out of range 36 - 44

~~~~~  
ID: EP-8 Seq. No.: 00006 A/S Pos.: 35 Date: 02/02/93

Replicate 1 Time: 09:44
Peak Area (A-s): 0.129 Peak Height (A): 0.223
Background Pk Area (A-s): 0.119 Background Pk Height (A): 0.194
Blank Corrected Pk Area (A-s): 0.125
Concentration (ug/L): 36.0

Check sample is out of range 36 - 44

~~~~~  
ID: EP-8      Seq. No.: 00007      A/S Pos.: 35      Date: 02/02/93

Replicate 1      Time: 09:57  
Peak Area (A-s): 0.136      Peak Height (A): 0.234  
Background Pk Area (A-s): 0.119      Background Pk Height (A): 0.199  
Blank Corrected Pk Area (A-s): 0.132  
Concentration (ug/L ): 38.3

Check sample is within range 36 - 44

~~~~~  
ID: PBS 1/26/93 Seq. No.: 00008 A/S Pos.: 1 Date: 02/02/93

Replicate 1 Time: 09:59
Peak Area (A-s): 0.001 Peak Height (A): 0.009
Background Pk Area (A-s): 0.011 Background Pk Height (A): 0.013
Blank Corrected Pk Area (A-s): -0.003
Concentration (ug/L): -0.7

~~~~~  
ID: PBS 1/26/93      Seq. No.: 00009      A/S Pos.: 1      Date: 02/02/93

Replicate 1      Time: 10:02  
Peak Area (A-s): 0.072      Peak Height (A): 0.130  
Background Pk Area (A-s): 0.062      Background Pk Height (A): 0.107  
Blank Corrected Pk Area (A-s): 0.068  
Concentration (ug/L ): 18.4

Recovery is 95.4%

~~~~~  
ID: LCSS 1/26 1-20 Seq. No.: 00010 A/S Pos.: 2 Date: 02/02/93

Replicate 1 Time: 10:05
Peak Area (A-s): 0.109 Peak Height (A): 0.211
Background Pk Area (A-s): 0.113 Background Pk Height (A): 0.177
Blank Corrected Pk Area (A-s): 0.105
Concentration (ug/L): 29.4 Corrected Conc (ug/L): 588.

~~~~~  
*117.6 mg/kg. 84.07. Rec.*

ID: LCSS 1/26 1-20      Seq. No.: 00011      A/S Pos.: 2      Date: 02/02/93

Replicate 1      Time: 10:08  
Peak Area (A-s): 0.161      Peak Height (A): 0.283  
Background Pk Area (A-s): 0.150      Background Pk Height (A): 0.254  
Blank Corrected Pk Area (A-s): 0.157  
Concentration (ug/L ):      47.3      Corrected Conc (ug/L ):      946.

Recovery is 89.4%

~~~~~  
ID: 0104-01 Seq. No.: 00012 A/S Pos.: 3 Date: 02/02/93

Replicate 1 Time: 10:11
Peak Area (A-s): 0.007 Peak Height (A): 0.015
Background Pk Area (A-s): 0.382 Background Pk Height (A): 0.164
Blank Corrected Pk Area (A-s): 0.003
Concentration (ug/L): 0.7

~~~~~  
ID: 0104-01      Seq. No.: 00013      A/S Pos.: 3      Date: 02/02/93

Replicate 1      Time: 10:14  
Peak Area (A-s): 0.074      Peak Height (A): 0.129  
Background Pk Area (A-s): 0.368      Background Pk Height (A): 0.137  
Blank Corrected Pk Area (A-s): 0.070  
Concentration (ug/L ):      18.8

Recovery is 90.1%

~~~~~  
ID: 0104-01 REP Seq. No.: 00014 A/S Pos.: 4 Date: 02/02/93

Replicate 1 Time: 10:17
Peak Area (A-s): 0.011 Peak Height (A): 0.020
Background Pk Area (A-s): 0.353 Background Pk Height (A): 0.150
Blank Corrected Pk Area (A-s): 0.007
Concentration (ug/L): 1.8

~~~~~  
ID: 0104-01 REP      Seq. No.: 00015      A/S Pos.: 4      Date: 02/02/93

Replicate 1      Time: 10:19  
Peak Area (A-s): 0.069      Peak Height (A): 0.132  
Background Pk Area (A-s): 0.355      Background Pk Height (A): 0.127  
Blank Corrected Pk Area (A-s): 0.065  
Concentration (ug/L ):      17.5

Recovery is 78.7%

~~~~~  
ID: 0104-01 SPK 50 Seq. No.: 00016 A/S Pos.: 5 Date: 02/02/93

Replicate 1 Time: 10:22
Peak Area (A-s): 0.112 Peak Height (A): 0.205
Background Pk Area (A-s): 0.115 Background Pk Height (A): 0.172
Blank Corrected Pk Area (A-s): 0.108
Concentration (ug/L): 30.3

75.8% Rec.

~~~~~  
ID: 0104-01 SPK 50      Seq. No.: 00017      A/S Pos.: 5      Date: 02/02/93

Replicate 1      Time: 10:25  
Peak Area (A-s): 0.163      Peak Height (A): 0.299  
Background Pk Area (A-s): 0.166      Background Pk Height (A): 0.271  
Blank Corrected Pk Area (A-s): 0.159  
Concentration (ug/L ): 48.1

Recovery is 88.6%

~~~~~  
ID: 0104-02 Seq. No.: 00018 A/S Pos.: 6 Date: 02/02/93

Replicate 1 Time: 10:28
Peak Area (A-s): 0.004 Peak Height (A): 0.010
Background Pk Area (A-s): 0.429 Background Pk Height (A): 0.208
Blank Corrected Pk Area (A-s): -0.001
Concentration (ug/L): -0.1

~~~~~  
ID: 0104-02      Seq. No.: 00019      A/S Pos.: 6      Date: 02/02/93

Replicate 1      Time: 10:31  
Peak Area (A-s): 0.044      Peak Height (A): 0.092  
Background Pk Area (A-s): 0.470      Background Pk Height (A): 0.214  
Blank Corrected Pk Area (A-s): 0.040  
Concentration (ug/L ): 10.4

Recovery is 52.8%

~~~~~  
ID: 0104-03 Seq. No.: 00020 A/S Pos.: 7 Date: 02/02/93

Replicate 1 Time: 10:33
Peak Area (A-s): -0.000 Peak Height (A): 0.010
Background Pk Area (A-s): 0.333 Background Pk Height (A): 0.162
Blank Corrected Pk Area (A-s): -0.004
Concentration (ug/L): -1.0

~~~~~  
ID: 0104-03      Seq. No.: 00021      A/S Pos.: 7      Date: 02/02/93

Replicate 1      Time: 10:36  
Peak Area (A-s): 0.031      Peak Height (A): 0.057  
Background Pk Area (A-s): 0.379      Background Pk Height (A): 0.170  
Blank Corrected Pk Area (A-s): 0.027  
Concentration (ug/L ): 7.0

Recovery is 40.2%

~~~~~  
ID: 0104-04 Seq. No.: 00022 A/S Pos.: 8 Date: 02/02/93

Replicate 1 Time: 10:39
Peak Area (A-s): 0.005 Peak Height (A): 0.013
Background Pk Area (A-s): 0.452 Background Pk Height (A): 0.210
Blank Corrected Pk Area (A-s): 0.001

Concentration (ug/L): 0.3

~~~~~  
ID: 0104-04 Seq. No.: 00023 A/S Pos.: 8 Date: 02/02/93

Replicate 1 Time: 10:42  
Peak Area (A-s): 0.040 Peak Height (A): 0.084  
Background Pk Area (A-s): 0.470 Background Pk Height (A): 0.204  
Blank Corrected Pk Area (A-s): 0.036  
Concentration (ug/L ): 9.3

Recovery is 45.2%

~~~~~  
ID: 0104-05 Seq. No.: 00024 A/S Pos.: 9 Date: 02/02/93

Replicate 1 Time: 10:45
Peak Area (A-s): 0.001 Peak Height (A): 0.008
Background Pk Area (A-s): 0.376 Background Pk Height (A): 0.169
Blank Corrected Pk Area (A-s): -0.003
Concentration (ug/L): -0.9

~~~~~  
ID: 0104-05 Seq. No.: 00025 A/S Pos.: 9 Date: 02/02/93

Replicate 1 Time: 10:47  
Peak Area (A-s): 0.065 Peak Height (A): 0.122  
Background Pk Area (A-s): 0.384 Background Pk Height (A): 0.149  
Blank Corrected Pk Area (A-s): 0.061  
Concentration (ug/L ): 16.3

Recovery is 85.8%

~~~~~  
ID: EP-8 Seq. No.: 00026 A/S Pos.: 35 Date: 02/02/93

Replicate 1 Time: 10:50
Peak Area (A-s): 0.133 Peak Height (A): 0.238
Background Pk Area (A-s): 0.130 Background Pk Height (A): 0.208
Blank Corrected Pk Area (A-s): 0.129
Concentration (ug/L): 37.4

Check sample is within range 36 - 44

Time: 10:53

Peak Height (A): 0.009

Background Pk Height (A): 0.013

Blank Corrected Fk Area (A-s): -0.004

Concentration (uo/L) : -0.9

ID:

Seq. No.: 00028

A/S Pos.: 10

Date: 02/02/93

Time: 10:56

Peak Height (A): 0.144

Background Pk Height (A): 0.115

Blank Corrected Fk Area (A-s): 0.068

Concentration (ug/L): 18.4

Recovery is 96.4%

~~~~~  
ID: IDL 1                      Seq. No.: 00029              A/S Pos.: 11              Date: 02/02/93

Replicate 1                      Time: 11:09  
Peak Area (A-s): 0.033              Peak Height (A): 0.067  
Background Pk Area (A-s): 0.046              Background Pk Height (A): 0.048  
Blank Corrected Pk Area (A-s): 0.029  
Concentration (ug/L ):              7.5

~~~~~  
ID: IDL 2 Seq. No.: 00030 A/S Pos.: 12 Date: 02/02/93

Replicate 1 Time: 11:12
Peak Area (A-s): 0.033 Peak Height (A): 0.068
Background Pk Area (A-s): 0.049 Background Pk Height (A): 0.051
Blank Corrected Pk Area (A-s): 0.029
Concentration (ug/L): 7.4

~~~~~  
ID: IDL 3                      Seq. No.: 00031              A/S Pos.: 13              Date: 02/02/93

Replicate 1                      Time: 11:15  
Peak Area (A-s): 0.034              Peak Height (A): 0.068  
Background Pk Area (A-s): 0.047              Background Pk Height (A): 0.050  
Blank Corrected Pk Area (A-s): 0.030  
Concentration (ug/L ):              7.9

~~~~~  
ID: IDL 4 Seq. No.: 00032 A/S Pos.: 14 Date: 02/02/93

Replicate 1 Time: 11:18
Peak Area (A-s): 0.035 Peak Height (A): 0.062
Background Pk Area (A-s): 0.041 Background Pk Height (A): 0.049
Blank Corrected Pk Area (A-s): 0.031
Concentration (ug/L): 7.9

~~~~~  
ID: IDL 5                      Seq. No.: 00033              A/S Pos.: 15              Date: 02/02/93

Replicate 1                      Time: 11:21  
Peak Area (A-s): 0.031              Peak Height (A): 0.064  
Background Pk Area (A-s): 0.038              Background Pk Height (A): 0.050  
Blank Corrected Pk Area (A-s): 0.027  
Concentration (ug/L ):              6.9

~~~~~  
ID: EP-8 Seq. No.: 00034 A/S Pos.: 35 Date: 02/02/93

Replicate 1 Time: 11:23
Peak Area (A-s): 0.129 Peak Height (A): 0.229
Background Pk Area (A-s): 0.131 Background Pk Height (A): 0.194
Blank Corrected Pk Area (A-s): 0.125
Concentration (ug/L): 35.9

~~~~~  
Check sample is out of range 36 - 44  
~~~~~

Replicate 1 Time: 11:33
Peak Area (A-s): 0.034 Peak Height (A): 0.061
Background Pk Area (A-s): 0.043 Background Pk Height (A): 0.052
Blank Corrected Pk Area (A-s): 0.030
Concentration (ug/L): 7.6

~~~~~  
ID: IDL 7 Seq. No.: 00037 A/S Pos.: 17 Date: 02/02/93

Replicate 1 Time: 11:36  
Peak Area (A-s): 0.037 Peak Height (A): 0.065  
Background Pk Area (A-s): 0.044 Background Pk Height (A): 0.053  
Blank Corrected Pk Area (A-s): 0.033  
Concentration (ug/L ): 8.6

$\bar{X} = 7.41$   
 $SP = .507$   
 $IPL = 1.58$



**SPEX**

# CERTIFICATE OF ANALYSIS

Catalog Number: SPIKE-1

Description: Spike Sample Standard

Lot Number: 4-257AS

*Rec 5/1/93***ICP CHECK:**

Element Labeled[ug/ml] Measured[ug/ml] Element Labeled[ug/ml] Measured[ug/ml]

|    |     |        |    |     |        |
|----|-----|--------|----|-----|--------|
| Al | 200 | 199.54 | Fe | 100 | 100.07 |
| Sb | 50  | 49.92  | Pb | 50  | 49.93  |
| As | 200 | 200.50 | Mn | 50  | 49.86  |
| Ba | 200 | 199.29 | Ni | 50  | 50.10  |
| Be | 5   | 4.99   | Se | 200 | 201.06 |
| Cd | 5   | 5.01   | Ag | 5   | 4.96   |
| Cr | 20  | 19.96  | Tl | 200 | 198.68 |
| Co | 50  | 50.09  | V  | 50  | 49.79  |
| Cu | 25  | 25.32  | Zn | 50  | 49.96  |

**Instrumental Analysis by Inductively Coupled Plasma Spectroscopy:**

The following SRM's were used in establishing the above results:

Al: NIST 3101, Sb: NIST 3102, As: NIST 3103, Ba: NIST 3104,  
Be: NIST 3105, Cd: NIST 3108, Cr: NIST 3112, Co: NIST 3113,  
Cu: NIST 3114, Fe: NIST 3126, Pb: NIST 3128, Mn: NIST 3132,  
Ni: NIST 3136, Se: NIST 3149, Ag: NIST 3151, Tl: NIST 3158,  
V: NIST 3165, Zn: NIST 3168.

Spex Reference Multi: Lot #1-79STDBD, 3-150AS

Balances are calibrated with NIST weight sets N.J. #76552, #76543, #82395, according to NIST circular 547 3.4.3.

SPEX plasma solution standards are guaranteed stable and accurate to  $\pm 0.5\%$  of labeled concentration for one year from date of purchase. This value is the sum of cumulative errors associated with analytical determinations, pipetting and diluting to final volume. For these solutions we use high purity acids, 18 megohm double deionized water and triple rinsed bottles. All glassware used is class A.

*N. Kochendakote*

MAY 93

Signed by: \_\_\_\_\_, Chemical Production Manager, Date: \_\_\_\_\_

**SPEX**

Printed on Recycled Paper

# HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-B

HAS Sample #93-0104-02

Date Sampled: 1/19/93

Date Prepared: 1/26/93

| ANALYTE   | EPA<br>METHOD | DATE<br>ANALYZE | RESULT | FINAL<br>VOL.(L) | (ml) (gm)<br>DRY WT. | METHO<br>DL | DET.LIMIT<br>(ug/kg) | RESULT<br>ug/kg | QC  |
|-----------|---------------|-----------------|--------|------------------|----------------------|-------------|----------------------|-----------------|-----|
| ALUMINUM  | 6010          | 1/28/93         | 15000  | 2                | 0.834                | 30.0        | 71900                | 36000000        | *95 |
| ANTIMONY  | 6010          | 1/28/93         | -15.7  | 2                | 0.834                | 50.0        | 120000               | <DL             | *95 |
| ARSENIC   | 7060          | 2/02/93         | 0.5    | 2                | 0.834                | 10.0        | 24000                | <DL             | *95 |
| BARIUM    | 6010          | 1/28/93         | 113    | 2                | 0.834                | 10.0        | 24000                | 271000          | *95 |
| BERYLLIUM | 6010          | 1/28/93         | 5.42   | 2                | 0.834                | 5.00        | 12000                | 13000           | *95 |
| CADMIUM   | 6010          | 1/28/93         | 0.0001 | 2                | 0.834                | 5.00        | 12000                | <DL             | *95 |
| CALCIUM   | 6010          | 1/28/93         | 91200  | 2                | 0.834                | 20.0        | 48000                | 219000000       | *95 |
| CHROMIUM  | 6010          | 1/28/93         | 19.5   | 2                | 0.834                | 10.0        | 24000                | 46800           | *95 |
| COBALT    | 6010          | 1/28/93         | -2.44  | 2                | 0.834                | 10.0        | 24000                | <DL             | *95 |
| COPPER    | 6010          | 1/28/93         | 6.03   | 2                | 0.834                | 10.0        | 24000                | <DL             | *95 |
| IRON      | 6010          | 1/28/93         | 11200  | 2                | 0.834                | 20.0        | 48000                | 26900000        | *95 |
| LEAD      | 7421          | 2/01/93         | 5      | 2                | 0.834                | 5.00        | 12000                | 12000           | *95 |
| MAGNESIUM | 6010          | 1/28/93         | 22700  | 2                | 0.834                | 40.0        | 95900                | 54400000        | *95 |
| MANGANESE | 6010          | 1/28/93         | 930    | 2                | 0.834                | 10.0        | 24000                | 2230000         | *95 |
| MERCURY   | 7471          | 1/27/93         | 0.064  | 0.1              | 0.192                | 0.20        | 104                  | <DL             | *95 |
| NICKEL    | 6010          | 1/28/93         | -2.93  | 2                | 0.834                | 40.0        | 95900                | <DL             | *95 |
| POTASSIUM | 6010          | 1/28/93         | 1770   | 2                | 0.834                | 3000        | 7190000              | <DL             | *95 |
| SELENIUM  | 7740          | 2/01/93         | 1.7    | 2                | 0.834                | 5.00        | 12000                | <DL             | *95 |
| SILVER    | 6010          | 1/28/93         | -1.14  | 2                | 0.834                | 10.0        | 24000                | <DL             | *95 |
| SODIUM    | 6010          | 1/28/93         | 626    | 2                | 0.834                | 50.0        | 120000               | 1500000         | *95 |
| THALLIUM  | 7841          | 2/02/93         | -0.1   | 0.2              | 0.834                | 10.0        | 2400                 | <DL             | *95 |
| VANADIUM  | 6010          | 1/28/93         | 3.12   | 2                | 0.834                | 20.0        | 48000                | <DL             | *95 |
| ZINC      | 6010          | 1/28/93         | 13.1   | 2                | 0.834                | 20.0        | 48000                | <DL             | *95 |

\*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON DRY WEIGHT.

THIS REPORT HAS BEEN AMENDED TO INCLUDE THE RESULTS OF THE ANALYSIS OF LEAD BY MSA.

# HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-C

HAS Sample #93-0104-03

Date Sampled: 1/19/93

Date Prepared: 1/26/93

| ANALYTE   | EPA<br>METHOD | DATE<br>ANALYZE | RESULT | FINAL<br>VOL.(L) | (ml) (gm)<br>DRY WT. | METHO<br>DL | DET.LIMIT<br>(ug/kg) | RESULT<br>ug/kg | QC  |
|-----------|---------------|-----------------|--------|------------------|----------------------|-------------|----------------------|-----------------|-----|
| ALUMINUM  | 6010          | 1/28/93         | 6860   | 2                | 0.797                | 30.0        | 75300                | 17200000        | *95 |
| ANTIMONY  | 6010          | 1/28/93         | -3.3   | 2                | 0.797                | 50.0        | 125000               | <DL             | *95 |
| ARSENIC   | 7060          | 2/02/93         | 2.6    | 2                | 0.797                | 10.0        | 25100                | <DL             | *95 |
| BARIUM    | 6010          | 1/28/93         | 56.2   | 2                | 0.797                | 10.0        | 25100                | 141000          | *95 |
| BERYLLIUM | 6010          | 1/28/93         | 2.85   | 2                | 0.797                | 5.00        | 12500                | <DL             | *95 |
| CADMIUM   | 6010          | 1/28/93         | 0.434  | 2                | 0.797                | 5.00        | 12500                | <DL             | *95 |
| CALCIUM   | 6010          | 1/28/93         | 91500  | 2                | 0.797                | 20.0        | 50200                | 230000000       | *95 |
| CHROMIUM  | 6010          | 1/28/93         | 241    | 2                | 0.797                | 10.0        | 25100                | 605000          | *95 |
| COBALT    | 6010          | 1/28/93         | 8.43   | 2                | 0.797                | 10.0        | 25100                | <DL             | *95 |
| COPPER    | 6010          | 1/28/93         | 32.5   | 2                | 0.797                | 10.0        | 25100                | 81600           | *95 |
| IRON      | 6010          | 1/28/93         | 56300  | 2                | 0.797                | 20.0        | 50200                | 141000000       | *95 |
| LEAD      | 7421          | 2/01/93         | 8      | 10               | 0.797                | 5.00        | 62700                | 100000          | *95 |
| MAGNESIUM | 6010          | 1/28/93         | 9850   | 2                | 0.797                | 40.0        | 100000               | 24700000        | *95 |
| MANGANESE | 6010          | 1/28/93         | 3890   | 2                | 0.797                | 10.0        | 25100                | 9760000         | *95 |
| MERCURY   | 7471          | 1/27/93         | 0.04   | 0.1              | 0.189                | 0.20        | 106                  | <DL             | *95 |
| NICKEL    | 6010          | 1/28/93         | -27    | 2                | 0.797                | 40.0        | 100000               | <DL             | *95 |
| POTASSIUM | 6010          | 1/28/93         | 1640   | 2                | 0.797                | 3000        | 7530000              | <DL             | *95 |
| SELENIUM  | 7740          | 2/01/93         | 0.4    | 2                | 0.797                | 5.00        | 12500                | <DL             | *95 |
| SILVER    | 6010          | 1/28/93         | -10.3  | 2                | 0.797                | 10.0        | 25100                | <DL             | *95 |
| SODIUM    | 6010          | 1/28/93         | 191    | 2                | 0.797                | 50.0        | 125000               | 479000          | *95 |
| THALLIUM  | 7841          | 2/02/93         | -1     | 0.2              | 0.797                | 10.0        | 2510                 | <DL             | *95 |
| VANADIUM  | 6010          | 1/28/93         | 81     | 2                | 0.797                | 20.0        | 50200                | 203000          | *95 |
| ZINC      | 6010          | 1/28/93         | 22.7   | 2                | 0.797                | 20.0        | 50200                | 57000           | *95 |

\*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON DRY WEIGHT.

THIS REPORT HAS BEEN AMENDED TO INCLUDE THE RESULTS OF THE ANALYSIS OF LEAD BY MSA.

# HUNTINGDON ANALYTICAL SERVICES

Sample ID: SS-D  
 HAS Sample #93-0104-04  
 Date Sampled: 1/19/93  
 Date Prepared: 1/26/93

| ANALYTE   | EPA<br>METHOD | DATE<br>ANALYZE | RESULT | FINAL<br>VOL.(L) | (ml) (gm)<br>DRY WT. | METHO<br>DL | DET.LIMIT<br>(ug/kg) | RESULT<br>ug/kg | QC  |
|-----------|---------------|-----------------|--------|------------------|----------------------|-------------|----------------------|-----------------|-----|
| ALUMINUM  | 6010          | 1/28/93         | 11600  | 2                | 0.827                | 30.0        | 72600                | 28100000        | *95 |
| ANTIMONY  | 6010          | 1/28/93         | 7.43   | 2                | 0.827                | 50.0        | 121000               | <DL             | *95 |
| ARSENIC   | 7060          | 2/02/93         | 3.1    | 2                | 0.827                | 10.0        | 24200                | <DL             | *95 |
| BARIUM    | 6010          | 1/28/93         | 90.8   | 2                | 0.827                | 10.0        | 24200                | 220000          | *95 |
| BERYLLIUM | 6010          | 1/28/93         | 4.47   | 2                | 0.827                | 5.00        | 12100                | <DL             | *95 |
| CADMIUM   | 6010          | 1/28/93         | -0.434 | 2                | 0.827                | 5.00        | 12100                | <DL             | *95 |
| CALCIUM   | 6010          | 1/28/93         | 76700  | 2                | 0.827                | 20.0        | 48400                | 185000000       | *95 |
| CHROMIUM  | 6010          | 1/28/93         | 54     | 2                | 0.827                | 10.0        | 24200                | 131000          | *95 |
| COBALT    | 6010          | 1/28/93         | 2.98   | 2                | 0.827                | 10.0        | 24200                | <DL             | *95 |
| COPPER    | 6010          | 1/28/93         | 15.5   | 2                | 0.827                | 10.0        | 24200                | 37500           | *95 |
| IRON      | 6010          | 1/28/93         | 37200  | 2                | 0.827                | 20.0        | 48400                | 90000000        | *95 |
| LEAD      | 7421          | 2/01/93         | 58     | 0.2              | 0.827                | 5.00        | 1210                 | 14000           | *95 |
| MAGNESIUM | 6010          | 1/28/93         | 17500  | 2                | 0.827                | 40.0        | 96700                | 42300000        | *95 |
| MANGANESE | 6010          | 1/28/93         | 1190   | 2                | 0.827                | 10.0        | 24200                | 2880000         | *95 |
| MERCURY   | 7471          | 1/27/93         | 0.064  | 0.1              | 0.198                | 0.20        | 101                  | <DL             | *95 |
| NICKEL    | 6010          | 1/28/93         | -16.6  | 2                | 0.827                | 40.0        | 96700                | <DL             | *95 |
| POTASSIUM | 6010          | 1/28/93         | 1670   | 2                | 0.827                | 3000        | 7260000              | <DL             | *95 |
| SELENIUM  | 7740          | 2/01/93         | -0.1   | 2                | 0.827                | 5.00        | 12100                | <DL             | *95 |
| SILVER    | 6010          | 1/28/93         | -4.56  | 2                | 0.827                | 10.0        | 24200                | <DL             | *95 |
| SODIUM    | 6010          | 1/28/93         | 329    | 2                | 0.827                | 50.0        | 121000               | 796000          | *95 |
| THALLIUM  | 7841          | 2/02/93         | 0.3    | 0.2              | 0.827                | 10.0        | 2420                 | <DL             | *95 |
| VANADIUM  | 6010          | 1/28/93         | 16.3   | 2                | 0.827                | 20.0        | 48400                | <DL             | *95 |
| ZINC      | 6010          | 1/28/93         | 8.54   | 2                | 0.827                | 20.0        | 48400                | <DL             | *95 |

\*THIS INDICATES A 95% CONFIDENCE LIMIT ACHIEVED WITH AN  
 EPA QUALITY CONTROL SOLUTION ANALYZED ALONG WITH YOUR  
 SAMPLE.

ALL SOIL/SLUDGE SAMPLE RESULTS ARE BASED UPON  
 DRY WEIGHT.

THIS REPORT HAS BEEN AMENDED TO INCLUDE THE RESULTS OF  
 THE ANALYSIS OF LEAD BY MSA.

Element File: PBMSA.GEL

Element: Pb

Print Data: Main+Suppl.

Print:

Remarks:

Matrix modifier: 4% NH<sub>4</sub>PO<sub>4</sub> .2% MoNO<sub>3</sub> Integration time: 3 s

Lamp: Temp. program:

Energy: Current:

Char. Mass: 20 ppb:

60 ppb: 100 ppb:

Analyst: ECG

Peak Storage: 1 Repl./Sample

INSTRUMENT: 4100 ZL

Wavelength: 283.3 Peak

Signal Type: Zeeman AA

Read Time: 3.0

Sample Replicates: 1

Standard Replicates: 1

Technique: HGA

Slit: 0.70 Low

Signal Measurement: Peak Area

Read Delay: 0.0

BOC Time: 2

Spike Replicates: Same as Sample

CALIBRATION:

| Solutions    | ID     | Conc | Location | Volume | Diluent | Modifier |
|--------------|--------|------|----------|--------|---------|----------|
|              |        |      |          |        | Volume  | #1 #2    |
| Calib. Blank |        |      | 37       | 20     |         | 5        |
| Standard 1   | 10 ppb | 10.  | 36       | 2      | 18      | 5        |
| Standard 2   | 30 ppb | 30.  | 36       | 6      | 14      | 5        |
| Standard 3   | 50 ppb | 50.  | 36       | 10     | 10      | 5        |
| Samples      |        |      |          | 20     |         | 5        |

Diluent Location: 0

Modifier #1 Location: 38

Calibration Units: ug/L

Calibration Type: Method of Add.

Modifier #2 Location:

Sample Units: ug/L

Furnace Time/Temperature Program:

| Step | Temp | Ramp | Hold | Gas Flow | Read | Gas Type |
|------|------|------|------|----------|------|----------|
| 1    | 120  | 1    | 40   | 250      |      | Norm     |
|      | 140  | 1    | 45   | 250      |      | Norm     |
|      | 800  | 7    | 20   | 250      |      | Norm     |
| 4    | 1500 | 0    | 3    | 0        | *    | Norm     |
| 5    | 2400 | 1    | 2    | 250      |      | Norm     |

Injection Temp: 120

Pipette Speed: 100%

Extraction System: On

SEQUENCE:

Step Action and Parameters

1 Pipet diluent + modifier 1 + spike + sample/std

2 Run HGA steps 1 to End

CHK DKS:

Recalibration Type: Completely Recalibrate

Locations: None

Do c. Above Calibration Action: Dilute & Reanalyze After 1 Rep

Alternate Sample Volumes (uL): 10.5.4

Run Alternate Volume Blanks: No

If %RSD > 20.0 and Peak Area > .05 then Retrv 1 times

Check %RSD on: Samples

Recovery Measurements:

10 uL of 40 ug/L Standard at Location 40 Gives 20.000 ug/L

Measure Recovery on Samples: All Samples

Add to QC Samples: No % Recovery Limits: 85 to 115

QC:

| # | A/S  | QC Sample | Conc. Limits | After | Periodic | At    | Count As   |
|---|------|-----------|--------------|-------|----------|-------|------------|
| 1 | Loc. | ID        | Lower        | Upper | Calib    | Check | End Sample |
| 1 | 35   | ICV       | 45           | 55    | X        |       |            |
| 2 | 37   | ICB       | -10          | 10    |          |       |            |
| 3 | 39   | CRA       | 0            | 10    |          |       |            |
| 4 | 35   | CCV       | 45           | 55    |          | X     | X          |
| 5 | 37   | CCB       | -10          | 10    |          |       |            |

QC Periodic QC Samples: Every 10

Out of Limit Action: Stop Element and Continue With Next

Matrix Check Calculations:

% Difference for Duols: No

Locations:

% Recovery for Spike: No

Locations:

Conc:

Element: File: PBMSA.GEL      Element: Pb      Wavelength: 283.3  
Date: 03/25/93      Time: 15:53      Slit: 0.70 L  
Data File: 32532.DAT      ID/Wt File: 0104.IDW      Lamp Current: 10  
Technique: HGA      Calib. Type: Method of Add.      Energy: 64  
Remark 1: Matrix modifier: 4% NH4PO4 .2% MoNO3      Integration time: 3 s  
Remark 2: Lamp:      Temp. program:  
Remark 3: Energy:      Current:  
Remark 4: Char. Mass: 20 pob:  
Remark 5: 60 pob:      100 pob:

Pb      ID: Calib. Blank      Seq. No.: 00001      A/S Pos.: 37      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 37  
Replicate 1 (Peak Stored)      Time: 15:56  
Peak Area (A-s): 0.000      Peak Height (A): 0.005  
Background Pk Area (A-s): 0.005      Background Pk Height (A): 0.006  
Blank Corrected Pk Area (A-s): 0.000

Auto-zero performed.

Pb      ID: ICV      Seq. No.: 00002      A/S Pos.: 35      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 35  
Replicate 1 (Peak Stored)      Time: 15:59  
Peak Area (A-s): -0.000      Peak Height (A): 0.003  
Background Pk Area (A-s): 0.004      Background Pk Height (A): 0.004  
Blank Corrected Pk Area (A-s): -0.001

Pb      ID: Calib. Blank      Seq. No.: 00003      A/S Pos.: 37      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 37  
Replicate 1 (Peak Stored)      Time: 16:02  
Peak Area (A-s): 0.000      Peak Height (A): 0.004  
Background Pk Area (A-s): 0.005      Background Pk Height (A): 0.006  
Blank Corrected Pk Area (A-s): 0.000

Auto-zero performed.

Pb      ID: ICV      Seq. No.: 00004      A/S Pos.: 35      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 35  
Replicate 1 (Peak Stored)      Time: 16:05  
Peak Area (A-s): 0.108      Peak Height (A): 0.207  
Background Pk Area (A-s): 0.045      Background Pk Height (A): 0.049  
Blank Corrected Pk Area (A-s): 0.107

Pb      ID: 0104-01      Seq. No.: 00005      A/S Pos.: 1      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 1  
Replicate 1 (Peak Stored)      Time: 16:08  
Peak Area (A-s): 0.927      Peak Height (A): 1.240  
Background Pk Area (A-s): 0.291      Background Pk Height (A): 0.497  
Blank Corrected Pk Area (A-s): 0.927

b. ID: 10 pob Seq. No.: 00006 A/S Pos.: 1 Date: 03/25/93

L dispensed: 18 from 0. 5 from 38. 2 from 36. 20 from 1  
eplicate 1 (Peak Stored) Time: 16:11  
ea. Area (A-s): 0.892 Peak Height (A): 1.167  
ackground Pk Area (A-s): 0.303 Background Pk Height (A): 0.418  
lank Corrected Pk Area (A-s): 0.892

b ID: 30 pob Seq. No.: 00007 A/S Pos.: 1 Date: 03/25/93

L dispensed: 14 from 0. 5 from 38. 6 from 36. 20 from 1  
eplicate 1 (Peak Stored) Time: 16:15  
Peak Area (A-s): 1.002 Peak Height (A): 1.200  
ackground Pk Area (A-s): 0.348 Background Pk Height (A): 0.472  
lank Corrected Pk Area (A-s): 1.002

b ID: 50 pob Seq. No.: 00008 A/S Pos.: 1 Date: 03/25/93

L dispensed: 10 from 0. 5 from 38. 10 from 36. 20 from 1  
eplicate 1 (Peak Stored) Time: 16:18  
Peak Area (A-s): 0.858 Peak Height (A): 1.156  
ackground Pk Area (A-s): 0.309 Background Pk Height (A): 0.391  
lank Corrected Pk Area (A-s): 0.858

Standard abs. & conc. values are not in the same order.

b ID: 0104-01 Seq. No.: 00005 A/S Pos.: 1 Date: 03/25/93

Concentration (ug/L ): -----

b ID: 0104-02 Seq. No.: 00009 A/S Pos.: 2 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 2  
Replicate 1 (Peak Stored) Time: 16:22  
Peak Area (A-s): 0.162 Peak Height (A): 0.133  
Background Pk Area (A-s): 0.084 Background Pk Height (A): 0.053  
Blank Corrected Pk Area (A-s): 0.162



Element File: PBMSA.GEL

Element: Pb

Analyst: ECG

Print Data: Main+Suppl.

Peak Storage: 1 Repl./Sample

Print:

Remarks:

Matrix modifier: 4% NH<sub>4</sub>PO<sub>4</sub> .2% MoNO<sub>3</sub> Integration time: 3 s

Lamp: Temp. program:

Energy: Current:

Clear. Mass: 20 ppb:

60 ppb: 100 ppb:

INSTRUMENT: 4100 ZL

Technique: HGA

Wavelength: 283.3 Peak

Slit: 0.70 Low

Signal Type: Zeeman AA

Signal Measurement: Peak Area

Read Time: 3.0

Read Delay: 0.0

BOC Time: 2

Sample Replicates: 1

Standard Replicates: 1

Spike Replicates: Same as Sample

#### CALIBRATION:

| Solutions    | ID     | Conc | Location | Volume | Diluent |  | Modifier |    |
|--------------|--------|------|----------|--------|---------|--|----------|----|
|              |        |      |          |        | Volume  |  | #1       | #2 |
| Calib. Blank |        |      | 37       | 20     |         |  | 5        |    |
| Standard 1   | 25 ppb | 25.  | 36       | 5      | 15      |  | 5        |    |
| Standard 2   | 50 ppb | 50.  | 36       | 10     | 5       |  | 5        |    |
| Samples      |        |      |          | 20     |         |  | 5        |    |

Diluent Location: 0

Modifier #1 Location: 38

Modifier #2 Location:

Calibration Units: ug/L

Sample Units: ug/L

Calibration Type: Method of Add.

#### Furnace Time/Temperature Program:

| Step | Temp | Ramp | Hold | Gas Flow | Read | Gas Type |
|------|------|------|------|----------|------|----------|
| 1    | 120  | 1    | 40   | 250      |      | Norm     |
| 2    | 140  | 1    | 45   | 250      |      | Norm     |
| 3    | 800  | 7    | 20   | 250      |      | Norm     |
| 4    | 1500 | 0    | 3    | 0        | *    | Norm     |
| 5    | 2400 | 1    | 2    | 250      |      | Norm     |

Injection Temp: 120

Pipette Speed: 100%

Extraction System: On

#### SEQUENCE:

Step Action and Parameters

1 Pipet diluent + modifier 1 + spike + sample/std

2 Run HGA steps 1 to End

CHECKS:  
Recalibration Type: Completely Recalibrate  
Locations: None

Conc.: Above Calibration Action: Dilute & Reanalyze After 1 Reo  
Alternate Sample Volumes (uL): 10.5.4  
Run Alternate Volume Blanks: No

If %RSD > 20.0 and Peak Area > .05 then Retrv 1 times  
Check %RSD on: Samples

Recovery Measurements:  
10 uL of 40 ug/L Standard at Location 40 Gives 20.000 ug/L  
Measure Recovery on Samples: All Samples  
Add to QC Samples: No % Recovery Limits: 85 to 115

QC:

| # | A/S  | QC Sample | Conc. | Limits | After | Periodic | At  | Count As |
|---|------|-----------|-------|--------|-------|----------|-----|----------|
|   | Loc. | ID        | Lower | Upper  | Calib | Check    | End | Sample   |
| 1 | 35   | ICV       | 45    | 55     | X     |          |     |          |
| 2 | 37   | ICB       | -10   | 10     |       |          |     |          |
| 3 | 39   | CRA       | 0     | 10     |       |          |     |          |
| 4 | 35   | CCV       | 45    | 55     |       | X        | X   |          |
| 5 | 37   | CCB       | -10   | 10     |       |          |     |          |

Run Periodic QC Samples: Every 10  
Out of Limit Action: Stop Element and Continue With Next

Matrix Check Calculations:  
% Difference for Duols: No Locations:  
% Recovery for Spike: No Locations: Conc:

~~~~~  
Pb ID: Calib. Blank Seq. No.: 00010 A/S Pos.: 37 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 37
Replicate 1 (Peak Stored) Time: 16:30
Peak Area (A-s): 0.001 Peak Height (A): 0.004
Background Pk Area (A-s): 0.008 Background Pk Height (A): 0.006
Blank Corrected Pk Area (A-s): 0.000

Auto-zero performed.

~~~~~  
Pb ID: ICV Seq. No.: 00011 A/S Pos.: 35 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 35  
Replicate 1 (Peak Stored) Time: 16:33  
Peak Area (A-s): 0.107 Peak Height (A): 0.165  
Background Pk Area (A-s): 0.057 Background Pk Height (A): 0.047  
Blank Corrected Pk Area (A-s): 0.106

~~~~~  
Pb ID: 0104-01 Seq. No.: 00012 A/S Pos.: 1 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 1
Replicate 1 (Peak Stored) Time: 16:36
Peak Area (A-s): 1.085 Peak Height (A): 0.954
Background Pk Area (A-s): 0.299 Background Pk Height (A): 0.255
Blank Corrected Pk Area (A-s): 1.085

Pb ID: 25 oob Seq. No.: 00013 A/S Pos.: 1 Date: 03/25/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 1
Replicate 1 (Peak Stored) Time: 16:39
Peak Area (A-s): 0.862 Peak Height (A): 1.220
Background Pk Area (A-s): 0.319 Background Pk Height (A): 0.483
Blank Corrected Pk Area (A-s): 0.861

Pb ID: 50 oob Seq. No.: 00014 A/S Pos.: 1 Date: 03/25/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 1
Replicate 1 (Peak Stored) Time: 16:42
Peak Area (A-s): 0.895 Peak Height (A): 1.193
Background Pk Area (A-s): 0.305 Background Pk Height (A): 0.408
Blank Corrected Pk Area (A-s): 0.895

Standard abs. & conc. values are not in the same order.

Pb ID: 0104-01 Seq. No.: 00012 A/S Pos.: 1 Date: 03/25/93

Concentration (ug/L): -----

~~~~~  
Pb ID: 0104-02 Seq. No.: 00015 A/S Pos.: 2 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 2  
Replicate 1 (Peak Stored) Time: 16:45  
Peak Area (A-s): 0.084 Peak Height (A): 0.130  
Background Pk Area (A-s): 0.072 Background Pk Height (A): 0.041  
Blank Corrected Pk Area (A-s): 0.083

Pb ID: 25 oob Seq. No.: 00016 A/S Pos.: 2 Date: 03/25/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 2  
Replicate 1 (Peak Stored) Time: 16:49  
Peak Area (A-s): 0.138 Peak Height (A): 0.191  
Background Pk Area (A-s): 0.106 Background Pk Height (A): 0.125  
Blank Corrected Pk Area (A-s): 0.137

Pb ID: 50 oob Seq. No.: 00017 A/S Pos.: 2 Date: 03/25/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 2  
Replicate 1 (Peak Stored) Time: 16:52  
Peak Area (A-s): 0.171 Peak Height (A): 0.293  
Background Pk Area (A-s): 0.101 Background Pk Height (A): 0.146  
Blank Corrected Pk Area (A-s): 0.170

The standard additions calibration curve may not be linear.

Pb ID: 0104-02 Seq. No.: 00015 A/S Pos.: 2 Date: 03/25/93

Concentration (ug/L ): 50.

Dc relation coefficient: 0.99031 Slope: 0.0017 Int: 0.087

Pb ID: 0104-03 Seq. No.: 00018 A/S Pos.: 3 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 3  
Replicate 1 (Peak Stored) Time: 16:55  
Peak Area (A-s): 0.447 Peak Height (A): 0.556  
Background Pk Area (A-s): 0.220 Background Pk Height (A): 0.124  
Blank Corrected Pk Area (A-s): 0.446

Pb ID: 25 oob Seq. No.: 00019 A/S Pos.: 3 Date: 03/25/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 3  
Replicate 1 (Peak Stored) Time: 16:58  
Peak Area (A-s): 0.477 Peak Height (A): 0.642  
Background Pk Area (A-s): 0.189 Background Pk Height (A): 0.208  
Blank Corrected Pk Area (A-s): 0.477

Pb ID: 50 oob Seq. No.: 00020 A/S Pos.: 3 Date: 03/25/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 3  
Replicate 1 (Peak Stored) Time: 17:01  
Peak Area (A-s): 0.464 Peak Height (A): 0.636  
Background Pk Area (A-s): 0.165 Background Pk Height (A): 0.248  
Blank Corrected Pk Area (A-s): 0.464

Standard abs. & conc. values are not in the same order.

Pb ID: 0104-03 Seq. No.: 00018 A/S Pos.: 3 Date: 03/25/93

Concentration (ug/L ): -----

Pb ID: ~~0104-03~~ Seq. No.: 00021 A/S Pos.: 4 Date: 03/25/93

uL dispensed: 5 from 38. 20 from 4  
Replicate 1 (Peak Stored) Time: 17:04  
Peak Area (A-s): 0.118 Peak Height (A): 0.198

3a koround Pk Area (A-s): 0.106      Background Pk Height (A): 0.063  
3l nk Corrected Pk Area (A-s): 0.117  
Pb ID: 25 pob      Seq. No.: 00022      A/S Pos.: 4      Date: 03/25/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 4  
Replicate 1 (Peak Stored)      Time: 17:07  
Peak Area (A-s): 0.175      Peak Height (A): 0.238  
Background Pk Area (A-s): 0.132      Background Pk Height (A): 0.130  
Blank Corrected Pk Area (A-s): 0.174

Pt ID: 50 pob      Seq. No.: 00023      A/S Pos.: 4      Date: 03/25/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 4  
Replicate 1 (Peak Stored)      Time: 17:11  
Peak Area (A-s): 0.221      Peak Height (A): 0.423  
Background Pk Area (A-s): 0.136      Background Pk Height (A): 0.146  
Blank Corrected Pk Area (A-s): 0.220

Pb ID: 0104-04      Seq. No.: 00021      A/S Pos.: 4      Date: 03/25/93

Concentration (ug/L ): 58  
Correlation coefficient: 0.99807      Slope: 0.0021      Int: 0.119

Pb ID: 0104-05      Seq. No.: 00024      A/S Pos.: 5      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 5  
Replicate 1 (Peak Stored)      Time: 17:14  
Peak Area (A-s): 0.908      Peak Height (A): 1.037  
Background Pk Area (A-s): 0.291      Background Pk Height (A): 0.284  
Blank Corrected Pk Area (A-s): 0.907

Pb ID: 25 pob      Seq. No.: 00025      A/S Pos.: 5      Date: 03/25/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 5  
Replicate 1 (Peak Stored)      Time: 17:17  
Peak Area (A-s): 0.977      Peak Height (A): 1.286  
Background Pk Area (A-s): 0.338      Background Pk Height (A): 0.515  
Blank Corrected Pk Area (A-s): 0.978

Pb ID: 50 pob      Seq. No.: 00026      A/S Pos.: 5      Date: 03/25/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 5  
Replicate 1 (Peak Stored)      Time: 17:20  
Peak Area (A-s): 0.922      Peak Height (A): 1.141  
Background Pk Area (A-s): 0.330      Background Pk Height (A): 0.435  
Blank Corrected Pk Area (A-s): 0.921

Standard abs. & conc. values are not in the same order.  
Pb ID: 0104-05      Seq. No.: 00024      A/S Pos.: 5      Date: 03/25/93

Concentration (ug/L ): -----

Pb ID: CCV      Seq. No.: 00027      A/S Pos.: 35      Date: 03/25/93

uL dispensed: 5 from 38. 20 from 35

Replicate 1 (Peak Stored)

Time: 17:23

Peak Area (A-s): 0.107

Peak Height (A): 0.158

Background Pk Area (A-s): 0.126

Background Pk Height (A): 0.069

Blank Corrected Pk Area (A-s): 0.107

~~~~~  
#b ID: Calib. Blank Seq. No.: 00028 A/S Pos.: 37 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 37

~~~~~  
#b ID: Calib. Blank Seq. No.: 00029 A/S Pos.: 37 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 37

Replicate 1 (Peak Stored) Time: 08:39  
Peak Area (A-s): 0.001 Peak Height (A): 0.005  
Background Pk Area (A-s): 0.016 Background Pk Height (A): 0.012  
Blank Corrected Pk Area (A-s): 0.000

Auto-zero performed.

~~~~~  
#b ID: ICV Seq. No.: 00030 A/S Pos.: 35 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 35

Replicate 1 (Peak Stored) Time: 08:42
Peak Area (A-s): 0.000 Peak Height (A): 0.003
Background Pk Area (A-s): 0.009 Background Pk Height (A): 0.008
Blank Corrected Pk Area (A-s): -0.001

~~~~~  
#b ID: EP-8 TV 50 Seq. No.: 00031 A/S Pos.: 1 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 1

Replicate 1 (Peak Stored) Time: 08:45  
Peak Area (A-s): 0.108 Peak Height (A): 0.161  
Background Pk Area (A-s): 0.088 Background Pk Height (A): 0.066  
Blank Corrected Pk Area (A-s): 0.107

#b ID: 25 ppb Seq. No.: 00032 A/S Pos.: 1 Date: 03/26/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 1

Replicate 1 (Peak Stored) Time: 08:48  
Peak Area (A-s): 0.143 Peak Height (A): 0.186  
Background Pk Area (A-s): 0.123 Background Pk Height (A): 0.079  
Blank Corrected Pk Area (A-s): 0.142

#b ID: 50 ppb Seq. No.: 00033 A/S Pos.: 1 Date: 03/26/93


uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 1

Replicate 1 (Peak Stored) Time: 08:51  
Peak Area (A-s): 0.193 Peak Height (A): 0.260  
Background Pk Area (A-s): 0.124 Background Pk Height (A): 0.087  
Blank Corrected Pk Area (A-s): 0.192

#b ID: EP-8 TV 50 Seq. No.: 00031 A/S Pos.: 1 Date: 03/26/93

Concentration (ug/L ): 62.

Correlation coefficient: 0.99565 Slope: 0.0017 Int: 0.105

~~~~~  
#b ID:  Seq. No.: 00034 A/S Pos.: 2 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 2

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.011
Background Pk Area (A-s): 0.152
Blank Corrected Pk Area (A-s): 0.011

Time: 08:54
Peak Height (A): 0.021
Background Pk Height (A): 0.097

Pt ID: 25 ppb Seq. No.: 00035 A/S Pos.: 2 Date: 03/26/93

ul dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 2

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.061
Background Pk Area (A-s): 0.148
Blank Corrected Pk Area (A-s): 0.061

Time: 08:58
Peak Height (A): 0.079
Background Pk Height (A): 0.081

Pb ID: 50 ppb Seq. No.: 00036 A/S Pos.: 2 Date: 03/26/93

ul dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 2

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.112
Background Pk Area (A-s): 0.153
Blank Corrected Pk Area (A-s): 0.112

Time: 09:01
Peak Height (A): 0.154
Background Pk Height (A): 0.098

P ID: 0104-02 1-10 Seq. No.: 00034 A/S Pos.: 2 Date: 03/26/93

Concentration (ug/L): 52. Corrected Conc (ug/L): 52.

Correlation coefficient: 0.99998 Slope: 0.0020 Int: 0.010

P ID: 0104-02 1-50 Seq. No.: 00037 A/S Pos.: 3 Date: 03/26/93

ul dispensed: 5 from 38. 20 from 3

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.005
Background Pk Area (A-s): 0.034
Blank Corrected Pk Area (A-s): 0.004

Time: 09:04
Peak Height (A): 0.016
Background Pk Height (A): 0.021

Pb ID: 25 ppb Seq. No.: 00038 A/S Pos.: 3 Date: 03/26/93

ul dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 3

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.054
Background Pk Area (A-s): 0.110
Blank Corrected Pk Area (A-s): 0.053

Time: 09:07
Peak Height (A): 0.076
Background Pk Height (A): 0.061

P ID: 50 ppb Seq. No.: 00039 A/S Pos.: 3 Date: 03/26/93

ul dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 3

Replicate 1 (Peak Stored)
Peak Area (A-s): 0.105
Background Pk Area (A-s): 0.111
Blank Corrected Pk Area (A-s): 0.104

Time: 09:10
Peak Height (A): 0.157
Background Pk Height (A): 0.074

P ID: 0104-02 1-50 Seq. No.: 00037 A/S Pos.: 3 Date: 03/26/93

Concentration (ug/L): 2. Corrected Conc (ug/L): 90.

Correlation coefficient: 0.99996 Slope: 0.0020 Int: 0.004


~~~~~  
Pb ID: 0104-03 1-10 Seq. No.: 00040 A/S Pos.: 4 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 4  
Replicate 1 (Peak Stored) Time: 09:13  
Peak Area (A-s): 0.075 Peak Height (A): 0.130  
Background Pk Area (A-s): 0.016 Background Pk Height (A): 0.026  
Blank Corrected Pk Area (A-s): 0.074

Pb ID: 25 ppb Seq. No.: 00041 A/S Pos.: 4 Date: 03/26/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 4  
Replicate 1 (Peak Stored) Time: 09:17  
Peak Area (A-s): 0.095 Peak Height (A): 0.201  
Background Pk Area (A-s): 0.131 Background Pk Height (A): 0.086  
Blank Corrected Pk Area (A-s): 0.094

Pb ID: 50 ppb Seq. No.: 00042 A/S Pos.: 4 Date: 03/26/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 4  
Replicate 1 (Peak Stored) Time: 09:20  
Peak Area (A-s): 0.136 Peak Height (A): 0.278  
Background Pk Area (A-s): 0.114 Background Pk Height (A): 0.093  
Blank Corrected Pk Area (A-s): 0.135

Th standard additions calibration curve may not be linear.

Pb ID: 0104-03 1-10 Seq. No.: 00040 A/S Pos.: 4 Date: 03/26/93

Co centration (ug/L ): 57. Corrected Conc (ug/L ): 570.

Correlation coefficient: 0.97998 Slope: 0.0012 Int: 0.070

~~~~~  
Pb ID: ~~0104-03 1-50~~ Seq. No.: 00043 A/S Pos.: 5 Date: 03/26/93

uL dispensed: 5 from 38. 20 from 5
Replicate 1 (Peak Stored) Time: 09:23
Peak Area (A-s): 0.013 Peak Height (A): 0.024
Background Pk Area (A-s): 0.047 Background Pk Height (A): 0.027
Blank Corrected Pk Area (A-s): 0.012

Pb ID: 25 ppb Seq. No.: 00044 A/S Pos.: 5 Date: 03/26/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 5
Replicate 1 (Peak Stored) Time: 09:26
Peak Area (A-s): 0.054 Peak Height (A): 0.097
Background Pk Area (A-s): 0.087 Background Pk Height (A): 0.054
Blank Corrected Pk Area (A-s): 0.053

Pb ID: 50 ppb Seq. No.: 00045 A/S Pos.: 5 Date: 03/26/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 5
Replicate 1 (Peak Stored) Time: 09:29
Peak Area (A-s): 0.094 Peak Height (A): 0.202
Background Pk Area (A-s): 0.098 Background Pk Height (A): 0.075
Blank Corrected Pk Area (A-s): 0.093

Pb ID: 0104-03 1-50 Seq. No.: 00043 A/S Pos.: 5 Date: 03/26/93

Concentration (ug/L): 8

Corrected Conc (ug/L): 378.

Correlation coefficient: 0.99999

Slope: 0.0016

Int: 0.012

Pb ID: CCV

Seq. No.: 00046

A/S Pos.: 35

Date: 03/26/93

uL dispensed: 5 from 38. 20 from 35

Replicate 1 (Peak Stored)

Time: 09:32

Peak Area (A-s): 0.000

Peak Height (A): 0.002

Background Pk Area (A-s): 0.009

Background Pk Height (A): 0.007

Blank Corrected Pk Area (A-s): -0.001

Pb ID: EP-8 TV 50

Seq. No.: 00047

A/S Pos.: 1

Date: 03/26/93

uL dispensed: 5 from 38. 20 from 1

Replicate 1 (Peak Stored)

Time: 09:45

Peak Area (A-s): 0.087

Peak Height (A): 0.184

Background Pk Area (A-s): 0.067

Background Pk Height (A): 0.054

Blank Corrected Pk Area (A-s): -0.086

Pb ID: 25 pob

Seq. No.: 00048

A/S Pos.: 1

Date: 03/26/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 1

Replicate 1 (Peak Stored)

Time: 09:48

Peak Area (A-s): 0.120

Peak Height (A): 0.264

Background Pk Area (A-s): 0.088

Background Pk Height (A): 0.079

Blank Corrected Pk Area (A-s): 0.120

Pb ID: 50 pob

Seq. No.: 00049

A/S Pos.: 1

Date: 03/26/93

uL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 1

Replicate 1 (Peak Stored)

Time: 09:51

Peak Area (A-s): 0.162

Peak Height (A): 0.366

Background Pk Area (A-s): 0.091

Background Pk Height (A): 0.098

Blank Corrected Pk Area (A-s): 0.161

Pb ID: EP-8 TV 50

Seq. No.: 00047

A/S Pos.: 1

Date: 03/26/93

Concentration (ug/L): 56.

Correlation coefficient: 0.99807

Slope: 0.0015

Int: 0.085

Pb ID: 0104-02 1-10

Seq. No.: 00050

A/S Pos.: 2

Date: 03/26/93

uL dispensed: 5 from 38. 20 from 2

Replicate 1 (Peak Stored)

Time: 09:54

Peak Area (A-s): 0.009

Peak Height (A): 0.016

Background Pk Area (A-s): 0.113

Background Pk Height (A): 0.060

Blank Corrected Pk Area (A-s): 0.008

Pb ID: 25 pob

Seq. No.: 00051

A/S Pos.: 2

Date: 03/26/93

uL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 2

Replicate 1 (Peak Stored)

Time: 09:57

Peak Area (A-s): 0.047

Peak Height (A): 0.096

Background Pk Area (A-s): 0.107

Background Pk Height (A): 0.069

Blank Corrected Pk Area (A-s): 0.047

Pb ID: 50 oob Seq. No.: 00052 A/S Pos.: 2 Date: 03/26/93

UL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 2

Replicate 1 (Peak Stored)

Time: 10:01

Peak Area (A-s): 0.930

Peak Height (A): 0.935

Background Pk Area (A-s): 0.277

Background Pk Height (A): 0.276

Blank Corrected Pk Area (A-s): 0.930

Standard additions calibration curve may not be linear.

Pt ID: 0104-02 1-10 Seq. No.: 00050 A/S Pos.: 2 Date: 03/26/93

Concentration (ug/L): -7.

Corrected Conc (ug/L): -72.

Correlation coefficient: 0.88405

Slope: 0.0184

Int: -0.133

Pt ID: 0104-02 1-50 Seq. No.: 00053 A/S Pos.: 3 Date: 03/26/93

UL dispensed: 5 from 38. 20 from 3

Replicate 1 (Peak Stored)

Time: 10:04

Peak Area (A-s): 0.086

Peak Height (A): 0.078

Background Pk Area (A-s): 0.057

Background Pk Height (A): 0.035

Blank Corrected Pk Area (A-s): 0.086

Pb ID: 25 oob Seq. No.: 00054 A/S Pos.: 3 Date: 03/26/93

UL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 3

Replicate 1 (Peak Stored)

Time: 10:07

Peak Area (A-s): 0.055

Peak Height (A): 0.101

Background Pk Area (A-s): 0.082

Background Pk Height (A): 0.053

Blank Corrected Pk Area (A-s): 0.054

P ID: 50 oob Seq. No.: 00055 A/S Pos.: 3 Date: 03/26/93

UL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 3

Replicate 1 (Peak Stored)

Time: 10:10

Peak Area (A-s): 0.082

Peak Height (A): 0.131

Background Pk Area (A-s): 0.080

Background Pk Height (A): 0.054

Blank Corrected Pk Area (A-s): 0.081

Standard abs. & conc. values are not in the same order.

Pt ID: 0104-02 1-50 Seq. No.: 00053 A/S Pos.: 3 Date: 03/26/93

Concentration (ug/L): -----

Pt ID: 0104-03 1-10 Seq. No.: 00056 A/S Pos.: 4 Date: 03/26/93

UL dispensed: 5 from 38. 20 from 4

Replicate 1 (Peak Stored)

Time: 10:13

Peak Area (A-s): 0.056

Peak Height (A): 0.083

Background Pk Area (A-s): 0.078

Background Pk Height (A): 0.047

Blank Corrected Pk Area (A-s): 0.055

Pb ID: 25 oob Seq. No.: 00057 A/S Pos.: 4 Date: 03/26/93

UL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 4

Replicate 1 (Peak Stored) Time: 10:16
Peak Area (A-s): 0.086 Peak Height (A): 0.166
Background Pk Area (A-s): 0.127 Background Pk Height (A): 0.081
Blank Corrected Pk Area (A-s): 0.085

Pb ID: 50 ppb Seq. No.: 00058 A/S Pos.: 4 Date: 03/26/93

µL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 4
Replicate 1 (Peak Stored) Time: 10:20
Peak Area (A-s): 0.128 Peak Height (A): 0.224
Background Pk Area (A-s): 0.093 Background Pk Height (A): 0.072
Blank Corrected Pk Area (A-s): 0.127

The standard additions calibration curve may not be linear.

Pt ID: 0104-03 1-10 Seq. No.: 00056 A/S Pos.: 4 Date: 03/26/93

Concentration (ug/L): 37. Corrected Conc (ug/L): 373.
Correlation coefficient: 0.99559 Slope: 0.0014 Int: 0.054

~~~~~  
Pt ID: 0104-03 1-50 Seq. No.: 00059 A/S Pos.: 5 Date: 03/26/93

µL dispensed: 5 from 38. 20 from 5  
Replicate 1 (Peak Stored) Time: 10:22  
Peak Area (A-s): 0.012 Peak Height (A): 0.025  
Background Pk Area (A-s): 0.043 Background Pk Height (A): 0.023  
Blank Corrected Pk Area (A-s): 0.011

Pb ID: 25 ppb Seq. No.: 00060 A/S Pos.: 5 Date: 03/26/93

µL dispensed: 15 from 0. 5 from 38. 5 from 36. 20 from 5  
Replicate 1 (Peak Stored) Time: 10:26  
Peak Area (A-s): 0.047 Peak Height (A): 0.103  
Background Pk Area (A-s): 0.073 Background Pk Height (A): 0.052  
Blank Corrected Pk Area (A-s): 0.047

Pt ID: 50 ppb Seq. No.: 00061 A/S Pos.: 5 Date: 03/26/93

µL dispensed: 5 from 0. 5 from 38. 10 from 36. 20 from 5  
Replicate 1 (Peak Stored) Time: 10:29  
Peak Area (A-s): 0.079 Peak Height (A): 0.132  
Background Pk Area (A-s): 0.062 Background Pk Height (A): 0.048  
Blank Corrected Pk Area (A-s): 0.078

Pt ID: 0104-03 1-50 Seq. No.: 00059 A/S Pos.: 5 Date: 03/26/93

Concentration (ug/L ): 9. Corrected Conc (ug/L ): 429.  
Correlation coefficient: 0.99923 Slope: 0.0013 Int: 0.012