



**Stantec Consulting Services Inc.**  
61 Commercial Street, Suite 100  
Rochester NY 14614-1009  
Tel: (585) 475-1440  
Fax: (585) 272-1814

June 10, 2015  
File: 190500751

Todd Caffoe, P.E  
New York State Department of Environmental Conservation  
Division of Environmental Remediation  
6274 East Avon-Lima Road  
Avon, NY 14414

**Reference:   Brownfield Cleanup Program  
                  Progress Report #23  
                  Site #C828184  
                  Former Carriage Factory  
                  33 Litchfield Street  
                  Rochester, Monroe County, New York**

Dear Todd,

On behalf of Carriage Factory Special Needs Apartments, LP (CFSNA), Stantec Consulting Services Inc. (Stantec) has prepared this Progress Report #23 for the Brownfield Cleanup Program (BCP) at the Former Carriage Factory located at 33 Litchfield Street in the City of Rochester, Monroe County, New York (Site). This report covers activities that took place since the submission of Progress Report #22 (dated March 10, 2015).

## **1. Actions Completed During The Reporting Period**

- On May 4 and 5, twelve groundwater monitoring wells (see Figure 1) were sampled as outlined in the Enhanced Reductive Dechlorination IRM Work Plan. The attached Table 1 summarizes the recorded groundwater field parameters for this event, and previous injection monitoring events. The field data indicate very anaerobic, reducing conditions [dissolved oxygen concentrations (DO) lower than 0.5 mg/L and oxidation-reduction potential (ORP) values generally lower than -150 mV] remain present in each of the wells that received injections of sodium lactate solution. Table 1 includes a chart showing the decreasing trend in average ORP values at these wells over time.

Exceptions to the decreasing ORP trend during this event include RW-4, RW-9, and RW-12. RW-4 was sampled by bailer due to poor groundwater recharge. Since the sample DO was not collected downhole, it is likely that this DO field measurement is biased high due to potential aeration during sample collection. A DO measurement collected during the initial attempt to use low-flow sampling, prior to volumetric bailer purging, was around 0.5 mg/L. This value is more in line with historical DO values for this well, and, therefore, is considered to be more representative of actual subsurface conditions within the vicinity of RW-4.



**Reference: Brownfield Cleanup Program  
Progress Report #23  
Site #C828184  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY**

RW-9 and RW-12 are both off-site monitoring wells and not located within the active treatment zone targeted by the sodium lactate injection activities. RW-9 is a side-gradient well where groundwater contaminant concentrations have been below groundwater standards for the last six sampling events. The DO and ORP values were both outside the performance metrics stated above, but still show anaerobic ( $\text{DO} < 2.0 \text{ mg/L}$ ) and reducing (negative ORP) conditions exist at this well location. RW-12 is located down gradient of the targeted treatment area and has consistently shown anaerobic and reducing conditions to be present in the vicinity of this well. The recorded ORP value ( $-113 \text{ mV}$ ) represents an increase since the previous sampling event, but desired geochemical conditions are still present at this location.

- On May 5, the quarterly elevator sump sample was collected, as required by MCDES to satisfy Sewer-Use Permit #996.
- On May 5, the purge water generated during the May 4 and 5 groundwater sampling event was sampled for discharge to the municipal sewer, pending approval by MCDES.
- On May 5, measurements were recorded from the three vacuum manometers installed in the fifth floor utility room to continuously monitor the vacuum induced by the three SSDS fans. The vacuum readings ranged from 2.0 to 2.2 inches  $\text{H}_2\text{O}$ . These readings were consistent with the readings recorded when the system was installed and monthly readings collected by Marc Banning, Property Manager for DePaul Properties, which ranged from 2.0 to 2.4 inches  $\text{H}_2\text{O}$  from January to March 2015.
- On May 22, the purge water generated during the May 2015 sampling event (approximately 27 gallons) was discharged to the municipal sewer, as approved by MCDES (see below).

## **2. Data Received or Generated in the Previous Reporting Period**

- Analytical results from the quarterly groundwater sampling event performed on May 4 and 5 are included in the attached Table 2, Figure 2 and in Appendix A. These results reflect groundwater conditions one year following injection of sodium lactate and indicate that the parent volatile organic compounds (VOCs) of tetrachloroethylene (PCE), and trichloroethylene (TCE) continue to degrade into the daughter compounds of the cis and trans isomers of 1,2-dichloroethylene (1,2-DCE) and vinyl chloride (VC), before proceeding to complete destruction. Of the twelve wells sampled, nine had both PCE and TCE concentrations below groundwater standards ( $5 \text{ ug/L}$ ). The parent compounds were also not detected at RW-6, but slightly elevated detection limits of  $10.0 \text{ ug/L}$  were due to elevated concentrations of cis-1,2-DCE which may have masked the potential presence of PCE and TCE between 5 and  $10 \text{ ug/L}$ . PCE concentrations above  $5 \text{ ug/L}$  were detected at B108-MW ( $14.4 \text{ ug/L}$ ) and RW-4 ( $18.7 \text{ ug/L}$ ), while RW-4 also had TCE concentrations greater than groundwater standards at  $8.94 \text{ ug/L}$ .



**Reference: Brownfield Cleanup Program  
Progress Report #23  
Site #C828184  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY**

In general, daughter compound concentrations have continued to decline at all wells from peak levels observed within the first 6 months of the remediation program. The highest groundwater concentrations of daughter compounds were observed at well RW-6 (373 ug/L cis-1,2-DCE, 6.18 ug/L trans-1,2-DCE, and 367 ug/L VC), but levels still represented a decline from the previous sampling event. These favorable results continue to indicate improving groundwater conditions as a result of the injection program. However, these data should be considered preliminary as they have not yet undergone data validation.

- Decreasing total organic carbon (TOC) concentrations at all well locations indicate that the injected sodium lactate continues to be consumed. Groundwater TOC concentrations ranged between 1.5 mg/L (B106-MW) and 6.2 mg/L (RW-12). These low levels may be an indication of the potential need to consider an additional, targeted sodium lactate injection event in the future.
- Gene-trac assay results from the groundwater samples obtained from B102-MW and B108-MW during the May sampling event are included in Appendix A. Results indicate that moderate population levels of dehalococcoides have been maintained, despite a slight decline since the last sampling event.
- Analytical results from the purge water sampled on May 5 are presented in Table 3. These results, which were well below discharge limits, were forwarded to Monroe County on May 19, and the water was approved for discharge to the municipal sewer on May 20.
- Analytical results from the sump water sampled on May 5 are also presented in Table 3. These results, which were well below discharge limits, were forwarded to Monroe County to fulfill the permit requirements.
- A Data Usability Summary Report (DUSR), dated May 22, 2015, was prepared for the analytical data from the October 2014 and February 2015 groundwater sampling events and is included in Appendix B. The attached Table 2 reflects the qualified data as a result of the DUSR findings.

### **3. Deliverables Completed and Submitted during the Previous Reporting Period**

- Progress Report No. 22 was submitted on March 10.

### **4. Actions Scheduled for the Next Reporting Period**

- Prepare a Data Usability Summary Report (DUSR) for the analytical data from the May 2015 groundwater sampling event.
- Submit an EDD for the May 2015 groundwater data.



Todd Caffoe, P.E.  
Page 4

**Reference: Brownfield Cleanup Program  
Progress Report #23  
Site #C828184  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY**

- It is proposed to conduct another round of quarterly groundwater sampling in early August 2014 to determine if contaminant levels continue to decline and assess both the need for a potential supplemental sodium lactate injection and the frequency and extent of the associated groundwater monitoring program. At this time, we propose to discontinue sampling at RW-9, as groundwater contaminant concentrations at this well have been below groundwater standards for the last six sampling events.

### Closing

If you have any questions or require further information, please contact us.

Regards,

**STANTEC CONSULTING SERVICES INC.**

Michael P. Storonsky  
Managing Principal  
Phone: 585-413-5266  
[mike.storonsky@stantec.com](mailto:mike.storonsky@stantec.com)

Peter Nielsen, P.E.  
Senior Associate  
Phone: 585-413-5280  
[peter.nielsen@stantec.com](mailto:peter.nielsen@stantec.com)

### Attachments:

Figure 1 – Site Map  
Figure 2 – Summary of CVOC Degradation Over Time – All Wells  
Table 1 – Summary of Groundwater Field Parameters  
Table 2 – Summary of Analytical Results in Groundwater (Preliminary)  
Table 3 – Summary of Analytical Results in Waste Water and Discharge Permit Samples  
Appendix A – Laboratory Analytical Reports  
Appendix B – Data Usability Summary Report (DUSR) – October 2014 and February 2015

ec:

James Mahoney (NYSDEC)  
Al Floro (Nixon Peabody)  
Justin Deming (NYSDOH)  
Jonathan Penna (Nixon Peabody)  
Stephanie Selmer (NYSDOH)  
Mark Gregor (City of Rochester)

James Whalen (CFSNA)  
Eleonora Bershadskaya (Goldman Sachs)  
Mark Fuller (CFSNA)  
Daniel Alger (Goldman Sachs)  
Gillian Conde (CFSNA)  
Linda Kaiser (Goldman Sachs)

Joy Cromwell (CFSNA)  
Patrick Miller (CPC)  
Chris Betts (Betts Housing)  
David Lent (IVI)  
Amy Reichhart (Nixon Peabody)



**Stantec Consulting Services Inc.**  
61 Commercial Street, Suite 100  
Rochester NY 14614-1009  
Tel: (585) 475-1440  
Fax: (585) 272-1814

## Figures



Copyright Reserved

The Contractor shall verify and be responsible for all dimensions. DO NOT scale the drawing - any errors or omissions shall be reported to Stantec without delay.

The Copyrights to all designs and drawings are the property of Stantec. Reproduction or use for any purpose other than that authorized by Stantec is forbidden.


## Consultants

### Legend



RW 2  
888 ug/l

LOCATIONS OF PHASE II ESA GROUNDWATER MONITORING WELLS



RI MONITORING WELL LOCATION

PROPERTY LINE

SANITARY SEWER (SA MH=MANHOLE)

## Notes

1. PLAN ADAPTED FROM BASE PLAN BY PARRONE ENGINEERING

|                 |      |       |          |
|-----------------|------|-------|----------|
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
| Revision        | By   | Appd. | YY.MM.DD |
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
|                 |      |       |          |
| Progress Report | RM   | MPS   | 15.06.10 |
| Issued          | By   | Appd. | YY.MM.DD |
|                 |      |       |          |
|                 |      |       |          |
| File Name:      |      |       |          |
|                 |      |       |          |
|                 | Dwn. | Chkd. | Dsgn.    |
| Permit-Seal     |      |       | YY.MM.DD |

Client/Project

CARRIAGE FACTORY SPECIAL NEEDS APARTMENTS, L.P.

**BROWNFIELD CLEANUP PROGRAM**  
**FORMER CARRIAGE FACTORY**  
33 LITCHFIELD STREET, ROCHESTER, NY

Title

## SITE PLAN

Project No.  
190500751

Scale  
AS SHOWN

Drawing No.


Shee

Revision

FIGURE 1

of


C



U:\190500751\drawing\CAD\Progress Reports\June 2015\Figure 1 Site Plan .dwg  
2015/06/05 5:52 PM By: Less. Andy

ORIGINAL SHEET - ANSI I

**Figure 2**  
**Summary of CVOC Degradation Over Time - All Wells**  
 Former Carriage Factory  
 33 Litchfield Street, Rochester, NY





**Stantec Consulting Services Inc.**  
61 Commercial Street, Suite 100  
Rochester NY 14614-1009  
Tel: (585) 475-1440  
Fax: (585) 272-1814

## Tables

Table 1  
Summary of Groundwater Field Parameters  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Sample Location               |       | B101-MW     |             | B102-MW     |             |             |             |             |             |             |  | B106-MW     |             |             |             |             |             |             |             |
|-------------------------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Purge Date                    |       | 21-May-13   | 22-May-13   | 27-Mar-14   | 28-May-14   | 2-Jul-14    | 6-Aug-14    | 28-Oct-14   | 3-Feb-15    | 4-May-15    |  | 23-May-13   | 26-Mar-14   | 28-May-14   | 2-Jul-14    | 7-Aug-14    | 28-Oct-14   | 3-Feb-15    | 5-May-15    |
| Purge Methodology             |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |  | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |
| Purge Method                  |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Sample Date                   |       | 21-May-13   | 22-May-13   | 27-Mar-14   | 28-May-14   | 2-Jul-14    | 6-Aug-14    | 28-Oct-14   | 3-Feb-15    | 4-May-15    |  | 23-May-13   | 26-Mar-14   | 28-May-14   | 2-Jul-14    | 7-Aug-14    | 28-Oct-14   | 3-Feb-15    | 5-May-15    |
| Sampling Method               |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Field Parameters              |       | Units       |             |             |             |             |             |             |             |             |  |             |             |             |             |             |             |             |             |
| Conductivity                  | mS/cm | 0.99        | 0.86        | 0.90        | 0.92        | 1.41        | 1.03        | 1.15        | 1.19        | 1.28        |  | 0.92        | 1.08        | 1.29        | 2.20        | 1.30        | 1.09        | 1.06        | 1.03        |
| Dissolved Oxygen              | mg/L  | 0.10        | 0.10        | 0.12        | 0.19        | 0.14        | 0.03        | 0.10        | 0.00        | 0.20        |  | 0.13        | 0.07        | 0.08        | 0.17        | 0.11        | 0.40        | 0.00        | 0.06        |
| Oxidation Reduction Potential | mV    | -25.0       | 13.3        | 73.6        | -49.7       | -271.6      | -284.0      | -118.9      | -154.7      | -233.3      |  | 17.8        | 90.8        | -96.3       | -231.4      | -274.4      | -138.8      | -172.9      | -241.4      |
| pH                            | S.U.  | 7.02        | 6.87        | 7.02        | 7.15        | 7.26        | 7.04        | 7.06        | 7.17        | 7.00        |  | 6.99        | 7.05        | 7.15        | 6.96        | 7.07        | 7.02        | 7.09        | 6.98        |
| Temperature                   | deg C | 13.4        | 20.5        | 3.7         | 18.4        | 16.2        | 20.4        | 15.9        | 7.7         | 10.9        |  | 16.1        | 3.0         | 18.3        | 15.7        | 16.5        | 15.4        | 16.2        | 16.7        |
| Turbidity                     | NTU   | 0.68        | 4.07        | 11.71       | 1.87        | 1.79        | 1.45        | 2.75        | 2.28        | 0.76        |  | 4.77        | 1.84        | 1.48        | 1.46        | 2.1         | 2.46        | 0.99        | 0.48        |
| Volume Purged                 | gal   | 0.8         | 1.2         | 0.5         | 2.6         | 2.0         | 2.0         | 0.7         | 0.5         | 1.8         |  | 1.1         | 0.7         | 1.8         | 1.5         | 1.7         | 1.4         | 1.1         | 1.7         |

| Sample Location               |       | B108-MW     |             |             |             |             |             |             |             |       | RW-1        |             |             |             |             |             |             |             |  |
|-------------------------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--|
| Purge Date                    |       | 23-May-13   | 26-Mar-14   | 28-May-14   | 2-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |       | 23-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  |
| Purge Methodology             |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |  |
| Purge Method                  |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |
| Sample Date                   |       | 23-May-13   | 26-Mar-14   | 28-May-14   | 2-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |       | 23-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  |
| Sampling Method               |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |
| Field Parameters              |       | Units       |             |             |             |             |             |             |             |       |             |             |             |             |             |             |             |             |  |
| Conductivity                  | mS/cm | 0.95        | 1.06        | 1.05        | 1.27        | 1.22        | 1.22        | 1.49        | 1.04        | 0.74  |             | 1.07        | 1.22        | 2.12        | 1.15        | 1.23        | 1.13        | 1.82        |  |
| Dissolved Oxygen              | mg/L  | 0.13        | 0.13        | 0.10        | 0.18        | 0.13        | 0.31        | 0.00        | 0.06        | 0.13  |             | 0.01        | 0.11        | 0.08        | 0.14        | 0.70        | 0.00        | 0.01        |  |
| Oxidation Reduction Potential | mV    | 29.1        | 137.1       | -69.9       | -216.0      | -293.4      | -354.1      | -327.4      | -241.5      | -94.3 |             | 179.0       | -147.8      | -252.9      | -313.0      | -297.2      | -321.0      | -266.7      |  |
| pH                            | S.U.  | 7.15        | 7.04        | 7.21        | 7.04        | 7.02        | 7.08        | 7.16        | 7.09        | 7.19  |             | 7.05        | 7.68        | 7.67        | 7.05        | 7.36        | 7.17        | 7.03        |  |
| Temperature                   | deg C | 13.6        | 10.6        | 19.5        | 16.1        | 15.4        | 16.0        | 16.7        | 16.1        | 12.5  |             | 8.6         | 18.8        | 16.5        | 15.0        | 15.3        | 15.2        | 15.3        |  |
| Turbidity                     | NTU   | 0.62        | 0.28        | 3.54        | 0.86        | 3.78        | 3.24        | 1.11        | 1.56        | 10.55 |             | 12.37       | 1.66        | 6.31        | 3.19        | 4.41        | 2.97        | 2.15        |  |
| Volume Purged                 | gal   | 0.5         | 0.7         | 1.8         | 1.1         | 1.55        | 1.7         | 0.7         | 1.8         | 0.7   |             | 0.7         | 1.5         | 1.4         | 1.8         | 0.9         | 1.2         | 2.3         |  |

| Sample Location               |       | RW-2        |             |             |             |             |             |             |             |  | RW-3        |             |             |             |             |             |             |             |  |  |  |
|-------------------------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--|--|--|
| Purge Date                    |       | 21-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  | 22-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 7-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  |  |  |
| Purge Methodology             |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |  | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |  |  |  |
| Purge Method                  |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |  |  |
| Sample Date                   |       | 21-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 8-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  | 22-May-13   | 26-Mar-14   | 29-May-14   | 1-Jul-14    | 7-Aug-14    | 29-Oct-14   | 3-Feb-15    | 5-May-15    |  |  |  |
| Sampling Method               |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |  |  |
| Field Parameters              |       | Units       |             |             |             |             |             |             |             |  |             |             |             |             |             |             |             |             |  |  |  |
| Conductivity                  | mS/cm | 0.85        | 1.08        | 2.34        | 1.70        | 1.68        | 1.27        | 1.27        | 1.03        |  | 0.87        | 1.09        | 1.79        | 1.31        | 1.00        | 1.05        | 1.23        | 1.22        |  |  |  |
| Dissolved Oxygen              | mg/L  | 0.28        | 0.03        | 0.20        | 0.11        | 0.16        | 0.65        | 0.11        | 0.08        |  | 0.15        | 0.06        | 0.08        | 0.06        | 0.23        | 0.37        | 0.00        | 0.10        |  |  |  |
| Oxidation Reduction Potential | mV    | -30.3       | 156.8       | -171.5      | -172.0      | -292.5      | -286.4      | -152.2      | -326.1      |  | 87.3        | 157.6       | -132.8      | -213.0      | -216.8      | -242.2      | -192.4      | -320.7      |  |  |  |
| pH                            | S.U.  | 7.36        | 7.11        | 6.94        | 7.56        | 6.93        | 7.52        | 7.61        | 7.09        |  | 7.39        | 7.07        | 7.45        | 7.67        | 7.35        | 7.71        | 7.48        | 7.20        |  |  |  |
| Temperature                   | deg C | 12.7        | 7.2         | 16.8        | 16.8        | 14.9        | 16.0        | 15.6        | 16.2        |  | 12.4        | 9.3         | 17.7        | 15.3        | 15          | 15.7        | 16.3        | 17.2        |  |  |  |
| Turbidity                     | NTU   | 5.23        | 3.81        | 7.53        | 2.34        | 1.71        | 3.71        | 2.92        | 1.45        |  | 0.88        | 1.29        | 1.24        | 1.72        | 1.62        | 2.42        | 2.62        | 0.48        |  |  |  |
| Volume Purged                 | gal   | 1.2         | 0.8         | 1.4         | 0.3         | 1.15        | 0.6         | 1.0         | 1.0         |  | 0.5         | 0.7         | 1.5         | 1.8         | 0.5         | 0.6         | 0.7         | 1.6         |  |  |  |

| Sample Location               |       | RW-4        |             |             |             |             |             |             |                         |      | RW-5        |             |             |             |             |             |             |             |  |
|-------------------------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------------------|------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|--|
| Purge Date                    |       | 22-May-13   | 26-Mar-14   | 29-May-14   | 2-Jul-14    | 6-Aug-14    | 29-Oct-14   | 4-Feb-15    | 4-May-15                |      | 21-May-13   | 27-Mar-14   | 29-May-14   | 2-Jul-14    | 7-Aug-14    | 28-Oct-14   | 3-Feb-15    | 4-May-15    |  |
| Purge Methodology             |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Volumetric <sup>d</sup> |      | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |  |
| Purge Method                  |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Bailer <sup>d</sup>     |      | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |
| Sample Date                   |       | 22-May-13   | 26-Mar-14   | 29-May-14   | 2-Jul-14    | 6-Aug-14    | 29-Oct-14   | 4-Feb-15    | 4-May-15                |      | 21-May-13   | 27-Mar-14   | 29-May-14   | 2-Jul-14    | 7-Aug-14    | 28-Oct-14   | 3-Feb-15    | 4-May-15    |  |
| Sampling Method               |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Bailer <sup>d</sup>     |      | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |  |
| Field Parameters              |       | Units       |             |             |             |             |             |             |                         |      |             |             |             |             |             |             |             |             |  |
| Conductivity                  | mS/cm | 0.91        | 0.88        | 0.89        | 1.94        | 1.67        | 1.00        | 1.48        | 1.77                    | 0.89 | 1.08        | 1.40        | 1.86        | 1.20        | 1.01        | 1.00        | 1.02        |             |  |
| Dissolved Oxygen              | mg/L  | 0.11        | 0.17        | 0.06        | 0.15        | 0.04        | 0.44        | 0.24        | 3.04                    | 0.28 | 0.00        | 0.06        | 0.19        | 0.08        | 0.43        | 0.00        | 0.39        |             |  |
| Oxidation Reduction Potential | mV    | 38.6        | 132.4       | 29.3        | -180.2      | -347        | -130.3      | -278.2      | -175.5                  | -2.3 | 74.7        | -95.6       | -137.8      | -170.0      | -164.1      | -269.2      | -230.4      |             |  |
| pH                            | S.U.  | 6.91        | 7.08        | 7.10        | 6.90        | 7.05        | 6.95        | 7.17        | 7.10                    | 7.07 | 7.29        | 7.27        | 7.03        | 7.07        | 7.23        | 7.31        | 7.19        |             |  |
| Temperature                   | deg C | 20.0        | 2.4         | 25.5        | 17.4        | 19.2        | 14.8        | 7.4         | 9.8                     | 16.2 | 5.7         | 22.8        | 17.3        | 19.9        | 17.5        | 5.2         | 11.6        |             |  |
| Turbidity                     | NTU   | 5.68        | 5.81        | 1.72        | 3.18        | 1.93        | 1.06        | 2.01        | 10.25                   | 2.98 | 1.22        | 7.10        | 1.88        | 3.89        | 1.77        | 3.60        | 6.21        |             |  |
| Volume Purged                 | gal   | 0.8         | 1.8         | 0.9         | 1.9         | 1.1         | 2.1         | 2.7         | 5.4                     | 1.1  | 3.2         | 0.5         | 1.2         | 1.5         | 0.8         | 1.4         | 0.9         |             |  |

See last page for Notes.

Table 1  
Summary of Groundwater Field Parameters  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Sample Location               |       | RW-6              |             |             |             |             |             |             |             | RW-7        |             |             |             |             |             |             |             |
|-------------------------------|-------|-------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Purge Date                    |       | 20-May-13         | 27-Mar-14   | 28-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    | 20-May-13   | 27-Mar-14   | 28-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    |
| Purge Methodology             |       | Low flow          | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |
| Purge Method                  |       | Peristaltic       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Sample Date                   |       | 20-May-13         | 27-Mar-14   | 28-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    | 20-May-13   | 27-Mar-14   | 28-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    |
| Sampling Method               |       | Peristaltic       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Field Parameters              |       | Units             |             |             |             |             |             |             |             |             |             |             |             |             |             |             |             |
| Conductivity                  | mS/cm | 0.93              | 1.07        | 1.72        | 1.34        | 1.30        | 1.21        | 1.08        | 1.01        | 1.02        | 1.21        | 1.30        | 1.17        | 1.07        | 0.96        | 1.16        | 1.08        |
| Dissolved Oxygen              | mg/L  | 0.08              | 0.01        | 0.07        | 0.10        | 0.14        | 0.42        | 0.28        | 0.08        | 0.08        | 0.38        | 0.31        | 0.13        | 0.11        | 0.44        | 0.39        | 0.07        |
| Oxidation Reduction Potential | mV    | -10.6             | 138.3       | -69.0       | -136.7      | -306.1      | -134.8      | -304.1      | -252.4      | 29.4        | 92.6        | -37.6       | -104.6      | -303.6      | -168.2      | -224.3      | -208.5      |
| pH                            | S.U.  | 7.13              | 7.33        | 7.03        | 6.91        | 7.00        | 7.06        | 7.22        | 7.14        | 7.06        | 7.27        | 7.08        | 6.99        | 7.07        | 7.11        | 7.12        | 6.99        |
| Temperature                   | deg C | 19.0              | 6.1         | 17.6        | 21.2        | 17.2        | 16.7        | 6.8         | 10.4        | 16.8        | 6.7         | 20.3        | 18.4        | 16.3        | 17.5        | 7.9         | 10.6        |
| Turbidity                     | NTU   | 7.08 <sup>a</sup> | 5.46        | 7.48        | 4.83        | 4.79        | 1.03        | 4.76        | 4.62        | 10.38       | 1.36        | 3.12        | 1.12        | 1.53        | 4.74        | 0.67        | 1.77        |
| Volume Purged                 | gal   | 1.3               | 1.1         | 1.2         | 0.7         | 1.0         | 0.7         | 1.2         | 2.0         | 1.2         | 0.9         | 1.8         | 1.2         | 1.5         | 1.3         | 2.0         | 1.8         |

| Sample Location               |       | RW-8        |             |             | RW-9        |             |             |             |             |             | RW-11             |             |  |
|-------------------------------|-------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------------|-------------|--|
| Purge Date                    |       | 20-May-13   | 21-May-13   | 27-Mar-14   | 29-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    | 22-May-13         | 27-Mar-14   |  |
| Purge Methodology             |       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow          | Low flow    |  |
| Purge Method                  |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic       | Peristaltic |  |
| Sample Date                   |       | 20-May-13   | 21-May-13   | 27-Mar-14   | 29-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14   | 4-Feb-15    | 4-May-15    | 22-May-13         | 27-Mar-14   |  |
| Sampling Method               |       | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic       | Peristaltic |  |
| Field Parameters              |       | Units       |             |             |             |             |             |             |             |             |                   |             |  |
| Conductivity                  | mS/cm | 1.04        | 0.94        | 1.05        | 0.68        | 0.74        | 0.85        | 0.98        | 1.03        | 0.97        | 0.79              | 0.82        |  |
| Dissolved Oxygen              | mg/L  | 1.06        | 2.48        | 2.45        | 5.52        | 2.37        | 2.43        | 0.50        | 0.45        | 0.61        | 2.36              | 1.62        |  |
| Oxidation Reduction Potential | mV    | 77.0        | 49.4        | 104.6       | 28.1        | 33.9        | 51.0        | 4.1         | -166.7      | -34.3       | 94.5              | 88.8        |  |
| pH                            | S.U.  | 7.05        | 7.13        | 7.29        | 7.44        | 7.12        | 7.06        | 7.04        | 7.12        | 6.99        | 7.15              | 7.33        |  |
| Temperature                   | deg C | 14.4        | 14.0        | 9.4         | 20.7        | 19.0        | 15.5        | 16.8        | 10.5        | 15.2        | 14.6              | 5.1         |  |
| Turbidity                     | NTU   | 2.54        | 0.33        | 0.50        | 3.62        | 1.80        | 1.06        | 1.61        | 0.71        | 2.88        | 0.11 <sup>b</sup> | 1.31        |  |
| Volume Purged                 | gal   | 1.0         | 0.8         | 1.2         | 0.7         | 0.35        | 0.7         | 2.9         | 1.5         | 1.6         | 0.4               | 0.7         |  |

| Sample Location               |       | RW-12          |             |             |             |             |             | RW-13       |             |             |
|-------------------------------|-------|----------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Purge Date                    |       | 20-May-13      | 28-May-14   | 2-Jul-14    | 7-Aug-14    | 29-Oct-14   | 4-Feb-15    | 4-May-15    | 20-May-13   | 27-Mar-14   |
| Purge Methodology             |       | Low flow       | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    | Low flow    |
| Purge Method                  |       | Peristaltic    | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Sample Date                   |       | 20-May-13      | 28-May-14   | 2-Jul-14    | 7-Aug-14    | 29-Oct-14   | 4-Feb-15    | 4-May-15    | 20-May-13   | 27-Mar-14   |
| Sampling Method               |       | Peristaltic    | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic | Peristaltic |
| Field Parameters              |       | Units          |             |             |             |             |             |             |             |             |
| Conductivity                  | mS/cm | 1.02           | 1.76        | 2.09        | 2.00        | 1.60        | 1.37        | 1.49        | 1.08        | 1.12        |
| Dissolved Oxygen              | mg/L  | 0.06           | 0.06        | 0.24        | 0.45        | 1.02        | 0.34        | 0.09        | 1.96        | 2.13        |
| Oxidation Reduction Potential | mV    | 20.0           | -149.5      | -204.6      | -159.7      | -44.7       | -284.1      | -113.1      | 48.6        | 101.8       |
| pH                            | S.U.  | 7.10           | 7.25        | 7.11        | 7.17        | 7.30        | 7.36        | 7.40        | 7.21        | 7.25        |
| Temperature                   | deg C | 16.0           | 24.1        | 17.4        | 18.1        | 14.8        | 6.8         | 12.4        | 17.2        | 6.0         |
| Turbidity                     | NTU   | — <sup>c</sup> | 1.10        | 5.55        | 2.82        | 2.45        | 1.40        | 0.61        | 5.10        | 1.86        |
| Volume Purged                 | gal   | 1.0            | 2.0         | 0.9         | 1.3         | 0.6         | 1.7         | 2.9         | 2.3         | 2.0         |

Notes:

- deg c
- degrees Celsius
- gal
- gallons
- mg/l
- milligrams per liter
- mS/cm
- milliSiemens per centimeter
- mV
- millivolts
- NTU
- nephelometric turbidity unit
- AU
- attenuation unit (equivalent to NTU)
- S.U.
- standard units
- <sup>a</sup>
- Sample turbidity measured approximately 10 minutes prior to sampling; subsequent measurements (-126 NTU) indicated that the turbidity meter was not functioning.
- <sup>b</sup>
- Sample turbidity measured approximately 5 minutes prior to sampling; subsequent measurement (-0.02 NTU) indicated that the turbidity meter was not functioning.
- <sup>c</sup>
- Turbidity meter was not functioning; groundwater was clear and did not have an odor.
- <sup>d</sup>
- Due to a large drop in water level, RW-4 was purged and sampled by bailer (during the May 2015 event); parameters provided were not measured downhole.

Parameter Average for All Wells  
Pre - Post Injection Comparison

|                  | Baseline | Month 1<br>Post Injection | Month 2<br>Post Injection | Month 3<br>Post Injection | Month 6<br>Post Injection | Month 9<br>Post Injection | Month 12<br>Post Injection |
|------------------|----------|---------------------------|---------------------------|---------------------------|---------------------------|---------------------------|----------------------------|
| Parameter        | Mar-14   | May-14                    | Jul-14                    | Aug-14                    | Oct-14                    | Feb-15                    | May-15                     |
| Conductivity     | 1.04     | 1.36                      | 1.60                      | 1.29                      | 1.15                      | 1.21                      | 1.23                       |
| Dissolved Oxygen | 0.55     | 0.57                      | 0.33                      | 0.34                      | 0.56                      | 0.15                      | 0.40                       |
| ORP              | 117.55   | -80.19                    | -173.91                   | -242.46                   | -172.97                   | -237.27                   | -220.33                    |
| pH               | 7.17     | 7.19                      | 7.11                      | 7.07                      | 7.20                      | 7.29                      | 7.09                       |
| Temperature      | 6.45     | 20.04                     | 17.28                     | 16.95                     | 16.03                     | 11.03                     | 13.53                      |
| Turbidity        | 3.76     | 3.46                      | 2.74                      | 2.49                      | 2.64                      | 2.17                      | 2.77                       |
| Volume Purged    | 1.15     | 1.48                      | 1.20                      | 1.32                      | 1.19                      | 1.31                      | 2.07                       |




Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area<br>Sample Location<br>Sample Date<br><br>Sample ID<br><br>Sampling Company<br>Laboratory<br>Laboratory Work Order<br>Laboratory Sample ID<br><br>Sample Type |          |                    | On-Site Parking Lot |                    |               |                  |                  |                    |                     |                      |                    |                    |                    |                    |                    |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------|--------------------|---------------------|--------------------|---------------|------------------|------------------|--------------------|---------------------|----------------------|--------------------|--------------------|--------------------|--------------------|--------------------|
|                                                                                                                                                                   |          |                    | B101MW              |                    |               |                  |                  |                    | B102MW              |                      |                    |                    |                    |                    |                    |
|                                                                                                                                                                   |          |                    | 21-May-13           | 21-May-13          | 22-May-13     | 27-Mar-14        | 27-Mar-14        | 28-May-14          | 2-Jul-14            | 6-Aug-14             | 28-Oct-14          | 3-Feb-15           | 3-Feb-15           | 4-May-15           | 4-May-15           |
|                                                                                                                                                                   |          |                    | LI-B101MW-GW1       | LI-B101MW-GW1DUP   | LI-B102MW-GW1 | LI-B102-MW       | LI-DUP-MW        | LI-B102-MW-PI1     | LI-B102-MW-PI2      | LI-B102-MW-PI3       | LI-B102-MW-PI6     | LI-B102-MW-PI9     | LI-DUP-PI9         | LI-B102-MW-PI12    | LI-DUP-PI12        |
|                                                                                                                                                                   |          |                    | STANTEC             | STANTEC            | STANTEC       | STANTEC          | STANTEC          | STANTEC            | STANTEC             | STANTEC              | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            |
|                                                                                                                                                                   |          |                    | CCGE                | CCGE               | CCGE          | PARAROCH         | PARAROCH         | PARAROCH           | PARAROCH            | PARAROCH             | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           |
|                                                                                                                                                                   |          |                    | E2314               | E2314              | E2342         | 141138           | 141138           | 142196             | 142794              | 143439               | 144730             | 150382             | 150382             | 151696             | 151696             |
| E2314-01                                                                                                                                                          | E2314-02 | E2342-04           | 141138-11           | 141138-14          | 142196-07     | 142794-09        | 143439-10        | 144730-10          | 150382-05           | 150382-13            | 151696-11          | 151696-10          |                    |                    |                    |
| Units                                                                                                                                                             | TOGS     |                    |                     |                    |               |                  |                  |                    |                     |                      |                    |                    |                    |                    |                    |
|                                                                                                                                                                   |          |                    | Field Duplicate     |                    |               |                  | Field Duplicate  |                    |                     |                      |                    | Field Duplicate    |                    | Field Duplicate    |                    |
| General Chemistry                                                                                                                                                 |          |                    |                     |                    |               |                  |                  |                    |                     |                      |                    |                    |                    |                    |                    |
| Total Organic Carbon                                                                                                                                              | µg/L     | n/v                | -                   | -                  | -             | 6000             | 4600             | 15200              | 146000              | 24600                | 7300               | 6500               | 6000               | 5400               | 5300               |
| Metals                                                                                                                                                            |          |                    |                     |                    |               |                  |                  |                    |                     |                      |                    |                    |                    |                    |                    |
| Aluminum                                                                                                                                                          | µg/L     | n/v                | 36.9                | 32.5               | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Antimony                                                                                                                                                          | µg/L     | 3 <sup>B</sup>     | 12.5 U              | 12.5 U             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Arsenic                                                                                                                                                           | µg/L     | 25 <sup>B</sup>    | 5.000 U             | 5.000 U            | -             | 10 U             | 10 U             | 10 U               | 10 U                | 10 UJ                | 10 U               | 5.98 J             | 6.89 J             | 7.92 J             | 10.4               |
| Barium                                                                                                                                                            | µg/L     | 1000 <sup>B</sup>  | 62                  | 69.6               | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Beryllium                                                                                                                                                         | µg/L     | 3 <sup>A</sup>     | 1.500 U             | 1.500 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Cadmium                                                                                                                                                           | µg/L     | 5 <sup>B</sup>     | 1.500 U             | 1.500 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Calcium                                                                                                                                                           | µg/L     | n/v                | 121000              | 132000             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Chromium (Total)                                                                                                                                                  | µg/L     | 50 <sup>B</sup>    | 2.500 U             | 2.500 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Cobalt                                                                                                                                                            | µg/L     | n/v                | 7.500 U             | 7.500 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Copper                                                                                                                                                            | µg/L     | 200 <sup>B</sup>   | 5.000 U             | 5.000 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Iron                                                                                                                                                              | µg/L     | 300. <sup>B</sup>  | 25.0 U              | 25.0 U             | -             | 100 U            | 100 U            | 4330 <sup>B</sup>  | 9940 <sup>B</sup>   | 6480 <sup>B</sup>    | 10700 <sup>B</sup> | 13900 <sup>B</sup> | 13600 <sup>B</sup> | 10000 <sup>B</sup> | 10100 <sup>B</sup> |
| Lead                                                                                                                                                              | µg/L     | 25 <sup>B</sup>    | 12.6                | 12.5               | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Magnesium                                                                                                                                                         | µg/L     | 35000 <sup>A</sup> | 30600               | 33100              | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Manganese                                                                                                                                                         | µg/L     | 300. <sup>B</sup>  | 5.42 J              | 5.53 J             | -             | 694 <sup>B</sup> | 675 <sup>B</sup> | 1070 <sup>B</sup>  | 2280 <sup>B</sup>   | 1200 <sup>B</sup>    | 1060 <sup>B</sup>  | 844 <sup>B</sup>   | 838 <sup>B</sup>   | 945 <sup>B</sup>   | 949 <sup>B</sup>   |
| Mercury                                                                                                                                                           | µg/L     | 0.7 <sup>B</sup>   | 0.200 U             | 0.200 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Nickel                                                                                                                                                            | µg/L     | 100 <sup>B</sup>   | 2.52 J              | 10.0 U             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Potassium                                                                                                                                                         | µg/L     | n/v                | 9810                | 11100              | -             | n/v              | 11100            | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Selenium                                                                                                                                                          | µg/L     | 10 <sup>B</sup>    | 5.92                | 4.23 J             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Silver                                                                                                                                                            | µg/L     | 50 <sup>B</sup>    | 2.500 U             | 2.500 U            | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Sodium                                                                                                                                                            | µg/L     | 20000 <sup>B</sup> | 24700 <sup>B</sup>  | 27600 <sup>B</sup> | -             | 18500            | 18100            | 41100 <sup>B</sup> | 169000 <sup>B</sup> | 83100 M <sup>B</sup> | 63800 <sup>B</sup> | 58000 <sup>B</sup> | 58900 <sup>B</sup> | 49800 <sup>B</sup> | 50300 <sup>B</sup> |
| Thallium                                                                                                                                                          | µg/L     | 0.5 <sup>A</sup>   | 10.0 U              | 10.0 U             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Vanadium                                                                                                                                                          | µg/L     | n/v                | 10.0 U              | 10.0 U             | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Zinc                                                                                                                                                              | µg/L     | 2000 <sup>A</sup>  | 12.4                | 10.5               | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Volatile Organic Compounds                                                                                                                                        |          |                    |                     |                    |               |                  |                  |                    |                     |                      |                    |                    |                    |                    |                    |
| Acetone                                                                                                                                                           | µg/L     | 50 <sup>A</sup>    | 25 U                | 25 U               | 25 U          | 10.0 U           | 10.0 U           | 10.0 U             | 6.54 J              | 10.0 U               | 10.0 U             | 10.0 UJ            | 10.0 UJ            | 10.0 U             | 10.0 U             |
| Benzene                                                                                                                                                           | µg/L     | 1 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 1 U              | 1 U              | 1 U                | 1 U                 | 1 U                  | 1 U                | 1 U                | 1 U                | 0.700 U            | 0.700 U            |
| Bromodichloromethane                                                                                                                                              | µg/L     | 50 <sup>A</sup>    | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Bromoform (Tribromomethane)                                                                                                                                       | µg/L     | 50 <sup>A</sup>    | 5 U                 | 5 U                | 5 U           | 5.00 U           | 5.00 U           | 5.00 U             | 5.00 U              | 5.00 U               | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             |
| Bromomethane (Methyl bromide)                                                                                                                                     | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 UJ             | 2.00 UJ              | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Butylbenzene, n-                                                                                                                                                  | µg/L     | 5.. <sup>B</sup>   | -                   | -                  | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Butylbenzene, sec- (2-Phenylbutane)                                                                                                                               | µg/L     | 5.. <sup>B</sup>   | -                   | -                  | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Butylbenzene, tert-                                                                                                                                               | µg/L     | 5.. <sup>B</sup>   | -                   | -                  | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Carbon Disulfide                                                                                                                                                  | µg/L     | 60 <sup>A</sup>    | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Carbon Tetrachloride (Tetrachloromethane)                                                                                                                         | µg/L     | 5 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Chlorobenzene (Monochlorobenzene)                                                                                                                                 | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Chlorobromomethane                                                                                                                                                | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 5.00 U           | 5.00 U           | 5.00 U             | 5.00 U              | 5.00 U               | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             |
| Chloroethane (Ethyl Chloride)                                                                                                                                     | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Chloroethyl Vinyl Ether, 2-                                                                                                                                       | µg/L     | n/v                | -                   | -                  | -             | -                | -                | -                  | -                   | -                    | -                  | -                  | -                  | -                  | -                  |
| Chloroform (Trichloromethane)                                                                                                                                     | µg/L     | 7 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Chloromethane                                                                                                                                                     | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Cyclohexane                                                                                                                                                       | µg/L     | n/v                | 5 U                 | 5 U                | 5 U           | 10.0 U           | 10.0 U           | 10.0 U             | 10.0 U              | 10.0 U               | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             |
| Dibromo-3-Chloropropane, 1,2- (DBCP)                                                                                                                              | µg/L     | 0.04 <sup>B</sup>  | 5 U                 | 5 U                | 5 U           | 10.0 U           | 10.0 U           | 10.0 U             | 10.0 U              | 10.0 U               | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             |
| Dibromochloromethane                                                                                                                                              | µg/L     | 50 <sup>A</sup>    | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,2-                                                                                                                                             | µg/L     | 3 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,3-                                                                                                                                             | µg/L     | 3 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,4-                                                                                                                                             | µg/L     | 3 <sup>B</sup>     | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 U              | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorodifluoromethane (Freon 12)                                                                                                                                | µg/L     | 5.. <sup>B</sup>   | 5 U                 | 5 U                | 5 U           | 2.00 U           | 2.00 U           | 2.00 U             | 2.00 UJ             | 2.00 U               | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| On-Site Parking Lot                      |       |                               |               |                  |                   |                   |                   |                   |                   |                   |                   |                      |                   |                   |                   |
|------------------------------------------|-------|-------------------------------|---------------|------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|----------------------|-------------------|-------------------|-------------------|
| Area                                     |       |                               | B101MW        |                  |                   |                   |                   |                   |                   |                   |                   |                      | B102MW            |                   |                   |
| Sample Location                          |       |                               | 21-May-13     | 21-May-13        | 22-May-13         | 27-Mar-14         | 27-Mar-14         | 28-May-14         | 2-Jul-14          | 6-Aug-14          | 28-Oct-14         | 3-Feb-15             | 3-Feb-15          | 4-May-15          | 4-May-15          |
| Sample Date                              |       |                               | LI-B101MW-GW1 | LI-B101MW-GW1DUP | LI-B102MW-GW1     | LI-B102-MW        | LI-DUP-MW         | LI-B102-MW-PI1    | LI-B102-MW-PI2    | LI-B102-MW-PI3    | LI-B102-MW-PI6    | LI-B102-MW-PI9       | LI-DUP-PI9        | LI-B102-MW-PI12   | LI-DUP-PI12       |
| Sample ID                                |       |                               |               |                  |                   |                   |                   |                   |                   |                   |                   |                      |                   |                   |                   |
| Sampling Company                         |       |                               | STANTEC       | STANTEC          | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC              | STANTEC           | STANTEC           | STANTEC           |
| Laboratory                               |       |                               | CCGE          | CCGE             | CCGE              | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH          | PARAROCH          | PARAROCH          |
| Laboratory Work Order                    |       |                               | E2314         | E2314            | E2342             | 141138            | 141138            | 142196            | 142794            | 143439            | 144730            | 150382               | 150382            | 151696            | 151696            |
| Laboratory Sample ID                     |       |                               | E2314-01      | E2314-02         | E2342-04          | 141138-11         | 141138-14         | 142196-07         | 142794-09         | 143439-10         | 144730-10         | 150382-05            | 150382-13         | 151696-11         | 151696-10         |
| Sample Type                              | Units | TOGS                          |               | Field Duplicate  |                   |                   | Field Duplicate   |                   |                   |                   |                   |                      | Field Duplicate   |                   | Field Duplicate   |
| Volatile Organic Compounds (cont'd)      |       |                               |               |                  |                   |                   |                   |                   |                   |                   |                   |                      |                   |                   |                   |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 7.5 <sup>B</sup>  | 4.45              | 4.44              | 4.61              | 7.04 <sup>B</sup> | 68.7 <sup>B</sup> | 7.01 <sup>B</sup> | 2.00 U               | 2.00 U            | 4.10              | 4.11              |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 UJ           | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Dioxane, 1,4-                            | µg/L  | n/v                           | 100 U         | 100 U            | 100 U R           | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U            | 20.0 U               | 20.0 U            | 20.0 U            | 20.0 U            |
| Ethylbenzene                             | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | 25 U          | 25 U             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Isopropylbenzene                         | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sup>n,B</sup>              | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Methyl Acetate                           | µg/L  | n/v                           | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | 25 U          | 25 U             | 25 U              | 10.0 UJ           | 10.0 UJ           | 10.0 U            | 27.8 J            | 10.0 U            | 10.0 UJ           | 10.0 UJ              | 10.0 UJ           | 10.0 U            | 10.0 U            |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | 25 U          | 25 U             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Methylcyclohexane                        | µg/L  | n/v                           | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Propylbenzene, n-                        | µg/L  | 5 <sup>n,B</sup>              | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Styrene                                  | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sup>n,B</sup>              | 1.6 J         | 1.2 J            | 20.9 <sup>B</sup> | 24.4 <sup>B</sup> | 25.4 <sup>B</sup> | 20.6 <sup>B</sup> | 26.4 <sup>B</sup> | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Toluene                                  | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U            | 5.00 U            |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Trichloroethene (TCE)                    | µg/L  | 5 <sup>n,B</sup>              | 0.51 J        | 5 U              | 14.9 <sup>B</sup> | 9.78 <sup>B</sup> | 10.2 <sup>B</sup> | 7.72 <sup>B</sup> | 15.3 <sup>B</sup> | 2.09              | 2.00 U            | 2.00 U               | 2.00 U            | 2.38              | 2.42              |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sup>n,B</sup>              | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sup>n,B</sup>              | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Vinyl Acetate                            | µg/L  | n/v                           | -             | -                | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 5 U           | 5 U              | 0.53 J            | 2.00 U            | 2.00 U            | 2.00 U            | 1.45 J            | 4.49 <sup>B</sup> | 20.8 <sup>B</sup> | 11.7 NJ <sup>B</sup> | 11.9 <sup>B</sup> | 11.0 <sup>B</sup> | 11.3 <sup>B</sup> |
| Xylene, m & p-                           | µg/L  | 5 <sup>n,B</sup>              | 10 U          | 10 U             | 10 U              | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Xylene, o-                               | µg/L  | 5 <sup>n,B</sup>              | 5 U           | 5 U              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U            | 2.00 U            |
| Total VOC                                | µg/L  | n/v                           | 2.11          | 1.2              | 43.83             | 38.63             | 40.04             | 32.93             | 84.53             | 75.28             | 27.81             | 11.7                 | 11.9              | 17.48             | 17.83             |
| VOC Tentatively Identified Compounds     |       |                               |               |                  |                   |                   |                   |                   |                   |                   |                   |                      |                   |                   |                   |
| Total VOC TICs                           | µg/L  | n/v                           | 2.5 U         | 2.5 U            | 2.5 U             | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -                 | -                 |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                      |       |                    | On-Site Parking Lot |                    |           |                    |                     |                     |                    |                     |                    |            |              |           |
|-------------------------------------------|-------|--------------------|---------------------|--------------------|-----------|--------------------|---------------------|---------------------|--------------------|---------------------|--------------------|------------|--------------|-----------|
| Sample Location                           |       |                    | RW-4                |                    |           |                    |                     |                     |                    |                     |                    | RW-11      |              |           |
| Sample Date                               |       |                    | 25-Apr-12           | 22-May-13          | 26-Mar-14 | 29-May-14          | 2-Jul-14            | 6-Aug-14            | 29-Oct-14          | 4-Feb-15            | 4-May-15           | 14-Jun-12  | 22-May-13    | 27-Mar-14 |
| Sample ID                                 |       |                    | RW-4                | LI-RW-4-GW1        | LI-RW-4   | LI-RW-4-PI1        | LI-RW-4-PI2         | LI-RW-4-PI3         | LI-RW4-PI6         | LI-RW-4-PI9         | LI-RW-4-PI12       | RW-11      | LI-RW-11-GW1 | LI-RW-11  |
| Sampling Company                          |       |                    | DECI                | STANTEC            | STANTEC   | STANTEC            | STANTEC             | STANTEC             | STANTEC            | STANTEC             | STANTEC            | DECI       | STANTEC      | STANTEC   |
| Laboratory                                |       |                    | PARAROCH            | CCGE               | PARAROCH  | PARAROCH           | PARAROCH            | PARAROCH            | PARAROCH           | PARAROCH            | PARAROCH           | PARAROCH   | CCGE         | PARAROCH  |
| Laboratory Work Order                     |       |                    | 12:1770             | E2342              | 141138    | 142196             | 142794              | 143439              | 144730             | 150382              | 151696             | 12:2523    | E2342        | 141138    |
| Laboratory Sample ID                      |       |                    | 12:1770-01          | E2342-03           | 141138-04 | 142196-13          | 142794-10           | 143439-04           | 144730-04          | 150382-11           | 151696-04          | 12:2523-03 | E2342-02     | 141138-09 |
| Sample Type                               | Units | TOGS               |                     |                    |           |                    |                     |                     |                    |                     |                    |            |              |           |
| General Chemistry                         |       |                    |                     |                    |           |                    |                     |                     |                    |                     |                    |            |              |           |
| Total Organic Carbon                      | µg/L  | n/v                | -                   | -                  | -         | 8200               | 339000              | 63000               | 6900               | 5900                | 5400               | -          | -            | -         |
| Metals                                    |       |                    |                     |                    |           |                    |                     |                     |                    |                     |                    |            |              |           |
| Aluminum                                  | µg/L  | n/v                | -                   | 43.8               | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Antimony                                  | µg/L  | 3 <sup>B</sup>     | -                   | 12.5 U             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Arsenic                                   | µg/L  | 25 <sup>B</sup>    | -                   | 5,000 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Barium                                    | µg/L  | 1000 <sup>B</sup>  | -                   | 151                | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Beryllium                                 | µg/L  | 3 <sup>A</sup>     | -                   | 1,500 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Cadmium                                   | µg/L  | 5 <sup>B</sup>     | -                   | 1,500 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Calcium                                   | µg/L  | n/v                | -                   | 141000             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Chromium (Total)                          | µg/L  | 50 <sup>B</sup>    | -                   | 2,500 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Cobalt                                    | µg/L  | n/v                | -                   | 7,500 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Copper                                    | µg/L  | 200 <sup>B</sup>   | -                   | 5,000 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Iron                                      | µg/L  | 300. <sup>B</sup>  | -                   | 11.7 J             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Lead                                      | µg/L  | 25 <sup>B</sup>    | -                   | 17                 | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Magnesium                                 | µg/L  | 35000 <sup>A</sup> | -                   | 29800              | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Manganese                                 | µg/L  | 300. <sup>B</sup>  | -                   | 667 J <sup>B</sup> | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Mercury                                   | µg/L  | 0.7 <sup>B</sup>   | -                   | 0.200 U            | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Nickel                                    | µg/L  | 100 <sup>B</sup>   | -                   | 6.32 J             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Potassium                                 | µg/L  | n/v                | -                   | 17800              | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Selenium                                  | µg/L  | 10 <sup>B</sup>    | -                   | 5.52               | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Silver                                    | µg/L  | 50 <sup>B</sup>    | -                   | 2,500 U N          | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Sodium                                    | µg/L  | 20000 <sup>B</sup> | -                   | 8750               | -         | 22300 <sup>B</sup> | 298000 <sup>B</sup> | 222000 <sup>B</sup> | 43500 <sup>B</sup> | 110000 <sup>B</sup> | 86900 <sup>B</sup> | -          | -            | -         |
| Thallium                                  | µg/L  | 0.5 <sup>A</sup>   | -                   | 10.0 U             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Vanadium                                  | µg/L  | n/v                | -                   | 10.0 U             | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Zinc                                      | µg/L  | 2000 <sup>A</sup>  | -                   | 18.2               | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Volatile Organic Compounds                |       |                    |                     |                    |           |                    |                     |                     |                    |                     |                    |            |              |           |
| Acetone                                   | µg/L  | 50 <sup>A</sup>    | 10.0 UJ             | 25 U               | 10.0 U    | 6.72 J             | 10.0 U              | 12.7 J              | 10.0 U             | 10.0 UJ             | 10.0 U             | -          | 25 U         | 10.0 U    |
| Benzene                                   | µg/L  | 1 <sup>B</sup>     | 0.700 UJ            | 5 U                | 1 U       | 1 U                | 1 U                 | 1 U                 | 1 U                | 1 U                 | 0.700 U            | -          | 5 U          | 1 U       |
| Bromodichloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Bromoform (Tribromomethane)               | µg/L  | 50 <sup>A</sup>    | 5.00 UJ             | 5 U                | 5.00 U    | 5.00 U             | 5.00 U              | 5.00 UJ             | 5.00 U             | 5.00 U              | 5.00 U             | 5.00 U     | 5 U          | 5.00 U    |
| Bromomethane (Methyl bromide)             | µg/L  | 5.. <sup>B</sup>   | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 UJ    | 5 U          | 2.00 U    |
| Butylbenzene, n-                          | µg/L  | 5.. <sup>B</sup>   | -                   | -                  | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Butylbenzene, sec- (2-Phenylbutane)       | µg/L  | 5.. <sup>B</sup>   | -                   | -                  | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Butylbenzene, tert-                       | µg/L  | 5.. <sup>B</sup>   | -                   | -                  | -         | -                  | -                   | -                   | -                  | -                   | -                  | -          | -            | -         |
| Carbon Disulfide                          | µg/L  | 60 <sup>A</sup>    | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 3.04                | 3.64                | 2.00 U             | 2.00 U              | 2.00 U             | -          | 5 U          | 2.00 U    |
| Carbon Tetrachloride (Tetrachloromethane) | µg/L  | 5 <sup>B</sup>     | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 UJ    | 5 U          | 2.00 U    |
| Chlorobenzene (Monochlorobenzene)         | µg/L  | 5.. <sup>B</sup>   | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Chlorobromomethane                        | µg/L  | 5.. <sup>B</sup>   | -                   | 5 U                | 5.00 U    | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 U              | 5.00 U             | -          | 5 U Q        | 5.00 U    |
| Chloroethane (Ethyl Chloride)             | µg/L  | 5.. <sup>B</sup>   | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Chloroethyl Vinyl Ether, 2-               | µg/L  | n/v                | 10.0 U R            | -                  | -         | -                  | -                   | -                   | -                  | -                   | -                  | 10.0 U R   | -            | -         |
| Chloroform (Trichloromethane)             | µg/L  | 7 <sup>B</sup>     | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 1.91 J              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Chloromethane                             | µg/L  | 5.. <sup>B</sup>   | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Cyclohexane                               | µg/L  | n/v                | -                   | 5 U                | 10.0 U    | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U              | 10.0 U             | -          | 5 U          | 10.0 U    |
| Dibromo-3-Chloropropane, 1,2- (DBCP)      | µg/L  | 0.04 <sup>B</sup>  | -                   | 5 U                | 10.0 U    | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U              | 10.0 U             | -          | 5 U          | 10.0 U    |
| Dibromochloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Dichlorobenzene, 1,2-                     | µg/L  | 3 <sup>B</sup>     | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Dichlorobenzene, 1,3-                     | µg/L  | 3 <sup>B</sup>     | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Dichlorobenzene, 1,4-                     | µg/L  | 3 <sup>B</sup>     | 2.00 UJ             | 5 U                | 2.00 U    | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | 2.00 U     | 5 U          | 2.00 U    |
| Dichlorodifluoromethane (Freon 12)        | µg/L  | 5.. <sup>B</sup>   | -                   | 5 U                | 2.00 U    | 2.00 U             | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             | -          | 5 U          | 2.00 U    |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                     |       |                               | On-Site Parking Lot |                   |                   |                   |                   |                   |                   |                      |                   |            |              |           |
|------------------------------------------|-------|-------------------------------|---------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|----------------------|-------------------|------------|--------------|-----------|
| Sample Location                          |       |                               |                     |                   |                   |                   |                   |                   |                   |                      |                   |            |              |           |
| Sample Date                              |       |                               | 25-Apr-12           | 22-May-13         | 26-Mar-14         | 29-May-14         | 2-Jul-14          | 6-Aug-14          | 29-Oct-14         | 4-Feb-15             | 4-May-15          | 14-Jun-12  | 22-May-13    | 27-Mar-14 |
| Sample ID                                |       |                               | RW-4                | LI-RW-4-GW1       | LI-RW-4           | LI-RW-4-PI1       | LI-RW-4-PI2       | LI-RW-4-PI3       | LI-RW4-PI6        | LI-RW-4-PI9          | LI-RW-4-PI12      | RW-11      | LI-RW-11-GW1 | LI-RW-11  |
| Sampling Company                         |       |                               | DECI                | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC              | STANTEC           | DECI       | STANTEC      | STANTEC   |
| Laboratory                               |       |                               | PARAROCH            | CCGE              | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH          | PARAROCH   | CCGE         | PARAROCH  |
| Laboratory Work Order                    |       |                               | 12:1770             | E2342             | 141138            | 142196            | 142794            | 143439            | 144730            | 150382               | 151696            | 12:2523    | E2342        | 141138    |
| Laboratory Sample ID                     |       |                               | 12:1770-01          | E2342-03          | 141138-04         | 142196-13         | 142794-10         | 143439-04         | 144730-04         | 150382-11            | 151696-04         | 12:2523-03 | E2342-02     | 141138-09 |
| Sample Type                              | Units | TOGS                          |                     |                   |                   |                   |                   |                   |                   |                      |                   |            |              |           |
| Volatile Organic Compounds (cont'd)      |       |                               |                     |                   |                   |                   |                   |                   |                   |                      |                   |            |              |           |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ    | 5 U          | 2.00 U    |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sup>B</sup>                | 23.1 J <sup>B</sup> | 14.9 <sup>B</sup> | 6.41 <sup>B</sup> | 9.56 <sup>B</sup> | 13.4 <sup>B</sup> | 87.9 <sup>B</sup> | 47.3 <sup>B</sup> | 23.7 <sup>B</sup>    | 14.8 <sup>B</sup> | 2.00 U     | 5 U          | 2.00 U    |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 1.11 J            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sup>p</sup> <sup>B</sup> | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sup>p</sup> <sup>B</sup> | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Dioxane, 1,4-                            | µg/L  | n/v                           | -                   | 100 U R           | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U            | 20.0 U               | 20.0 U            | -          | 100 U R      | 20.0 U R  |
| Ethylbenzene                             | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | 5.00 UJ             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | -          | 25 U         | 5.00 U    |
| Isopropylbenzene                         | µg/L  | 5 <sup>B</sup>                | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sup>B</sup>                | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Methyl Acetate                           | µg/L  | n/v                           | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5.93              | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | 10.0 UJ             | 25 U              | 10.0 UJ           | 10.0 U            | 20.8 J            | 15.8              | 10.0 UJ           | 10.0 UJ              | 10.0 U            | -          | 25 U         | 10.0 UJ   |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | 5.00 UJ             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | -          | 25 U         | 5.00 U    |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Methylcyclohexane                        | µg/L  | n/v                           | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sup>B</sup>                | 5.00 UJ             | 5 U               | 4.35 JB           | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 U     | 5 U          | 5.00 U    |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Propylbenzene, n-                        | µg/L  | 5 <sup>B</sup>                | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Styrene                                  | µg/L  | 5 <sup>B</sup>                | 5.00 UJ             | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | -          | 5 U          | 5.00 U    |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sup>B</sup>                | 62.6 J <sup>B</sup> | 55.8 <sup>B</sup> | 62.7 <sup>B</sup> | 76.0 <sup>B</sup> | 73.0 <sup>B</sup> | 54.5 <sup>B</sup> | 10.3 <sup>B</sup> | 9.17 <sup>B</sup>    | 18.7 <sup>B</sup> | 2.00 U     | 1.3 J        | 1.11 J    |
| Toluene                                  | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sup>B</sup>                | -                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | -          | 5 U          | 5.00 U    |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sup>B</sup>                | -                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            | -          | 5 U          | 5.00 U    |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 U     | 5 U          | 2.00 U    |
| Trichloroethene (TCE)                    | µg/L  | 5 <sup>B</sup>                | 21.4 J <sup>B</sup> | 19.8 <sup>B</sup> | 10.3 <sup>B</sup> | 18.0 <sup>B</sup> | 20.4 <sup>B</sup> | 34.3 <sup>B</sup> | 13.7 <sup>B</sup> | 5.85 <sup>B</sup>    | 8.94 <sup>B</sup> | 2.00 U     | 5 U          | 2.00 U    |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ    | 5 U          | 2.00 U    |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sup>B</sup>                | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sup>B</sup>                | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sup>B</sup>                | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Vinyl Acetate                            | µg/L  | n/v                           | 5.00 UJ             | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | -            | -         |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 3.86 J <sup>B</sup> | 1.8 J             | 1.72 J            | 2.00 U            | 3.07 <sup>B</sup> | 2.00 U            | 28.4 <sup>B</sup> | 4.58 NJ <sup>B</sup> | 5.15 <sup>B</sup> | 2.00 U     | 5 U          | 2.00 U    |
| Xylene, m & p-                           | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 10 U              | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 10 U         | 2.00 U    |
| Xylene, o-                               | µg/L  | 5 <sup>B</sup>                | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            | -          | 5 U          | 2.00 U    |
| Total VOC                                | µg/L  | n/v                           | 110.96              | 92.3              | 85.48             | 110.28            | 135.62            | 209.95            | 105.63            | 43.3                 | 47.59             | ND         | 1.3          | 1.11      |
| VOC Tentatively Identified Compounds     |       |                               |                     |                   |                   |                   |                   |                   |                   |                      |                   |            |              |           |
| Total VOC TICs                           | µg/L  | n/v                           | -                   | 2.5 U             | -                 | -                 | -                 | -                 | -                 | -                    | -                 | -          | 2.5 U        | -         |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                      |       |                    | B106MW        |            |                     |                     |                     |                    |                    |                    | On-Site Building   |                    |                     |                     | B108MW                |                     |                    |                     |                    |  |  |
|-------------------------------------------|-------|--------------------|---------------|------------|---------------------|---------------------|---------------------|--------------------|--------------------|--------------------|--------------------|--------------------|---------------------|---------------------|-----------------------|---------------------|--------------------|---------------------|--------------------|--|--|
| Sample Location                           |       |                    | 23-May-13     | 26-Mar-14  | 28-May-14           | 2-Jul-14            | 7-Aug-14            | 28-Oct-14          | 3-Feb-15           | 5-May-15           | 23-May-13          | 26-Mar-14          | 28-May-14           | 28-May-14           | 2-Jul-14              | 8-Aug-14            | 29-Oct-14          | 3-Feb-15            | 5-May-15           |  |  |
| Sample Date                               |       |                    | LI-B106MW-GW1 | LI-B106-MW | LI-B106-MW-PI1      | LI-B106-MW-PI2      | LI-B106-MW-PI3      | LI-B106-MW-PI6     | LI-B106-MW-PI9     | LI-B106-MW-PI12    | LI-B108MW-GW1      | LI-B108-MW         | LI-B108-MW-PI1      | LI-MW-DUP-PI1       | LI-B108-MW-PI2        | LI-B108-MW-PI3      | LI-B108-MW-PI6     | LI-B108-MW-PI9      | LI-B108-MW-PI12    |  |  |
| Sample ID                                 |       |                    | LI-B106MW-GW1 | LI-B106-MW | LI-B106-MW-PI1      | LI-B106-MW-PI2      | LI-B106-MW-PI3      | LI-B106-MW-PI6     | LI-B106-MW-PI9     | LI-B106-MW-PI12    | LI-B108MW-GW1      | LI-B108-MW         | LI-B108-MW-PI1      | LI-MW-DUP-PI1       | LI-B108-MW-PI2        | LI-B108-MW-PI3      | LI-B108-MW-PI6     | LI-B108-MW-PI9      | LI-B108-MW-PI12    |  |  |
| Sampling Company                          |       |                    | STANTEC       | STANTEC    | STANTEC             | STANTEC             | STANTEC             | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC             | STANTEC             | STANTEC               | STANTEC             | STANTEC            | STANTEC             | STANTEC            |  |  |
| Laboratory                                |       |                    | CCGE          | PARAROCH   | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH           | PARAROCH           | PARAROCH           | CCGE               | PARAROCH           | PARAROCH            | PARAROCH            | PARAROCH              | PARAROCH            | PARAROCH           | PARAROCH            | PARAROCH           |  |  |
| Laboratory Work Order                     |       |                    | E2363         | 141138     | 142196              | 142794              | 143439              | 144730             | 150382             | 151696             | E2363              | 141138             | 142196              | 142196              | 142794                | 143439              | 144730             | 150382              | 151696             |  |  |
| Laboratory Sample ID                      |       |                    | E2363-03      | 141138-12  | 142196-06           | 142794-11           | 143439-11           | 144730-11          | 150382-06          | 151696-12          | E2363-02           | 141138-13          | 142196-04           | 142196-05           | 142794-12             | 143439-12           | 144730-12          | 150382-07           | 151696-13          |  |  |
| Sample Type                               | Units | TOGS               |               |            |                     |                     |                     |                    |                    |                    |                    |                    |                     | Field Duplicate     |                       |                     |                    |                     |                    |  |  |
| General Chemistry                         |       |                    |               |            |                     |                     |                     |                    |                    |                    |                    |                    |                     |                     |                       |                     |                    |                     |                    |  |  |
| Total Organic Carbon                      | µg/L  | n/v                | -             | -          | 188000              | 514000              | 77600               | 4000 J-            | 3100 J+            | 1500               | -                  | 3300               | 60300               | 60200               | 86100                 | 72200               | 45000              | 18100 J             | 1700               |  |  |
| Metals                                    |       |                    |               |            |                     |                     |                     |                    |                    |                    |                    |                    |                     |                     |                       |                     |                    |                     |                    |  |  |
| Aluminum                                  | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 66                 | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Antimony                                  | µg/L  | 3 <sup>B</sup>     | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 12.5 U             | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Arsenic                                   | µg/L  | 25 <sup>B</sup>    | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 6.2                | 10 U               | 10 U                | 10 U                | 10 U                  | 10 U                | 10 U               | 5.92 J              | 10.0 U             |  |  |
| Barium                                    | µg/L  | 1000 <sup>B</sup>  | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 54.9               | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Beryllium                                 | µg/L  | 3 <sup>A</sup>     | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 1.500 U            | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Cadmium                                   | µg/L  | 5 <sup>B</sup>     | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 0.7 J              | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Calcium                                   | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 97000              | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Chromium (Total)                          | µg/L  | 50 <sup>B</sup>    | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 2.500 U            | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Cobalt                                    | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 7.500 U            | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Copper                                    | µg/L  | 200 <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 4.16 J             | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Iron                                      | µg/L  | 300. <sup>B</sup>  | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 45.3               | 100 U              | 1400 <sup>B</sup>   | 978 <sup>B</sup>    | 3520 <sup>B</sup>     | 2480 <sup>B</sup>   | 2350 <sup>B</sup>  | 2660 <sup>B</sup>   | 999 <sup>B</sup>   |  |  |
| Lead                                      | µg/L  | 25 <sup>B</sup>    | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 4.9                | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Magnesium                                 | µg/L  | 35000 <sup>A</sup> | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 23200              | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Manganese                                 | µg/L  | 300. <sup>B</sup>  | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 46.4 J             | 187                | 184                 | 179                 | 217                   | 158                 | 106                | 87.6                | 81.8               |  |  |
| Mercury                                   | µg/L  | 0.7 <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 0.200 U            | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Nickel                                    | µg/L  | 100 <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 2.1 J              | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Potassium                                 | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 10500              | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Selenium                                  | µg/L  | 10 <sup>B</sup>    | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 5.03               | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Silver                                    | µg/L  | 50 <sup>B</sup>    | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 2.500 U            | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Sodium                                    | µg/L  | 20000 <sup>B</sup> | -             | -          | 162000 <sup>B</sup> | 375000 <sup>B</sup> | 185000 <sup>B</sup> | 59200 <sup>B</sup> | 50200 <sup>B</sup> | 40100 <sup>B</sup> | 26300 <sup>B</sup> | 33000 <sup>B</sup> | 103000 <sup>B</sup> | 101000 <sup>B</sup> | 100000 M <sup>B</sup> | 115000 <sup>B</sup> | 82900 <sup>B</sup> | 130000 <sup>B</sup> | 42400 <sup>B</sup> |  |  |
| Thallium                                  | µg/L  | 0.5 <sup>A</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 10.0 U             | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Vanadium                                  | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 10.0 U             | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Zinc                                      | µg/L  | 2000 <sup>A</sup>  | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | 8.94 J             | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Volatile Organic Compounds                |       |                    |               |            |                     |                     |                     |                    |                    |                    |                    |                    |                     |                     |                       |                     |                    |                     |                    |  |  |
| Acetone                                   | µg/L  | 50 <sup>A</sup>    | 25 U          | 10.0 U     | 10.0 U              | 12.9                | 10.0 U              | 10.0 U             | 10.0 UJ            | 10.0 U             | 25 U               | 10.0 U             | 10.0 U              | 10.0 U              | 6.04 J                | 8.49 J              | 10.0 U             | 6.51 J              | 10.0 U             |  |  |
| Benzene                                   | µg/L  | 1 <sup>B</sup>     | 5 U           | 1 U        | 1 U                 | 0.842 J             | 0.391 J             | 1 U                | 1 U                | 0.700 U            | 5 U                | 1 U                | 1 U                 | 1 U                 | 1 U                   | 1 U                 | 1 U                | 1 U                 | 0.700 U            |  |  |
| Bromodichloromethane                      | µg/L  | 50 <sup>A</sup>    | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Bromoform (Tribromomethane)               | µg/L  | 50 <sup>A</sup>    | 5 U           | 5.00 U     | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 U             | 5.00 U             | 5 U                | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 U                | 5.00 U              | 5.00 U             | 5.00 U              | 5.00 U             |  |  |
| Bromomethane (Methyl bromide)             | µg/L  | 5.. <sup>B</sup>   | 5 U           | 2.00 U     | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 UJ               | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Butylbenzene, n-                          | µg/L  | 5.. <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | -                  | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Butylbenzene, sec- (2-Phenylbutane)       | µg/L  | 5.. <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | -                  | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Butylbenzene, tert-                       | µg/L  | 5.. <sup>B</sup>   | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | -                  | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Carbon Disulfide                          | µg/L  | 60 <sup>A</sup>    | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Carbon Tetrachloride (Tetrachloromethane) | µg/L  | 5 <sup>B</sup>     | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Chlorobenzene (Monochlorobenzene)         | µg/L  | 5.. <sup>B</sup>   | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Chlorobromomethane                        | µg/L  | 5.. <sup>B</sup>   | 5 U           | 5.00 U     | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 U             | 5.00 U             | 5 U                | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 U                | 5.00 U              | 5.00 U             | 5.00 U              | 5.00 U             |  |  |
| Chloroethane (Ethyl Chloride)             | µg/L  | 5.. <sup>B</sup>   | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Chloroethyl Vinyl Ether, 2-               | µg/L  | n/v                | -             | -          | -                   | -                   | -                   | -                  | -                  | -                  | -                  | -                  | -                   | -                   | -                     | -                   | -                  | -                   | -                  |  |  |
| Chloroform (Trichloromethane)             | µg/L  | 7 <sup>B</sup>     | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Chloromethane                             | µg/L  | 5.. <sup>B</sup>   | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Cyclohexane                               | µg/L  | n/v                | 0.69 J        | 10.0 U     | 10.0 U              | 15.8                | 7.47 J              | 10.0 U             | 10.0 U             | 10.0 U             | 5 U                | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U                | 10.0 U              | 10.0 U             | 10.0 U              | 10.0 U             |  |  |
| Dibromo-3-Chloropropane, 1,2- (DBCP)      | µg/L  | 0.04 <sup>B</sup>  | 5 U           | 10.0 U     | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U             | 10.0 U             | 5 U                | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U                | 10.0 U              | 10.0 U             | 10.0 U              | 10.0 U             |  |  |
| Dibromochloromethane                      | µg/L  | 50 <sup>A</sup>    | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 10.4               | 2.00 U              | 2.00 U             |  |  |
| Dichlorobenzene, 1,2-                     | µg/L  | 3 <sup>B</sup>     | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Dichlorobenzene, 1,3-                     | µg/L  | 3 <sup>B</sup>     | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Dichlorobenzene, 1,4-                     | µg/L  | 3 <sup>B</sup>     | 5 U           | 2.00 U     | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U                | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |
| Dichlorodifluoromethane (Freon 12)        | µg/L  | 5.. <sup>B</sup>   | 5 U           | 2.00 U     | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 UJ               | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U             |  |  |

See last page for notes.



Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                     |       |                               | B106MW             |                   |                   |                    |                   |                   |                   |                   |                   |                   | On-Site Building  |                   |                   |                   | B108MW            |                   |                   |  |  |  |
|------------------------------------------|-------|-------------------------------|--------------------|-------------------|-------------------|--------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|--|--|--|
| Sample Location                          |       |                               |                    |                   |                   |                    |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |  |  |  |
| Sample Date                              |       |                               | 23-May-13          | 26-Mar-14         | 28-May-14         | 2-Jul-14           | 7-Aug-14          | 28-Oct-14         | 3-Feb-15          | 5-May-15          | 23-May-13         | 26-Mar-14         | 28-May-14         | 28-May-14         | 2-Jul-14          | 8-Aug-14          | 29-Oct-14         | 3-Feb-15          | 5-May-15          |  |  |  |
| Sample ID                                |       |                               | LI-B106MW-GW1      | LI-B106-MW        | LI-B106-MW-PI1    | LI-B106-MW-PI2     | LI-B106-MW-PI3    | LI-B106-MW-PI6    | LI-B106-MW-PI9    | LI-B106-MW-PI12   | LI-B108MW-GW1     | LI-B108-MW        | LI-B108-MW-PI1    | LI-MW-DUP-PI1     | LI-B108-MW-PI2    | LI-B108-MW-PI3    | LI-B108-MW-PI6    | LI-B108-MW-PI9    | LI-B108-MW-PI12   |  |  |  |
| Sampling Company                         |       |                               | STANTEC            | STANTEC           | STANTEC           | STANTEC            | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           |  |  |  |
| Laboratory                               |       |                               | CCGE               | PARAROCH          | PARAROCH          | PARAROCH           | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | CCGE              | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          |  |  |  |
| Laboratory Work Order                    |       |                               | E2363              | 141138            | 142196            | 142794             | 143439            | 144730            | 150382            | 151696            | E2363             | 141138            | 142196            | 142196            | 142794            | 143439            | 144730            | 150382            | 151696            |  |  |  |
| Laboratory Sample ID                     |       |                               | E2363-03           | 141138-12         | 142196-06         | 142794-11          | 143439-11         | 144730-11         | 150382-06         | 151696-12         | E2363-02          | 141138-13         | 142196-04         | 142196-05         | 142794-12         | 143439-12         | 144730-12         | 150382-07         | 151696-13         |  |  |  |
| Sample Type                              | Units | TOGS                          |                    |                   |                   |                    |                   |                   |                   |                   |                   |                   |                   | Field Duplicate   |                   |                   |                   |                   |                   |  |  |  |
| Volatile Organic Compounds (cont'd)      |       |                               |                    |                   |                   |                    |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |  |  |  |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 16.9 <sup>B</sup>  | 6.89 <sup>B</sup> | 8.67 <sup>B</sup> | 28.4 <sup>B</sup>  | 16.3 <sup>B</sup> | 40.4 <sup>B</sup> | 26.0 <sup>B</sup> | 10.2 <sup>B</sup> | 5.7 <sup>B</sup>  | 2.00 U            | 11.0 <sup>B</sup> | 10.9 <sup>B</sup> | 23.2 <sup>B</sup> | 4.99              | 1.96 J            | 1.87 J            | 1.49 J            |  |  |  |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 1.4 J              | 2.00 U            | 2.00 U            | 3.84               | 1.61 J            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.04              | 1.37 J            | 2.00 U            | 1.42 J            | 2.00 U            |  |  |  |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Dioxane, 1,4-                            | µg/L  | n/v                           | 100 U R            | 20.0 U R          | 20.0 U R          | 20.0 U R           | 20.0 U R          | 20.0 U            | 20.0 U            | 20.0 U            | 100 U R           | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U            | 20.0 U            | 20.0 U            |  |  |  |
| Ethylbenzene                             | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 1.79 J             | 1.20 J            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | 25 U               | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Isopropylbenzene                         | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Methyl Acetate                           | µg/L  | n/v                           | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | 25 U               | 10.0 UJ           | 10.7              | 151 J <sup>A</sup> | 31.6              | 10.0 UJ           | 10.0 UJ           | 10.0 U            | 25 U              | 10.0 UJ           | 13.9              | 12.8              | 41.5 J            | 63.8 <sup>A</sup> | 22.9 J            | 23.0 J            | 10.0 U            |  |  |  |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | 25 U               | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Methylcyclohexane                        | µg/L  | n/v                           | 0.77 J             | 2.00 U            | 2.03              | 11.7               | 6.30              | 1.21 J            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5 U               | 4.93 JB           | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Propylbenzene, n-                        | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Styrene                                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 14.8 <sup>B</sup>  | 21.7 <sup>B</sup> | 9.51 <sup>B</sup> | 11.7 <sup>B</sup>  | 7.73 <sup>B</sup> | 2.00 U            | 2.00 U            | 2.00 U            | 15.9 <sup>B</sup> | 6.45 <sup>B</sup> | 10.1 <sup>B</sup> | 9.75 <sup>B</sup> | 10.7 <sup>B</sup> | 9.63 <sup>B</sup> | 10.4 <sup>B</sup> | 6.73 <sup>B</sup> | 14.4 <sup>B</sup> |  |  |  |
| Toluene                                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 5.00 U            | 5.00 U            | 5.00 U             | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            |  |  |  |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Trichloroethene (TCE)                    | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 12 <sup>B</sup>    | 8.27 <sup>B</sup> | 5.11 <sup>B</sup> | 9.44 <sup>B</sup>  | 16.6 <sup>B</sup> | 2.23              | 2.00 U            | 2.12              | 8.5 <sup>B</sup>  | 1.05 J            | 4.17              | 4.15              | 4.21              | 1.65 J            | 4.04              | 2.93              | 2.72              |  |  |  |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Vinyl Acetate                            | µg/L  | n/v                           | -                  | -                 | -                 | -                  | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 2.1 J <sup>B</sup> | 2.00 U            | 2.84 <sup>B</sup> | 15.2 <sup>B</sup>  | 7.60 <sup>B</sup> | 15.2 <sup>B</sup> | 12.8 <sup>B</sup> | 4.89 <sup>B</sup> | 5 U               | 2.00 U            | 2.75 <sup>B</sup> | 2.61 <sup>B</sup> | 10.2 <sup>B</sup> | 14.6 <sup>B</sup> | 4.23 <sup>B</sup> | 2.00 U            | 2.00 U            |  |  |  |
| Xylene, m & p-                           | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 10 U               | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 10 U              | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Xylene, o-                               | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5 U                | 2.00 U            | 2.00 U            | 2.00 U             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            |  |  |  |
| Total VOC                                | µg/L  | n/v                           | 48.66              | 36.86             | 38.86             | 262.612            | 96.801            | 59.04             | 38.8              | 17.21             | 30.1              | 12.43             | 41.92             | 40.21             | 97.89             | 104.53            | 53.93             | 42.46             | 18.61             |  |  |  |
| VOC Tentatively Identified Compounds     |       |                               |                    |                   |                   |                    |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |                   |  |  |  |
| Total VOC TICs                           | µg/L  | n/v                           | 2.5 U              | -                 | -                 | -                  | -                 | -                 | -                 | -                 | 2.5 U             | -                 | -                 | -                 | -                 | -                 | -                 | -                 | -                 |  |  |  |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area<br>Sample Location<br>Sample Date    |       |                    | On-Site Building |             |           |                     |                     |                     |                     |                    |                     |            |                    |                    |                     |                     |                     |                     |                     |                    |
|-------------------------------------------|-------|--------------------|------------------|-------------|-----------|---------------------|---------------------|---------------------|---------------------|--------------------|---------------------|------------|--------------------|--------------------|---------------------|---------------------|---------------------|---------------------|---------------------|--------------------|
|                                           |       |                    | RW-1             |             |           |                     |                     |                     |                     |                    |                     |            | RW-2               |                    |                     |                     |                     |                     |                     |                    |
| Sample ID                                 |       |                    | 23-Mar-12        | 23-May-13   | 26-Mar-14 | 29-May-14           | 1-Jul-14            | 8-Aug-14            | 29-Oct-14           | 3-Feb-15           | 5-May-15            | 23-Mar-12  | 21-May-13          | 26-Mar-14          | 29-May-14           | 1-Jul-14            | 8-Aug-14            | 29-Oct-14           | 3-Feb-15            | 5-May-15           |
| Sample ID                                 |       |                    | RW-1             | LI-RW-1-GW1 | LI-RW-1   | LI-RW-1-PI1         | LI-RW-1-PI2         | LI-RW-1-PI3         | LI-RW1-PI6          | LI-RW-1-PI9        | LI-RW-1-PI12        | RW-2       | LI-RW-2-GW1        | LI-RW-2            | LI-RW-2-PI1         | LI-RW-2-PI2         | LI-RW-2-PI3         | LI-RW2-PI6          | LI-RW-2-PI9         | LI-RW-2-PI12       |
| Sampling Company                          |       |                    | DECI             | STANTEC     | STANTEC   | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC            | STANTEC             | DECI       | STANTEC            | STANTEC            | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC            |
| Laboratory                                |       |                    | PARAROCH         | CCGE        | PARAROCH  | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH           | PARAROCH            | PARAROCH   | CCGE               | PARAROCH           | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH           |
| Laboratory Work Order                     |       |                    | 12:1239          | E2363       | 141138    | 142196              | 142794              | 143439              | 144730              | 150382             | 151696              | 12:1239    | E2314              | 141138             | 142196              | 142794              | 143439              | 144730              | 150382              | 151696             |
| Laboratory Sample ID                      |       |                    | 12:1239-01       | E2363-01    | 141138-01 | 142196-09           | 142794-08           | 143439-01           | 144730-01           | 150382-01          | 151696-01           | 12:1239-02 | E2314-03           | 141138-02          | 142196-10           | 142794-07           | 143439-02           | 144730-02           | 150382-02           | 151696-02          |
| Sample Type                               | Units | TOGS               |                  |             |           |                     |                     |                     |                     |                    |                     |            |                    |                    |                     |                     |                     |                     |                     |                    |
| General Chemistry                         |       |                    |                  |             |           |                     |                     |                     |                     |                    |                     |            |                    |                    |                     |                     |                     |                     |                     |                    |
| Total Organic Carbon                      | µg/L  | n/v                | -                | -           | -         | 1060000             | 415000              | 43500               | 103000              | 9900               | 4500                | -          | -                  | 3200               | 553000              | 150000              | 259000              | 23900               | 9800                | 2700               |
| Metals                                    |       |                    |                  |             |           |                     |                     |                     |                     |                    |                     |            |                    |                    |                     |                     |                     |                     |                     |                    |
| Aluminum                                  | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 64.5               | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Antimony                                  | µg/L  | 3 <sup>B</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 12.5 U             | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Arsenic                                   | µg/L  | 25 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 5.000 U            | 10 U               | 10 U                | 10 U                | 10 U                | 10 U                | 10 U                | 10.0 U             |
| Barium                                    | µg/L  | 1000 <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 59.7 N             | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Beryllium                                 | µg/L  | 3 <sup>A</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 1.500 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Cadmium                                   | µg/L  | 5 <sup>B</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 1.500 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Calcium                                   | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 87300              | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Chromium (Total)                          | µg/L  | 50 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 2.500 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Cobalt                                    | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 7.500 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Copper                                    | µg/L  | 200 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 5.000 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Iron                                      | µg/L  | 300. <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 169                | 300                | 2220 <sup>B</sup>   | 1210 <sup>B</sup>   | 937 <sup>B</sup>    | 1430 <sup>B</sup>   | 498 <sup>B</sup>    | 1850 <sup>B</sup>  |
| Lead                                      | µg/L  | 25 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 9.61               | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Magnesium                                 | µg/L  | 35000 <sup>A</sup> | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 29500              | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Manganese                                 | µg/L  | 300. <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 305 J <sup>B</sup> | 120                | 233                 | 60.8                | 108                 | 187                 | 47.5                | 66.3               |
| Mercury                                   | µg/L  | 0.7 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 0.200 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Nickel                                    | µg/L  | 100 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 10.0 U             | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Potassium                                 | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 22600              | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Selenium                                  | µg/L  | 10 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 5.000 U N          | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Silver                                    | µg/L  | 50 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 2.500 U            | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Sodium                                    | µg/L  | 20000 <sup>B</sup> | -                | -           | -         | 146000 <sup>B</sup> | 331000 <sup>B</sup> | 137000 <sup>B</sup> | 146000 <sup>B</sup> | 85700 <sup>B</sup> | 175000 <sup>B</sup> | -          | 35600 <sup>B</sup> | 39100 <sup>B</sup> | 370000 <sup>B</sup> | 290000 <sup>B</sup> | 197000 <sup>B</sup> | 152000 <sup>B</sup> | 129000 <sup>B</sup> | 60600 <sup>B</sup> |
| Thallium                                  | µg/L  | 0.5 <sup>A</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 10.0 U             | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Vanadium                                  | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 10.0 U             | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Zinc                                      | µg/L  | 2000 <sup>A</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | 14.6               | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Volatile Organic Compounds                |       |                    |                  |             |           |                     |                     |                     |                     |                    |                     |            |                    |                    |                     |                     |                     |                     |                     |                    |
| Acetone                                   | µg/L  | 50 <sup>A</sup>    | 10.0 U           | 25 U        | 10.0 U    | 10.0 U              | 10.0 U              | 10.0 U              | 15.2                | 10.0 UJ            | 10.0 U              | 10.0 U     | 160 <sup>A</sup>   | 10.0 U             | 32.4                | 19.4                | 9.47 J              | 10.0 U              | 10.0 UJ             | 10.0 U             |
| Benzene                                   | µg/L  | 1 <sup>B</sup>     | 0.700 U          | 0.49 NJ     | 1 U       | 1 U                 | 1 U                 | 1 U                 | 0.561 J             | 1 U                | 0.700 U             | 0.700 U    | 5 U                | 1 U                | 1 U                 | 1 U                 | 1 U                 | 1 U                 | 1 U                 | 0.700 U            |
| Bromodichloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Bromoform (Tribromomethane)               | µg/L  | 50 <sup>A</sup>    | 5.00 U           | 5 U         | 5.00 U    | 5.00 U              | 5.00 U              | 5.00 UJ             | 5.00 U              | 5.00 U             | 5.00 U              | 5.00 U     | 5 U                | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 UJ             | 5.00 U              | 5.00 U              | 5.00 U             |
| Bromomethane (Methyl bromide)             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Butylbenzene, n-                          | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | -                  | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Butylbenzene, sec- (2-Phenylbutane)       | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | -                  | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Butylbenzene, tert-                       | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | -          | -                  | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Carbon Disulfide                          | µg/L  | 60 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Carbon Tetrachloride (Tetrachloromethane) | µg/L  | 5 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Chlorobenzene (Monochlorobenzene)         | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Chlorobromomethane                        | µg/L  | 5.. <sup>B</sup>   | -                | 5 U         | 5.00 U    | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 U              | -          | 5 U                | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             |
| Chloroethane (Ethyl Chloride)             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Chloroethyl Vinyl Ether, 2-               | µg/L  | n/v                | 10.0 U R         | -           | -         | -                   | -                   | -                   | -                   | -                  | -                   | 10.0 U R   | -                  | -                  | -                   | -                   | -                   | -                   | -                   | -                  |
| Chloroform (Trichloromethane)             | µg/L  | 7 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 0.67 J             | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Chloromethane                             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Cyclohexane                               | µg/L  | n/v                | -                | 4.5 NJ      | 10.0 U    | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U              | -          | 5 U                | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             |
| Dibromo-3-Chloropropane, 1,2- (DBCP)      | µg/L  | 0.04 <sup>B</sup>  | -                | 5 U         | 10.0 U    | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U              | -          | 5 U                | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             |
| Dibromochloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U             |
| Dichlorobenzene, 1,2-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Dichlorobenzene, 1,3-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Dichlorobenzene, 1,4-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | 2.00 U     | 5 U                | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |
| Dichlorodifluoromethane (Freon 12)        | µg/L  | 5.. <sup>B</sup>   | -                | 5 U         | 2.00 U    | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U              | -          | 5 U                | 2.00 U             | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area<br>Sample Location<br>Sample Date<br><br>Sample ID<br>Sampling Company<br>Laboratory<br>Laboratory Work Order<br>Laboratory Sample ID<br><br>Sample Type |      |                               | On-Site Building                    |                   |                   |                   |                      |                      |                     |             |                   |                   |                    |                   |                     |                   |                   |                      |             |              |  |  |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------|------|-------------------------------|-------------------------------------|-------------------|-------------------|-------------------|----------------------|----------------------|---------------------|-------------|-------------------|-------------------|--------------------|-------------------|---------------------|-------------------|-------------------|----------------------|-------------|--------------|--|--|
|                                                                                                                                                               |      |                               | RW-1                                |                   |                   |                   |                      |                      |                     |             |                   |                   | RW-2               |                   |                     |                   |                   |                      |             |              |  |  |
|                                                                                                                                                               |      |                               | 23-Mar-12                           | 23-May-13         | 26-Mar-14         | 29-May-14         | 1-Jul-14             | 8-Aug-14             | 29-Oct-14           | 3-Feb-15    | 5-May-15          | 23-Mar-12         | 21-May-13          | 26-Mar-14         | 29-May-14           | 1-Jul-14          | 8-Aug-14          | 29-Oct-14            | 3-Feb-15    | 5-May-15     |  |  |
|                                                                                                                                                               |      |                               | RW-1                                | LI-RW-1-GW1       | LI-RW-1           | LI-RW-1-PI1       | LI-RW-1-PI2          | LI-RW-1-PI3          | LI-RW1-PI6          | LI-RW-1-PI9 | LI-RW-1-PI12      | RW-2              | LI-RW-2-GW1        | LI-RW-2           | LI-RW-2-PI1         | LI-RW-2-PI2       | LI-RW-2-PI3       | LI-RW2-PI6           | LI-RW-2-PI9 | LI-RW-2-PI12 |  |  |
|                                                                                                                                                               |      |                               | DECI                                | STANTEC           | STANTEC           | STANTEC           | STANTEC              | STANTEC              | STANTEC             | STANTEC     | STANTEC           | DECI              | STANTEC            | STANTEC           | STANTEC             | STANTEC           | STANTEC           | STANTEC              | STANTEC     | STANTEC      |  |  |
|                                                                                                                                                               |      |                               | PARAROCH                            | CCGE              | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH             | PARAROCH            | PARAROCH    | PARAROCH          | PARAROCH          | CCGE               | PARAROCH          | PARAROCH            | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH    | PARAROCH     |  |  |
|                                                                                                                                                               |      |                               | 12:1239                             | E2363             | 141138            | 142196            | 142794               | 143439               | 144730              | 150382      | 151696            | 12:1239           | E2314              | 141138            | 142196              | 142794            | 143439            | 144730               | 150382      | 151696       |  |  |
|                                                                                                                                                               |      |                               | 12:1239-01                          | E2363-01          | 141138-01         | 142196-09         | 142794-08            | 143439-01            | 144730-01           | 150382-01   | 151696-01         | 12:1239-02        | E2314-03           | 141138-02         | 142196-10           | 142794-07         | 143439-02         | 144730-02            | 150382-02   | 151696-02    |  |  |
|                                                                                                                                                               |      |                               | Units                               | TOGS              |                   |                   |                      |                      |                     |             |                   |                   |                    |                   |                     |                   |                   |                      |             |              |  |  |
|                                                                                                                                                               |      |                               | Volatile Organic Compounds (cont'd) |                   |                   |                   |                      |                      |                     |             |                   |                   |                    |                   |                     |                   |                   |                      |             |              |  |  |
| Dichloroethane, 1,1-                                                                                                                                          | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dichloroethane, 1,2-                                                                                                                                          | µg/L | 0.6 <sup>B</sup>              | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dichloroethene, 1,1-                                                                                                                                          | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 1 J                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dichloroethene, cis-1,2-                                                                                                                                      | µg/L | 5 <sup>B</sup>                | 6.88 <sup>B</sup>                   | 14.5 <sup>B</sup> | 5.57 <sup>B</sup> | 4.53              | 4.71                 | 8.12 <sup>B</sup>    | 2.00 U              | 2.00 U      | 2.00 U            | 26.6 <sup>B</sup> | 360 D <sup>B</sup> | 38.8 <sup>B</sup> | 55.7 <sup>B</sup>   | 51.3 <sup>B</sup> | 23.6 <sup>B</sup> | 87.7 <sup>B</sup>    | 4.37        | 2.00 U       |  |  |
| Dichloroethene, trans-1,2-                                                                                                                                    | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 4.2 J             | 2.00 U            | 2.00 U            | 1.03 J               | 2.00 U               | 1.34 J              | 2.00 U      | 2.00 U            | 2.43              | 11.4 <sup>B</sup>  | 2.39              | 3.06                | 2.50              | 3.57              | 12.8 <sup>B</sup>    | 2.00 U      | 1.17 J       |  |  |
| Dichloropropane, 1,2-                                                                                                                                         | µg/L | 1 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dichloropropene, cis-1,3-                                                                                                                                     | µg/L | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dichloropropene, trans-1,3-                                                                                                                                   | µg/L | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Dioxane, 1,4-                                                                                                                                                 | µg/L | n/v                           | -                                   | 100 U R           | 20.0 U R          | 20.0 U R          | 20.0 U R             | 20.0 U R             | 20.0 U              | 20.0 U      | 20.0 U            | -                 | 100 U              | 20.0 U R          | 20.0 U R            | 20.0 U R          | 20.0 U R          | 20.0 U               | 20.0 U      | 20.0 U       |  |  |
| Ethylbenzene                                                                                                                                                  | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Ethylene Dibromide (Dibromoethane, 1,2-)                                                                                                                      | µg/L | 0.0006 <sup>B</sup>           | -                                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Hexanone, 2- (Methyl Butyl Ketone)                                                                                                                            | µg/L | 50 <sup>A</sup>               | 5.00 U                              | 25 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | 5.00 U            | 25 U               | 5.00 U            | 5.00 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U      | 5.00 U       |  |  |
| Isopropylbenzene                                                                                                                                              | µg/L | 5 <sup>B</sup>                | -                                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Isopropyltoluene, p- (Cymene)                                                                                                                                 | µg/L | 5 <sup>B</sup>                | -                                   | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Methyl Acetate                                                                                                                                                | µg/L | n/v                           | -                                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 3.03 NJ              | 2.00 U      | 2.00 U       |  |  |
| Methyl Ethyl Ketone (MEK)                                                                                                                                     | µg/L | 50 <sup>A</sup>               | 10.0 U                              | 25 U              | 10.0 UJ           | 6.42 J            | 87.3 J <sup>A</sup>  | 9.42 NJ              | 57.3 J <sup>A</sup> | 10.0 UJ     | 10.0 U            | 10.0 U            | 110 <sup>A</sup>   | 10.0 UJ           | 175 NJ <sup>A</sup> | 29.3 J            | 38.1              | 10.2 J               | 10.0 UJ     | 10.0 U       |  |  |
| Methyl Isobutyl Ketone (MIBK)                                                                                                                                 | µg/L | n/v                           | 5.00 U                              | 25 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | 5.00 U            | 25 U               | 5.00 U            | 5.00 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U      | 5.00 U       |  |  |
| Methyl tert-butyl ether (MTBE)                                                                                                                                | µg/L | 10 <sup>A</sup>               | -                                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 2.4 J              | 1.08 J            | 1.61 NJ             | 2.00 U            | 1.92 J            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Methylcyclohexane                                                                                                                                             | µg/L | n/v                           | -                                   | 3.1 J             | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Methylene Chloride (Dichloromethane)                                                                                                                          | µg/L | 5 <sup>B</sup>                | 5.00 U                              | 5 U               | 2.84 JB           | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | 5.00 U            | 5 U                | 3.76 JB           | 5.00 U              | 5.00 U            | 5.00 U            | 7.55 <sup>B</sup>    | 5.00 U      | 5.00 U       |  |  |
| Naphthalene                                                                                                                                                   | µg/L | 10 <sup>B</sup>               | -                                   | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Propylbenzene, n-                                                                                                                                             | µg/L | 5 <sup>B</sup>                | -                                   | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Styrene                                                                                                                                                       | µg/L | 5 <sup>B</sup>                | 5.00 U                              | 5 U               | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | 5.00 U            | 5 U                | 5.00 U            | 5.00 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U      | 5.00 U       |  |  |
| Tetrachloroethane, 1,1,2,2-                                                                                                                                   | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Tetrachloroethene (PCE)                                                                                                                                       | µg/L | 5 <sup>B</sup>                | 6.72 <sup>B</sup>                   | 3.6 J             | 5.35 <sup>B</sup> | 10.1 <sup>B</sup> | 6.14 <sup>B</sup>    | 2.65                 | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 110 <sup>B</sup>   | 4.44              | 3.08                | 1.42 J            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Toluene                                                                                                                                                       | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Trichlorobenzene, 1,2,3-                                                                                                                                      | µg/L | 5 <sup>B</sup>                | -                                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | -                 | 5 U                | 5.00 U            | 5.00 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U      | 5.00 U       |  |  |
| Trichlorobenzene, 1,2,4-                                                                                                                                      | µg/L | 5 <sup>B</sup>                | -                                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U               | 5.00 U              | 5.00 U      | 5.00 U            | -                 | 5 U                | 5.00 U            | 5.00 U              | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U      | 5.00 U       |  |  |
| Trichloroethane, 1,1,1-                                                                                                                                       | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Trichloroethane, 1,1,2-                                                                                                                                       | µg/L | 1 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Trichloroethene (TCE)                                                                                                                                         | µg/L | 5 <sup>B</sup>                | 7.15 <sup>B</sup>                   | 8.1 <sup>B</sup>  | 4.02              | 6.09 <sup>B</sup> | 4.52                 | 5.49 <sup>B</sup>    | 2.00 U              | 2.00 U      | 2.00 U            | 9.19 <sup>B</sup> | 76.4 <sup>B</sup>  | 27.6 <sup>B</sup> | 21.5 <sup>B</sup>   | 6.31 <sup>B</sup> | 2.39              | 1.05 J               | 2.00 U      | 2.00 U       |  |  |
| Trichlorofluoromethane (Freon 11)                                                                                                                             | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Trichlorotrifluoroethane (Freon 113)                                                                                                                          | µg/L | 5 <sup>B</sup>                | -                                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | -                 | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Trimethylbenzene, 1,2,4-                                                                                                                                      | µg/L | 5 <sup>B</sup>                | -                                   | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Trimethylbenzene, 1,3,5-                                                                                                                                      | µg/L | 5 <sup>B</sup>                | -                                   | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Vinyl Acetate                                                                                                                                                 | µg/L | n/v                           | 5.00 U                              | -                 | -                 | -                 | -                    | -                    | -                   | -           | -                 | 5.00 U            | -                  | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |
| Vinyl chloride                                                                                                                                                | µg/L | 2 <sup>B</sup>                | 3.99 <sup>B</sup>                   | 7.7 <sup>B</sup>  | 2.00 U            | 1.45 NJ           | 4.61 NJ <sup>B</sup> | 5.29 NJ <sup>B</sup> | 2.00 U              | 2.00 U      | 2.98 <sup>B</sup> | 2.00 U            | 5.9 <sup>B</sup>   | 1.24 J            | 1.64 NJ             | 7.48 <sup>B</sup> | 56.4 <sup>B</sup> | 23.9 NJ <sup>B</sup> | 1.17 NJ     | 1.26 J       |  |  |
| Xylene, m & p-                                                                                                                                                | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 10 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 10 U               | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Xylene, o-                                                                                                                                                    | µg/L | 5 <sup>B</sup>                | 2.00 U                              | 5 U               | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U               | 2.00 U              | 2.00 U      | 2.00 U            | 2.00 U            | 5 U                | 2.00 U            | 2.00 U              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U      | 2.00 U       |  |  |
| Total VOC                                                                                                                                                     | µg/L | n/v                           | 24.74                               | 46.19             | 17.78             | 28.59             | 108.31               | 30.97                | 74.40               | 1.43        | 2.98              | 38.22             | 837.77             | 79.31             | 293.99              | 117.71            | 135.45            | 146.23               | 5.54        | 2.43         |  |  |
| VOC Tentatively Identified Compounds                                                                                                                          |      |                               |                                     |                   |                   |                   |                      |                      |                     |             |                   |                   |                    |                   |                     |                   |                   |                      |             |              |  |  |
| Total VOC TICs                                                                                                                                                | µg/L | n/v                           | -                                   | 4,900 J           | -                 | -                 | -                    | -                    | -                   | -           | -                 | -                 | 770,000 J          | -                 | -                   | -                 | -                 | -                    | -           | -            |  |  |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                      |       |                    | On-Site Building |             |           |                     |                     |                     |                     |                     |                    | Off-Site Locations  |             |                    |                     |                     |                     |                    |                    |                    |
|-------------------------------------------|-------|--------------------|------------------|-------------|-----------|---------------------|---------------------|---------------------|---------------------|---------------------|--------------------|---------------------|-------------|--------------------|---------------------|---------------------|---------------------|--------------------|--------------------|--------------------|
| Sample Location                           |       |                    | RW-3             |             |           |                     |                     |                     |                     |                     |                    | RW-5                |             |                    |                     |                     |                     |                    |                    |                    |
| Sample Date                               |       |                    | 23-Mar-12        | 22-May-13   | 26-Mar-14 | 29-May-14           | 1-Jul-14            | 7-Aug-14            | 29-Oct-14           | 3-Feb-15            | 5-May-15           | 25-Apr-12           | 21-May-13   | 27-Mar-14          | 29-May-14           | 2-Jul-14            | 7-Aug-14            | 28-Oct-14          | 3-Feb-15           | 4-May-15           |
| Sample ID                                 |       |                    | RW-3             | LI-RW-3-GW1 | LI-RW-3   | LI-RW-3-PI1         | LI-RW-3-PI2         | LI-RW-3-PI3         | LI-RW3-PI6          | LI-RW-3-PI9         | LI-RW-3-PI12       | RW-5                | LI-RW-5-GW1 | LI-RW-5            | LI-RW-5-PI1         | LI-RW-5-PI2         | LI-RW-5-PI3         | LI-RW5-PI6         | LI-RW-5-PI9        | LI-RW-5-PI12       |
| Sampling Company                          |       |                    | DECI             | STANTEC     | STANTEC   | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC            | DECI                | STANTEC     | STANTEC            | STANTEC             | STANTEC             | STANTEC             | STANTEC            | STANTEC            | STANTEC            |
| Laboratory                                |       |                    | PARAROCH         | CCGE        | PARAROCH  | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH           | PARAROCH            | CCGE        | PARAROCH           | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH           | PARAROCH           | PARAROCH           |
| Laboratory Work Order                     |       |                    | 12:1239          | E2342       | 141138    | 142196              | 142794              | 143439              | 144730              | 150382              | 151696             | 12:1770             | E2314       | 141138             | 142196              | 142794              | 143439              | 144730             | 150382             | 151696             |
| Laboratory Sample ID                      |       |                    | 12:1239-03       | E2342-01    | 141138-03 | 142196-11           | 142794-06           | 143439-03           | 144730-03           | 150382-03           | 151696-03          | 12:1770-02          | E2314-06    | 141138-05          | 142196-14           | 142794-13           | 143439-05           | 144730-05          | 150382-04          | 151696-05          |
| Sample Type                               | Units | TOGS               |                  |             |           |                     |                     |                     |                     |                     |                    |                     |             |                    |                     |                     |                     |                    |                    |                    |
| General Chemistry                         |       |                    |                  |             |           |                     |                     |                     |                     |                     |                    |                     |             |                    |                     |                     |                     |                    |                    |                    |
| Total Organic Carbon                      | µg/L  | n/v                | -                | -           | -         | 229000              | 87900               | 12700               | 11000               | 10300               | 6100               | -                   | -           | 3300               | 141000              | 299000              | 86700               | 8700               | 4600 J+            | 2200               |
| Metals                                    |       |                    |                  |             |           |                     |                     |                     |                     |                     |                    |                     |             |                    |                     |                     |                     |                    |                    |                    |
| Aluminum                                  | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Antimony                                  | µg/L  | 3 <sup>B</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Arsenic                                   | µg/L  | 25 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | 10 U               | 10 U                | 10 U                | 10 U                | 10 U               | 10 U               | 10.0 U             |
| Barium                                    | µg/L  | 1000 <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Beryllium                                 | µg/L  | 3 <sup>A</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Cadmium                                   | µg/L  | 5 <sup>B</sup>     | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Calcium                                   | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Chromium (Total)                          | µg/L  | 50 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Cobalt                                    | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Copper                                    | µg/L  | 200 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Iron                                      | µg/L  | 300. <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | 100 U              | 2500 <sup>B</sup>   | 6250 <sup>B</sup>   | 6000 <sup>B</sup>   | 4420 <sup>B</sup>  | 4760 <sup>B</sup>  | 9910 <sup>B</sup>  |
| Lead                                      | µg/L  | 25 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Magnesium                                 | µg/L  | 35000 <sup>A</sup> | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Manganese                                 | µg/L  | 300. <sup>B</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | 69.2               | 69.1                | 102                 | 60.4 B              | 47.8               | 25.7               | 29.8               |
| Mercury                                   | µg/L  | 0.7 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Nickel                                    | µg/L  | 100 <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Potassium                                 | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Selenium                                  | µg/L  | 10 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Silver                                    | µg/L  | 50 <sup>B</sup>    | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Sodium                                    | µg/L  | 20000 <sup>B</sup> | -                | -           | -         | 252000 <sup>B</sup> | 199000 <sup>B</sup> | 103000 <sup>B</sup> | 125000 <sup>B</sup> | 120000 <sup>B</sup> | 85300 <sup>B</sup> | -                   | -           | 39500 <sup>B</sup> | 242000 <sup>B</sup> | 312000 <sup>B</sup> | 164000 <sup>B</sup> | 85200 <sup>B</sup> | 66600 <sup>B</sup> | 44600 <sup>B</sup> |
| Thallium                                  | µg/L  | 0.5 <sup>A</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Vanadium                                  | µg/L  | n/v                | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Zinc                                      | µg/L  | 2000 <sup>A</sup>  | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Volatile Organic Compounds                |       |                    |                  |             |           |                     |                     |                     |                     |                     |                    |                     |             |                    |                     |                     |                     |                    |                    |                    |
| Acetone                                   | µg/L  | 50 <sup>A</sup>    | 10.0 U           | 25 U        | 10.0 U    | 132 <sup>A</sup>    | 43.2 J              | 47.6 J              | 10.0 U              | 10.0 UJ             | 10.0 U             | 10.0 UJ             | 2.6 J       | 10.0 U             | 10.0 U              | 7.44 J              | 10.0 U              | 10.0 U             | 10.0 UJ            | 10.0 U             |
| Benzene                                   | µg/L  | 1 <sup>B</sup>     | 0.700 U          | 5 U         | 1 U       | 5 U                 | 1 U                 | 1 U                 | 1 U                 | 1 U                 | 0.700 U            | 1.13 J <sup>B</sup> | 5 U         | 1 U                | 0.737 J             | 0.358 J             | 1 U                 | 0.507 J            | 1 U                | 0.700 U            |
| Bromodichloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Bromoform (Tribromomethane)               | µg/L  | 50 <sup>A</sup>    | 5.00 U           | 5 U         | 5.00 U    | 25.0 U              | 25.0 U              | 5.00 UJ             | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 UJ             | 5 U         | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 UJ             | 5.00 U             | 5.00 U             | 5.00 U             |
| Bromomethane (Methyl bromide)             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Butylbenzene, n-                          | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Butylbenzene, sec- (2-Phenylbutane)       | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Butylbenzene, tert-                       | µg/L  | 5.. <sup>B</sup>   | -                | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | -                   | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Carbon Disulfide                          | µg/L  | 60 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Carbon Tetrachloride (Tetrachloromethane) | µg/L  | 5 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Chlorobenzene (Monochlorobenzene)         | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Chlorobromomethane                        | µg/L  | 5.. <sup>B</sup>   | -                | 5 U         | 5.00 U    | 25.0 U              | 25.0 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             | -                   | 5 U         | 5.00 U             | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U             | 5.00 U             | 5.00 U             |
| Chloroethane (Ethyl Chloride)             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Chloroethyl Vinyl Ether, 2-               | µg/L  | n/v                | 10.0 U R         | -           | -         | -                   | -                   | -                   | -                   | -                   | -                  | 10.0 U R            | -           | -                  | -                   | -                   | -                   | -                  | -                  | -                  |
| Chloroform (Trichloromethane)             | µg/L  | 7 <sup>B</sup>     | 3.78             | 3.9 J       | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Chloromethane                             | µg/L  | 5.. <sup>B</sup>   | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Cyclohexane                               | µg/L  | n/v                | -                | 5 U         | 10.0 U    | 50.0 U              | 50.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | -                   | 5 U         | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U             | 10.0 U             |
| Dibromo-3-Chloropropane, 1,2- (DBCP)      | µg/L  | 0.04 <sup>B</sup>  | -                | 5 U         | 10.0 U    | 50.0 U              | 50.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | -                   | 5 U         | 10.0 U             | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U             | 10.0 U             | 10.0 U             |
| Dibromochloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 UJ             | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,2-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,3-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorobenzene, 1,4-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U           | 5 U         | 2.00 U    | 10.0 U              | 10.0 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 UJ             | 5 U         | 2.00 U             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |
| Dichlorodifluoromethane (Freon 12)        | µg/L  | 5.. <sup>B</sup>   | -                | 5 U         | 2.00 U    | 10.0 U              | 10.0 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U             | -                   | 5 U         | 2.00 U             | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U             | 2.00 U             | 2.00 U             |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                     |       |                               | On-Site Building  |                    |                   |                   |                     |                      |                   |                      |                   | Off-Site Locations  |                   |                   |                   |                   |                   |                   |                      |              |  |
|------------------------------------------|-------|-------------------------------|-------------------|--------------------|-------------------|-------------------|---------------------|----------------------|-------------------|----------------------|-------------------|---------------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|----------------------|--------------|--|
| Sample Location                          |       |                               | RW-3              |                    |                   |                   |                     |                      |                   |                      |                   | RW-5                |                   |                   |                   |                   |                   |                   |                      |              |  |
| Sample Date                              |       |                               | 23-Mar-12         | 22-May-13          | 26-Mar-14         | 29-May-14         | 1-Jul-14            | 7-Aug-14             | 29-Oct-14         | 3-Feb-15             | 5-May-15          | 25-Apr-12           | 21-May-13         | 27-Mar-14         | 29-May-14         | 2-Jul-14          | 7-Aug-14          | 28-Oct-14         | 3-Feb-15             | 4-May-15     |  |
| Sample ID                                |       |                               | RW-3              | LI-RW-3-GW1        | LI-RW-3           | LI-RW-3-PI1       | LI-RW-3-PI2         | LI-RW-3-PI3          | LI-RW3-PI6        | LI-RW-3-PI9          | LI-RW-3-PI12      | RW-5                | LI-RW-5-GW1       | LI-RW-5           | LI-RW-5-PI1       | LI-RW-5-PI2       | LI-RW-5-PI3       | LI-RW5-PI6        | LI-RW-5-PI9          | LI-RW-5-PI12 |  |
| Sampling Company                         |       |                               | DECI              | STANTEC            | STANTEC           | STANTEC           | STANTEC             | STANTEC              | STANTEC           | STANTEC              | STANTEC           | DECI                | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC              | STANTEC      |  |
| Laboratory                               |       |                               | PARAROCH          | CCGE               | PARAROCH          | PARAROCH          | PARAROCH            | PARAROCH             | PARAROCH          | PARAROCH             | PARAROCH          | PARAROCH            | CCGE              | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH     |  |
| Laboratory Work Order                    |       |                               | 12:1239           | E2342              | 141138            | 142196            | 142794              | 143439               | 144730            | 150382               | 151696            | 12:1770             | E2314             | 141138            | 142196            | 142794            | 143439            | 144730            | 150382               | 151696       |  |
| Laboratory Sample ID                     |       |                               | 12:1239-03        | E2342-01           | 141138-03         | 142196-11         | 142794-06           | 143439-03            | 144730-03         | 150382-03            | 151696-03         | 12:1770-02          | E2314-06          | 141138-05         | 142196-14         | 142794-13         | 143439-05         | 144730-05         | 150382-04            | 151696-05    |  |
| Sample Type                              | Units | TOGS                          |                   |                    |                   |                   |                     |                      |                   |                      |                   |                     |                   |                   |                   |                   |                   |                   |                      |              |  |
| Volatile Organic Compounds (cont'd)      |       |                               |                   |                    |                   |                   |                     |                      |                   |                      |                   |                     |                   |                   |                   |                   |                   |                   |                      |              |  |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sup>n,B</sup>              | 81.8 <sup>B</sup> | 130 <sup>B</sup>   | 3.77              | 30.1 <sup>B</sup> | 90.5 <sup>B</sup>   | 143 <sup>B</sup>     | 3.35              | 1.40 J               | 1.23 J            | 49.5 J <sup>B</sup> | 18.2 <sup>B</sup> | 7.64 <sup>B</sup> | 32.7 <sup>B</sup> | 45.7 <sup>B</sup> | 46.0 <sup>B</sup> | 132 <sup>B</sup>  | 8.81 <sup>B</sup>    | 4.52         |  |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sup>n,B</sup>              | 10.2 <sup>B</sup> | 18.8 <sup>B</sup>  | 2.00 U            | 10.0 U            | 7.12 J <sup>B</sup> | 3.16                 | 4.47              | 6.02 <sup>B</sup>    | 3.63              | 5.63 J <sup>B</sup> | 2.2 J             | 1.10 J            | 2.92              | 1.89 J            | 1.32 J            | 3.78              | 2.00 U               | 2.00 U       |  |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Dioxane, 1,4-                            | µg/L  | n/v                           | -                 | 100 U R            | 20.0 U R          | 100 U R           | 100 U R             | 20.0 U R             | 20.0 U            | 20.0 U               | 20.0 U            | -                   | 100 U             | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U            | 20.0 U               | 20.0 U       |  |
| Ethylbenzene                             | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | -                 | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | 5.00 U            | 25 U               | 5.00 U            | 25.0 U            | 25.0 U              | 5.00 U               | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 UJ             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Isopropylbenzene                         | µg/L  | 5 <sup>n,B</sup>              | -                 | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sup>n,B</sup>              | -                 | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Methyl Acetate                           | µg/L  | n/v                           | -                 | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.87              | 2.00 U               | 2.00 U            | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | 10.0 U            | 25 U               | 10.0 UJ           | 40.4 <sup>A</sup> | 139 J <sup>A</sup>  | 60.0 <sup>A</sup>    | 10.0 UJ           | 10.0 UJ              | 10.0 U            | 10.0 UJ             | 25 U              | 10.0 UJ           | 10.0 U            | 43.1 J            | 10.8              | 10.0 UJ           | 10.0 UJ              | 10.0 U       |  |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | 5.00 U            | 25 U               | 5.00 U            | 25.0 U            | 25.0 U              | 5.00 U               | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 UJ             | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | -                 | 7.1                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 3.43                 | 2.12              | -                   | 1.3 J             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Methylcyclohexane                        | µg/L  | n/v                           | -                 | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sup>n,B</sup>              | 5.00 U            | 5 U                | 4.04 JB           | 25.0 U            | 25.0 U              | 5.00 U               | 6.12 <sup>B</sup> | 5.00 U               | 5.00 U            | 5.00 UJ             | 5 U               | 4.53 JB           | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -                 | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Propylbenzene, n-                        | µg/L  | 5 <sup>n,B</sup>              | -                 | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Styrene                                  | µg/L  | 5 <sup>n,B</sup>              | 5.00 U            | 5 U                | 5.00 U            | 25.0 U            | 25.0 U              | 5.00 U               | 5.00 U            | 5.00 U               | 5.00 U            | 5.00 UJ             | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sup>n,B</sup>              | 2.81              | 7.8 <sup>B</sup>   | 2.36              | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 12.2 J <sup>B</sup> | 5.6 <sup>B</sup>  | 2.75              | 11.2 <sup>B</sup> | 2.44              | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Toluene                                  | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sup>n,B</sup>              | -                 | 5 U                | 5.00 U            | 25.0 U            | 25.0 U              | 5.00 U               | 5.00 U            | 5.00 U               | 5.00 U            | -                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sup>n,B</sup>              | -                 | 5 U                | 5.00 U            | 25.0 U            | 25.0 U              | 5.00 U               | 5.00 U            | 5.00 U               | 5.00 U            | -                   | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U       |  |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Trichloroethene (TCE)                    | µg/L  | 5 <sup>n,B</sup>              | 125 <sup>B</sup>  | 320 D <sup>B</sup> | 10.5 <sup>B</sup> | 83.9 <sup>B</sup> | 36.6 <sup>B</sup>   | 2.00 U               | 2.00 U            | 2.00 U               | 1.04 J            | 48.5 J <sup>B</sup> | 25.2 <sup>B</sup> | 6.65 <sup>B</sup> | 40.0 <sup>B</sup> | 14.2 <sup>B</sup> | 1.10 J            | 2.76              | 2.00 U               | 2.00 U       |  |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sup>n,B</sup>              | -                 | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | -                   | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sup>n,B</sup>              | -                 | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sup>n,B</sup>              | -                 | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Vinyl Acetate                            | µg/L  | n/v                           | 5.00 U            | -                  | -                 | -                 | -                   | -                    | -                 | -                    | -                 | 5.00 UJ             | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 2.00 U            | 3 J <sup>B</sup>   | 2.00 U            | 10.0 U            | 18.1 <sup>B</sup>   | 10.1 NJ <sup>B</sup> | 22.5 <sup>B</sup> | 4.14 NJ <sup>B</sup> | 4.41 <sup>B</sup> | 2.93 J <sup>B</sup> | 0.6 J             | 2.00 U            | 2.00 U            | 1.28 NJ           | 3.76 <sup>B</sup> | 12.8 <sup>B</sup> | 2.30 NJ <sup>B</sup> | 1.61 J       |  |
| Xylene, m & p-                           | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 10 U               | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 10 U              | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Xylene, o-                               | µg/L  | 5 <sup>n,B</sup>              | 2.00 U            | 5 U                | 2.00 U            | 10.0 U            | 10.0 U              | 2.00 U               | 2.00 U            | 2.00 U               | 2.00 U            | 2.00 UJ             | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U       |  |
| Total VOC                                | µg/L  | n/v                           | 223.59            | 490.6              | 20.67             | 650               | 334.52              | 263.86               | 39.31             | 14.99                | 12.43             | 119.89              | 55.7              | 22.67             | 87.557            | 116.408           | 62.98             | 151.85            | 11.11                | 6.13         |  |
| VOC Tentatively Identified Compounds     |       |                               |                   |                    |                   |                   |                     |                      |                   |                      |                   |                     |                   |                   |                   |                   |                   |                   |                      |              |  |
| Total VOC TICs                           | µg/L  | n/v                           | -                 | 2.5 U              | -                 | -                 | -                   | -                    | -                 | -                    | -                 | -                   | 5.500 J           | -                 | -                 | -                 | -                 | -                 | -                    | -            |  |

See last page for notes.



Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                     |       |                               | Off-Site Locations  |                   |                    |                   |                     |                   |                     |                     |                   |                     |                     |            |                  |           |             |             |             |                      |                   |                   |
|------------------------------------------|-------|-------------------------------|---------------------|-------------------|--------------------|-------------------|---------------------|-------------------|---------------------|---------------------|-------------------|---------------------|---------------------|------------|------------------|-----------|-------------|-------------|-------------|----------------------|-------------------|-------------------|
| Sample Location                          |       |                               |                     |                   |                    |                   |                     |                   |                     |                     |                   |                     |                     |            |                  |           |             |             |             |                      |                   |                   |
| Sample Date                              |       |                               | 25-Apr-12           | 4-May-12          | 20-May-13          | 27-Mar-14         | 28-May-14           | 1-Jul-14          | 7-Aug-14            | 7-Aug-14            | 28-Oct-14         | 4-Feb-15            | 4-May-15            | 12-Jun-12  | 20-May-13        | 27-Mar-14 | 28-May-14   | 1-Jul-14    | 7-Aug-14    | 28-Oct-14            | 4-Feb-15          | 4-May-15          |
| Sample ID                                |       |                               | RW-6                | RW-6              | LI-RW-6-GW1        | LI-RW-6           | LI-RW-6-PI1         | LI-RW-6-PI2       | LI-RW-6-PI3         | LI-FD-PI3           | LI-RW6-PI6        | LI-RW-6-PI9         | LI-RW-6-PI12        | RW-7       | LI-RW-7-GW1      | LI-RW-7   | LI-RW-7-PI1 | LI-RW-7-PI2 | LI-RW-7-PI3 | LI-RW7-PI6           | LI-RW-7-PI9       | LI-RW-7-PI12      |
| Sampling Company                         |       |                               | DECI                | DECI              | STANTEC            | STANTEC           | STANTEC             | STANTEC           | STANTEC             | STANTEC             | STANTEC           | STANTEC             | STANTEC             | DECI       | STANTEC          | STANTEC   | STANTEC     | STANTEC     | STANTEC     | STANTEC              | STANTEC           | STANTEC           |
| Laboratory                               |       |                               | PARAROCH            | PARAROCH          | CCGE               | PARAROCH          | PARAROCH            | PARAROCH          | PARAROCH            | PARAROCH            | PARAROCH          | PARAROCH            | PARAROCH            | PARAROCH   | CCGE             | PARAROCH  | PARAROCH    | PARAROCH    | PARAROCH    | PARAROCH             | PARAROCH          | PARAROCH          |
| Laboratory Work Order                    |       |                               | 12:1770             | 12:1927           | E2301              | 141138            | 142196              | 142794            | 143439              | 143439              | 144730            | 150382              | 151696              | 12:2486    | E2301            | 141138    | 142196      | 142794      | 143439      | 144730               | 150382            | 151696            |
| Laboratory Sample ID                     |       |                               | 12:1770-03          | 12:1927-01        | E2301-01           | 141138-06         | 142196-02           | 142794-03         | 143439-06           | 143439-13           | 144730-06         | 150382-09           | 151696-06           | 12:2486-02 | E2301-02         | 141138-07 | 142196-01   | 142794-02   | 143439-07   | 144730-07            | 150382-10         | 151696-07         |
| Sample Type                              | Units | TOGS                          | Field Duplicate     |                   |                    |                   |                     |                   |                     |                     |                   |                     |                     |            |                  |           |             |             |             |                      |                   |                   |
| Volatile Organic Compounds (cont'd)      |       |                               |                     |                   |                    |                   |                     |                   |                     |                     |                   |                     |                     |            |                  |           |             |             |             |                      |                   |                   |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 59.8 J <sup>B</sup> | 63.1 <sup>B</sup> | 47.3 <sup>B</sup>  | 81.9 <sup>B</sup> | 670 <sup>B</sup>    | 86.7 <sup>B</sup> | 3980 <sup>B</sup>   | 4070 <sup>B</sup>   | 2730 <sup>B</sup> | 687 <sup>B</sup>    | 373 <sup>B</sup>    | 4.28       | 8.2 <sup>B</sup> | 2.35      | 2.65        | 2.43        | 2.96        | 4.44                 | 1.33 J            | 2.00 U            |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 1.1 J              | 40.0 U            | 76.1 <sup>B</sup>   | 3.31              | 76.6 J <sup>B</sup> | 77.6 J <sup>B</sup> | 100 U             | 6.64 J <sup>B</sup> | 6.18 J <sup>B</sup> | 2.00 U     | 0.92 J           | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Dioxane, 1,4-                            | µg/L  | n/v                           | -                   | -                 | 100 U R            | 400 U R           | 400 U R             | 20.0 U R          | 1000 U R            | 1000 U R            | 1000 U            | 100 U               | 100 U               | -          | 100 U R          | 20.0 U R  | 20.0 U R    | 20.0 U R    | 20.0 U R    | 20.0 U               | 20.0 U            | 20.0 U            |
| Ethylbenzene                             | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | -                   | -                 | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | 5.00 UJ             | 50.0 U            | 25 U               | 100 U             | 100 U               | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | -          | 25 U             | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Isopropylbenzene                         | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Methyl Acetate                           | µg/L  | n/v                           | -                   | -                 | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | 10.0 UJ             | 100 U             | 25 U               | 200 UJ            | 200 U               | 13.3 J            | 500 U               | 500 U               | 500 UJ            | 50.0 UJ             | 50.0 U              | -          | 25 U             | 10.0 UJ   | 10.0 U      | 10.0 UJ     | 10.0 U      | 10.0 UJ              | 10.0 UJ           | 10.0 U            |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | 5.00 UJ             | 50.0 U            | 25 U               | 100 U             | 100 U               | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | -          | 25 U             | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | -                   | -                 | 2.1 J              | 40.0 U            | 40.0 U              | 1.03 J            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 1.8 J            | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Methylcyclohexane                        | µg/L  | n/v                           | -                   | -                 | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5.00 UJ             | 50.0 UJ           | 5 U                | 100 U             | 56.8 J <sup>B</sup> | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | 5.00 U     | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -                   | -                 | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Propylbenzene, n-                        | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Styrene                                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 5.00 UJ             | 50.0 U            | 5 U                | 100 U             | 100 U               | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | -          | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 881 J <sup>B</sup>  | 732 <sup>B</sup>  | 880 D <sup>B</sup> | 3380 <sup>B</sup> | 84.6 <sup>B</sup>   | 3.26              | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 0.76 J           | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Toluene                                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | 5 U                | 100 U             | 100 U               | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | -          | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | 5 U                | 100 U             | 100 U               | 5.00 U            | 250 U               | 250 U               | 250 U             | 25.0 U              | 25.0 U              | -          | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U      | 5.00 U               | 5.00 U            | 5.00 U            |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Trichloroethene (TCE)                    | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 112 J <sup>B</sup>  | 93.2 <sup>B</sup> | 140 <sup>B</sup>   | 283 <sup>B</sup>  | 752 <sup>B</sup>    | 35.8 <sup>B</sup> | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5.8 <sup>B</sup> | 2.85      | 2.99        | 3.05        | 3.12        | 2.00 U               | 2.00 U            | 1.29 J            |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | -                   | -                 | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Vinyl Acetate                            | µg/L  | n/v                           | 5.00 UJ             | 50.0 U            | -                  | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | -                | -         | -           | -           | -           | -                    | -                 | -                 |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 2.00 UJ             | 20.0 U            | 0.52 NJ            | 40.0 U            | 40.0 U              | 2.00 U            | 115 <sup>B</sup>    | 116 <sup>B</sup>    | 868 <sup>B</sup>  | 455 <sup>B</sup>    | 367 <sup>B</sup>    | 2.00 U     | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 4.58 NJ <sup>B</sup> | 5.43 <sup>B</sup> | 2.87 <sup>B</sup> |
| Xylene, m & p-                           | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 U            | 10 U               | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 10 U             | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Xylene, o-                               | µg/L  | 5 <sub>n</sub> <sup>B</sup>   | 2.00 UJ             | 20.0 UJ           | 5 U                | 40.0 U            | 40.0 U              | 2.00 U            | 100 U               | 100 U               | 100 U             | 10.0 U              | 10.0 U              | -          | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U      | 2.00 U               | 2.00 U            | 2.00 U            |
| Total VOC                                | µg/L  | n/v                           | 1052.8              | 888.3             | 1075.22            | 3744.9            | 1639.5              | 143.4             | 4171.6              | 4263.6              | 3598              | 1148.64             | 746.18              | 4.28       | 17.48            | 5.2       | 5.64        | 5.48        | 6.08        | 9.02                 | 6.76              | 4.16              |
| VOC Tentatively Identified Compounds     |       |                               |                     |                   |                    |                   |                     |                   |                     |                     |                   |                     |                     |            |                  |           |             |             |             |                      |                   |                   |
| Total VOC TICs                           | µg/L  | n/v                           | -                   | -                 | 5.800 J            | -                 | -                   | -                 | -                   | -                   | -                 | -                   | -                   | -          | 2.5 U            | -         | -           | -           | -           | -                    | -                 | -                 |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                      |       |                    | Off-Site Locations |             |            |             |                    |                    |                    |                    |                    |                    |                    |                    |                    |            |              |                     |                     |                     |                     |                     |                     |
|-------------------------------------------|-------|--------------------|--------------------|-------------|------------|-------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|------------|--------------|---------------------|---------------------|---------------------|---------------------|---------------------|---------------------|
| Sample Location                           |       |                    | RW-8               |             | RW-9       |             |                    |                    |                    |                    | RW-12              |                    |                    |                    |                    |            | RW-12        |                     |                     |                     |                     |                     |                     |
| Sample Date                               |       |                    | 14-Jun-12          | 20-May-13   | 8-Jun-12   | 21-May-13   | 27-Mar-14          | 29-May-14          | 1-Jul-14           | 1-Jul-14           | 7-Aug-14           | 28-Oct-14          | 28-Oct-14          | 4-Feb-15           | 4-May-15           | 8-Jun-12   | 20-May-13    | 28-May-14           | 2-Jul-14            | 7-Aug-14            | 29-Oct-14           | 4-Feb-15            | 4-May-15            |
| Sample ID                                 |       |                    | RW-8               | LI-RW-8-GW1 | RW-9       | LI-RW-9-GW1 | LI-RW-9            | LI-RW-9-PI1        | LI-RW-9-PI2        | LI-RW-DUP-PI2      | LI-RW-9-PI3        | LI-RW9-PI6         | LI-DUP-PI6         | LI-RW-9-PI9        | LI-RW-9-PI12       | RW-12      | LI-RW-12-GW1 | LI-RW-12-PI1        | LI-RW-12-PI2        | LI-RW-12-PI3        | LI-RW12-PI6         | LI-RW-12-PI9        | LI-RW-12-PI12       |
| Sampling Company                          |       |                    | DECI               | STANTEC     | DECI       | STANTEC     | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | STANTEC            | DECI       | STANTEC      | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC             | STANTEC             |
| Laboratory                                |       |                    | PARAROCH           | CCGE        | PARAROCH   | CCGE        | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH           | PARAROCH   | CCGE         | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            | PARAROCH            |
| Laboratory Work Order                     |       |                    | 12:2523            | E2301       | 12:2431    | E2314       | 141138             | 142196             | 142794             | 142794             | 143439             | 144730             | 144730             | 150382             | 151696             | 12:2431    | E2301        | 142196              | 142794              | 143439              | 144730              | 150382              | 151696              |
| Laboratory Sample ID                      |       |                    | 12:2523-01         | E2301-03    | 12:2431-01 | E2314-07    | 141138-08          | 142196-12          | 142794-04          | 142794-05          | 143439-08          | 144730-08          | 144730-13          | 150382-12          | 151696-08          | 12:2431-02 | E2301-04     | 142196-03           | 142794-14           | 143439-09           | 144730-09           | 150382-08           | 151696-09           |
| Sample Type                               | Units | TOGS               |                    |             |            |             |                    |                    |                    | Field Duplicate    |                    |                    | Field Duplicate    |                    |                    |            |              |                     |                     |                     |                     |                     |                     |
| General Chemistry                         |       |                    |                    |             |            |             |                    |                    |                    |                    |                    |                    |                    |                    |                    |            |              |                     |                     |                     |                     |                     |                     |
| Total Organic Carbon                      | µg/L  | n/v                | -                  | -           | -          | -           | 2000               | 2000               | 2500               | 2100               | 2100               | 2000 J-            | 2000 J-            | 2400 J+            | 1700               | -          | -            | 103000              | 186000              | 44800               | 5700                | 33900               | 6200                |
| Metals                                    |       |                    |                    |             |            |             |                    |                    |                    |                    |                    |                    |                    |                    |                    |            |              |                     |                     |                     |                     |                     |                     |
| Aluminum                                  | µg/L  | n/v                | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Antimony                                  | µg/L  | 3 <sup>B</sup>     | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Arsenic                                   | µg/L  | 25 <sup>B</sup>    | -                  | -           | -          | -           | 10 U               | 10 U               | 10 U               | 10 U               | 10 U               | 10 U               | 10 U               | 10 U               | 10.0 U             | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Barium                                    | µg/L  | 1000 <sup>B</sup>  | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Beryllium                                 | µg/L  | 3 <sup>A</sup>     | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Cadmium                                   | µg/L  | 5 <sup>B</sup>     | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Calcium                                   | µg/L  | n/v                | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Chromium (Total)                          | µg/L  | 50 <sup>B</sup>    | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Cobalt                                    | µg/L  | n/v                | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Copper                                    | µg/L  | 200 <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Iron                                      | µg/L  | 300. <sup>B</sup>  | -                  | -           | -          | -           | 100 U              | 91.9 J             | 129                | 91.0 J             | 86.4 J             | 100 U              | 100 U              | 68.3 J             | 76.5 J             | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Lead                                      | µg/L  | 25 <sup>B</sup>    | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Magnesium                                 | µg/L  | 35000 <sup>A</sup> | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Manganese                                 | µg/L  | 300. <sup>B</sup>  | -                  | -           | -          | -           | 15 U               | 19.8               | 98.1               | 94.4               | 220                | 153                | 161                | 284                | 214                | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Mercury                                   | µg/L  | 0.7 <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Nickel                                    | µg/L  | 100 <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Potassium                                 | µg/L  | n/v                | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Selenium                                  | µg/L  | 10 <sup>B</sup>    | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Silver                                    | µg/L  | 50 <sup>B</sup>    | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Sodium                                    | µg/L  | 20000 <sup>B</sup> | -                  | -           | -          | -           | 38100 <sup>B</sup> | 25200 <sup>B</sup> | 29000 <sup>B</sup> | 28800 <sup>B</sup> | 27700 <sup>B</sup> | 39100 <sup>B</sup> | 38600 <sup>B</sup> | 41600 <sup>B</sup> | 32000 <sup>B</sup> | -          | -            | 200000 <sup>B</sup> | 255000 <sup>B</sup> | 282000 <sup>B</sup> | 193000 <sup>B</sup> | 167000 <sup>B</sup> | 213000 <sup>B</sup> |
| Thallium                                  | µg/L  | 0.5 <sup>A</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Vanadium                                  | µg/L  | n/v                | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Zinc                                      | µg/L  | 2000 <sup>A</sup>  | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Volatile Organic Compounds                |       |                    |                    |             |            |             |                    |                    |                    |                    |                    |                    |                    |                    |                    |            |              |                     |                     |                     |                     |                     |                     |
| Acetone                                   | µg/L  | 50 <sup>A</sup>    | -                  | 25 U        | -          | 25 U        | 10.0 U             | 6.70 J             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 UJ            | 10.0 U             | -          | 25 U         | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 UJ             | 10.0 U              |
| Benzene                                   | µg/L  | 1 <sup>B</sup>     | -                  | 5 U         | -          | 5 U         | 1 U                | 1 U                | 1 U                | 1 U                | 1 U                | 1 U                | 1 U                | 1 U                | 0.700 U            | -          | 5 U          | 1 U                 | 1 U                 | 1 U                 | 1 U                 | 1 U                 | 0.700 U             |
| Bromodichloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Bromoform (Tribromomethane)               | µg/L  | 50 <sup>A</sup>    | 5.00 U             | 5 U         | 5.00 U     | 5 U         | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 UJ            | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U     | 5 U          | 5.00 U              | 5.00 U              | 5.00 UJ             | 5.00 U              | 5.00 U              | 5.00 U              |
| Bromomethane (Methyl bromide)             | µg/L  | 5.. <sup>B</sup>   | 2.00 UJ            | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 UJ            | 2.00 UJ            | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Butylbenzene, n-                          | µg/L  | 5.. <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Butylbenzene, sec- (2-Phenylbutane)       | µg/L  | 5.. <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Butylbenzene, tert-                       | µg/L  | 5.. <sup>B</sup>   | -                  | -           | -          | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -          | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Carbon Disulfide                          | µg/L  | 60 <sup>A</sup>    | -                  | 5 U         | -          | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | -          | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Carbon Tetrachloride (Tetrachloromethane) | µg/L  | 5 <sup>B</sup>     | 2.00 UJ            | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Chlorobenzene (Monochlorobenzene)         | µg/L  | 5.. <sup>B</sup>   | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Chlorobromomethane                        | µg/L  | 5.. <sup>B</sup>   | -                  | 5 U         | -          | 5 U         | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | 5.00 U             | -          | 5 U          | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              | 5.00 U              |
| Chloroethane (Ethyl Chloride)             | µg/L  | 5.. <sup>B</sup>   | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Chloroethyl Vinyl Ether, 2-               | µg/L  | n/v                | 10.0 U R           | -           | 10.0 U R   | -           | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | -                  | 10.0 U R   | -            | -                   | -                   | -                   | -                   | -                   | -                   |
| Chloroform (Trichloromethane)             | µg/L  | 7 <sup>B</sup>     | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Chloromethane                             | µg/L  | 5.. <sup>B</sup>   | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Cyclohexane                               | µg/L  | n/v                | -                  | 5 UJ        | -          | 5 U         | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | -          | 5 UJ         | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              |
| Dibromo-3-Chloropropane, 1,2- (DBCP)      | µg/L  | 0.04 <sup>B</sup>  | -                  | 5 U         | -          | 5 U         | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | 10.0 U             | -          | 5 U          | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              | 10.0 U              |
| Dibromochloromethane                      | µg/L  | 50 <sup>A</sup>    | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 UJ            | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U              |
| Dichlorobenzene, 1,2-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Dichlorobenzene, 1,3-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Dichlorobenzene, 1,4-                     | µg/L  | 3 <sup>B</sup>     | 2.00 U             | 5 U         | 2.00 U     | 5 U         | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U     | 5 U          | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |
| Dichlorodifluoromethane (Freon 12)        | µg/L  | 5.. <sup>B</sup>   | -                  | 5 U         | -          | 5 U         | 2.00 U             | 2.00 U             | 2.00 UJ            | 2.00 UJ            | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | 2.00 U             | -          | 5 U          | 2.00 U              | 2.00 UJ             | 2.00 U              | 2.00 U              | 2.00 U              | 2.00 U              |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area                                     |       |                               | Off-Site Locations |                   |                   |                  |           |             |             |                 |             |            |                 |             |              |                   |                   |                   |                   |                   |                   |                      |                   |  |  |
|------------------------------------------|-------|-------------------------------|--------------------|-------------------|-------------------|------------------|-----------|-------------|-------------|-----------------|-------------|------------|-----------------|-------------|--------------|-------------------|-------------------|-------------------|-------------------|-------------------|-------------------|----------------------|-------------------|--|--|
| Sample Location                          |       |                               | RW-8               |                   |                   |                  | RW-9      |             |             |                 | RW-12       |            |                 |             |              |                   |                   |                   |                   |                   |                   |                      |                   |  |  |
| Sample Date                              |       |                               | 14-Jun-12          | 20-May-13         | 8-Jun-12          | 21-May-13        | 27-Mar-14 | 29-May-14   | 1-Jul-14    | 1-Jul-14        | 7-Aug-14    | 28-Oct-14  | 28-Oct-14       | 4-Feb-15    | 4-May-15     | 8-Jun-12          | 20-May-13         | 28-May-14         | 2-Jul-14          | 7-Aug-14          | 29-Oct-14         | 4-Feb-15             | 4-May-15          |  |  |
| Sample ID                                |       |                               | RW-8               | LI-RW-8-GW1       | RW-9              | LI-RW-9-GW1      | LI-RW-9   | LI-RW-9-PI1 | LI-RW-9-PI2 | LI-RW-DUP-PI2   | LI-RW-9-PI3 | LI-RW9-PI6 | LI-DUP-PI6      | LI-RW-9-PI9 | LI-RW-9-PI12 | RW-12             | LI-RW-12-GW1      | LI-RW-12-PI1      | LI-RW-12-PI2      | LI-RW-12-PI3      | LI-RW12-PI6       | LI-RW-12-PI9         | LI-RW-12-PI12     |  |  |
| Sampling Company                         |       |                               | DECI               | STANTEC           | DECI              | STANTEC          | STANTEC   | STANTEC     | STANTEC     | STANTEC         | STANTEC     | STANTEC    | STANTEC         | STANTEC     | STANTEC      | DECI              | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC           | STANTEC              | STANTEC           |  |  |
| Laboratory                               |       |                               | PARAROCH           | CCGE              | PARAROCH          | CCGE             | PARAROCH  | PARAROCH    | PARAROCH    | PARAROCH        | PARAROCH    | PARAROCH   | PARAROCH        | PARAROCH    | PARAROCH     | PARAROCH          | CCGE              | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH          | PARAROCH             | PARAROCH          |  |  |
| Laboratory Work Order                    |       |                               | 12:2523            | E2301             | 12:2431           | E2314            | 141138    | 142196      | 142794      | 142794          | 143439      | 144730     | 144730          | 150382      | 151696       | 12:2431           | E2301             | 142196            | 142794            | 143439            | 144730            | 150382               | 151696            |  |  |
| Laboratory Sample ID                     |       |                               | 12:2523-01         | E2301-03          | 12:2431-01        | E2314-07         | 141138-08 | 142196-12   | 142794-04   | 142794-05       | 143439-08   | 144730-08  | 144730-13       | 150382-12   | 151696-08    | 12:2431-02        | E2301-04          | 142196-03         | 142794-14         | 143439-09         | 144730-09         | 150382-08            | 151696-09         |  |  |
| Sample Type                              | Units | TOGS                          |                    |                   |                   |                  |           |             |             | Field Duplicate |             |            | Field Duplicate |             |              |                   |                   |                   |                   |                   |                   |                      |                   |  |  |
| Volatile Organic Compounds (cont'd)      |       |                               |                    |                   |                   |                  |           |             |             |                 |             |            |                 |             |              |                   |                   |                   |                   |                   |                   |                      |                   |  |  |
| Dichloroethane, 1,1-                     | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloroethane, 1,2-                     | µg/L  | 0.6 <sup>B</sup>              | 2.00 UJ            | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloroethene, 1,1-                     | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloroethene, cis-1,2-                 | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 6.50 <sup>B</sup>  | 17.8 <sup>B</sup> | 2.00 U            | 1.2 J            | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 1.35 J     | 1.37 J          | 1.66 J      | 2.00 U       | 24.5 <sup>B</sup> | 26.5 <sup>B</sup> | 79.5 <sup>B</sup> | 118 <sup>B</sup>  | 36.2 <sup>B</sup> | 6.01 <sup>B</sup> | 20.9 <sup>B</sup>    | 5.24 <sup>B</sup> |  |  |
| Dichloroethene, trans-1,2-               | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 2.1 J             | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 1.2 J             | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloropropane, 1,2-                    | µg/L  | 1 <sup>B</sup>                | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloropropene, cis-1,3-                | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dichloropropene, trans-1,3-              | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Dioxane, 1,4-                            | µg/L  | n/v                           | -                  | 100 U R           | -                 | 100 U            | 20.0 U R  | 20.0 U R    | 20.0 U R    | 20.0 U R        | 20.0 U R    | 20.0 U     | 20.0 U          | 20.0 U      | 20.0 U       | -                 | 100 U R           | 20.0 U R          | 20.0 U R          | 20.0 U R          | 20.0 U            | 20.0 U               | 20.0 U            |  |  |
| Ethylbenzene                             | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Ethylene Dibromide (Dibromoethane, 1,2-) | µg/L  | 0.0006 <sup>B</sup>           | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Hexanone, 2- (Methyl Butyl Ketone)       | µg/L  | 50 <sup>A</sup>               | -                  | 25 U              | -                 | 25 U             | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | -                 | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Isopropylbenzene                         | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Isopropyltoluene, p- (Cymene)            | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Methyl Acetate                           | µg/L  | n/v                           | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Methyl Ethyl Ketone (MEK)                | µg/L  | 50 <sup>A</sup>               | -                  | 25 U              | -                 | 25 U             | 10.0 UJ   | 10.0 U      | 10.0 UJ     | 10.0 UJ         | 10.0 UJ     | 10.0 UJ    | 10.0 UJ         | 10.0 UJ     | 10.0 UJ      | -                 | 25 U              | 10.0 U            | 10.0 UJ           | 10.0 U            | 10.0 UJ           | 10.0 UJ              | 10.0 U            |  |  |
| Methyl Isobutyl Ketone (MIBK)            | µg/L  | n/v                           | -                  | 25 U              | -                 | 25 U             | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | -                 | 25 U              | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Methyl tert-butyl ether (MTBE)           | µg/L  | 10 <sup>A</sup>               | -                  | 3.3 J             | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 0.85 J            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Methylcyclohexane                        | µg/L  | n/v                           | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Methylene Chloride (Dichloromethane)     | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 5.00 U             | 5 U               | 5.00 U            | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | 5.00 U            | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Naphthalene                              | µg/L  | 10 <sup>B</sup>               | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Propylbenzene, n-                        | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Styrene                                  | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | -                 | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Tetrachloroethane, 1,1,2,2-              | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Tetrachloroethene (PCE)                  | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 4.3 J             | 11.3 <sup>B</sup> | 8.5 <sup>B</sup> | 3.04      | 3.58        | 4.10        | 4.11            | 3.20        | 3.28       | 3.40            | 2.67        | 2.13         | 2.71              | 4.9 J             | 5.52 <sup>B</sup> | 4.37              | 2.78              | 4.74              | 7.82 <sup>B</sup>    | 2.79              |  |  |
| Toluene                                  | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Trichlorobenzene, 1,2,3-                 | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | -                 | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Trichlorobenzene, 1,2,4-                 | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 5.00 U    | 5.00 U      | 5.00 U      | 5.00 U          | 5.00 U      | 5.00 U     | 5.00 U          | 5.00 U      | 5.00 U       | -                 | 5 U               | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U            | 5.00 U               | 5.00 U            |  |  |
| Trichloroethane, 1,1,1-                  | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Trichloroethane, 1,1,2-                  | µg/L  | 1 <sup>B</sup>                | 2.00 U             | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Trichloroethene (TCE)                    | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 7.59 <sup>B</sup>  | 20.7 <sup>B</sup> | 2.00 U            | 1.5 J            | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.45        | 2.00 U       | 6.80 <sup>B</sup> | 15 <sup>B</sup>   | 25.1 <sup>B</sup> | 29.8 <sup>B</sup> | 4.38              | 7.10 <sup>B</sup> | 14.5 <sup>B</sup>    | 4.85              |  |  |
| Trichlorofluoromethane (Freon 11)        | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | 2.00 UJ            | 5 U               | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Trichlorotrifluoroethane (Freon 113)     | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 5 U               | -                 | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | -                 | 5 U               | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U            | 2.00 U               | 2.00 U            |  |  |
| Trimethylbenzene, 1,2,4-                 | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Trimethylbenzene, 1,3,5-                 | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Vinyl Acetate                            | µg/L  | n/v                           | -                  | -                 | -                 | -                | -         | -           | -           | -               | -           | -          | -               | -           | -            | -                 | -                 | -                 | -                 | -                 | -                 | -                    | -                 |  |  |
| Vinyl chloride                           | µg/L  | 2 <sup>B</sup>                | 2.00 U             | 0.63 NJ           | 2.00 U            | 5 U              | 2.00 U    | 2.00 U      | 2.00 U      | 2.00 U          | 2.00 U      | 2.00 U     | 2.00 U          | 2.00 U      | 2.00 U       | 2.00 U            | 0.55 J            | 2.00 U            | 1.17 J            | 2.27 <sup>B</sup> | 2.00 U            | 2.28 NJ <sup>B</sup> | 2.00 U            |  |  |
| Xylene, m & p-                           | µg/L  | 5 <sup>n</sup> , <sup>B</sup> | -                  | 10 U              | -                 | 10 U             | 2.00 U    | 2.00 U      | 2.00 U      |                 |             |            |                 |             |              |                   |                   |                   |                   |                   |                   |                      |                   |  |  |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area<br>Sample Location<br>Sample Date<br><br>Sample ID<br><br>Sampling Company<br>Laboratory<br>Laboratory Work Order<br>Laboratory Sample ID<br><br>Sample Type |       |                    | Off-Site Locations<br>RW-13       |                           |                                 | QA/QC<br>Trip Blank               |                           |                           |                                 |                                 |                                 |                                 |                                 |                                 |                                 |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|--------------------|-----------------------------------|---------------------------|---------------------------------|-----------------------------------|---------------------------|---------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
|                                                                                                                                                                   |       |                    | 12-Jun-12                         | 20-May-13                 | 27-Mar-14                       | 12-Jun-12                         | 20-May-13                 | 21-May-13                 | 27-Mar-14                       | 29-May-14                       | 1-Jul-14                        | 8-Aug-14                        | 28-Oct-14                       | 3-Feb-15                        | 4-May-15                        |
|                                                                                                                                                                   |       |                    | RW-13                             | LI-RW-13-GW1              | LI-RW-13                        | Trip Blank 7346                   | Trip Blank                | Trip Blank                | Trip Blank                      | LI-Trip Blank-PI1               | LI-TRIPBLANK-PI2                | Trip Blank (T-532)              | Trip Blank (T-570)              | LI-TRIPBLANK-PI9 (T-586)        | Trip Blank (T-614)              |
|                                                                                                                                                                   |       |                    | DECI                              | STANTEC                   | STANTEC                         | DECI                              | STANTEC                   | STANTEC                   | STANTEC                         | STANTEC                         | STANTEC                         | STANTEC                         | STANTEC                         | STANTEC                         | STANTEC                         |
|                                                                                                                                                                   |       |                    | PARAROCH<br>12:2486<br>12:2486-01 | CCGE<br>E2301<br>E2301-05 | PARAROCH<br>141138<br>141138-10 | PARAROCH<br>12:2486<br>12:2486-03 | CCGE<br>E2301<br>E2301-07 | CCGE<br>E2314<br>E2314-08 | PARAROCH<br>141138<br>141138-15 | PARAROCH<br>142196<br>142196-08 | PARAROCH<br>142794<br>142794-01 | PARAROCH<br>143439<br>143439-14 | PARAROCH<br>144730<br>144730-14 | PARAROCH<br>150382<br>150382-14 | PARAROCH<br>151696<br>151696-14 |
|                                                                                                                                                                   | Units | TOGS               |                                   |                           |                                 | Trip Blank                        | Trip Blank                | Trip Blank                | Trip Blank                      | Trip Blank                      | Trip Blank                      | Trip Blank                      | Trip Blank                      | Trip Blank                      | Trip Blank                      |
| General Chemistry                                                                                                                                                 |       |                    |                                   |                           |                                 |                                   |                           |                           |                                 |                                 |                                 |                                 |                                 |                                 |                                 |
| Total Organic Carbon                                                                                                                                              | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Metals                                                                                                                                                            |       |                    |                                   |                           |                                 |                                   |                           |                           |                                 |                                 |                                 |                                 |                                 |                                 |                                 |
| Aluminum                                                                                                                                                          | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Antimony                                                                                                                                                          | µg/L  | 3 <sup>B</sup>     | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Arsenic                                                                                                                                                           | µg/L  | 25 <sup>B</sup>    | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Barium                                                                                                                                                            | µg/L  | 1000 <sup>B</sup>  | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Beryllium                                                                                                                                                         | µg/L  | 3 <sup>A</sup>     | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Cadmium                                                                                                                                                           | µg/L  | 5 <sup>B</sup>     | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Calcium                                                                                                                                                           | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Chromium (Total)                                                                                                                                                  | µg/L  | 50 <sup>B</sup>    | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Cobalt                                                                                                                                                            | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Copper                                                                                                                                                            | µg/L  | 200 <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Iron                                                                                                                                                              | µg/L  | 300. <sup>B</sup>  | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Lead                                                                                                                                                              | µg/L  | 25 <sup>B</sup>    | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Magnesium                                                                                                                                                         | µg/L  | 35000 <sup>A</sup> | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Manganese                                                                                                                                                         | µg/L  | 300. <sup>B</sup>  | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Mercury                                                                                                                                                           | µg/L  | 0.7 <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Nickel                                                                                                                                                            | µg/L  | 100 <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Potassium                                                                                                                                                         | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Selenium                                                                                                                                                          | µg/L  | 10 <sup>B</sup>    | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Silver                                                                                                                                                            | µg/L  | 50 <sup>B</sup>    | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Sodium                                                                                                                                                            | µg/L  | 20000 <sup>B</sup> | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Thallium                                                                                                                                                          | µg/L  | 0.5 <sup>A</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Vanadium                                                                                                                                                          | µg/L  | n/v                | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Zinc                                                                                                                                                              | µg/L  | 2000 <sup>A</sup>  | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Volatile Organic Compounds                                                                                                                                        |       |                    |                                   |                           |                                 |                                   |                           |                           |                                 |                                 |                                 |                                 |                                 |                                 |                                 |
| Acetone                                                                                                                                                           | µg/L  | 50 <sup>A</sup>    | -                                 | 25 U                      | 10.0 U                          | -                                 | 25 U                      | 25 U                      | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 UJ                         | 10.0 U                          |
| Benzene                                                                                                                                                           | µg/L  | 1 <sup>B</sup>     | -                                 | 5 U                       | 1 U                             | -                                 | 5 U                       | 5 U                       | 1 U                             | 1 U                             | 1 U                             | 1 U                             | 1 U                             | 1 U                             | 0.700 U                         |
| Bromodichloromethane                                                                                                                                              | µg/L  | 50 <sup>A</sup>    | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Bromoform (Tribromomethane)                                                                                                                                       | µg/L  | 50 <sup>A</sup>    | 5.00 U                            | 5 U                       | 5.00 U                          | 5.00 U                            | 5 U                       | 5 U                       | 5.00 U                          | 5.00 U                          | 5.00 U                          | 5.00 UJ                         | 5.00 U                          | 5.00 U                          | 5.00 U                          |
| Bromomethane (Methyl bromide)                                                                                                                                     | µg/L  | 5.. <sup>B</sup>   | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 UJ                         | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Butylbenzene, n-                                                                                                                                                  | µg/L  | 5.. <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Butylbenzene, sec- (2-Phenylbutane)                                                                                                                               | µg/L  | 5.. <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Butylbenzene, tert-                                                                                                                                               | µg/L  | 5.. <sup>B</sup>   | -                                 | -                         | -                               | -                                 | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Carbon Disulfide                                                                                                                                                  | µg/L  | 60 <sup>A</sup>    | -                                 | 5 U                       | 2.00 U                          | -                                 | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Carbon Tetrachloride (Tetrachloromethane)                                                                                                                         | µg/L  | 5 <sup>B</sup>     | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Chlorobenzene (Monochlorobenzene)                                                                                                                                 | µg/L  | 5.. <sup>B</sup>   | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Chlorobromomethane                                                                                                                                                | µg/L  | 5.. <sup>B</sup>   | -                                 | 5 U                       | 5.00 U                          | -                                 | 5 U                       | 5 U                       | 5.00 U                          | 5.00 U                          | 5.00 U                          | 5.00 U                          | 5.00 U                          | 5.00 U                          | 5.00 U                          |
| Chloroethane (Ethyl Chloride)                                                                                                                                     | µg/L  | 5.. <sup>B</sup>   | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Chloroethyl Vinyl Ether, 2-                                                                                                                                       | µg/L  | n/v                | 10.0 U R                          | -                         | -                               | 10.0 U R                          | -                         | -                         | -                               | -                               | -                               | -                               | -                               | -                               | -                               |
| Chloroform (Trichloromethane)                                                                                                                                     | µg/L  | 7 <sup>B</sup>     | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Chloromethane                                                                                                                                                     | µg/L  | 5.. <sup>B</sup>   | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Cyclohexane                                                                                                                                                       | µg/L  | n/v                | -                                 | 5 UJ                      | 10.0 U                          | -                                 | 5 UJ                      | 5 U                       | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          |
| Dibromo-3-Chloropropane, 1,2- (DBCP)                                                                                                                              | µg/L  | 0.04 <sup>B</sup>  | -                                 | 5 U                       | 10.0 U                          | -                                 | 5 U                       | 5 U                       | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          | 10.0 U                          |
| Dibromochloromethane                                                                                                                                              | µg/L  | 50 <sup>A</sup>    | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 UJ                         | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Dichlorobenzene, 1,2-                                                                                                                                             | µg/L  | 3 <sup>B</sup>     | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Dichlorobenzene, 1,3-                                                                                                                                             | µg/L  | 3 <sup>B</sup>     | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Dichlorobenzene, 1,4-                                                                                                                                             | µg/L  | 3 <sup>B</sup>     | 2.00 U                            | 5 U                       | 2.00 U                          | 2.00 U                            | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |
| Dichlorodifluoromethane (Freon 12)                                                                                                                                | µg/L  | 5.. <sup>B</sup>   | -                                 | 5 U                       | 2.00 U                          | -                                 | 5 U                       | 5 U                       | 2.00 U                          | 2.00 U                          | 2.00 UJ                         | 2.00 U                          | 2.00 U                          | 2.00 U                          | 2.00 U                          |

See last page for notes.

Table 2  
Summary of Analytical Results in Groundwater  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Area<br>Sample Location<br>Sample Date<br><br>Sample ID<br><br>Sampling Company<br>Laboratory<br>Laboratory Work Order<br>Laboratory Sample ID<br><br>Sample Type |       |                               | Off-Site Locations<br>RW-13               |                                      |                                            | QA/QC<br>Trip Blank                       |                                      |                                      |                                  |                                            |                                            |                                            |                                            |                                            |                                            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|-------------------------------|-------------------------------------------|--------------------------------------|--------------------------------------------|-------------------------------------------|--------------------------------------|--------------------------------------|----------------------------------|--------------------------------------------|--------------------------------------------|--------------------------------------------|--------------------------------------------|--------------------------------------------|--------------------------------------------|
|                                                                                                                                                                   |       |                               | 12-Jun-12                                 | 20-May-13                            | 27-Mar-14                                  | 12-Jun-12                                 | 20-May-13                            | 21-May-13                            | 27-Mar-14                        | 29-May-14                                  | 1-Jul-14                                   | 8-Aug-14                                   | 28-Oct-14                                  | 3-Feb-15                                   | 4-May-15                                   |
|                                                                                                                                                                   |       |                               | RW-13                                     | LI-RW-13-GW1                         | LI-RW-13                                   | Trip Blank 7346                           | Trip Blank                           | Trip Blank                           | Trip Blank                       | LI-Trip Blank-PI1                          | LI-TRIPBLANK-PI2                           | Trip Blank (T-532)                         | Trip Blank (T-570)                         | LI-TRIPBLANK-PI9 (T-586)                   | Trip Blank (T-614)                         |
|                                                                                                                                                                   |       |                               | DECI<br>PARAROCH<br>12:2486<br>12:2486-01 | STANTEC<br>CCGE<br>E2301<br>E2301-05 | STANTEC<br>PARAROCH<br>141138<br>141138-10 | DECI<br>PARAROCH<br>12:2486<br>12:2486-03 | STANTEC<br>CCGE<br>E2301<br>E2301-07 | STANTEC<br>CCGE<br>E2314<br>E2314-08 | STANTEC<br>PARAROCH<br>141138-15 | STANTEC<br>PARAROCH<br>142196<br>142196-08 | STANTEC<br>PARAROCH<br>142794<br>142794-01 | STANTEC<br>PARAROCH<br>143439<br>143439-14 | STANTEC<br>PARAROCH<br>144730<br>144730-14 | STANTEC<br>PARAROCH<br>150382<br>150382-14 | STANTEC<br>PARAROCH<br>151696<br>151696-14 |
|                                                                                                                                                                   | Units | TOGS                          |                                           |                                      |                                            | Trip Blank                                | Trip Blank                           | Trip Blank                           | Trip Blank                       | Trip Blank                                 | Trip Blank                                 | Trip Blank                                 | Trip Blank                                 | Trip Blank                                 | Trip Blank                                 |
| Volatile Organic Compounds (cont'd)                                                                                                                               |       |                               |                                           |                                      |                                            |                                           |                                      |                                      |                                  |                                            |                                            |                                            |                                            |                                            |                                            |
| Dichloroethane, 1,1-                                                                                                                                              | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloroethane, 1,2-                                                                                                                                              | µg/L  | 0.6 <sup>B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloroethene, 1,1-                                                                                                                                              | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloroethene, cis-1,2-                                                                                                                                          | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloroethene, trans-1,2-                                                                                                                                        | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloropropane, 1,2-                                                                                                                                             | µg/L  | 1 <sup>B</sup>                | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloropropene, cis-1,3-                                                                                                                                         | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dichloropropene, trans-1,3-                                                                                                                                       | µg/L  | 0.4 <sub>p</sub> <sup>B</sup> | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Dioxane, 1,4-                                                                                                                                                     | µg/L  | n/v                           | -                                         | 100 U R                              | 20.0 U R                                   | -                                         | 100 U R                              | 100 U                                | 20.0 U R                         | 20.0 U R                                   | 20.0 U R                                   | 20.0 U R                                   | 20.0 U                                     | 20.0 U                                     | 20.0 U                                     |
| Ethylbenzene                                                                                                                                                      | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Ethylene Dibromide (Dibromoethane, 1,2-)                                                                                                                          | µg/L  | 0.0006 <sup>B</sup>           | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Hexanone, 2- (Methyl Butyl Ketone)                                                                                                                                | µg/L  | 50 <sup>A</sup>               | -                                         | 25 U                                 | 5.00 U                                     | -                                         | 25 U                                 | 25 U                                 | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Isopropylbenzene                                                                                                                                                  | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Isopropyltoluene, p- (Cymene)                                                                                                                                     | µg/L  | 5 <sup>n,B</sup>              | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Methyl Acetate                                                                                                                                                    | µg/L  | n/v                           | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Methyl Ethyl Ketone (MEK)                                                                                                                                         | µg/L  | 50 <sup>A</sup>               | -                                         | 25 U                                 | 10.0 UJ                                    | -                                         | 25 U                                 | 25 U                                 | 10.0 UJ                          | 10.0 U                                     | 10.0 UJ                                    | 10.0 U                                     | 10.0 UJ                                    | 10.0 UJ                                    | 10.0 U                                     |
| Methyl Isobutyl Ketone (MIBK)                                                                                                                                     | µg/L  | n/v                           | -                                         | 25 U                                 | 5.00 U                                     | -                                         | 25 U                                 | 25 U                                 | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Methyl tert-butyl ether (MTBE)                                                                                                                                    | µg/L  | 10 <sup>A</sup>               | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Methylcyclohexane                                                                                                                                                 | µg/L  | n/v                           | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Methylene Chloride (Dichloromethane)                                                                                                                              | µg/L  | 5 <sup>n,B</sup>              | 5.00 U                                    | 5 U                                  | 5.00 U                                     | 5.00 U                                    | 5 U                                  | 3.4 J                                | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Naphthalene                                                                                                                                                       | µg/L  | 10 <sup>B</sup>               | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Propylbenzene, n-                                                                                                                                                 | µg/L  | 5 <sup>n,B</sup>              | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Styrene                                                                                                                                                           | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 5.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Tetrachloroethane, 1,1,2,2-                                                                                                                                       | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Tetrachloroethene (PCE)                                                                                                                                           | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 2.8 J                                | 2.00                                       | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Toluene                                                                                                                                                           | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trichlorobenzene, 1,2,3-                                                                                                                                          | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 5.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Trichlorobenzene, 1,2,4-                                                                                                                                          | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 5.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 5.00 U                           | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     | 5.00 U                                     |
| Trichloroethane, 1,1,1-                                                                                                                                           | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trichloroethane, 1,1,2-                                                                                                                                           | µg/L  | 1 <sup>B</sup>                | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trichloroethene (TCE)                                                                                                                                             | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 0.99 J                               | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trichlorofluoromethane (Freon 11)                                                                                                                                 | µg/L  | 5 <sup>n,B</sup>              | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trichlorotrifluoroethane (Freon 113)                                                                                                                              | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Trimethylbenzene, 1,2,4-                                                                                                                                          | µg/L  | 5 <sup>n,B</sup>              | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Trimethylbenzene, 1,3,5-                                                                                                                                          | µg/L  | 5 <sup>n,B</sup>              | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Vinyl Acetate                                                                                                                                                     | µg/L  | n/v                           | -                                         | -                                    | -                                          | -                                         | -                                    | -                                    | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |
| Vinyl chloride                                                                                                                                                    | µg/L  | 2 <sup>B</sup>                | 2.00 U                                    | 5 U                                  | 2.00 U                                     | 2.00 U                                    | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Xylene, m & p-                                                                                                                                                    | µg/L  | 5 <sup>n,B</sup>              | -                                         | 10 U                                 | 2.00 U                                     | -                                         | 10 U                                 | 10 U                                 | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Xylene, o-                                                                                                                                                        | µg/L  | 5 <sup>n,B</sup>              | -                                         | 5 U                                  | 2.00 U                                     | -                                         | 5 U                                  | 5 U                                  | 2.00 U                           | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     | 2.00 U                                     |
| Total VOC                                                                                                                                                         | µg/L  | n/v                           | ND                                        | 3.79                                 | 2.00                                       | ND                                        | ND                                   | 3.4                                  | ND                               | ND                                         | ND                                         | ND                                         | ND                                         | ND                                         | ND                                         |
| VOC Tentatively Identified Compounds                                                                                                                              |       |                               |                                           |                                      |                                            |                                           |                                      |                                      |                                  |                                            |                                            |                                            |                                            |                                            |                                            |
| Total VOC TICs                                                                                                                                                    | µg/L  | n/v                           | -                                         | 2.5 U                                | -                                          | -                                         | 2.5 U                                | 2.5 U                                | -                                | -                                          | -                                          | -                                          | -                                          | -                                          | -                                          |

See last page for notes.



Table 3  
Summary of Analytical Results in Waste Water and Discharge Permit Samples  
Former Carriage Factory  
33 Litchfield Street, Rochester, NY

| Sample Date                           |       |                  | 14-Jun-12        |                 | 2-Aug-12      |                |               | 4-Dec-12        | 22-May-13 | 24-May-13 | 7-Jun-13  | 17-Jun-13 | 5-Jul-13  | 12-Aug-13 | 4-Mar-14  | 28-Mar-14 | 10-Apr-14 | 21-Apr-14 | 2-May-14  | 12-May-14 | 27-May-14    | 10-Jul-14 | 15-Sep-14 | 4-Nov-14  | 4-Feb-15  | 17-Feb-15 | 5-May-15  |           |
|---------------------------------------|-------|------------------|------------------|-----------------|---------------|----------------|---------------|-----------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|--------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Sample Location                       |       |                  | DRILL WATER      | PUMP WATER      |               | WATER TANK     |               | Discharge Water | LI-EL-WI  | LI-WW1    | LI-EL-W2  | LI-WW2    | LI-WW3    | LI-WW4    | LI-EL-W3  | LI-EL-W4  | LI-EL-W5  | LI-EL-W6  | LI-EL-W7  | LI-EL-W8  | LI-EL-W9     | LI-EL-W10 | LI-EL-W11 | LI-EL-W12 | LI-EL-W13 | LI-EL-W14 | LI-EL-W15 | LI-EL-W16 |
| Sample ID                             |       |                  | DRILL WATER (DW) | PUMP WATER (PW) | 02 WATER TANK | 03- WATER TANK | 04-WATER TANK | S-1             | LI-EL-WI  | LI-WW1    | LI-EL-W2  | LI-WW2    | LI-WW3    | LI - WW4  | LI-EL-W3  | LI-EL-W4  | LI-EL-W5  | LI-EL-W6  | LI-EL-W7  | LI-EL-W8  | LI - EL - W9 | LI-EL-W10 | LI-EL-W11 | LI-EL-W12 | LI-EL-W13 | LI-EL-W14 | LI-EL-W15 | LI-EL-W16 |
| Sampling Company                      |       | County of        | DECI             | DECI            | DECI          | DECI           | DECI          | STANTEC         | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC      | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   | STANTEC   |           |
| Laboratory                            |       | Monroe           | PARAROCH         | PARAROCH        | PARAROCH      | PARAROCH       | PARAROCH      | PARAROCH        | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH     | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  | PARAROCH  |           |
| Laboratory Work Order                 |       | Sewer Use Permit | 12:2523          | 12:2523         | 12:3240       | 12:3240        | 12:3240       | 12:4966         | 131862    | 131903    | 132076    | 132221    | 132504    | 133068    | 140757    | 141139    | 141378    | 141521    | 141728    | 141895    | 142140       | 142903    | 144025    | 144818    | 150381    | 150502    | 151695    |           |
| Laboratory Sample ID                  |       |                  | 12:2523-02       | 12:2523-04      | 12:3240-03    | 12:3240-04     | 12:3240-05    | 12:4966-01      | 131862-01 | 131903-01 | 132076-01 | 132221-01 | 132504-01 | 133068-01 | 140757-01 | 141139-01 | 141378-01 | 141521-01 | 141728-01 | 141895-01 | 142140-01    | 142903-01 | 144025-01 | 144818-01 | 150381-01 | 150502-01 | 151695-01 |           |
| Sample Type                           | Units | Enclosure        |                  |                 |               |                |               |                 |           |           |           |           |           |           |           |           |           |           |           |           |              |           |           |           |           |           |           |           |
| Semi-Volatile Organic Compounds       |       |                  |                  |                 |               |                |               |                 |           |           |           |           |           |           |           |           |           |           |           |           |              |           |           |           |           |           |           |           |
| Acenaphthene                          | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Acenaphthylene                        | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Aniline                               | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Anthracene                            | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzidine                             | µg/L  | n/v              | -                | -               | 25.0 U        | -              | -             | 100 U           | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | 20.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzo(a)anthracene                    | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzo(a)pyrene                        | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzo(b)fluoranthene                  | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzo(g,h,i)perylene                  | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzo(k)fluoranthene                  | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzoic acid                          | µg/L  | n/v              | -                | -               | 25.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Benzyl Alcohol                        | µg/L  | n/v              | -                | -               | 25.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Bis(2-Chloroethoxy)methane            | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Bis(2-Chloroethyl)ether               | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Bis(2-Chloroisopropyl)ether           | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Bis(2-Ethylhexyl)phthalate (DEHP)     | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Bromophenyl Phenyl Ether, 4-          | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Butyl Benzyl Phthalate                | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chloro-3-methyl phenol, 4-            | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 25.0 U          | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chloroaniline, 4-                     | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chloronaphthalene, 2-                 | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chlorophenol, 2- (ortho-Chlorophenol) | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 10.0 U          | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chlorophenyl Phenyl Ether, 4-         | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Chrysene                              | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Cresol, m- (Methylphenol, 3-)         | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Cresol, o- (Methylphenol, 2-)         | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Dibenzo(a,h)anthracene                | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 250 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Dibenzofuran                          | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | -               | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Dibutyl Phthalate (DBP)               | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 647             | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Dichlorobenzene, 1,2-                 | µg/L  | n/v              | -                | -               | 10.0 U        | -              | -             | 100 U           | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | 10.0 U    | -         | -            | -         | -         | -         | -         | -         | -         |           |
| Dichlorobenzene, 1,3-                 | µg/L  | n/v              | -                | -               |               |                |               |                 |           |           |           |           |           |           |           |           |           |           |           |           |              |           |           |           |           |           |           |           |





**Stantec Consulting Services Inc.**  
61 Commercial Street, Suite 100  
Rochester NY 14614-1009  
Tel: (585) 475-1440  
Fax: (585) 272-1814

## **Appendix A**

### **Laboratory Analytical Reports**



**PARADIGM**  
ENVIRONMENTAL SERVICES, INC.

*Analytical Report For*

**Stantec**

*For Lab Project ID*

**151695**

*Referencing*

**Carriage Factory**

*Prepared*

**Wednesday, May 13, 2015**

Any noncompliant QC parameters or other notes impacting data interpretation are flagged or documented on the final report or are noted below.

A handwritten signature in black ink, reading "K. R. Hansen". The signature is written in a cursive, flowing style. The first name "K" is large and prominent, followed by "R." and "Hansen". The signature is positioned above a horizontal line.

Certifies that this report has been approved by the Technical Director or Designee

179 Lake Avenue • Rochester, NY 14608 • (585) 647-2530 • Fax (585) 647-3311 • ELAP ID# 10958



Lab Project ID: 151695

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-EL-W15

Lab Sample ID: 151695-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Metals**

| Analyte | Result    | Units | Qualifier | Date Analyzed  |
|---------|-----------|-------|-----------|----------------|
| Cadmium | < 0.00500 | mg/L  |           | 5/8/2015 14:23 |
| Copper  | < 0.0250  | mg/L  |           | 5/8/2015 14:23 |
| Lead    | < 0.0100  | mg/L  |           | 5/8/2015 14:23 |
| Zinc    | < 0.0600  | mg/L  |           | 5/8/2015 14:23 |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

**Volatile Organics (Halogenated)**

| Analyte                   | Result      | Units | Qualifier | Date Analyzed  |
|---------------------------|-------------|-------|-----------|----------------|
| 1,1,1-Trichloroethane     | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,1,2,2-Tetrachloroethane | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,1,2-Trichloroethane     | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,1-Dichloroethane        | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,1-Dichloroethene        | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,2-Dichlorobenzene       | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,2-Dichloroethane        | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,2-Dichloropropane       | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,3-Dichlorobenzene       | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| 1,4-Dichlorobenzene       | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Bromodichloromethane      | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Bromoform                 | < 5.00      | ug/L  |           | 5/9/2015 02:59 |
| Bromomethane              | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Carbon Tetrachloride      | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Chlorobenzene             | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Chloroethane              | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Chloroform                | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| Chloromethane             | < 2.00      | ug/L  |           | 5/9/2015 02:59 |
| cis-1,2-Dichloroethene    | <b>19.2</b> | ug/L  |           | 5/9/2015 02:59 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Page 2 of 9

Report Prepared Wednesday, May 13, 2015



**Lab Project ID:** 151695

**Client:** Stantec

**Project Reference:** Carriage Factory

**Sample Identifier:** LI-EL-W15

**Lab Sample ID:** 151695-01

**Date Sampled:** 5/5/2015

**Matrix:** Groundwater

**Date Received:** 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 02:59 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 02:59 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 02:59 |
| Tetrachloroethene         | <b>2.25</b> | ug/L | 5/9/2015 02:59 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 02:59 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 02:59 |
| Trichloroethene           | <b>2.62</b> | ug/L | 5/9/2015 02:59 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 02:59 |
| Vinyl chloride            | <b>10.3</b> | ug/L | 5/9/2015 02:59 |

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |
|-----------------------|-------------------------|---------------|-----------------|----------------------|
| 1,2-Dichloroethane-d4 | <b>107</b>              | 82.3 - 115    |                 | 5/9/2015 02:59       |
| 4-Bromofluorobenzene  | <b>92.8</b>             | 85.5 - 111    |                 | 5/9/2015 02:59       |
| Pentafluorobenzene    | <b>101</b>              | 91.2 - 107    |                 | 5/9/2015 02:59       |
| Toluene-D8            | <b>97.5</b>             | 90.9 - 108    |                 | 5/9/2015 02:59       |

**Method Reference(s):** EPA 8260C

EPA 5030

**Data File:** x22671.D



Lab Project ID: 151695

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-EL-W16

Lab Sample ID: 151695-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Metals**

| Analyte | Result       | Units | Qualifier | Date Analyzed  |
|---------|--------------|-------|-----------|----------------|
| Cadmium | < 0.00500    | mg/L  |           | 5/8/2015 14:28 |
| Copper  | < 0.0250     | mg/L  |           | 5/8/2015 14:28 |
| Lead    | < 0.0100     | mg/L  |           | 5/8/2015 14:28 |
| Zinc    | <b>0.177</b> | mg/L  |           | 5/8/2015 14:28 |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

**Volatile Organics (Halogenated)**

| Analyte                   | Result | Units | Qualifier | Date Analyzed  |
|---------------------------|--------|-------|-----------|----------------|
| 1,1,1-Trichloroethane     | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,1,2,2-Tetrachloroethane | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,1,2-Trichloroethane     | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,1-Dichloroethane        | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,1-Dichloroethene        | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,2-Dichlorobenzene       | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,2-Dichloroethane        | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,2-Dichloropropane       | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,3-Dichlorobenzene       | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| 1,4-Dichlorobenzene       | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Bromodichloromethane      | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Bromoform                 | < 5.00 | ug/L  |           | 5/9/2015 03:23 |
| Bromomethane              | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Carbon Tetrachloride      | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Chlorobenzene             | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Chloroethane              | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Chloroform                | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| Chloromethane             | < 2.00 | ug/L  |           | 5/9/2015 03:23 |
| cis-1,2-Dichloroethene    | < 2.00 | ug/L  |           | 5/9/2015 03:23 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Page 4 of 9

Report Prepared Wednesday, May 13, 2015



Lab Project ID: 151695

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-EL-W16

Lab Sample ID: 151695-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |        |      |                |
|---------------------------|--------|------|----------------|
| cis-1,3-Dichloropropene   | < 2.00 | ug/L | 5/9/2015 03:23 |
| Dibromochloromethane      | < 2.00 | ug/L | 5/9/2015 03:23 |
| Methylene chloride        | < 5.00 | ug/L | 5/9/2015 03:23 |
| Tetrachloroethene         | < 2.00 | ug/L | 5/9/2015 03:23 |
| trans-1,2-Dichloroethene  | < 2.00 | ug/L | 5/9/2015 03:23 |
| trans-1,3-Dichloropropene | < 2.00 | ug/L | 5/9/2015 03:23 |
| Trichloroethene           | < 2.00 | ug/L | 5/9/2015 03:23 |
| Trichlorofluoromethane    | < 2.00 | ug/L | 5/9/2015 03:23 |
| Vinyl chloride            | < 2.00 | ug/L | 5/9/2015 03:23 |

| Surrogate             | Percent Recovery | Limits     | Outliers | Date Analyzed  |
|-----------------------|------------------|------------|----------|----------------|
| 1,2-Dichloroethane-d4 | 105              | 82.3 - 115 |          | 5/9/2015 03:23 |
| 4-Bromofluorobenzene  | 91.1             | 85.5 - 111 |          | 5/9/2015 03:23 |
| Pentafluorobenzene    | 97.8             | 91.2 - 107 |          | 5/9/2015 03:23 |
| Toluene-D8            | 95.5             | 90.9 - 108 |          | 5/9/2015 03:23 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22672.D



## Analytical Report Appendix

The reported results relate only to the samples as they have been received by the laboratory.

Each page of this document is part of a multipage report. This document may not be reproduced except in its entirety, without the prior consent of Paradigm Environmental Services, Inc.

All soil/sludge samples have been reported on a dry weight basis, unless qualified "reported as received". Other solids are reported as received.

Low level Volatiles blank reports for soil/solid matrix are based on a nominal 5 gram weight. Sample results and reporting limits are based on actual weight, which may be more or less than 5 grams.

The Chain of Custody provides additional information, including compliance with sample condition requirements upon receipt. Sample condition requirements are defined under the 2003 NELAC Standard, sections 5.5.8.3.1 and 5.5.8.3.2.

NYSDOH ELAP does not certify for all parameters. Paradigm Environmental Services or the indicated subcontracted laboratory does hold certification for all analytes where certification is offered by ELAP unless otherwise specified. Aliquots separated for certain tests, such as TCLP, are indicated on the Chain of Custody and final reports with an "A" suffix.

Data qualifiers are used, when necessary, to provide additional information about the data. This information may be communicated as a flag or as text at the bottom of the report. Please refer to the following list of analyte-specific, frequently used data flags and their meaning:

*"<" = Analyzed for but not detected at or above the quantitation limit.*

*"E" = Result has been estimated, calibration limit exceeded.*

*"Z" = See case narrative.*

*"D" = Sample, Laboratory Control Sample, or Matrix Spike Duplicate results above Relative Percent Difference limit.*

*"M" = Matrix spike recoveries outside QC limits. Matrix bias indicated.*

*"B" = Method blank contained trace levels of analyte. Refer to included method blank report.*

*"J" = Result estimated between the quantitation limit and half the quantitation limit.*

*"L" = Laboratory Control Sample recovery outside accepted QC limits.*

*"P" = Concentration differs by more than 40% between the primary and secondary analytical columns.*  
*"NC" = Not calculable. Applicable to RPD if sample or duplicate result is non-detect or estimated (see primary report for data flags). Applicable to MS if sample is greater or equal to ten times the spike added. Applicable to sample surrogates or MS if sample dilution is 10x or higher.*

*"\*" = Indicates any recoveries outside associated acceptance windows. Surrogate outliers in samples are presumed matrix effects. LCS demonstrates method compliance unless otherwise noted.*

*"(1)" = Indicates data from primary column used for QC calculation.*

# GENERAL TERMS AND CONDITIONS

## LABORATORY SERVICES

These Terms and Conditions embody the whole agreement of the parties in the absence of a signed and executed contract between the Laboratory (LAB) and Client. They shall supersede all previous communications, representations, or agreements, either verbal or written, between the parties. The LAB specifically rejects all additional, inconsistent, or conflicting terms, whether printed or otherwise set forth in any purchase order or other communication from the Client to the LAB. The invalidity or unenforceability in whole or in part of any provision, term or condition hereof shall not affect in any way the validity or enforceability of the remainder of the Terms and Conditions. No waiver by LAB of any provision, term, or condition hereof or of any breach by or obligation of the Client hereunder shall constitute a waiver of such provision, term, or condition on any other occasion or a waiver of any other breach by or obligation of the Client. This agreement shall be administered and interpreted under the laws of the state which services are procured.

### Warranty.

Recognizing that the nature of many samples is unknown and that some may contain potentially hazardous components, LAB warrants only that it will perform testing services, obtain findings, and prepare reports in accordance with generally accepted analytical laboratory principles and practices at the time of performance of services. LAB makes no other warranty, express or implied.

### Scope and Compensation.

LAB agrees to perform the services described in the chain of custody to which these terms and conditions are attached. Unless the parties agree in writing to the contrary, the duties of LAB shall not be construed to exceed the services specifically described. LAB will use LAB default method for all tests unless specified otherwise on the Work Order.

Payment terms are net 30 days from the date of invoice. All overdue payments are subject to an interest charge of one and one-half percent (1-1/2%) per month or a portion thereof. Client shall also be responsible for costs of collection, including payment of reasonable attorney fees if such expense is incurred. The prices, unless stated, do not include any sale, use or other taxes. Such taxes will be added to invoice prices when required.

### Prices.

Compensation for services performed will be based on the current Lab Analytical Fee Schedule or on quotations agreed to in writing by the parties. Turnaround time based charges are determined from the time of resolution of all work order questions. Testimony, court appearances or data compilation for legal action will be charged separately. Evaluation and reporting of initial screening runs may incur additional fees.

### Limitations of Liability.

In the event of any error, omission, or other professional negligence, the sole and exclusive responsibility of LAB shall be to re-perform the deficient work at its own expense and LAB shall have no other liability whatsoever. All claims shall be deemed waived unless made in writing and received by LAB within ninety (90) days following completion of services.

LAB shall have no liability, obligation, or responsibility of any kind for losses, costs, expenses, or other damages (including but not limited to any special, direct, incidental or consequential damages) with respect to LAB's services or results.

All results provided by LAB are strictly for the use of its clients and LAB is in no way responsible for the use of such results by clients or third parties. All reports should be considered in their entirety, and LAB is not responsible for the separation, detachment, or other use of any portion of these reports. Client may not assign the lab report without the written consent of the LAB.

Client covenants and agrees, at its/his/her sole expense, to indemnify, protect, defend, and save harmless the LAB from and against any and all damages, losses, liabilities, obligations, penalties, claims, litigation, demands, defenses, judgments, suits, actions, proceedings, costs, disbursements and/or expenses (including, without limitation attorneys' and experts' fees and disbursements) of any kind whatsoever which may at any time be imposed upon, incurred by or asserted or awarded against client relating to, resulting from or arising out of (a) the breach of this agreement by this client, (b) the negligence of the client in handling, delivering or disclosing any hazardous substance, (c) the violation of the Client of any applicable law, (d) non-compliance by the Client with any environmental permit or (e) a material misrepresentation in disclosing the materials to be tested.

### Hazard Disclosure.

Client represents and warrants that any sample delivered to LAB will be preceded or accompanied by complete written disclosure of the presence of any hazardous substances known or suspected by Client. Client further warrants that any sample containing any hazardous substance that is to be delivered to LAB will be packaged, labeled, transported, and delivered properly and in accordance with applicable laws.

### Sample Handling.

Prior to LAB's acceptance of any sample (or after any revocation of acceptance), the entire risk of loss or of damage to such sample remains with Client. Samples are accepted when receipt is acknowledged on chain of custody documentation. In no event will LAB have any responsibility for the action or inaction of any carrier shipping or delivering any sample to or from LAB premises. Client authorizes LAB to proceed with the analysis of samples as received by the laboratory, recognizing that any samples not in compliance with all current DOH-ELAP-NELAP requirements for containers, preservation or holding time will be noted as such on the final report.

Disposal of hazardous waste samples is the responsibility of the Client. If the Client does not wish such samples returned, LAB may add storage and disposal fees to the final invoice. Maximum storage time for samples is 30 days after completion of analysis unless modified by applicable state or federal laws. Client will be required to give the LAB written instructions concerning disposal of these samples.

LAB reserves the absolute right, exercisable at any time, to refuse to receive delivery of, refuse to accept, or revoke acceptance of any sample, which, in the sole judgment of LAB (a) is of unsuitable volume, (b) may be or become unsuitable for or may pose a risk in handling, transport, or processing for any health, safety, environmental or other reason whether or not due to the presence in the sample of any hazardous substance, and whether or not such presence has been disclosed to LAB by Client or (c) if the condition or sample date make the sample unsuitable for analysis.

### Legal Responsibility.

LAB is solely responsible for performance of this contract, and no affiliated company, director, officer, employee, or agent shall have any legal responsibility hereunder, whether in contract or tort including negligence.

### Assignment.

LAB may assign its performance obligations under this contract to other parties, as it deems necessary. LAB shall disclose to Client any assignee (subcontractor) by ELAP ID # on the submitted final report.

### Force Majeure.

LAB shall have no responsibility or liability to the Client for any failure or delay in performance by LAB, which results in whole or in part from any cause or circumstance beyond the reasonable control of LAB. Such causes and circumstances shall include, but not limited to, acts of God, acts or orders of any government authority, strikes or other labor disputes, natural disasters, accidents, wars, civil disturbances, difficulties or delays in transportation, mail or delivery services, inability to obtain sufficient services or supplies from LAB's usual suppliers, or any other cause beyond LAB's reasonable control.

### Law.

This contract shall be continued under the laws of the State of New York without regard to its conflicts of laws provision.

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Page 7 of 9

Report Prepared Wednesday, May 13, 2015

Page 8 of 9

PARADIGM LAB  
SAMPLE  
NUMBER

1

2072



# Chain of Custody Supplement

Client: Stantec Completed by: mollyvail  
 Lab Project ID: 151695 Date: 5/5/15

## Sample Condition Requirements

Per NELAC/ELAP 210/241/242/243/244

| NELAC compliance with the sample condition requirements upon receipt |                                     |                          |                                            |
|----------------------------------------------------------------------|-------------------------------------|--------------------------|--------------------------------------------|
| Condition                                                            | Yes                                 | No                       | N/A                                        |
| Container Type                                                       | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Transferred to method-compliant container                            | <input type="checkbox"/>            | <input type="checkbox"/> | <input checked="" type="checkbox"/>        |
| Headspace (<1 mL)                                                    | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Preservation                                                         | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Chlorine Absent (<0.10 ppm per test strip)                           | <input type="checkbox"/>            | <input type="checkbox"/> | <input checked="" type="checkbox"/>        |
| Comments                                                             |                                     |                          |                                            |
| Holding Time                                                         | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Temperature                                                          | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> metals |
| Comments                                                             | 15°C iced started in field 5/5/15   |                          |                                            |
| Sufficient Sample Quantity                                           | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

---

Sample Identifier: LI-RW-1-PI12

Lab Sample ID: 151696-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 175                  | mg/L                |                         | 5/8/2015 14:36              |

Method Reference(s): EPA 6010C  
EPA 3005  
Preparation Date: 5/7/2015  
Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-1-PI12

Lab Sample ID: 151696-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 18:13 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 18:13 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 18:13 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 18:13 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 18:13 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 18:13 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 18:13 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 18:13 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-1-PI12

Lab Sample ID: 151696-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/7/2015 18:13 |
| Chloroform                | < 2.00      | ug/L | 5/7/2015 18:13 |
| Chloromethane             | < 2.00      | ug/L | 5/7/2015 18:13 |
| cis-1,2-Dichloroethene    | < 2.00      | ug/L | 5/7/2015 18:13 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/7/2015 18:13 |
| Cyclohexane               | < 10.0      | ug/L | 5/7/2015 18:13 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/7/2015 18:13 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/7/2015 18:13 |
| Ethylbenzene              | < 2.00      | ug/L | 5/7/2015 18:13 |
| Freon 113                 | < 2.00      | ug/L | 5/7/2015 18:13 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/7/2015 18:13 |
| m,p-Xylene                | < 2.00      | ug/L | 5/7/2015 18:13 |
| Methyl acetate            | < 2.00      | ug/L | 5/7/2015 18:13 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/7/2015 18:13 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/7/2015 18:13 |
| Methylene chloride        | < 5.00      | ug/L | 5/7/2015 18:13 |
| o-Xylene                  | < 2.00      | ug/L | 5/7/2015 18:13 |
| Styrene                   | < 5.00      | ug/L | 5/7/2015 18:13 |
| Tetrachloroethene         | < 2.00      | ug/L | 5/7/2015 18:13 |
| Toluene                   | < 2.00      | ug/L | 5/7/2015 18:13 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/7/2015 18:13 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/7/2015 18:13 |
| Trichloroethene           | < 2.00      | ug/L | 5/7/2015 18:13 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/7/2015 18:13 |
| Vinyl chloride            | <b>2.98</b> | ug/L | 5/7/2015 18:13 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-1-PI12

Lab Sample ID: 151696-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |
|-----------------------|-------------------------|---------------|-----------------|----------------------|
| 1,2-Dichloroethane-d4 | <b>96.4</b>             | 82.3 - 115    |                 | 5/7/2015 18:13       |
| 4-Bromofluorobenzene  | <b>89.0</b>             | 85.5 - 111    |                 | 5/7/2015 18:13       |
| Pentafluorobenzene    | <b>94.9</b>             | 91.2 - 107    |                 | 5/7/2015 18:13       |
| Toluene-D8            | <b>97.3</b>             | 90.9 - 108    |                 | 5/7/2015 18:13       |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22615.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-1-PI12

Lab Sample ID: 151696-01

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 4.5           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-2-PI12

Lab Sample ID: 151696-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | < 0.0100             | mg/L                |                         | 5/8/2015 15:58              |
| Iron                  | <b>1.85</b>          | mg/L                |                         | 5/8/2015 15:58              |
| Manganese             | <b>0.0663</b>        | mg/L                |                         | 5/8/2015 15:58              |
| Sodium                | <b>60.6</b>          | mg/L                |                         | 5/8/2015 15:58              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-2-PI12

Lab Sample ID: 151696-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 18:36 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 18:36 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 18:36 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 18:36 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 18:36 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 18:36 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 18:36 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 18:36 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-2-PI12

Lab Sample ID: 151696-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |   |                |
|---------------------------|-------------|------|---|----------------|
| Chloroethane              | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Chloroform                | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Chloromethane             | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| cis-1,2-Dichloroethene    | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Cyclohexane               | < 10.0      | ug/L |   | 5/7/2015 18:36 |
| Dibromochloromethane      | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Dichlorodifluoromethane   | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Ethylbenzene              | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Freon 113                 | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Isopropylbenzene          | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| m,p-Xylene                | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Methyl acetate            | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Methylcyclohexane         | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Methylene chloride        | < 5.00      | ug/L |   | 5/7/2015 18:36 |
| o-Xylene                  | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Styrene                   | < 5.00      | ug/L |   | 5/7/2015 18:36 |
| Tetrachloroethene         | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Toluene                   | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| trans-1,2-Dichloroethene  | <b>1.17</b> | ug/L | J | 5/7/2015 18:36 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Trichloroethene           | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Trichlorofluoromethane    | < 2.00      | ug/L |   | 5/7/2015 18:36 |
| Vinyl chloride            | <b>1.26</b> | ug/L | J | 5/7/2015 18:36 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-2-PI12

Lab Sample ID: 151696-02

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>96.0</b>             | 82.3 - 115    |                 | 5/7/2015             | 18:36 |
| 4-Bromofluorobenzene  | <b>89.5</b>             | 85.5 - 111    |                 | 5/7/2015             | 18:36 |
| Pentafluorobenzene    | <b>95.2</b>             | 91.2 - 107    |                 | 5/7/2015             | 18:36 |
| Toluene-D8            | <b>96.5</b>             | 90.9 - 108    |                 | 5/7/2015             | 18:36 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22616.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



**Lab Project ID:** 151696

**Client:** **Stantec**

**Project Reference:** Carriage Factory

---

---

**Sample Identifier:** LI-RW-2-PI12

**Lab Sample ID:** 151696-02

**Date Sampled:** 5/5/2015

**Matrix:** Groundwater

**Date Received:** 5/5/2015

---

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 2.7           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

*Report Prepared Friday, May 15, 2015*



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-3-PI12

Lab Sample ID: 151696-03

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 85.3                 | mg/L                |                         | 5/8/2015 16:02              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-3-PI12

Lab Sample ID: 151696-03

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 19:00 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 19:00 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 19:00 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 19:00 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 19:00 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 19:00 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 19:00 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 19:00 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-3-PI12

Lab Sample ID: 151696-03

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |   |                |
|---------------------------|-------------|------|---|----------------|
| Chloroethane              | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Chloroform                | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Chloromethane             | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| cis-1,2-Dichloroethene    | <b>1.23</b> | ug/L | J | 5/7/2015 19:00 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Cyclohexane               | < 10.0      | ug/L |   | 5/7/2015 19:00 |
| Dibromochloromethane      | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Dichlorodifluoromethane   | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Ethylbenzene              | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Freon 113                 | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Isopropylbenzene          | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| m,p-Xylene                | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Methyl acetate            | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Methyl tert-butyl Ether   | <b>2.12</b> | ug/L |   | 5/7/2015 19:00 |
| Methylcyclohexane         | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Methylene chloride        | < 5.00      | ug/L |   | 5/7/2015 19:00 |
| o-Xylene                  | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Styrene                   | < 5.00      | ug/L |   | 5/7/2015 19:00 |
| Tetrachloroethene         | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Toluene                   | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| trans-1,2-Dichloroethene  | <b>3.63</b> | ug/L |   | 5/7/2015 19:00 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Trichloroethene           | <b>1.04</b> | ug/L | J | 5/7/2015 19:00 |
| Trichlorofluoromethane    | < 2.00      | ug/L |   | 5/7/2015 19:00 |
| Vinyl chloride            | <b>4.41</b> | ug/L |   | 5/7/2015 19:00 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-3-PI12

Lab Sample ID: 151696-03

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | 97.1                    | 82.3 - 115    |                 | 5/7/2015             | 19:00 |
| 4-Bromofluorobenzene  | 86.8                    | 85.5 - 111    |                 | 5/7/2015             | 19:00 |
| Pentafluorobenzene    | 97.5                    | 91.2 - 107    |                 | 5/7/2015             | 19:00 |
| Toluene-D8            | 96.4                    | 90.9 - 108    |                 | 5/7/2015             | 19:00 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22617.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-3-PI12

Lab Sample ID: 151696-03

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Total Organic Carbon  | 6.1                  | mg/L                |                         | 5/12/2015                   |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-4-PI12

Lab Sample ID: 151696-04

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 86.9                 | mg/L                |                         | 5/8/2015 16:07              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-4-PI12

Lab Sample ID: 151696-04

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 19:23 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 19:23 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 19:23 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 19:23 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 19:23 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 19:23 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 19:23 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 19:23 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-4-PI12

Lab Sample ID: 151696-04

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/7/2015 19:23 |
| Chloroform                | < 2.00      | ug/L | 5/7/2015 19:23 |
| Chloromethane             | < 2.00      | ug/L | 5/7/2015 19:23 |
| cis-1,2-Dichloroethene    | <b>14.8</b> | ug/L | 5/7/2015 19:23 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/7/2015 19:23 |
| Cyclohexane               | < 10.0      | ug/L | 5/7/2015 19:23 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/7/2015 19:23 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/7/2015 19:23 |
| Ethylbenzene              | < 2.00      | ug/L | 5/7/2015 19:23 |
| Freon 113                 | < 2.00      | ug/L | 5/7/2015 19:23 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/7/2015 19:23 |
| m,p-Xylene                | < 2.00      | ug/L | 5/7/2015 19:23 |
| Methyl acetate            | < 2.00      | ug/L | 5/7/2015 19:23 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/7/2015 19:23 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/7/2015 19:23 |
| Methylene chloride        | < 5.00      | ug/L | 5/7/2015 19:23 |
| o-Xylene                  | < 2.00      | ug/L | 5/7/2015 19:23 |
| Styrene                   | < 5.00      | ug/L | 5/7/2015 19:23 |
| Tetrachloroethene         | <b>18.7</b> | ug/L | 5/7/2015 19:23 |
| Toluene                   | < 2.00      | ug/L | 5/7/2015 19:23 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/7/2015 19:23 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/7/2015 19:23 |
| Trichloroethene           | <b>8.94</b> | ug/L | 5/7/2015 19:23 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/7/2015 19:23 |
| Vinyl chloride            | <b>5.15</b> | ug/L | 5/7/2015 19:23 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-4-PI12

Lab Sample ID: 151696-04

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>96.9</b>             | 82.3 - 115    |                 | 5/7/2015             | 19:23 |
| 4-Bromofluorobenzene  | <b>88.8</b>             | 85.5 - 111    |                 | 5/7/2015             | 19:23 |
| Pentafluorobenzene    | <b>95.7</b>             | 91.2 - 107    |                 | 5/7/2015             | 19:23 |
| Toluene-D8            | <b>97.7</b>             | 90.9 - 108    |                 | 5/7/2015             | 19:23 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22618.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-4-PI12

Lab Sample ID: 151696-04

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 5.4           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

---

Sample Identifier: LI-RW-5-PI12

Lab Sample ID: 151696-05

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | < 0.0100             | mg/L                |                         | 5/8/2015 16:11              |
| Iron                  | <b>9.91</b>          | mg/L                |                         | 5/8/2015 16:11              |
| Manganese             | <b>0.0298</b>        | mg/L                |                         | 5/8/2015 16:11              |
| Sodium                | <b>44.6</b>          | mg/L                |                         | 5/8/2015 16:11              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-5-PI12

Lab Sample ID: 151696-05

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 19:47 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 19:47 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 19:47 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 19:47 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 19:47 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 19:47 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 19:47 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 19:47 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-5-PI12

Lab Sample ID: 151696-05

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                  |
|---------------------------|-------------|------|------------------|
| Chloroethane              | < 2.00      | ug/L | 5/7/2015 19:47   |
| Chloroform                | < 2.00      | ug/L | 5/7/2015 19:47   |
| Chloromethane             | < 2.00      | ug/L | 5/7/2015 19:47   |
| cis-1,2-Dichloroethene    | <b>4.52</b> | ug/L | 5/7/2015 19:47   |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/7/2015 19:47   |
| Cyclohexane               | < 10.0      | ug/L | 5/7/2015 19:47   |
| Dibromochloromethane      | < 2.00      | ug/L | 5/7/2015 19:47   |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/7/2015 19:47   |
| Ethylbenzene              | < 2.00      | ug/L | 5/7/2015 19:47   |
| Freon 113                 | < 2.00      | ug/L | 5/7/2015 19:47   |
| Isopropylbenzene          | < 2.00      | ug/L | 5/7/2015 19:47   |
| m,p-Xylene                | < 2.00      | ug/L | 5/7/2015 19:47   |
| Methyl acetate            | < 2.00      | ug/L | 5/7/2015 19:47   |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/7/2015 19:47   |
| Methylcyclohexane         | < 2.00      | ug/L | 5/7/2015 19:47   |
| Methylene chloride        | < 5.00      | ug/L | 5/7/2015 19:47   |
| o-Xylene                  | < 2.00      | ug/L | 5/7/2015 19:47   |
| Styrene                   | < 5.00      | ug/L | 5/7/2015 19:47   |
| Tetrachloroethene         | < 2.00      | ug/L | 5/7/2015 19:47   |
| Toluene                   | < 2.00      | ug/L | 5/7/2015 19:47   |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/7/2015 19:47   |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/7/2015 19:47   |
| Trichloroethene           | < 2.00      | ug/L | 5/7/2015 19:47   |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/7/2015 19:47   |
| Vinyl chloride            | <b>1.61</b> | ug/L | J 5/7/2015 19:47 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-5-PI12

Lab Sample ID: 151696-05

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>98.2</b>             | 82.3 - 115    |                 | 5/7/2015             | 19:47 |
| 4-Bromofluorobenzene  | <b>89.3</b>             | 85.5 - 111    |                 | 5/7/2015             | 19:47 |
| Pentafluorobenzene    | <b>93.3</b>             | 91.2 - 107    |                 | 5/7/2015             | 19:47 |
| Toluene-D8            | <b>96.5</b>             | 90.9 - 108    |                 | 5/7/2015             | 19:47 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22619.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-5-PI12

Lab Sample ID: 151696-05

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 2.2           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-6-PI12

Lab Sample ID: 151696-06

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | < 0.0100             | mg/L                |                         | 5/8/2015 16:15              |
| Iron                  | <b>0.864</b>         | mg/L                |                         | 5/8/2015 16:15              |
| Manganese             | <b>0.0309</b>        | mg/L                |                         | 5/8/2015 16:15              |
| Sodium                | <b>68.8</b>          | mg/L                |                         | 5/8/2015 16:15              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-6-PI12

Lab Sample ID: 151696-06

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result | Units | Qualifier | Date Analyzed  |
|-----------------------------|--------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,1,2,2-Tetrachloroethane   | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,1,2-Trichloroethane       | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,1-Dichloroethane          | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,1-Dichloroethene          | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2,3-Trichlorobenzene      | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2,4-Trichlorobenzene      | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2-Dibromo-3-Chloropropane | < 50.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2-Dibromoethane           | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2-Dichlorobenzene         | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2-Dichloroethane          | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,2-Dichloropropane         | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,3-Dichlorobenzene         | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,4-Dichlorobenzene         | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| 1,4-dioxane                 | < 100  | ug/L  |           | 5/7/2015 20:32 |
| 2-Butanone                  | < 50.0 | ug/L  |           | 5/7/2015 20:32 |
| 2-Hexanone                  | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| 4-Methyl-2-pentanone        | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| Acetone                     | < 50.0 | ug/L  |           | 5/7/2015 20:32 |
| Benzene                     | < 3.50 | ug/L  |           | 5/7/2015 20:32 |
| Bromochloromethane          | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| Bromodichloromethane        | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| Bromoform                   | < 25.0 | ug/L  |           | 5/7/2015 20:32 |
| Bromomethane                | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| Carbon disulfide            | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| Carbon Tetrachloride        | < 10.0 | ug/L  |           | 5/7/2015 20:32 |
| Chlorobenzene               | < 10.0 | ug/L  |           | 5/7/2015 20:32 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-6-PI12

Lab Sample ID: 151696-06

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                  |
|---------------------------|-------------|------|------------------|
| Chloroethane              | < 10.0      | ug/L | 5/7/2015 20:32   |
| Chloroform                | < 10.0      | ug/L | 5/7/2015 20:32   |
| Chloromethane             | < 10.0      | ug/L | 5/7/2015 20:32   |
| cis-1,2-Dichloroethene    | <b>373</b>  | ug/L | 5/7/2015 20:32   |
| cis-1,3-Dichloropropene   | < 10.0      | ug/L | 5/7/2015 20:32   |
| Cyclohexane               | < 50.0      | ug/L | 5/7/2015 20:32   |
| Dibromochloromethane      | < 10.0      | ug/L | 5/7/2015 20:32   |
| Dichlorodifluoromethane   | < 10.0      | ug/L | 5/7/2015 20:32   |
| Ethylbenzene              | < 10.0      | ug/L | 5/7/2015 20:32   |
| Freon 113                 | < 10.0      | ug/L | 5/7/2015 20:32   |
| Isopropylbenzene          | < 10.0      | ug/L | 5/7/2015 20:32   |
| m,p-Xylene                | < 10.0      | ug/L | 5/7/2015 20:32   |
| Methyl acetate            | < 10.0      | ug/L | 5/7/2015 20:32   |
| Methyl tert-butyl Ether   | < 10.0      | ug/L | 5/7/2015 20:32   |
| Methylcyclohexane         | < 10.0      | ug/L | 5/7/2015 20:32   |
| Methylene chloride        | < 25.0      | ug/L | 5/7/2015 20:32   |
| o-Xylene                  | < 10.0      | ug/L | 5/7/2015 20:32   |
| Styrene                   | < 25.0      | ug/L | 5/7/2015 20:32   |
| Tetrachloroethene         | < 10.0      | ug/L | 5/7/2015 20:32   |
| Toluene                   | < 10.0      | ug/L | 5/7/2015 20:32   |
| trans-1,2-Dichloroethene  | <b>6.18</b> | ug/L | J 5/7/2015 20:32 |
| trans-1,3-Dichloropropene | < 10.0      | ug/L | 5/7/2015 20:32   |
| Trichloroethene           | < 10.0      | ug/L | 5/7/2015 20:32   |
| Trichlorofluoromethane    | < 10.0      | ug/L | 5/7/2015 20:32   |
| Vinyl chloride            | <b>367</b>  | ug/L | 5/7/2015 20:32   |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-6-PI12

Lab Sample ID: 151696-06

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>98.7</b>             | 82.3 - 115    |                 | 5/7/2015             | 20:32 |
| 4-Bromofluorobenzene  | <b>87.5</b>             | 85.5 - 111    |                 | 5/7/2015             | 20:32 |
| Pentafluorobenzene    | <b>95.0</b>             | 91.2 - 107    |                 | 5/7/2015             | 20:32 |
| Toluene-D8            | <b>94.5</b>             | 90.9 - 108    |                 | 5/7/2015             | 20:32 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22621.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-6-PI12

Lab Sample ID: 151696-06

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 3.0           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-7-PI12

Lab Sample ID: 151696-07

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 49.9                 | mg/L                |                         | 5/8/2015 16:20              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-7-PI12

Lab Sample ID: 151696-07

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 20:09 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 20:09 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 20:09 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 20:09 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 20:09 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 20:09 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 20:09 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 20:09 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-7-PI12

Lab Sample ID: 151696-07

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                  |
|---------------------------|-------------|------|------------------|
| Chloroethane              | < 2.00      | ug/L | 5/7/2015 20:09   |
| Chloroform                | < 2.00      | ug/L | 5/7/2015 20:09   |
| Chloromethane             | < 2.00      | ug/L | 5/7/2015 20:09   |
| cis-1,2-Dichloroethene    | < 2.00      | ug/L | 5/7/2015 20:09   |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/7/2015 20:09   |
| Cyclohexane               | < 10.0      | ug/L | 5/7/2015 20:09   |
| Dibromochloromethane      | < 2.00      | ug/L | 5/7/2015 20:09   |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/7/2015 20:09   |
| Ethylbenzene              | < 2.00      | ug/L | 5/7/2015 20:09   |
| Freon 113                 | < 2.00      | ug/L | 5/7/2015 20:09   |
| Isopropylbenzene          | < 2.00      | ug/L | 5/7/2015 20:09   |
| m,p-Xylene                | < 2.00      | ug/L | 5/7/2015 20:09   |
| Methyl acetate            | < 2.00      | ug/L | 5/7/2015 20:09   |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/7/2015 20:09   |
| Methylcyclohexane         | < 2.00      | ug/L | 5/7/2015 20:09   |
| Methylene chloride        | < 5.00      | ug/L | 5/7/2015 20:09   |
| o-Xylene                  | < 2.00      | ug/L | 5/7/2015 20:09   |
| Styrene                   | < 5.00      | ug/L | 5/7/2015 20:09   |
| Tetrachloroethene         | < 2.00      | ug/L | 5/7/2015 20:09   |
| Toluene                   | < 2.00      | ug/L | 5/7/2015 20:09   |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/7/2015 20:09   |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/7/2015 20:09   |
| Trichloroethene           | <b>1.29</b> | ug/L | J 5/7/2015 20:09 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/7/2015 20:09   |
| Vinyl chloride            | <b>2.87</b> | ug/L | 5/7/2015 20:09   |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-7-PI12

Lab Sample ID: 151696-07

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>98.8</b>             | 82.3 - 115    |                 | 5/7/2015             | 20:09 |
| 4-Bromofluorobenzene  | <b>87.9</b>             | 85.5 - 111    |                 | 5/7/2015             | 20:09 |
| Pentafluorobenzene    | <b>92.8</b>             | 91.2 - 107    |                 | 5/7/2015             | 20:09 |
| Toluene-D8            | <b>96.0</b>             | 90.9 - 108    |                 | 5/7/2015             | 20:09 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22620.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-7-PI12

Lab Sample ID: 151696-07

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 3.1           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-9-PI12

Lab Sample ID: 151696-08

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | < 0.0100             | mg/L                |                         | 5/8/2015 16:24              |
| Iron                  | <b>0.0765</b>        | mg/L                | J                       | 5/8/2015 16:24              |
| Manganese             | <b>0.214</b>         | mg/L                |                         | 5/8/2015 16:24              |
| Sodium                | <b>32.0</b>          | mg/L                |                         | 5/8/2015 16:24              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-9-PI12

Lab Sample ID: 151696-08

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  | M         | 5/9/2015 01:49 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/9/2015 01:49 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/9/2015 01:49 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/9/2015 01:49 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/9/2015 01:49 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/9/2015 01:49 |
| Acetone                     | < 10.0  | ug/L  |           | 5/9/2015 01:49 |
| Benzene                     | < 0.700 | ug/L  |           | 5/9/2015 01:49 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/9/2015 01:49 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| Bromoform                   | < 5.00  | ug/L  | M         | 5/9/2015 01:49 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/9/2015 01:49 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/9/2015 01:49 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-9-PI12

Lab Sample ID: 151696-08

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/9/2015 01:49 |
| Chloroform                | < 2.00      | ug/L | 5/9/2015 01:49 |
| Chloromethane             | < 2.00      | ug/L | 5/9/2015 01:49 |
| cis-1,2-Dichloroethene    | < 2.00      | ug/L | 5/9/2015 01:49 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 01:49 |
| Cyclohexane               | < 10.0      | ug/L | 5/9/2015 01:49 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 01:49 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/9/2015 01:49 |
| Ethylbenzene              | < 2.00      | ug/L | 5/9/2015 01:49 |
| Freon 113                 | < 2.00      | ug/L | 5/9/2015 01:49 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/9/2015 01:49 |
| m,p-Xylene                | < 2.00      | ug/L | 5/9/2015 01:49 |
| Methyl acetate            | < 2.00      | ug/L | 5/9/2015 01:49 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/9/2015 01:49 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/9/2015 01:49 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 01:49 |
| o-Xylene                  | < 2.00      | ug/L | 5/9/2015 01:49 |
| Styrene                   | < 5.00      | ug/L | 5/9/2015 01:49 |
| Tetrachloroethene         | <b>2.13</b> | ug/L | 5/9/2015 01:49 |
| Toluene                   | < 2.00      | ug/L | 5/9/2015 01:49 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 01:49 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 01:49 |
| Trichloroethene           | < 2.00      | ug/L | 5/9/2015 01:49 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 01:49 |
| Vinyl chloride            | < 2.00      | ug/L | 5/9/2015 01:49 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-9-PI12

Lab Sample ID: 151696-08

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>106</b>              | 82.3 - 115    |                 | 5/9/2015             | 01:49 |
| 4-Bromofluorobenzene  | <b>87.8</b>             | 85.5 - 111    |                 | 5/9/2015             | 01:49 |
| Pentafluorobenzene    | <b>96.8</b>             | 91.2 - 107    |                 | 5/9/2015             | 01:49 |
| Toluene-D8            | <b>94.0</b>             | 90.9 - 108    |                 | 5/9/2015             | 01:49 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22668.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-9-PI12

Lab Sample ID: 151696-08

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 1.7           | mg/L         |                  | 5/12/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-12-PI12

Lab Sample ID: 151696-09

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 213                  | mg/L                |                         | 5/8/2015 16:45              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-12-PI12

Lab Sample ID: 151696-09

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/9/2015 01:26 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/9/2015 01:26 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/9/2015 01:26 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| Acetone                     | < 10.0  | ug/L  |           | 5/9/2015 01:26 |
| Benzene                     | < 0.700 | ug/L  |           | 5/9/2015 01:26 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/9/2015 01:26 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/9/2015 01:26 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/9/2015 01:26 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-12-PI12

Lab Sample ID: 151696-09

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/9/2015 01:26 |
| Chloroform                | < 2.00      | ug/L | 5/9/2015 01:26 |
| Chloromethane             | < 2.00      | ug/L | 5/9/2015 01:26 |
| cis-1,2-Dichloroethene    | <b>5.24</b> | ug/L | 5/9/2015 01:26 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 01:26 |
| Cyclohexane               | < 10.0      | ug/L | 5/9/2015 01:26 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 01:26 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/9/2015 01:26 |
| Ethylbenzene              | < 2.00      | ug/L | 5/9/2015 01:26 |
| Freon 113                 | < 2.00      | ug/L | 5/9/2015 01:26 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/9/2015 01:26 |
| m,p-Xylene                | < 2.00      | ug/L | 5/9/2015 01:26 |
| Methyl acetate            | < 2.00      | ug/L | 5/9/2015 01:26 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/9/2015 01:26 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/9/2015 01:26 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 01:26 |
| o-Xylene                  | < 2.00      | ug/L | 5/9/2015 01:26 |
| Styrene                   | < 5.00      | ug/L | 5/9/2015 01:26 |
| Tetrachloroethene         | <b>2.79</b> | ug/L | 5/9/2015 01:26 |
| Toluene                   | < 2.00      | ug/L | 5/9/2015 01:26 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 01:26 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 01:26 |
| Trichloroethene           | <b>4.85</b> | ug/L | 5/9/2015 01:26 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 01:26 |
| Vinyl chloride            | < 2.00      | ug/L | 5/9/2015 01:26 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-RW-12-PI12

Lab Sample ID: 151696-09

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>106</b>              | 82.3 - 115    |                 | 5/9/2015             | 01:26 |
| 4-Bromofluorobenzene  | <b>89.6</b>             | 85.5 - 111    |                 | 5/9/2015             | 01:26 |
| Pentafluorobenzene    | <b>101</b>              | 91.2 - 107    |                 | 5/9/2015             | 01:26 |
| Toluene-D8            | <b>98.3</b>             | 90.9 - 108    |                 | 5/9/2015             | 01:26 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22667.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-RW-12-PI12

Lab Sample ID: 151696-09

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Total Organic Carbon  | 6.2                  | mg/L                |                         | 5/13/2015                   |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

---

Sample Identifier: LI-DUP-PI12

Lab Sample ID: 151696-10

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | <b>0.0104</b>        | mg/L                |                         | 5/8/2015 16:49              |
| Iron                  | <b>10.1</b>          | mg/L                |                         | 5/8/2015 16:49              |
| Manganese             | <b>0.949</b>         | mg/L                |                         | 5/8/2015 16:49              |
| Sodium                | <b>50.3</b>          | mg/L                |                         | 5/8/2015 16:49              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-DUP-PI12

Lab Sample ID: 151696-10

Matrix: Groundwater

Date Sampled: 5/4/2015

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/9/2015 01:03 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/9/2015 01:03 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/9/2015 01:03 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| Acetone                     | < 10.0  | ug/L  |           | 5/9/2015 01:03 |
| Benzene                     | < 0.700 | ug/L  |           | 5/9/2015 01:03 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/9/2015 01:03 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/9/2015 01:03 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/9/2015 01:03 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-DUP-PI12

Lab Sample ID: 151696-10

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/9/2015 01:03 |
| Chloroform                | < 2.00      | ug/L | 5/9/2015 01:03 |
| Chloromethane             | < 2.00      | ug/L | 5/9/2015 01:03 |
| cis-1,2-Dichloroethene    | <b>4.11</b> | ug/L | 5/9/2015 01:03 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 01:03 |
| Cyclohexane               | < 10.0      | ug/L | 5/9/2015 01:03 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 01:03 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/9/2015 01:03 |
| Ethylbenzene              | < 2.00      | ug/L | 5/9/2015 01:03 |
| Freon 113                 | < 2.00      | ug/L | 5/9/2015 01:03 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/9/2015 01:03 |
| m,p-Xylene                | < 2.00      | ug/L | 5/9/2015 01:03 |
| Methyl acetate            | < 2.00      | ug/L | 5/9/2015 01:03 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/9/2015 01:03 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/9/2015 01:03 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 01:03 |
| o-Xylene                  | < 2.00      | ug/L | 5/9/2015 01:03 |
| Styrene                   | < 5.00      | ug/L | 5/9/2015 01:03 |
| Tetrachloroethene         | < 2.00      | ug/L | 5/9/2015 01:03 |
| Toluene                   | < 2.00      | ug/L | 5/9/2015 01:03 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 01:03 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 01:03 |
| Trichloroethene           | <b>2.42</b> | ug/L | 5/9/2015 01:03 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 01:03 |
| Vinyl chloride            | <b>11.3</b> | ug/L | 5/9/2015 01:03 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-DUP-PI12

Lab Sample ID: 151696-10

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>105</b>              | 82.3 - 115    |                 | 5/9/2015             | 01:03 |
| 4-Bromofluorobenzene  | <b>92.3</b>             | 85.5 - 111    |                 | 5/9/2015             | 01:03 |
| Pentafluorobenzene    | <b>99.0</b>             | 91.2 - 107    |                 | 5/9/2015             | 01:03 |
| Toluene-D8            | <b>98.0</b>             | 90.9 - 108    |                 | 5/9/2015             | 01:03 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22666.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-DUP-PI12

Lab Sample ID: 151696-10

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 5.3           | mg/L         |                  | 5/13/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-B102-MW-PI12

Lab Sample ID: 151696-11

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | 0.00792              | mg/L                | J                       | 5/8/2015 16:53              |
| Iron                  | 10.0                 | mg/L                |                         | 5/8/2015 16:53              |
| Manganese             | 0.945                | mg/L                |                         | 5/8/2015 16:53              |
| Sodium                | 49.8                 | mg/L                |                         | 5/8/2015 16:53              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B102-MW-PI12

Lab Sample ID: 151696-11

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/9/2015 00:39 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/9/2015 00:39 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/9/2015 00:39 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| Acetone                     | < 10.0  | ug/L  |           | 5/9/2015 00:39 |
| Benzene                     | < 0.700 | ug/L  |           | 5/9/2015 00:39 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/9/2015 00:39 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/9/2015 00:39 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/9/2015 00:39 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B102-MW-PI12

Lab Sample ID: 151696-11

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/9/2015 00:39 |
| Chloroform                | < 2.00      | ug/L | 5/9/2015 00:39 |
| Chloromethane             | < 2.00      | ug/L | 5/9/2015 00:39 |
| cis-1,2-Dichloroethene    | <b>4.10</b> | ug/L | 5/9/2015 00:39 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 00:39 |
| Cyclohexane               | < 10.0      | ug/L | 5/9/2015 00:39 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 00:39 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/9/2015 00:39 |
| Ethylbenzene              | < 2.00      | ug/L | 5/9/2015 00:39 |
| Freon 113                 | < 2.00      | ug/L | 5/9/2015 00:39 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/9/2015 00:39 |
| m,p-Xylene                | < 2.00      | ug/L | 5/9/2015 00:39 |
| Methyl acetate            | < 2.00      | ug/L | 5/9/2015 00:39 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/9/2015 00:39 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/9/2015 00:39 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 00:39 |
| o-Xylene                  | < 2.00      | ug/L | 5/9/2015 00:39 |
| Styrene                   | < 5.00      | ug/L | 5/9/2015 00:39 |
| Tetrachloroethene         | < 2.00      | ug/L | 5/9/2015 00:39 |
| Toluene                   | < 2.00      | ug/L | 5/9/2015 00:39 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 00:39 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 00:39 |
| Trichloroethene           | <b>2.38</b> | ug/L | 5/9/2015 00:39 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 00:39 |
| Vinyl chloride            | <b>11.0</b> | ug/L | 5/9/2015 00:39 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B102-MW-PI12

Lab Sample ID: 151696-11

Date Sampled: 5/4/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>106</b>              | 82.3 - 115    |                 | 5/9/2015             | 00:39 |
| 4-Bromofluorobenzene  | <b>89.6</b>             | 85.5 - 111    |                 | 5/9/2015             | 00:39 |
| Pentafluorobenzene    | <b>100</b>              | 91.2 - 107    |                 | 5/9/2015             | 00:39 |
| Toluene-D8            | <b>96.2</b>             | 90.9 - 108    |                 | 5/9/2015             | 00:39 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22665.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



**PARADIGM**  
ENVIRONMENTAL SERVICES, INC.

**Lab Project ID:** 151696

**Client:** **Stantec**

**Project Reference:** Carriage Factory

---

**Sample Identifier:** LI-B102-MW-PI12

**Lab Sample ID:** 151696-11

**Date Sampled:** 5/4/2015

**Matrix:** Groundwater

**Date Received:** 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | <b>5.4</b>    | mg/L         |                  | 5/13/2015            |

**Method Reference(s):** SM 5310 C  
**Subcontractor ELAP ID:** 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

*Report Prepared Friday, May 15, 2015*



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-B106-MW-PI12

Lab Sample ID: 151696-12

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Sodium                | 40.1                 | mg/L                |                         | 5/8/2015 16:58              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B106-MW-PI12

Lab Sample ID: 151696-12

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/9/2015 00:16 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/9/2015 00:16 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/9/2015 00:16 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| Acetone                     | < 10.0  | ug/L  |           | 5/9/2015 00:16 |
| Benzene                     | < 0.700 | ug/L  |           | 5/9/2015 00:16 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/9/2015 00:16 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/9/2015 00:16 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/9/2015 00:16 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B106-MW-PI12

Lab Sample ID: 151696-12

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |                |
|---------------------------|-------------|------|----------------|
| Chloroethane              | < 2.00      | ug/L | 5/9/2015 00:16 |
| Chloroform                | < 2.00      | ug/L | 5/9/2015 00:16 |
| Chloromethane             | < 2.00      | ug/L | 5/9/2015 00:16 |
| cis-1,2-Dichloroethene    | <b>10.2</b> | ug/L | 5/9/2015 00:16 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L | 5/9/2015 00:16 |
| Cyclohexane               | < 10.0      | ug/L | 5/9/2015 00:16 |
| Dibromochloromethane      | < 2.00      | ug/L | 5/9/2015 00:16 |
| Dichlorodifluoromethane   | < 2.00      | ug/L | 5/9/2015 00:16 |
| Ethylbenzene              | < 2.00      | ug/L | 5/9/2015 00:16 |
| Freon 113                 | < 2.00      | ug/L | 5/9/2015 00:16 |
| Isopropylbenzene          | < 2.00      | ug/L | 5/9/2015 00:16 |
| m,p-Xylene                | < 2.00      | ug/L | 5/9/2015 00:16 |
| Methyl acetate            | < 2.00      | ug/L | 5/9/2015 00:16 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L | 5/9/2015 00:16 |
| Methylcyclohexane         | < 2.00      | ug/L | 5/9/2015 00:16 |
| Methylene chloride        | < 5.00      | ug/L | 5/9/2015 00:16 |
| o-Xylene                  | < 2.00      | ug/L | 5/9/2015 00:16 |
| Styrene                   | < 5.00      | ug/L | 5/9/2015 00:16 |
| Tetrachloroethene         | < 2.00      | ug/L | 5/9/2015 00:16 |
| Toluene                   | < 2.00      | ug/L | 5/9/2015 00:16 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L | 5/9/2015 00:16 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L | 5/9/2015 00:16 |
| Trichloroethene           | <b>2.12</b> | ug/L | 5/9/2015 00:16 |
| Trichlorofluoromethane    | < 2.00      | ug/L | 5/9/2015 00:16 |
| Vinyl chloride            | <b>4.89</b> | ug/L | 5/9/2015 00:16 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B106-MW-PI12

Lab Sample ID: 151696-12

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>109</b>              | 82.3 - 115    |                 | 5/9/2015             | 00:16 |
| 4-Bromofluorobenzene  | <b>90.0</b>             | 85.5 - 111    |                 | 5/9/2015             | 00:16 |
| Pentafluorobenzene    | <b>100</b>              | 91.2 - 107    |                 | 5/9/2015             | 00:16 |
| Toluene-D8            | <b>97.0</b>             | 90.9 - 108    |                 | 5/9/2015             | 00:16 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22664.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



**PARADIGM**  
ENVIRONMENTAL SERVICES, INC.

**Lab Project ID:** 151696

**Client:** **Stantec**

**Project Reference:** Carriage Factory

---

**Sample Identifier:** LI-B106-MW-PI12

**Lab Sample ID:** 151696-12

**Date Sampled:** 5/5/2015

**Matrix:** Groundwater

**Date Received:** 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 1.5           | mg/L         |                  | 5/13/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

*Report Prepared Friday, May 15, 2015*



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-B108-MW-PI12

Lab Sample ID: 151696-13

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Metals**

| <b><u>Analyte</u></b> | <b><u>Result</u></b> | <b><u>Units</u></b> | <b><u>Qualifier</u></b> | <b><u>Date Analyzed</u></b> |
|-----------------------|----------------------|---------------------|-------------------------|-----------------------------|
| Arsenic               | < 0.0100             | mg/L                |                         | 5/8/2015 17:02              |
| Iron                  | <b>0.999</b>         | mg/L                |                         | 5/8/2015 17:02              |
| Manganese             | <b>0.0818</b>        | mg/L                |                         | 5/8/2015 17:02              |
| Sodium                | <b>42.4</b>          | mg/L                |                         | 5/8/2015 17:02              |

Method Reference(s): EPA 6010C

EPA 3005

Preparation Date: 5/7/2015

Data File: 050815a

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B108-MW-PI12

Lab Sample ID: 151696-13

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/8/2015 23:52 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/8/2015 23:52 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/8/2015 23:52 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| Acetone                     | < 10.0  | ug/L  |           | 5/8/2015 23:52 |
| Benzene                     | < 0.700 | ug/L  |           | 5/8/2015 23:52 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/8/2015 23:52 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/8/2015 23:52 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/8/2015 23:52 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B108-MW-PI12

Lab Sample ID: 151696-13

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

|                           |             |      |   |                |
|---------------------------|-------------|------|---|----------------|
| Chloroethane              | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Chloroform                | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Chloromethane             | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| cis-1,2-Dichloroethene    | <b>1.49</b> | ug/L | J | 5/8/2015 23:52 |
| cis-1,3-Dichloropropene   | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Cyclohexane               | < 10.0      | ug/L |   | 5/8/2015 23:52 |
| Dibromochloromethane      | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Dichlorodifluoromethane   | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Ethylbenzene              | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Freon 113                 | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Isopropylbenzene          | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| m,p-Xylene                | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Methyl acetate            | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Methyl tert-butyl Ether   | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Methylcyclohexane         | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Methylene chloride        | < 5.00      | ug/L |   | 5/8/2015 23:52 |
| o-Xylene                  | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Styrene                   | < 5.00      | ug/L |   | 5/8/2015 23:52 |
| Tetrachloroethene         | <b>14.4</b> | ug/L |   | 5/8/2015 23:52 |
| Toluene                   | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| trans-1,2-Dichloroethene  | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| trans-1,3-Dichloropropene | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Trichloroethene           | <b>2.72</b> | ug/L |   | 5/8/2015 23:52 |
| Trichlorofluoromethane    | < 2.00      | ug/L |   | 5/8/2015 23:52 |
| Vinyl chloride            | < 2.00      | ug/L |   | 5/8/2015 23:52 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: LI-B108-MW-PI12

Lab Sample ID: 151696-13

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | <b>104</b>              | 82.3 - 115    |                 | 5/8/2015             | 23:52 |
| 4-Bromofluorobenzene  | <b>93.0</b>             | 85.5 - 111    |                 | 5/8/2015             | 23:52 |
| Pentafluorobenzene    | <b>99.5</b>             | 91.2 - 107    |                 | 5/8/2015             | 23:52 |
| Toluene-D8            | <b>97.6</b>             | 90.9 - 108    |                 | 5/8/2015             | 23:52 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22663.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

---

Sample Identifier: LI-B108-MW-PI12

Lab Sample ID: 151696-13

Date Sampled: 5/5/2015

Matrix: Groundwater

Date Received: 5/5/2015

---

**Total Organic Carbon**

| <u>Analyte</u>       | <u>Result</u> | <u>Units</u> | <u>Qualifier</u> | <u>Date Analyzed</u> |
|----------------------|---------------|--------------|------------------|----------------------|
| Total Organic Carbon | 1.7           | mg/L         |                  | 5/13/2015            |

Method Reference(s): SM 5310 C  
Subcontractor ELAP ID: 10709

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: Trip Blank (T-614)

Lab Sample ID: 151696-14

Date Sampled: 5/4/2015

Matrix: Water

Date Received: 5/5/2015

**Volatile Organics**

| Analyte                     | Result  | Units | Qualifier | Date Analyzed  |
|-----------------------------|---------|-------|-----------|----------------|
| 1,1,1-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,1,2,2-Tetrachloroethane   | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,1,2-Trichloroethane       | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,1-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,1-Dichloroethene          | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2,3-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2,4-Trichlorobenzene      | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2-Dibromo-3-Chloropropane | < 10.0  | ug/L  |           | 5/7/2015 17:50 |
| 1,2-Dibromoethane           | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2-Dichloroethane          | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,2-Dichloropropane         | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,3-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,4-Dichlorobenzene         | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| 1,4-dioxane                 | < 20.0  | ug/L  |           | 5/7/2015 17:50 |
| 2-Butanone                  | < 10.0  | ug/L  |           | 5/7/2015 17:50 |
| 2-Hexanone                  | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| 4-Methyl-2-pentanone        | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| Acetone                     | < 10.0  | ug/L  |           | 5/7/2015 17:50 |
| Benzene                     | < 0.700 | ug/L  |           | 5/7/2015 17:50 |
| Bromochloromethane          | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| Bromodichloromethane        | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| Bromoform                   | < 5.00  | ug/L  |           | 5/7/2015 17:50 |
| Bromomethane                | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| Carbon disulfide            | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| Carbon Tetrachloride        | < 2.00  | ug/L  |           | 5/7/2015 17:50 |
| Chlorobenzene               | < 2.00  | ug/L  |           | 5/7/2015 17:50 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: Trip Blank (T-614)

Lab Sample ID: 151696-14

Date Sampled: 5/4/2015

Matrix: Water

Date Received: 5/5/2015

|                           |        |      |                |
|---------------------------|--------|------|----------------|
| Chloroethane              | < 2.00 | ug/L | 5/7/2015 17:50 |
| Chloroform                | < 2.00 | ug/L | 5/7/2015 17:50 |
| Chloromethane             | < 2.00 | ug/L | 5/7/2015 17:50 |
| cis-1,2-Dichloroethene    | < 2.00 | ug/L | 5/7/2015 17:50 |
| cis-1,3-Dichloropropene   | < 2.00 | ug/L | 5/7/2015 17:50 |
| Cyclohexane               | < 10.0 | ug/L | 5/7/2015 17:50 |
| Dibromochloromethane      | < 2.00 | ug/L | 5/7/2015 17:50 |
| Dichlorodifluoromethane   | < 2.00 | ug/L | 5/7/2015 17:50 |
| Ethylbenzene              | < 2.00 | ug/L | 5/7/2015 17:50 |
| Freon 113                 | < 2.00 | ug/L | 5/7/2015 17:50 |
| Isopropylbenzene          | < 2.00 | ug/L | 5/7/2015 17:50 |
| m,p-Xylene                | < 2.00 | ug/L | 5/7/2015 17:50 |
| Methyl acetate            | < 2.00 | ug/L | 5/7/2015 17:50 |
| Methyl tert-butyl Ether   | < 2.00 | ug/L | 5/7/2015 17:50 |
| Methylcyclohexane         | < 2.00 | ug/L | 5/7/2015 17:50 |
| Methylene chloride        | < 5.00 | ug/L | 5/7/2015 17:50 |
| o-Xylene                  | < 2.00 | ug/L | 5/7/2015 17:50 |
| Styrene                   | < 5.00 | ug/L | 5/7/2015 17:50 |
| Tetrachloroethene         | < 2.00 | ug/L | 5/7/2015 17:50 |
| Toluene                   | < 2.00 | ug/L | 5/7/2015 17:50 |
| trans-1,2-Dichloroethene  | < 2.00 | ug/L | 5/7/2015 17:50 |
| trans-1,3-Dichloropropene | < 2.00 | ug/L | 5/7/2015 17:50 |
| Trichloroethene           | < 2.00 | ug/L | 5/7/2015 17:50 |
| Trichlorofluoromethane    | < 2.00 | ug/L | 5/7/2015 17:50 |
| Vinyl chloride            | < 2.00 | ug/L | 5/7/2015 17:50 |

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



Lab Project ID: 151696

Client: **Stantec**

Project Reference: Carriage Factory

Sample Identifier: Trip Blank (T-614)

Lab Sample ID: 151696-14

Date Sampled: 5/4/2015

Matrix: Water

Date Received: 5/5/2015

| <b>Surrogate</b>      | <b>Percent Recovery</b> | <b>Limits</b> | <b>Outliers</b> | <b>Date Analyzed</b> |       |
|-----------------------|-------------------------|---------------|-----------------|----------------------|-------|
| 1,2-Dichloroethane-d4 | 95.5                    | 82.3 - 115    |                 | 5/7/2015             | 17:50 |
| 4-Bromofluorobenzene  | 90.6                    | 85.5 - 111    |                 | 5/7/2015             | 17:50 |
| Pentafluorobenzene    | 93.9                    | 91.2 - 107    |                 | 5/7/2015             | 17:50 |
| Toluene-D8            | 95.2                    | 90.9 - 108    |                 | 5/7/2015             | 17:50 |

Method Reference(s): EPA 8260C

EPA 5030

Data File: x22614.D

This report is part of a multipage document and should only be evaluated in its entirety. The Chain of Custody provides additional sample information, including compliance with the sample condition requirements upon receipt.

Report Prepared Friday, May 15, 2015



## Analytical Report Appendix

The reported results relate only to the samples as they have been received by the laboratory.

Each page of this document is part of a multipage report. This document may not be reproduced except in its entirety, without the prior consent of Paradigm Environmental Services, Inc.

All soil/sludge samples have been reported on a dry weight basis, unless qualified "reported as received". Other solids are reported as received.

Low level Volatiles blank reports for soil/solid matrix are based on a nominal 5 gram weight. Sample results and reporting limits are based on actual weight, which may be more or less than 5 grams.

The Chain of Custody provides additional information, including compliance with sample condition requirements upon receipt. Sample condition requirements are defined under the 2003 NELAC Standard, sections 5.5.8.3.1 and 5.5.8.3.2.

NYSDOH ELAP does not certify for all parameters. Paradigm Environmental Services or the indicated subcontracted laboratory does hold certification for all analytes where certification is offered by ELAP unless otherwise specified. Aliquots separated for certain tests, such as TCLP, are indicated on the Chain of Custody and final reports with an "A" suffix.

Data qualifiers are used, when necessary, to provide additional information about the data. This information may be communicated as a flag or as text at the bottom of the report. Please refer to the following list of analyte-specific, frequently used data flags and their meaning:

*"<" = Analyzed for but not detected at or above the quantitation limit.*

*"E" = Result has been estimated, calibration limit exceeded.*

*"Z" = See case narrative.*

*"D" = Sample, Laboratory Control Sample, or Matrix Spike Duplicate results above Relative Percent Difference limit.*

*"M" = Matrix spike recoveries outside QC limits. Matrix bias indicated.*

*"B" = Method blank contained trace levels of analyte. Refer to included method blank report.*

*"J" = Result estimated between the quantitation limit and half the quantitation limit.*

*"L" = Laboratory Control Sample recovery outside accepted QC limits.*

*"P" = Concentration differs by more than 40% between the primary and secondary analytical columns.*

*"NC" = Not calculable. Applicable to RPD if sample or duplicate result is non-detect or estimated (see primary report for data flags). Applicable to MS if sample is greater or equal to ten times the spike added. Applicable to sample surrogates or MS if sample dilution is 10x or higher.*

*"\*" = Indicates any recoveries outside associated acceptance windows. Surrogate outliers in samples are presumed matrix effects. LCS demonstrates method compliance unless otherwise noted.*

*"(1)" = Indicates data from primary column used for QC calculation.*

# GENERAL TERMS AND CONDITIONS

## LABORATORY SERVICES

These Terms and Conditions embody the whole agreement of the parties in the absence of a signed and executed contract between the Laboratory (LAB) and Client. They shall supersede all previous communications, representations, or agreements, either verbal or written, between the parties. The LAB specifically rejects all additional, inconsistent, or conflicting terms, whether printed or otherwise set forth in any purchase order or other communication from the Client to the LAB. The invalidity or unenforceability in whole or in part of any provision, term, or condition hereof shall not affect in any way the validity or enforceability of the remainder of the Terms and Conditions. No waiver by LAB of any provision, term, or condition hereof or of any breach by or obligation of the Client hereunder shall constitute a waiver of such provision, term, or condition on any other occasion or a waiver of any other breach by or obligation of the Client. This agreement shall be administered and interpreted under the laws of the state which services are procured.

### **Warranty.**

Recognizing that the nature of many samples is unknown and that some may contain potentially hazardous components, LAB warrants only that it will perform testing services, obtain findings, and prepare reports in accordance with generally accepted analytical laboratory principles and practices at the time of performance of services. LAB makes no other warranty, express or implied.

### **Scope and Compensation.**

LAB agrees to perform the services described in the chain of custody to which these terms and conditions are attached. Unless the parties agree in writing to the contrary, the duties of LAB shall not be construed to exceed the services specifically described. LAB will use LAB default method for all tests unless specified otherwise on the Work Order.

Payment terms are net 30 days from the date of invoice. All overdue payments are subject to an interest charge of one and one-half percent (1-1/2%) per month or a portion thereof. Client shall also be responsible for costs of collection, including payment of reasonable attorney fees if such expense is incurred. The prices, unless stated, do not include any sale, use or other taxes. Such taxes will be added to invoice prices when required.

### **Prices.**

Compensation for services performed will be based on the current Lab Analytical Fee Schedule or on quotations agreed to in writing by the parties. Turnaround time based charges are determined from the time of resolution of all work order questions. Testimony, court appearances or data compilation for legal action will be charged separately. Evaluation and reporting of initial screening runs may incur additional fees.

### **Limitations of Liability.**

In the event of any error, omission, or other professional negligence, the sole and exclusive responsibility of LAB shall be to re-perform the deficient work at its own expense and LAB shall have no other liability whatsoever. All claims shall be deemed waived unless made in writing and received by LAB within ninety (90) days following completion of services.

LAB shall have no liability, obligation, or responsibility of any kind for losses, costs, expenses, or other damages (including but not limited to any special, direct, incidental or consequential damages) with respect to LAB's services or results.

All results provided by LAB are strictly for the use of its clients and LAB is in no way responsible for the use of such results by clients or third parties. All reports should be considered in their entirety, and LAB is not responsible for the separation, detachment, or other use of any portion of these reports. Client may not assign the lab report without the written consent of the LAB.

Client covenants and agrees, at its/his/her sole expense, to indemnify, protect, defend, and save harmless the LAB from and against any and all damages, losses, liabilities, obligations, penalties, claims, litigation, demands, defenses, judgments, suits, actions, proceedings, costs, disbursements and/or expenses (including, without limitation attorneys' and experts' fees and disbursements) of any kind whatsoever which may at any time be imposed upon, incurred by or asserted or awarded against client relating to, resulting from or arising out of (a) the breach of this agreement by this client, (b) the negligence of the client in handling, delivering or disclosing any hazardous substance, (c) the violation of the Client of any applicable law, (d) non-compliance by the Client with any environmental permit or (e) a material misrepresentation in disclosing the materials to be tested.

### **Hazard Disclosure.**

Client represents and warrants that any sample delivered to LAB will be preceded or accompanied by complete written disclosure of the presence of any hazardous substances known or suspected by Client. Client further warrants that any sample containing any hazardous substance that is to be delivered to LAB will be packaged, labeled, transported, and delivered properly and in accordance with applicable laws.

### **Sample Handling.**

Prior to LAB's acceptance of any sample (or after any revocation of acceptance), the entire risk of loss or of damage to such sample remains with Client. Samples are accepted when receipt is acknowledged on chain of custody documentation. In no event will LAB have any responsibility for the action or inaction of any carrier shipping or delivering any sample to or from LAB premises. Client authorizes LAB to proceed with the analysis of samples as received by the laboratory, recognizing that any samples not in compliance with all current DOH-ELAP-NELAP requirements for containers, preservation or holding time will be noted as such on the final report.

Disposal of hazardous waste samples is the responsibility of the Client. If the Client does not wish such samples returned, LAB may add storage and disposal fees to the final invoice. Maximum storage time for samples is 30 days after completion of analysis unless modified by applicable state or federal laws. Client will be required to give the LAB written instructions concerning disposal of these samples.

LAB reserves the absolute right, exercisable at any time, to refuse to receive delivery of, refuse to accept, or revoke acceptance of any sample, which, in the sole judgment of LAB (a) is of unsuitable volume, (b) may be or become unsuitable for or may pose a risk in handling, transport, or processing for any health, safety, environmental or other reason whether or not due to the presence in the sample of any hazardous substance, and whether or not such presence has been disclosed to LAB by Client or (c) if the condition or sample date make the sample unsuitable for analysis.

### **Legal Responsibility.**

LAB is solely responsible for performance of this contract, and no affiliated company, director, officer, employee, or agent shall have any legal responsibility hereunder, whether in contract or tort including negligence.

### **Assignment.**

LAB may assign its performance obligations under this contract to other parties, as it deems necessary. LAB shall disclose to Client any assignee (subcontractor) by ELAP ID # on the submitted final report.


### **Force Majeure.**

LAB shall have no responsibility or liability to the Client for any failure or delay in performance by LAB, which results in whole or in part from any cause or circumstance beyond the reasonable control of LAB. Such causes and circumstances shall include, but not limited to, acts of God, acts or orders of any government authority, strikes or other labor disputes, natural disasters, accidents, wars, civil disturbances, difficulties or delays in transportation, mail or delivery services, inability to obtain sufficient services or supplies from LAB's usual suppliers, or any other cause beyond LAB's reasonable control.

### **Law.**

This contract shall be continued under the laws of the State of New York without regard to its conflicts of laws provision.

## CHAIN OF CUSTODY



| REPORT TO:        |                   | INVOICE TO:                                                                                                                                                                         |               | LAB PROJECT ID            |  |
|-------------------|-------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|---------------------------|--|
| CLIENT:           | Starck            | CLIENT:                                                                                                                                                                             | Starck        | 151696                    |  |
| ADDRESS:          | 61 Commercial St. | ADDRESS:                                                                                                                                                                            |               | Quotation #:              |  |
| CITY:             | Rochester         | CITY:                                                                                                                                                                               |               | Email:                    |  |
| STATE:            | NY                | STATE:                                                                                                                                                                              |               | mike.storonsky@starck.com |  |
| ZIP:              | 14614             | ZIP:                                                                                                                                                                                |               |                           |  |
| PHONE:            | 413-5266          | PHONE:                                                                                                                                                                              | 978-5248      |                           |  |
| ATTN:             | Mike Storonsky    | ATTN:                                                                                                                                                                               | Ben Haravitch |                           |  |
| PROJECT REFERENCE |                   | Carriage Factory                                                                                                                                                                    |               |                           |  |
| Matrix Codes:     |                   | WA - Water<br>WG - Groundwater<br>DW - Drinking Water<br>WW - Wastewater<br>SO - Soil<br>SL - Sludge<br>SD - Solid<br>PT - Paint<br>WP - Wipe<br>CK - Caulk<br>OL - Oil<br>AR - Air |               |                           |  |
| Matrix Codes:     |                   | AQ - Aqueous Liquid<br>NQ - Non-Aqueous Liquid                                                                                                                                      |               |                           |  |

| DATE COLLECTED | TIME COLLECTED | C O M P O S I T E | G R A B | SAMPLE IDENTIFIER     | M C A D R E I S | C O N T A I N E R S | REQUESTED ANALYSIS | REMARKS | PARADIGM LAB SAMPLE NUMBER |
|----------------|----------------|-------------------|---------|-----------------------|-----------------|---------------------|--------------------|---------|----------------------------|
| 1 5/5/15       | 0930           |                   | X       | LI-RW-1-PI12          | WGS             | 5                   |                    |         | 01                         |
| 2 5/5/15       | 1020           |                   |         | LI-RW-2-PI12          |                 | 5                   |                    |         | 02                         |
| 3 5/5/15       | 1140           |                   |         | LI-RW-3-PI12          |                 | 5                   |                    |         | 03                         |
| 4 5/4/15       | 0855           |                   |         | LI-RW-4-PI12          |                 | 5                   |                    |         | 04                         |
| 5 5/4/15       | 1125           |                   |         | LI-RW-5-PI12          |                 | 5                   |                    |         | 05                         |
| 6 5/4/15       | 1415           |                   |         | LI-RW-6-PI12          |                 | 5                   |                    |         | 06                         |
| 7 5/4/15       | 1525           |                   |         | LI-RW-7-PI12          |                 | 5                   |                    |         | 07                         |
| 8 5/4/15       | 1640           |                   |         | LI-RW-9-PI12 (MS/MSD) |                 | 15                  |                    |         | 08                         |
| 9 5/4/15       | 1310           |                   |         | LI-RW-12-PI12         |                 | 5                   |                    |         | 09                         |
| 10 5/4/15      | 1020           |                   |         | LI-DUP-PI12           |                 | 5                   |                    |         | 10                         |

| Turnaround Time                                                       | Report Supplements |
|-----------------------------------------------------------------------|--------------------|
| Availability contingent upon lab approval; additional fees may apply. |                    |
| Standard 5 day                                                        | Batch QC           |
| Rush 3 day                                                            | Category A         |
| Rush 2 day                                                            | Category B         |
| Rush 1 day                                                            | Other              |
| Other                                                                 | Other EDD          |
| please indicate:                                                      | please indicate:   |
| today                                                                 | Starck             |

Sampled By: James Best Date/Time: 5/5/15 16:15  
 Relinquished By: James Best Date/Time: 5/5/15 16:20  
 Received By: James Best Date/Time: 5/5/15 16:20  
 Received @ Lab By: James Best Date/Time: 5/5/15 16:53


P.I.F.



60cc each 5/5/15 16:47  
 Custody Seals w/ 1 Samples delivered by client. on 5/5/15

1083

# CHAIN OF CUSTODY



|                                   |                           |                                                                                                                                                                                                                             |  |                            |  |
|-----------------------------------|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|----------------------------|--|
| <b>REPORT TO:</b>                 |                           | <b>INVOICE TO:</b>                                                                                                                                                                                                          |  | <b>LAB PROJECT ID</b>      |  |
| CLIENT: <u>Stattac</u>            | CLIENT: <u>Stattac</u>    | LAB PROJECT ID: <u>151696</u>                                                                                                                                                                                               |  | Quotation #: <u>151696</u> |  |
| ADDRESS: <u>61 Commercial St.</u> | ADDRESS: <u>Stattac</u>   | Email: <u>mike.stormsky@stattac.com</u>                                                                                                                                                                                     |  |                            |  |
| CITY: <u>Rochester</u>            | CITY: <u>Stattac</u>      |                                                                                                                                                                                                                             |  |                            |  |
| STATE: <u>NY</u>                  | STATE: <u>NY</u>          |                                                                                                                                                                                                                             |  |                            |  |
| ZIP: <u>14614</u>                 | ZIP: <u>14614</u>         |                                                                                                                                                                                                                             |  |                            |  |
| PHONE: <u>413-5266</u>            | PHONE: <u>978-5248</u>    |                                                                                                                                                                                                                             |  |                            |  |
| ATTN: <u>Mike Stormsky</u>        | ATTN: <u>Ben Havavick</u> |                                                                                                                                                                                                                             |  |                            |  |
| <b>PROJECT REFERENCE</b>          |                           | <b>REQUESTED ANALYSIS</b>                                                                                                                                                                                                   |  |                            |  |
| <u>Carriage Factory</u>           |                           | Matrix Codes: WA - Water, WG - Groundwater, SO - Soil, SD - Solid, PT - Paint, WP - Wipe, CK - Caulk, OL - Oil, AR - Air<br>AQ - Aqueous Liquid, NQ - Non-Aqueous Liquid, DW - Drinking Water, WW - Wastewater, SL - Sludge |  |                            |  |

| DATE COLLECTED | TIME COLLECTED         | C O M P O S I T E | G R A B | SAMPLE IDENTIFIER  | M C A O T D R E I S | C O N T A I N E R N U M B E R S | REMARKS | PARADIGM LAB SAMPLE NUMBER |
|----------------|------------------------|-------------------|---------|--------------------|---------------------|---------------------------------|---------|----------------------------|
| 1 5/4/15       | 10:15                  |                   | X       | LI-B102-MW-PI12    | WG                  | 5                               | X X X X | 1 1                        |
| 2 5/5/15       | 12:45                  |                   | X       | LI-B106-MW-PI12    | WG                  | 5                               | X X X   | 1 2                        |
| 3 5/5/15       | 12:45                  |                   | X       | LI-B108-MW-PI12    | WG                  | 5                               | X X X X | 1 3                        |
| 4 5/4/15       | —                      |                   |         | Trip Blank (T-614) | WA                  | 1                               | X       | 1 4                        |
| 5              | per Trip Blank (T-614) |                   |         |                    |                     |                                 |         |                            |
| 6              | per 5/5/15             |                   |         |                    |                     |                                 |         |                            |
| 7              |                        |                   |         |                    |                     |                                 |         |                            |
| 8              |                        |                   |         |                    |                     |                                 |         |                            |
| 9              |                        |                   |         |                    |                     |                                 |         |                            |
| 10             |                        |                   |         |                    |                     |                                 |         |                            |

|                                                                       |                                     |                                |                                     |
|-----------------------------------------------------------------------|-------------------------------------|--------------------------------|-------------------------------------|
| <b>Turnaround Time</b>                                                |                                     | <b>Report Supplements</b>      |                                     |
| Availability contingent upon lab approval; additional fees may apply. |                                     |                                |                                     |
| Standard 5 day                                                        | <input type="checkbox"/>            | Batch QC                       | <input type="checkbox"/>            |
| Rush 3 day                                                            | <input type="checkbox"/>            | Category A                     | <input type="checkbox"/>            |
| Rush 2 day                                                            | <input type="checkbox"/>            | Category B                     | <input checked="" type="checkbox"/> |
| Rush 1 day                                                            | <input type="checkbox"/>            |                                |                                     |
| Other                                                                 | <input checked="" type="checkbox"/> | Other EDD                      | <input type="checkbox"/>            |
| please indicate: <u>10 day</u>                                        |                                     | please indicate: <u>static</u> |                                     |

|                                     |                                |                               |
|-------------------------------------|--------------------------------|-------------------------------|
| Sampled By: <u>Anna Best</u>        | Date/Time: <u>5/5/15 10:15</u> | Total Cost: <u>          </u> |
| Relinquished By: <u>Anna Best</u>   | Date/Time: <u>5/5/15 10:20</u> |                               |
| Received By: <u>Anna Best</u>       | Date/Time: <u>5/5/15 10:20</u> |                               |
| Received @ Lab By: <u>Anna Best</u> | Date/Time: <u>5/5/15 16:53</u> |                               |




## Chain of Custody Supplement

Client: Startec Completed by: Glenn Pezzullo  
 Lab Project ID: 151696 Date: 5/5/15

### Sample Condition Requirements

Per NELAC/ELAP 210/241/242/243/244

| NELAC compliance with the sample condition requirements upon receipt |                                     |                          |                                            |
|----------------------------------------------------------------------|-------------------------------------|--------------------------|--------------------------------------------|
| Condition                                                            | Yes                                 | No                       | N/A                                        |
| Container Type                                                       | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Transferred to method-compliant container                            | <input type="checkbox"/>            | <input type="checkbox"/> | <input checked="" type="checkbox"/>        |
| Headspace (<1 mL)                                                    | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> metals |
| Comments                                                             |                                     |                          |                                            |
| Preservation                                                         | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Chlorine Absent (<0.10 ppm per test strip)                           | <input type="checkbox"/>            | <input type="checkbox"/> | <input checked="" type="checkbox"/>        |
| Comments                                                             |                                     |                          |                                            |
| Holding Time                                                         | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |
| Temperature                                                          | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> metals |
| Comments                                                             | 6°C iced                            |                          |                                            |
| Sufficient Sample Quantity                                           | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/>                   |
| Comments                                                             |                                     |                          |                                            |



150502018 179 Lake Avenue, Rochester, NY 14608 Office (585) 647-2530 Fax (585) 647-3311

# CHAIN OF CUSTODY

ADIRONDACK: ELAP ID: 10709

|                                                                                         |                        |                   |          |
|-----------------------------------------------------------------------------------------|------------------------|-------------------|----------|
| REPORT TO:                                                                              |                        | INVOICE TO:       |          |
| COMPANY: Paradigm Environmental                                                         | ADDRESS:               | COMPANY: Same     | ADDRESS: |
| CITY:                                                                                   | STATE:                 | CITY:             | STATE:   |
| PHONE:                                                                                  | FAX:                   | PHONE:            | FAX:     |
| ATTN: Kate Hansen                                                                       | ATTN: Meredith Dillman | LAB PROJECT #:    |          |
| COMMENTS: Please email results to khansen@paradigmenv.com and reporting@paradigmenv.com |                        | CLIENT PROJECT #: |          |

| DATE    | TIME   | COMMENTS | GRADES    | SAMPLE LOCATION/FIELD ID | MATERIALS | CONTAMINANTS | TESTS | REMARKS                                                              | PARADIGM LAB SAMPLE NUMBER |
|---------|--------|----------|-----------|--------------------------|-----------|--------------|-------|----------------------------------------------------------------------|----------------------------|
| 15/5/15 | 09:30  | X        | 151696-01 | Ground water             | 2         | X            | TOC   | Report 5/14/15. ASP Cat B Package due 5/27. SDG closed. SW-846 HT'S. | 001                        |
| 2       | 10:20  |          |           |                          |           |              |       | LI-RW-2-PI12                                                         | 002                        |
| 3       | 11:40  |          |           |                          |           |              |       | LI-RW-3-PI12                                                         | 003                        |
| 4       | 5/4/15 | 08:55    |           |                          |           |              |       | LI-RW-4-PI12                                                         | 004                        |
| 5       | 11:25  |          |           |                          |           |              |       | LI-RW-5-PI12                                                         | 005                        |
| 6       | 14:15  |          |           |                          |           |              |       | LI-RW-6-PI12                                                         | 006                        |
| 7       | 15:25  |          |           |                          |           |              |       | LI-RW-7-PI12                                                         | 007                        |
| 8       | 16:40  |          |           |                          |           |              |       | LI-RW-9-PI12                                                         | 008                        |
| 9       | 13:10  |          |           |                          |           |              |       | LI-RW-12-PI12                                                        | 009                        |
| 10      | 10:20  |          |           |                          |           |              |       | LI-RW-08-PI12                                                        | 010                        |

LAB USE ONLY BELOW THIS LINE  
Sample Condition: Per NELAC/ELAP 210/241/242/243/244  
Receipt Parameter: NELAC Compliance  
Comments: Container Type: Y ☐ N ☐  
Preservation: Y ☐ N ☐  
Holding Time: Y ☐ N ☐  
Temperature: 40C Y ☐ N ☐  
Comments: Y ☐ N ☐

|                                |                            |                   |             |
|--------------------------------|----------------------------|-------------------|-------------|
| Client                         | Sampled By: [Signature]    | Date/Time: 5/6/15 | Total Cost: |
| Relinquished By: [Signature]   | Date/Time: 16:00           |                   |             |
| Received By: [Signature]       | Date/Time: 5/7/15 10:13 AM | P.L.F.            |             |
| Received @ Lab By: [Signature] | Date/Time: [Blank]         |                   |             |



24

ADIRONDACK. EI AP ID-10709

REQUESTED ANALYSIS

LAB USE ONLY BELOW THIS LINE

| Y                        | Z                        |
|--------------------------|--------------------------|
| <input type="checkbox"/> | <input type="checkbox"/> |

## Certificate of Analysis: Gene-Trac® *Dehalococcoides* Assay

**Customer:** Mike Storonsky, Stantec

**Project:** Carriage Factory

**Customer Reference:** 190500751

**SiREM Reference:** S-3572

**Report Date:** 21-May-15

**Data Files:** iQ5-DHC-QPCR-1227  
iQ5-DB-DHC-QPCR-0581

**Table 1a: Test Results**

| Customer Sample ID | SiREM Sample ID | Sample Collection Date | Sample Matrix | Percent Dhc * | <i>Dehalococcoides</i> Enumeration/Liter ** |
|--------------------|-----------------|------------------------|---------------|---------------|---------------------------------------------|
| LI-B102-MW-PI12    | DHC-11718       | 4-May-15               | Groundwater   | 0.04 - 0.1 %  | $6 \times 10^5$                             |
| LI-B108-MW-PI12    | DHC-11719       | 5-May-15               | Groundwater   | 0.02 - 0.05 % | $3 \times 10^5$                             |

**Notes:**

\* Percent *Dehalococcoides* (Dhc) in microbial population. This value is calculated by dividing the number of Dhc 16S ribosomal ribonucleic acid (rRNA) gene copies by the total number of bacteria as estimated by the mass of DNA extracted from the sample. Range represents normal variation in Dhc enumeration.

\*\* Based on quantification of Dhc 16S rRNA gene copies. Dhc are generally reported to contain one 16S rRNA gene copy per cell; therefore, this number is often interpreted to represent the number of Dhc cells present in the sample.

J The associated value is an estimated quantity between the method detection limit and quantitation limit.

U Not detected, associated value is the quantification limit.

B Analyte was detected in the method blank within an order of magnitude of the test sample

NA Not applicable as *Dehalococcoides* not detected and/or quantifiable DNA not extracted from the sample.

I Sample inhibited the test reaction based on inability to PCR amplify extracted DNA with universal primers.

E Extracted genomic DNA was not detected in sample.

**Analyst:**   
**Ben Reside**  
Laboratory Technician

**Approved:**   
**Ximena Druar, B.Sc.**  
Genetic Testing Coordinator

## Certificate of Analysis: Gene-Trac® VC, Vinyl Chloride Reductase (*vcrA*) Assay

**Customer:** Mike Storonsky, Stantec

**Project:** Carriage Factory

**Customer Reference:** 190500751

**SiREM Reference:** S-3572

**Report Date:** 21-May-15

**Data Files:** MyiQ-VC-QPCR-0763  
MyiQ-DB-VC-QPCR-0474

**Table 1b: Test Results**

| Customer Sample ID | SiREM Sample ID | Sample Collection Date | Sample Matrix | Percent <i>vcrA</i> * | Vinyl Chloride Reductase ( <i>vcrA</i> ) Gene Copies/Liter |
|--------------------|-----------------|------------------------|---------------|-----------------------|------------------------------------------------------------|
| LI-B102-MW-PI12    | VCR-5854        | 4-May-15               | Groundwater   | 0.05 - 0.2 %          | $1 \times 10^6$                                            |
| LI-B108-MW-PI12    | VCR-5855        | 5-May-15               | Groundwater   | 0.02 - 0.06 %         | $4 \times 10^5$                                            |

**Notes:**

\* Percent *vcrA* in microbial population. This value is calculated by dividing the number of vinyl chloride reductase A (*vcrA*) gene copies quantified by the total number of bacteria estimated to be in the sample based on the mass of DNA extracted from the sample. Range represents normal variation in enumeration of *vcrA*.

J The associated value is an estimated quantity between the method detection limit and quantitation limit.

U Not detected, associated value is the quantification limit.

B Analyte was detected in the method blank within an order of magnitude of the test sample.

NA Not applicable as *vcrA* not detected and/or quantifiable DNA not extracted from the sample.

I Sample inhibited the test reaction based on inability to PCR amplify extracted DNA with universal primers.

C Correction factor applied to correct for non-specific PCR amplification products, value is an estimated quantity.

**Analyst:**

  
**Ben Reside**  
Laboratory Technician

**Approved:**

  
**Ximena Druar, B.Sc.**  
Genetic Testing Coordinator

**Table 2: Detailed Test Parameters, Gene-Trac Test Reference S-3572**

|                                           |                 |                 |
|-------------------------------------------|-----------------|-----------------|
| Customer Sample ID                        | LI-B102-MW-PI12 | LI-B108-MW-PI12 |
| SiREM Dhc Sample ID                       | DHC-11718       | DHC-11719       |
| SiREM <i>vcrA</i> Sample ID               | VCR-5854        | VCR-5855        |
| Date Received *                           | 7-May-15        | 7-May-15        |
| Sample Temperature                        | 4.4 °C          | 4.4 °C          |
| Filtration Date *                         | 12-May-15       | 12-May-15       |
| Volume Used for DNA Extraction            | 500 mL          | 500 mL          |
| DNA Extraction Date                       | 14-May-15       | 14-May-15       |
| DNA Concentration in Sample (extractable) | 3549 ng/L       | 3625 ng/L       |
| PCR Amplifiable DNA                       | Detected        | Detected        |
| Dhc qPCR Date Analyzed                    | 14-May-15       | 14-May-15       |
| <i>vcrA</i> qPCR Date Analyzed            | 14-May-15       | 14-May-15       |
| Laboratory Controls (see Tables 3 & 4)    | Passed          | Passed          |
| Comments                                  | --              | --              |

**Notes:**

Refer to Tables 3 & 4 for detailed results of controls.

°C = degrees Celsius

Dhc = *Dehalococcoides*

DNA = Deoxyribonucleic acid

mL = milliliters

ng/L = nanograms per liter

NA = Not Applicable

ND = Not Detected

PCR = polymerase chain reaction

qPCR = quantitative PCR

*vcrA* = vinyl chloride reductase

\*Samples are stabilized by freezing at -80 °C upon sample reception (field filters) or in-lab filtration (groundwater). Hold time not exceeded if sampling date is within 14 days of date received.

**Table 3: Gene-Trac Dhc Control Results, Test Reference S-3572**

| Laboratory Control                         | Analysis Date | Control Description                    | Spiked<br>Dhc 16S rRNA Gene<br>Copies per Liter | Recovered<br>Dhc 16S rRNA Gene<br>Copies per Liter | Comments |
|--------------------------------------------|---------------|----------------------------------------|-------------------------------------------------|----------------------------------------------------|----------|
| <b>Positive Control Low Concentration</b>  | 14-May-15     | qPCR with KB1 genomic DNA (CSLD-0865)  | $4.2 \times 10^6$                               | $4.4 \times 10^6$                                  | --       |
| <b>Positive Control High Concentration</b> | 14-May-15     | qPCR with KB1 genomic DNA (CSHD-0865)  | $1.0 \times 10^{10}$                            | $1.3 \times 10^{10}$                               | --       |
| <b>DNA Extraction Blank</b>                | 14-May-15     | DNA extraction sterile water (FB-2402) | 0                                               | $2.6 \times 10^3$ U                                | --       |
| <b>Negative Control</b>                    | 14-May-15     | Tris Reagent Blank (TBD-0824)          | 0                                               | $2.6 \times 10^3$ U                                | --       |

**Notes:**

Dhc = *Dehalococcoides*

DNA = Deoxyribonucleic acid

qPCR = quantitative PCR

16S rRNA = 16S ribosomal ribonucleic acid

U Not detected, associated value is the quantification limit.

**Table 4: Gene-Trac VC Control Results, Test Reference S-3572**

| Laboratory Control                         | Analysis Date | Control Description                    | Spiked <i>vcrA</i> reductase Gene Copies per Liter | Recovered <i>vcrA</i> reductase Gene Copies per Liter | Comments |
|--------------------------------------------|---------------|----------------------------------------|----------------------------------------------------|-------------------------------------------------------|----------|
| <b>Positive Control Low Concentration</b>  | 14-May-15     | qPCR with KB1 genomic DNA (CSLV-0631)  | $1.1 \times 10^8$                                  | $1.6 \times 10^8$                                     | --       |
| <b>Positive Control High Concentration</b> | 14-May-15     | qPCR with KB1 genomic DNA (CSHV-0631)  | $2.9 \times 10^{10}$                               | $2.0 \times 10^{10}$                                  | --       |
| <b>DNA Extraction Blank</b>                | 14-May-15     | DNA extraction sterile water (FB-2402) | 0                                                  | $2.6 \times 10^3$ U                                   | --       |
| <b>Negative Control</b>                    | 14-May-15     | Tris Reagent Blank (TBV-0602)          | 0                                                  | $2.6 \times 10^3$ U                                   | --       |

**Notes:**

DNA = Deoxyribonucleic acid

qPCR = quantitative PCR

16S rRNA = 16S ribosomal ribonucleic acid

U Not detected, associated value is the quantification limit.

*vcrA* = vinyl chloride reductase





**Stantec Consulting Services Inc.**  
61 Commercial Street, Suite 100  
Rochester NY 14614-1009  
Tel: (585) 475-1440  
Fax: (585) 272-1814

## **Appendix B**

### **Data Usability Summary Report (DUSR) – October 2014 and February 2015**

# Data Validation Services

120 Cobble Creek Road P.O. Box 208

North Creek, NY 12853

Phone 518-251-4429

harry@frontiernet.net

May 22, 2015

Laura Best  
Stantec  
61 Commercial St.  
Rochester, NY 14614

RE: Data Usability Summary Report (DUSR)  
Validation of the 33 Litchfield Old Carriage Factory Remediation Site Analytical Data  
Paradigm SDG Nos. 144730 and 150382

Dear Ms. Best:

Review has been completed for the data packages generated by Paradigm Environmental Services, Inc that pertain to samples collected between 10/28/14 and 02/04/15 at the 33 Litchfield Carriage Factory site. Twenty-four aqueous samples and two field duplicates were analyzed for TCL volatiles, total sodium, and TOC. Six of the samples and the field duplicates were also processed for total arsenic, iron, and manganese. Matrix spikes and trip blanks were also processed. Analytical methodologies are those of the USEPA SW846 6010, 8260, and Standard Methods SM5310C. TOC results were subcontracted to Adirondack Environmental Services (AES).

The data package submitted by the laboratory contains full deliverables for validation, but this usability report is generated from review of the QC summary form information, with full review of sample raw data and limited review of associated QC raw data. The reported QC summary forms and sample raw data have been reviewed for application of validation qualifiers, in accordance with the project QAPP, with guidance from the USEPA national and regional validation documents, and in consideration for the specific requirements of the analytical methodology. The following items were reviewed:

- \* Data Completeness
- \* Case Narrative
- \* Custody Documentation/Sample Receipt
- \* Holding Times
- \* Surrogate and Internal Standard Recoveries
- \* Calibration/Trip/Method Blanks
- \* Laboratory Control Sample (LCS)
- \* Blind Field Duplicate Correlations
- \* Instrumental Tunes
- \* Calibration Standards
- \* ICP Serial Dilution Evaluations
- \* ICP Interference Check Samples
- \* Method Compliance
- \* Sample Result Verification

The data review includes evaluation of the specific items noted in The NYS DER-10 Appendix B section 2.0 (c). The items listed above that show deficiencies are discussed within the text of this narrative. The laboratory QC forms illustrating the excursions can be found within the laboratory data package.

**In summary**, most results are usable either as reported or with minor qualification or edit. However, the following item of concern is noted:

- the reporting limits for benzene in aqueous samples are edited upward (by 43%) to reflect processing

The accuracy and precision evaluations, data completeness, representativeness, and the analytical method comparability are acceptable.

Copies of the client sample identifications are attached to this text, and should be reviewed in conjunction with this report. Also attached are the client EDD files, with recommended qualifiers/edits applied in red.

#### **Chain-of-Custody/Sample Receipt**

The laboratory chains of custody do not have sufficient fields for relinquish entries, and therefore the final receipt entries are not preceded by relinquish entries.

The custody forms should have fields to indicate preservation. The supplement to the chain-of-custody notes “yes” to indicate that the samples were preserved. Volatiles were processed within the holding time for unpreserved samples. Neither the volatile nor the metals preparation/analysis logs include the pH of the samples.

#### **General**

The data deliverables are not in compliance with NYS category B. The Paradigm sample report forms do not include the required information such as volume of sample, dilution factor, preparation date etc.

The Paradigm data packages do not provide a method summary nor are the specific methods mentioned anywhere in the data packages.

The method 415.1 was requested on the chain-of-custody for the TOC analyses. The laboratory utilized method SM 5310C.

#### **Field Duplicate Correlations**

Aqueous field duplicates were collected at locations LI-RW-9-PI6 and LI-B102-MW-PI9, and show acceptable correlations.

#### **Volatile Analyses by EPA8260**

A reporting limit for benzene in aqueous samples of 0.7 ug/L (prior to the application of the dilution factor) was reported by the laboratory. However, the lowest initial calibration standard concentration is 1ug/L, and detection of that analyte at a lower concentration is not assured. Per the analytical protocol, the reporting limit in the samples has been elevated accordingly. Additionally detections of benzene below the established linear range have been qualified as estimated in value.

The following detections are qualified tentative in identification and estimated in value due to poor mass spectral identification:

- vinyl chloride in LI-RW2-PI6, LI-RW7-PI6, LI-RW-2-PI9, LI-RW-3-PI9, LI-RW-5-PI9, LI-RW-5-PI9, LI-RW-12-PI9, and LI-B102-MW-PI9
- methyl acetate in LI-RW2-PI6

The detection of LI-RW-1-PI9 are qualified tentative in identification and estimated in value due to poor mass spectral identification:

The matrix spikes of LI-B102-MW-PI6 and LI-B108-MW-PI9 show acceptable recoveries and correlations, with the exception of one slightly low recovery for bromomethane. No qualification to the data is made.

Calibration standards show acceptable responses, with the following exceptions, results for which are qualified as estimated in the indicated samples:

- 2-butanone (low RRFs) in all project samples
- acetone (low RRFs) in the samples reported in SDG 150382

The date of analysis for Trip Blank T-570 entered onto the tune QC summary Form 5 is incorrect.

#### **Arsenic, Iron Manganese, and Sodium by EPA 6010**

Matrix spikes/duplicate evaluations for arsenic, iron, manganese, and sodium on LI-B102-MW-PI6 and LI-B108-MW-PI9 show acceptable recoveries and correlations.

An ICP serial dilution evaluation was performed on LI-B102-MW-PI6, with correlations falling within validation guidelines.

The required low level calibration standard analyses were not reported with the data reported in SDG 144730. Therefore accuracy of reporting limits and low level detections has not been fully evaluated.

The blanks show no contamination affecting sample reported results.

#### **TOC by SM5310C**

Review was conducted for method compliance, holding times, transcription, calculations, standard and blank acceptability, accuracy and precision, etc., as applicable to each procedure.


Matrix spike/duplicate evaluations were performed for TOC on LI-B102-MW-PI6 and LI-B108-MW-PI9. That for the former shows acceptable recoveries and correlations. The matrix spike of the latter produced a low recovery (50%), with a marginally elevated correlation. Therefore, the result TOC in LI-B108-MW-PI9 has been qualified as estimated.

Due to outlying recoveries of associated low concentration standards, the following TOC results are qualified as estimated:

- LI-B106-MW-PI6, LI-RW9-PI6, and LI-DUP-PI6 have a low bias (51% recovery)
- LI-RW5-PI9, LI-RW7-PI9, LI-RW9-PI9, and LI-B106-MW-PI9 have a high bias (186%)

Please do not hesitate to contact me if questions or comments arise during your review of this report.

Very truly yours,

A handwritten signature in black ink, appearing to read 'Judy Harry', with a long, sweeping horizontal stroke extending to the right.

Judy Harry

Att: Validation Qualifier Definitions  
Client and Laboratory Sample IDs  
Qualified Client EQuis EDDs

## VALIDATION DATA QUALIFIER DEFINITIONS

|             |                                                                                                                                                                                                                                      |
|-------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>U</b>    | The analyte was analyzed for, but was not detected above the level of the associated reported quantitation limit.                                                                                                                    |
| <b>J</b>    | The analyte was positively identified; the associated numerical value is an approximate concentration of the analyte in the sample.                                                                                                  |
| <b>J-</b>   | The analyte was positively identified; the associated numerical value is an estimated quantity that may be biased low.                                                                                                               |
| <b>J+</b>   | The analyte was positively identified; the associated numerical value is an estimated quantity that may be biased high.                                                                                                              |
| <b>UJ</b>   | The analyte was analyzed for, but was not detected. The associated reported quantitation limit is approximate and may be inaccurate or imprecise.                                                                                    |
| <b>NJ</b>   | The detection is tentative in identification and estimated in value. Although there is presumptive evidence of the analyte, the result should be used with caution as a potential false positive and/or elevated quantitative value. |
| <b>R</b>    | The data are unusable. The sample results are rejected due to serious deficiencies in meeting Quality Control limits. The analyte may or may not be present.                                                                         |
| <b>EMPC</b> | The results do not meet all criteria for a confirmed identification. The quantitative value represents the Estimated Maximum Possible Concentration of the analyte in the sample.                                                    |

## **CLIENT and LABORATORY SAMPLE IDs**

[illegible]

BATCH COMPLETE: 2/19/2015  
DATE: 2/19/2015  
DATE DUE: 3/4/2015  
PROTOCOL: SW846

[illegible]