REMEDIAL INVESTIGATION WORK PLAN

BROWNFIELD CLEANUP PROGRAM For **OCTAVUS STORAGE NY, LLC** 3678 Southwestern Boulevard, Orchard Park, NY 14127 BCP # C915342



Prepared For: **OCTAVUS STORAGE NY, LLC** 1810 West Kennedy Boulevard, Tampa, FL 33606 HEI Project No: e1728

Prepared By:

Hazard Evaluations, Inc. 3636 North Buffalo Road Orchard Park, New York 14127 Orchard Park, New York 14127 (716) 667-3130

Wittman GeoSciences, PLLC 3636 North Buffalo Road 716-574-1513

April 25, 2019





CERTIFICATIONI, Michele M. Wittman, P.G., certify that I am currently a NYS registered professional geologist and that this Remedial Investigation Work Plan was prepared in accordance with applicable statutes and regulations and in substantial conformance with the DER Technical Guidance for Site Investigation and Remediation (DER-10).

Michele M. Wittman, P.G.
Wittman GeoSciences, PLLC
New York License Number 000726

Date

New York License Number 000726

04/25/19





TABLE OF CONTENTS

	<u>Page</u>
1.0 INTRODUCTION	1
1.1 Project Background	
1.2 Site Background	
1.3 Summary of Environmental Conditions	2
1.4 Site Conditions	2
2.0 PROJECT OBJECTIVES	2
2.1 Regulatory Criteria	3
2.2 Project Organization	3
3.0 INVESTIGATION SCOPE OF WORK	4
3.1 Introduction	4
3.2 Building Survey	4
3.3 Field Investigation Activities	4
3.3.1 Surface Soil Investigation	4
3.3.2 Industrial Pond Samples	5
3.3.3 Test Pit Excavations	6
3.3.4 Subsurface Soil Investigation	6
3.3.5 Subsurface Soil Sample Laboratory Analysis	7
3.3.6 Monitoring Well Installation	7
3.3.7 Monitoring Well Development	8
3.3.8 Groundwater Sampling	8
3.3.9 Field Specific Quality Assurance/Quality Control Sam	pling9
3.4 Investigation- Derived Waste Management	10
3.5 Site Mapping	10
3.6 Personnel Decontamination	
3.7 Decontamination of Equipment	
3.8 Disposal of Contaminated Materials	
4.0 REMEDIAL INVESTIGATION/INTERIM REMEDIAL MEA ALTERNATAIVES ANALYSIS REPORT	
5.0 ADDITIONAL PROJECT DOCUMENTS	12
5.1 Quality Assurance Project Plan	12
5.2 Health and Safety Plan	12
6.0 PROJECT SCHEDULE	12





TABLE OF CONTENTS (Cont'd)

FIGURES

IIOOKEO	
Figure 1	Site Location Map
Figure 2	Site Limits
Figure 3	Previous Sampling Locations
Figure 4	Nearby Wetland and Floodplain Locations
Figure 5	Proposed Sampling Locations
Figure 6	Draft Concept Site Plan
Figure 7	Project Schedule

TABLES

Table 1 Analytical Testing Program Summary

APPENDIX

Appendix A	Historical Investigation Information
Appendix B	Quality Assurance Project Plan
Appendix C	Health & Safety Plan





1.0 INTRODUCTION

1.1 **Project Background**

This Remedial Investigation (RI) presents the proposed scope of work (Work Plan) at the future Octavus Storage NY, LLC (Octavus) facility located at 3678 Southwestern Boulevard in Orchard Park, New York (Site), as shown on Figure 1 and Figure 2. The Applicant, Octavus Storage NY, LLC, is requesting entry into the Brownfield Cleanup Program (BCP) as a Volunteer. The BCP Application has been submitted concurrently with this RI Work Plan.

The RI will be completed by Hazard Evaluations, Inc. (HEI) and Wittman GeoSciences, PLLC (WGS) on behalf of Octavus Storage NY, LLC. The work will be completed in general accordance with New York State Department of Environmental Conservation (NYSDEC) DER-10 guidelines. The work plan provides details on the site investigation to be undertaken. The Site investigation will be focused on subsurface conditions in areas of former incineration waste disposal and industrial scrubber discharge retention pond. Following RI and required remedial action, Octavus Storage NY, LLC will redevelop the property as a self-storage unit facility.

1.2 Site Background

The Site is addressed as 3678 Southwestern Boulevard located in Orchard Park, Erie County, New York. The Site consists of one parcel, totaling approximately 7.34 acres of land, and is bound to the north and east by Rout 219 Expressway, to the south by Southwestern Boulevard, and to the west by vacant land. The property is located within a commercial area of Orchard Park, utilized for light industrial, commercial and residential purposes.

The Site is improved with a two-story, approximate 4,756 square foot building. The parcel includes an asphalt driveway, concrete loading ramp, perimeter fence around portions, and wooded area in the northern portion which contains a former industrial discharge retention pond.

The subject site was initially utilized agriculturally, from at least 1926 through at least 1966. In 1967, the property was purchased by the County of Erie (Torrax) to be utilized for the disposal of municipal solid waste through incineration. The facility was identified as under construction in 1971, and operated by the Erie County Torrax Facility from 1972 to 1979. The building was used for an experimental pilot program for the incineration of refuse. During its operational phase, the building housed a high temperature incinerator capable of using 13,000 cu ft/hr natural gas. The maximum waste to be burned was 8,400 lb/hr. In 1974, the experimental program was halted and all of the major process equipment was removed. During the period of operation, the northern portion of the site was used as an ash landfill and scrubber discharge retention pond.

The building was then occupied between 1974 and 1977 by the Erie County Sewerage and Management Department as a vehicle repair and/or storage facility, and





has been vacant since that time. The property was purchased in 1993 by Torrax Building Joint Venture and remained unoccupied, but was used for occasional storage of equipment. The building is currently vacant.

1.3 <u>Summary of Environmental Conditions</u>

A Phase II Environmental Site Assessment (ESA) was completed by Great Lakes Environmental & Safety Consultants, Inc. (Great Lakes) in April 2018. The purpose of the investigation was to determine the presence of contaminants in the soil, and identify potential contamination source areas and mitigation pathways. The work included completion of 3 soil borings and 5 surface soil samples. Great Lakes indicated each boring location was composited for laboratory analysis from, and five (5) surface soil samples (depth not identified) were also collected for analysis. A copy of the Phase II report is included in Appendix A.

The Phase II testing identified metals in the fill areas associated with on-site ash disposal located north of the building within a wooded area. Elevated levels of barium and copper were present at concentrations exceeding the commercial use soil cleanup objective (CUSCO). In addition, chromium and lead were detected at concentrations exceeding restricted residential use soil cleanup objectives (RRUSCO). Groundwater samples were not collected during the Phase II ESA.

1.4 <u>Site Conditions</u>

Three soil borings were completed to a depth of approximately 12 feet below grade. However, Great Lakes did not provide information on subsurface conditions, nor were boring logs included within their report. Additionally, no information was provided on groundwater depth, if encountered.

Based on a review of the site topographic conditions as depicted on the USGS 7.5 minute Topographic Quadrangle Map of Buffalo SE, New York, shallow regional groundwater flow is expected in a northwesterly direction toward an unnamed tributary of South Branch Smokes Creek located approximately 500 feet northwest of the Site. This tributary flows into South Branch Smokes Creek and toward Lake Erie located approximately four miles west of the Site.

2.0 PROJECT OBJECTIVES

The Site has not been comprehensively characterized; therefore, Octavus Storage NY, LLC intends to further investigate the soil/fill, former industrial scrubber discharge retention pond and groundwater (if encountered) at the Site. Areas of concern (AOCs) identified during previous investigation work at the Site include an incinerator ash/cinder disposal area and a former industrial process pond" identified on Figure 5. Data collected during the RI will be used to identify potential health risks and to evaluate remedial alternatives. The objectives of the RI include the following:

 Define the nature and extent of on-site contamination in both soil and groundwater.





- Identify on-site source areas of contamination, if any.
- Collect data of sufficient quantity and quality to evaluate potential threats to the public health and environment.
- Collect data of sufficient quantity and quality to evaluate remedial alternatives.
- After completion of RI investigation work and laboratory data review, a
 determination will be made for the need of a soil vapor intrusion (SVI).

2.1 Regulatory Criteria

NYSDEC has applicable standards, criteria and guidance (SCG) values that will be used for this project. These goals are applicable when considering remedial alternatives. For purposes of the RI the following SCG will be utilized:

- 6 NYCRR Part 375-3 Brownfield Cleanup Program dated December 14, 2006.
- NYSDEC Policy CP-51/Soil Cleanup Guidance dated October 21, 2010.
- NYSDEC "DER-10 Technical guidance for Investigation and Remediation", dated May 2010.
- NYSDEC Division of Water Technical and Operational Guidance Series (1.1.1) document "Ambient Water Quality Standards and Guidance Values and Groundwater Effluent Limitations" dated June 1998, amended January 1999 Errata Sheet, April 2000 Addendum and June 2004 Addendum.
- State Department of Health (NYSDOH) "Guidance for Evaluating Soil Vapor Intrusion in the State of New York" dated October 2006.

Soil and groundwater samples will be collected in general accordance with NYSDEC and United Sates Environmental Protection Agency (USEPA) sample collection and handling methodologies. Samples selected for laboratory analysis will be submitted to a NYSDOH Environmental Laboratory Accreditation Program (ELAP) Contract Laboratory Protocol (CLP) certified laboratory, with a Category B deliverables package. Additionally, a Data Usability Summary Report (DUSR) will be prepared by a third-party data validator.

Sampling data will be used to evaluate remedial alternatives to meet the objectives identified above. Two data confidence levels will be considered, including field screening data and analytical level data. Field screening will include photoionization detector (PID), groundwater elevation measurement, and field groundwater characteristic measurements (pH, temperature, specific conductivity, turbidity). Analytical level data will be associated with select soil and groundwater samples submitted for chemical analysis to an independent laboratory.

2.2 **Project Organization**

HEI/WGS will establish a project team for successful completion of the project. The project team has not been finalized and subcontractors will be determined. Once the team has been finalized, appropriate resumes and information will be provided to NYSDEC. The anticipated project team is listed below:





Company	Name	Role						
Octavus Storage NY, LLC	Joe McNeil	Applicant and Property Owner						
Wittman GeoSciences, PLLC	Michele Wittman, PE	Project Director						
Hazard Evaluations	Mark Hanna, CHMM	Environmental Health & Safety Manager						
Carmina, Wood, Morris, DPC	R. Christopher Wood, PE	Project Engineer						
Hazard Evaluations	Eric Betzold	Project Scientist/Site Safety Officer						
Alpha Analytical	Candace Fox	Analytical Laboratory Subcontractor						
TBD	TBD	Test pit/monitoring well						
Data Validation Services	Judy Harry	Data Usability Subcontractor						

Michele Wittman – Michele will be the Project Director for the work and will be responsible for completion of each task, including coordination and supervision of field activities, adherence to work plan, schedule and budget. Additionally, Ms. Wittman will be responsible for development of the work plan, coordination of subcontractors, field project oversight and report preparations.

3.0 INVESTIGATION SCOPE OF WORK

3.1 Introduction

The proposed RI scope of work will include investigation for potential site contaminants in the soil/fill, former industrial process pond, and groundwater at the Site. The scope of work includes thirty-five (35) test pits, installation of five (5) monitoring wells, eight (8) surface soil samples, and location of one (1) surface water sample and three (3) sediment samples from the pond area. Proposed sampling locations are included on Figure 5 and summary of proposed analytical testing is presented on Table 1.

3.2 Building Survey

A building survey will be completed to identify any chemicals or wastes present on-site. An inventory of identified chemicals and/or wastes will be recorded at the time of the building survey.

3.3 Field Investigation Activities

Prior to intrusive activities, HEI/WGS and appropriate subcontractors will contact Dig Safely New York a minimum of three business days prior to the commencement of the field work. Investigative procedures are described below:

3.3.1 Surface Soil Investigation





Eight (8) surface soil samples will be collected from distributed locations as shown on Figure 5. Surface soil samples will be collected from 0 to 2 inches below grade. If topsoil is present, the topsoil will be removed, and the immediate 2 inches will be sampled. The sampled depth interval will be screened for VOCs using a calibrated organic vapor meter equipped with a photoionization detector (PID). Organic vapor meter results and soil descriptions will be recorded on the field soil boring logs.

Surface soil samples will be selected for analysis for the following as shown on Table 1 and briefly summarized below:

- 8 samples for Target Compound List (TCL) VOCs
- 8 samples for TCL semi-volatile organic compounds (SVOCs)
- 8 samples for Target Analyte List (TAL) metals
- 4 samples for Polychlorinated bi-phenyls (PCBs)
- 4 samples for pesticides and herbicides
- 4 samples for 1,4-dioxane
- 4 samples for per- and polyfluoroalkyl substances (PFAS)

Actual sample locations will be selected in the field based on utility locations, field observations, screening results, and engineering judgement. Subsurface soil samples will be collected using dedicated stainless-steel sampling tools. Select representative soil samples will be placed in pre-cleaned laboratory-provide sample bottles, labeled and cooled in 4°C in the field, and transported under chain-of-custody to a NYSDEC ELAP certified analytical laboratory.

3.3.2 Industrial Pond Samples

A former industrial scrubber discharge retention pond is located in the northern portion of the site. The pond is man-made and was constructed strictly for industrial use purposes and was not meant to support aquatic life. Further, during site development plans, the pond water will be removed and the pond backfilled with NYSDEC-approved fill material, and finished with either a storage building or parking area. As the pond was not designed to support aquatic life, the proposed surface water and sediment sampling to be done is to assess potential impacts within the sediment material, prior to pond backfilling activities, and will be used to assess disposal options for the pond water and soil to be removed. The following samples will be collected from the pond:

- One **surface water sample** will be collected from the pond, and will be analyzed for VOCs, SVOCs, total Metals, PCBs, pesticides/herbicides, 1,4-dioxane, and PFSA. As the pond is one continuous waster feature, only one surface water sample is needed.
- Three (3) **sediment samples** will be collected from depths of 0 to 6 inches, within the pond. One sediment sample will be taken adjacent to the inflow discharge pipe, one sediment sample will be taken adjacent to the outflow discharge pipe, and one sediment sample will be taken from a





low point within the pond. The proposed sediment sample locations are included in Figure 5. The three sediment samples will be analyzed for the following parameters, as shown on Table 1.

- TCL VOCs
- TCL SVOCs
- TAL metals
- PCBs
- pesticides and herbicides
- 1.4-dioxane
- PFAS

3.3.3 Test Pit Excavations

Test pit excavations will be completed throughout the site, but will focus on the incinerator waste disposal area in the northern wooded area. Thirty five (35) test pits will be completed to further investigate the vertical and horizontal extent of the incinerator, as well as general site soil conditions. The proposed locations are included on Figure 5. Test pits will be completed with a tracked excavator capable of reaching a minimum of 15 feet below grade. The depth of the test pit will extend 2 to 3 feet into underlying native soils.

In accordance with NYSDEC DER-10 Technical Guidance to Site Investigation and Remediation (DER-10), a Community Air Monitoring Plan (CAMP) must be implemented during all ground intrusion activities. A CAMP will be completed to achieve real-time monitoring of volatile organic compounds (VOCs) and particulates (dust). A copy of the CAMP is provided in Attachment A of the Health & Safety Plan (HASP), included in Appendix C.

Soil samples will be selected for analytical analysis based on in field screening results, visual and olfactory observations. Soil samples will be taken from different test pit and soil boring locations, distributed throughout the site. HEI will collect representative samples from each of the identified fill types, as well as the underlying native soils, for appropriate laboratory analysis.

3.3.4 Subsurface Soil Investigation

Seven (7) soil borings will be completed throughout the site, as shown on Figure 5. Five of the seven soil borings will be converted to monitoring well locations. The soil borings and monitoring well locations will be advanced using a drill rig or direct push drill rig capable of advancing hollow-stem augers for installing 2-inch monitoring wells. Monitoring wells are expected to be completed to greater depths of up to 16 to 20 feet below grade.

Discrete subsurface soil samples will be field screened in approximate two-foot depth intervals for VOCs with a calibrated organic vapor meter equipped with a photoionization detector (PID). Organic vapor meter results and soil descriptions will be recorded on the field soil boring logs.





Soil samples will be selected for analysis based in field screening results, visual and olfactory observations. During initial investigations, incinerator waste was present in the northern portion of the site. Remaining soil conditions are currently not known. WGS/HEI will collect representative samples from each of the identified fill types, as well as the underlying native soils, for geologic classification. Currently, WGS/HEI does not plan on further soil sample analysis from soil boring locations.

3.3.5 Subsurface Soil Sample Laboratory Analysis

The sample interval from test pits and/or soil borings identified as the most impacted (i.e., highest PID reading, visual/olfactory evidence of odors, staining, or product) will be selected for analysis. Should fill material be encountered, a discrete sample will be collected from each type of fill soil. In the event that visual and olfactory evidence of impact is not encountered in the soil samples, a representative sample will be collected from the native soils directly below the fill/native interface for analysis. Additionally, attempt will be made to collected soil samples at vertical variations within the native soil.

Subsurface soil samples will be selected for analysis for the following as shown on Table 1 and briefly summarized below:

- 10 samples for Target Compound List (TCL) VOCs
- 10 samples for TCL semi-volatile organic compounds (SVOCs)
- 15 samples for Target Analyte List (TAL) metals
- 4 samples for polychlorinated bi-phenyls (PCBs)
- 4 samples for pesticides and herbicides
- 4 samples for 1,4-dioxane
- 4 samples for per- and polyfluoroalkyl substances (PFAS)

Actual sample numbers and locations will be selected in the field based on utility locations, field observations, screening results, and engineering judgment. Additional soil samples will be selected, as needed, based on observed subsurface conditions. Subsurface soil samples will be collected using dedicated stainless steel sampling tools. Select representative soil samples will be placed in pre-cleaned laboratory-provided sample bottles, labeled and cooled to 4°C in the field, and transported under chain-of-custody to a NYSDOH ELAP certified analytical laboratory.

3.3.6 Monitoring Well Installation

Five (5) of the seven (7) soil boring locations will be converted to monitoring wells using a drill rig capable of advancing a borehole, using a direct push method via Geoprobe drill rig, to allow for installation of 2-inch diameter wells. The wells will be utilized for measurement of groundwater depth and collection of groundwater samples. The five proposed locations are included on Figure 5.





After completion of the soil borings to depths expected to range from 16 to 20 feet below grade, a 2-inch diameter, schedule 40 PVC monitoring well will be installed at each location. Groundwater samples will be evaluated as per NYSDEC DER-10 Technical Guidance for Site Investigation and Remediation (DER-10) and, if necessary, soil borings/wells will be advanced to the required depths of groundwater. An approximate 10 foot length of 0.010-inch machine slotted well screen will be installed at each location attached to the riser. The well screen depth will be backfilled with silica sand filter pack (estimated at size #0) from the base to approximately 2 feet above the well screen. A bentonite seal will be placed above the sand and hydrated to limit potential for down-hole contamination. The top of the well riser will be flush with the ground surface and completed with a locking J-plug. The well will be finished with protective casing.

Groundwater samples will be collected from each of the monitoring wells using low flow sampling techniques. The total depth of the wells is expected to range from 16 to 20 feet below grade.

3.3.7 Monitoring Well Development

After a minimum of 24-hours from installation, the monitoring wells will be developed using dedicated disposable polyethylene bailers via purge methodology. Field parameters, including pH, temperature, turbidity, and specific conductance will be measured periodically until they become relatively stable (approximately 10% fluctuation or less). A minimum of three well volumes will be removed from each monitoring well, unless dry well conditions are encountered. Development water will be containerized in 55-gallon drums and sampled for future off-site disposal.

3.3.8 Groundwater Sampling

Prior to sample collection, static groundwater levels will be measured at each of the monitoring wells. The wells will be purged and field measurements of pH, specific conductivity, temperature and turbidity will be recorded and monitored for stabilization prior to sampling. Groundwater samples will be collected using low flow sampling techniques. If insufficient groundwater, new dedicated disposable bailers may be used to collect the groundwater samples. Purge water will be containerized in 55-gallon drums and sampled for future off-site disposal.

The five (5) groundwater samples will be analyzed for the following parameters as summarized on Table 1:

- Target Compound List (TCL) VOCs
- TCL semi-volatile organic compounds (SVOCs)
- Target Analyte List (TAL) metals (dissolved phase only)
- Polychlorinated bi-phenyls (PCBs)
- Pesticides and herbicides
- 1,4-dioxane (3 samples only)
- PFAS (3 samples only)





Groundwater samples will be placed in pre-cleaned laboratory-provided sample bottles, labeled and preserved in accordance with USEPA SW-846 methodology, and transported under chain-of-custody to a NYSDOH ELAP certified analytical laboratory.

Per- and polyfluoroalkyl substances (PFAS) Sampling Protocol

Groundwater sample collection procedure for PFAS will be done in accordance with NYSDEC protocol, as included in the Quality Assurance/Quality Control (QAPP) found in Appendix B. At this time acceptable materials for sampling include stainless steel, high density polyethylene (HDPE), PVC, silicone, acetate and polypropylene. Equipment blanks should be generated at least daily. Additional materials may be acceptable if pre-approved by NYSDEC. Requests to use alternate equipment should include clean equipment blanks. NOTE: Grunfos pumps and bladder pumps are known to contain PFC materials (e.g. Teflon™ washers for Grunfos pumps and LDPE bladders for bladder pumps). All sampling equipment components and sample containers should not come in contact with aluminum foil, low density polyethylene (LDPE), glass or polytetrafluoroethylene (PTFE, Teflon™) materials including sample bottle cap liners with a PTFE layer. Standard two step decontamination using detergent and clean water rinse will be performed for equipment that does come in contact with PFC materials. Clothing that contains PTFE material (including GORE-TEX®) or that have been waterproofed with PFC materials must be avoided. Many food and drink packaging materials and "plumbers thread seal tape" contain PFCs.

All clothing worn by sampling personnel must have been laundered multiple times. The sampler must wear nitrile gloves while filling and sealing the sample bottles.

Pre-cleaned sample bottles with closures, coolers, ice, sample labels and a chain of custody form will be provided by the laboratory.

- 1. Fill two pre-cleaned 500 mL HDPE or polypropylene bottle with the sample.
- 2. Cap the bottles with an acceptable cap and liner closure system.
- 3. Label the sample bottles.
- 4. Fill out the chain of custody.
- 5. Place in a cooler maintained at 4 ± 2° Celsius.

Groundwater samples will be analyzed for PFAS via EPA method 537, to achieve reporting limits of 2 ng/l (parts per trillion (ppt)).

3.3.9 Field Specific Quality Assurance/Quality Control Sampling

Field-specific quality assurance/quality control samples will be collected and analyzed, as summarized on Table 1 to support third-party data usability assessment effort. Site-specific QA/QC samples will include blind duplicate, matrix spike/matrix spike duplicate, rinsate blank, and trip blank.





3.4 <u>Investigation- Derived Waste Management</u>

During the completion of soil borings, removed materials will be placed into the borehole. The excess soil cuttings that cannot be replaced into the borehole will placed next to the soil boring location. Based on analytical testing results, the excess soil may be utilized on-site, or disposed off-site. Development/purge water generated during well development and/or sampling activities will be containerized in 55-gallon drums for testing and future off-site disposal.

3.5 Site Mapping

A base map will be prepared by a New York State-licensed surveyor. The surveyor will establish the horizontal location and vertical elevations. The map will include the RI investigation/sampling locations. Soil/fill boring locations will be field located and incorporated within the survey. Elevations of the ground surface and top of PVC riser will be measured for each monitoring well.

3.6 Personnel Decontamination

The degree of decontamination is a function of both the particular task and the physical environment in which it takes place. Decontamination procedures will remain flexible, thereby allowing the decontamination crew to respond appropriately to changing conditions at the Site. On-site sampling activities will be carried out in such a manner as to avoid gross contamination of site workers, personal protective equipment, machinery and equipment.

Between sampling locations (or sometimes between samples at one sampling location), and upon the completion of the daily field activities, site workers will proceed to the Contaminated Reduction Zone (CRZ) or mobile reduction zone area. Equipment (e.g., sampling tubes, shovels, tools, etc.) will be decontaminated in this area. Prior to leaving the Site for breaks, at the end of the work shift, or when PPE has been grossly contaminated, disposable boot covers, gloves, and suits will be removed and placed in a drum designated for the disposal of these materials. After removing PPE, each Site worker will wash with soap and fresh water prior to donning new PPE or leaving the Site for the day. All wash water and rinse water will be collected and disposed of in accordance with appropriate regulations.

3.7 Decontamination of Equipment

Equipment decontamination efforts will be conducted in the CRZ or mobile reduction zone areas. Gross contamination will first be removed with plastic scrapers or other appropriate tools. The equipment will be decontaminated at a temporary equipment decontamination pad in the CRZ via hand washing or pressure washing. Downhole tools and augers can be hand washed or pressure washed.

The decontamination of the direct push drilling rig, excavator, or other heavy equipment will be undertaken as necessary. Initially, scraping of the equipment will remove heavily caked materials prior to washing. Washing will then be accomplished by pressure washing. Water generated during decontamination activities may be





allowed to discharge in the ground at the test pit location or may be collected, stored and profiled for future off-site disposal.

3.8 <u>Disposal of Contaminated Materials</u>

Potentially contaminated materials (gloves, clothing, sample sleeves etc.) will be bagged and segregated for proper disposal. Investigation-derived waste will be managed in accordance with NYSDEC guidance regulations. All fluids collected during groundwater sampling and decontamination will be containerized and managed appropriately subsequent to field activities and decontamination procedures.

4.0 REMEDIAL INVESTIGATIONALTERNATAIVES ANALYSIS REPORT

Upon completion of the RI tasks, a RI report will be generated in general requirements as identified in DER-10 Section 3.14. The report will include the following information.

- Background and Site information.
- Description of investigation areas.
- Identify and characterize the sources of contamination.
- Comparison with cleanup levels during the alternatives analysis report (AAR).
- Describe the amount, concentration, environmental fate and transport (if necessary), location and other significant characteristics of the contaminants present.
- Define hydraulic factors, as needed.
- Provide a qualitative human exposure assessment.
- Identify actual or potential adverse impacts to fish and wildlife resources

An independent data validation expert will complete a third-party data view of the analytical data generated during the RI work. A Data Usability Summary Report (DUSR) will be prepared, with appropriate data qualifiers added to the results.

The RI report will also include an alternatives analysis report to evaluate a remedial approach. The AAR will evaluate the need for further remedial activities.

Remedial action objectives will be evaluated and developed to assure the selected remedy is protective of human health and the environment under the proposed future Site usage. Proposed soil cleanup objectives will be based on proposed future usage. Should further remedial requirements be identified, a list of potentially applicable remedial technologies will be developed and evaluated. Criteria to be evaluated for the remedy and protectiveness to public health and the environment include:

- Overall protection of the public health and the environment
- Standards, criteria and guidance (SCG)
- Long-term effectiveness and permanence
- Reduction of toxicity, mobility or volume of contamination through treatment
- Short-term impact and effectiveness





- Implementability
- Cost effectiveness
- Land use

The results of the AAR will identify a remedial alternative to be recommended for the Site, which will include a discussion on the reasons for the selection. Community acceptance and comments will be evaluated within the alternative selection.

5.0 ADDITIONAL PROJECT DOCUMENTS

Various supporting documents have been prepared associated with the RI work plan and included in the appendix as listed below.

5.1 **Quality Assurance Project Plan**

The Quality Assurance Project Plan (QAPP) was generated in general accordance with Section 2.4 in DER-10. The QAPP describes the quality assurance/quality control (QA/QC) protocols and guidance associated with the RI Work Plan to ensure the suitability and verifiable data result from the sampling and analysis. The QAPP also provides procedures to be used during sampling of various media, field activities, and analytical laboratory testing. The QAPP is included in Appendix B.

5.2 **Health and Safety Plan**

A Site-specific Health and Safety Plan (HASP) has been prepared for this project and included in Appendix C. The HASP will be enforced by HEI, and will apply to all Site visitors and subcontractors associated with the RI field activities. The HASP covers the on-Site investigation and interim remedial work. Subcontractors will be required to develop and implement their health and safety plan which reflects the requirements of this document.

The HASP will include a Community Air Monitoring Plan (CAMP) to describe particulate and volatile organic vapor monitoring to protect nearby community during the investigative activities.


6.0 PROJECT SCHEDULE

Figure 6 presents the tentative schedule for planned activities. A certificate of completion (COC) is anticipated by December 2019.






FIGURES



THIS DRAWING IS FOR ILLUSTRATIVE AND INFORMATIONAL PURPOSES ONLY AND WAS ADAPTED FROM USGS, BUFFALO SE, NEW YORK 2016 QUADRANGLE.

HAZARD EVALUATIONS, INC. Phase I/II Audits – Site Investigations – Facility Inspections SITE LOCATION 3678 SOUTHWESTERN BOULEVARD ORCHARD PARK, NEW YORK OCTAVUS STORAGE NY, LLC TAMPA, FLORIDA DRAWN BY: SS SCALE: NOT TO SCALE PROJECT: e1728 CHECKED BY: MW DATE: 09/18 FIGURE NO: 1




HAZARD EVALUATIONS, INC.

Phase I/II Audits – Site Investigations – Facility Inspections

SITE LIMITS

3678 SOUTHWESTERN BOULEVARD ORCHARD PARK, NEW YORK

DRAWN BY: SS	SCALE: not to scale	PROJECT: e1728
CHECKED BY: MMW	DATE: 09/18	FIGURE NO: 2



KEY

Soil Boring Location (04/2018)

Surface Soil Location (04/2018)

exceeds Unrestriced Use SCO
exceeds Restricted Residential Use SCO
exceeds Commercial Use SCO


HAZARD EVALUATIONS, INC.


Phase I/II Audits – Site Investigations – Facility Inspections

SOIL BORING LOCATION PLAN

3678 SOUTHWESTERN BOULEVARD ORCHARD PARK, NEW YORK

DRAWN BY: SS	SCALE: not to scale	PROJECT: e1/28
CHECKED BY: MMW	DATE: 09/18	FIGURE NO: 3








Figure 7 BCP Project Schedule Octavus Storage NY, LLC

3678 Southwestern Boulevard, Orchard Park, NY

	2	2018															20	19												
(F) 1	November	De	cember		Jan	uary	February	I	March		Apı	ril		May	J	une		Jul	y		August	S	eptembe	er	0	ctober		November	D	December
Task	5 12 19 26	3 10	17 24	31	7 14	21 28	4 11 18 2	5 4 1	11 18 2	5 1	8 15	5 22 2	29 6	5 13 20 27	3 10	17 24	1	8 15	5 22 2	29 5	12 19 26	2 9	9 16 2	3 30	7 1	4 21	28	4 11 18 25	2 9	16 23 3
Application and Work Plan																														
Submittal of BCP App and RI/AAR Work Plan																														
NYSDEC Review for Completeness and Revisions																														
Application Deemed complete																														
30 day public notice for Application																														
BCP Agreement																														
NYSDEC Review of RI/AAR Work Plan																														
30 day public notice for RI/AAR Work Plan																														
RI/AAR Comments and revisions																														
Acceptance of BCP and Work Plan approval										Ü																				
Remedial Investigation																														
Test Pits and Soil Borings																														
Groundwater Sampling																														
Analaytical Testing																														
DUSR Preparation																														
RI/AAR Report																														
NYSDEC Review and Decision Document																														
Remedial Action Work Plan																														
Complete and Submitt RAWP																														
NYSDEC Review of RAWP																														
45 day public notice																														
RAWP Approval																														
Remedial Work and Site Development																														
Remedial and Construction Activities																														
Reporting																														
Final Engineering Report																														
NYSDEC Review																														
Certificate of Completion													Ī																	

Milestone Date
Task by HEI/WGS
NYSEC Review
Public Comment
Laboratory analysis/DUSR by Subcontractor
Owner/Applicant

TABLES

TABLE 1 Analytical Testing Program Summary Octavus Storage NY, LLC

3678 Southwester Boulevard, Orchard Park, NY NYSDEC Brownfield Cleanup Program

Location	Number of Proposed Locations	Matrix	TCL VOCS	TCL SVOCs	TAL METALS Total	TAL METALS dissolved	PCBs	Pest/ Herbs	1,4-dioxane	PFAS
Industrial Pond Sedim	ent Samples									
Sediment	3	Soil	3	3	3		3	3	3	3
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1	-	1	1	1	1
Total			7	7	7	0	7	7	7	7
Industrial Pond Surfa	ce Water Samp	les								
Water	1	Water	1	1	1	1	1	1	1	1
Rinsate		Water	1	1	1	1	1	1	1	1
Total			2	2	2	2	2	2	2	2
Test Pit Locations - Su	bsurface Samp	les								
Test Pits	17	Soil	10	10	15	-	4	4	4	4
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1	-	1	1	1	1
Total			14	14	19	0	8	8	8	8
Surface Soil Samples										
Soil	5	Soil	8	8	8	-	4	4	4	4
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1		1	1	1	1
Total			12	12	12	0	8	8	8	8
Monitoring Wells										
Monitoring Well	5	Groundwater	5	5	5	5	5	5	3	3
Duplicate		Groundwater	1	1	1	1	1	1	1	1
MS/MSD		Groundwater	2	2	2	2	2	2	2	2
Rinsate		Water	1	1	1	1	1	1	1	1
Trip Blank		Water	1	-	-	-	-	-	-	-
Total			10	9	9	9	9	9	7	7
	то	TAL SAMPLES	VOCs 45	SVOCs 44	METALS 49	METALS 11	PCBs 34	Pest / Herbs 34	1,4-dioxane 32	PFAS 32

Notes:

TCL VOCs - Target Compound List Volatile Organic Compounds. TCL SVOCs - Target Compound List Semi-volatile Organic Compounds.

TAL Metals - Target Analyte List Metals.
TCL PCBs - Target Compound List Polychlorinated Biphenyls.

PFAS - Polyfluoroalkyl Substances

APPENDIX A HISTORICAL INVESTIGATION INFORMATION

PHASE II ENVIRONMENTAL SITE ASSESSMENT

- for -

Octavus Storage Joe McNeil 1810 W. Kennedy Blvd Tampa, FL 33606

- location -

3678 Southwestern Boulevard Orchard Park, New York 14127

April 2018

Prepared By:

GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS, INC.

50 Ridge Road
Buffalo, New York 14218
(716) 827-0700

<u>Tabl</u>	e of Co	<u>ontents</u>	<u>Page</u>
1.0	Intro	duction	1
	1.1	General	1
	1.2	Project Understanding	1
	1.3	Site Investigation Objectives	1
	1.4	Report Organization	1
2.0	Site 1	Investigation	2
	2.1	General	2
	2.2	Soil Borings	2
	2.3	PID Analysis	2
3.0	Anal	3	
	3.1	Subsurface Soil & Groundwater	3
4.0	Conc	4	
	4.1	Conclusions	4
	4.2	Recommendations	4

Appendices

Appendix A – Site Map Appendix B – Analytical Data Report

Λ 4•	1
COCTION	
Section	
~ • • • • • • • • • • • • • • • • • • •	

GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS

Introduction

1.0 Introduction

1.1 General

This report has been prepared by Great Lakes Environmental & Safety Consultants, Inc. ("Great Lakes") to document the findings of a Phase II Environmental Site Assessment ("ESA") performed on April 6, 2018 for a parcel located at 3678 Southwestern Boulevard in Orchard Park, New York (the "Site"). The Site is comprised of a degraded building, paved areas, overgrown wooded areas, and an industrial process pond, occupying approximately 7.33 acres. The Site was previously a waste burning demonstration project, process waste materials are found at the North end of the site in a blackened field.

This report presents and summarizes the methodology and findings of the ESA for the site. The scope of work associated with the ESA is based on discussions with site representatives as well as review of documents related to previous assessments conducted at the site.

1.2 Project Understanding

Based on the review and evaluation of background information pertaining to the Site and neighboring properties, strategically selected areas were investigated and evaluated during the ESA.

1.3 Site Investigation Objectives

The purpose and objectives of the ESA include the following:

- Determine the presence of contaminants in the soil and the extent if any;
- Identify potential contamination source areas and migration pathways; and
- Develop and evaluate potential remedial measures, as necessary.

1.4 Report Organization

This report presents the findings from data obtained during the Phase II ESA. Section 2 discusses the investigation activities that were performed at the Site. Section 3 provides an overview of the analytical data obtained during the Phase II ESA. Section 4 presents our conclusions regarding the interpretation and findings of the data obtained during the Phase II ESA.

GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS

Site Investigation

2.0 Site Investigation

2.1 General

This section presents the approach and methodology used in performing the Phase II ESA at the Site. To meet the objectives of the project, various field activities were conducted at the Site on April 6, 2018 that included the following:

- Soil borings;
- Monitoring of all recovered soil from borings with photo ionization detector ("PID") for elevated levels of volatile organic chemicals ("VOCs");
- Sampling and analysis of boring soils; and
- Sampling and analysis of surface samples of the blackened field area to the North of the Site.

The following subsections briefly describe the implementation of the above noted field activities.

2.2 Soil Borings

On April 6, 2018 soil borings were performed at the Site with oversight by a Great Lakes environmental professional. The borings were selected in a manner to most effectively test for intrusion of contaminants from the previous industrial activities. A site map indicating the boring locations is presented in Appendix A.

The borings were advanced to a goal depth of twelve (12) feet below grade. Variations in Site geology caused variations in the actual final depth of the borings. The borings were performed using a drill rig equipped with a hydraulic powered, vibratory press.

2.3 PID Analysis

During drilling activities, the exposed soil was visually inspected for evidence of contamination and screened for total volatile organic vapors using a field calibrated Industrial Scientific IBRID MX6 PID. Each boring location was composited and sent in for analysis.

Section 3 Anaytical Results

GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS

3.0 Analytical Results

This section presents the results of the environmental media sampling performed as part of the Phase II ESA.

3.1 Subsurface Soil & Ground Water

All recovered soil from the Site was screened for total volatile organic compounds (VOCs) using an Industrial Scientific IBRID MX6 PID. Headspace VOCs were collected in the field, any elevated readings would indicate a requirement for further analysis. These readings were taken from all borings, at four-foot intervals.

The boring soils and the selected surface soil samples were analyzed for the STARS volatile organic compounds and STARS semi-volatile organic compounds, polychlorinated biphenyls, and total metals in accordance with USEPA Testing Methods. The complete analytical data report is presented in Appendix C.



GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS

4.0 Conclusions and Recommendations

This section summarizes the interpretation of the field data and associated findings obtained during the Phase II ESA.

4.1 Conclusions

The major findings associated with the ESA conducted at the study area include the following:

Soils


Soil borings were conducted, and soil samples were visually inspected for evidence of contamination and screened for total volatile organic vapors using a field calibrated Industrial Scientific IBRID MX6 PID. All the boring locations (B-1, B-2, B-3) and the surface samples (SS-1, SS-2, SS-3, SS-4, SS-5) were submitted for laboratory analysis. The analytical results were then compared with the standards set forth by the New York State Department of Environmental Conservation ("NYSDEC") in 6 NYCRR Subpart 375-6 – Remedial Program Soil Cleanup Objectives ("Subpart 375-6"). Neither B-1, B-2, and B-3 composite boring soil samples returned any levels of VOC's, SVOC's, PCBs, or Metals that exceeded NYSDEC's cleanup objectives when compared to Table 375-6.8(a) Unrestricted Use Soil Cleanup Objectives. The surface samples however did have exceedances for Commercial usages when comparing to the cleanup objectives for the metals Barium and Copper.

4.2 Recommendations

Based on the information obtained during the Phase II ESA, and through consulting Great Lakes contacts within the New York State Department of Environmental Conservation, further investigation of the northern blackened area and of the industrial settling pond is required. Commercial development could not proceed on this site without remediation, so investigation of these areas is mandatory to assess the impacts of this dumped waste on the site and adjacent areas. Once impacts to this northern area could be assessed, remedial activities could proceed to clean up the site to the commercial cleanup objectives. Depending on the extent of this metal contamination, this site may be accepted into New York State's Brownfield tax credits. Further investigation is required, and Great Lakes would be pleased to provide a quote to investigate the extent of the metal contamination on the site.

Appendix A Site Map

GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS





GREAT LAKES ENVIRONMENTAL & SAFETY CONSULTANTS



ANALYTICAL REPORT

Lab Number: L1812163

Client: Great Lake Environmental & Safety Cons.

50 Ridge Road Buffalo, NY 14218

ATTN: Mark Mol

Phone: (716) 827-0700

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified Report Date: 04/16/18

The original project report/data package is held by Alpha Analytical. This report/data package is paginated and should be reproduced only in its entirety. Alpha Analytical holds no responsibility for results and/or data that are not consistent with the original.

Certifications & Approvals: MA (M-MA086), NH NELAP (2064), NJ NELAP (MA935), CT (PH-0574), IL (200077), ME (MA00086), MD (348), NY (11148), NC (25700/666), PA (68-03671), RI (LAO00065), TX (T104704476), VT (VT-0935), VA (460195), USDA (Permit #P330-14-00197).

Eight Walkup Drive, Westborough, MA 01581-1019 508-898-9220 (Fax) 508-898-9193 800-624-9220 - www.alphalab.com



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163 **Report Date:** 04/16/18

Alpha Sample ID	Client ID	Matrix	Sample Location	Collection Date/Time	Receive Date
L1812163-01	SS-1	SOIL	ORCHARD PARK, NY	04/06/18 08:17	04/09/18
L1812163-02	SS-2	SOIL	ORCHARD PARK, NY	04/06/18 08:20	04/09/18
L1812163-03	SS-3	SOIL	ORCHARD PARK, NY	04/06/18 08:22	04/09/18
L1812163-04	SS-4	SOIL	ORCHARD PARK, NY	04/06/18 08:25	04/09/18
L1812163-05	SS-5	SOIL	ORCHARD PARK, NY	04/06/18 08:29	04/09/18
L1812163-06	B-1	SOIL	ORCHARD PARK, NY	04/06/18 09:30	04/09/18
L1812163-07	B-2	SOIL	ORCHARD PARK, NY	04/06/18 09:50	04/09/18
L1812163-08	B-3	SOIL	ORCHARD PARK, NY	04/06/18 10:10	04/09/18



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified **Report Date:** 04/16/18

Case Narrative

The samples were received in accordance with the Chain of Custody and no significant deviations were encountered during the preparation or analysis unless otherwise noted. Sample Receipt, Container Information, and the Chain of Custody are located at the back of the report.

Results contained within this report relate only to the samples submitted under this Alpha Lab Number and meet NELAP requirements for all NELAP accredited parameters unless otherwise noted in the following narrative. The data presented in this report is organized by parameter (i.e. VOC, SVOC, etc.). Sample specific Quality Control data (i.e. Surrogate Spike Recovery) is reported at the end of the target analyte list for each individual sample, followed by the Laboratory Batch Quality Control at the end of each parameter. Tentatively Identified Compounds (TICs), if requested, are reported for compounds identified to be present and are not part of the method/program Target Compound List, even if only a subset of the TCL are being reported. If a sample was re-analyzed or re-extracted due to a required quality control corrective action and if both sets of data are reported, the Laboratory ID of the re-analysis or re-extraction is designated with an "R" or "RE", respectively. When multiple Batch Quality Control elements are reported (e.g. more than one LCS), the associated samples for each element are noted in the grey shaded header line of each data table. Any Laboratory Batch, Sample Specific % recovery or RPD value that is outside the listed Acceptance Criteria is bolded in the report. All specific QC information is also incorporated in the Data Usability format of our Data Merger tool where it can be reviewed along with any associated usability implications. Soil/sediments, solids and tissues are reported on a dry weight basis unless otherwise noted. Definitions of all data qualifiers and acronyms used in this report are provided in the Glossary located at the back of the report.

In reference to questions H (CAM) or 4 (RCP) when "NO" is checked, the performance criteria for CAM and RCP methods allow for some quality control failures to occur and still be within method compliance. In these instances the specific failure is not narrated but noted in the associated QC table. The information is also incorporated in the Data Usability format of our Data Merger tool where it can be reviewed along with any associated usability implications.

Please see the associated ADEx data file for a comparison of laboratory reporting limits that were achieved with the regulatory Numerical Standards requested on the Chain of Custody.

HOLD POLICY

For samples submitted on hold, Alpha's policy is to hold samples (with the exception of Air canisters) free of charge for 21 calendar days from the date the project is completed. After 21 calendar days, we will dispose of all samples submitted including those put on hold unless you have contacted your Client Service Representative and made arrangements for Alpha to continue to hold the samples. Air canisters will be disposed after 3 business days from the date the project is completed.

Please contact Client Services at 800-624-9220 with any questions.



Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

Case Narrative (continued)

Report Submission

All non-detect (ND) or estimated concentrations (J-qualified) have been quantitated to the limit noted in the MDL column.

Volatile Organics

Any reported concentrations that are below 200 ug/kg may be low due to the sample not being collected according to 5035-L/5035A-L low-level specifications.

Total Metals

L1812163-01 through -08: The sample has elevated detection limits for all elements, with the exception of mercury, due to the dilution required by matrix interferences encountered during analysis.

I, the undersigned, attest under the pains and penalties of perjury that, to the best of my knowledge and belief and based upon my personal inquiry of those responsible for providing the information contained in this analytical report, such information is accurate and complete. This certificate of analysis is not complete unless this page accompanies any and all pages of this report.

- Amita Naik

Authorized Signature:

Title: Technical Director/Representative Date: 04/16/18

Nails

ORGANICS



VOLATILES



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1

Sample Location: ORCHARD PARK, NY Date Received: Field Prep:

Lab Number:

Report Date:

04/09/18 Not Specified

L1812163

04/16/18

Sample Depth:

Matrix: Soil 1,8260C Analytical Method: Analytical Date: 04/11/18 12:31

JC Analyst: 96% Percent Solids:

Volatile Organics by GC/MS - Westborough Lab Methylene chloride ND ug/kg 5.2 1,1-Dichloroethane ND ug/kg 0.78 Chloroform ND ug/kg 0.78 Carbon tetrachloride ND ug/kg 0.52 1,2-Dichloropropane ND ug/kg 1.8 Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 0.52 1,2-Dichloroethane ND ug/kg 0.52 Bromodichloromethane ND ug/kg 0.52	0.86 0.14	1
1,1-Dichloroethane ND ug/kg 0.78 Chloroform ND ug/kg 0.78 Carbon tetrachloride ND ug/kg 0.52 1,2-Dichloropropane ND ug/kg 1.8 Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 0.52 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52		1
Chloroform ND ug/kg 0.78 Carbon tetrachloride ND ug/kg 0.52 1,2-Dichloropropane ND ug/kg 1.8 Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 0.52 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.14	
Carbon tetrachloride ND ug/kg 0.52 1,2-Dichloropropane ND ug/kg 1.8 Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52		1
1,2-Dichloropropane ND ug/kg 1.8 Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.19	1
Dibromochloromethane ND ug/kg 0.52 1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.18	1
1,1,2-Trichloroethane ND ug/kg 0.78 Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.12	1
Tetrachloroethene ND ug/kg 0.52 Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.09	1
Chlorobenzene ND ug/kg 0.52 Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.16	1
Trichlorofluoromethane ND ug/kg 2.6 1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.16	1
1,2-Dichloroethane ND ug/kg 0.52 1,1,1-Trichloroethane ND ug/kg 0.52	0.18	1
1,1,1-Trichloroethane ND ug/kg 0.52	0.22	1
ÿ ÿ	0.13	1
Promodichloromethons ND value 0.52	0.18	1
Bromodichloromethane ND ug/kg 0.52	0.16	1
trans-1,3-Dichloropropene ND ug/kg 0.52	0.11	1
cis-1,3-Dichloropropene ND ug/kg 0.52	0.12	1
1,3-Dichloropropene, Total ND ug/kg 0.52	0.11	1
Bromoform ND ug/kg 2.1	0.12	1
1,1,2,2-Tetrachloroethane ND ug/kg 0.52	0.16	1
Benzene ND ug/kg 0.52	0.10	1
Toluene ND ug/kg 0.78	0.10	1
Ethylbenzene ND ug/kg 0.52	0.09	1
Chloromethane ND ug/kg 2.6	0.23	1
Bromomethane ND ug/kg 1.0	0.18	1
Vinyl chloride ND ug/kg 1.0	0.16	1
Chloroethane ND ug/kg 1.0	0.16	1
1,1-Dichloroethene ND ug/kg 0.52	0.19	1
trans-1,2-Dichloroethene ND ug/kg 0.78		
Trichloroethene ND ug/kg 0.52	0.13	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	tborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	2.6	0.10	1
1,3-Dichlorobenzene	ND		ug/kg	2.6	0.11	1
1,4-Dichlorobenzene	ND		ug/kg	2.6	0.10	1
Methyl tert butyl ether	ND		ug/kg	1.0	0.08	1
p/m-Xylene	ND		ug/kg	1.0	0.18	1
o-Xylene	ND		ug/kg	1.0	0.18	1
Xylenes, Total	ND		ug/kg	1.0	0.18	1
cis-1,2-Dichloroethene	ND		ug/kg	0.52	0.18	1
1,2-Dichloroethene, Total	ND		ug/kg	0.52	0.13	1
Styrene	ND		ug/kg	1.0	0.21	1
Dichlorodifluoromethane	ND		ug/kg	5.2	0.26	1
Acetone	3.2	J	ug/kg	5.2	1.2	1
Carbon disulfide	ND		ug/kg	5.2	0.58	1
2-Butanone	ND		ug/kg	5.2	0.36	1
4-Methyl-2-pentanone	ND		ug/kg	5.2	0.13	1
2-Hexanone	ND		ug/kg	5.2	0.35	1
Bromochloromethane	ND		ug/kg	2.6	0.19	1
1,2-Dibromoethane	ND		ug/kg	2.1	0.10	1
n-Butylbenzene	ND		ug/kg	0.52	0.12	1
sec-Butylbenzene	ND		ug/kg	0.52	0.11	1
tert-Butylbenzene	ND		ug/kg	2.6	0.13	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	2.6	0.21	1
Isopropylbenzene	ND		ug/kg	0.52	0.10	1
p-Isopropyltoluene	ND		ug/kg	0.52	0.10	1
Naphthalene	ND		ug/kg	2.6	0.07	1
n-Propylbenzene	ND		ug/kg	0.52	0.11	1
1,2,3-Trichlorobenzene	ND		ug/kg	2.6	0.13	1
1,2,4-Trichlorobenzene	ND		ug/kg	2.6	0.11	1
1,3,5-Trimethylbenzene	ND		ug/kg	2.6	0.08	1
1,2,4-Trimethylbenzene	ND		ug/kg	2.6	0.10	1
Methyl Acetate	ND		ug/kg	10	0.24	1
Cyclohexane	ND		ug/kg	10	0.23	1
1,4-Dioxane	ND		ug/kg	21	7.5	1
Freon-113	ND		ug/kg	10	0.27	1
Methyl cyclohexane	ND		ug/kg	2.1	0.12	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	101	70-130	
Toluene-d8	97	70-130	
4-Bromofluorobenzene	105	70-130	
Dibromofluoromethane	99	70-130	



L1812163

04/09/18

Not Specified

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Date Collected:

Report Date: 04/16/18

Lab Number:

Date Received:

Field Prep:

SAMPLE RESULTS

Lab ID: 04/06/18 08:20 L1812163-02

Client ID: SS-2

Sample Location: ORCHARD PARK, NY

Sample Depth:

Matrix: Soil Analytical Method: 1,8260C Analytical Date: 04/11/18 12:58

JC Analyst: 96% Percent Solids:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Volatile Organics by GC/MS - West	borough Lab						
Methylene chloride	ND		ug/kg	8.8	1.4	1	
1,1-Dichloroethane	ND		ug/kg	1.3	0.24	1	
Chloroform	ND		ug/kg	1.3	0.33	1	
Carbon tetrachloride	ND		ug/kg	0.88	0.30	1	
1,2-Dichloropropane	ND		ug/kg	3.1	0.20	1	
Dibromochloromethane	ND		ug/kg	0.88	0.16	1	
1,1,2-Trichloroethane	ND		ug/kg	1.3	0.28	1	
Tetrachloroethene	ND		ug/kg	0.88	0.27	1	
Chlorobenzene	ND		ug/kg	0.88	0.31	1	
Trichlorofluoromethane	ND		ug/kg	4.4	0.37	1	
1,2-Dichloroethane	ND		ug/kg	0.88	0.22	1	
1,1,1-Trichloroethane	ND		ug/kg	0.88	0.31	1	
Bromodichloromethane	ND		ug/kg	0.88	0.27	1	
trans-1,3-Dichloropropene	ND		ug/kg	0.88	0.18	1	
cis-1,3-Dichloropropene	ND		ug/kg	0.88	0.20	1	
1,3-Dichloropropene, Total	ND		ug/kg	0.88	0.18	1	
Bromoform	ND		ug/kg	3.5	0.21	1	
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.88	0.26	1	
Benzene	ND		ug/kg	0.88	0.17	1	
Toluene	ND		ug/kg	1.3	0.17	1	
Ethylbenzene	ND		ug/kg	0.88	0.15	1	
Chloromethane	ND		ug/kg	4.4	0.38	1	
Bromomethane	ND		ug/kg	1.8	0.30	1	
Vinyl chloride	ND		ug/kg	1.8	0.28	1	
Chloroethane	ND		ug/kg	1.8	0.28	1	
1,1-Dichloroethene	ND		ug/kg	0.88	0.33	1	
trans-1,2-Dichloroethene	ND		ug/kg	1.3	0.21	1	
Trichloroethene	ND		ug/kg	0.88	0.27	1	



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	stborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	4.4	0.16	1
1,3-Dichlorobenzene	ND		ug/kg	4.4	0.19	1
1,4-Dichlorobenzene	ND		ug/kg	4.4	0.16	1
Methyl tert butyl ether	ND		ug/kg	1.8	0.14	1
p/m-Xylene	ND		ug/kg	1.8	0.31	1
o-Xylene	ND		ug/kg	1.8	0.30	1
Xylenes, Total	ND		ug/kg	1.8	0.30	1
cis-1,2-Dichloroethene	ND		ug/kg	0.88	0.30	1
1,2-Dichloroethene, Total	ND		ug/kg	0.88	0.21	1
Styrene	ND		ug/kg	1.8	0.35	1
Dichlorodifluoromethane	ND		ug/kg	8.8	0.44	1
Acetone	ND		ug/kg	8.8	2.0	1
Carbon disulfide	ND		ug/kg	8.8	0.97	1
2-Butanone	ND		ug/kg	8.8	0.61	1
4-Methyl-2-pentanone	ND		ug/kg	8.8	0.22	1
2-Hexanone	ND		ug/kg	8.8	0.59	1
Bromochloromethane	ND		ug/kg	4.4	0.32	1
1,2-Dibromoethane	ND		ug/kg	3.5	0.18	1
n-Butylbenzene	ND		ug/kg	0.88	0.20	1
sec-Butylbenzene	ND		ug/kg	0.88	0.19	1
tert-Butylbenzene	ND		ug/kg	4.4	0.22	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.4	0.35	1
Isopropylbenzene	ND		ug/kg	0.88	0.17	1
p-Isopropyltoluene	ND		ug/kg	0.88	0.18	1
Naphthalene	ND		ug/kg	4.4	0.12	1
n-Propylbenzene	ND		ug/kg	0.88	0.19	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.4	0.22	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.4	0.19	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.4	0.14	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.4	0.16	1
Methyl Acetate	ND		ug/kg	18	0.41	1
Cyclohexane	ND		ug/kg	18	0.38	1
1,4-Dioxane	ND		ug/kg	35	13.	1
Freon-113	ND		ug/kg	18	0.45	1
Methyl cyclohexane	ND		ug/kg	3.5	0.21	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	100	70-130	
Toluene-d8	97	70-130	
4-Bromofluorobenzene	103	70-130	
Dibromofluoromethane	97	70-130	



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:22

Lab ID: L1812163-03

Client ID: SS-3

Sample Location: ORCHARD PARK, NY

Sample Depth:

Matrix: Soil
Analytical Method: 1,8260C
Analytical Date: 04/11/18 13:25

Analyst: PK Percent Solids: 93%

Date Collected:	04/06/18 08:22
Date Received:	04/09/18
Field Pren:	Not Specified

Lab Number:

Report Date:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	stborough Lab					
Methylene chloride	ND		ug/kg	9.1	1.5	1
1,1-Dichloroethane	ND		ug/kg	1.4	0.24	1
Chloroform	ND		ug/kg	1.4	0.34	1
Carbon tetrachloride	ND		ug/kg	0.91	0.31	1
1,2-Dichloropropane	ND		ug/kg	3.2	0.21	1
Dibromochloromethane	ND		ug/kg	0.91	0.16	1
1,1,2-Trichloroethane	ND		ug/kg	1.4	0.28	1
Tetrachloroethene	ND		ug/kg	0.91	0.27	1
Chlorobenzene	ND		ug/kg	0.91	0.32	1
Trichlorofluoromethane	ND		ug/kg	4.5	0.38	1
1,2-Dichloroethane	ND		ug/kg	0.91	0.22	1
1,1,1-Trichloroethane	ND		ug/kg	0.91	0.32	1
Bromodichloromethane	ND		ug/kg	0.91	0.28	1
trans-1,3-Dichloropropene	ND		ug/kg	0.91	0.19	1
cis-1,3-Dichloropropene	ND		ug/kg	0.91	0.21	1
1,3-Dichloropropene, Total	ND		ug/kg	0.91	0.19	1
Bromoform	ND		ug/kg	3.6	0.22	1
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.91	0.27	1
Benzene	ND		ug/kg	0.91	0.18	1
Toluene	ND		ug/kg	1.4	0.18	1
Ethylbenzene	ND		ug/kg	0.91	0.15	1
Chloromethane	ND		ug/kg	4.5	0.40	1
Bromomethane	ND		ug/kg	1.8	0.31	1
Vinyl chloride	ND		ug/kg	1.8	0.29	1
Chloroethane	ND		ug/kg	1.8	0.29	1
1,1-Dichloroethene	ND		ug/kg	0.91	0.34	1
trans-1,2-Dichloroethene	ND		ug/kg	1.4	0.22	1
Trichloroethene	ND		ug/kg	0.91	0.27	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03 Date Collected: 04/06/18 08:22

Client ID: SS-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	stborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	4.5	0.16	1
1,3-Dichlorobenzene	ND		ug/kg	4.5	0.20	1
1,4-Dichlorobenzene	ND		ug/kg	4.5	0.16	1
Methyl tert butyl ether	ND		ug/kg	1.8	0.14	1
p/m-Xylene	ND		ug/kg	1.8	0.32	1
o-Xylene	ND		ug/kg	1.8	0.31	1
Xylenes, Total	ND		ug/kg	1.8	0.31	1
cis-1,2-Dichloroethene	ND		ug/kg	0.91	0.31	1
1,2-Dichloroethene, Total	ND		ug/kg	0.91	0.22	1
Styrene	ND		ug/kg	1.8	0.36	1
Dichlorodifluoromethane	ND		ug/kg	9.1	0.45	1
Acetone	ND		ug/kg	9.1	2.1	1
Carbon disulfide	ND		ug/kg	9.1	1.0	1
2-Butanone	ND		ug/kg	9.1	0.63	1
4-Methyl-2-pentanone	ND		ug/kg	9.1	0.22	1
2-Hexanone	ND		ug/kg	9.1	0.60	1
Bromochloromethane	ND		ug/kg	4.5	0.32	1
1,2-Dibromoethane	ND		ug/kg	3.6	0.18	1
n-Butylbenzene	ND		ug/kg	0.91	0.21	1
sec-Butylbenzene	ND		ug/kg	0.91	0.20	1
tert-Butylbenzene	ND		ug/kg	4.5	0.22	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.5	0.36	1
Isopropylbenzene	ND		ug/kg	0.91	0.18	1
p-Isopropyltoluene	ND		ug/kg	0.91	0.18	1
Naphthalene	ND		ug/kg	4.5	0.12	1
n-Propylbenzene	ND		ug/kg	0.91	0.20	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.5	0.23	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.5	0.20	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.5	0.15	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.5	0.17	1
Methyl Acetate	ND		ug/kg	18	0.42	1
Cyclohexane	ND		ug/kg	18	0.39	1
1,4-Dioxane	ND		ug/kg	36	13.	1
Freon-113	ND		ug/kg	18	0.47	1
Methyl cyclohexane	ND		ug/kg	3.6	0.22	1



04/06/18 08:22

Date Collected:

Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

L1812163-03

Client ID: SS-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Lab ID:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Acceptance Qualifier Criteria
1,2-Dichloroethane-d4	99	70-130
Toluene-d8	96	70-130
4-Bromofluorobenzene	105	70-130
Dibromofluoromethane	97	70-130



L1812163

04/16/18

Not Specified

Dilution Factor

04/09/18

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Result

Date Collected: 04/06/18 08:25

Lab Number:

Report Date:

Date Received:

Field Prep:

RL

MDL

Lab ID: L1812163-04

Client ID: SS-4

Sample Location: ORCHARD PARK, NY

Sample Depth:

Parameter

Matrix: Soil
Analytical Method: 1,8260C
Analytical Date: 04/11/18 13:52

Analyst: PK Percent Solids: 96%

raiaillelei	Nesuit	Qualifiei	Ullita	NL.	IVIDE	Dilution i actor	
Volatile Organics by GC/MS - Wes	stborough Lab						
Methylene chloride	ND		ug/kg	8.6	1.4	1	
1,1-Dichloroethane	ND		ug/kg	1.3	0.23	1	
Chloroform	ND		ug/kg	1.3	0.32	1	
Carbon tetrachloride	ND		ug/kg	0.86	0.30	1	
1,2-Dichloropropane	ND		ug/kg	3.0	0.20	1	
Dibromochloromethane	ND		ug/kg	0.86	0.15	1	
1,1,2-Trichloroethane	ND		ug/kg	1.3	0.27	1	
Tetrachloroethene	ND		ug/kg	0.86	0.26	1	
Chlorobenzene	ND		ug/kg	0.86	0.30	1	
Trichlorofluoromethane	ND		ug/kg	4.3	0.36	1	
1,2-Dichloroethane	ND		ug/kg	0.86	0.21	1	
1,1,1-Trichloroethane	ND		ug/kg	0.86	0.30	1	
Bromodichloromethane	ND		ug/kg	0.86	0.26	1	
trans-1,3-Dichloropropene	ND		ug/kg	0.86	0.18	1	
cis-1,3-Dichloropropene	ND		ug/kg	0.86	0.20	1	
1,3-Dichloropropene, Total	ND		ug/kg	0.86	0.18	1	
Bromoform	ND		ug/kg	3.4	0.20	1	
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.86	0.26	1	
Benzene	ND		ug/kg	0.86	0.16	1	
Toluene	ND		ug/kg	1.3	0.17	1	
Ethylbenzene	ND		ug/kg	0.86	0.14	1	
Chloromethane	ND		ug/kg	4.3	0.37	1	
Bromomethane	ND		ug/kg	1.7	0.29	1	
Vinyl chloride	ND		ug/kg	1.7	0.27	1	
Chloroethane	ND		ug/kg	1.7	0.27	1	
1,1-Dichloroethene	ND		ug/kg	0.86	0.32	1	
trans-1,2-Dichloroethene	ND		ug/kg	1.3	0.21	1	
Trichloroethene	ND		ug/kg	0.86	0.26	1	

Qualifier

Units



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04 Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	stborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	4.3	0.16	1
1,3-Dichlorobenzene	ND		ug/kg	4.3	0.19	1
1,4-Dichlorobenzene	ND		ug/kg	4.3	0.16	1
Methyl tert butyl ether	ND		ug/kg	1.7	0.13	1
p/m-Xylene	ND		ug/kg	1.7	0.30	1
o-Xylene	ND		ug/kg	1.7	0.29	1
Xylenes, Total	ND		ug/kg	1.7	0.29	1
cis-1,2-Dichloroethene	ND		ug/kg	0.86	0.29	1
1,2-Dichloroethene, Total	ND		ug/kg	0.86	0.21	1
Styrene	ND		ug/kg	1.7	0.34	1
Dichlorodifluoromethane	ND		ug/kg	8.6	0.43	1
Acetone	ND		ug/kg	8.6	2.0	1
Carbon disulfide	ND		ug/kg	8.6	0.94	1
2-Butanone	ND		ug/kg	8.6	0.59	1
4-Methyl-2-pentanone	ND		ug/kg	8.6	0.21	1
2-Hexanone	ND		ug/kg	8.6	0.57	1
Bromochloromethane	ND		ug/kg	4.3	0.31	1
1,2-Dibromoethane	ND		ug/kg	3.4	0.17	1
n-Butylbenzene	ND		ug/kg	0.86	0.20	1
sec-Butylbenzene	ND		ug/kg	0.86	0.19	1
tert-Butylbenzene	ND		ug/kg	4.3	0.21	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.3	0.34	1
Isopropylbenzene	ND		ug/kg	0.86	0.17	1
p-lsopropyltoluene	ND		ug/kg	0.86	0.17	1
Naphthalene	ND		ug/kg	4.3	0.12	1
n-Propylbenzene	ND		ug/kg	0.86	0.18	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.3	0.22	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.3	0.18	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.3	0.14	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.3	0.16	1
Methyl Acetate	ND		ug/kg	17	0.40	1
Cyclohexane	ND		ug/kg	17	0.37	1
1,4-Dioxane	ND		ug/kg	34	12.	1
Freon-113	ND		ug/kg	17	0.44	1
Methyl cyclohexane	ND		ug/kg	3.4	0.20	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04 Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	99	70-130	
Toluene-d8	97	70-130	
4-Bromofluorobenzene	105	70-130	
Dibromofluoromethane	98	70-130	



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

L1812163-05

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:29

Lab Number:

Report Date:

SAMI LE NEGOLIO

SS-5 Date Received: 04/09/18
ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Sample Location:

Lab ID:

Client ID:

Matrix: Soil
Analytical Method: 1,8260C
Analytical Date: 04/11/18 14:18

Analyst: PK Percent Solids: 96%

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Volatile Organics by GC/MS - We	stborough Lab						
Methylene chloride	ND		ug/kg	9.2	1.5	1	
1,1-Dichloroethane	ND		ug/kg	1.4	0.25	1	
Chloroform	ND		ug/kg	1.4	0.34	1	
Carbon tetrachloride	ND		ug/kg	0.92	0.32	1	
1,2-Dichloropropane	ND		ug/kg	3.2	0.21	1	
Dibromochloromethane	ND		ug/kg	0.92	0.16	1	
1,1,2-Trichloroethane	ND		ug/kg	1.4	0.29	1	
Tetrachloroethene	ND		ug/kg	0.92	0.28	1	
Chlorobenzene	ND		ug/kg	0.92	0.32	1	
Trichlorofluoromethane	ND		ug/kg	4.6	0.38	1	
1,2-Dichloroethane	ND		ug/kg	0.92	0.22	1	
1,1,1-Trichloroethane	ND		ug/kg	0.92	0.32	1	
Bromodichloromethane	ND		ug/kg	0.92	0.28	1	
trans-1,3-Dichloropropene	ND		ug/kg	0.92	0.19	1	
cis-1,3-Dichloropropene	ND		ug/kg	0.92	0.21	1	
1,3-Dichloropropene, Total	ND		ug/kg	0.92	0.19	1	
Bromoform	ND		ug/kg	3.7	0.22	1	
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.92	0.27	1	
Benzene	ND		ug/kg	0.92	0.18	1	
Toluene	ND		ug/kg	1.4	0.18	1	
Ethylbenzene	ND		ug/kg	0.92	0.16	1	
Chloromethane	ND		ug/kg	4.6	0.40	1	
Bromomethane	ND		ug/kg	1.8	0.31	1	
Vinyl chloride	ND		ug/kg	1.8	0.29	1	
Chloroethane	ND		ug/kg	1.8	0.29	1	
1,1-Dichloroethene	ND		ug/kg	0.92	0.34	1	
trans-1,2-Dichloroethene	ND		ug/kg	1.4	0.22	1	
Trichloroethene	ND		ug/kg	0.92	0.28	1	



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05 Date Collected: 04/06/18 08:29

Client ID: SS-5 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	stborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	4.6	0.17	1
1,3-Dichlorobenzene	ND		ug/kg	4.6	0.20	1
1,4-Dichlorobenzene	ND		ug/kg	4.6	0.17	1
Methyl tert butyl ether	ND		ug/kg	1.8	0.14	1
p/m-Xylene	ND		ug/kg	1.8	0.32	1
o-Xylene	ND		ug/kg	1.8	0.31	1
Xylenes, Total	ND		ug/kg	1.8	0.31	1
cis-1,2-Dichloroethene	ND		ug/kg	0.92	0.31	1
1,2-Dichloroethene, Total	ND		ug/kg	0.92	0.22	1
Styrene	ND		ug/kg	1.8	0.37	1
Dichlorodifluoromethane	ND		ug/kg	9.2	0.46	1
Acetone	ND		ug/kg	9.2	2.1	1
Carbon disulfide	ND		ug/kg	9.2	1.0	1
2-Butanone	ND		ug/kg	9.2	0.63	1
4-Methyl-2-pentanone	ND		ug/kg	9.2	0.22	1
2-Hexanone	ND		ug/kg	9.2	0.61	1
Bromochloromethane	ND		ug/kg	4.6	0.33	1
1,2-Dibromoethane	ND		ug/kg	3.7	0.18	1
n-Butylbenzene	ND		ug/kg	0.92	0.21	1
sec-Butylbenzene	ND		ug/kg	0.92	0.20	1
tert-Butylbenzene	ND		ug/kg	4.6	0.23	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.6	0.36	1
Isopropylbenzene	ND		ug/kg	0.92	0.18	1
p-Isopropyltoluene	ND		ug/kg	0.92	0.18	1
Naphthalene	ND		ug/kg	4.6	0.13	1
n-Propylbenzene	ND		ug/kg	0.92	0.20	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.6	0.23	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.6	0.20	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.6	0.15	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.6	0.17	1
Methyl Acetate	ND		ug/kg	18	0.42	1
Cyclohexane	ND		ug/kg	18	0.40	1
1,4-Dioxane	ND		ug/kg	37	13.	1
Freon-113	ND		ug/kg	18	0.47	1
Methyl cyclohexane	ND		ug/kg	3.7	0.22	1



Project Name: Lab Number: SOUTHWESTERN PHASE 2 L1812163

Project Number: Report Date: Not Specified 04/16/18

SAMPLE RESULTS

Lab ID: Date Collected: 04/06/18 08:29 L1812163-05

Date Received: Client ID: SS-5 04/09/18 Sample Location: Field Prep: ORCHARD PARK, NY Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL **Dilution Factor**

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	100	70-130	
Toluene-d8	97	70-130	
4-Bromofluorobenzene	105	70-130	
Dibromofluoromethane	98	70-130	



L1812163

Project Name: SOUTHWESTERN PHASE 2

ORCHARD PARK, NY

L1812163-06

B-1

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 09:30

Report Date: 04/16/18

Lab Number:

Date Received: 04/09/18
Field Prep: Not Specified

Sample Depth:

Sample Location:

Lab ID:

Client ID:

Matrix: Soil
Analytical Method: 1,8260C
Analytical Date: 04/11/18 14:45

Analyst: PK Percent Solids: 93%

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Volatile Organics by GC/MS - Wes	tborough Lab						
Methylene chloride	ND		ug/kg	8.0	1.3	1	
1,1-Dichloroethane	ND		ug/kg	1.2	0.22	1	
Chloroform	ND		ug/kg	1.2	0.30	1	
Carbon tetrachloride	ND		ug/kg	0.80	0.28	1	
1,2-Dichloropropane	ND		ug/kg	2.8	0.18	1	
Dibromochloromethane	ND		ug/kg	0.80	0.14	1	
1,1,2-Trichloroethane	ND		ug/kg	1.2	0.25	1	
Tetrachloroethene	ND		ug/kg	0.80	0.24	1	
Chlorobenzene	ND		ug/kg	0.80	0.28	1	
Trichlorofluoromethane	ND		ug/kg	4.0	0.33	1	
1,2-Dichloroethane	ND		ug/kg	0.80	0.20	1	
1,1,1-Trichloroethane	ND		ug/kg	0.80	0.28	1	
Bromodichloromethane	ND		ug/kg	0.80	0.25	1	
trans-1,3-Dichloropropene	ND		ug/kg	0.80	0.17	1	
cis-1,3-Dichloropropene	ND		ug/kg	0.80	0.18	1	
1,3-Dichloropropene, Total	ND		ug/kg	0.80	0.17	1	
Bromoform	ND		ug/kg	3.2	0.19	1	
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.80	0.24	1	
Benzene	ND		ug/kg	0.80	0.15	1	
Toluene	ND		ug/kg	1.2	0.16	1	
Ethylbenzene	ND		ug/kg	0.80	0.14	1	
Chloromethane	ND		ug/kg	4.0	0.35	1	
Bromomethane	ND		ug/kg	1.6	0.27	1	
Vinyl chloride	ND		ug/kg	1.6	0.25	1	
Chloroethane	ND		ug/kg	1.6	0.25	1	
1,1-Dichloroethene	ND		ug/kg	0.80	0.30	1	
trans-1,2-Dichloroethene	ND		ug/kg	1.2	0.19	1	
Trichloroethene	ND		ug/kg	0.80	0.24	1	

Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Wes	tborough Lab					
1,2-Dichlorobenzene	ND		ug/kg	4.0	0.15	1
1,3-Dichlorobenzene	ND		ug/kg	4.0	0.17	1
1,4-Dichlorobenzene	ND		ug/kg	4.0	0.15	1
Methyl tert butyl ether	ND		ug/kg	1.6	0.12	1
p/m-Xylene	ND		ug/kg	1.6	0.28	1
o-Xylene	ND		ug/kg	1.6	0.27	1
Xylenes, Total	ND		ug/kg	1.6	0.27	1
cis-1,2-Dichloroethene	ND		ug/kg	0.80	0.27	1
1,2-Dichloroethene, Total	ND		ug/kg	0.80	0.19	1
Styrene	ND		ug/kg	1.6	0.32	1
Dichlorodifluoromethane	ND		ug/kg	8.0	0.40	1
Acetone	ND		ug/kg	8.0	1.8	1
Carbon disulfide	1.1	J	ug/kg	8.0	0.88	1
2-Butanone	ND		ug/kg	8.0	0.55	1
4-Methyl-2-pentanone	ND		ug/kg	8.0	0.20	1
2-Hexanone	ND		ug/kg	8.0	0.53	1
Bromochloromethane	ND		ug/kg	4.0	0.29	1
1,2-Dibromoethane	ND		ug/kg	3.2	0.16	1
n-Butylbenzene	ND		ug/kg	0.80	0.18	1
sec-Butylbenzene	ND		ug/kg	0.80	0.17	1
tert-Butylbenzene	ND		ug/kg	4.0	0.20	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.0	0.32	1
Isopropylbenzene	ND		ug/kg	0.80	0.16	1
p-Isopropyltoluene	ND		ug/kg	0.80	0.16	1
Naphthalene	ND		ug/kg	4.0	0.11	1
n-Propylbenzene	ND		ug/kg	0.80	0.17	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.0	0.20	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.0	0.17	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.0	0.13	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.0	0.15	1
Methyl Acetate	ND		ug/kg	16	0.37	1
Cyclohexane	ND		ug/kg	16	0.35	1
1,4-Dioxane	ND		ug/kg	32	12.	1
Freon-113	ND		ug/kg	16	0.41	1
Methyl cyclohexane	0.21	J	ug/kg	3.2	0.19	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: B-1 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	99	70-130	
Toluene-d8	99	70-130	
4-Bromofluorobenzene	116	70-130	
Dibromofluoromethane	98	70-130	



L1812163

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Report Date: 04/16/18

Lab Number:

Lab ID: L1812163-07

Client ID: B-2

Sample Location: ORCHARD PARK, NY

Sample Depth:

Matrix: Soil 1,8260C Analytical Method: Analytical Date: 04/11/18 15:12

Analyst: PΚ 87% Percent Solids:

Date Collected:	04/06/18 09:50
Date Received:	04/09/18
Field Prep:	Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - We	estborough Lab					
Methylene chloride	ND		ug/kg	8.3	1.4	1
1,1-Dichloroethane	ND		ug/kg	1.2	0.22	1
Chloroform	ND		ug/kg	1.2	0.31	1
Carbon tetrachloride	ND		ug/kg	0.83	0.29	1
1,2-Dichloropropane	ND		ug/kg	2.9	0.19	1
Dibromochloromethane	ND		ug/kg	0.83	0.15	1
1,1,2-Trichloroethane	ND		ug/kg	1.2	0.26	1
Tetrachloroethene	ND		ug/kg	0.83	0.25	1
Chlorobenzene	ND		ug/kg	0.83	0.29	1
Trichlorofluoromethane	ND		ug/kg	4.2	0.35	1
1,2-Dichloroethane	ND		ug/kg	0.83	0.20	1
1,1,1-Trichloroethane	ND		ug/kg	0.83	0.29	1
Bromodichloromethane	ND		ug/kg	0.83	0.26	1
trans-1,3-Dichloropropene	ND		ug/kg	0.83	0.17	1
cis-1,3-Dichloropropene	ND		ug/kg	0.83	0.19	1
1,3-Dichloropropene, Total	ND		ug/kg	0.83	0.17	1
Bromoform	ND		ug/kg	3.3	0.20	1
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.83	0.25	1
Benzene	ND		ug/kg	0.83	0.16	1
Toluene	ND		ug/kg	1.2	0.16	1
Ethylbenzene	ND		ug/kg	0.83	0.14	1
Chloromethane	ND		ug/kg	4.2	0.36	1
Bromomethane	ND		ug/kg	1.7	0.28	1
Vinyl chloride	ND		ug/kg	1.7	0.26	1
Chloroethane	ND		ug/kg	1.7	0.26	1
1,1-Dichloroethene	ND		ug/kg	0.83	0.31	1
trans-1,2-Dichloroethene	ND		ug/kg	1.2	0.20	1
Trichloroethene	ND		ug/kg	0.83	0.25	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07 Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Volatile Organics by GC/MS - Westborough Lab 1,2-Dichlorobenzene 1,3-Dichlorobenzene 1,4-Dichlorobenzene Methyl tert butyl ether p/m-Xylene o-Xylene Xylenes, Total	ND ND ND ND ND ND ND ND ND	ug/kg ug/kg ug/kg ug/kg ug/kg	4.2 4.2 4.2 1.7	0.15 0.18 0.15 0.13	1 1 1
1,3-Dichlorobenzene 1,4-Dichlorobenzene Methyl tert butyl ether p/m-Xylene o-Xylene	ND ND ND ND ND ND ND ND	ug/kg ug/kg ug/kg	4.2 4.2 1.7	0.18 0.15	1
1,3-Dichlorobenzene 1,4-Dichlorobenzene Methyl tert butyl ether p/m-Xylene o-Xylene	ND ND ND ND ND ND ND ND	ug/kg ug/kg ug/kg	4.2 4.2 1.7	0.18 0.15	1
1,4-Dichlorobenzene Methyl tert butyl ether p/m-Xylene o-Xylene	ND ND ND ND	ug/kg ug/kg	4.2 1.7	0.15	1
Methyl tert butyl ether p/m-Xylene o-Xylene	ND ND ND	ug/kg	1.7		
p/m-Xylene o-Xylene	ND ND ND			0.13	1
o-Xylene	ND ND	ug/kg		0.29	1
·	ND	//	1.7	0.29	1
Aylenes, rotal		ug/kg	1.7	0.28	1
air 1.2 Diablaraethana	ND	ug/kg	0.83		1
cis-1,2-Dichloroethene	ND	ug/kg	0.83	0.28	1
1,2-Dichloroethene, Total		ug/kg			
Styrene	ND	ug/kg	1.7	0.33	1
Dichlorodifluoromethane	ND	ug/kg	8.3	0.42	1
Acetone	ND	ug/kg	8.3	1.9	1
Carbon disulfide	ND	ug/kg	8.3	0.92	<u> </u>
2-Butanone	ND	ug/kg	8.3	0.57	1
4-Methyl-2-pentanone	ND	ug/kg	8.3	0.20	1
2-Hexanone	ND	ug/kg	8.3	0.55	1
Bromochloromethane	ND	ug/kg	4.2	0.30	1
1,2-Dibromoethane	ND	ug/kg	3.3	0.16	1
n-Butylbenzene	ND	ug/kg	0.83	0.19	1
sec-Butylbenzene	ND	ug/kg	0.83	0.18	1
tert-Butylbenzene	ND	ug/kg	4.2	0.20	1
1,2-Dibromo-3-chloropropane	ND	ug/kg	4.2	0.33	1
Isopropylbenzene	ND	ug/kg	0.83	0.16	1
p-Isopropyltoluene	ND	ug/kg	0.83	0.17	1
Naphthalene	ND	ug/kg	4.2	0.11	1
n-Propylbenzene	ND	ug/kg	0.83	0.18	1
1,2,3-Trichlorobenzene	ND	ug/kg	4.2	0.21	1
1,2,4-Trichlorobenzene	ND	ug/kg	4.2	0.18	1
1,3,5-Trimethylbenzene	ND	ug/kg	4.2	0.13	1
1,2,4-Trimethylbenzene	ND	ug/kg	4.2	0.15	1
Methyl Acetate	ND	ug/kg	17	0.38	1
Cyclohexane	ND	ug/kg	17	0.36	1
1,4-Dioxane	ND	ug/kg	33	12.	1
Freon-113	ND	ug/kg	17	0.43	1
Methyl cyclohexane	ND	ug/kg	3.3	0.20	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07 Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL Dilution Factor

Surrogate	% Recovery	Qualifier	Acceptance Criteria	
1,2-Dichloroethane-d4	100		70-130	
Toluene-d8	98		70-130	
4-Bromofluorobenzene	108		70-130	
Dibromofluoromethane	99		70-130	



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 10:10

Lab Number:

Report Date:

L1812163-08 Date Received: 04/09/18 B-3 ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Sample Location:

Lab ID:

Client ID:

Matrix: Soil Analytical Method: 1,8260C Analytical Date: 04/11/18 15:39

Analyst: PΚ 91% Percent Solids:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Westborou	gh Lab					
Methylene chloride	ND		ug/kg	8.5	1.4	1
1,1-Dichloroethane	ND		ug/kg	1.3	0.23	1
Chloroform	ND		ug/kg	1.3	0.31	1
Carbon tetrachloride	ND		ug/kg	0.85	0.29	1
1,2-Dichloropropane	ND		ug/kg	3.0	0.19	1
Dibromochloromethane	ND		ug/kg	0.85	0.15	1
1,1,2-Trichloroethane	ND		ug/kg	1.3	0.26	1
Tetrachloroethene	ND		ug/kg	0.85	0.26	1
Chlorobenzene	ND		ug/kg	0.85	0.30	1
Trichlorofluoromethane	ND		ug/kg	4.2	0.35	1
1,2-Dichloroethane	ND		ug/kg	0.85	0.21	1
1,1,1-Trichloroethane	ND		ug/kg	0.85	0.30	1
Bromodichloromethane	ND		ug/kg	0.85	0.26	1
trans-1,3-Dichloropropene	ND		ug/kg	0.85	0.18	1
cis-1,3-Dichloropropene	ND		ug/kg	0.85	0.20	1
1,3-Dichloropropene, Total	ND		ug/kg	0.85	0.18	1
Bromoform	ND		ug/kg	3.4	0.20	1
1,1,2,2-Tetrachloroethane	ND		ug/kg	0.85	0.25	1
Benzene	ND		ug/kg	0.85	0.16	1
Toluene	ND		ug/kg	1.3	0.16	1
Ethylbenzene	ND		ug/kg	0.85	0.14	1
Chloromethane	ND		ug/kg	4.2	0.37	1
Bromomethane	ND		ug/kg	1.7	0.29	1
Vinyl chloride	ND		ug/kg	1.7	0.27	1
Chloroethane	ND		ug/kg	1.7	0.27	1
1,1-Dichloroethene	ND		ug/kg	0.85	0.32	1
trans-1,2-Dichloroethene	ND		ug/kg	1.3	0.20	1
Trichloroethene	ND		ug/kg	0.85	0.26	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08 Date Collected: 04/06/18 10:10

Client ID: B-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Volatile Organics by GC/MS - Westborough Lab						
1,2-Dichlorobenzene	ND		ug/kg	4.2	0.15	1
1,3-Dichlorobenzene	ND		ug/kg	4.2	0.18	1
1,4-Dichlorobenzene	ND		ug/kg	4.2	0.15	1
Methyl tert butyl ether	ND		ug/kg	1.7	0.13	1
p/m-Xylene	ND		ug/kg	1.7	0.30	1
o-Xylene	ND		ug/kg	1.7	0.29	1
Xylenes, Total	ND		ug/kg	1.7	0.29	1
cis-1,2-Dichloroethene	ND		ug/kg	0.85	0.29	1
1,2-Dichloroethene, Total	ND		ug/kg	0.85	0.20	1
Styrene	ND		ug/kg	1.7	0.34	1
Dichlorodifluoromethane	ND		ug/kg	8.5	0.42	1
Acetone	ND		ug/kg	8.5	1.9	1
Carbon disulfide	ND		ug/kg	8.5	0.93	1
2-Butanone	ND		ug/kg	8.5	0.58	1
4-Methyl-2-pentanone	ND		ug/kg	8.5	0.21	1
2-Hexanone	ND		ug/kg	8.5	0.56	1
Bromochloromethane	ND		ug/kg	4.2	0.30	1
1,2-Dibromoethane	ND		ug/kg	3.4	0.17	1
n-Butylbenzene	ND		ug/kg	0.85	0.19	1
sec-Butylbenzene	ND		ug/kg	0.85	0.18	1
tert-Butylbenzene	ND		ug/kg	4.2	0.21	1
1,2-Dibromo-3-chloropropane	ND		ug/kg	4.2	0.34	1
Isopropylbenzene	ND		ug/kg	0.85	0.16	1
p-Isopropyltoluene	ND		ug/kg	0.85	0.17	1
Naphthalene	ND		ug/kg	4.2	0.12	1
n-Propylbenzene	ND		ug/kg	0.85	0.18	1
1,2,3-Trichlorobenzene	ND		ug/kg	4.2	0.21	1
1,2,4-Trichlorobenzene	ND		ug/kg	4.2	0.18	1
1,3,5-Trimethylbenzene	ND		ug/kg	4.2	0.14	1
1,2,4-Trimethylbenzene	ND		ug/kg	4.2	0.16	1
Methyl Acetate	ND		ug/kg	17	0.39	1
Cyclohexane	ND		ug/kg	17	0.37	1
1,4-Dioxane	ND		ug/kg	34	12.	1
Freon-113	ND		ug/kg	17	0.44	1
Methyl cyclohexane	0.58	J	ug/kg	3.4	0.20	1



Project Name: Lab Number: SOUTHWESTERN PHASE 2 L1812163

Project Number: Report Date: Not Specified 04/16/18

SAMPLE RESULTS

Lab ID: Date Collected: 04/06/18 10:10 L1812163-08

Date Received: Client ID: 04/09/18 B-3 Sample Location: Field Prep: ORCHARD PARK, NY Not Specified

Sample Depth:

Parameter Result Qualifier Units RL MDL **Dilution Factor**

Surrogate	% Recovery	Acceptance Qualifier Criteria	
1,2-Dichloroethane-d4	100	70-130	
Toluene-d8	98	70-130	
4-Bromofluorobenzene	108	70-130	
Dibromofluoromethane	98	70-130	



Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: 1,8260C Analytical Date: 04/11/18 12:05

Analyst: MV

Methylene chloride ND ug/kg 10 1.6 1,1-Dichloroethane ND ug/kg 1.5 0.27 Chloroform ND ug/kg 1.5 0.37 Carbon tetrachloride ND ug/kg 1.0 0.34 1,2-Dichloropropane ND ug/kg 1.0 0.18 1,2-Trichloroethane ND ug/kg 1.0 0.18 1,1,2-Trichloroethane ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichloroffluoromethane ND ug/kg 1.0 0.35 Trichloroffluoromethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.35 Bromotichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total	Parameter	Result	Qualifier U	Units		RL	MDL	
1,1-Dichloroethane ND ug/kg 1.5 0.27 Chloroform ND ug/kg 1.5 0.37 Carbon tetrachloride ND ug/kg 1.0 0.34 1,2-Dichloropropane ND ug/kg 3.5 0.23 Dibromochloromethane ND ug/kg 1.0 0.18 1,1,2-Trichloroethane ND ug/kg 1.5 0.31 Tetrachloroethene ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.24 1,1,2,2-Tetrachloroethane ND<	Volatile Organics by GC/MS	- Westborough Lab	for sample(s): (01-08	Batch:	WG1105645-5	
Chloroform ND ug/kg 1.5 0.37 Carbon tetrachloride ND ug/kg 1.0 0.34 1,2-Dichloropropane ND ug/kg 3.5 0.23 Dibromochloromethane ND ug/kg 1.0 0.18 1,1,2-Trichloroethane ND ug/kg 1.0 0.31 Tetrachloroethene ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 1.0 0.25 1,2-Dichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.24 1,2,2-Tetrachloroethane N	Methylene chloride	ND	1	ug/kg		10	1.6	
Carbon tetrachloride ND ug/kg 1.0 0.34 1,2-Dichloropropane ND ug/kg 3.5 0.23 Dibromochloromethane ND ug/kg 1.0 0.18 1,1,2-Trichloroethane ND ug/kg 1.5 0.31 Tetrachloroethane ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.30 Benzene ND	1,1-Dichloroethane	ND	I	ug/kg		1.5	0.27	
1,2-Dichloropropane ND	Chloroform	ND		ug/kg		1.5	0.37	
Dibromochloromethane ND ug/kg 1.0 0.18 1,1,2-Trichloroethane ND ug/kg 1.5 0.31 Tetrachloroethene ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.5 0.20 Ethylbenzene ND </td <td>Carbon tetrachloride</td> <td>ND</td> <td></td> <td>ug/kg</td> <td></td> <td>1.0</td> <td>0.34</td> <td></td>	Carbon tetrachloride	ND		ug/kg		1.0	0.34	
1,1,2-Trichloroethane ND ug/kg 1.5 0.31 Tetrachloroethene ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.0 0.17 Chloromethane ND	1,2-Dichloropropane	ND		ug/kg		3.5	0.23	
Tetrachloroethene ND ug/kg 1.0 0.30 Chlorobenzene ND ug/kg 1.0 0.35 Trichloroftuoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 1.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg	Dibromochloromethane	ND		ug/kg		1.0	0.18	
Chlorobenzene ND ug/kg 1.0 0.35 Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.32 Chloroethane ND ug/kg	1,1,2-Trichloroethane	ND		ug/kg		1.5	0.31	
Trichlorofluoromethane ND ug/kg 5.0 0.42 1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene, Total ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 4.0 0.24 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/	Tetrachloroethene	ND		ug/kg		1.0	0.30	
1,2-Dichloroethane ND ug/kg 1.0 0.25 1,1,1-Trichloroethane ND ug/kg 1.0 0.35 Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 Ti-Dichloroethene ND ug/kg	Chlorobenzene	ND		ug/kg		1.0	0.35	
1,1,1-Trichloroethane	Trichlorofluoromethane	ND		ug/kg		5.0	0.42	
Bromodichloromethane ND ug/kg 1.0 0.31 trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg	1,2-Dichloroethane	ND		ug/kg		1.0	0.25	
trans-1,3-Dichloropropene ND ug/kg 1.0 0.21 cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.0 0.19 Ethylbenzene ND ug/kg 1.0 0.17 Chloroethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg <td< td=""><td>1,1,1-Trichloroethane</td><td>ND</td><td></td><td>ug/kg</td><td></td><td>1.0</td><td>0.35</td><td></td></td<>	1,1,1-Trichloroethane	ND		ug/kg		1.0	0.35	
cis-1,3-Dichloropropene ND ug/kg 1.0 0.23 1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Bromodichloromethane	ND		ug/kg		1.0	0.31	
1,3-Dichloropropene, Total ND ug/kg 1.0 0.21 Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	trans-1,3-Dichloropropene	ND		ug/kg		1.0	0.21	
Bromoform ND ug/kg 4.0 0.24 1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	cis-1,3-Dichloropropene	ND		ug/kg		1.0	0.23	
1,1,2,2-Tetrachloroethane ND ug/kg 1.0 0.30 Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	1,3-Dichloropropene, Total	ND		ug/kg		1.0	0.21	
Benzene ND ug/kg 1.0 0.19 Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Bromoform	ND		ug/kg		4.0	0.24	
Toluene ND ug/kg 1.5 0.20 Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	1,1,2,2-Tetrachloroethane	ND		ug/kg		1.0	0.30	
Ethylbenzene ND ug/kg 1.0 0.17 Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Benzene	ND	ı	ug/kg		1.0	0.19	
Chloromethane ND ug/kg 5.0 0.44 Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Toluene	ND		ug/kg		1.5	0.20	
Bromomethane ND ug/kg 2.0 0.34 Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Ethylbenzene	ND	l	ug/kg		1.0	0.17	
Vinyl chloride ND ug/kg 2.0 0.32 Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Chloromethane	ND	l	ug/kg		5.0	0.44	
Chloroethane ND ug/kg 2.0 0.32 1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Bromomethane	ND		ug/kg		2.0	0.34	
1,1-Dichloroethene ND ug/kg 1.0 0.37 trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Vinyl chloride	ND	ı	ug/kg		2.0	0.32	
trans-1,2-Dichloroethene ND ug/kg 1.5 0.24 Trichloroethene ND ug/kg 1.0 0.30	Chloroethane	ND		ug/kg		2.0	0.32	
Trichloroethene ND ug/kg 1.0 0.30	1,1-Dichloroethene	ND		ug/kg		1.0	0.37	
	trans-1,2-Dichloroethene	ND		ug/kg		1.5	0.24	
1,2-Dichlorobenzene ND ug/kg 5.0 0.18	Trichloroethene	ND		ug/kg		1.0	0.30	
	1,2-Dichlorobenzene	ND		ug/kg		5.0	0.18	



Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: 1,8260C Analytical Date: 04/11/18 12:05

Analyst: MV

/olatile Organics by GC/MS - V 1,3-Dichlorobenzene 1,4-Dichlorobenzene Methyl tert butyl ether	Nestborough Lab ND ND ND ND	ug/l	kg	5.0	WG1105645-5 0.22
1,4-Dichlorobenzene	ND ND	ug/			0.22
	ND		kg	- 0	
Methyl tert butyl ether		ua/	J	5.0	0.18
,	ND	ug/	kg	2.0	0.15
p/m-Xylene		ug/	kg	2.0	0.35
o-Xylene	ND	ug/	kg	2.0	0.34
Xylenes, Total	ND	ug/	kg	2.0	0.34
cis-1,2-Dichloroethene	ND	ug/	kg	1.0	0.34
1,2-Dichloroethene, Total	ND	ug/	kg	1.0	0.24
Styrene	ND	ug/	kg	2.0	0.40
Dichlorodifluoromethane	ND	ug/	kg	10	0.50
Acetone	ND	ug/	kg	10	2.3
Carbon disulfide	ND	ug/	kg	10	1.1
2-Butanone	ND	ug/	kg	10	0.69
4-Methyl-2-pentanone	ND	ug/	kg	10	0.24
2-Hexanone	ND	ug/	kg	10	0.67
Bromochloromethane	ND	ug/	kg	5.0	0.36
1,2-Dibromoethane	ND	ug/	kg	4.0	0.20
n-Butylbenzene	ND	ug/	kg	1.0	0.23
sec-Butylbenzene	ND	ug/	kg	1.0	0.22
tert-Butylbenzene	ND	ug/	kg	5.0	0.25
1,2-Dibromo-3-chloropropane	ND	ug/	kg	5.0	0.40
Isopropylbenzene	ND	ug/	kg	1.0	0.19
p-Isopropyltoluene	ND	ug/	kg	1.0	0.20
Naphthalene	ND	ug/	kg	5.0	0.14
n-Propylbenzene	ND	ug/	kg	1.0	0.22
1,2,3-Trichlorobenzene	ND	ug/	kg	5.0	0.25
1,2,4-Trichlorobenzene	ND	ug/	kg	5.0	0.22
1,3,5-Trimethylbenzene	ND	ug/	kg	5.0	0.16
1,2,4-Trimethylbenzene	ND	ug/	kg	5.0	0.19



Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: 1,8260C

Analytical Date: 04/11/18 12:05

Analyst: MV

Parameter	Result C	Qualifier Units	RL	MDL	
Volatile Organics by GC/MS -	Westborough Lab fo	or sample(s): 0	01-08 Batch:	WG1105645-5	
Methyl Acetate	ND	ug/kg	20	0.46	
Cyclohexane	ND	ug/kg	20	0.43	
1,4-Dioxane	ND	ug/kg	40	14.	
Freon-113	ND	ug/kg	20	0.51	
Methyl cyclohexane	ND	ug/kg	4.0	0.24	

	Acceptance					
Surrogate	%Recovery Qu	ualifier Criteria				
1,2-Dichloroethane-d4	100	70-130				
Toluene-d8	96	70-130				
4-Bromofluorobenzene	104	70-130				
Dibromofluoromethane	99	70-130				



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

arameter	LCS %Recovery	Qual	LCSD %Recovery	Qual	%Recovery Limits	RPD	Qual	RPD Limits	
olatile Organics by GC/MS - West	borough Lab Associated	sample(s):	01-08 Batch: V	VG1105645-3	WG1105645-4				
Methylene chloride	76		75		70-130	1		30	
1,1-Dichloroethane	104		105		70-130	1		30	
Chloroform	97		98		70-130	1		30	
Carbon tetrachloride	97		96		70-130	1		30	
1,2-Dichloropropane	103		105		70-130	2		30	
Dibromochloromethane	93		93		70-130	0		30	
1,1,2-Trichloroethane	99		99		70-130	0		30	
Tetrachloroethene	96		94		70-130	2		30	
Chlorobenzene	96		95		70-130	1		30	
Trichlorofluoromethane	110		113		70-139	3		30	
1,2-Dichloroethane	99		101		70-130	2		30	
1,1,1-Trichloroethane	100		100		70-130	0		30	
Bromodichloromethane	100		102		70-130	2		30	
trans-1,3-Dichloropropene	99		97		70-130	2		30	
cis-1,3-Dichloropropene	102		102		70-130	0		30	
Bromoform	94		95		70-130	1		30	
1,1,2,2-Tetrachloroethane	98		101		70-130	3		30	
Benzene	101		101		70-130	0		30	
Toluene	96		93		70-130	3		30	
Ethylbenzene	96		95		70-130	1		30	
Chloromethane	102		96		52-130	6		30	
Bromomethane	136		140		57-147	3		30	
Vinyl chloride	115		113		67-130	2		30	



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Parameter	LCS %Recovery	Qual	LCSD %Recovery	' Qual	%Recovery Limits	RPD	RPD Qual Limits
/olatile Organics by GC/MS - Westborough	Lab Associated	sample(s):	01-08 Batch:	WG1105645-3	WG1105645-4		
Chloroethane	107		107		50-151	0	30
1,1-Dichloroethene	101		101		65-135	0	30
trans-1,2-Dichloroethene	103		103		70-130	0	30
Trichloroethene	101		102		70-130	1	30
1,2-Dichlorobenzene	93		94		70-130	1	30
1,3-Dichlorobenzene	94		94		70-130	0	30
1,4-Dichlorobenzene	93		95		70-130	2	30
Methyl tert butyl ether	102		103		66-130	1	30
p/m-Xylene	96		95		70-130	1	30
o-Xylene	95		95		70-130	0	30
cis-1,2-Dichloroethene	103		102		70-130	1	30
Styrene	93		92		70-130	1	30
Dichlorodifluoromethane	100		95		30-146	5	30
Acetone	206	Q	224	Q	54-140	8	30
Carbon disulfide	104		91		59-130	13	30
2-Butanone	129		139	Q	70-130	7	30
4-Methyl-2-pentanone	99		99		70-130	0	30
2-Hexanone	102		109		70-130	7	30
Bromochloromethane	100		101		70-130	1	30
1,2-Dibromoethane	99		99		70-130	0	30
n-Butylbenzene	98		98		70-130	0	30
sec-Butylbenzene	98		97		70-130	1	30
tert-Butylbenzene	96		96		70-130	0	30



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Parameter	LCS %Recovery	Qual	_	CSD covery	Qual	%Recovery Limits	RPD	Qual	RPD Limits
/olatile Organics by GC/MS - Westborough L	ab Associated	sample(s):	01-08	Batch:	WG1105645-3	WG1105645-4			
1,2-Dibromo-3-chloropropane	92			97		68-130	5		30
Isopropylbenzene	97			97		70-130	0		30
p-Isopropyltoluene	95			95		70-130	0		30
Naphthalene	89			92		70-130	3		30
n-Propylbenzene	98			99		70-130	1		30
1,2,3-Trichlorobenzene	94			95		70-130	1		30
1,2,4-Trichlorobenzene	96			97		70-130	1		30
1,3,5-Trimethylbenzene	96			95		70-130	1		30
1,2,4-Trimethylbenzene	94			94		70-130	0		30
Methyl Acetate	81			83		51-146	2		30
Cyclohexane	102			100		59-142	2		30
1,4-Dioxane	96			101		65-136	5		30
Freon-113	105			104		50-139	1		30
Methyl cyclohexane	105			104		70-130	1		30

Surrogate	LCS	LCSD	Acceptance
	%Recovery Qual	%Recovery Qual	Criteria
1,2-Dichloroethane-d4	98	99	70-130
Toluene-d8	98	96	70-130
4-Bromofluorobenzene	104	104	70-130
Dibromofluoromethane	99	99	70-130



SEMIVOLATILES



L1812163

Project Name: SOUTHWESTERN PHASE 2

L1812163-01

04/12/18 18:31

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:17

Report Date: 04/16/18

Lab Number:

SS-1 Date Received: Client ID: 04/09/18

ORCHARD PARK, NY Sample Location: Field Prep: Not Specified

Sample Depth:

Analytical Date:

Lab ID:

Extraction Method: EPA 3546 Matrix: Soil **Extraction Date:** 04/10/18 17:49 Analytical Method: 1,8270D

Analyst: SZ 96% Percent Solids:

		Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - Westbook	rough Lab					
Acenaphthene	ND		ug/kg	140	18.	1
1,2,4-Trichlorobenzene	ND		ug/kg	170	20.	1
Hexachlorobenzene	ND		ug/kg	100	19.	1
Bis(2-chloroethyl)ether	ND		ug/kg	150	23.	1
2-Chloronaphthalene	ND		ug/kg	170	17.	1
1,2-Dichlorobenzene	ND		ug/kg	170	31.	1
1,3-Dichlorobenzene	ND		ug/kg	170	29.	1
1,4-Dichlorobenzene	ND		ug/kg	170	30.	1
3,3'-Dichlorobenzidine	ND		ug/kg	170	46.	1
2,4-Dinitrotoluene	ND		ug/kg	170	34.	1
2,6-Dinitrotoluene	ND		ug/kg	170	29.	1
Fluoranthene	ND		ug/kg	100	20.	1
4-Chlorophenyl phenyl ether	ND		ug/kg	170	18.	1
4-Bromophenyl phenyl ether	ND		ug/kg	170	26.	1
Bis(2-chloroisopropyl)ether	ND		ug/kg	200	29.	1
Bis(2-chloroethoxy)methane	ND		ug/kg	180	17.	1
Hexachlorobutadiene	ND		ug/kg	170	25.	1
Hexachlorocyclopentadiene	ND		ug/kg	490	160	1
Hexachloroethane	ND		ug/kg	140	28.	1
Isophorone	ND		ug/kg	150	22.	1
Naphthalene	ND		ug/kg	170	21.	1
Nitrobenzene	ND		ug/kg	150	25.	1
NDPA/DPA	ND		ug/kg	140	20.	1
n-Nitrosodi-n-propylamine	ND		ug/kg	170	26.	1
Bis(2-ethylhexyl)phthalate	ND		ug/kg	170	59.	1
Butyl benzyl phthalate	ND		ug/kg	170	43.	1
Di-n-butylphthalate	ND		ug/kg	170	32.	1
Di-n-octylphthalate	ND		ug/kg	170	58.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - We	estborough Lab					
Diethyl phthalate	ND		ug/kg	170	16.	1
Dimethyl phthalate	ND		ug/kg	170	36.	1
Benzo(a)anthracene	ND		ug/kg	100	19.	1
Benzo(a)pyrene	ND		ug/kg	140	42.	1
Benzo(b)fluoranthene	ND		ug/kg	100	29.	1
Benzo(k)fluoranthene	ND		ug/kg	100	27.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	26.	1
Anthracene	ND		ug/kg	100	33.	1
Benzo(ghi)perylene	ND		ug/kg	140	20.	1
Fluorene	ND		ug/kg	170	17.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	24.	1
Pyrene	ND		ug/kg	100	17.	1
Biphenyl	ND		ug/kg	390	40.	1
4-Chloroaniline	ND		ug/kg	170	31.	1
2-Nitroaniline	ND		ug/kg	170	33.	1
3-Nitroaniline	ND		ug/kg	170	32.	1
4-Nitroaniline	ND		ug/kg	170	71.	1
Dibenzofuran	ND		ug/kg	170	16.	1
2-Methylnaphthalene	ND		ug/kg	200	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	170	18.	1
Acetophenone	ND		ug/kg	170	21.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	32.	1
p-Chloro-m-cresol	ND		ug/kg	170	26.	1
2-Chlorophenol	ND		ug/kg	170	20.	1
2,4-Dichlorophenol	ND		ug/kg	150	28.	1
2,4-Dimethylphenol	ND		ug/kg	170	56.	1
2-Nitrophenol	ND		ug/kg	370	64.	1
4-Nitrophenol	ND		ug/kg	240	70.	1
2,4-Dinitrophenol	ND		ug/kg	820	80.	1
4,6-Dinitro-o-cresol	ND		ug/kg	440	82.	1
Pentachlorophenol	ND		ug/kg	140	38.	1
Phenol	ND		ug/kg	170	26.	1
2-Methylphenol	ND		ug/kg	170	26.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	250	27.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	170	33.	1
Benzoic Acid	ND		ug/kg	560	170	1
Benzyl Alcohol	ND		ug/kg	170	52.	1
Carbazole	ND		ua/ka	170	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria	
2-Fluorophenol	79	25-120	
Phenol-d6	81	10-120	
Nitrobenzene-d5	71	23-120	
2-Fluorobiphenyl	72	30-120	
2,4,6-Tribromophenol	89	10-136	
4-Terphenyl-d14	76	18-120	



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

L1812163-02

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:20

Lab Number:

Report Date:

Client ID: SS-2 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Lab ID:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8270D Extraction Date: 04/10/18 17:49
Analytical Date: 04/12/18 18:58

Analyst: SZ Percent Solids: 96%

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - Westbord	ough Lab					
Acenaphthene	ND		ug/kg	140	18.	1
1,2,4-Trichlorobenzene	ND		ug/kg	170	20.	1
Hexachlorobenzene	ND		ug/kg	100	19.	1
Bis(2-chloroethyl)ether	ND		ug/kg	150	23.	1
2-Chloronaphthalene	ND		ug/kg	170	17.	1
1,2-Dichlorobenzene	ND		ug/kg	170	31.	1
1,3-Dichlorobenzene	ND		ug/kg	170	29.	1
1,4-Dichlorobenzene	ND		ug/kg	170	30.	1
3,3'-Dichlorobenzidine	ND		ug/kg	170	46.	1
2,4-Dinitrotoluene	ND		ug/kg	170	34.	1
2,6-Dinitrotoluene	ND		ug/kg	170	29.	1
Fluoranthene	24	J	ug/kg	100	20.	1
4-Chlorophenyl phenyl ether	ND		ug/kg	170	18.	1
4-Bromophenyl phenyl ether	ND		ug/kg	170	26.	1
Bis(2-chloroisopropyl)ether	ND		ug/kg	200	29.	1
Bis(2-chloroethoxy)methane	ND		ug/kg	180	17.	1
Hexachlorobutadiene	ND		ug/kg	170	25.	1
Hexachlorocyclopentadiene	ND		ug/kg	490	160	1
Hexachloroethane	ND		ug/kg	140	28.	1
Isophorone	ND		ug/kg	150	22.	1
Naphthalene	ND		ug/kg	170	21.	1
Nitrobenzene	ND		ug/kg	150	25.	1
NDPA/DPA	ND		ug/kg	140	19.	1
n-Nitrosodi-n-propylamine	ND		ug/kg	170	26.	1
Bis(2-ethylhexyl)phthalate	ND		ug/kg	170	59.	1
Butyl benzyl phthalate	ND		ug/kg	170	43.	1
Di-n-butylphthalate	ND		ug/kg	170	32.	1
Di-n-octylphthalate	ND		ug/kg	170	58.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - \	Westborough Lab					
Diethyl phthalate	ND		ug/kg	170	16.	1
Dimethyl phthalate	ND		ug/kg	170	36.	1
Benzo(a)anthracene	ND		ug/kg	100	19.	1
Benzo(a)pyrene	ND		ug/kg	140	42.	1
Benzo(b)fluoranthene	ND		ug/kg	100	29.	1
Benzo(k)fluoranthene	ND		ug/kg	100	27.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	26.	1
Anthracene	ND		ug/kg	100	33.	1
Benzo(ghi)perylene	ND		ug/kg	140	20.	1
Fluorene	ND		ug/kg	170	17.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	24.	1
Pyrene	18	J	ug/kg	100	17.	1
Biphenyl	ND		ug/kg	390	40.	1
4-Chloroaniline	ND		ug/kg	170	31.	1
2-Nitroaniline	ND		ug/kg	170	33.	1
3-Nitroaniline	ND		ug/kg	170	32.	1
4-Nitroaniline	ND		ug/kg	170	71.	1
Dibenzofuran	ND		ug/kg	170	16.	1
2-Methylnaphthalene	ND		ug/kg	200	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	170	18.	1
Acetophenone	ND		ug/kg	170	21.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	32.	1
p-Chloro-m-cresol	ND		ug/kg	170	26.	1
2-Chlorophenol	ND		ug/kg	170	20.	1
2,4-Dichlorophenol	ND		ug/kg	150	28.	1
2,4-Dimethylphenol	ND		ug/kg	170	56.	1
2-Nitrophenol	ND		ug/kg	370	64.	1
4-Nitrophenol	ND		ug/kg	240	70.	1
2,4-Dinitrophenol	ND		ug/kg	820	80.	1
4,6-Dinitro-o-cresol	ND		ug/kg	440	82.	1
Pentachlorophenol	ND		ug/kg	140	38.	1
Phenol	ND		ug/kg	170	26.	1
2-Methylphenol	ND		ug/kg	170	26.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	250	27.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	170	33.	1
Benzoic Acid	ND		ug/kg	550	170	1
Benzyl Alcohol	ND		ug/kg	170	52.	1
Carbazole	ND		ua/ka	170	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	67	25-120
Phenol-d6	71	10-120
Nitrobenzene-d5	68	23-120
2-Fluorobiphenyl	73	30-120
2,4,6-Tribromophenol	94	10-136
4-Terphenyl-d14	84	18-120

L1812163

04/09/18

Dilution Factor

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:22

Report Date: 04/16/18

Lab ID: L1812163-03

Client ID: SS-3

Sample Location: ORCHARD PARK, NY

Not Specified

Lab Number:

Date Received:

Field Prep:

RL

Sample Depth:

Parameter

Matrix: Soil 1,8270D Analytical Method: Analytical Date: 04/12/18 19:24

Analyst: SZ 93% Percent Solids:

Extraction Method: EPA 3546 **Extraction Date:** 04/10/18 17:49

MDL

Parameter	Result	Qualifier	Units	KL	MDL	Dilution Factor	
Semivolatile Organics by GC/MS - V	Vestborough Lab						
Acenaphthene	ND		ug/kg	140	18.	1	
1,2,4-Trichlorobenzene	ND		ug/kg	180	20.	1	
Hexachlorobenzene	ND		ug/kg	100	20.	1	
Bis(2-chloroethyl)ether	ND		ug/kg	160	24.	1	
2-Chloronaphthalene	ND		ug/kg	180	18.	1	
1,2-Dichlorobenzene	ND		ug/kg	180	32.	1	
1,3-Dichlorobenzene	ND		ug/kg	180	30.	1	
1,4-Dichlorobenzene	ND		ug/kg	180	31.	1	
3,3'-Dichlorobenzidine	ND		ug/kg	180	47.	1	
2,4-Dinitrotoluene	ND		ug/kg	180	35.	1	
2,6-Dinitrotoluene	ND		ug/kg	180	30.	1	
Fluoranthene	23	J	ug/kg	100	20.	1	
4-Chlorophenyl phenyl ether	ND		ug/kg	180	19.	1	
4-Bromophenyl phenyl ether	ND		ug/kg	180	27.	1	
Bis(2-chloroisopropyl)ether	ND		ug/kg	210	30.	1	
Bis(2-chloroethoxy)methane	ND		ug/kg	190	18.	1	
Hexachlorobutadiene	ND		ug/kg	180	26.	1	
Hexachlorocyclopentadiene	ND		ug/kg	500	160	1	
Hexachloroethane	ND		ug/kg	140	28.	1	
Isophorone	ND		ug/kg	160	23.	1	
Naphthalene	ND		ug/kg	180	22.	1	
Nitrobenzene	ND		ug/kg	160	26.	1	
NDPA/DPA	ND		ug/kg	140	20.	1	
n-Nitrosodi-n-propylamine	ND		ug/kg	180	27.	1	
Bis(2-ethylhexyl)phthalate	ND		ug/kg	180	61.	1	
Butyl benzyl phthalate	ND		ug/kg	180	44.	1	
Di-n-butylphthalate	ND		ug/kg	180	33.	1	
Di-n-octylphthalate	ND		ug/kg	180	60.	1	

Qualifier

Result

Units



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03 Date Collected: 04/06/18 08:22

Client ID: SS-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS -	Westborough Lab					
Diethyl phthalate	ND		ug/kg	180	16.	1
Dimethyl phthalate	ND		ug/kg	180	37.	1
Benzo(a)anthracene	ND		ug/kg	100	20.	1
Benzo(a)pyrene	ND		ug/kg	140	43.	1
Benzo(b)fluoranthene	ND		ug/kg	100	30.	1
Benzo(k)fluoranthene	ND		ug/kg	100	28.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	27.	1
Anthracene	ND		ug/kg	100	34.	1
Benzo(ghi)perylene	ND		ug/kg	140	21.	1
Fluorene	ND		ug/kg	180	17.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	25.	1
Pyrene	18	J	ug/kg	100	18.	1
Biphenyl	ND		ug/kg	400	41.	1
4-Chloroaniline	ND		ug/kg	180	32.	1
2-Nitroaniline	ND		ug/kg	180	34.	1
3-Nitroaniline	ND		ug/kg	180	33.	1
4-Nitroaniline	ND		ug/kg	180	73.	1
Dibenzofuran	ND		ug/kg	180	17.	1
2-Methylnaphthalene	ND		ug/kg	210	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	180	18.	1
Acetophenone	ND		ug/kg	180	22.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	33.	1
p-Chloro-m-cresol	ND		ug/kg	180	26.	1
2-Chlorophenol	ND		ug/kg	180	21.	1
2,4-Dichlorophenol	ND		ug/kg	160	28.	1
2,4-Dimethylphenol	ND		ug/kg	180	58.	1
2-Nitrophenol	ND		ug/kg	380	66.	1
4-Nitrophenol	ND		ug/kg	250	72.	1
2,4-Dinitrophenol	ND		ug/kg	850	82.	1
4,6-Dinitro-o-cresol	ND		ug/kg	460	85.	1
Pentachlorophenol	ND		ug/kg	140	39.	1
Phenol	ND		ug/kg	180	27.	1
2-Methylphenol	ND		ug/kg	180	27.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	250	28.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03 Date Collected: 04/06/18 08:22

Client ID: SS-3 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	180	34.	1
Benzoic Acid	280	J	ug/kg	570	180	1
Benzyl Alcohol	ND		ug/kg	180	54.	1
Carbazole	ND		ua/ka	180	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	71	25-120
Phenol-d6	73	10-120
Nitrobenzene-d5	69	23-120
2-Fluorobiphenyl	66	30-120
2,4,6-Tribromophenol	80	10-136
4-Terphenyl-d14	65	18-120



L1812163

04/09/18

Not Specified

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:25

04/16/18

Report Date:

Lab Number:

Date Received:

Field Prep:

Lab ID: L1812163-04

SS-4 Client ID:

Sample Location: ORCHARD PARK, NY

Sample Depth:

Matrix: Soil Analytical Method: 1,8270D Analytical Date: 04/12/18 19:49

Analyst: SZ 96% Percent Solids:

Extraction Method: EPA 3546 **Extraction Date:** 04/10/18 17:49

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Semivolatile Organics by GC/MS -	Westborough Lab						
Acenaphthene	ND		ug/kg	140	18.	1	
1,2,4-Trichlorobenzene	ND		ug/kg	170	20.	1	
Hexachlorobenzene	ND		ug/kg	100	19.	1	
Bis(2-chloroethyl)ether	ND		ug/kg	160	24.	1	
2-Chloronaphthalene	ND		ug/kg	170	17.	1	
1,2-Dichlorobenzene	ND		ug/kg	170	31.	1	
1,3-Dichlorobenzene	ND		ug/kg	170	30.	1	
1,4-Dichlorobenzene	ND		ug/kg	170	30.	1	
3,3'-Dichlorobenzidine	ND		ug/kg	170	46.	1	
2,4-Dinitrotoluene	ND		ug/kg	170	35.	1	
2,6-Dinitrotoluene	ND		ug/kg	170	30.	1	
Fluoranthene	ND		ug/kg	100	20.	1	
4-Chlorophenyl phenyl ether	ND		ug/kg	170	18.	1	
4-Bromophenyl phenyl ether	ND		ug/kg	170	26.	1	
Bis(2-chloroisopropyl)ether	ND		ug/kg	210	30.	1	
Bis(2-chloroethoxy)methane	ND		ug/kg	190	17.	1	
Hexachlorobutadiene	ND		ug/kg	170	25.	1	
Hexachlorocyclopentadiene	ND		ug/kg	500	160	1	
Hexachloroethane	ND		ug/kg	140	28.	1	
Isophorone	ND		ug/kg	160	22.	1	
Naphthalene	ND		ug/kg	170	21.	1	
Nitrobenzene	ND		ug/kg	160	26.	1	
NDPA/DPA	ND		ug/kg	140	20.	1	
n-Nitrosodi-n-propylamine	ND		ug/kg	170	27.	1	
Bis(2-ethylhexyl)phthalate	ND		ug/kg	170	60.	1	
Butyl benzyl phthalate	ND		ug/kg	170	44.	1	
Di-n-butylphthalate	ND		ug/kg	170	33.	1	
Di-n-octylphthalate	ND		ug/kg	170	59.	1	



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04 Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - W	/estborough Lab					
Disthict white state	ND		4	470	40	4
Diethyl phthalate	ND		ug/kg	170	16.	1
Dimethyl phthalate	ND		ug/kg	170	36.	1
Benzo(a)anthracene	ND		ug/kg	100	20.	1
Benzo(a)pyrene	ND		ug/kg	140	42.	1
Benzo(b)fluoranthene	ND		ug/kg	100	29.	1
Benzo(k)fluoranthene	ND		ug/kg	100	28.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	27.	1
Anthracene	ND		ug/kg	100	34.	1
Benzo(ghi)perylene	ND		ug/kg	140	20.	1
Fluorene	ND		ug/kg	170	17.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	24.	1
Pyrene	ND		ug/kg	100	17.	1
Biphenyl	ND		ug/kg	400	40.	1
4-Chloroaniline	ND		ug/kg	170	32.	1
2-Nitroaniline	ND		ug/kg	170	33.	1
3-Nitroaniline	ND		ug/kg	170	33.	1
4-Nitroaniline	ND		ug/kg	170	72.	1
Dibenzofuran	ND		ug/kg	170	16.	1
2-Methylnaphthalene	ND		ug/kg	210	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	170	18.	1
Acetophenone	ND		ug/kg	170	21.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	33.	1
p-Chloro-m-cresol	ND		ug/kg	170	26.	1
2-Chlorophenol	ND		ug/kg	170	20.	1
2,4-Dichlorophenol	ND		ug/kg	160	28.	1
2,4-Dimethylphenol	ND		ug/kg	170	57.	1
2-Nitrophenol	ND		ug/kg	370	65.	1
4-Nitrophenol	ND		ug/kg	240	71.	1
2,4-Dinitrophenol	ND		ug/kg	830	81.	1
4,6-Dinitro-o-cresol	ND		ug/kg	450	83.	1
Pentachlorophenol	ND		ug/kg	140	38.	1
Phenol	ND		ug/kg	170	26.	1
2-Methylphenol	ND		ug/kg	170	27.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	250	27.	1
	110		ug/kg	200	۷۱.	·



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04 Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - We	stborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	170	33.	1
Benzoic Acid	250	J	ug/kg	560	180	1
Benzyl Alcohol	ND		ug/kg	170	53.	1
Carbazole	ND		ug/kg	170	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	77	25-120
Phenol-d6	80	10-120
Nitrobenzene-d5	72	23-120
2-Fluorobiphenyl	73	30-120
2,4,6-Tribromophenol	89	10-136
4-Terphenyl-d14	74	18-120



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

L1812163-05

04/12/18 20:15

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 08:29

Lab Number:

Report Date:

Date Received: SS-5 04/09/18 ORCHARD PARK, NY Field Prep: Not Specified

Sample Location:

Sample Depth:

Analytical Date:

Lab ID:

Client ID:

Extraction Method: EPA 3546 Matrix: Soil **Extraction Date:** 04/10/18 17:49 Analytical Method: 1,8270D

Analyst: SZ 96% Percent Solids:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Semivolatile Organics by GC/MS - Wes	stborough Lab						
Acenaphthene	ND		ug/kg	140	18.	1	
1,2,4-Trichlorobenzene	ND		ug/kg	170	20.	1	
Hexachlorobenzene	ND		ug/kg	100	19.	1	
Bis(2-chloroethyl)ether	ND		ug/kg	150	23.	1	
2-Chloronaphthalene	ND		ug/kg	170	17.	1	
1,2-Dichlorobenzene	ND		ug/kg	170	31.	1	
1,3-Dichlorobenzene	ND		ug/kg	170	29.	1	
1,4-Dichlorobenzene	ND		ug/kg	170	30.	1	
3,3'-Dichlorobenzidine	ND		ug/kg	170	45.	1	
2,4-Dinitrotoluene	ND		ug/kg	170	34.	1	
2,6-Dinitrotoluene	ND		ug/kg	170	29.	1	
Fluoranthene	29	J	ug/kg	100	20.	1	
4-Chlorophenyl phenyl ether	ND		ug/kg	170	18.	1	
4-Bromophenyl phenyl ether	ND		ug/kg	170	26.	1	
Bis(2-chloroisopropyl)ether	ND		ug/kg	200	29.	1	
Bis(2-chloroethoxy)methane	ND		ug/kg	180	17.	1	
Hexachlorobutadiene	ND		ug/kg	170	25.	1	
Hexachlorocyclopentadiene	ND		ug/kg	490	150	1	
Hexachloroethane	ND		ug/kg	140	28.	1	
Isophorone	ND		ug/kg	150	22.	1	
Naphthalene	ND		ug/kg	170	21.	1	
Nitrobenzene	ND		ug/kg	150	25.	1	
NDPA/DPA	ND		ug/kg	140	19.	1	
n-Nitrosodi-n-propylamine	ND		ug/kg	170	26.	1	
Bis(2-ethylhexyl)phthalate	ND		ug/kg	170	59.	1	
Butyl benzyl phthalate	ND		ug/kg	170	43.	1	
Di-n-butylphthalate	ND		ug/kg	170	32.	1	
Di-n-octylphthalate	ND		ug/kg	170	58.	1	



Project Name: Lab Number: SOUTHWESTERN PHASE 2 L1812163

Project Number: Report Date: Not Specified 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05 Date Collected: 04/06/18 08:29

Client ID: SS-5 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS -	· Westborough Lab					
District white state	ND			470	40	4
Diethyl phthalate	ND		ug/kg	170	16.	1
Dimethyl phthalate	ND		ug/kg	170	36.	1
Benzo(a)anthracene	ND		ug/kg	100	19.	1
Benzo(a)pyrene	ND		ug/kg	140	42.	1
Benzo(b)fluoranthene	ND		ug/kg	100	29.	1
Benzo(k)fluoranthene	ND		ug/kg	100	27.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	26.	1
Anthracene	ND		ug/kg	100	33.	1
Benzo(ghi)perylene	ND		ug/kg	140	20.	1
Fluorene	ND		ug/kg	170	16.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	24.	1
Pyrene	24	J	ug/kg	100	17.	1
Biphenyl	ND		ug/kg	390	40.	1
4-Chloroaniline	ND		ug/kg	170	31.	1
2-Nitroaniline	ND		ug/kg	170	33.	1
3-Nitroaniline	ND		ug/kg	170	32.	1
4-Nitroaniline	ND		ug/kg	170	71.	1
Dibenzofuran	ND		ug/kg	170	16.	1
2-Methylnaphthalene	ND		ug/kg	200	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	170	18.	1
Acetophenone	ND		ug/kg	170	21.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	32.	1
p-Chloro-m-cresol	ND		ug/kg	170	25.	1
2-Chlorophenol	ND		ug/kg	170	20.	1
2,4-Dichlorophenol	ND		ug/kg	150	27.	1
2,4-Dimethylphenol	ND		ug/kg	170	56.	1
2-Nitrophenol	ND		ug/kg	370	64.	1
4-Nitrophenol	ND		ug/kg	240	70.	1
2,4-Dinitrophenol	ND		ug/kg	820	80.	1
4,6-Dinitro-o-cresol	ND		ug/kg	440	82.	1
Pentachlorophenol	ND		ug/kg	140	38.	1
Phenol	ND		ug/kg	170	26.	1
2-Methylphenol	ND		ug/kg	170	26.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	240	27.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05 Date Collected: 04/06/18 08:29

Client ID: SS-5 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Semivolatile Organics by GC/M	S - Westborough Lab						
2,4,5-Trichlorophenol	ND		ug/kg	170	33.	1	
Benzoic Acid	220	J	ug/kg	550	170	1	
Benzyl Alcohol	ND		ug/kg	170	52.	1	
Carbazole	ND		ua/ka	170	16.	1	

Surrogate	% Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	69	25-120
Phenol-d6	73	10-120
Nitrobenzene-d5	68	23-120
2-Fluorobiphenyl	66	30-120
2,4,6-Tribromophenol	80	10-136
4-Terphenyl-d14	68	18-120

L1812163

04/09/18

Not Specified

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 09:30

Report Date: 04/16/18

Lab Number:

Date Received:

Field Prep:

Lab ID: L1812163-06 Date

Client ID: B-1

Sample Location: ORCHARD PARK, NY

Sample Depth:

Matrix: Soil
Analytical Method: 1,8270D
Analytical Date: 04/12/18 20:41

Analyst: SZ Percent Solids: 93% Extraction Method: EPA 3546
Extraction Date: 04/10/18 17:49

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - \	Westborough Lab					
Acenaphthene	ND		ug/kg	140	18.	1
1,2,4-Trichlorobenzene	ND		ug/kg	170	20.	1
Hexachlorobenzene	ND		ug/kg	100	19.	1
Bis(2-chloroethyl)ether	ND		ug/kg	160	24.	1
2-Chloronaphthalene	ND		ug/kg	170	17.	1
1,2-Dichlorobenzene	ND		ug/kg	170	31.	1
1,3-Dichlorobenzene	ND		ug/kg	170	30.	1
1,4-Dichlorobenzene	ND		ug/kg	170	30.	1
3,3'-Dichlorobenzidine	ND		ug/kg	170	46.	1
2,4-Dinitrotoluene	ND		ug/kg	170	35.	1
2,6-Dinitrotoluene	ND		ug/kg	170	30.	1
Fluoranthene	ND		ug/kg	100	20.	1
4-Chlorophenyl phenyl ether	ND		ug/kg	170	19.	1
4-Bromophenyl phenyl ether	ND		ug/kg	170	26.	1
Bis(2-chloroisopropyl)ether	ND		ug/kg	210	30.	1
Bis(2-chloroethoxy)methane	ND		ug/kg	190	17.	1
Hexachlorobutadiene	ND		ug/kg	170	25.	1
Hexachlorocyclopentadiene	ND		ug/kg	500	160	1
Hexachloroethane	ND		ug/kg	140	28.	1
Isophorone	ND		ug/kg	160	22.	1
Naphthalene	ND		ug/kg	170	21.	1
Nitrobenzene	ND		ug/kg	160	26.	1
NDPA/DPA	ND		ug/kg	140	20.	1
n-Nitrosodi-n-propylamine	ND		ug/kg	170	27.	1
Bis(2-ethylhexyl)phthalate	ND		ug/kg	170	60.	1
Butyl benzyl phthalate	ND		ug/kg	170	44.	1
Di-n-butylphthalate	ND		ug/kg	170	33.	1
Di-n-octylphthalate	ND		ug/kg	170	59.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS - W	estborough Lab					
Diethyl phthalate	ND		ug/kg	170	16.	1
Dimethyl phthalate	ND		ug/kg	170	36.	1
Benzo(a)anthracene	ND		ug/kg	100	20.	1
Benzo(a)pyrene	ND		ug/kg	140	42.	1
Benzo(b)fluoranthene	ND		ug/kg	100	29.	1
Benzo(k)fluoranthene	ND		ug/kg	100	28.	1
Chrysene	ND		ug/kg	100	18.	1
Acenaphthylene	ND		ug/kg	140	27.	1
Anthracene	ND		ug/kg	100	34.	1
Benzo(ghi)perylene	ND		ug/kg	140	20.	1
Fluorene	ND		ug/kg	170	17.	1
Phenanthrene	ND		ug/kg	100	21.	1
Dibenzo(a,h)anthracene	ND		ug/kg	100	20.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	24.	1
Pyrene	ND		ug/kg	100	17.	1
Biphenyl	ND		ug/kg	400	40.	1
4-Chloroaniline	ND		ug/kg	170	32.	1
2-Nitroaniline	ND		ug/kg	170	34.	1
3-Nitroaniline	ND		ug/kg	170	33.	1
4-Nitroaniline	ND		ug/kg	170	72.	1
Dibenzofuran	ND		ug/kg	170	16.	1
2-Methylnaphthalene	ND		ug/kg	210	21.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	170	18.	1
Acetophenone	ND		ug/kg	170	22.	1
2,4,6-Trichlorophenol	ND		ug/kg	100	33.	1
p-Chloro-m-cresol	ND		ug/kg	170	26.	1
2-Chlorophenol	ND		ug/kg	170	20.	1
2,4-Dichlorophenol	ND		ug/kg	160	28.	1
2,4-Dimethylphenol	ND		ug/kg	170	57.	1
2-Nitrophenol	ND		ug/kg	380	65.	1
4-Nitrophenol	ND		ug/kg	240	71.	1
2,4-Dinitrophenol	ND		ug/kg	840	81.	1
4,6-Dinitro-o-cresol	ND		ug/kg	450	84.	1
Pentachlorophenol	ND		ug/kg	140	38.	1
Phenol	ND		ug/kg	170	26.	1
2-Methylphenol	ND		ug/kg	170	27.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	250	27.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: B-1 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	170	33.	1
Benzoic Acid	ND		ug/kg	560	180	1
Benzyl Alcohol	ND		ug/kg	170	53.	1
Carbazole	ND		ua/ka	170	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria	
2-Fluorophenol	74	25-120	
Phenol-d6	80	10-120	
Nitrobenzene-d5	75	23-120	
2-Fluorobiphenyl	69	30-120	
2,4,6-Tribromophenol	55	10-136	
4-Terphenyl-d14	71	18-120	



L1812163

04/09/18

Not Specified

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 09:50

04/16/18

Lab Number:

Report Date:

Date Received:

Lab ID: L1812163-07

Client ID: B-2

ORCHARD PARK, NY Sample Location:

Field Prep:

Sample Depth:

Matrix: Soil Analytical Method: 1,8270D Analytical Date: 04/12/18 21:07

Analyst: SZ 87% Percent Solids:

Extraction Method: EPA 3546 **Extraction Date:** 04/10/18 17:49

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Semivolatile Organics by GC/MS - V	Vestborough Lab						
Acenaphthene	ND		ug/kg	150	20.	1	
1,2,4-Trichlorobenzene	ND		ug/kg	190	22.	1	
Hexachlorobenzene	ND		ug/kg	110	21.	1	
Bis(2-chloroethyl)ether	ND		ug/kg	170	26.	1	
2-Chloronaphthalene	ND		ug/kg	190	19.	1	
1,2-Dichlorobenzene	ND		ug/kg	190	34.	1	
1,3-Dichlorobenzene	ND		ug/kg	190	33.	1	
1,4-Dichlorobenzene	ND		ug/kg	190	33.	1	
3,3'-Dichlorobenzidine	ND		ug/kg	190	51.	1	
2,4-Dinitrotoluene	ND		ug/kg	190	38.	1	
2,6-Dinitrotoluene	ND		ug/kg	190	33.	1	
Fluoranthene	ND		ug/kg	110	22.	1	
4-Chlorophenyl phenyl ether	ND		ug/kg	190	20.	1	
4-Bromophenyl phenyl ether	ND		ug/kg	190	29.	1	
Bis(2-chloroisopropyl)ether	ND		ug/kg	230	32.	1	
Bis(2-chloroethoxy)methane	ND		ug/kg	200	19.	1	
Hexachlorobutadiene	ND		ug/kg	190	28.	1	
Hexachlorocyclopentadiene	ND		ug/kg	540	170	1	
Hexachloroethane	ND		ug/kg	150	31.	1	
Isophorone	ND		ug/kg	170	25.	1	
Naphthalene	ND		ug/kg	190	23.	1	
Nitrobenzene	ND		ug/kg	170	28.	1	
NDPA/DPA	ND		ug/kg	150	22.	1	
n-Nitrosodi-n-propylamine	ND		ug/kg	190	29.	1	
Bis(2-ethylhexyl)phthalate	ND		ug/kg	190	66.	1	
Butyl benzyl phthalate	ND		ug/kg	190	48.	1	
Di-n-butylphthalate	ND		ug/kg	190	36.	1	
Di-n-octylphthalate	ND		ug/kg	190	65.	1	

Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07 Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS	- Westborough Lab					
Diethyl phthalate	ND		ug/kg	190	18.	1
Dimethyl phthalate	ND		ug/kg	190	40.	1
Benzo(a)anthracene	ND		ug/kg	110	21.	1
Benzo(a)pyrene	ND		ug/kg	150	46.	1
Benzo(b)fluoranthene	ND		ug/kg	110	32.	1
Benzo(k)fluoranthene	ND		ug/kg	110	30.	1
Chrysene	ND		ug/kg	110	20.	1
Acenaphthylene	ND		ug/kg	150	29.	1
Anthracene	ND		ug/kg	110	37.	1
Benzo(ghi)perylene	ND		ug/kg	150	22.	1
Fluorene	ND		ug/kg	190	18.	1
Phenanthrene	ND		ug/kg	110	23.	1
Dibenzo(a,h)anthracene	ND		ug/kg	110	22.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	150	26.	1
Pyrene	ND		ug/kg	110	19.	1
Biphenyl	ND		ug/kg	430	44.	1
4-Chloroaniline	ND		ug/kg	190	35.	1
2-Nitroaniline	ND		ug/kg	190	37.	1
3-Nitroaniline	ND		ug/kg	190	36.	1
4-Nitroaniline	ND		ug/kg	190	79.	1
Dibenzofuran	ND		ug/kg	190	18.	1
2-Methylnaphthalene	ND		ug/kg	230	23.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	190	20.	1
Acetophenone	ND		ug/kg	190	24.	1
2,4,6-Trichlorophenol	ND		ug/kg	110	36.	1
p-Chloro-m-cresol	ND		ug/kg	190	28.	1
2-Chlorophenol	ND		ug/kg	190	22.	1
2,4-Dichlorophenol	ND		ug/kg	170	31.	1
2,4-Dimethylphenol	ND		ug/kg	190	63.	1
2-Nitrophenol	ND		ug/kg	410	72.	1
4-Nitrophenol	ND		ug/kg	270	78.	1
2,4-Dinitrophenol	ND		ug/kg	910	89.	1
4,6-Dinitro-o-cresol	ND		ug/kg	490	91.	1
Pentachlorophenol	ND		ug/kg	150	42.	1
Phenol	ND		ug/kg	190	29.	1
2-Methylphenol	ND		ug/kg	190	29.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	270	30.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07 Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	190	36.	1
Benzoic Acid	ND		ug/kg	620	190	1
Benzyl Alcohol	ND		ug/kg	190	58.	1
Carbazole	ND		ua/ka	190	18.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	80	25-120
Phenol-d6	84	10-120
Nitrobenzene-d5	76	23-120
2-Fluorobiphenyl	72	30-120
2,4,6-Tribromophenol	80	10-136
4-Terphenyl-d14	74	18-120



L1812163

04/16/18

Project Name: SOUTHWESTERN PHASE 2

B-3

L1812163-08

ORCHARD PARK, NY

Project Number: Not Specified

SAMPLE RESULTS

Date Collected: 04/06/18 10:10

Lab Number:

Report Date:

Date Received: 04/09/18

Field Prep: Not Specified

Sample Depth:

Sample Location:

Lab ID:

Client ID:

Matrix: Soil
Analytical Method: 1,8270D
Analytical Date: 04/12/18 21:33

Analyst: SZ Percent Solids: 91% Extraction Method: EPA 3546
Extraction Date: 04/10/18 17:49

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	
Semivolatile Organics by GC/MS -	Westborough Lab						
Acenaphthene	ND		ug/kg	140	18.	1	
1,2,4-Trichlorobenzene	ND		ug/kg	180	20.	1	
Hexachlorobenzene	ND		ug/kg	110	20.	1	
Bis(2-chloroethyl)ether	ND		ug/kg	160	24.	1	
2-Chloronaphthalene	ND		ug/kg	180	18.	1	
1,2-Dichlorobenzene	ND		ug/kg	180	32.	1	
1,3-Dichlorobenzene	ND		ug/kg	180	31.	1	
1,4-Dichlorobenzene	ND		ug/kg	180	31.	1	
3,3'-Dichlorobenzidine	ND		ug/kg	180	48.	1	
2,4-Dinitrotoluene	ND		ug/kg	180	36.	1	
2,6-Dinitrotoluene	ND		ug/kg	180	31.	1	
Fluoranthene	ND		ug/kg	110	20.	1	
4-Chlorophenyl phenyl ether	ND		ug/kg	180	19.	1	
4-Bromophenyl phenyl ether	ND		ug/kg	180	27.	1	
Bis(2-chloroisopropyl)ether	ND		ug/kg	210	30.	1	
Bis(2-chloroethoxy)methane	ND		ug/kg	190	18.	1	
Hexachlorobutadiene	ND		ug/kg	180	26.	1	
Hexachlorocyclopentadiene	ND		ug/kg	510	160	1	
Hexachloroethane	ND		ug/kg	140	29.	1	
Isophorone	ND		ug/kg	160	23.	1	
Naphthalene	ND		ug/kg	180	22.	1	
Nitrobenzene	ND		ug/kg	160	26.	1	
NDPA/DPA	ND		ug/kg	140	20.	1	
n-Nitrosodi-n-propylamine	ND		ug/kg	180	28.	1	
Bis(2-ethylhexyl)phthalate	ND		ug/kg	180	62.	1	
Butyl benzyl phthalate	ND		ug/kg	180	45.	1	
Di-n-butylphthalate	ND		ug/kg	180	34.	1	
Di-n-octylphthalate	ND		ug/kg	180	61.	1	



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08 Date Collected: 04/06/18 10:10

Client ID: B-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/MS	- Westborough Lab					
Diethyl phthalate	ND		ug/kg	180	16.	1
Dimethyl phthalate	ND		ug/kg	180	38.	1
Benzo(a)anthracene	ND		ug/kg	110	20.	1
Benzo(a)pyrene	ND		ug/kg	140	44.	1
Benzo(b)fluoranthene	ND		ug/kg	110	30.	1
Benzo(k)fluoranthene	ND		ug/kg	110	29.	1
Chrysene	ND		ug/kg	110	18.	1
Acenaphthylene	ND		ug/kg	140	28.	1
Anthracene	ND		ug/kg	110	35.	1
Benzo(ghi)perylene	ND		ug/kg	140	21.	1
Fluorene	ND		ug/kg	180	17.	1
Phenanthrene	ND		ug/kg	110	22.	1
Dibenzo(a,h)anthracene	ND		ug/kg	110	21.	1
Indeno(1,2,3-cd)pyrene	ND		ug/kg	140	25.	1
Pyrene	ND		ug/kg	110	18.	1
Biphenyl	ND		ug/kg	410	41.	1
4-Chloroaniline	ND		ug/kg	180	32.	1
2-Nitroaniline	ND		ug/kg	180	34.	1
3-Nitroaniline	ND		ug/kg	180	34.	1
4-Nitroaniline	ND		ug/kg	180	74.	1
Dibenzofuran	ND		ug/kg	180	17.	1
2-Methylnaphthalene	ND		ug/kg	210	22.	1
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	180	19.	1
Acetophenone	ND		ug/kg	180	22.	1
2,4,6-Trichlorophenol	ND		ug/kg	110	34.	1
p-Chloro-m-cresol	ND		ug/kg	180	27.	1
2-Chlorophenol	ND		ug/kg	180	21.	1
2,4-Dichlorophenol	ND		ug/kg	160	29.	1
2,4-Dimethylphenol	ND		ug/kg	180	59.	1
2-Nitrophenol	ND		ug/kg	390	67.	1
4-Nitrophenol	ND		ug/kg	250	73.	1
2,4-Dinitrophenol	ND		ug/kg	860	83.	1
4,6-Dinitro-o-cresol	ND		ug/kg	460	86.	1
Pentachlorophenol	ND		ug/kg	140	39.	1
Phenol	ND		ug/kg	180	27.	1
2-Methylphenol	ND		ug/kg	180	28.	1
3-Methylphenol/4-Methylphenol	ND		ug/kg	260	28.	1



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08 Date Collected: 04/06/18 10:10

Client ID: B-3 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor
Semivolatile Organics by GC/M	S - Westborough Lab					
2,4,5-Trichlorophenol	ND		ug/kg	180	34.	1
Benzoic Acid	ND		ug/kg	580	180	1
Benzyl Alcohol	ND		ug/kg	180	55.	1
Carbazole	ND		ua/ka	180	17.	1

Surrogate	% Recovery	Acceptance Qualifier Criteria	
2-Fluorophenol	81	25-120	
Phenol-d6	85	10-120	
Nitrobenzene-d5	80	23-120	
2-Fluorobiphenyl	73	30-120	
2,4,6-Tribromophenol	62	10-136	
4-Terphenyl-d14	77	18-120	



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: 1,8270D Analytical Date: 04/11/18 22:56

Analyst: TT

Extraction Method: EPA 3546
Extraction Date: 04/10/18 17:49

arameter	Result	Qualifier	Units	RL		MDL
emivolatile Organics by GC/MS -	Westborough	Lab for s	sample(s):	01-08	Batch:	WG1105283-1
Acenaphthene	ND		ug/kg	130		17.
1,2,4-Trichlorobenzene	ND		ug/kg	160		18.
Hexachlorobenzene	ND		ug/kg	97		18.
Bis(2-chloroethyl)ether	ND		ug/kg	150		22.
2-Chloronaphthalene	ND		ug/kg	160		16.
1,2-Dichlorobenzene	ND		ug/kg	160		29.
1,3-Dichlorobenzene	ND		ug/kg	160		28.
1,4-Dichlorobenzene	ND		ug/kg	160		28.
3,3'-Dichlorobenzidine	ND		ug/kg	160		43.
2,4-Dinitrotoluene	ND		ug/kg	160		32.
2,6-Dinitrotoluene	ND		ug/kg	160		28.
Fluoranthene	ND		ug/kg	97		19.
4-Chlorophenyl phenyl ether	ND		ug/kg	160		17.
4-Bromophenyl phenyl ether	ND		ug/kg	160		25.
Bis(2-chloroisopropyl)ether	ND		ug/kg	190		28.
Bis(2-chloroethoxy)methane	ND		ug/kg	180		16.
Hexachlorobutadiene	ND		ug/kg	160		24.
Hexachlorocyclopentadiene	ND		ug/kg	460		150
Hexachloroethane	ND		ug/kg	130		26.
Isophorone	ND		ug/kg	150		21.
Naphthalene	ND		ug/kg	160		20.
Nitrobenzene	ND		ug/kg	150		24.
NDPA/DPA	ND		ug/kg	130		18.
n-Nitrosodi-n-propylamine	ND		ug/kg	160		25.
Bis(2-ethylhexyl)phthalate	ND		ug/kg	160		56.
Butyl benzyl phthalate	ND		ug/kg	160		41.
Di-n-butylphthalate	ND		ug/kg	160		31.
Di-n-octylphthalate	ND		ug/kg	160		55.
Diethyl phthalate	ND		ug/kg	160		15.



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: 1,8270D Analytical Date: 04/11/18 22:56

Analyst: TT

Extraction Method: EPA 3546
Extraction Date: 04/10/18 17:49

arameter	Result	Qualifier	Units	RL		MDL
Semivolatile Organics by GC/MS	- Westborough	Lab for s	ample(s):	01-08	Batch:	WG1105283-1
Dimethyl phthalate	ND		ug/kg	160		34.
Benzo(a)anthracene	ND		ug/kg	97		18.
Benzo(a)pyrene	ND		ug/kg	130		40.
Benzo(b)fluoranthene	ND		ug/kg	97		27.
Benzo(k)fluoranthene	ND		ug/kg	97		26.
Chrysene	ND		ug/kg	97		17.
Acenaphthylene	ND		ug/kg	130		25.
Anthracene	ND		ug/kg	97		32.
Benzo(ghi)perylene	ND		ug/kg	130		19.
Fluorene	ND		ug/kg	160		16.
Phenanthrene	ND		ug/kg	97		20.
Dibenzo(a,h)anthracene	ND		ug/kg	97		19.
Indeno(1,2,3-cd)pyrene	ND		ug/kg	130		23.
Pyrene	ND		ug/kg	97		16.
Biphenyl	ND		ug/kg	370		38.
4-Chloroaniline	ND		ug/kg	160		30.
2-Nitroaniline	ND		ug/kg	160		31.
3-Nitroaniline	ND		ug/kg	160		31.
4-Nitroaniline	ND		ug/kg	160		67.
Dibenzofuran	ND		ug/kg	160		15.
2-Methylnaphthalene	ND		ug/kg	190		20.
1,2,4,5-Tetrachlorobenzene	ND		ug/kg	160		17.
Acetophenone	ND		ug/kg	160		20.
2,4,6-Trichlorophenol	ND		ug/kg	97		31.
p-Chloro-m-cresol	ND		ug/kg	160		24.
2-Chlorophenol	ND		ug/kg	160		19.
2,4-Dichlorophenol	ND		ug/kg	150		26.
2,4-Dimethylphenol	ND		ug/kg	160		54.
2-Nitrophenol	ND		ug/kg	350		61.



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date: 04/16/18

Method Blank Analysis Batch Quality Control

Analytical Method: Analytical Date: 1,8270D 04/11/18 22:56

Analyst:

TT

Extraction Method: EPA 3546

Extraction Date:

04/10/18 17:49

Parameter	Result	Qualifier	Units	RL		MDL	
Semivolatile Organics by GC/MS -	Westboroug	h Lab for s	ample(s):	01-08	Batch:	WG1105283-1	
4-Nitrophenol	ND		ug/kg	230		66.	
2,4-Dinitrophenol	ND		ug/kg	780		76.	-
4,6-Dinitro-o-cresol	ND		ug/kg	420		78.	-
Pentachlorophenol	ND		ug/kg	130		36.	
Phenol	ND		ug/kg	160		24.	
2-Methylphenol	ND		ug/kg	160		25.	
3-Methylphenol/4-Methylphenol	ND		ug/kg	230		25.	
2,4,5-Trichlorophenol	ND		ug/kg	160		31.	
Benzoic Acid	ND		ug/kg	530		160	
Benzyl Alcohol	ND		ug/kg	160		50.	
Carbazole	ND		ug/kg	160		16.	

Tentatively Identified Compounds		
No Tentatively Identified Compounds	ND	ug/kg

Surrogate	%Recovery	Acceptance Qualifier Criteria
2-Fluorophenol	84	25-120
Phenol-d6	88	10-120
Nitrobenzene-d5	72	23-120
2-Fluorobiphenyl	71	30-120
2,4,6-Tribromophenol	82	10-136
4-Terphenyl-d14	76	18-120



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Parameter	LCS %Recovery	Qual	LCSD %Recovery		Recovery Limits	RPD	Qual	RPD Limits
Semivolatile Organics by GC/MS -	Westborough Lab Assoc	iated sample(s):	01-08 Bat	tch: WG110528	3-2 WG11052	83-3		
Acenaphthene	69		67		31-137	3		50
1,2,4-Trichlorobenzene	72		71		38-107	1		50
Hexachlorobenzene	70		69		40-140	1		50
Bis(2-chloroethyl)ether	73		72		40-140	1		50
2-Chloronaphthalene	71		69		40-140	3		50
1,2-Dichlorobenzene	66		67		40-140	2		50
1,3-Dichlorobenzene	66		66		40-140	0		50
1,4-Dichlorobenzene	65		67		28-104	3		50
3,3'-Dichlorobenzidine	50		48		40-140	4		50
2,4-Dinitrotoluene	74		72		40-132	3		50
2,6-Dinitrotoluene	82		82		40-140	0		50
Fluoranthene	68		66		40-140	3		50
4-Chlorophenyl phenyl ether	68		67		40-140	1		50
4-Bromophenyl phenyl ether	66		65		40-140	2		50
Bis(2-chloroisopropyl)ether	82		80		40-140	2		50
Bis(2-chloroethoxy)methane	78		77		40-117	1		50
Hexachlorobutadiene	69		70		40-140	1		50
Hexachlorocyclopentadiene	56		56		40-140	0		50
Hexachloroethane	71		72		40-140	1		50
Isophorone	87		85		40-140	2		50
Naphthalene	67		67		40-140	0		50
Nitrobenzene	77		76		40-140	1		50
NDPA/DPA	69		67		36-157	3		50



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Parameter	LCS %Recovery	Qual	LCSI %Recov		%Recovery Qual Limits	RPD	RPD Qual Limits	
Semivolatile Organics by GC/MS - Westborou	ugh Lab Assoc	iated sample(s):	01-08	Batch:	WG1105283-2 WG110	05283-3		
n-Nitrosodi-n-propylamine	86		85		32-121	1	50	
Bis(2-ethylhexyl)phthalate	82		81		40-140	1	50	
Butyl benzyl phthalate	78		76		40-140	3	50	
Di-n-butylphthalate	74		72		40-140	3	50	
Di-n-octylphthalate	79		77		40-140	3	50	
Diethyl phthalate	73		71		40-140	3	50	
Dimethyl phthalate	76		73		40-140	4	50	
Benzo(a)anthracene	72		70		40-140	3	50	
Benzo(a)pyrene	73		71		40-140	3	50	
Benzo(b)fluoranthene	72		70		40-140	3	50	
Benzo(k)fluoranthene	69		68		40-140	1	50	
Chrysene	71		70		40-140	1	50	
Acenaphthylene	76		75		40-140	1	50	
Anthracene	67		67		40-140	0	50	
Benzo(ghi)perylene	72		72		40-140	0	50	
Fluorene	69		67		40-140	3	50	
Phenanthrene	66		65		40-140	2	50	
Dibenzo(a,h)anthracene	71		71		40-140	0	50	
Indeno(1,2,3-cd)pyrene	74		72		40-140	3	50	
Pyrene	67		66		35-142	2	50	
Biphenyl	70		70		54-104	0	50	
4-Chloroaniline	88		87		40-140	1	50	
2-Nitroaniline	77		77		47-134	0	50	



Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Parameter	LCS %Recovery	Qual	LCSI %Recov		Qual	%Recovery Limits	RPD	Qual	RPD Limits
Semivolatile Organics by GC/MS - Westb	orough Lab Assoc	iated sample(s):	01-08	Batch:	WG110528	3-2 WG1105	283-3		
3-Nitroaniline	58		57			26-129	2		50
4-Nitroaniline	66		66			41-125	0		50
Dibenzofuran	66		66			40-140	0		50
2-Methylnaphthalene	71		69			40-140	3		50
1,2,4,5-Tetrachlorobenzene	70		70			40-117	0		50
Acetophenone	79		79			14-144	0		50
2,4,6-Trichlorophenol	76		75			30-130	1		50
p-Chloro-m-cresol	81		80			26-103	1		50
2-Chlorophenol	71		72			25-102	1		50
2,4-Dichlorophenol	76		75			30-130	1		50
2,4-Dimethylphenol	81		78			30-130	4		50
2-Nitrophenol	77		77			30-130	0		50
4-Nitrophenol	74		74			11-114	0		50
2,4-Dinitrophenol	44		42			4-130	5		50
4,6-Dinitro-o-cresol	70		67			10-130	4		50
Pentachlorophenol	54		51			17-109	6		50
Phenol	70		70			26-90	0		50
2-Methylphenol	79		80			30-130.	1		50
3-Methylphenol/4-Methylphenol	84		83			30-130	1		50
2,4,5-Trichlorophenol	77		78			30-130	1		50
Benzoic Acid	34		29			10-110	16		50
Benzyl Alcohol	78		78			40-140	0		50
Carbazole	67		66			54-128	2		50



Lab Control Sample Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Lab Number:

L1812163

Project Number: Not Specified

Report Date:

04/16/18

LCS LCSD %Recovery RPD Parameter %Recovery Qual %Recovery Qual Limits RPD Qual Limits

Semivolatile Organics by GC/MS - Westborough Lab Associated sample(s): 01-08 Batch: WG1105283-2 WG1105283-3

Surrogate	LCS %Recovery Qual	LCSD %Recovery Qual	Acceptance Criteria
2-Fluorophenol	71	69	25-120
Phenol-d6	80	78	10-120
Nitrobenzene-d5	72	71	23-120
2-Fluorobiphenyl	65	64	30-120
2,4,6-Tribromophenol	77	74	10-136
4-Terphenyl-d14	64	61	18-120

PCBS



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01 Date Collected: 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 15:36 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 96% Cleanup Method: EPA 3660B
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC - V	Vestborough Lab						
Aroclor 1016	ND		ug/kg	33.5	3.80	1	А
Aroclor 1221	ND		ug/kg	33.5	5.10	1	Α
Aroclor 1232	ND		ug/kg	33.5	3.29	1	Α
Aroclor 1242	ND		ug/kg	33.5	4.10	1	Α
Aroclor 1248	ND		ug/kg	33.5	3.76	1	Α
Aroclor 1254	ND		ug/kg	33.5	2.73	1	А
Aroclor 1260	ND		ug/kg	33.5	3.50	1	А
Aroclor 1262	ND		ug/kg	33.5	2.75	1	Α
Aroclor 1268	ND		ug/kg	33.5	2.37	1	Α
PCBs, Total	ND		ug/kg	33.5	2.37	1	Α

Surrogate	% Recovery	Qualifier	Acceptance Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	84		30-150	Α
Decachlorobiphenyl	83		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	88		30-150	В
Decachlorobiphenyl	94		30-150	В

Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 15:49 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 96% Cleanup Date: 04/11/18
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC - We	estborough Lab						
Aroclor 1016	ND		ug/kg	33.3	3.78	1	Α
Aroclor 1221	ND		ug/kg	33.3	5.07	1	Α
Aroclor 1232	ND		ug/kg	33.3	3.28	1	Α
Aroclor 1242	ND		ug/kg	33.3	4.08	1	А
Aroclor 1248	ND		ug/kg	33.3	3.74	1	Α
Aroclor 1254	ND		ug/kg	33.3	2.72	1	Α
Aroclor 1260	ND		ug/kg	33.3	3.48	1	Α
Aroclor 1262	ND		ug/kg	33.3	2.74	1	Α
Aroclor 1268	ND		ug/kg	33.3	2.36	1	Α
PCBs, Total	ND		ug/kg	33.3	2.36	1	Α

Our manufacture	0/ B	0	Acceptance	
Surrogate	% Recovery	Qualifier	Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	81		30-150	Α
Decachlorobiphenyl	80		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	85		30-150	В
Decachlorobiphenyl	90		30-150	В



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03 Date Collected: 04/06/18 08:22

Client ID: SS-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 16:03 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 93% Cleanup Method: EPA 3660B
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC - W	estborough Lab						
Aroclor 1016	ND		ug/kg	34.1	3.86	1	Α
Aroclor 1221	ND		ug/kg	34.1	5.19	1	Α
Aroclor 1232	ND		ug/kg	34.1	3.35	1	Α
Aroclor 1242	ND		ug/kg	34.1	4.17	1	А
Aroclor 1248	ND		ug/kg	34.1	3.82	1	Α
Aroclor 1254	ND		ug/kg	34.1	2.78	1	Α
Aroclor 1260	ND		ug/kg	34.1	3.56	1	Α
Aroclor 1262	ND		ug/kg	34.1	2.80	1	Α
Aroclor 1268	ND		ug/kg	34.1	2.41	1	Α
PCBs, Total	ND		ug/kg	34.1	2.41	1	А

	24.5	o ""	Acceptance	
Surrogate	% Recovery	Qualifier	Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	95		30-150	Α
Decachlorobiphenyl	95		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	103		30-150	В
Decachlorobiphenyl	109		30-150	В



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 16:17 Cleanup Method: EPA 3665A

Analytical Date: 04/11/18 16:17 Cleanup Method: EPA 3665A
Analyst: WR
Percent Solids: 96% Cleanup Date: 04/11/18
Cleanup Method: EPA 3660B
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC - We	stborough Lab						
Aroclor 1016	ND		ug/kg	33.2	3.76	1	А
Aroclor 1221	ND		ug/kg	33.2	5.05	1	Α
Aroclor 1232	ND		ug/kg	33.2	3.26	1	Α
Aroclor 1242	ND		ug/kg	33.2	4.06	1	Α
Aroclor 1248	ND		ug/kg	33.2	3.72	1	Α
Aroclor 1254	ND		ug/kg	33.2	2.70	1	Α
Aroclor 1260	ND		ug/kg	33.2	3.46	1	Α
Aroclor 1262	ND		ug/kg	33.2	2.72	1	Α
Aroclor 1268	ND		ug/kg	33.2	2.35	1	Α
PCBs, Total	ND		ug/kg	33.2	2.35	1	Α

Surrogate	% Recovery	Qualifier	Acceptance Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	78		30-150	Α
Decachlorobiphenyl	77		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	81		30-150	В
Decachlorobiphenyl	84		30-150	В

Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05 Date Collected: 04/06/18 08:29

Client ID: SS-5 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Total Top. Not opcomed

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 16:31 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 96% Cleanup Date: 04/11/18
Cleanup Date: 04/11/18
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC	- Westborough Lab						
Aroclor 1016	ND		ug/kg	34.2	3.88	1	Α
Aroclor 1221	ND		ug/kg	34.2	5.21	1	Α
Aroclor 1232	ND		ug/kg	34.2	3.37	1	Α
Aroclor 1242	ND		ug/kg	34.2	4.19	1	Α
Aroclor 1248	ND		ug/kg	34.2	3.84	1	Α
Aroclor 1254	ND		ug/kg	34.2	2.79	1	Α
Aroclor 1260	ND		ug/kg	34.2	3.57	1	Α
Aroclor 1262	ND		ug/kg	34.2	2.81	1	Α
Aroclor 1268	ND		ug/kg	34.2	2.42	1	Α
PCBs, Total	ND		ug/kg	34.2	2.42	1	Α

Surremente	9/ Bassyami	Ouglifier	Acceptance	Caluman
Surrogate	% Recovery	Qualifier	Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	88		30-150	Α
Decachlorobiphenyl	90		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	95		30-150	В
Decachlorobiphenyl	101		30-150	В



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: B-1 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 16:45 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 93% Cleanup Date: 04/11/18
Cleanup Date: 04/11/18
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC	- Westborough Lab						
Aroclor 1016	ND		ug/kg	34.8	3.94	1	Α
Aroclor 1221	ND		ug/kg	34.8	5.29	1	Α
Aroclor 1232	ND		ug/kg	34.8	3.42	1	Α
Aroclor 1242	ND		ug/kg	34.8	4.26	1	Α
Aroclor 1248	ND		ug/kg	34.8	3.90	1	А
Aroclor 1254	ND		ug/kg	34.8	2.84	1	Α
Aroclor 1260	ND		ug/kg	34.8	3.63	1	Α
Aroclor 1262	ND		ug/kg	34.8	2.86	1	Α
Aroclor 1268	ND		ug/kg	34.8	2.46	1	Α
PCBs, Total	ND		ug/kg	34.8	2.46	1	Α

Surrogate	% Recovery	Qualifier	Acceptance Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	61		30-150	Α
Decachlorobiphenyl	55		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	62		30-150	В
Decachlorobiphenyl	63		30-150	В

Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 16:59 Cleanup Method: EPA 3665A

Analytical Date: 04/11/18 16:59

Analyst: WR

Percent Solids: 87%

Cleanup Method: EPA 3665A

Cleanup Date: 04/11/18

Cleanup Method: EPA 3660B

Cleanup Method: EPA 3660B

Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC - V	Westborough Lab						
Aroclor 1016	ND		ug/kg	38.1	4.32	1	Α
Aroclor 1221	ND		ug/kg	38.1	5.80	1	Α
Aroclor 1232	ND		ug/kg	38.1	3.75	1	Α
Aroclor 1242	ND		ug/kg	38.1	4.67	1	Α
Aroclor 1248	ND		ug/kg	38.1	4.28	1	Α
Aroclor 1254	ND		ug/kg	38.1	3.11	1	Α
Aroclor 1260	ND		ug/kg	38.1	3.98	1	Α
Aroclor 1262	ND		ug/kg	38.1	3.13	1	Α
Aroclor 1268	ND		ug/kg	38.1	2.70	1	Α
PCBs, Total	ND		ug/kg	38.1	2.70	1	Α

			Acceptance	
Surrogate	% Recovery	Qualifier	Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	66		30-150	Α
Decachlorobiphenyl	54		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	70		30-150	В
Decachlorobiphenyl	64		30-150	В



Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08 Date Collected: 04/06/18 10:10

Client ID: B-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Extraction Method: EPA 3546
Analytical Method: 1,8082A Extraction Date: 04/10/18 21:36
Analytical Date: 04/11/18 17:13 Cleanup Method: EPA 3665A

Analyst: WR Cleanup Date: 04/11/18
Percent Solids: 91% Cleanup Method: EPA 3660B
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Column
Polychlorinated Biphenyls by GC -	· Westborough Lab						
Aroclor 1016	ND		ug/kg	36.0	4.08	1	Α
Aroclor 1221	ND		ug/kg	36.0	5.48	1	Α
Aroclor 1232	ND		ug/kg	36.0	3.54	1	Α
Aroclor 1242	ND		ug/kg	36.0	4.41	1	Α
Aroclor 1248	ND		ug/kg	36.0	4.04	1	Α
Aroclor 1254	ND		ug/kg	36.0	2.94	1	Α
Aroclor 1260	ND		ug/kg	36.0	3.76	1	Α
Aroclor 1262	ND		ug/kg	36.0	2.96	1	Α
Aroclor 1268	ND		ug/kg	36.0	2.55	1	Α
PCBs, Total	ND		ug/kg	36.0	2.55	1	Α

Surrogate	% Recovery	Qualifier	Acceptance Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	85		30-150	Α
Decachlorobiphenyl	78		30-150	Α
2,4,5,6-Tetrachloro-m-xylene	86		30-150	В
Decachlorobiphenyl	94		30-150	В

L1812163

Lab Number:

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified Report Date: 04/16/18

Method Blank Analysis
Batch Quality Control

Analytical Method: 1,8082A Analytical Date: 04/11/18 12:25

Analyst: WR

Extraction Method: EPA 3546
Extraction Date: 04/10/18 21:36
Cleanup Method: EPA 3665A
Cleanup Date: 04/11/18
Cleanup Method: EPA 3660B
Cleanup Date: 04/11/18

Parameter	Result	Qualifier	Units	RL		MDL	Column
Polychlorinated Biphenyls by GC - V	Vestborough	n Lab for s	ample(s):	01-08	Batch:	WG110	05337-1
Aroclor 1016	ND		ug/kg	31.4		3.57	А
Aroclor 1221	ND		ug/kg	31.4		4.79	А
Aroclor 1232	ND		ug/kg	31.4		3.09	Α
Aroclor 1242	ND		ug/kg	31.4		3.85	Α
Aroclor 1248	ND		ug/kg	31.4		3.53	Α
Aroclor 1254	ND		ug/kg	31.4		2.57	Α
Aroclor 1260	ND		ug/kg	31.4		3.28	Α
Aroclor 1262	ND		ug/kg	31.4		2.58	Α
Aroclor 1268	ND		ug/kg	31.4		2.23	Α
PCBs, Total	ND		ug/kg	31.4		2.23	Α

		Acceptance	e
Surrogate	%Recovery Qualifie	r Criteria	Column
2,4,5,6-Tetrachloro-m-xylene	102	30-150	Α
Decachlorobiphenyl	85	30-150	Α
2,4,5,6-Tetrachloro-m-xylene	110	30-150	В
Decachlorobiphenyl	77	30-150	В



Lab Control Sample Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

	LCS		LCSD		%Recovery				
Parameter	%Recovery	Qual	%Recovery	Qual Limits		RPD	Qual	Limits	Column
Polychlorinated Biphenyls by GC - Westbor	ough Lab Associa	ated sample(s)	: 01-08 Batch:	WG1105	337-2 WG110533	37-3			
Aroclor 1016	107		104		40-140	3		50	Α
Aroclor 1260	114		108		40-140	5		50	А

Surrogate	LCS %Recovery Qu	LCSD al %Recovery Qual	Acceptance Criteria Column
2,4,5,6-Tetrachloro-m-xylene	100	100	30-150 A
Decachlorobiphenyl	88	85	30-150 A
2,4,5,6-Tetrachloro-m-xylene	100	105	30-150 B
Decachlorobiphenyl	73	75	30-150 B



METALS



Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-01

Client ID: SS-1 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Percent Solids: 96%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 27600 mg/kg 8.09 2.18 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC 5.59 4.04 0.307 2 1,6010C LC Antimony, Total mg/kg 04/10/18 19:50 04/11/18 19:52 EPA 3050B 2 Arsenic, Total 1.30 mg/kg 0.809 0.168 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC 2 Barium, Total 701 0.809 0.141 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC mg/kg 0.380 J 0.404 0.027 2 1,6010C LC Beryllium, Total mg/kg 04/10/18 19:50 04/11/18 19:52 EPA 3050B 2 1.80 0.079 1,6010C LC Cadmium, Total mg/kg 0.809 04/10/18 19:50 04/11/18 19:52 EPA 3050B Calcium, Total 36300 8.09 2.83 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C mg/kg LC Chromium, Total 2 1,6010C LC 272 0.809 0.078 04/10/18 19:50 04/11/18 19:52 EPA 3050B mg/kg 2 1,6010C LC Cobalt, Total 12.1 mg/kg 1.62 0.134 04/10/18 19:50 04/11/18 19:52 EPA 3050B 2 1,6010C Copper, Total 961 mg/kg 0.809 0.209 04/10/18 19:50 04/11/18 19:52 EPA 3050B LC 7.31 20 1,6010C LC Iron, Total 166000 40.4 04/10/18 19:50 04/12/18 16:44 EPA 3050B mg/kg 2 1,6010C Lead, Total 206 mg/kg 4.04 0.217 04/10/18 19:50 04/11/18 19:52 EPA 3050B LC Magnesium, Total 5070 8.09 1.25 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC mg/kg 0.809 0.129 2 1,6010C LC Manganese, Total 1170 mg/kg 04/10/18 19:50 04/11/18 19:52 EPA 3050B Mercury, Total ND mg/kg 0.065 0.014 1 04/11/18 07:00 04/11/18 18:21 EPA 7471B 1,7471B EΑ Nickel, Total 28.7 2.02 0.196 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC mg/kg 4450 2 1,6010C LC Potassium, Total mg/kg 202 11.6 04/10/18 19:50 04/11/18 19:52 EPA 3050B Selenium, Total ND mg/kg 1.62 0.209 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC Silver, Total 1.24 mg/kg 0.809 0.229 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC Sodium, Total 32400 mg/kg 162 2.55 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC Thallium, Total 0.348 J mg/kg 1.62 0.255 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC Vanadium, Total 14.1 0.809 2 04/10/18 19:50 04/11/18 19:52 EPA 3050B 1,6010C LC mg/kg 0.164 2 1,6010C LC 733 4.04 0.237 Zinc, Total mg/kg 04/10/18 19:50 04/11/18 19:52 EPA 3050B



Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02

Client ID: SS-2 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil
Percent Solids: 96%

Percent Solids: Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 27600 mg/kg 8.32 2.25 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC 8.51 2 1,6010C LC Antimony, Total mg/kg 4.16 0.316 04/10/18 19:50 04/11/18 19:57 EPA 3050B Arsenic, Total 1.02 mg/kg 0.832 0.173 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC Barium, Total 702 0.832 0.145 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC mg/kg 0.416 0.028 2 1,6010C LC Beryllium, Total mg/kg 0.416 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1.89 0.082 2 1,6010C LC Cadmium, Total mg/kg 0.832 04/10/18 19:50 04/11/18 19:57 EPA 3050B Calcium, Total 36100 8.32 2.91 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C mg/kg LC 2 1,6010C LC Chromium, Total 286 0.832 0.080 04/10/18 19:50 04/11/18 19:57 EPA 3050B mg/kg 2 9.03 1,6010C LC Cobalt, Total mg/kg 1.66 0.138 04/10/18 19:50 04/11/18 19:57 EPA 3050B 2 1,6010C Copper, Total 1010 mg/kg 0.832 0.215 04/10/18 19:50 04/11/18 19:57 EPA 3050B LC 7.52 20 1,6010C LC Iron, Total 163000 41.6 04/10/18 19:50 04/12/18 16:49 EPA 3050B mg/kg 446 0.223 2 1,6010C Lead, Total mg/kg 4.16 04/10/18 19:50 04/11/18 19:57 EPA 3050B LC Magnesium, Total 5150 8.32 1.28 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC mg/kg 0.832 0.132 2 1,6010C LC Manganese, Total 1190 mg/kg 04/10/18 19:50 04/11/18 19:57 EPA 3050B Mercury, Total ND mg/kg 0.066 0.014 1 04/11/18 07:00 04/11/18 18:23 EPA 7471B 1,7471B EΑ Nickel, Total 21.8 2.08 0.201 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC mg/kg 4410 208 12.0 2 1,6010C LC Potassium, Total mg/kg 04/10/18 19:50 04/11/18 19:57 EPA 3050B Selenium, Total ND mg/kg 1.66 0.215 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC Silver, Total 1.45 mg/kg 0.832 0.236 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC Sodium, Total 33900 mg/kg 166 2.62 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC Thallium, Total 0.374 J mg/kg 1.66 0.262 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC Vanadium, Total 14.4 0.832 0.169 2 04/10/18 19:50 04/11/18 19:57 EPA 3050B 1,6010C LC mg/kg 2 1,6010C LC 740 0.244 Zinc, Total mg/kg 4.16 04/10/18 19:50 04/11/18 19:57 EPA 3050B



Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03

Client ID: SS-3 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil
Percent Solids: 93%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 24600 mg/kg 8.32 2.25 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC 4.97 2 1,6010C LC Antimony, Total mg/kg 4.16 0.316 04/10/18 19:50 04/11/18 20:02 EPA 3050B Arsenic, Total 1.93 mg/kg 0.832 0.173 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC Barium, Total 617 0.832 0.145 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC mg/kg 0.333 J 0.028 2 1,6010C LC Beryllium, Total mg/kg 0.416 04/10/18 19:50 04/11/18 20:02 EPA 3050B 0.082 2 1,6010C LC Cadmium, Total 2.14 mg/kg 0.832 04/10/18 19:50 04/11/18 20:02 EPA 3050B Calcium, Total 32300 8.32 2.91 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C mg/kg LC Chromium, Total 2 1,6010C LC 308 0.832 0.080 04/10/18 19:50 04/11/18 20:02 EPA 3050B mg/kg 2 1,6010C LC Cobalt, Total 14.3 mg/kg 1.66 0.138 04/10/18 19:50 04/11/18 20:02 EPA 3050B 2 1,6010C Copper, Total 1110 mg/kg 0.832 0.215 04/10/18 19:50 04/11/18 20:02 EPA 3050B LC 7.52 20 1,6010C LC Iron, Total 168000 41.6 04/10/18 19:50 04/12/18 16:54 EPA 3050B mg/kg 0.223 2 Lead, Total 184 mg/kg 4.16 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC Magnesium, Total 4680 8.32 1.28 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC mg/kg 946 0.832 0.132 2 1,6010C LC Manganese, Total mg/kg 04/10/18 19:50 04/11/18 20:02 EPA 3050B Mercury, Total ND mg/kg 0.068 0.014 1 04/11/18 07:00 04/11/18 18:25 EPA 7471B 1,7471B EΑ Nickel, Total 45.0 2.08 0.201 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC mg/kg 3350 208 12.0 2 1,6010C LC Potassium, Total mg/kg 04/10/18 19:50 04/11/18 20:02 EPA 3050B Selenium, Total ND mg/kg 1.66 0.215 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC Silver, Total 1.80 mg/kg 0.832 0.236 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC Sodium, Total 27700 mg/kg 166 2.62 2 04/10/18 19:50 04/11/18 20:02 EPA 3050B 1,6010C LC



1,6010C

1,6010C

1,6010C

LC

LC

LC

Thallium, Total

Vanadium, Total

Zinc, Total

ND

12.8

546

mg/kg

mg/kg

mg/kg

1.66

0.832

4.16

0.262

0.169

0.244

2

2

2

04/10/18 19:50 04/11/18 20:02 EPA 3050B

04/10/18 19:50 04/11/18 20:02 EPA 3050B

04/10/18 19:50 04/11/18 20:02 EPA 3050B

Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04

Client ID: SS-4 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil

96% Percent Solids: Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 26600 mg/kg 8.16 2.20 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC 5.77 4.08 2 1,6010C LC Antimony, Total mg/kg 0.310 04/10/18 19:50 04/11/18 20:07 EPA 3050B Arsenic, Total 1.28 mg/kg 0.816 0.170 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC Barium, Total 608 0.816 0.142 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC mg/kg 0.432 0.027 2 1,6010C LC Beryllium, Total mg/kg 0.408 04/10/18 19:50 04/11/18 20:07 EPA 3050B 2.48 0.080 2 1,6010C LC Cadmium, Total mg/kg 0.816 04/10/18 19:50 04/11/18 20:07 EPA 3050B Calcium, Total 34700 8.16 2.86 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C mg/kg LC 2 1,6010C LC Chromium, Total 340 0.816 0.078 04/10/18 19:50 04/11/18 20:07 EPA 3050B mg/kg 2 1,6010C LC Cobalt, Total 11.5 mg/kg 1.63 0.135 04/10/18 19:50 04/11/18 20:07 EPA 3050B 2 1,6010C Copper, Total 1130 mg/kg 0.816 0.210 04/10/18 19:50 04/11/18 20:07 EPA 3050B LC 7.37 20 1,6010C LC Iron, Total 198000 40.8 04/10/18 19:50 04/12/18 16:58 EPA 3050B mg/kg 308 2 1,6010C Lead, Total mg/kg 4.08 0.219 04/10/18 19:50 04/11/18 20:07 EPA 3050B LC Magnesium, Total 5260 8.16 1.26 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC mg/kg 986 0.816 0.130 2 1,6010C LC Manganese, Total mg/kg 04/10/18 19:50 04/11/18 20:07 EPA 3050B Mercury, Total ND mg/kg 0.066 0.014 1 04/11/18 07:00 04/11/18 18:27 EPA 7471B 1,7471B EΑ Nickel, Total 28.5 2.04 0.198 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC mg/kg 4090 2 1,6010C LC Potassium, Total mg/kg 204 11.8 04/10/18 19:50 04/11/18 20:07 EPA 3050B Selenium, Total 3.35 J mg/kg 4.08 0.526 5 04/10/18 19:50 04/13/18 12:10 EPA 3050B 1,6010C LC Silver, Total 1.62 mg/kg 0.816 0.231 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC Sodium, Total 31200 mg/kg 163 2.57 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC Thallium, Total 0.604 J mg/kg 1.63 0.257 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC Vanadium, Total 17.6 0.816 2 04/10/18 19:50 04/11/18 20:07 EPA 3050B 1,6010C LC mg/kg 0.166

2

04/10/18 19:50 04/11/18 20:07 EPA 3050B

0.239

4.08

mg/kg



1,6010C

LC

Zinc, Total

885

Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05

Client ID: SS-5 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Percent Solids: 96%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 26200 mg/kg 8.05 2.17 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC 4.74 4.02 2 1,6010C LC Antimony, Total mg/kg 0.306 04/10/18 19:50 04/11/18 20:12 EPA 3050B Arsenic, Total 1.33 mg/kg 0.805 0.167 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC 2 Barium, Total 779 0.805 0.140 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC mg/kg 0.346 J 0.402 0.027 2 1,6010C LC Beryllium, Total mg/kg 04/10/18 19:50 04/11/18 20:12 EPA 3050B 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 2.45 0.079 1,6010C LC Cadmium, Total mg/kg 0.805 Calcium, Total 36100 8.05 2.82 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C mg/kg LC 0.077 2 1,6010C LC Chromium, Total 322 0.805 04/10/18 19:50 04/11/18 20:12 EPA 3050B mg/kg 2 12.7 1,6010C LC Cobalt, Total mg/kg 1.61 0.134 04/10/18 19:50 04/11/18 20:12 EPA 3050B 2 1,6010C Copper, Total 1050 mg/kg 0.805 0.208 04/10/18 19:50 04/11/18 20:12 EPA 3050B LC 185000 7.26 20 1,6010C LC Iron, Total 40.2 04/10/18 19:50 04/12/18 17:03 EPA 3050B mg/kg 2 1,6010C Lead, Total 214 mg/kg 4.02 0.216 04/10/18 19:50 04/11/18 20:12 EPA 3050B LC Magnesium, Total 5310 8.05 1.24 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC mg/kg 0.805 0.128 2 1,6010C LC Manganese, Total 1110 mg/kg 04/10/18 19:50 04/11/18 20:12 EPA 3050B J Mercury, Total 0.022 mg/kg 0.066 0.014 1 04/11/18 07:00 04/11/18 18:32 EPA 7471B 1,7471B EΑ Nickel, Total 29.2 2.01 0.195 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC mg/kg 4090 2 1,6010C LC Potassium, Total mg/kg 201 11.6 04/10/18 19:50 04/11/18 20:12 EPA 3050B Selenium, Total 3.10 J mg/kg 4.02 0.519 5 04/10/18 19:50 04/13/18 12:14 EPA 3050B 1,6010C LC Silver, Total 2.20 mg/kg 0.805 0.228 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC Sodium, Total 33800 mg/kg 161 2.53 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC Thallium, Total 0.628 J mg/kg 1.61 0.253 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC Vanadium, Total 17.5 0.805 2 04/10/18 19:50 04/11/18 20:12 EPA 3050B 1,6010C LC mg/kg 0.163 2 1,6010C LC 919 4.02 0.236 Zinc, Total mg/kg 04/10/18 19:50 04/11/18 20:12 EPA 3050B



04/06/18 09:30

Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06

Client ID: B-1 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil
Percent Solids: 93%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 5430 mg/kg 8.26 2.23 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC J 2 1,6010C LC Antimony, Total 0.719 mg/kg 4.13 0.314 04/10/18 19:50 04/11/18 20:17 EPA 3050B Arsenic, Total 5.19 mg/kg 0.826 0.172 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC 2 Barium, Total 25.5 0.826 0.144 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC mg/kg 0.289 J 0.027 2 1,6010C LC Beryllium, Total mg/kg 0.413 04/10/18 19:50 04/11/18 20:17 EPA 3050B J 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 0.081 1,6010C LC Cadmium, Total 0.289 mg/kg 0.826 Calcium, Total 41400 8.26 2.89 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C mg/kg LC Chromium, Total 2 1,6010C LC 10.1 0.826 0.079 04/10/18 19:50 04/11/18 20:17 EPA 3050B mg/kg 2 1,6010C LC Cobalt, Total 6.57 mg/kg 1.65 0.137 04/10/18 19:50 04/11/18 20:17 EPA 3050B 2 1,6010C Copper, Total 24.4 mg/kg 0.826 0.213 04/10/18 19:50 04/11/18 20:17 EPA 3050B LC 14600 2 1,6010C LC Iron, Total 0.746 04/10/18 19:50 04/11/18 20:17 EPA 3050B mg/kg 4.13 8.31 2 1,6010C Lead, Total mg/kg 4.13 0.222 04/10/18 19:50 04/11/18 20:17 EPA 3050B LC Magnesium, Total 10800 8.26 1.27 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC mg/kg 0.826 0.131 2 1,6010C LC Manganese, Total 310 mg/kg 04/10/18 19:50 04/11/18 20:17 EPA 3050B J Mercury, Total 0.026 mg/kg 0.068 0.014 1 04/11/18 07:00 04/11/18 18:34 EPA 7471B 1,7471B EΑ Nickel, Total 18.5 2.07 0.200 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC mg/kg 460 2 1,6010C LC Potassium, Total mg/kg 207 11.9 04/10/18 19:50 04/11/18 20:17 EPA 3050B Selenium, Total ND mg/kg 1.65 0.213 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC Silver, Total ND mg/kg 0.826 0.234 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC J Sodium, Total 124 mg/kg 165 2.60 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC Thallium, Total ND mg/kg 1.65 0.260 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC Vanadium, Total 0.826 2 04/10/18 19:50 04/11/18 20:17 EPA 3050B 1,6010C LC 11.0 mg/kg 0.168 2 1,6010C LC 55.4 0.242 Zinc, Total mg/kg 4.13 04/10/18 19:50 04/11/18 20:17 EPA 3050B



04/06/18 09:50

Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07

Client ID: B-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil Percent Solids: 87%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 9260 mg/kg 8.72 2.35 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC 0.933 J 4.36 0.331 2 1,6010C LC Antimony, Total mg/kg 04/10/18 19:50 04/11/18 20:22 EPA 3050B 2 Arsenic, Total 2.82 mg/kg 0.872 0.181 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC 2 Barium, Total 36.7 0.872 0.152 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC mg/kg 0.410 J 0.436 0.029 2 1,6010C LC Beryllium, Total mg/kg 04/10/18 19:50 04/11/18 20:22 EPA 3050B J 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 0.086 1,6010C LC Cadmium, Total 0.288 mg/kg 0.872 04/10/18 19:50 04/12/18 17:08 EPA 3050B Calcium, Total 88500 87.2 30.5 20 1,6010C mg/kg LC 2 1,6010C LC Chromium, Total 17.0 0.872 0.084 04/10/18 19:50 04/11/18 20:22 EPA 3050B mg/kg 2 1,6010C LC Cobalt, Total 10.8 mg/kg 1.74 0.145 04/10/18 19:50 04/11/18 20:22 EPA 3050B 2 1,6010C Copper, Total 48.5 mg/kg 0.872 0.225 04/10/18 19:50 04/11/18 20:22 EPA 3050B LC 0.788 2 1,6010C LC Iron, Total 23000 4.36 04/10/18 19:50 04/11/18 20:22 EPA 3050B mg/kg 8.04 2 1,6010C Lead, Total mg/kg 4.36 0.234 04/10/18 19:50 04/11/18 20:22 EPA 3050B LC Magnesium, Total 5920 8.72 1.34 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC mg/kg 358 0.872 2 1,6010C LC Manganese, Total mg/kg 0.139 04/10/18 19:50 04/11/18 20:22 EPA 3050B J Mercury, Total 0.039 mg/kg 0.072 0.015 1 04/11/18 07:00 04/11/18 18:36 EPA 7471B 1,7471B EΑ Nickel, Total 35.4 2.18 0.211 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC mg/kg 560 2 1,6010C LC Potassium, Total mg/kg 218 12.6 04/10/18 19:50 04/11/18 20:22 EPA 3050B Selenium, Total ND mg/kg 1.74 0.225 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC Silver, Total ND mg/kg 0.872 0.247 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC J Sodium, Total 127 mg/kg 174 2.75 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC Thallium, Total ND mg/kg 1.74 0.275 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC Vanadium, Total 13.9 0.872 2 04/10/18 19:50 04/11/18 20:22 EPA 3050B 1,6010C LC mg/kg 0.177 2 1,6010C LC 50.6 4.36 0.256 Zinc, Total mg/kg 04/10/18 19:50 04/11/18 20:22 EPA 3050B



04/06/18 10:10

Date Collected:

Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08

Client ID: B-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil
Percent Solids: 91%

Dilution Date Date Prep **Analytical** Method Qualifier Factor **Prepared** Analyzed Method **Parameter** Result Units RL MDL Analyst Total Metals - Mansfield Lab Aluminum, Total 5630 mg/kg 8.48 2.29 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC J 4.24 0.322 2 1,6010C LC Antimony, Total 1.69 mg/kg 04/10/18 19:50 04/11/18 20:50 EPA 3050B Arsenic, Total 11.3 mg/kg 0.848 0.176 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC 2 Barium, Total 21.1 0.848 0.148 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC mg/kg 0.271 J 0.424 0.028 2 1,6010C LC Beryllium, Total mg/kg 04/10/18 19:50 04/11/18 20:50 EPA 3050B J 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 0.083 1,6010C LC Cadmium, Total 0.271 mg/kg 0.848 04/10/18 19:50 04/12/18 17:42 EPA 3050B Calcium, Total 104000 84.8 29.7 20 1,6010C mg/kg LC Chromium, Total 2 1,6010C LC 11.1 0.848 0.081 04/10/18 19:50 04/11/18 20:50 EPA 3050B mg/kg 2 26.4 1,6010C LC Cobalt, Total mg/kg 1.70 0.141 04/10/18 19:50 04/11/18 20:50 EPA 3050B 2 1,6010C Copper, Total 48.3 mg/kg 0.848 0.219 04/10/18 19:50 04/11/18 20:50 EPA 3050B LC 0.766 2 1,6010C LC Iron, Total 19100 4.24 04/10/18 19:50 04/11/18 20:50 EPA 3050B mg/kg 2 1,6010C Lead, Total 26.1 mg/kg 4.24 0.227 04/10/18 19:50 04/11/18 20:50 EPA 3050B LC Magnesium, Total 5880 8.48 1.31 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC mg/kg 300 0.848 0.135 2 1,6010C LC Manganese, Total mg/kg 04/10/18 19:50 04/11/18 20:50 EPA 3050B J Mercury, Total 0.035 mg/kg 0.070 0.015 1 04/11/18 07:00 04/11/18 18:38 EPA 7471B 1,7471B EΑ Nickel, Total 45.7 2.12 0.205 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC mg/kg 523 2 1,6010C LC Potassium, Total mg/kg 212 12.2 04/10/18 19:50 04/11/18 20:50 EPA 3050B Selenium, Total 0.568 J mg/kg 1.70 0.219 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC Silver, Total ND mg/kg 0.848 0.240 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC J Sodium, Total 116 mg/kg 170 2.67 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC Thallium, Total ND mg/kg 1.70 0.267 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC Vanadium, Total 0.848 2 04/10/18 19:50 04/11/18 20:50 EPA 3050B 1,6010C LC 11.6 mg/kg 0.172 2 1,6010C LC 32.5 4.24 0.248 Zinc, Total mg/kg 04/10/18 19:50 04/11/18 20:50 EPA 3050B



L1812163

04/16/18

Lab Number:

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified **Report Date:**

Method Blank Analysis

Batch Quality Control

Parameter	Result Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	
Total Metals - Mansfield	Lab for sample(s):	01-08 B	atch: W	G110528	35-1				
Aluminum, Total	ND	mg/kg	4.00	1.08	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Antimony, Total	ND	mg/kg	2.00	0.152	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Arsenic, Total	ND	mg/kg	0.400	0.083	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Barium, Total	ND	mg/kg	0.400	0.070	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Beryllium, Total	ND	mg/kg	0.200	0.013	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Cadmium, Total	ND	mg/kg	0.400	0.039	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Calcium, Total	ND	mg/kg	4.00	1.40	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Chromium, Total	ND	mg/kg	0.400	0.038	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Cobalt, Total	ND	mg/kg	0.800	0.066	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Copper, Total	ND	mg/kg	0.400	0.103	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Iron, Total	ND	mg/kg	2.00	0.361	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Lead, Total	ND	mg/kg	2.00	0.107	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Magnesium, Total	ND	mg/kg	4.00	0.616	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Manganese, Total	ND	mg/kg	0.400	0.064	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Nickel, Total	ND	mg/kg	1.00	0.097	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Potassium, Total	ND	mg/kg	100	5.76	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Selenium, Total	ND	mg/kg	0.800	0.103	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Silver, Total	ND	mg/kg	0.400	0.113	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Sodium, Total	ND	mg/kg	80.0	1.26	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Thallium, Total	ND	mg/kg	0.800	0.126	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Vanadium, Total	ND	mg/kg	0.400	0.081	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC
Zinc, Total	ND	mg/kg	2.00	0.117	1	04/10/18 19:50	04/11/18 18:46	1,6010C	LC

Prep Information

Digestion Method: EPA 3050B

Parameter	Result Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	
Total Metals - Mar	nsfield Lab for sample(s):	01-08 B	atch: W	G11053	88-1				
Mercury, Total	ND	mg/kg	0.083	0.018	1	04/11/18 07:00	04/11/18 18:10	1,7471B	EA



Project Name: SOUTHWESTERN PHASE 2 **Lab Number:** L1812163

Project Number: Not Specified Report Date: 04/16/18

Method Blank Analysis Batch Quality Control

Prep Information

Digestion Method: EPA 7471B



Lab Control Sample Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Report Date: 04/16/18

arameter	LCS %Recovery	y Qual	LCSD %Recover	У Qual	%Recovery Limits	RPD	Qual	RPD Limits
otal Metals - Mansfield Lab Associated sample	e(s): 01-08 E	Batch: WG110)5285-2 SR	M Lot Number:	D098-540			
Aluminum, Total	59		-		47-153	-		
Antimony, Total	148		-		6-194	-		
Arsenic, Total	90		-		83-117	-		
Barium, Total	88		-		82-118	-		
Beryllium, Total	84		-		83-117	-		
Cadmium, Total	99		-		82-117	-		
Calcium, Total	91		-		81-118	-		
Chromium, Total	88		-		83-119	-		
Cobalt, Total	101		-		84-116	-		
Copper, Total	90		-		84-116	-		
Iron, Total	75		-		60-140	-		
Lead, Total	86		-		82-117	-		
Magnesium, Total	77		-		76-124	-		
Manganese, Total	89		-		82-118	-		
Nickel, Total	99		-		82-117	-		
Potassium, Total	70		-		69-131	-		
Selenium, Total	95		-		78-121	-		
Silver, Total	89		-		80-120	-		
Sodium, Total	93		-		74-126	-		
Thallium, Total	85		-		80-119	-		
Vanadium, Total	85		-		79-121	-		

Lab Control Sample Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Parameter	LCS %Recov	_ 	SD % covery	%Recovery Limits	RPD	RPD Limits
Total Metals - Mansfield Lab	Associated sample(s): 01-08	Batch: WG1105285-2	SRM Lot Number: D09	98-540		
Zinc, Total	87		-	81-119	-	
Total Metals - Mansfield Lab	Associated sample(s): 01-08	Batch: WG1105388-2	SRM Lot Number: D09	98-540		
Mercury, Total	112		-	50-149	-	



Matrix Spike Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number: L1812163

Report Date: 04/16/18

<u>Parameter</u>	Native Sample	MS Added	MS Found	MS %Recovery	Qual	MSD Found	MSD %Recovery Q	Recovery ual Limits	RPD Qual	RPD Limits
Total Metals - Mansfield L	ab Associated sar	mple(s): 01-08	QC Ba	tch ID: WG110	5285-3	QC San	nple: L1810464-05	Client ID: MS	S Sample	
Aluminum, Total	3720	183	3830	60	Q	-	-	75-125	-	20
Antimony, Total	0.758J	45.7	43.4	95		-	-	75-125	-	20
Arsenic, Total	2.40	11	13.1	98		-	-	75-125	-	20
Barium, Total	22.8	183	192	93		-	-	75-125	-	20
Beryllium, Total	0.238J	4.57	4.48	98		-	-	75-125	-	20
Cadmium, Total	0.141J	4.66	4.72	101		-	-	75-125	-	20
Calcium, Total	1540	913	2340	88		-	-	75-125	-	20
Chromium, Total	9.51	18.3	26.4	92		-	-	75-125	-	20
Cobalt, Total	2.69	45.7	44.6	92		-	-	75-125	-	20
Copper, Total	13.2	22.8	35.8	99		-	-	75-125	-	20
Iron, Total	7020	91.3	6930	0	Q	-	-	75-125	-	20
Lead, Total	30.4	46.6	70.7	86		-	-	75-125	-	20
Magnesium, Total	2050	913	2330	31	Q	-	-	75-125	-	20
Manganese, Total	95.0	45.7	116	46	Q	-	-	75-125	-	20
Nickel, Total	5.61	45.7	46.8	90		-	-	75-125	-	20
Potassium, Total	969.	913	1700	80		-	-	75-125	-	20
Selenium, Total	ND	11	10.5	96		-	-	75-125	-	20
Silver, Total	ND	27.4	26.6	97		-	-	75-125	-	20
Sodium, Total	169.J	913	1070	117		-	-	75-125	-	20
Thallium, Total	ND	11	9.77	89		-	-	75-125	-	20
Vanadium, Total	16.3	45.7	58.2	92		-	-	75-125	-	20

Matrix Spike Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Parameter	Native Sample	MS Added	MS Found	MS %Recovery	MSD Found	MSD %Recovery	Recovery Limits	RPD	RPD Limits
Total Metals - Mansfield Lab	Associated sam	ple(s): 01-08	QC Bat	ch ID: WG1105285-3	QC Sam	ple: L1810464-05	Client ID: MS	Sample	
Zinc, Total	136.	45.7	170	74 Q	-	-	75-125	-	20
Total Metals - Mansfield Lab	Associated sam	ple(s): 01-08	QC Bat	tch ID: WG1105388-3	3 WG11053	88-4 QC Sample:	L1812198-11	Client ID:	MS Sample
Mercury, Total	0.188	0.14	0.366	127 Q	0.290	72	Q 80-120	23	Q 20



Lab Duplicate Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified Lab Number: L1812163

04/16/18 Report Date:

arameter	Native Sample	Duplicate Sample	Units	RPD	Qual	RPD Limits
otal Metals - Mansfield Lab Associated sample(s): 01-0	08 QC Batch ID: \	WG1105285-4 QC Sample:	L1810464-05	Client ID:	DUP Sam	ple
Aluminum, Total	3720	3980	mg/kg	7		20
Antimony, Total	0.758J	0.618J	mg/kg	NC		20
Arsenic, Total	2.40	2.79	mg/kg	15		20
Barium, Total	22.8	19.0	mg/kg	18		20
Beryllium, Total	0.238J	0.218J	mg/kg	NC		20
Cadmium, Total	0.141J	0.118J	mg/kg	NC		20
Calcium, Total	1540	1330	mg/kg	15		20
Chromium, Total	9.51	11.5	mg/kg	19		20
Cobalt, Total	2.69	3.94	mg/kg	38	Q	20
Copper, Total	13.2	18.4	mg/kg	33	Q	20
Iron, Total	7020	8410	mg/kg	18		20
Lead, Total	30.4	42.0	mg/kg	32	Q	20
Magnesium, Total	2050	2250	mg/kg	9		20
Manganese, Total	95.0	97.8	mg/kg	3		20
Nickel, Total	5.61	6.03	mg/kg	7		20
Potassium, Total	969.	1310	mg/kg	30	Q	20
Selenium, Total	ND	ND	mg/kg	NC		20
Silver, Total	ND	ND	mg/kg	NC		20
Sodium, Total	169.J	150J	mg/kg	NC		20



Lab Duplicate Analysis Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Parameter	Native Sample	Duplica	Duplicate Sample		RPD	R	PD Limits
Total Metals - Mansfield Lab Associated sample(s): 01-0	08 QC Batch ID:	WG1105285-4	QC Sample:	L1810464-05	Client ID:	DUP Sample	!
Thallium, Total	ND		ND	mg/kg	NC		20
Vanadium, Total	16.3		17.1	mg/kg	5		20
Zinc, Total	136.		182	mg/kg	29	Q	20



INORGANICS & MISCELLANEOUS



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Report Date: **Project Number:** 04/16/18 Not Specified

SAMPLE RESULTS

Lab ID: Date Collected: L1812163-01 04/06/18 08:17

Client ID: SS-1 Date Received: 04/09/18

Not Specified Sample Location: ORCHARD PARK, NY Field Prep:

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry -	Westborough Lab									
Solids, Total	96.4		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-02 Date Collected: 04/06/18 08:20

Client ID: SS-2 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry -	Westborough Lab)								
Solids, Total	95.8		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-03 Date Collected: 04/06/18 08:22

Client ID: SS-3 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry - V	Westborough Lab									
Solids, Total	93.2		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-04 Date Collected: 04/06/18 08:25

Client ID: SS-4 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry -	Westborough Lab)								
Solids, Total	95.5		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-05 Date Collected: 04/06/18 08:29

Client ID: SS-5 Date Received: 04/09/18

Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry - \	Westborough Lab)								
Solids, Total	95.6		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-06 Date Collected: 04/06/18 09:30

Client ID: B-1 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry - \	Westborough Lab)								
Solids, Total	93.0		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-07 Date Collected: 04/06/18 09:50

Client ID: B-2 Date Received: 04/09/18
Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Sample Depth:

Matrix: Soil

Parameter	Result	Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry -	Westborough Lab									
Solids, Total	87.0		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Project Name: SOUTHWESTERN PHASE 2 Lab Number: L1812163

Project Number: Not Specified Report Date: 04/16/18

SAMPLE RESULTS

Lab ID: L1812163-08 Date Collected: 04/06/18 10:10

Client ID: B-3 Date Received: 04/09/18 Sample Location: ORCHARD PARK, NY Field Prep: Not Specified

Campio 200anom Cittoria in Cit

Sample Depth:

Matrix: Soil

Parameter	Result (Qualifier	Units	RL	MDL	Dilution Factor	Date Prepared	Date Analyzed	Analytical Method	Analyst
General Chemistry - \	Westborough Lab									
Solids, Total	90.7		%	0.100	NA	1	-	04/10/18 10:13	121,2540G	RI



Lab Duplicate Analysis
Batch Quality Control

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Lab Number:

L1812163

Report Date:

04/16/18

Parameter	Native Sam	ple D	uplicate Sample	Units	RPD	Qual	RPD Limits
General Chemistry - Westborough Lab	Associated sample(s): 01-08	QC Batch ID:	WG1105098-1	QC Sample:	L1812163-01	Client ID:	SS-1
Solids, Total	96.4		97.0	%	1		20



Project Name: **SOUTHWESTERN PHASE 2**

Project Number: Not Specified

Lab Number: L1812163 **Report Date:** 04/16/18

Sample Receipt and Container Information

Were project specific reporting limits specified?

YES

Cooler Information

Custody Seal Cooler

Α Absent

Container Info	Container Information		Initial	Final	Temp			Frozen	
Container ID	Container Type	Cooler	рН	pН	deg C	Pres	Seal	Date/Time	Analysis(*)
L1812163-01A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG- TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL- TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE- TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE- TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA- TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-01B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-01C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-01X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-01Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-01Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-02A	Metals Only-Glass 60mL/2oz unpreserved	А	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-02B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-02C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-02X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-02Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-02Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-03A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-03B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-03C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-03X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)



Lab Number: L1812163

Report Date: 04/16/18

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Container Info	Container Information		Initial	Final	Temp			Frozen	
Container ID	Container Type	Cooler		рН	deg C	Pres	Seal	Date/Time	Analysis(*)
L1812163-03Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-03Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-04A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-04B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-04C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-04X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-04Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-04Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-05A	Metals Only-Glass 60mL/2oz unpreserved	A	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-05B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-05C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-05X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-05Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-05Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-06A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-06B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-06C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-06X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-06Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-06Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)



Lab Number: L1812163

Report Date: 04/16/18

Project Name: SOUTHWESTERN PHASE 2

Project Number: Not Specified

Container Information			Initial	Final	Temp			Frozen	
Container ID	Container Type	Cooler	рH	рН	deg C	Pres	Seal	Date/Time	Analysis(*)
L1812163-07A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-07B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-07C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-07X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-07Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-07Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-08A	Metals Only-Glass 60mL/2oz unpreserved	Α	NA		3.1	Y	Absent		BE-TI(180),AS-TI(180),BA-TI(180),AG-TI(180),AL-TI(180),CR-TI(180),NI-TI(180),TL-TI(180),CU-TI(180),PB-TI(180),SB-TI(180),SE-TI(180),ZN-TI(180),CO-TI(180),V-TI(180),FE-TI(180),HG-T(28),MG-TI(180),MN-TI(180),CA-TI(180),CD-TI(180),K-TI(180),NA-TI(180)
L1812163-08B	Vial Large Septa unpreserved (4oz)	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-08C	Glass 250ml/8oz unpreserved	Α	NA		3.1	Υ	Absent		NYTCL-8270(14),TS(7),NYTCL-8082(14)
L1812163-08X	Vial MeOH preserved split	Α	NA		3.1	Υ	Absent		NYTCL-8260-R2(14)
L1812163-08Y	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)
L1812163-08Z	Vial Water preserved split	Α	NA		3.1	Υ	Absent	10-APR-18 13:21	NYTCL-8260-R2(14)



Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

GLOSSARY

Acronyms

EDL - Estimated Detection Limit: This value represents the level to which target analyte concentrations are reported as estimated

values, when those target analyte concentrations are quantified below the reporting limit (RL). The EDL includes any adjustments from dilutions, concentrations or moisture content, where applicable. The use of EDLs is specific to the analysis

of PAHs using Solid-Phase Microextraction (SPME).

EPA - Environmental Protection Agency.

LCS - Laboratory Control Sample: A sample matrix, free from the analytes of interest, spiked with verified known amounts of

analytes or a material containing known and verified amounts of analytes.

LCSD - Laboratory Control Sample Duplicate: Refer to LCS.

LFB - Laboratory Fortified Blank: A sample matrix, free from the analytes of interest, spiked with verified known amounts of

analytes or a material containing known and verified amounts of analytes.

MDL - Method Detection Limit: This value represents the level to which target analyte concentrations are reported as estimated values, when those target analyte concentrations are quantified below the reporting limit (RL). The MDL includes any

adjustments from dilutions, concentrations or moisture content, where applicable.

MS - Matrix Spike Sample: A sample prepared by adding a known mass of target analyte to a specified amount of matrix sample for

which an independent estimate of target analyte concentration is available.

MSD - Matrix Spike Sample Duplicate: Refer to MS.

NA - Not Applicable.

NC - Not Calculated: Term is utilized when one or more of the results utilized in the calculation are non-detect at the parameter's

reporting unit.

NDPA/DPA - N-Nitrosodiphenylamine/Diphenylamine.

NI - Not Ignitable.

NP - Non-Plastic: Term is utilized for the analysis of Atterberg Limits in soil.

RL - Reporting Limit: The value at which an instrument can accurately measure an analyte at a specific concentration. The RL

includes any adjustments from dilutions, concentrations or moisture content, where applicable.

RPD - Relative Percent Difference: The results from matrix and/or matrix spike duplicates are primarily designed to assess the precision of analytical results in a given matrix and are expressed as relative percent difference (RPD). Values which are less

precision of analytical results in a given matrix and are expressed as relative percent difference (RPD). Values which are less than five times the reporting limit for any individual parameter are evaluated by utilizing the absolute difference between the

values; although the RPD value will be provided in the report.

SRM - Standard Reference Material: A reference sample of a known or certified value that is of the same or similar matrix as the

associated field samples.

STLP - Semi-dynamic Tank Leaching Procedure per EPA Method 1315.

TIC - Tentatively Identified Compound: A compound that has been identified to be present and is not part of the target compound

list (TCL) for the method and/or program. All TICs are qualitatively identified and reported as estimated concentrations.

Footnotes

- The reference for this analyte should be considered modified since this analyte is absent from the target analyte list of the original method.

Terms

Analytical Method: Both the document from which the method originates and the analytical reference method. (Example: EPA 8260B is shown as 1,8260B.) The codes for the reference method documents are provided in the References section of the Addendum.

Final pH: As it pertains to Sample Receipt & Container Information section of the report, Final pH reflects pH of container determined after adjustment at the laboratory, if applicable. If no adjustment required, value reflects Initial pH.

Frozen Date/Time: With respect to Volatile Organics in soil, Frozen Date/Time reflects the date/time at which associated Reagent Water-preserved vials were initially frozen. Note: If frozen date/time is beyond 48 hours from sample collection, value will be reflected in 'bold'.

Initial pH: As it pertains to Sample Receipt & Container Information section of the report, Initial pH reflects pH of container determined upon receipt, if applicable.

Total: With respect to Organic analyses, a 'Total' result is defined as the summation of results for individual isomers or Aroclors. If a 'Total' result is requested, the results of its individual components will also be reported. This is applicable to 'Total' results for methods 8260, 8081 and 8082.

Data Qualifiers

A - Spectra identified as "Aldol Condensation Product".

B - The analyte was detected above the reporting limit in the associated method blank. Flag only applies to associated field samples that have detectable concentrations of the analyte at less than ten times (10x) the concentration found in the blank. For MCP-related

Report Format: DU Report with 'J' Qualifiers



Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

Data Qualifiers

projects, flag only applies to associated field samples that have detectable concentrations of the analyte at less than ten times (10x) the concentration found in the blank. For DOD-related projects, flag only applies to associated field samples that have detectable concentrations of the analyte at less than ten times (10x) the concentration found in the blank AND the analyte was detected above one-half the reporting limit (or above the reporting limit for common lab contaminants) in the associated method blank. For NJ-Air-related projects, flag only applies to associated field samples that have detectable concentrations of the analyte above the reporting limit. For NJ-related projects (excluding Air), flag only applies to associated field samples that have detectable concentrations of the analyte, which was detected above the reporting limit in the associated method blank or above five times the reporting limit for common lab contaminants (Phthalates, Acetone, Methylene Chloride, 2-Butanone).

- Co-elution: The target analyte co-elutes with a known lab standard (i.e. surrogate, internal standards, etc.) for co-extracted analyses.
- Concentration of analyte was quantified from diluted analysis. Flag only applies to field samples that have detectable concentrations
 of the analyte.
- E Concentration of analyte exceeds the range of the calibration curve and/or linear range of the instrument.
- G The concentration may be biased high due to matrix interferences (i.e, co-elution) with non-target compound(s). The result should be considered estimated.
- H The analysis of pH was performed beyond the regulatory-required holding time of 15 minutes from the time of sample collection.
- I The lower value for the two columns has been reported due to obvious interference.
- M Reporting Limit (RL) exceeds the MCP CAM Reporting Limit for this analyte.
- NJ Presumptive evidence of compound. This represents an estimated concentration for Tentatively Identified Compounds (TICs), where the identification is based on a mass spectral library search.
- P The RPD between the results for the two columns exceeds the method-specified criteria.
- Q The quality control sample exceeds the associated acceptance criteria. For DOD-related projects, LCS and/or Continuing Calibration Standard exceedences are also qualified on all associated sample results. Note: This flag is not applicable for matrix spike recoveries when the sample concentration is greater than 4x the spike added or for batch duplicate RPD when the sample concentrations are less than 5x the RL. (Metals only.)
- **R** Analytical results are from sample re-analysis.
- **RE** Analytical results are from sample re-extraction.
- S Analytical results are from modified screening analysis.
- Estimated value. The Target analyte concentration is below the quantitation limit (RL), but above the Method Detection Limit (MDL) or Estimated Detection Limit (EDL) for SPME-related analyses. This represents an estimated concentration for Tentatively Identified Compounds (TICs).
- ND Not detected at the method detection limit (MDL) for the sample, or estimated detection limit (EDL) for SPME-related analyses.

Report Format: DU Report with 'J' Qualifiers



Project Name:SOUTHWESTERN PHASE 2Lab Number:L1812163Project Number:Not SpecifiedReport Date:04/16/18

REFERENCES

Test Methods for Evaluating Solid Waste: Physical/Chemical Methods. EPA SW-846. Third Edition. Updates I - IV, 2007.

121 Standard Methods for the Examination of Water and Wastewater. APHA-AWWA-WEF. Standard Methods Online.

LIMITATION OF LIABILITIES

Alpha Analytical performs services with reasonable care and diligence normal to the analytical testing laboratory industry. In the event of an error, the sole and exclusive responsibility of Alpha Analytical shall be to re-perform the work at it's own expense. In no event shall Alpha Analytical be held liable for any incidental, consequential or special damages, including but not limited to, damages in any way connected with the use of, interpretation of, information or analysis provided by Alpha Analytical.

We strongly urge our clients to comply with EPA protocol regarding sample volume, preservation, cooling, containers, sampling procedures, holding time and splitting of samples in the field.



Alpha Analytical, Inc. Facility: Company-wide

Department: Quality Assurance

Title: Certificate/Approval Program Summary

ID No.:17873

Revision 11

Page 1 of 1

Published Date: 1/8/2018 4:15:49 PM

Certification Information

The following analytes are not included in our Primary NELAP Scope of Accreditation:

Westborough Facility

EPA 624: m/p-xylene, o-xylene

EPA 8260C: NPW: 1,2,4,5-Tetramethylbenzene; 4-Ethyltoluene, Azobenzene; SCM: lodomethane (methyl iodide), Methyl methacrylate, 1,2,4,5-Tetramethylbenzene; 4-Ethyltoluene.

EPA 8270D: NPW: Dimethylnaphthalene,1,4-Diphenylhydrazine; SCM: Dimethylnaphthalene,1,4-Diphenylhydrazine.

EPA 300: DW: Bromide EPA 6860: SCM: Perchlorate

EPA 9010: NPW and SCM: Amenable Cyanide Distillation

SM4500: NPW: Amenable Cyanide, Dissolved Oxygen; SCM: Total Phosphorus, TKN, NO2, NO3.

Mansfield Facility

SM 2540D: TSS

EPA 8082A: NPW: PCB: 1, 5, 31, 87,101, 110, 141, 151, 153, 180, 183, 187.

EPA TO-15: Halothane, 2,4,4-Trimethyl-2-pentene, 2,4,4-Trimethyl-1-pentene, Thiophene, 2-Methylthiophene,

3-Methylthiophene, 2-Ethylthiophene, 1,2,3-Trimethylbenzene, Indan, Indene, 1,2,4,5-Tetramethylbenzene, Benzothiophene, 1-Methylnaphthalene.

Biological Tissue Matrix: EPA 3050B

The following analytes are included in our Massachusetts DEP Scope of Accreditation

Westborough Facility:

Drinking Water

EPA 300.0: Chloride, Nitrate-N, Fluoride, Sulfate; EPA 353.2: Nitrate-N, Nitrite-N; SM4500NO3-F: Nitrate-N, Nitrite-N; SM4500F-C, SM4500CN-CE, EPA 180.1, SM2130B, SM4500CI-D, SM2320B, SM2540C, SM4500H-B

EPA 332: Perchlorate; EPA 524.2: THMs and VOCs; EPA 504.1: EDB, DBCP.

Microbiology: SM9215B; SM9223-P/A, SM9223B-Colilert-QT,SM9222D.

Non-Potable Water

SM4500H,B, EPA 120.1, SM2510B, SM2540C, SM2320B, SM4500CL-E, SM4500F-BC, SM4500NH3-BH: Ammonia-N and Kjeldahl-N, EPA 350.1: Ammonia-N, LACHAT 10-107-06-1-B: Ammonia-N, EPA 351.1, SM4500NO3-F, EPA 353.2: Nitrate-N, EPA 351.1, SM4500P-B, E, EA 351.1, SM4500P-B, EA 351.1, SM500P-B, EA 351.1, SM500P-B, EA 351.1, SM500P-B, EA SM4500SO4-E, SM5220D, EPA 410.4, SM5210B, SM5310C, SM4500CL-D, EPA 1664, EPA 420.1, SM4500-CN-CE, SM2540D.

EPA 624: Volatile Halocarbons & Aromatics, EPA 608: Chlordane, Toxaphene, Aldrin, alpha-BHC, beta-BHC, gamma-BHC, delta-BHC, Dieldrin, DDD, DDE, DDT, Endosulfan I, Endosulfan II, Endosulfan sulfate, Endrin, Endrin Aldehyde, Heptachlor, Heptachlor Epoxide, PCBs

EPA 625: SVOC (Acid/Base/Neutral Extractables), EPA 600/4-81-045: PCB-Oil.

Microbiology: SM9223B-Colilert-QT; Enterolert-QT, SM9221E, SM9222D.

Mansfield Facility:

Drinking Water

EPA 200.7: Al, Ba, Be, Cd, Cr, Cu, Mn, Ni, Na, Ag, Ca, Zn. EPA 200.8: Al, Sb, As, Ba, Be, Cd, Cr, Cu, Pb, Mn, Ni, Se, Ag, TL, Zn. EPA 245.1 Hg. EPA 522.

Non-Potable Water

EPA 200.7: Al, Sb, As, Be, Cd, Ca, Cr, Co, Cu, Fe, Pb, Mg, Mn, Mo, Ni, K, Se, Ag, Na, Sr, TL, Ti, V, Zn.

EPA 200.8: Al, Sb, As, Be, Cd, Cr, Cu, Pb, Mn, Ni, Se, Ag, TL, Zn.

EPA 245.1 Hg.

SM2340B

For a complete listing of analytes and methods, please contact your Alpha Project Manager.

Pre-Qualtrax Document ID: 08-113 Document Type: Form

	NEWYORK	Service Centers			Page	4			SERVICE	4500		HOLDE,		
ALPHA	NEW YORK CHAIN OF	Mahwah, NJ 07430: 35 White Albany, NY 12205: 14 Walker				1	Date Rec'd In Lab					ALPHA Job#		
	CUSTODY	Tonawanda, NY 14150: 275 C		05				111	Lau	111	3/16		L1812163	
Westborough, MA 0158		Project Information	No.			118	Delive	erable	5				Billing Information	
8 Walkup Dr. TEL: 508-898-9220	320 Forbes Blvd TEL: 508-822-9300	Project Name:	Southwester	n Phase 2				ASP.	-A		ASP-B		Same as Client Info	
FAX: 508-898-9193	FAX 508-822-3268		chard Park	L, NY				EQu	S (1 File)		EQuIS (4 File)	PO# 0 / 01	F
Client Information		Project #					☐ Other						Octovus Storage	- 1
Control of the contro	ikes Env & Safety	(Use Project name as I	Project #)				Regu	latory	Requirem	ent			Disposal Site Information	
ddress: 50 Ridge		Project Manager: Ma						NYT	ogs		NY Part 3	375	Please identify below location of	1
datess. 50 mag	711000	ALPHAQuote #:						AWQ	Standards		NY CP-5	1	applicable disposal facilities.	
hone: 716-827	-0700	Turn-Around Time	A THE REAL PROPERTY.	30/200	B. St.	N. E-20		NYR	estricted Us	е 🗌	Other		Disposal Facility:	-
11000	-0700	Standa	rd 🔀	Due Date:			X	NY U	nrestricted U	lse			☐ NJ ☐ NY	
ax: mail: www.gre	atlakesenvironmental.co		Charles of the Control of the Contro	# of Days:				NYC	Sewer Disch	narge			Other: NA	
	e been previously analyz						ANAL	YSIS					Sample Filtration	
	ific requirements/com												Done	
Atter project spec								SCB	TALS				Lab to do Preservation Lab to do	
Please specify Met	als or TAL.						8260	8270/PCB	TAL METALS				(Please Specify below)	
ALPHA Lab ID	s	ample ID		ection	Sample Matrix	Sampler's Initials	1		=				Sample Specific Comments	- 10
(Lab Use Only)		versitassautes .	Date	Time	-				-	+	+-+	+	Gample Specific Comments	-
12/63-01	55-1		4/6/18	0817	Soil	ММ	X	X	X	+-	-	_		$^{+}$
02	55-2		4/6/18	0830	Soil	MM	X	X	X	+	+ +	+		$^{+}$
05	55-3		4/6/18	0833	Soil	MM	X	X	X	+	-	_		_
04	55-4		4/4/18	0835	Soil	MM	Х	X	X	+	-	_		+
05	55-5		4/6/18	0829	Soil	MM	X	X	X	_	+	-		+
16	B-1		4-6-18	0930	Soil	ww	Х	X	X	+	+-+	_		+
07	B-2		4-6-18	0950	Soil	mm	Х	Х	X	+	-		-	+
86	8-3		4-6-18	1010	Soil	MM	X	X	X	-	-	_	-	+
					Soil		X	X	X	_	+	_		5
					Soil-3		X -	X	X	_	-	_		1
Preservative Code: \(= \text{None} \) \(= \text{HCI} \)	Container Code P = Plastic A = Amber Glass	Westboro: Certification Mansfield: Certification			Cor	ntainer Type	A	A	A				Please print clearly, legit and completely. Sample:	
= HNO ₃ = H ₂ SO ₄ = NaOH	V = Vial G = Glass B = Bacteria Cup					Preservative	A	A	A				not be logged in and turnaround time clock wi start until any ambiguitie	
= NaOH = MeOH S = NaHSO ₄	MeOH C = Cube Relinquished By:			4-9-18	/Time //:10 //:40	Jan.	Recai		y: AAC	4/	Date/T	-	resolved, BY EXECUTIN	IG T

APPENDIX B QUALITY ASSURANCE PROJECT PLAN

QUALITY ASSURANCE PROJECT PLAN

BROWNFIELD CLEANUP PROGRAM For **OCTAVUS STORAGE NY, LLC.** 3678 Southwestern Boulevard, Orchard Park, NY 14127 BCP # C915342



Prepared For: **OCTAVUS STORAGE NY, LLC** 1810 West Kennedy Boulevard, Tampa, FL 33606 HEI Project No: e1728

Prepared By:

Hazard Evaluations, Inc. 3636 North Buffalo Road Orchard Park, New York 14127 Orchard Park, New York 14127 (716) 667-3130

Wittman GeoSciences, PLLC 3636 North Buffalo Road 716-574-1513

March 12, 2019



TABLE OF CONTENTS

		<u>Page</u>
1.0	INTRODUCTION	1
1.1	Project Scope	1
1.2		
1.3	PROJECT SUB-CONTRACTORS	2
2.0	FIELD INVESTIGATION PROCEDURES	2
2.1	AIR MONITORING	2
2.2	SOIL SCREENING AND LOGGING	
2.3	SOIL SAMPLE COLLECTION	3
2.4	Soil Borings	4
2.5	TEST PITS	
2.6	MONITORING WELL INSTALLATION	
2.7		
	.7.1 Monitoring Well Development	
	.7.2 Groundwater Sampling	
2.8		
2.9	= *** ::=:::= = = * * : :: ::::::::::::	
2.10		
2.1′		
3.0	SAMPLE HANDLING AND MANAGEMENT	9
3.1	SAMPLE LABEL AND IDENTIFICATION	9
3.2	CHAIN OF CUSTODY	10
4.0	FIELD DOCUMENTATION	10
5.0	ANALYTICAL LABORATORY QA/QC PROTOCOLS	11
5.1	ANALYTICAL METHODS, PROCEDURES AND CALIBRATION	11
5	.1.2 Analytical Methods	
5	.1.3 Laboratory Instrumentation & Equipment	11
5.	.1.3 Field Equipment	
5.2		
5.3	CORRECTIVE ACTIONS	13
6.0	DATA USABII ITY	13



TABLE OF CONTENTS (CONT'D)

LIST OF TABLES

Table 1 Analytical Testing Program Summary

Table 2 Sample Container, Volume, Preservation and Holding Time Requirements

LIST OF FIGURES

Figure 1 Proposed Sampling Locations

Figure 2 Proposed Overburden Monitoring Well Installation Diagram

LIST OF ATTACHMENTS

Attachment 1 Resumes
Attachment 2 Field Forms

Attachment 3 Emergent Contaminant Sampling and Laboratory Analysis



1.0 INTRODUCTION

This Quality Assurance Project Plan (QAPP) has been cooperatively developed by Hazard Evaluations Inc. (HEI) and Wittman GeoSciences, PLLC, (WGS) and was prepared for Octavus Storage NY, LLC (Octavus) located at 3678 Southwestern Boulevard in Orchard Park, New York. The QAPP was prepared in general accordance with the requirements of Section 2.4 of the New York State Department of Environmental Conservation (NYSDEC) DER-10, Technical Guidance for Site Investigation and Remediation, dated May 2010 (DER-10).

The QAPP is designed to produce data of the quality necessary to achieve the project objectives. The objective of the quality assurance/quality control (QA/QC) protocol and procedures is to ensure the information, data, and decisions associated with the project are technically sound and properly documented.

1.1 **Project Scope**

This QAPP presents the project scope, objectives, organization, planned activities, data quality objectives, QA/QC procedures and sampling procedures. This project involves test borings, test pits, monitoring well installation, monitoring well development, and subsurface soil and groundwater sample collection. Proposed sampling locations are included on Figure 1 and a summary of the anticipated number of samples and analytical testing is included on Table 1. The project goal associated with the Remedial Investigation (RI) includes the following:

- Define the nature and extent of on-site contamination in both soil and groundwater;
- Identify on-site source areas of contamination, if any;
- Collect data of sufficient quantity and quality to evaluate potential threats to the public health and environment; and
- Collect data of sufficient quantity and quality to evaluate remedial alternatives.

1.2 **Project Organization**

The general responsibilities of key project personnel are listed below and resumes are included in Attachment 1.

Project Manager: Ms. Michele Wittman, WGS Principal, has responsibility for

overall program/project management and coordination with

NYSDEC and subcontractors.

Engineering: R. Christopher Wood, PE, is responsible for engineering aspects

and responsibilities.

QA Officer: Mark Hanna, CHMM, HEI Principal, will serve as Quality

Assurance Officer (QAO), and will be responsible for laboratory and data validation, subcontractor procurement and assignment,



as well as data usability reports. The QAO may conduct audits of the operations at the site to ensure that work is being performed in accordance with the QAAP.

Field Team: Eric Betzold will have overall responsibility for on-site

implementation of the Site Investigation project activities. The technical team will consist of experienced professionals (i.e.; engineers, geologists, scientists) to gather and analyze data, prepare project documentation and collection of various soil and

groundwater samples.

1.3 Project Sub-Contractors

Subcontractor specialists will be contracted for services relating to drilling and monitoring well installation, laboratory/analytical services, data validation services, field surveying, and waste transportation and disposal. The subcontractors will be approved by NYSDEC prior to beginning site work:

Laboratory Analysis: Alpha Analytical - A laboratory certified under the New

York State Department of Health (NYSDOH)

Environmental Laboratory Approval Program (ELAP)

Data Validation: Data Validation Services

Exploration Services: To be determined. Surveying: To be determined.

2.0 FIELD INVESTIGATION PROCEDURES

Field sampling at Octavus (Site) has been designed to obtain representative samples of various environmental media to assess the impacts that the Site may have to human health and the environment. The field investigation procedures include sampling for subsurface soils, groundwater, air and vapor samples.

Proposed sampling locations are included within the RI Work Plan. Environmental sampling and other field activities will be performed in general accordance with the appropriate techniques presented in the following guidance document:

DRAFT DER-10: Technical Guidance for Site Investigations and Remediation; NYSDEC Division of Environmental Remediation, May 2010.

Field activities are described in the following sections and in the RI Work Plan.

2.1 <u>Air Monitoring</u>

Air monitoring/screening of volatile organic compounds (VOC) for health and safety concerns will be performed with a portable organic vapor meter (OVM) equipped with a photoionization detector (PID) that is using a 10.6 electron volt (eV)



bulb. Monitoring will be done during invasive activities, which include soil borings, monitoring well installation, well development, and sampling. Detections of volatile compounds above background levels during air monitoring will require that the work be stopped until air monitoring levels decrease to acceptable background levels or until health and safety protocol are upgraded and approved by NYSDEC. On-site personnel will be outfitted in modified Level D personal protective equipment (i.e., hardhat, safety glasses, work boots and gloves).

2.2 Soil Screening and Logging

Subsurface soil samples will be collected from direct push macro-core samplers in general accordance with American Society for Testing and Material (ASTM) D6282-98: Standard Guide for Direct Push Soil Samples for Environmental Site Characteristics. Subsurface soil sampling from split-spoon samples advanced ahead of hollow steam augers will be completed in general accordance with ASTM D1586-99. A soil boring log will be prepared for each location to include date, boring location, drill rig type, blow counts, sample identification, sample depth interval, percent recovery, OVM reading, stratigraphic boundaries, and well installation information.

Subsurface soil will be sampled by opening the split-spoon sampler (borings) or slicing the core vertically down the middle with a sharp blade. Soil samples will be visually examined for evidence of suspect contamination (e.g., staining, odor) and field screened with a calibrated OVM. Prior to screening, the soil samples will be allowed to equilibrate to ambient temperature. Different portions of the soil samples will be placed within sealable plastic bags, the OVM sampling port will be placed within a corner of the bag, and the peak reading will be recorded on the boring log. Portions of the soil samples may be placed in containers for future analytical testing.

2.3 Soil Sample Collection

Soil samples selected for VOC analysis will be collected using an Encore or Terracore sampling kit. Headspace will be limited by compacting the soil into the container. Samples for VOC analysis will be placed into the appropriate container immediately after opening of the sampler, prior to making any field measurements or sample homogenization.

The soil samples will be homogenized using a "coning and quartering" procedure. The soil will be removed from the sampling equipment, transferred to a clean surface (e.g., metal foil, steel pan, bowl, etc.), and thoroughly mixed to provide a more homogeneous sample to the lab. An aliquot of the sample will then be transferred to the required sample containers and sealed with the appropriate cap.

Due to emergent contaminant sampling requirements, at this time acceptable materials for sampling include stainless steel, high density polyethylene (HDPE), PVC, silicone, acetate and polypropylene. Additional materials may be acceptable if pre-approved by NYSDEC. All sampling equipment components and sample containers should not come in contact with aluminum foil, low density polyethylene



(LDPE), glass or polytetrafluoroethylene (PTFE, Teflon™) materials including sample bottle cap liners with a PTFE layer.

Clothing that contains PTFE material (including GORE-TEX®) or that have been waterproofed with PFC materials must be avoided. Many food and drink packaging materials and "plumbers thread seal tape" contain PFCs.

All clothing worn by sampling personnel must have been laundered multiple times. The sampler must wear nitrile gloves while filling and sealing the sample bottles.

2.4 Soil Borings

Soil borings will be completed to define subsurface conditions using either direct push subsurface investigation techniques or rotary drilling with continuous split-spoon sampling and hollow-stem augers. Drilling cuttings will be visually inspected, screened for the presence of VOCs using an OVM, and managed consistent with DER-10 requirements. Soil samples for laboratory analysis will be selected in the field based on visual/olfactory observations and OVM screening results.

The drill rig, soil probe rig, tools, augers, etc., will be decontaminated between holes at an on-site temporary decontamination pad or area. Decontamination will be accomplished using steam cleaning or high pressure wash equipment. Direct push sampling equipment and split-spoon sampling devices will be cleaned manually with non-phosphate detergent (i.e., Alconox) wash and potable water, followed by a potable water rinse or a second steam cleaning, followed by a distilled/deionized water rinse. All equipment will be cleaned prior to leaving the Site.

2.5 Test Pits

Test pits will be completed using a track-mounted excavator and bucket to provide a detailed visual examination of near surface soil, and fill materials present on-site. Samples will be collected from the walls and/or floor through the use of a bucket, steel trowels and bowls. The samples will be placed directly into the appropriate containers, inspected and screened for the presence of VOCs using an OVM, and managed consistent with DER-10 requirements. Soil samples for laboratory analysis will be selected in the field based on visual and olfactory observations and OVM screening results.

Sampling equipment will be cleaned manually with non-phosphate detergent (i.e., Alconox) wash and potable water, followed by a potable water rinse or a second steam cleaning, followed by a distilled/deionized water rinse. All equipment will be cleaned prior to leaving the Site.

2.6 Monitoring Well Installation

Monitoring wells will be constructed of two-inch ID flush coupled with Schedule 40, polyvinyl chloride (PVC) riser and screen. The actual installation depth and screen depth will be selected based on groundwater depth, observation of



subsurface materials, and headspace screening test results. In general, the screen will consist of a maximum ten foot length of 0.010-inch machine slotted well screen. A schematic of the well construction detail is provided as Figure 2.

Following the placement of the assembled screen and riser, the borehole will be backfilled. The well screen depth will be backfilled with silica sand filter pack (estimated at size #0) from the base to a minimum of one (1) foot above the well screen. A minimum one-foot layer of bentonite pellets will be placed above the sand filter and allowed to hydrate. Then, a mixture of cement/bentonite water will be placed above the bentonite seal, and the monitoring well will be completed by placing a locking steel casing or road box over the riser. Concrete will be then placed in the borehole around the protective casing and sloped away from the casing.

2.7 <u>Monitoring Well Development and Sampling</u>

2.7.1 Monitoring Well Development

Monitoring wells will be developed by utilizing either a dedicated tubing or new dedicated disposable bailer, depending on the field conditions. Fluids will not be added during the development process. New dedicated well development equipment will be utilized prior to development of each well. The well development procedure is listed below.

- Unlock the well cover and survey the ambient air and air directly at the top of the well using the OVM.
- Take a pre-development static water level measurement.
- Sound the bottom of the well and agitate/loosen accumulated sediment.
- Calculate the water volume in the well.
- Obtain initial field water quality measurements, including pH, specific conductance, turbidity, and temperature using a Horiba U-22 water quality meter (or equivalent).
- Alternate water agitation methods, such as moving a bailer or pump tubing up and down inside screened interval coupled with water removal methods (pumping or bailing) in order to suspend and remove solids/sediment from the wells.
- Record water quality meter measurements every one to three gallons of water removed and record the quantity of water removed.
- Development can cease when at least five well volumes have been removed, or the following water quality criteria are met:
 - Water is clear and free of sediment and turbidity is less than 50 nephelometric turbity units (NTUs);
 - pH is +/- 0.1 standard unit between readings;
 - Specific conductivities is +/-3% between readings; and
 - Temperature is +/-10% between readings.
- Record post-development water level readings and development information on well development logs.



After the water level has returned to its pre-purge level (or within a maximum of two hours, if the well has recharged sufficiently to allow sampling), samples will be collected from the middle of the screened portion of the well for overburden wells. If the water level is slow to recharge and does not reach its pre-purge level within two hours, then samples can be collected after sufficient water has recharged. The degree of recharge will be indicated in field notes, with time and depth-to-water noted.

2.7.2 Groundwater Sampling

Groundwater samples will be collected by utilizing low-flow sampling techniques with dedicated tubing or by conventional methods using a new dedicated disposable bailer. A peristaltic pump and new disposable high density polyethylene (HDPE) tubing will be used at each location. Tubing and sampling equipment will be cleaned upon arrival at the Site. The well will be sampled after removal of three well volumes or well purging.

A Well Data Sheet should be completed during groundwater sampling. Each well to be sampled will have designated pre-labeled and certified clean sample bottles. The following steps describe the groundwater sample procedure:

- Unlock and remove the well cap. Test the air at the wellhead with the OVM.
- Measure the static water level and determine the total well volume.
- Slowly lower the dedicated bailer or tubing into the well. Purge the well
 a minimum of three well volumes. If the well goes dry during bailing,
 allow for full recovery and sample. If recovery takes longer than 20
 minutes, proceed to next well but return to sample within 24 hours.
- Fill the appropriate sample bottles. Two or three (depending on laboratory-specific requirements) 40-ml glass vials (with Teflon septa) will be used to collect samples for VOCs. Samples will be collected in the following order: volatile organic compounds, semi-volatile organic compounds, PCBs/pesticides/herbicides, and metals. If the well should go dry during sampling, the well will be re-sampled the next day. The second attempt to sample the well will proceed with the same sample order.
- Sampling preservatives will be added to the provided jars by the laboratory. The following parameters require additional special handling:
 - VOC samples must be free of air bubbles. When the container is determined to be bubble free, the sample containers should be immediately chilled.
 - Metals analysis should be preserved with nitric acid to a pH less than two.
- Record pertinent information in the field logbook and well data sheet.
- Lock the well, inspect the well site, and note any maintenance required.
- Purged water will be containerized for future disposal.



Per- and polyfluoroalkyl substances (PFAS) Sampling Protocol

Groundwater sample collection procedure for PFAS will be done in accordance with NYSDEC protocol, as included in the Quality Assurance/Quality Control (QAPP) found in Appendix B. At this time acceptable materials for sampling include stainless steel, high density polyethylene (HDPE), PVC, silicone, acetate and polypropylene. Equipment blanks should be generated at least daily. Additional materials may be acceptable if pre-approved by NYSDEC. Requests to use alternate equipment should include clean equipment blanks. NOTE: Grunfos pumps and bladder pumps are known to contain PFC materials (e.g. Teflon™ washers for Grunfos pumps and LDPE bladders for bladder pumps). All sampling equipment components and sample containers should not come in contact with aluminum foil, low density polyethylene (LDPE), glass or polytetrafluoroethylene (PTFE, Teflon™) materials including sample bottle cap liners with a PTFE layer. Standard two step decontamination using detergent and clean water rinse will be performed for equipment that does come in contact with PFC materials. Clothing that contains PTFE material (including GORE-TEX®) or that have been waterproofed with PFC materials must be avoided. Many food and drink packaging materials and "plumbers thread seal tape" contain PFCs.

All clothing worn by sampling personnel must have been laundered multiple times. The sampler must wear nitrile gloves while filling and sealing the sample bottles.

Pre-cleaned sample bottles with closures, coolers, ice, sample labels and a chain of custody form will be provided by the laboratory.

- 1. Fill two pre-cleaned 500 mL HDPE or polypropylene bottle with the sample.
- 2. Cap the bottles with an acceptable cap and liner closure system.
- 3. Label the sample bottles.
- 4. Fill out the chain of custody.
- 5. Place in a cooler maintained at $4 \pm 2^{\circ}$ Celsius.

2.8 Background Samples

Due to the historically undeveloped surrounding parcels and limited industrial usage on-site, soil and groundwater samples have not been pre-designed as likely to characterize Site background conditions.

2.9 **Equipment Decontamination**

In order to reduce the potential for cross-contamination of samples collected during the project, sampling equipment will be decontaminated to ensure that data is acceptable. It is anticipated that most of the materials used in sample collection will be disposable one-time use materials, such as sampling containers, bailers, tubing, gloves, etc.



Non-dedicated material such as split-spoon samples, stainless steel mixing bowls, drill rig, water-level indicator, etc., will be decontaminated by the following methods:

- Steam clean the equipment within a dedicated decontamination area; or
- Decontaminate the equipment by scrubbing/washing with a laboratory-grade detergent (e.g. Alconox) to remove visible contamination, followed by potable (tap) water and analyte-free water rinses. Tap water may be used from any treated municipal water system.

The effectiveness of the equipment decontamination of non-dedicated sampling equipment will be evaluated via analytical testing of rinsate blanks. Decontamination liquids, disposable equipment, and personal protective equipment (PPE) will be containerized for future disposal.

2.10 Storage and Disposal of Investigation-Derived Waste

The sampling methods and equipment have been selected to limit the need for decontamination and the volume of waste material to be generated. Investigation-derived material (e.g., drill cuttings and purge water) generated will be presumed to be non-hazardous waste and will be disposed at the boring or well from which the material was derived. Excess auger cuttings will be drummed and stored on Site for future disposal. Monitoring well development/purge water will be containerized in 55-gallon drums for testing and future off-site disposal.

Personal protective equipment and disposable sampling equipment will be placed in plastic garbage bags for disposal as a non-hazardous waste.

Decontamination water used in steam cleaning and/or spoon washing and rinse-water, including detergent, may be generated during Site work. Tap and analyte-free water used for rinsing will be allowed to percolate back into the ground, or will be disposed into a sanitary system. Non-phosphate detergent and rinse-water will be disposed into a sanitary system.

2.11 Survey/Site Mapping

A topographic base map will be prepared by a New York State licensed surveyor. The base map will include property lines, buildings, fence lines, and other key Site features. This will allow measurement of the actual exploration locations and elevations. The surveyor will establish the horizontal location and vertical elevations. The map will include the RI investigation/sampling locations. The monitoring well vertical measurements will include the ground surface at exploration locations, plus the top of casing and top of riser at monitoring well locations. The top of riser will serve as the water level monitoring point. Soil/fill boring locations will be field located and incorporated within the survey.



3.0 SAMPLE HANDLING AND MANAGEMENT

Various environmental samples will be collected during the RI investigation work. The procedures below will assist in documentation and tracing of the various samples. During sampling, field personnel will wear disposable latex or nitrile gloves. Gloves will be changed and discarded between sampling locations.

Laboratory analysis samples will be placed in new laboratory-grade containers. Appropriate sample preservatives will be added to the sample containers by the laboratory prior to delivery to the project site. The specific volume and preservation of samples, if any, is summarized on Table 2. Samples will be kept in coolers, on ice, and shipped to the laboratory within 48-hours from sample collection.

3.1 Sample Label and Identification

Each field and QC sample will be identified by a self-adhesive, non-removable label placed on the sample containers. The label information will include, at a minimum, client name, Site location, data and time of collection, sample identification number, sampler's name, and notes, as needed, recorded in waterproof ink. All sample bottles within each shipping container will be individually labeled with the laboratory-provided label.

Each sample will be labeled with a unique identification using the following test location designations:

Designation	Media Type	Sample Location	Example
SB	Soil	Soil boring number with sample depth interval (x-x')	SB1 (8-10')
TP	Soil	Test pit number with sample depth interval (x-x')	TP3 (4-8')
MW	Groundwater	Monitoring well with well number	MW2
EX	Soil	Excavation confirmation sample with sample depth interval	EX3 (1-2')
TB	Trip blank	None – include day/month/year	TB1 – 10/25/16
RB	Rinsate blank	Any – rinsate of sampling equipment; include day/month/year	RB2 – 10/25/16
MS/MSD	Matrix spike/ matrix spike duplicate	Any – identify original sample location	SB1 MS MW2 MSD

Quality control (QC) field duplicate samples will be submitted blind to the laboratory. The sample identifications (of the original sample and its field duplicate) will be marked in the project-specific field book and on the copy of the chain-of-custody kept by the sampler and copied to the project manager.



3.2 Chain of Custody

A chain-of-custody form will trace the path of sample containers from the project Site to the laboratory. An example chain-of-custody is included in Attachment 2. The chain-of-custody documentation will accompany the samples from their inception until analysis. Pertinent field information will be included on the chain-of-custody, including client name, project name/location, sampler name, sample identification number, date, time, media, sample type (i.e., grab/composite), number of containers, analysis required, and preservation.

Samples will be packaged into coolers used for shipment. The cooler will be packed with ice (or equivalent) to maintain sample temperature at 4°C. The chain-of-custody forms will be signed and placed in a sealed plastic bag in the cooler. The cooler will be sealed and the custody seal placed over the cooler opening, designed to break if opened or disturbed. The custody seal will be signed and dated and shipping tape will be wrapped around the cooler and over the custody seal. Sample receipt personnel at the laboratory will document whether the custody seals remained intact upon arrival and lab personnel will sign the chain-of-custody form.

4.0 FIELD DOCUMENTATION

Daily field activities will be recorded in a bound field notebook. The field notebook will include the following daily information for Site activities:

- Date, time of arrival, time of departure, weather conditions;
- Field staff, sub-contractors or other personnel on Site;
- Description of field activities and location of work area;
- Equipment used on Site (such as drill rig, operator);
- Field observations and descriptions, such as soil descriptions, well/piezometer installation information, evidence of contamination, staining, odors, etc.;
- Field measurements (OVM, water quality readings) and calibration;
- Sampling locations, depths, identification numbers, time, etc.;
- Sampling location measurements;
- Chain-of-custody information; and
- Modifications to scope of work or issues encountered.

Field notes may be transferred to soil boring logs or monitoring well forms as part of the RI. Typical forms to be utilized during the field investigation are presented in Attachment 2 and include:

- Daily Field Report;
- Soil Boring Log or Test Pit Log;
- Monitoring Well Installation Log;
- Well Development Data Sheet;
- Chain of Custody;
- Building Inventory; and
- SVI Sampling Data Sheet.



5.0 ANALYTICAL LABORATORY QA/QC PROTOCOLS

This section describes the analytical methods, principles and procedures that will be used to generate quality data. These protocols include laboratory calibration, field equipment calibration, QC sample collection and analysis, quantitative evaluation of data quality protocols and data qualification, if necessary.

5.1 Analytical Methods, Procedures and Calibration

Chemical analysis for samples collected during the field work will be completed by a laboratory capable of performing project specific analysis as included in this QAAP.

5.1.2 Analytical Methods

Sample analytical analysis will be consistent with the NYSDEC Analytical Services Protocols (ASP) Category B requirements. Specific methods and references for each parameter, including sample preservation and holding times, are shown on Table 2. Quantification and detections limits for all analysis are those specified under the appropriate test methods.

NYSDEC has provided specific laboratory requirements associated with 1,4-dioxane and PFSA analytical methods, which are included in Attachment 3.

5.1.3 Laboratory Instrumentation & Equipment

Laboratory instruments and equipment will be calibrated following USEPA's SW-846 analytical methods protocol and laboratory requirements.

5.1.3 Field Equipment

Various field equipment will be used during the project. Calibration of the field equipment will be completed in accordance with manufacturer's specifications, prior to the start of each day.

Organic Vapor Meter – Real-time monitoring for VOCs will be done with an organic vapor meter (OVM) equipped with a photoionization detector (PID) to evaluate the nature and extent of potential petroleum or solvent impacts at the Site. The OVM will be calibrated on a daily basis in accordance with manufacturer's specifications.

Particulate Monitoring Equipment – Particulate air monitoring will be completed during soil excavation activities as noted in the Community Air Monitoring Plan (CAMP). Measurements will be collected along the upwind perimeter of the excavation areas to assess the amount of particulates naturally occurring in the air. The particulate meter will be regularly calibrated in accordance with the manufacturer's specifications.



Additional Field Equipment – Additional field equipment will be used as part of the project including an electric static water level indicator and Horiba U-22 water quality meter that measures pH, specific conductivity, temperature, dissolved oxygen, oxygen reduction potential, and turbidity. The meters will be calibrated in accordance with the manufacturer's specifications.

5.2 **Quality Control Samples**

Analytical methods, summarized on Table 2, to be utilized for laboratory sample analysis address the quality control to be used and the frequency of replicates, blanks and calibration standards for laboratory analytical equipment. Several types of field QC samples will be collected and submitted for laboratory analysis including trip blanks, sample duplicate, matrix spike and matrix spike duplicate.

Trip blanks – A trip blank sample monitors for potential impacts due to handling, transport, and cross-contamination from other samples during storage or laboratory contamination. The trip blanks, for aqueous VOCs only, will consist of analyte–free, reagent-grade water in VOC sampling containers to be used for the project. Trip blanks will be prepared at the laboratory, sealed, transported to the Site and returned without being opened to assess contamination that may have occurred during transport. Trip blanks will be submitted at a rate of one per cooler when aqueous VOCs are shipped to the laboratory.

Blind duplicates – Blind duplicate samples are used to monitor field and laboratory precision, as well as matrix heterogeneity. The samples are separate aliquots of the same sample, collected from the same location, at the same time, in the same manner as the first, and placed into a separate container. Each duplicate sample will be analyzed for the same parameters as the original sample collected that day. Blind duplicates will be collected at a frequency of one per 20 environmental samples of a given matrices (i.e. soil or groundwater).

Matrix spike/matrix spike duplicate (MS/MSD) – MS/MSDs are used to monitor precision and accuracy of the analytical method on various matrices. The samples are spiked with known quantities of target analytes at the laboratory. The MS/MSD will be collected at a frequency of one pair per 20 environmental samples of a given matrices (i.e. soil or groundwater).

Rinsate Blanks – The rinsate blank is used to indicate potential contamination from sample instruments used to collect and/or transfer samples. The rinsate blank will be generated by passing distilled water through and over cleaned sampling equipment. Rinsate blank samples will not be performed when dedicated disposal equipment is used. The rinsate blank will be collected at a frequency of one per 20 environmental samples of a given matrices (i.e. soil or



groundwater).

5.3 Corrective Actions

If instrument performance or data fall outside acceptable limits, then corrective actions will be taken to resolve problems and restore proper functioning of the analytical system. Actions may include recalibration or standardization of instruments, acquiring new standards, replacing equipment, repairing equipment, and reanalyzing samples or redoing sections of work. Subcontractors providing analytical services should perform their own internal laboratory audits and calibration procedures with data review conducted at a frequency so that errors and problems are detected early, thus avoiding the prospect of redoing large segments of work.

6.0 DATA USABILITY

The main objective of the Data Usability Summary Report (DUSR) is to determine whether the data presented meets the project-specific needs for data quality and data use. Data validation will be performed and a DUSR will be prepared to meet the NYSDEC requirements for analytical data generated during the RI. The DUSR will be completed in general accordance with Appendix 2B of DER-10. The findings of the DUSR will be incorporated in the RI report. Waste characterization samples will not be validated.



TABLES

TABLE 1 Analytical Testing Program Summary Octavus Storage NY, LLC

3678 Southwester Boulevard, Orchard Park, NY NYSDEC Brownfield Cleanup Program

Location	Number of Proposed Locations	Matrix	TCL VOCS	TCL SVOCs	TAL METALS Total	TAL METALS dissolved	PCBs	Pest/ Herbs	1,4-dioxane	PFAS
Industrial Pond Sedim	ent Samples									
Sediment	3	Soil	3	3	3	-	3	3	3	3
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1	-	1	1	1	1
Total			7	7	7	0	7	7	7	7
Industrial Pond Surfa	ce Water Samp	les								
Water	1	Water	1	1	1	1	1	1	1	1
Rinsate		Water	1	1	1	1	1	1	1	1
Total			2	2	2	2	2	2	2	2
Test Pit Locations - Su	bsurface Samp	les								
Test Pits	17	Soil	10	10	15	-	4	4	4	4
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1	-	1	1	1	1
Total			14	14	19	0	8	8	8	8
Surface Soil Samples										
Soil	5	Soil	8	8	8	-	4	4	4	4
Duplicate		Soil	1	1	1	-	1	1	1	1
MS/MSD		Soil	2	2	2	-	2	2	2	2
Rinsate		Water	1	1	1		1	1	1	1
Total			12	12	12	0	8	8	8	8
Monitoring Wells										
Monitoring Well	5	Groundwater	5	5	5	5	5	5	3	3
Duplicate		Groundwater	1	1	1	1	1	1	1	1
MS/MSD		Groundwater	2	2	2	2	2	2	2	2
Rinsate		Water	1	1	1	1	1	1	1	1
Trip Blank		Water	1	-	-	-	-	-	-	-
Total			10	9	9	9	9	9	7	7
	то	TAL SAMPLES	VOCs 45	SVOCs 44	METALS 49	METALS 11	PCBs 34	Pest / Herbs 34	1,4-dioxane 32	PFAS 32

Notes:

TCL VOCs - Target Compound List Volatile Organic Compounds. TCL SVOCs - Target Compound List Semi-volatile Organic Compounds.

TAL Metals - Target Analyte List Metals.
TCL PCBs - Target Compound List Polychlorinated Biphenyls.

PFAS - Polyfluoroalkyl Substances


TABLE 2 Sample Container, Volume, Preserving and Holding Time Requirements Octavus Storage NY, LLC 3678 Southwestern Boulevarad, Orchard Park, NY

NYSDEC Brownfield Cleanup Program

			Quantity/		
PARAMETER DESCRIPTION	MATRIX	METHOD NO.	Bottle Type	Preservation	Holding Time
l Samples					
			Encore or Terracore		Freeze within 48 hours
Volatiles, TCL list	Soil	5035/3035A/8260	Samplers	Freeze withint 48 hours	14 days
Semi-Volatiles, TCL list	Soil	8270	(1) 4oz glass jar	Cool, 4 C	14 days
Metals, TAL (no CN)	Soil	6010/7000	(1) 4oz glass jar	none	180 days, Mercury 28 days
PCBs	Soil	8082	(1) 4oz glass jar	Cool, 4 C	365 days/40 days from extraction
Pesticides	Soil	8081	(1) 4oz glass jar	Cool, 4 C	14 days/40 days from extraction
Herbicides	Soil	8151	(1) 4oz glass jar	Cool, 4 C	14 days/40 days from extraction
1,4-Dioxane	Soil	8270D-SIM	(1) 8oz glass jar	Cool, 4 C	14 days
Polyfluoroalkyl Substances (PFAS)	Soil	Method 537	(1) 8oz plastic	Cool, 4 C	28 days
onitoring Wells					
Volatiles, TCL list	Water	8260	(3) 40ml vial	Cool, 4 C, HCL	14 days
Semi-Volatiles, TCL list	Water	8270	(2) 1 liter amber	Cool, 4 C	7 days
PCBs	Water	8082	(2) 1 liter amber	Cool, 4 C	7 days/40 days from extraction
Pesticides	Water	8081	(2) 500ml amber	Cool, 4 C	7 days/40 days from extraction
Herbicides	Water	8151	(2) 1 liter amber	Cool, 4 C	7 days/40 days from extraction
Metals, TAL	Water	6010	(1) 250ml plastic	HNO3	180 days
Mercury, Total	Water	7000	(1) 250ml plastic	HNO3	28 days
Metals, TAL (dissolved) field filtered	Water	6010	(1) 250ml plastic	HNO3	180 days
Mercury, Dissolved	Water	7000	(1) 250ml plastic	HNO3	28 days
1,4-Dioxane	Water	8270D-SIM	(2) 500ml amber	Cool, 4 C	7 days
			2 plastic/1 plastic/1 H20		
Polyfluoroalkyl Substances (PFAS)	Water	Method 537	plastic*	Cool, 4 C	14 days
r Samples - Sub-slab; Ambient Air; So	oil Vapor		•	· · · · · · · · · · · · · · · · · · ·	·
Volatiles	Air	TO-15	1 liter canister	none	30 days

^{*} Three bottles in kit provided by lab. These bottles will be kept separate from other bottles that may have teflon lids

FIGURES



Wittm	an C	SeoScier	nces, PLLC		Hole No.:			
locatio	n·				Date finished:		Sheet 1 o	f 1
Locatio								
Project				Drilling Co.			Weather:	
Proj. M	lgr:			Driller:				
				Drill Rig:		T =		
Depth			Sample	Well Con		Field Analytical	Well	Groundwater and Other
(ft.)	No.	Depth (ft.)	Blows /6"	Deta	ails	Readings	Details	Observations
				1" well completed w/	flush road box			
				Cement/bentonite	miv (1' 2')			
	0	4.0		Cement/bentonite	5 HHX (1 - 2)			
4 _	2	4-8		5	(0) (1)			
				Bentonite pelle				
				1" sch. 40 PVC	riser (0'-5')			
_8 _	3	8-12						
				#0 sand (4	'-15')			
_ 12	4	12-15	NI/A NA/all a a marala ta al co/	1" sch. 40 PVC (.01	0 slot screen).			
			N/A: Well completed w/ geoprobe drill rig					
				Bottom of screen	15 feet bg.			
				Bottom of borehole	e 15 feet bg.			
_ 16 _								
_ 10 _								
– 24 <i>–</i>								
30					1	<u>L</u>		
,	-0-	lit Conner	T- 0	aalby Tuba	 	Backfill W	ell Key:	Cement/
		ock Core:		nelby Tube: Weight of Hammer	1	Grout		Bentonite
		STM D15		TT OIGHT OF TIGHTING	<u> </u>	Sand		Bentonite

Attachment 1

Resumes

3636 N. Buffalo Road, Orchard Park, NY 14127 michelewittmangeo@gmail.com 716-574-1513

Education

B.A., 1994, Geology, State University of New York at Buffalo

B.S., 1994, Social Sciences-Environmental Studies, State University of New York at Buffalo

Professional Registrations

2018, Professional Geologist, New York, #000726

2002, Professional Geologist, Washington, #1772

Affiliations and Certifications

New York State Council of Professional Geologists, Member Buffalo Association of Professional Geologists, Member Air and Waste Management Association of Western New York, Member OSHA 40 Hour 29 CFR 1910. (HAZWOPER) Certification Ms. Wittman is a Professional Geologist with over 24 years of professional experience in conducting a variety of environmental projects for both private and public clients. Clients have included industry, governmental agencies, developers, legal firms, financial institutions, and engineering firms. Project work has included conducting and managing Phase I and Phase II Environmental Site Assessments throughout New York and surrounding states, Brownfield Cleanup Program project investigations and site remediation, hydrogeologic investigations, remedial option evaluation and cost estimating, and remediation of soil and groundwater.

Wittman's responsibilities have ranged from Ms. supervising field and technical activities, completion of field work including soil classification, well installation, collection environmental laboratory samples, excavation oversight; training staff, data analysis, report preparation and review, and client contact. Additionally, responsible for developing and maintaining client relationships, account and project management, bidding, contracting and scheduling and financial management including budgets, proposals, profit/loss assessment. Ms. Wittman has also acted as business manager which included business development and client management, marketing generation of materials; supervising administration staff, and office management.

Ms. Wittman also previously held the position as Assistant Vice President and Environmental Risk Analysis Officer at an international financial institution. During her tenure at this position, Ms. Wittman reviewed hundreds of environmental reports and provided remedial cost estimates to evaluate the potential risk and future losses.

Areas of Specialization

- ✓ Brownfield Cleanup Program
- ✓ Remedial Investigations
- ✓ Feasibility Studies
- ✓ Hydrogeologic Investigations
- ✓ Petroleum and Chemical Bulk Storage
- ✓ Environmental Site Assessments
- ✓ Geologic Evaluations
- ✓ Soil Testing
- ✓ Budgeting & Cost Controls
- ✓ Subcontractor/Crew Management

3636 N. Buffalo Road, Orchard Park, NY 14127 michelewittmangeo@gmail.com 716-574-1513

Environmental Project Highlights

Phase I Environmental Site Assessments – Various

Ms. Wittman has performed, completed, managed or reviewed over 1,500 Phase I Environmental Site Assessments (ESAs) from areas throughout the United States, with focus in the north east and Western New York area. Site assessments have ranged from small, vacant properties, apartment complexes, office buildings, commercial buildings, shopping plaza, automotive dealerships, gasoline stations, as well as small to large 1,000,000-square foot manufacturing and industrial facilities. Ms. Wittman has completed all aspects of Phase I ESAs including site visits, historical review, municipal agency review, database evaluations, and report preparation.

Phase II Environmental Site Assessments - Various

Ms. Wittman has been involved with hundreds of Phase II ESAs at various commercial, manufacturing, industrial and gasoline station properties. Work has included completion of soil borings and/or test pits, installation of groundwater monitoring wells, collection of soil and/or groundwater samples, and vapor intrusion sampling. Ms. Wittman completes data evaluation including with a final report with conclusions and recommendations, if appropriate.

Brownfield Cleanup Program - Commercial Facility, Cheektowaga, New York

Project Manager and Geologist for the investigation and remediation of a former gasoline station for future commercial and residential usage. Site work involved remedial investigation and an interim remedial measure including removal of two underground storage tanks and petroleum impacted soil. Limited remedial work is required and the facility is anticipated to receive certificate of completion in Fall 2018.

Brownfield Cleanup Program - Commercial Facility, Buffalo, New York

Project Manager and Geologist for the investigation and remediation of a former industrial facility for future proposed commercial and residential usage. Site work included remedial investigation, which identified high concentrations of PCBs and hazardous concentrations of lead within site soil. Remedial design included removal and out-of-state disposal of PCB soils and on-site stabilization and off-site disposal of lead impacted soils. Remedial work also included excavation of underground storage tank, and removal of impacted soil and concrete within the building. Additionally, significant asbestos abatement was completed. Remedial design included installation of a sub-slab vapor mitigation system. Facility received certificate of completion in less than nine months from work plan approval.

Remedial Cleanup - Commercial Facility, Amherst, New York

Project Manager and Geologist for the remedial oversight during new building construction, which resulted in identifying former oil/water separator pits, hydraulic lifts, and underground tanks. Each underground structure was evaluated upon discovery, removed, and appropriate samples collected for laboratory analysis. NYSDEC oversight was present during the construction process and one NY Spill was assigned to the site. Upon completion of the project, a final report was done to summarize the findings and the NY Spill was closed.

3636 N. Buffalo Road, Orchard Park, NY 14127 michelewittmangeo@gmail.com 716-574-1513

Remedial Action Plan Evaluation – Former Bulk Petroleum Terminal, Rochester, New York

Developed Remedial Action Plan for former terminal property that underwent extensive subsurface investigations resulting in over 70 borings and 80 soil sample analyses. Initial remedial estimates (by others) included significant soil excavation and remedial costs. Our evaluation included comparison to NYSDEC CP-51 soil guidance for assessment of potential remediation. As such, based on minimal groundwater contamination and identification of significant impacts at greater depths, and negotiation with NYSDEC, no soil remediation was needed.

Management of Environmental Conditions – Retail Gasoline Chain, Western New York

Evaluated environmental concerns associated with 75 different retail gasoline stations. Reviewed regulatory information, previous reports, and data analysis to assess current environmental status. Developed a summary of findings and recommendation of action for each property. Further evaluations included Phase II investigation and continued monitoring of remedial efforts. Developed remedial cost estimate ranges for locations current undergoing remedial work.

Voluntary Cleanup Program - Commercial Facility, Hamburg, New York

Completed a Phase I ESA and identified historical dry cleaner. Conducted investigation and identified contamination beneath the building floor slab and behind the building (i.e. back door). Interim remedial measures (IRM) included soil removal, resulting in approximately 200 tons of soil that was disposed at a hazardous waste landfill. A soil vapor intrusion study was done and identified the presence of compounds To achieve site closure, negotiated a remedial solution that included confirmation sampling of soils around the building structure and installation of a subslab depressurization/vent system.

Contract to Closure, Remedial Activities, Commercial Facility, Rochester, New York.

Two former gasoline stations were located at adjoining properties. Our client wanted to develop the Site for commercial use. Completed a Site Investigation and identified subsurface soil contamination, groundwater contamination and separate phase product. Developed a Remedial Work Plan that included removal of separate phase product and implementation of in-situ chemical oxidation via hydrogen peroxide injections to further reduce contaminants in soil and groundwater. Remedial action also included asbestos abatement and building. The Site received a "no further action" letter and has been developed as a retail bank.

True Bethel Baptist Church – Technical Consultant

Senior Project Manager on the NYSDEC first ever Technical Assistance Grant (TAG) to a community group impacted by a brownfield site. Reviewed site technical documents, attended public meetings and interacted on behalf of the community with NYSDEC and its representatives and contractors on the Site.



Mr. Hanna has over 34 years of experience in environmental pollution control and health/safety services. As principal for Hazard Evaluations, Inc., Mr. Hanna is responsible for all technical services. He specializes in hazardous materials/wastes management, site assessment and remediation, industrial compliance auditing, chemical exposure assessment, safety program development and implementation, and Process Safety Management and Risk Management Planning programs.

Mr. Hanna's career has included over 40 federal/state Superfund projects and over 1,500 due diligence projects. His industrial experience focuses on air, water, waste and chemical management compliance aspects at metal working, wood working, foundry, electroplating, printing and food production facilities.

Education

B.A., 1975, Biology, S.U.C. at Oswego, N.Y.

M.S., 1977, Natural Sciences (Toxicology Concentration), S.U.N.Y. at Buffalo, N.Y.

MEPC, 1982, Pollution Control, Pennsylvania State University

M.S., 1983 Forest Hydrology (Hydrogeology Minor), Pennsylvania State University

Professional Registrations

1985, Certified Hazardous Materials Manager, Senior Level

1989-1998, Registered Environmental Professional

1997, Certified Hazardous Materials Manager, Master Level

Key Skills

- Industrial Emission Permits and Controls
- Hazardous/Solid Waste Management
- Industrial Wastewater Pretreatment and Discharge Permits
- Waste Reduction and Pollution Prevention Programs
- Petroleum and Chemical Bulk Storage
- Industrial Stormwater Management
- Environmental Site Assessments
- Environmental Compliance Assessment
- Industrial Risk Management Program and Audit
- Remedial Investigations
- Brownfield Cleanup Program
- Budgeting & Cost Controls

Affiliations and Certifications

Academy of Hazardous Materials Management, Member Erie County Local Emergency Planning Committee, Member New York Water Environment Association, Member International Institute of Ammonia Refrigeration, Member OSHA 40 Hour 29 CFR 1910. (HAZWOPER) Certification



Environmental Project Highlights

- Performed site characterization for subsurface TCE contamination from historical improper disposal via septic system. Developed Interim Remedial Measures and Remedial Alternatives Reports and Work Plan for this Voluntary Brownfield Cleanup. Installed two banks of piezometers to allow both extraction of contaminated groundwater and injection of Potassium permanganate using continuously operating metering pumps. Recovered over 60 gallons of free product and significantly reduced contamination in groundwater in one year.
- Project Manager for the remediation of numerous (85+) underground petroleum storage tank sites located throughout Western New York. The primary method of remediation has been excavation/removal with appropriate management of tank contents and/or residues, cleaning and scrapping of the tanks and piping, and site restoration. Where petroleum releases were detected, excavation/removal of contaminated soil/fill was completed the majority of the time, with soil management including off-site disposal or on-site bio-treatment. In several cases, on-site vapor extraction systems or chemical oxidation systems with groundwater monitoring have been installed as the recommended remedial method.
- Project Manager for industrial site restoration project which involved the characterization of Leadcontaminated kiln brick surfaces. Appropriate characterization allowed demolition debris from kiln to
 be disposed of in-place on-site as solid waste material as authorized by NYSDEC. Area was then
 backfilled with structural flowable fill to allow reuse of floor space for manufacturing.
- Completed investigation and remediation (excavate and remove) of subsurface Lead contamination at an historical industrial site in Buffalo (NY).
- Project Manager for non-hazardous aspects of site remediation at former Frontier Chemical-Pendleton Site. Remedial tasks included sampling/analysis of wastes, emptying, cleaning and scrapping of bulk storage tanks and collecting/disposing of various on-site residuals.
- Project Manager for the installation of groundwater monitoring wells at AL Tech Specialty Steel's solid
 waste management unit located in Watervliet, NY. Prepared Closure Plan and Bid Specifications for
 the related RCRA surface impoundment. Addressed technical impact of surface run-off from adjacent
 landfill, steep terrain and on-site source for cover material. Prepared response package required by
 NYSDEC regarding the basis of design and construction practices completed during closure.
- Project Manager for the remediation of a cutting oil spill at a Lockport, NY machine shop. Cleanup
 activities included an underground storage tank removal, scarification of surface soils and inoculation
 of contaminated soils with petroleum biodegrading bacteria. Responsibilities included coordination of
 subcontractors, soil sampling, and preparation of report certifying contamination removal.
- Project Manager for industrial site restoration project for solid waste materials abandoned on-site in the on-site production of flowable fill as authorized by a NYSDEC Beneficial Use Determination. Flowable fill produced was used as structural fill to backfill subfloor tanks and large vaults to grade within the facility to allow reuse of the floor space. Tasks included CBS-registered process tank fluid removal and management, basement vault water management, chemical lab packing and disposal, PCBs-contaminated concrete characterization and disposal, UST closure and soil management, scrap and demolition debris management, and subsequent SEQR filing and Phase I Environmental Site Assessment.



Regulatory Compliance Project Highlights

- Project Manager for the development of numerous Process Safety Management and/or Risk Management Plan programs utilizing anhydrous ammonia for refrigeration, including Sorrento Lactalis, Inc.'s South Park (Buffalo, NY), Goshen, NY, Nampa, ID and San Jose, CA facilities, Upstate Niagara Cooperative, Inc.'s Culture (West Seneca, NY), Dale Road (Cheektowaga, NY) and Fulton (Rochester, NY) facilities, as well as Rosina Foods, Inc. (West Seneca, NY), Steuben Foods, Inc. (Elma, NY), Elmhurst Dairy, Inc. (Jamaica, NY), and Sodus Cold Storage, Inc. (Sodus, NY). Responsibilities included coordinating written program preparation, Process Hazard Analysis development, preparing release scenarios, evaluating and upgrading SOPs, developing MOC methods, etc.
- Provided consulting services to over 75 facilities nationwide regarding SARA Title III reporting requirements. Services included regulations and process reviews, mass balance calculations, purchasing and process data evaluation, database development and USEPA Tier Two and Form R preparation.
- Project Manager for numerous environmental compliance audits including, Mod-Pac Corp., Buffalo, NY (commercial printing), Sahlen Packing Co., Inc., Buffalo, NY (meat packing), Upstate Niagara Cooperative, Inc., Buffalo, NY (dairy products), MoldTech, Inc., Lancaster, NY (plastics), Sorrento Lactalis, Inc., Buffalo, NY (cheese manufacturing), Chautauqua Hardware Corp., Jamestown, NY (brass hardware), Thomson Professional Publishing, Webster, NY (printed media), Buffalo China, Inc., Buffalo, NY (lead glazed china), Brainerd Manufacturing Co., East Rochester, NY (electroplating and finishing), Falconer Die Casting Co., Inc., Lakewood, NY (aluminum and zinc casting), and Jensen Fittings Corp., North Tonawanda, NY (stainless pipe fittings). These audits emphasized the inspection of all manufacturing operations, hazardous materials and hazardous waste handling, wastewater treatment operations, air emissions and facility records to evaluate current practices with regard to RCRA, SARA, New York State Parts 200 (air), 360 (solid waste) and 370 (hazardous waste) regulations, USEPA Categorical Pretreatment Standards, UIC NESHAP & CFATS regulations, New York State SPDES regulations, and local sewer authority and fire and building department codes.
- Oversaw the modification of an industrial wastewater pre-treatment system for Whiting Door Manufacturing. Evaluated plant manufacturing wastewater sources, modified existing pretreatment system, developed wastewater pretreatment schedule, and completed wastewater discharge monitoring. Developed a Toxic Organics Management Plan to reduce cost of wastewater monitoring. Evaluated and assisted with the revision of municipal Industrial User Permit.
- Project Manager for Title V Clean Air Act permit development for Whiting Door Manufacturing Corp., Dinaire, Inc., Metalico Aluminum Recovery, Inc. and Flexo Transparent, Inc. Continued services include annual emission statements, 12-month rolling emissions determinations and semi-annual compliance reporting.
- Project Manager for Clean Air Act and/or NYSDEC Part 228 determinations and State Air Facility Permit or Air Facility Registration development for numerous industrial clients including Niagara Ceramics Corporation, Buffalo Metal Casting Co., Inc., ITT Standard/XYLEM, Metalico Rochester, Inc., Ulrich Planfiling Equipment Corp., United Silicone, Inc., U.S. Chrome Corp., Metalico Aluminum Recovery, Inc., Truck-Lite Co., Inc., Jensen Fittings Corp., API Delavan, Inc., Tapecon Inc., Dura-Plating, Inc., Buffalo China, Inc., Forsyth Industries, Inc., Jamestown Laminating Co., Classic Brass Inc., Ivaco Steel Processing (New York), LLC, Innovative Tool & Machine Co., Inc., and Whiting Door Manufacturing, Inc.

Attachment 2

SOP-Technical Guides



Sampling for 1,4-Dioxane and Per- and Polyfluoroalkyl Substances (PFAS) Under DEC's Part 375 Remedial Programs

Objective

The Department of Environmental Conservation (DEC) is requiring sampling of all environmental media and subsequent analysis for the emerging contaminants 1,4-Dioxane and PFAS as part of all remedial programs implemented under 6 NYCRR Part 375, as further described in the guidance below.

Sample Planning

The number of samples required for emerging contaminant analyses is to be the same number of samples where "full TAL/TCL sampling" would typically be required in an investigation or remedial action compliance program.

Upon a new site being brought into any program (e.g., SSF, BCP), PFAS and 1,4-dioxane will be incorporated into the investigation of potentially affected media, including soil, groundwater, surface water, and sediment as an addition to the standard "full TAL/TCL sampling." Biota sampling may be necessary based upon the potential for biota to be affected as determined pursuant to a Fish and Wildlife Impact analysis. Soil vapor sampling for PFAS and 1,4-dioxane is not required.

Upon an emerging contaminant being identified as a contaminant of concern (COC) for a site, those compounds must be assessed as part of the remedy selection process in accordance with Part 375 and DER-10 and included as part of the monitoring program upon entering the site management phase.

Soil imported to a site for use in a soil cap, soil cover, or as backfill must be sampled for 1,4-dioxane and PFAS contamination in general conformance with DER-10, section 5.4(e). Assessment of the soil data will be made on a site-specific basis to determine appropriateness for use.

The work plan should explicitly describe analysis and reporting requirements, including laboratory analytical procedures for modified methods discussed below.

Analysis and Reporting

Labs should provide a full category B deliverable, and a DUSR should be prepared by an independent 3rd party data validator. QA/QC samples should be collected as required in DER-10, Section 2.3(c). The electronic data submission should meet the requirements provided at: https://www.dec.ny.gov/chemical/62440.html.

<u>PFAS analysis and reporting:</u> DEC has developed a *PFAS Analyte List* (below) for remedial programs. It is expected that reported results for PFAS will include, at a minimum, all the compounds listed. If lab and/or matrix specific issues are encountered for any compounds, the DEC project manager, in consultation with the DEC remedial program chemist, will make case-by-case decisions as to whether certain analytes may be temporarily or permanently discontinued from analysis at each site.

Currently, ELAP does not offer certification for PFAS compounds in matrices other than finished drinking water. However, laboratories analyzing environmental samples (e.g., soil, sediments, and groundwater) are required by DER to hold ELAP certification for PFOA and PFOS in drinking water by EPA Method 537 or ISO 25101.



Modified EPA Method 537 is the preferred method to use for environmental samples due to its ability to achieve very low detection limits. Reporting limits for PFAS in groundwater and soil are to be 2 ng/L (ppt) and 1 ug/kg (ppb), respectively. If contract labs or work plans submitted by responsible parties indicate that they are not able to achieve these reporting limits for the entire list of 21 PFAS, site-specific decisions will need to be made by the DEC project manager in consultation with the DEC remedial program chemist. Note: Reporting limits for PFOA and PFOS in groundwater should not exceed 2 ng/L.

Additional laboratory methods for analysis of PFAS may be warranted at a site. These methods include Synthetic Precipitation Leaching Procedure (SPLP) by EPA Method 1312 and Total Oxidizable Precursor Assay (TOP Assay).

SPLP is a technique for determining the potential for chemicals in soil to leach to groundwater and may be helpful in determining the need for addressing PFAS-containing soils or other solid material as part of the remedy. SPLP sampling need not be considered if there are no elevated PFAS levels in groundwater. If elevated levels of PFAS are detected in water, and PFAS are also seen in soil, then an SPLP test should be considered to better understand the relationship between the PFAS in the two media.

The TOP Assay can assist in determining the potential PFAS risk at a site. For example, some polyfluoroalkyl substances may transform to form perfluoroalkyl substances, resulting in an increase in perfluoroalkyl substance concentrations as contaminated groundwater moves away from the site. To conceptualize the amount and type of oxidizable perfluoroalkyl substances which could be liberated in the environment, a "TOP Assay" analysis can be performed, which approximates the maximum concentration of perfluoroalkyl substances that could be generated if all polyfluoroalkyl substances were oxidized.

PFAS-containing materials can be made up of per- and polyfluoroalkyl substances that are not analyzable by routine analytical methodology (LC-MS/MS). The TOP assay converts, through oxidation, polyfluoroalkyl substances (precursors) into perfluoroalkyl substances that can be detected by current analytical methodology. Please note that analysis of highly contaminated samples, such as those from an AFFF site, can result in incomplete oxidation of the samples and an underestimation of the total perfluoroalkyl substances. Please consult with a DEC remedial program chemist for assistance interpreting the results.

1,4-Dioxane analysis and reporting: The reporting limit for 1,4-dioxane in groundwater should be no higher than 0.35 μg/L (ppb) and no higher than 0.1 mg/kg (ppm) in soil. Although ELAP offers certification for both EPA Method 8260 SIM and EPA Method 8270 SIM in waters, DER is advising the use of Method 8270 SIM because it provides a more robust extraction procedure, uses a larger sample volume, and is less vulnerable to interference from chlorinated solvents. The analysis currently performed for SVOCs in soil is adequate for evaluation of 1,4-dioxane in soil, which already has an established SCO.



Refinement of sample analyses

As with other contaminants that are analyzed for at a site, the emerging contaminant analyte list may be refined for future sampling events based on investigative findings. Initially, however, sampling using this PFAS Analyte List and 1,4-dioxane is needed to understand the nature of contamination.

PFAS Analyte List

Group	Chemical Name	Abbreviation	CAS Number
	Perfluorobutanesulfonic acid	PFBS	375-73-5
	Perfluorohexanesulfonic acid	PFHxS	355-46-4
Perfluoroalkyl sulfonates	Perfluoroheptanesulfonic acid	PFHpS	375-92-8
dimoniation	Perfluorooctanessulfonic acid	PFOS	1763-23-1
	Perfluorodecanesulfonic acid	PFDS	335-77-3
	Perfluorobutanoic acid	PFBA	375-22-4
	Perfluoropentanoic acid	PFPeA	2706-90-3
	Perfluorohexanoic acid	PFHxA	307-24-4
	Perfluoroheptanoic acid	PFHpA	375-85-9
5 (1 11 1	Perfluorooctanoic acid	PFOA	335-67-1
Perfluoroalkyl carboxylates	Perfluorononanoic acid	PFNA	375-95-1
darboxylatoo	Perfluorodecanoic acid	PFDA	335-76-2
	Perfluoroundecanoic acid	PFUA/PFUdA	2058-94-8
	Perfluorododecanoic acid	PFDoA	307-55-1
	Perfluorotridecanoic acid	PFTriA/PFTrDA	72629-94-8
	Perfluorotetradecanoic acid	PFTA/PFTeDA	376-06-7
Fluorinated Telomer	d Telomer 6:2 Fluorotelomer sulfonate		27619-97-2
Sulfonates 8:2 Fluorotelomer sulfonate		8:2 FTS	39108-34-4
Perfluorooctane- sulfonamides	Perfluroroctanesulfonamide	FOSA	754-91-6
Perfluorooctane-	N-methyl perfluorooctanesulfonamidoacetic acid	N-MeFOSAA	2355-31-9
sulfonamidoacetic acids	N-ethyl perfluorooctanesulfonamidoacetic acid	N-EtFOSAA	2991-50-6



Laboratory Guidance for Analysis of Per- and Polyfluoroalkyl Substances (PFAS) in Non-Potable Water and Solids

The Division of Environmental Remediation (DER) developed the following guidance for laboratories submitting PFAS data to DER. If laboratories cannot comply with any of the following requirements, they must contact Dana Maikels at dana.maikels@dec.ny.gov prior to analysis of samples.

- Standards containing both branched and linear isomers must be used when standards are commercially available. Currently, quantitative standards are available for PFHxS, PFOS, NMeFOSAA, and NEtFOSAA. All isomer peaks present in the standard must be integrated and the areas summed. Samples must be integrated in the same manner as the standards.
 - Since a quantitative standard does not exist for branched isomers of PFOA, the instrument must be calibrated using just the linear isomer and a technical (qualitative) PFOA standard must be used to identify the retention time of the branched PFOA isomers in the sample. The total response of PFOA branched and linear isomers must be integrated in the samples and quantitated using the calibration curve of the linear standard.
- Quantifier and qualifier ions must be monitored for all target analytes (PFPeA and PFBA are an
 exception). The ratio of quantifier ion response to qualifier ion response must be calculated for
 each target analyte and the ratio compared to standards. Lab derived criteria can be used to
 determine if the ratios are acceptable.
- 3. The ion transitions below must be used for the following PFASs:

PFOA	413 > 369
PFOS	499 > 80
PFHxS	399 > 80
PFBS	299 > 80
6:2 FTS	427 > 407
8:2 FTS	527 > 507
NEtFOSAA	584 > 419
NMeFOSAA	570 > 419

- 4. For all target analyte ions used for quantification, signal to noise ratio must be 3:1 or greater.
- 5. For water samples, the entire sample bottle must be extracted, and the sample bottle rinsed with appropriate solvent to remove any residual PFAS.
- 6. Detections below the reporting limit should be reported and qualified with a J qualifier.

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:**23528** Revision 12 Published Date: 2/22/2019 3:48:15 PM

Page 1 of 28

Determination of Selected Perfluorinated Alkyl Substances by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry Isotope Dilution (LC/MS/MS)

Reference: EPA Method 537, Version 1.1, September

2009, EPA Document #: EPA/600/R-08/09

EPA Method 537.1, Version 1, November 2018, EPA

Document #: EPA/600/R-18/352

Department of Defense, Quality Systems Manual for Environmental Laboratories, Version 5.2, .2019

1. Scope and Application

Matrices: Drinking water, Non-potable Water, and Soil Matrices

Definitions: Refer to Alpha Analytical Quality Manual.

- 1.1 This is a liquid chromatography/tandem mass spectrometry (LC/MS/MS) method for the determination of selected perfluorinated alkyl substances (PFAS) in Non-Drinking Water and soil Matrices. Accuracy and precision data have been generated in reagent water, and finished ground and surface waters for the compounds listed in Table 1.
- 1.2 The data report packages present the documentation of any method modification related to the samples tested. Depending upon the nature of the modification and the extent of intended use, the laboratory may be required to demonstrate that the modifications will produce equivalent results for the matrix. Approval of all method modifications is by one or more of the following laboratory personnel before performing the modification: Area Supervisor, Department Supervisor, Laboratory Director, or Quality Assurance Officer.
- 1.3 This method is restricted to use by or under the supervision of analysts experienced in the operation of the LC/MS/MS and in the interpretation of LC/MS/MS data. Each analyst must demonstrate the ability to generate acceptable results with this method by performing an initial demonstration of capability.

2. Summary of Method

2.1 A 250-mL water sample is fortified with extracted internal standards (EIS) and passed through a solid phase extraction (WAX) cartridge containing a mixed mode, Weak Anion Exchange, reversed phase, water-wettable polymer to extract the method analytes and isotopically-labeled compounds. The compounds are eluted from the solid phase in two fractions with methanol followed by a small amount of 2% ammonium hydroxide in methanol solution. The extract is concentrated with nitrogen in a heated water bath, and then adjusted to a 1-mL volume with 80:20% (vol/vol) methanol:water. A 3 µl injection is made into an LC equipped with a C18 column that is interfaced to an MS/MS. The analytes are separated and identified by comparing the acquired mass spectra and retention times to reference spectra and retention times for calibration standards acquired under identical LC/MS/MS conditions. The concentration of each analyte is determined by using the isotope dilution technique. Extracted Internal Standards (EIS) analytes are used to monitor the extraction efficiency of the method analytes.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 2 of 28

2.2 Method Modifications from Reference

None.

Table 1

Parameter	Acronym	CAS					
PERFLUOROALKYL ETHER CARBOXYLIC ACIDS (PFECAs)							
Tetrafluoro-2-(heptafluoropropoxy)propanoic acid	HFPO-DA	62037-80-3					
4,8-dioxa-3H-perfluorononanoic acid	ADONA	919005-14-4					
PERFLUOROALKYLCARBOXILIC ACIDS (PFCAs)							
Perfluorobutanoic acid	PFBA	375-22-4					
Perfluoropentanoic acid	PFPeA	2706-90-3					
Perfluorohexanoic acid	PFHxA *	307-24-4					
Perfluoroheptanoic acid	PFHpA *	375-85-9					
Perfluorooctanoic acid	PFOA *	335-67-1					
Perfluorononanoic acid	PFNA *	375-95-1					
Perfluorodecanoic acid	PFDA *	335-76-2					
Perfluoroundecanoic acid	PFUnA *	2058-94-8					
Perfluorododecanoic acid	PFDoA *	307-55-1					
Perfluorotridecanoic acid	PFTrDA *	72629-94-8					
Perfluorotetradecanoic acid	PFTA *	376-06-7					
Perfluorohexadecanoic acid	PFHxDA	67905-19-5					
Perfluorooctadecanoic acid	PFODA	16517-11-6					
PERFLUOROALKYLSULFONATES (PFASs)	·	•					
Perfluorobutanesulfonic acid	PFBS *	375-73-5					
Perfluoropentanesulfonic acid	PFPeS	2706-91-4					
Perfluorohexanesulfonic acid	PFHxS *	355-46-4					
Perfluoroheptanesulfonic acid	PFHpS	375-92-8					
Perfluorooctanesulfonic acid	PFOS *	1763-23-1					
Perfluorononanesulfonic acid	PFNS	68259-12-1					
Perfluorodecanesulfonic acid	PFDS	335-77-3					
Perfluorododecanesulfonic acid	PFDoS	79780-39-5					

^{*} also reportable via the standard 537 method

Alpha Analytical, Inc.ID No.:23528Facility: Mansfield, MARevision 12Department: SemivolatilesPublished Date: 2/22/2019 3:48:15 PM

Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 3 of 28

Table 1 Cont.

Acronym	CAS						
CHLORO-PERFLUOROALKYLSULFONATE							
11CI- PF3OUdS	763051-92-9						
9CI-PF3ONS	756426-58-1						
•							
PFOSA	754-91-6						
NMeFOSA	31506-32-8						
NEtFOSA	4151-50-2						
•							
4:2FTS	27619-93-8						
6:2FTS	27619-97-2						
8:2FTS	39108-34-4						
10:2FTS	120226-60-0						
S							
NMeFOSAA *	2355-31-9						
NEtFOSAA *	2991-50-6						
ANOLS (FOSEs)							
NMeFOSE	24448-09-7						
NEtFOSE	1691-99-2						
	11CI- PF3OUdS 9CI-PF3ONS PFOSA NMeFOSA NEtFOSA 4:2FTS 6:2FTS 8:2FTS 10:2FTS 0S NMeFOSAA* NEtFOSAA* NEtFOSAA*						

^{*} also reportable via the standard 537 method

3. Reporting Limits

The reporting limit for PFAS's is 2 ng/L for aqueous samples (20 ng/L for HFPO-DA) and 1 ng/g (10 ng/g for HFPO-DA) for soil samples.

4. Interferences

- 4.1 PFAS standards, extracts and samples should not come in contact with any glass containers or pipettes as these analytes can potentially adsorb to glass surfaces. PFAS analyte and EIS standards commercially purchased in glass ampoules are acceptable; however, all subsequent transfers or dilutions performed by the analyst must be prepared and stored in polypropylene containers.
- 4.2 Method interferences may be caused by contaminants in solvents, reagents (including reagent water), sample bottles and caps, and other sample processing hardware that lead to discrete artifacts and/or elevated baselines in the chromatograms. The method analytes in this method can also be found in many common laboratory supplies and equipment, such

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Published Date: 2/22/2019 3:48:15 PM

Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 4 of 28

as PTFE (polytetrafluoroethylene) products, LC solvent lines, methanol, aluminum foil, SPE sample transfer lines, etc. All items such as these must be routinely demonstrated to be free from interferences (less than 1/3 the RL for each method analyte) under the conditions of the analysis by analyzing laboratory reagent blanks as described in Section 9.2. Subtracting blank values from sample results is not permitted.

- 4.3 Matrix interferences may be caused by contaminants that are co-extracted from the sample. The extent of matrix interferences will vary considerably from source to source, depending upon the nature of the water. Humic and/or fulvic material can be co-extracted during SPE and high levels can cause enhancement and/or suppression in the electrospray ionization source or low recoveries on the SPE sorbent. Total organic carbon (TOC) is a good indicator of humic content of the sample.
- **4.4** SPE cartridges can be a source of interferences. The analysis of field and laboratory reagent blanks can provide important information regarding the presence or absence of such interferences. Brands and lots of SPE devices should be tested to ensure that contamination does not preclude analyte identification and quantitation.

5. Health and Safety

- 5.1 The toxicity or carcinogenicity of each reagent and standard used in this method is not fully established; however, each chemical compound should be treated as a potential health hazard. From this viewpoint, exposure to these chemicals must be reduced to the lowest possible level by whatever means available. A reference file of material safety data sheets is available to all personnel involved in the chemical analysis. Additional references to laboratory safety are available in the Chemical Hygiene Plan.
- 5.2 All personnel handling environmental samples known to contain or to have been in contact with municipal waste must follow safety practices for handling known disease causative agents.
- 5.3 PFOA has been described as "likely to be carcinogenic to humans." Pure standard materials and stock standard solutions of these method analytes should be handled with suitable protection to skin and eyes, and care should be taken not to breathe the vapors or ingest the materials.

6. Sample Collection, Preservation, Shipping and Handling

6.1 Sample Collection for Aqueous Samples

- **6.1.1** Samples must be collected in two (2) 250-mL high density polyethylene (HDPE) container with an unlined plastic screw cap.
- 6.1.2 The sample handler must wash their hands before sampling and wear nitrile gloves while filling and sealing the sample bottles. PFAS contamination during sampling can occur from a number of common sources, such as food packaging and certain foods and beverages. Proper hand washing and wearing nitrile gloves will aid in minimizing this type of accidental contamination of the samples.
- **6.1.3** Open the tap and allow the system to flush until the water temperature has stabilized (approximately 3 to 5 min). Collect samples from the flowing system.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 5 of 28

6.1.4 Fill sample bottles. Samples do not need to be collected headspace free.

- **6.1.5** After collecting the sample and cap the bottle. Keep the sample sealed from time of collection until extraction.
- 6.1.6 Field Reagent Blank (FRB)
 - 6.1.6.1 A FRB must be handled along with each sample set. The sample set is composed of samples collected from the same sample site and at the same time. At the laboratory, fill the field blank sample bottle with reagent water and preservatives, seal, and ship to the sampling site along with the sample bottles. For each FRB shipped, an empty sample bottle (no preservatives) must also be shipped. At the sampling site, the sampler must open the shipped FRB and pour the reagent water into the empty shipped sample bottle, seal and label this bottle as the FRB. The FRB is shipped back to the laboratory along with the samples and analyzed to ensure that PFAS's were not introduced into the sample during sample collection/handling.

The reagent water used for the FRBs must be initially analyzed for method analytes as a MB and must meet the MB criteria in Section 9.2.1 prior to use. This requirement will ensure samples are not being discarded due to contaminated reagent water rather than contamination during sampling.

6.2 Sample Collection for Soil and Sediment samples.

Grab samples are collected in polypropylene containers. Sample containers and contact surfaces containing PTFE shall be avoided.

6.3 Sample Preservation

Not applicable.

6.4 Sample Shipping

Samples must be chilled during shipment and must not exceed 10 °C during the first 48 hours after collection. Sample temperature must be confirmed to be at or below 10 °C when the samples are received at the laboratory. Samples stored in the lab must be held at or below 6 °C until extraction, but should not be frozen.

NOTE: Samples that are significantly above 10° C, at the time of collection, may need to be iced or refrigerated for a period of time, in order to chill them prior to shipping. This will allow them to be shipped with sufficient ice to meet the above requirements.

6.5 Sample Handling

- 6.5.1 Holding Times
 - 6.5.1.1 Water samples should be extracted as soon as possible but must be extracted within 14 days. Soil samples should be extracted within 28 days. Extracts are stored at < 10 ° C and analyzed within 28 days after extraction.</p>

7. Equipment and Supplies

Printouts of this document may be out of date and should be considered uncontrolled. To accomplish work, the published version of the document should be viewed online.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 6 of 28

7.1 SAMPLE CONTAINERS – 250-mL high density polyethylene (HDPE) bottles fitted with unlined screw caps. Sample bottles must be discarded after use.

- 7.2 POLYPROPYLENE BOTTLES 4-mL narrow-mouth polypropylene bottles.
- **7.3** CENTRIFUGE TUBES 50-mL conical polypropylene tubes with polypropylene screw caps for storing standard solutions and for collection of the extracts.
- **7.4** AUTOSAMPLER VIALS Polypropylene 0.7-mL autosampler vials with polypropylene caps.
 - 7.4.1 NOTE: Polypropylene vials and caps are necessary to prevent contamination of the sample from PTFE coated septa. However, polypropylene caps do not reseal, so evaporation occurs after injection. Thus, multiple injections from the same vial are not possible.
- **7.5** POLYPROPYLENE GRADUATED CYLINDERS Suggested sizes include 25, 50, 100 and 1000-mL cylinders.
- **7.6** Auto Pipets Suggested sizes include 5, 10, 25, 50, 100, 250, 500, 1000, 5000 and 10,000-μls.
- 7.7 PLASTIC PIPETS Polypropylene or polyethylene disposable pipets.
- 7.8 ANALYTICAL BALANCE Capable of weighing to the nearest 0.0001 g.
- 7.9 SOLID PHASE EXTRACTION (SPE) APPARATUS FOR USING CARTRIDGES
 - **7.9.1** SPE CARTRIDGES 0.5 g SPE cartridges containing a reverse phase copolymer characterized by a weak anion exchanger (WAX) sorbent phase.
 - 7.9.2 VACUUM EXTRACTION MANIFOLD A manual vacuum manifold with large volume sampler for cartridge extractions, or an automatic/robotic sample preparation system designed for use with SPE cartridges, may be used if all QC requirements discussed in Section 9 are met. Extraction and/or elution steps may not be changed or omitted to accommodate the use of an automated system. Care must be taken with automated SPE systems to ensure the PTFE commonly used in these systems does not contribute to unacceptable analyte concentrations in the MB (Sect. 9.2.1).
 - 7.9.3 SAMPLE DELIVERY SYSTEM Use of a polypropylene transfer tube system, which transfers the sample directly from the sample container to the SPE cartridge, is recommended, but not mandatory. Standard extraction manifolds come equipped with PTFE transfer tube systems. These can be replaced with 1/8" O.D. x 1/16" I.D. polypropylene or polyethylene tubing cut to an appropriate length to ensure no sample contamination from the sample transfer lines. Other types of non-PTFE tubing may be used provided it meets the MB (Sect. 9.2.1) and LCS (Sect. 9.3) QC requirements. The PTFE transfer tubes may be used, but an MB must be run on each PFTE transfer tube and the QC requirements in Section 13.2.2 must be met. In the case of automated SPE, the removal of PTFE lines may not be feasible; therefore, MBs will need to be rotated among the ports and must meet the QC requirements of Sections 13.2.2 and 9.2.1.
- **7.10** Extract Clean-up Cartridge 250 mg 6ml SPE Cartridge containing graphitized polymer carbon

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 7 of 28

7.11 EXTRACT CONCENTRATION SYSTEM – Extracts are concentrated by evaporation with nitrogen using a water bath set no higher than 65 °C.

- **7.12** LABORATORY OR ASPIRATOR VACUUM SYSTEM Sufficient capacity to maintain a vacuum of approximately 10 to 15 inches of mercury for extraction cartridges.
- 7.13 LIQUID CHROMATOGRAPHY (LC)/TANDEM MASS SPECTROMETER (MS/MS) WITH DATA SYSTEM
 - 7.13.1 LC SYSTEM Instrument capable of reproducibly injecting up to 10-µL aliquots, and performing binary linear gradients at a constant flow rate near the flow rate used for development of this method (0.4 mL/min). The LC must be capable of pumping the water/methanol mobile phase without the use of a degasser which pulls vacuum on the mobile phase bottle (other types of degassers are acceptable). Degassers which pull vacuum on the mobile phase bottle will volatilize the ammonium acetate mobile phase causing the analyte peaks to shift to earlier retention times over the course of the analysis batch. The usage of a column heater is optional.

NOTE: During the course of method development, it was discovered that while idle for more than one day, PFAS's built up in the PTFE solvent transfer lines. To prevent long delays in purging high levels of PFAS's from the LC solvent lines, they were replaced with PEEK tubing and the PTFE solvent frits were replaced with stainless steel frits. It is not possible to remove all PFAS background contamination, but these measures help to minimize their background levels.

- 7.13.2 LC/TANDEM MASS SPECTROMETER The LC/MS/MS must be capable of negative ion electrospray ionization (ESI) near the suggested LC flow rate of 0.4 mL/min. The system must be capable of performing MS/MS to produce unique product ions for the method analytes within specified retention time segments. A minimum of 10 scans across the chromatographic peak is required to ensure adequate precision.
- 7.13.3 DATA SYSTEM An interfaced data system is required to acquire, store, reduce, and output mass spectral data. The computer software should have the capability of processing stored LC/MS/MS data by recognizing an LC peak within any given retention time window. The software must allow integration of the ion abundance of any specific ion within specified time or scan number limits. The software must be able to calculate relative response factors, construct linear regressions or quadratic calibration curves, and calculate analyte concentrations.
- 7.13.4 ANALYTICAL COLUMN An LC BEH C_{18} column (2.1 x 50 mm) packed with 1.7 μ m d_p C_{18} solid phase particles was used. Any column that provides adequate resolution, peak shape, capacity, accuracy, and precision (Sect. 9) may be used.

8. Reagents and Standards

- 8.1 GASES, REAGENTS, AND SOLVENTS Reagent grade or better chemicals should be used.
 - 8.1.1 REAGENT WATER Purified water which does not contain any measurable quantities of any method analytes or interfering compounds greater than 1/3 the RL for each method analyte of interest. Prior to daily use, at least 3 L of reagent water should be flushed from the purification system to rinse out any build-up of analytes in the system's tubing.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 8 of 28

8.1.2 METHANOL (CH₃OH, CAS#: 67-56-1) – High purity, demonstrated to be free of analytes and interferences.

- **8.1.3** AMMONIUM ACETATE (NH₄C₂H₃O₂, CAS#: 631-61-8) High purity, demonstrated to be free of analytes and interferences.
- **8.1.4** ACETIC ACID (H₃CCOOH, CAS#: 64-19-7) High purity, demonstrated to be free of analytes and interferences.
- **8.1.5** 1M AMMONIUM ACETATE/REAGENT WATER High purity, demonstrated to be free of analytes and interferences.
- 8.1.6 2mM AMMONIUM ACETATE/METHANOL:WATER (5:95) To prepare, mix 2 ml of 1M AMMONIUM ACETATE,1 ml ACETIC ACID and 50 ml METHANOL into I Liter of REAGENT WATER.
- **8.1.7** Methanol/Water (80:20) To prepare a 1 Liter bottle, mix 200 ml of REAGENT WATER with 800 ml of METHANOL.
- **8.1.8** AMMONIUM HYDROXIDE (NH₃, CAS#: 1336-21-6) High purity, demonstrated to be free of analytes and interferences.
- **8.1.9** Sodium Acetate (NaOOCCH₃, CAS#: 127-09-3) High purity, demonstrated to be free of analytes and interferences.
- 8.1.10 25 mM Sodium Acetate Buffer To prepare 250mls, dissolve .625 grams of sodium acetate into 100 mls of reagent water. Add 4 mls Acetic Acid and adjust the final volume to 250 mls with reagent water.
- 8.1.11 NITROGEN Used for the following purposes: Nitrogen aids in aerosol generation of the ESI liquid spray and is used as collision gas in some MS/MS instruments. The nitrogen used should meet or exceed instrument manufacturer's specifications. In addition, Nitrogen is used to concentrate sample extracts (Ultra High Purity or equivalent).
- 8.1.12 ARGON Used as collision gas in MS/MS instruments. Argon should meet or exceed instrument manufacturer's specifications. Nitrogen gas may be used as the collision gas provided sufficient sensitivity (product ion formation) is achieved.
- 8.2 STANDARD SOLUTIONS When a compound purity is assayed to be 96% or greater, the weight can be used without correction to calculate the concentration of the stock standard. PFAS analyte and IS standards commercially purchased in glass ampoules are acceptable; however, all subsequent transfers or dilutions performed by the analyst must be prepared and stored in polypropylene containers. Standards for sample fortification generally should be prepared in the smallest volume that can be accurately measured to minimize the addition of excess organic solvent to aqueous samples.

NOTE: Stock standards and diluted stock standards are stored at ≤4 °C.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 9 of 28

8.2.1 ISOTOPE DILUTION Extracted Internal Standard (ID EIS) STOCK SOLUTIONS
 ID EIS stock standard solutions are stable for at least 6 months when stored at 4 °C. The stock solution is purchased at a concentration of 1000 ng/mL.

8.2.2 ISOTOPE DILUTION Extracted Internal Standard PRIMARY DILUTION STANDARD (ID EIS PDS) – Prepare the ID EIS PDS at a concentration of 500 ng/mL. The ID PDS is prepared in 80:20% (vol/vol) methanol:water. The ID PDS is stable for 6 months when stored at ≤4 °C.

Table 2

Isotope Labeled	Conc. of EIS	Vol. of EIS Stock	Final Vol. of EIS	Final Conc. of
Standard	Stock (ng/mL)	(mL)	PDS (mL)	EIS PDS (ng/mL)
M4PFBA	1000	1.0	2.0	500
M5PFPeA	1000	1.0	2.0	500
M5PFHxA	1000	1.0	2.0	500
M4PFHpA	1000	1.0	2.0	500
M8PFOA	1000	1.0	2.0	500
M9PFNA	1000	1.0	2.0	500
M6PFDA	1000	1.0	2.0	500
M7PFUdA	1000	1.0	2.0	500
MPFDoA	1000	1.0	2.0	500
M2PFTeDA	1000	1.0	2.0	500
M2PFHxDA	50,000	.02	2.0	500
d3-N-MeFOSA	50,000	.02	2.0	500
d5-N-EtFOSA	50,000	.02	2.0	500
d7-N-MeFOSE	50,000	.02	2.0	500
d9-N-EtFOSE	50,000	.02	2.0	500
M8FOSA	1000	1.0	2.0	500
d3-N-MeFOSAA	1000	1.0	2.0	500
d5-N-EtFOSAA	1000	1.0	2.0	500
M3PFBS	929	1.0	2.0	464.5
M3PFHxS	946	1.0	2.0	473
M8PFOS	957	1.0	2.0	478.5
M2-4:2FTS	935	1.0	2.0	467.5
M2-6:2FTS	949	1.0	2.0	474.5
M2-8:2FTS	958	1.0	2.0	479
M3HFPO-DA	50,000	.4	2.0	10,000

- 8.2.3 ANALYTE STOCK STANDARD SOLUTION Analyte stock standards are stable for at least 6 months when stored at 4 °C. When using these stock standards to prepare a PDS, care must be taken to ensure that these standards are at room temperature and adequately vortexed.
- 8.2.4 Analyte Secondary Spiking Standard Prepare the spiking solution of additional add on components for project specific requirements only. ANALYTE PRIMARY SPIKING STANDARD Prepare the spiking standard at a concentration of 500 ng/mL in methanol. The spiking standard is stable for at least two months when stored in polypropylene centrifuge tubes at room temperature.

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 10 of 28

Table 3

Analyte	Conc. of IS	Vol. of IS	Final Vol. of IS PDS	Final Conc. of IS
	Stock (ng/mL)	Stock (mL)	(mL)	PDS (ng/mL)
PFBA	2000	1	4	500
PFPeA	2000	1	4	500
PFHxA	2000	1	4	500
PFHpA	2000	1	4	500
PFOA	2000	1	4	500
PFNA	2000	1	4	500
PFDA	2000	1	4	500
PFUdA	2000	1	4	500
PFDoA	2000	1	4	500
PFTrDA	2000	1	4	500
PFTeDA	2000	1	4	500
FOSA	2000	1	4	500
N-MeFOSAA	2000	1	4	500
N-EtFOSAA	2000	1	4	500
L-PFBS	1770	1	4	442.5
L-PFPeS	1880	1	4	470
L-PFHxSK	1480	1	4	370
Br-PFHxSK	344	1	4	86
L-PFHpS	1900	1	4	475
L-PFOSK	1460	1	4	365
Br-PFOSK	391	1	4	97.75
L-PFNS	1920	1	4	480
L-PFDS	1930	1	4	482.5
4:2FTS	1870	1	4	467.5
6:2FTS	1900	1	4	475
8:2FTS	1920	1	4	480

8.2.5 Analyte Secondary Spiking Standard Prepare the spiking solution of additional add on components for project specific requirements only.

Table 4

Analyte	Conc. of IS	Vol. of IS Stock	Final Vol. of IS PDS	Final Conc. of IS
	Stock (ng/mL)	(mL)	(mL)	PDS (ng/mL)
ADONA	2000	1	4	500
PFHxDA	2000	1	4	500
PFODA	2000	1	4	500
HFPO-DA	100,000	.4	4	10,000
9CIPF3ONS	50,000	0.04	4	500
11CIPF3OUdS	50,000	0.04	4	500

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Page 11 of 28

8.2.6 LOW, MEDIUM AND HIGH LEVEL LCS – The LCS's will be prepared at the following concentrations and rotated per batch; 2 ng/L, 40 ng/L, 500 ng/l for drinking waters. The analyte PDS contains all the method analytes of interest at various concentrations in methanol. The analyte PDS has been shown to be stable for six months when stored at ≤4 °C.

- 8.2.7 Isotope Dilution Labeled Recovery Stock Solutions (ID REC) ID REC Stock solutions are stable for at least 6 months when stored at 4 °C. The stock solution is purchased at a concentration of 1000 ng/mL.
- 8.2.8 Isotope Dilution Labeled Recovery Primary Dilution Standard (ID REC PDS) Prepare the ID REC PDS at a concentration of 500 ng/mL. The ID REC PDS is prepared in 80:20% (vol/vol) methanol:water. The ID REC PDS is stable for at least six months when stored in polypropylene centrifuge tubes at ≤4 °C.

Table 5

Analyte	Conc. of REC Stock (ng/mL)	Vol. of REC Stock (mL)	Final Vol. of REC PDS (mL)	Final Conc. of REC PDS (ng/mL)
M2PFOA	2000	1	4	500
		1	7	
M2PFDA	2000	11	4	500
M3PFBA	2000	1	4	500
M4PFOS	2000	1	4	500

8.2.9 CALIBRATION STANDARDS (CAL) -

Current Concentrations (ng/mL): 0.5, 1.0, 5.0, 10.0, 50.0, 125, 150, 250, 500

Prepare the CAL standards over the concentration range of interest from dilutions of the analyte PDS in methanol containing 20% reagent water. 20 μ l of the EIS PDS and REC PDS are added to the CAL standards to give a constant concentration of 10 ng/ml. The lowest concentration CAL standard must be at or below the RL (2 ng/L), which may depend on system sensitivity. The CAL standards may also be used as CCVs (Sect. 9.8). To make calibration stock standards:

Table 6

Calibration Standard Concentration	Final Aqueous Cal STD Level Concentration	Final Soil Cal STD Level Concentration	24 compound stock added (ul)	PFHxDA Stock added (ul)	500 ng/ml PFHxDA dilution added (ul)	PFODA Stock added (ul)	500 ng/ml PFODA dilution added (ul)	ADONA, HFPO-DA, 11CI- PF3OUdS, 9CI- PF3ONS Stock added (ul)	500 ng/ml ADONA dilution added (ul)	Final Volume in MeOH/H ₂ O (82:20)
.5 ng/ml	2 ng/L	.25 ng/g	6.25		25		25		25	25 mls
1 ng/ml	4 ng/L	.5 ng/g	5		20		20		20	10 mls
5 ng/ml	20 ng/L	1 ng/g	25		100		100		100	10 mls
10 ng/ml	40 ng/L	5 ng/g	125	5		5		5		25 mls

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

ID No.:23528

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 12 of 28

50 ng/ml	200 ng/L	25 ng/g	250	10	10	10	10 mls
125 ng/ml	500 ng/L	62.5 ng/g	625	25	25	25	10 mls
150 ng/ml	600 ng/L	75 ng/g	750	30	30	30	10 mls
250 ng/ml	1000 ng/L	125 ng/g	625				5 mls
500 ng/ml	2000 ng/L	250 ng/g	1250			0	5 mls

9. Quality Control

The laboratory must maintain records to document the quality of data that is generated. Ongoing data quality checks are compared with established performance criteria to determine if the results of analyses meet the performance characteristics of the method.

9.1 MINIMUM REPORTING LIMIT (MRL) CONFIRMATION

9.1.1 Fortify, extract, and analyze seven replicate LCSs at 2 ng/l. Calculate the mean measured concentration (*Mean*) and standard deviation for these replicates. Determine the Half Range for the prediction interval of results (*HR_{PIR}*) using the equation below

$$HR_{PIR} = 3.963s$$

Where:

s = the standard deviation 3.963 = a constant value for seven replicates.

9.1.2 Confirm that the upper and lower limits for the Prediction Interval of Result ($PIR = Mean \pm HR_{PIR}$) meet the upper and lower recovery limits as shown below

The Upper PIR Limit must be ≤150% recovery.

Mean + HR
$$_{PIR}$$
 x 100% ≤ 150%
Fortified Concentration

The Lower PIR Limit must be $\geq 50\%$ recovery.

Mean − HR
$$_{PIR}$$
 x 100% ≥ 50%
Fortified Concentration

9.1.3 The RL is validated if both the Upper and Lower PIR Limits meet the criteria described above. If these criteria are not met, the RL has been set too low and must be determined again at a higher concentration.

9.2 Blank(s)

9.2.1 METHOD BLANK (MB) - A Method Blank (MB) is required with each extraction batch to confirm that potential background contaminants are not interfering with the identification or quantitation of method analytes. Prep and analyze a MB for every 20 samples. If the MB produces a peak within the retention time window of any analyte that would prevent the determination of that analyte, determine the source of contamination and eliminate the interference before processing samples. Background contamination must be reduced to an acceptable level before proceeding. Background from method analytes or other contaminants that

ID No.:23528 Revision 12 Published Date: 2/22/2019 3:48:15 PM

Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 13 of 28

> interfere with the measurement of method analytes must be below the RL. If the method analytes are detected in the MB at concentrations equal to or greater than this level, then all data for the problem analyte(s) must be considered invalid for all samples in the extraction batch. Because background contamination is a significant problem for several method analytes, it is highly recommended that the analyst maintain a historical record of MB data.

9.2.2 FIELD REAGENT BLANK (FRB) - The purpose of the FRB is to ensure that PFAS's measured in the Field Samples were not inadvertently introduced into the sample during sample collection/handling. Analysis of the FRB is required only if a Field Sample contains a method analyte or analytes at or above the RL. The FRB is processed, extracted and analyzed in exactly the same manner as a Field Sample.

9.3 Laboratory Control Sample (LCS) and Laboratory Control Sample Duplicates (LCSD)

An LCS is required with each extraction batch. The fortified concentration of the 9.3.1 LCS may be rotated between low, medium, and high concentrations from batch to batch. Default limits of 50-150% of the true value may be used for analytes until sufficient replicates have been analyzed to generate proper control limits. Calculate the percent recovery (%R) for each analyte using the equation

$$%R = A \times 100$$
B

Where:

A = measured concentration in the fortified sample B =fortification concentration.

9.3.2 Where applicable, LCSD's are to be extracted and analyzed. The concentration and analyte recovery criteria for the LCSD must be the same as the batch LCS The RSD's must fall within ≤30% of the true value for medium and high level replicates, and ≤50% for low level replicates. Calculate the relative percent difference (RPD) for duplicate MSs (MS and MSD) using the equation

$$RPD = |LCS - LCSD| \times 100$$

$$(LCS + LCSD)/2$$

9.3.3 If the LCS and or LCSD results do not meet these criteria for method analytes, then all data for the problem analyte(s) must be considered invalid for all samples in the extraction batch.

9.4 Labeled Recovery Standards (REC)

The analyst must monitor the peak areas of the REC(s) in all injections during each analysis day. **9.5** Extracted Internal Standards (EIS)

The EIS standard is fortified into all samples, CCVs, MBs, LCSs, MSs, MSDs, FD, and FRB prior to extraction. It is also added to the CAL standards. The EIS is a means of assessing method performance from extraction to final

Printouts of this document may be out of date and should be considered uncontrolled. To accomplish work, the published version of the document should be viewed online.

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528 Revision 12 Published Date: 2/22/2019 3:48:15 PM

Page 14 of 28

chromatographic measurement. Calculate the recovery (%R) for the EIS using the following equation

 $%R = (A/B) \times 100$

Where:

A = calculated EIS concentration for the QC or Field Sample

B = fortified concentration of the EIS.

9.5.2 Default limits of 50-150% may be used for analytes until sufficient replicates have been analyzed to generate proper control limits. A low or high percent recovery for a sample, blank, or CCV does not require discarding the analytical data but it may indicate a potential problem with future analytical data. When EIS recovery from a sample, blank, or CCV are outside control limits, check 1) calculations to locate possible errors, 2) standard solutions for degradation, 3) contamination, and 4) instrument performance. For CCVs and QC elements spiked with all target analytes, if the recovery of the corresponding target analytes meet the acceptance criteria for the EIS in question, the data can be used but all potential biases in the recovery of the EIS must be documented in the sample report. If the associated target analytes do not meet the acceptance criteria, the data must be reanalyzed.

9.6 Matrix Spike (MS)

- 9.6.1 Analysis of an MS is required in each extraction batch and is used to determine that the sample matrix does not adversely affect method accuracy. Assessment of method precision is accomplished by analysis of a Field Duplicate (FD) (Sect. 9.6); however, infrequent occurrence of method analytes would hinder this assessment. If the occurrence of method analytes in the samples is infrequent, or if historical trends are unavailable, a second MS, or MSD, must be prepared, extracted, and analyzed from a duplicate of the Field Sample. Extraction batches that contain MSDs will not require the extraction of a field sample duplicate. If a variety of different sample matrices are analyzed regularly, for example, drinking water from groundwater and surface water sources, method performance should be established for each. Over time, MS data should be documented by the laboratory for all routine sample sources.
- 9.6.2 Within each extraction batch, a minimum of one Field Sample is fortified as an MS for every 20 Field Samples analyzed. The MS is prepared by spiking a sample with an appropriate amount of the Analyte Stock Standard (Sect. 8.2.3). Use historical data and rotate through the low, mid and high concentrations when selecting a fortifying concentration. Calculate the percent recovery (%R) for each analyte using the equation

$$%R = (A - B) \times 100$$

C

Where:

A = measured concentration in the fortified sample

B = measured concentration in the unfortified sample

C = fortification concentration.

9.6.3 Analyte recoveries may exhibit matrix bias. For samples fortified at or above their native concentration, recoveries should range between 50-150%. If the accuracy of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the LCS, the recovery is judged to be

Document Type: SOP-Technical

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Revision 12 Published Date: 2/22/2019 3:48:15 PM Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 15 of 28

> matrix biased. The result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.

ID No.:23528

9.7 Laboratory Duplicate

- 9.7.1 FIELD DUPLICATE OR LABORATORY FORTIFIED SAMPLE MATRIX DUPLICATE (FD or MSD) - Within each extraction batch (not to exceed 20 Field Samples), a minimum of one FD or MSD must be analyzed. Duplicates check the precision associated with sample collection, preservation, storage, and laboratory procedures. If method analytes are not routinely observed in Field Samples, an MSD should be analyzed rather than an FD.
- 9.7.2 Calculate the relative percent difference (RPD) for duplicate measurements (FD1 and FD2) using the equation

RPD =
$$|FD1 - FD2|$$
 x 100
(FD1 + FD2) / 2

- 9.7.3 RPDs for FDs should be ≤30%. Greater variability may be observed when FDs have analyte concentrations that are within a factor of 2 of the RL. At these concentrations, FDs should have RPDs that are ≤50%. If the RPD of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the CCV, the recovery is judged to be matrix biased. The result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.
- 9.7.4 If an MSD is analyzed instead of a FD, calculate the relative percent difference (RPD) for duplicate MSs (MS and MSD) using the equation

$$RPD = \underline{|MS - MSD|} x 100$$

$$(MS + MSD) / 2$$

9.7.5 RPDs for duplicate MSs should be ≤30% for samples fortified at or above their native concentration. Greater variability may be observed when MSs are fortified at analyte concentrations that are within a factor of 2 of the RL. MSs fortified at these concentrations should have RPDs that are ≤50% for samples fortified at or above their native concentration. If the RPD of any analyte falls outside the designated range, and the laboratory performance for that analyte is shown to be in control in the LCSD where applicable, the result is judged to be matrix biased. If no LCSD is present, the associated MS and MSD are to be re-analyzed to determine if any analytical has occurred. If the resulting RPDs are still outside control limits, the result for that analyte in the unfortified sample is labeled suspect/matrix to inform the data user that the results are suspect due to matrix effects.

9.8 Initial Calibration Verification (ICV)

As part of the IDC (Sect. 13.2), and after each ICAL, analyze a QCS sample 9.8.1 from a source different from the source of the CAL standards. If a second vendor is not available, then a different lot of the standard should be used. The QCS should be prepared and analyzed just like a CCV. Acceptance criteria for the QCS are identical to the CCVs; the calculated amount for each analyte must be ±

Revision 12 Published Date: 2/22/2019 3:48:15 PM Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 16 of 28


> 30% of the expected value. If measured analyte concentrations are not of acceptable accuracy, check the entire analytical procedureto locate and correct the problem.

9.9 Continuing Calibration Verification (CCV)

CCV Standards are analyzed at the beginning of each analysis batch, after every 10 Field Samples, and at the end of the analysis batch. See Section 10.7 for concentration requirements and acceptance criteria.

9.10 Method-specific Quality Control Samples

9.10.1 PEAK ASYMMETRY FACTOR - A peak asymmetry factor must be calculated using the equation below during the IDL and every time a calibration curve is generated. The peak asymmetry factor for the first two eluting peaks in a midlevel CAL standard (if only two analytes are being analyzed, both must be evaluated) must fall in the range of 0.8 to 1.5. Modifying the standard or extract composition to more aqueous content to prevent poor shape is not permitted. See guidance in Section 10.6.4.1 if the calculated peak asymmetry factors do not meet the criteria.



ID No.:23528

 $A_c = b/a$

Where:

 A_s = peak asymmetry factor

- b = width of the back half of the peak measured (at 10% peak height) from the trailing edge of the peak to a line dropped perpendicularly from the peak apex
- a = the width of the front half of the peak measured (at 10% peak height) from the leading edge of the peak to a line dropped perpendicularly from the apex.

9.11 Method Sequence

- **CCV-LOW**
- MB
- LCS
- **LCSD**
- MS
- Duplicate or MSD
- Field Samples (1-10)
- **CCV-MID**
- Field Samples (11-20)
- **CCV-LOW**

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Department: Semivolatiles Published Date: 2/22/2019 3:48:15 PM
Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 17 of 28

10. Procedure

10.1 Equipment Set-up

10.1.1 This procedure may be performed manually or in an automated mode using a robotic or automatic sample preparation device. If an automated system is used to prepare samples, follow the manufacturer's operating instructions, but all extraction and elution steps must be the same as in the manual procedure. Extraction and/or elution steps may not be changed or omitted to accommodate the use of an automated system. If an automated system is used, the MBs should be rotated among the ports to ensure that all the valves and tubing meet the MB requirements (Sect. 9.2).

ID No.:23528

Revision 12

- 10.1.2 Some of the PFAS's adsorb to surfaces, including polypropylene. Therefore, the aqueous sample bottles must be rinsed with the elution solvent (Sect 10.3.4) whether extractions are performed manually or by automation. The bottle rinse is passed through the cartridge to elute the method analytes and is then collected (Sect. 10.3.4).
- **10.1.3 NOTE:** The SPE cartridges and sample bottles described in this section are designed as single use items and should be discarded after use. They may not be refurbished for reuse in subsequent analyses.

10.2 Sample Preparation and Extraction of Aqueous Samples

10.2.1 Samples are preserved, collected and stored as presented in Section 6.

The entire sample that is received must be sent through the SPE cartridge. In addition, the bottle must be solvent rinsed and this rinse must be sent through the SPE cartridge as well. The method blank (MB) and laboratory control sample (LCS) must be extracted in exactly the same manner (i.e., must include the bottle solvent rinse). It should be noted that a water rinse alone is not sufficient. This does not apply to samples with high concentrations of PFAS that are prepared using serial dilution and not SPE.

- 10.2.2 Determine sample volume. Weigh all samples to the nearest 1g. If visible sediment is present, centrifuge and decant into a new 250mL HDPE bottle and record the weight of the new container.
 - NOTE: Some of the PFAS's adsorb to surfaces, thus the sample volume may **NOT** be transferred to a graduated cylinder for volume measurement.
- 10.2.3 The MB, LCS and FRB may be prepared by measuring 250 mL of reagent water with a polypropylene graduated cylinder or filling a 250-mL sample bottle to near the top.
- 10.2.4 Adjust the QC and sample pH to 3 by adding acetic acid in water dropwise
- 10.2.5 Add 20 μL of the EIS PDS (Sect. 8.2.2) to each sample and QC, cap and invert to mix.
- 10.2.6 If the sample is an LCS, LCSD, MS, or MSD, add the necessary amount of analyte PDS (Sect. 8.2.3). Cap and invert each sample to mix.

10.3 Cartridge SPE Procedure

Alpha Analytical, Inc.
Facility: Mansfield, MA
Department: Semivolatiles
Title: PFAS by SPE and LC/MS/MS Isotope Dilution

Revision 12 Published Date: 2/22/2019 3:48:15 PM Page 18 of 28

ID No.:23528

10.3.1 CARTRIDGE CLEAN-UP AND CONDITIONING – DO NOT allow cartridge packing material to go dry during any of the conditioning steps. Rinse each cartridge with 3 X 5 mL of 2% ammonium hydroxide in methanol, followed by 5mls of methanol. Next, rinse each cartridge with 5 mls of the 25 mM acetate buffer, followed by 15 mL of reagent water, without allowing the water to drop below the top edge of the packing. If the cartridge goes dry during the conditioning phase, the conditioning must be started over. Add 4-5 mL of reagent water to each cartridge, attach the sample transfer tubes (Sect. 7.9.3), turn on the vacuum, and begin adding sample to the cartridge.

- 10.3.2 SAMPLE EXTRACTON Adjust the vacuum so that the approximate flow rate is approximately 4 mL/min. Do not allow the cartridge to go dry before all the sample has passed through.
- 10.3.3 SAMPLE BOTTLE AND CARTRIDGE RINSE After the entire sample has passed through the cartridge, rinse the sample bottles with 4 ml reagent water followed by 4 ml 25 mM acetate buffer at pH 4 and draw the aliquot through the sample transfer tubes and the cartridges. Draw air or nitrogen through the cartridge for 5-10 min at high vacuum (10-15 in. Hg). NOTE: If empty plastic reservoirs are used in place of the sample transfer tubes to pass the samples through the cartridges, these reservoirs must be treated like the transfer tubes. After the entire sample has passed through the cartridge, the reservoirs must be rinsed to waste with reagent water.
- 10.3.4 SAMPLE BOTTLE AND CARTRIDGE ELUTION, Fraction 1 Turn off and release the vacuum. Lift the extraction manifold top and insert a rack with collection tubes into the extraction tank to collect the extracts as they are eluted from the cartridges. Rinse the sample bottles with 12 mls of methanol and draw the aliquot through the sample transfer tubes and cartridges. Use a low vacuum such that the solvent exits the cartridge in a dropwise fashion.

SAMPLE BOTTLE AND CARTRIDGE ELUTION, Fraction 2 In a separate collection vial, rinse the sample bottles with 12 mL of 2% ammonium hydroxide in methanol and elute the analytes from the cartridges by pulling the 4 mL of methanol through the sample transfer tubes and the cartridges. Use a low vacuum such that the solvent exits the cartridge in a dropwise fashion. To the final extract, add 50 ul of acetic acid.

NOTE: If empty plastic reservoirs are used in place of the sample transfer tubes to pass the samples through the cartridges, these reservoirs must be treated like the transfer tubes. After the reservoirs have been rinsed in Section 10.3.3, the elution solvent used to rinse the sample bottles must be swirled down the sides of the reservoirs while eluting the cartridge to ensure that any method analytes on the surface of the reservoirs are transferred to the extract.

CLEAN-UP CARTRIDGE ELUTION, Elute the clean-up cartridge with 8 additional mls of methanol and draw the aliquot through the cartridge. Use a low vacuum such that the solvent exits the cartridge in a dropwise fashion.

10.3.5 Fractions 1 and 2 are to be combined during the concentration stage (section 10.6)

10.4 Sample Prep and Extraction Protocol for Soils

Alpha Analytical, Inc.
Facility: Mansfield, MA
Department: Semivolatiles
Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:**23528**Revision 12
Published Date: 2/22/2019 3:48:15 PM
tion Page 19 of 28

- 10.4.1 Homogenize and weigh 2 grams of sample (measured to the nearest hundredth of a gram) into a 50 ml polypropylene centrifuge tube. For laboratory control blanks and spikes, 2 grams of clean sand is used.
- 10.4.2 Add 20 µL of the EIS PDS (Sect. 8.2.2) to each sample and QC.
- **10.4.3** If the sample is an LCS, LCSD, MS, or MSD, add the necessary amount of analyte PDS (Sect. 8.2.3). Cap and invert each sample to mix.
- **10.4.4** To all samples, add 10 mls of methanol, cap, vortex for 25 seconds at 3000RPM and mix for 30 minutes using a shaker table of tumbler at 120RPM.
- **10.4.5** Following mixing, sonicate each sample for 30 minutes and let samples sit overnight (at least 2 hours is required for RUSH samples).
- 10.4.6 Centrifuge each sample at 3500RPM for 10 minutes.
- 10.4.7 Remove supernatant, and reserve for clean-up.

10.5 Extract Clean-up

- 10.5.1 CARTRIDGE CLEAN-UP AND CONDITIONING –. Rinse each cartridge with 15 mL of methanol and discard. If the cartridge goes dry during the conditioning phase, the conditioning must be started over. Attach the sample transfer tubes (Sect. 7.9.3), turn on the vacuum, and begin adding sample to the cartridge.
- **10.5.2** Adjust the vacuum so that the approximate flow rate is 1-2 mL/min. Do not allow the cartridge to go dry before all the sample has passed through.
- 10.5.3 SAMPLE BOTTLE AND CARTRIDGE RINSE After the entire sample has passed through the cartridge, rinse the sample collection vial with two 1-mL aliquots of methanol and draw each aliquot through the cartridges. Draw air or nitrogen through the cartridge for 5 min at high vacuum (10-15 in. Hg).
- **10.5.4** If extracts are not to be immediately evaporated, cover collection tubes and store at ambient temperature till concentration.

10.6 Extract Concentration

10.6.1 Concentrate the extract to dryness under a gentle stream of nitrogen in a heated water bath (60-65 °C) to remove all the water/methanol mix. Add the appropriate amount of 80:20% (vol/vol) methanol:water solution and 20 μl of the ID REC PDS (Sect. 8.2.7) to the collection vial to bring the volume to 1 mL and vortex. Transfer two aliquots with a plastic pipet (Sect. 7.6) into 2 polypropylene autosampler vials.

NOTE: It is recommended that the entire 1-mL aliquot not be transferred to the autosampler vial because the polypropylene autosampler caps do not reseal after injection. Therefore, do not store the extracts in the autosampler vials as evaporation losses can occur occasionally in these autosampler vials. Extracts can be split between 2 X 700 μ l vials (Sect. 7.4).

10.7 Sample Volume Determination

Alpha Analytical, Inc.

Facility: Mansfield, MA

Department: Semivolatiles

Revision 12

Published Date: 2/22/2019 3:48:15 PM

Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 20 of 28

10.7.1 If the level of the sample was marked on the sample bottle, use a graduated cylinder to measure the volume of water required to fill the original sample bottle to the mark made prior to extraction. Determine to the nearest 10 mL.

- 10.7.2 If using weight to determine volume, weigh the empty bottle to the nearest 10 g and determine the sample weight by subtraction of the empty bottle weight from the original sample weight (Sect. 10.2.2). Assume a sample density of 1.0 g/mL. In either case, the sample volume will be used in the final calculations of the analyte concentration (Sect. 11.2).
- 10.8 Initial Calibration Demonstration and documentation of acceptable initial calibration is required before any samples are analyzed. After the initial calibration is successful, a CCV is required at the beginning and end of each period in which analyses are performed, and after every tenth Field Sample.

10.8.1 ESI-MS/MS TUNE

- **10.8.1.1** Calibrate the mass scale of the MS with the calibration compounds and procedures prescribed by the manufacturer.
- 10.8.1.2 Optimize the [M-H]- for each method analyte by infusing approximately 0.5-1.0 µg/mL of each analyte (prepared in the initial mobile phase conditions) directly into the MS at the chosen LC mobile phase flow rate (approximately 0.4 mL/min). This tune can be done on a mix of the method analytes. The MS parameters (voltages, temperatures, gas flows, etc.) are varied until optimal analyte responses are determined. The method analytes may have different optima requiring some compromise between the optima.
- 10.8.1.3 Optimize the product ion for each analyte by infusing approximately 0.5-1.0 μg/mL of each analyte (prepared in the initial mobile phase conditions) directly into the MS at the chosen LC mobile phase flow rate (approximately 0.4 mL/min). This tune can be done on a mix of the method analytes. The MS/MS parameters (collision gas pressure, collision energy, etc.) are varied until optimal analyte responses are determined. Typically, the carboxylic acids have very similar MS/MS conditions and the sulfonic acids have similar MS/MS conditions.
- 10.8.2 Establish LC operating parameters that optimize resolution and peak shape. Modifying the standard or extract composition to more aqueous content to prevent poor shape is not permitted.

Cautions: LC system components, as well as the mobile phase constituents, contain many of the method analytes in this method. Thus, these PFAS's will build up on the head of the LC column during mobile phase equilibration. To minimize the background PFAS peaks and to keep background levels constant, the time the LC column sits at initial conditions must be kept constant and as short as possible (while ensuring reproducible retention times). In addition, prior to daily use, flush the column with 100% methanol for at least 20 min before initiating a sequence. It may be necessary on some systems to flush other LC components such as wash syringes, sample needles or any other system components before daily use.

10.8.3 Inject a mid-level CAL standard under LC/MS conditions to obtain the retention times of each method analyte. If analyzing for PFTA, ensure that the LC

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:**23528** Revision 12 Published Date: 2/22/2019 3:48:15 PM

Page 21 of 28

conditions are adequate to prevent co-elution of PFTA and the mobile phase interferants. These interferants have the same precursor and products ions as PFTA, and under faster LC conditions may co-elute with PFTA. Divide the chromatogram into retention time windows each of which contains one or more chromatographic peaks. During MS/MS analysis, fragment a small number of selected precursor ions ([M-H]-) for the analytes in each window and choose the most abundant product ion. For maximum sensitivity, small mass windows of ±0.5 daltons around the product ion mass were used for quantitation.

- 10.8.4 Inject a mid-level CAL standard under optimized LC/MS/MS conditions to ensure that each method analyte is observed in its MS/MS window and that there are at least 10 scans across the peak for optimum precision.
 - 10.8.4.1 If broad, split or fronting peaks are observed for the first two eluting chromatographic peaks (if only two analytes are being analyzed, both must be evaluated), change the initial mobile phase conditions to higher aqueous content until the peak asymmetry ratio for each peak is 0.8 1.5. The peak asymmetry factor is calculated as described in Section 9.9.1 on a mid-level CAL standard. The peak asymmetry factor must meet the above criteria for the first two eluting peaks during the IDL and every time a new calibration curve is generated. Modifying the standard or extract composition to more aqueous content to prevent poor shape is not permitted.

NOTE: PFHxS, PFOS, NMeFOSAA, and NEtFOSAA have multiple chromatographic peaks using the LC conditions in Table 5 due to chromatographic resolution of the linear and branched isomers of these compounds. Most PFAS's are produced by two different processes. One process gives rise to linear PFAS's only while the other process produces both linear and branched isomers. Thus, both branched and linear PFAS's can potentially be found in the environment. For the aforementioned compounds that give rise to more than one peak, all the chromatographic peaks observed in the standard must be integrated and the areas totaled. Chromatographic peaks in a sample must be integrated in the same way as the CAL standard.

- 10.8.5 Prepare a set of CAL standards as described in Section 8.2.5. The lowest concentration CAL standard must be at or below the RL (2 ng/L), which may depend on system sensitivity.
- 10.8.6 The LC/MS/MS system is calibrated using the IS technique. Use the LC/MS/MS data system software to generate a linear regression or quadratic calibration curve for each of the analytes. This curve must always be forced through zero and may be concentration weighted, if necessary. Forcing zero allows for a better estimate of the background levels of method analytes. A minimum of 5 levels are required for a linear calibration model and a minimum of 6 levels are required for a quadratic calibration model.
- 10.8.7 CALIBRATION ACCEPTANCE CRITERIA A linear fit is acceptable if the coefficient of determination (r²) is greater than 0.99. When quantitated using the initial calibration curve, each calibration point, except the lowest point, for each analyte should calculate to be within 70-130% of its true value. The lowest CAL point should calculate to be within 50-150% of its true value. If these criteria cannot be met, the analyst will have difficulty meeting ongoing QC criteria. It is

Alpha Analytical, Inc.
Facility: Mansfield, MA
Department: Semiyolatiles

Department: Semivolatiles Published Date: 2/22/2019 3:48:15 PM
Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 22 of 28

recommended that corrective action is taken to reanalyze the CAL standards, restrict the range of calibration, or select an alternate method of calibration (forcing the curve through zero is still required).

ID No.:23528

Revision 12

- 10.8.7.1 CAUTION: When acquiring MS/MS data, LC operating conditions must be carefully reproduced for each analysis to provide reproducible retention times. If this is not done, the correct ions will not be monitored at the appropriate times. As a precautionary measure, the chromatographic peaks in each window must not elute too close to the edge of the segment time window.
- 10.9 CONTINUING CALIBRATION CHECK (CCV) Minimum daily calibration verification is as follows. Verify the initial calibration at the beginning and end of each group of analyses, and after every tenth sample during analyses. In this context, a "sample" is considered to be a Field Sample. MBs, CCVs, LCSs, MSs, FDs FRBs and MSDs are not counted as samples. The beginning CCV of each analysis batch must be at or below the RL in order to verify instrument sensitivity prior to any analyses. If standards have been prepared such that all low CAL points are not in the same CAL solution, it may be necessary to analyze two CAL standards to meet this requirement. Alternatively, the analyte concentrations in the analyte PDS may be customized to meet these criteria. Subsequent CCVs should alternate between a medium and Low concentration CAL standard.
 - **10.9.1** Inject an aliquot of the appropriate concentration CAL standard and analyze with the same conditions used during the initial calibration.
 - 10.9.2 Calculate the concentration of each analyte and EIS in the CCV. The calculated amount for each analyte for medium level CCVs must be within ± 30% of the true value with an allowance of 10% of the reported analytes to be greater than 30%, but less than 40%. The calculated amount for each EIS must be within ± 50% of the true value. The calculated amount for the lowest calibration point for each analyte must be within ± 50%. If these conditions do not exist, then all data for the problem analyte must be considered invalid, and remedial action should be taken (Sect. 10.7.4) which may require recalibration. Any Field or QC Samples that have been analyzed since the last acceptable calibration verification should be reanalyzed after adequate calibration has been restored, with the following exception. If the CCV fails because the calculated concentration is greater than 130% (150% for the low-level CCV) for a particular method analyte, and Field Sample extracts show no detection for that method analyte, non-detects may be reported without re-analysis.
 - 10.9.3 REMEDIAL ACTION Failure to meet CCV QC performance criteria may require remedial action. Major maintenance, such as cleaning the electrospray probe, atmospheric pressure ionization source, cleaning the mass analyzer, replacing the LC column, etc., requires recalibration (Sect 10.6) and verification of sensitivity by analyzing a CCV at or below the RL (Sect 10.7).

Alpha Analytical, Inc.
Facility: Mansfield, MA
Department: Semivolatiles
Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528 Revision 12 Published Date: 2/22/2019 3:48:15 PM

Page 23 of 28

10.10 EXTRACT ANALYSIS

- 10.10.1 Establish operating conditions equivalent to those summarized in Tables 6-8 of Section 16. Instrument conditions and columns should be optimized prior to the initiation of the IDC.
- 10.10.2 Establish an appropriate retention time window for each analyte. This should be based on measurements of actual retention time variation for each method analyte in CAL standard solutions analyzed on the LC over the course of time. A value of plus or minus three times the standard deviation of the retention time obtained for each method analyte while establishing the initial calibration and completing the IDC can be used to calculate a suggested window size. However, the experience of the analyst should weigh heavily on the determination of the appropriate retention window size.
- 10.10.3 Calibrate the system by either the analysis of a calibration curve (Sect. 10.6) or by confirming the initial calibration is still valid by analyzing a CCV as described in Section 10.7. If establishing an initial calibration, complete the IDC as described in Section 13.2.
- 10.10.4 Begin analyzing Field Samples, including QC samples, at their appropriate frequency by injecting the same size aliquots under the same conditions used to analyze the CAL standards.
- 10.10.5 At the conclusion of data acquisition, use the same software that was used in the calibration procedure to identify peaks of interest in predetermined retention time windows. Use the data system software to examine the ion abundances of the peaks in the chromatogram. Identify an analyte by comparison of its retention time with that of the corresponding method analyte peak in a reference standard.
- 10.10.6 The analyst must not extrapolate beyond the established calibration range. If an analyte peak area exceeds the range of the initial calibration curve, the sample should be re-extracted with a reduced sample volume in order to bring the out of range target analytes into the calibration range. If a smaller sample size would not be representative of the entire sample, the following options are recommended. Re-extract an additional aliquot of sufficient size to insure that it is representative of the entire sample. Spike it with a higher concentration of internal standard. Prior to LC/MS analysis, dilute the sample so that it has a concentration of internal standard equivalent to that present in the calibration standard. Then, analyze the diluted extract.

11. Data Evaluation, Calculations and Reporting

- 11.1 Complete chromatographic resolution is not necessary for accurate and precise measurements of analyte concentrations using MS/MS. In validating this method, concentrations were calculated by measuring the product ions listed in Table 7.
- 11.2 Calculate analyte concentrations using the multipoint calibration established in Section 10.6. Do not use daily calibration verification data to quantitate analytes in samples. Adjust final analyte concentrations to reflect the actual sample volume determined in Section 10.6 where:

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles

Revision 12 Published Date: 2/22/2019 3:48:15 PM Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 24 of 28

Cex = (Area of target analyte * Concentration of Labeled analog) / (area of labeled analog *

ID No.:23528

 $C_s = (C_{ex} / sample volume in ml) * 1000$

 C_{ex} = The concentration of the analyte in the extract

CF = calibration factor from calibration.

- 11.3 Prior to reporting the data, the chromatogram should be reviewed for any incorrect peak identification or poor integration.
- 11.4 PFHxS, PFOS, PFOA, NMeFOSAA, and NEtFOSAA have multiple chromatographic peaks using the LC conditions in Table 5 due to the linear and branch isomers of these compounds (Sect. 10.6.4.1). The areas of all the linear and branched isomer peaks observed in the CAL standards for each of these analytes must be summed and the concentrations reported as a total for each of these analytes.
- 11.5 Calculations must utilize all available digits of precision, but final reported concentrations should be rounded to an appropriate number of significant figures (one digit of uncertainty), typically two, and not more than three significant figures.

12. Contingencies for Handling Out-of-Control Data or Unacceptable Data

- 12.1 Section 9.0 outlines sample batch QC acceptance criteria. If non-compliant organic compound results are to be reported, the Organic Section Head and/or the Laboratory Director, and the Operations Manager must approve the reporting of these results. The laboratory Project Manager shall be notified, and may choose to relay the non-compliance to the client, for approval, or other corrective action, such as re-sampling and re-analysis. The analyst, Data Reviewer, or Department Supervisor performing the secondary review initiates the project narrative, and the narrative must clearly document the non-compliance and provide a reason for acceptance of these results.
- 12.2 All results for the organic compounds of interest are reportable without qualification if extraction and analytical holding times are met, preservation requirements (including cooler temperatures) are met, all QC criteria are met, and matrix interference is not suspected during extraction or analysis of the samples. If any of the below QC parameters are not met, all associated samples must be evaluated for re-extraction and/or re-analysis.

13. Method Performance

13.1 Detection Limit Study (DL) / Limit of Detection Study (LOD) / Limit of Quantitation (LOQ)

13.1.1 The laboratory follows the procedure to determine the DL, LOD, and/or LOQ as outlined in Alpha SOP ID 1732. These studies performed by the laboratory are maintained on file for review.

Alpha Analytical, Inc. Facility: Mansfield, MA

Revision 12 Department: Semivolatiles Published Date: 2/22/2019 3:48:15 PM Title: PFAS by SPE and LC/MS/MS Isotope Dilution Page 25 of 28

13.2 Demonstration of Capability Studies

13.2.1 The IDC must be successfully performed prior to analyzing any Field Samples. Prior to conducting the IDC, the analyst must first generate an acceptable Initial Calibration following the procedure outlined in Section 10.6.

ID No.:23528

- 13.2.2 INITIAL DEMONSTRATION OF LOW SYSTEM BACKGROUND Any time a new lot of SPE cartridges, solvents, centrifuge tubes, disposable pipets, and autosampler vials are used, it must be demonstrated that an MB is reasonably free of contamination and that the criteria in Section 9.2.1 are met. If an automated extraction system is used, an MB should be extracted on each port to ensure that all the valves and tubing are free from potential PFAS contamination.
- 13.2.3 INITIAL DEMONSTRATION OF PRECISION (IDP) Prepare, extract, and analyze four to seven replicate LCSs fortified near the midrange of the initial calibration curve according to the procedure described in Section 10. Sample preservatives as described in Section 6.2.1 must be added to these samples. The relative standard deviation (RSD) of the results of the replicate analyses must be less than 20%.
- 13.2.4 INITIAL DEMONSTRATION OF ACCURACY (IDA) Using the same set of replicate data generated for Section 13.2.3, calculate average recovery. The average recovery of the replicate values must be within ± 30% of the true value.
- 13.2.5 INITIAL DEMONSTRATION OF PEAK ASYMMETRY FACTOR Peak asymmetry factors must be calculated using the equation in Section 9.10.1 for the first two eluting peaks (if only two analytes are being analyzed, both must be evaluated) in a mid-level CAL standard. The peak asymmetry factors must fall in the range of 0.8 to 1.5. See guidance in Section 10.6.4.1 if the calculated peak asymmetry factors do not meet the criteria.
- 13.2.6 Refer to Alpha SOP ID 1739 for further information regarding IDC/DOC Generation.
- 13.2.7 The analyst must make a continuing, annual, demonstration of the ability to generate acceptable accuracy and precision with this method.

14. Pollution Prevention and Waste Management

- 14.1 Refer to Alpha's Chemical Hygiene Plan and Hazardous Waste Management and Disposal SOP for further pollution prevention and waste management information.
- 14.2 This method utilizes SPE to extract analytes from water. It requires the use of very small volumes of organic solvent and very small quantities of pure analytes, thereby minimizing the potential hazards to both the analyst and the environment as compared to the use of large volumes of organic solvents in conventional liquid-liquid extractions.
- 14.3 The analytical procedures described in this method generate relatively small amounts of waste since only small amounts of reagents and solvents are used. The matrices of concern are finished drinking water or source water. However, laboratory waste management practices must be conducted consistent with all applicable rules and regulations, and that laboratories protect the air, water, and land by minimizing and controlling all releases from fume hoods and bench operations. Also, compliance is required with any sewage discharge permits and regulations, particularly the hazardous waste identification rules and land disposal restrictions.

Alpha Analytical, Inc. Facility: Mansfield, MA Department: Semivolatiles Published Date: 2/22/2019 3:48:15 PM Title: PFAS by SPE and LC/MS/MS Isotope Dilution

15. Referenced Documents

Chemical Hygiene Plan - ID 2124

SOP ID 1732 Detection Limit (DL), Limit of Detection (LOD) & Limit of Quantitation (LOQ) SOP

ID No.:23528

Revision 12

Page 26 of 28

SOP ID 1739 Demonstration of Capability (DOC) Generation SOP

SOP ID 1728 Hazardous Waste Management and Disposal SOP

16. Attachments

Table 7: LC Method Conditions

Time (min)	2 mM Ammonium Acetate (5:95 MeOH/H ₂ O)	100% Methanol
Initial	100.0	0.0
1.0	100.0	0.0
2.2	85.0	15.0
11	20.0	80.0
11.4	0.0	100.0
12.4	100.0	0.00
15.5	100.0	0.0
Motoro Aguita III	DI C @ DEUC 21 v E0 mm packed wit	h 17 um DEU C

Waters Aquity UPLC ® BEHC₁₈ 2.1 x 50 mm packed with 1.7 µm BEH C₁₈ stationary phase Flow rate of 0.4 mL/min 2-5 µL injection

Table 8: ESI-MS Method Conditions

ESI Conditions					
Polarity	Negative ion				
Capillary needle voltage	.5 kV				
Cone Gas Flow	25 L/hr				
Nitrogen desolvation gas	1000 L/hr				
Desolvation gas temp.	500 °C				

Table 9: Method Analyte Source, Retention Times (RTs), and EIS References

#	Analyte	Transition	RT	IS	Туре
1	МЗРВА	216>171	2.65		REC
2	PFBA	213 > 169	2.65	2: M4PFBA	
3	M4PFBA	217 > 172	2.65	1: M3PBA	EIS
4	PFPeA	263 > 219	5.67	4: M5PFPEA	
5	M5PFPEA	268 > 223	5.66	1: M3PBA	EIS
6	PFBS	299 > 80	6.35	6: M3PFBS	
7	M3PFBS	302 > 80	6.35	29:M4PFOS	EIS
8	FtS 4:2	327 > 307	7.47	9: M2-4:2FTS	

Alpha Analytical, Inc.
Facility: Mansfield, MA
Department: Semivolatiles
Title: PFAS by SPE and LC/MS/MS Isotope Dilution

ID No.:23528 Revision 12 Published Date: 2/22/2019 3:48:15 PM Page 27 of 28

#	Analyte	Transition	RT	IS	Туре
9	M2-4:2FTS	329 > 81	7.47	29:M4PFOS	EIS
10	PFHxA	303 > 269	7.57	10: M5PFHxA	
11	M5PFHxA	318 > 273	7.57	19:M2PFOA	EIS
12	PFPeS	349 > 80	7.88	18: M3PFHxS	
13	PFHpA	363 > 319	8.80	14: M4PFHpA	
14	M4PFHpA	367 > 322	8.80	19:M2PFOA	EIS
15	L-PFHxS	399 > 80	8.94	18: M3PFHxS	
16	br-PFHxS	399 > 80	8.72	18: M3PFHxS	
17	PFHxS Total	399 > 80	8.94	18: M3PFHxS	
18	M3PFHxS	402 > 80	8.94	29:M4PFOS	EIS
19	MPFOA	415 > 370	9.7		REC
20	PFOA	413 > 369	9.7	23: M8PFOA	
21	br-PFOA	413 > 369	9.48	23: M8PFOA	
22	PFOA Total	413 > 369	9.7	23: M8PFOA	
23	M8PFOA	421 > 376	9.7	19: M2PFOA	EIS
24	FtS 6:2	427 > 407	9.66	25: M2-6:2FTS	
25	M2-6:2FTS	429 > 409	9.66	29:M4PFOS	EIS
26	PFHpS	449 > 80	9.78	33: M8PFOS	
27	PFNA	463 > 419	10.41	33: M8PFOS	
28	M9PFNA	472 > 427	10.41	19: M2PFOA	EIS
29	M4PFOS	501 > 80	10.45		REC
30	PFOS	499 > 80	10.45	33: M8PFOS	
31	br-PFOS	499 > 80	10.27	33: M8PFOS	
32	PFOS Total	499 > 80	10.45	33: M8PFOS	
33	M8PFOS	507 > 80	10.45	29: M4PFOS	EIS
34	FtS 8:2	527 > 507	10.99	38: M2-8:2FTS	
35	M2-8:2FTS	529 > 509	10.99	29:M4PFOS	EIS
36	M2PFDA	515 > 470	11.00		REC
37	PFDA	513 > 469	11.00	38: M6PFDA	
38	M6PFDA	519 > 474	11.00	36: M2PFDA	EIS
39	PFNS	549 > 80	11.02	33:M8PFOS	
40	NMeFOSAA	570 > 419	11.41	41: D3-NMeFOSAA	
41	d3-NMeFOSAA	573 > 419	11.41	36: M2PFDA	EIS
42	PFOSA	498 > 78	11.48	29: M8FOSA	
43	M8FOSA	506 > 78	11.48	19: M2PFOA	EIS
44	PFUnDA	563 > 519	11.51	41: M7-PFUDA	
45	M7-PFUDA	570 > 525	11.51	36: M2PFDA	EIS
46	PFDS	599 > 80	11.51	33:M8PFOS	
47	NEtFOSAA	584 > 419	11.68	48: d5-NEtFOSAA	

Facility: Mansfield, MA
Department: Semivolatiles
Title: PFAS by SPE and LC/MS/MS Isotope Dilution

Alpha Analytical, Inc.

ID No.:23528 Revision 12 Published Date: 2/22/2019 3:48:15 PM Page 28 of 28

#	Analyte	Transition	RT	IS	Туре
48	d5-NEtFOSAA	589 > 419	11.68	36: M2PFDA	EIS
49	PFDoA	613 > 569	11.96	50: MPFDOA	
50	MPFDOA	615 > 570	11.96	36: M2PFDA	EIS
51	PFTriA	663 > 619	12.34	50: MPFDOA	
52	PFTeA	713 > 669	12.6	53: M2PFTEDA	
53	M2PFTEDA	715 > 670	12.6	36: M2PFDA	EIS
54	M3HFPO-DA	329>285	7.97	19: M2PFOA	EIS
55	HFPO-DA	332>287	7.97	54: M3HFPO-DA	
56	ADONA	377>251	8.00	23: M8PFOA	
57	PFHxDA	813>769	13.20	59: M2PFHxDA	
58	PFODA	913>869	13.50	59: M2PFHxDA	
59	M2PFHxDA	815>770	13.20	36:M2PFDA	EIS
60	NEtFOSA	526>169	11.00	61: NMeFOSA	
61	NMeFOSA	512>169	10.50	63: d3-NMeFOSA	
62	d3-NMeFOSA	515>169	10.50	29: M4PFOS	EIS
63	d5-NEtFOSA	531>169	11.00	29: M4PFOS	EIS
64	NMeFOSE	556>122	11.25	66: d7-NMeFOSE	
65	NEtFOSE	570>136	10.75	67: d9-NEtFOSE	
66	d7-NMeFOSE	563>126	11.25	29: M4PFOS	EIS
67	d9-NEtFOSE	579>142	10.75	29: M4PFOS	EIS
68	FtS 10:2	627>607	11.50	25: M2-6:2FTS	
69	PFDoS	699>99	12.50	33: M8PFOS	

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164** Revision 15 Published Date:3/4/2019 10:48:12 AM

Page 1 of 18

1,4-Dioxane

By Gas Chromatography / Mass Spectrometry in Selected Ion Mode (GC/MS-SIM) with Isotope Dilution Modification

References:

Method 8270D, Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC/MS), Revision 4, February 2007, Test Methods for Evaluating Solid Waste, SW-846.

EPA 8000C, SW-846,Test Methods for Evaluating Solid Waste: Physical/Chemical Methods, EPA SW-846. Update III, March 2003.

1,4-Dioxane Analytical Notes, Appendix II-B-4, WSC-CAM-II-B, Revision 1, July 2010, (Massachusetts Department of Environmental Protection Bureau of Waste Site Cleanup).

1. Scope and Application

Matrices: Aqueous, soil/sediment and non-aqueous waste matrices

Definitions: Refer to Alpha Analytical Quality Manual.

This method is applicable to the quantification of 1,4-Dioxane extracted from aqueous samples in methylene chloride and analyzed by GC/MS-SIM. The extraction method is listed below and should be referenced for more details. Detection limits will vary with instrument calibration range, and volume of sample analyzed. 1,4-Dioxane detected over the calibration ranges of the instrument it is being analyzed on will be diluted and re-analyzed for accurate quantification.

The following extraction method applies:

- Extraction of Water Samples by Separatory Funnel (SOP 2165)
- Microscale Solvent Extraction (MSE) (SOP 2172)
- Organic Waste Dilution Extraction (SOP 2265)

The data report packages present the documentation of any method modification related to the samples tested. Depending upon the nature of the modification and the extent of intended use, the laboratory may be required to demonstrate that the modifications will produce equivalent results for the matrix. Approval of all method modifications is by one or more of the following laboratory personnel before performing the modification: Area Supervisor, Department Supervisor, Laboratory Director, or Quality Assurance Officer.

This method is restricted to use by or under the supervision of analysts experienced in the operation of the GC/MS-SIM and in the interpretation of GC/MS-SIM data. Each analyst must demonstrate the ability to generate acceptable results with this method by performing an initial demonstration of capability.

2. Summary of Method

Generally 1000 or 500ml of aqueous samples are serially extracted with methylene chloride in a 2 Liter Separatory Funnel (Method 3510) at a neutral pH. The extract is concentrated in MeCl2 to a 10 or 5ml final volume respectively depending on the volume of sample extracted. Soil/sediment samples are extracted by Microscale Solvent Extraction (MSE Method 3570). Approximately 5g of sample is extracted and concentrated to a 4mL final volume. Non-aqueous waste samples are

Alpha Analytical, Inc. Facility: Mansfield

Revision 15 Department: GC/MS-Semivolatiles Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM Page 2 of 18

extracted by Organic Waste Dilution Extraction (Method 3580). Approximately 1g is diluted to 10mL final volume.

ID No.:2164

Analytes are introduced into the GC/MS using a large volume injector and injecting 3ul of the calibration standards, quality control samples, and sample extracts into the GC equipped with a narrow-bore capillary column. The GC column is temperature programmed to separate the analytes, which are then detected with a mass spectrometer (MS) in selective ion mode. Identification of target analytes is accomplished by comparing their mass spectra with the electron impact spectra of the calibration standards. Concentrations are determined using mean relative response factors from a multi-level calibration curve. Response factors for target analytes and surrogate compounds are determined relative to the internal standards.

Isotope dilution quantification is achieved by spiking 1,4-Dioxane-d8 at extraction which is then in turn used as both an internal standard (IS) and surrogate. For quantification, the 1,4-Dioxane-d8 IS quantifies 1,4-Dioxane in samples. An additional IS added prior to analysis, 1,4-Dichlorobenzened4, quantifies 1,4-Dioxane-d8 as a surrogate.

2.1 Method Modifications from Reference

SIM option with Isotope dilution is utilized to increase sensitivity for this analyte. Note that while this method may be used for the evaluation of 1,4-Dioxane in soil/sediment samples, the MA DEP CAM does not approve of this method for determining presumptive certainty for MA DEP cleanup sites. Method 8260 should be employed in these cases.

3. Reporting Limits

Concentrations for 1,4-Dioxane can be detected in water samples in the range of 150ng/L to 100,000ng/L. Soil/sediment samples can be detected in the range of 8ug/Kg to 8,000ug/Kg. Standard reporting limits for aqueous samples are 150ng/L and for soil/sediment samples 8ug/kg respectively.

4. Interferences

- 4.1 Phthalate esters can be a major source of contamination if any material containing plasticizers (phthalates) comes in contact with the sample during the extraction process. Use of plastic or any material containing plasticizers (phthalates) should be avoided during extraction or analysis.
- 4.2 The injection port of the gas chromatograph can become contaminated with high boiling compounds resulting in the loss of sensitivity. It may be necessary to replace the injection port liner routinely to prevent this loss of sensitivity. Clipping off approximately four inches of the column at the injection end may also increase sensitivity. Low instrument response can be detected during the daily tuning procedure by including pentachlorophenol and benzidine in the daily tuning mix.
- 4.3 Raw GC/MS data from all blanks, samples, and spikes must be evaluated for interferences or carryover. Contamination by carryover can occur whenever high-concentration and lowconcentration samples are sequentially analyzed.
- 4.4 Solvents, reagents and glassware may introduce interferences. These must be demonstrated to be free of interferences by the analysis of a method blank. See the SOP Reagent, Solvent and Standard Control (G-008) and Laboratory Glassware Cleaning (G-002), for additional details.

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Alpha Analytical, Inc. Facility: Mansfield Department: GC/MS-Semivolatiles

Revision 15 Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM Page 3 of 18

4.5 It should be noted that there are some chromatographic consequences observed due to the acetone used in the MSE extraction process. This is characterized by a shift in retention time for 1,4-Dioxane and the 1,4-Dioxane-d8 surrogate (~0.4minutes), as well as the presence of a large peak (likely acetone) in all chromatograms.

ID No.:2164

5. Health and Safety

The toxicity or carcinogenicity of each reagent and standard used in this method is not fully established; however, each chemical compound should be treated as a potential health hazard. From this viewpoint, exposure to these chemicals must be reduced to the lowest possible level by whatever means available. A reference file of material safety data sheets is available to all personnel involved in the chemical analysis. Additional references to laboratory safety are available in the Chemical Hygiene Plan.

All personnel handling environmental samples known to contain or to have been in contact with municipal waste must follow safety practices for handling known disease causative agents.

6. Sample Collection, Preservation, Shipping and Handling

6.1 Sample Collection

Solid samples; A minimum of 100grams of sample must be collected in a glass jar with a Teflon lined screw cap.

Water samples: A minimum of 0.5 to 2 liters of sample must be collected in amber glass bottles.

6.2 Sample Preservation

Solid samples: The sample must be refrigerated and maintained at 4±2 °C until extraction and analysis. Sediment samples can be frozen at -20±5 °C until extraction to extend hold time. The extracts must be refrigerated and maintained at 4±2 °C until analysis.

Water samples: The samples must not be preserved except by refrigeration at 4±2 °C until extraction and analysis. The extracts must be refrigerated and maintained at 4±2 °C until analysis.

Non-aqueous waste samples: Concentrated sample extracts must be stored in contaminantfree containers and preserved in a refrigerator when not used for more than four hours.

6.3 Sample Shipping

No special shipping requirements.

6.4 Sample Handling

Solid samples: All solid samples must be extracted within 14 days from the date of collection. Frozen sample hold times are monitored up to 14 days from the date removed from freezer. The extracts must be refrigerated and maintained at 4±2 °C until analysis. Sample extracts must be analyzed within 40 days from date of extraction.

Water samples: All water samples must be extracted within 7 days from the date of collection. Sample extracts must be analyzed within 40 days from date of extraction.

Non-aqueous waste samples: Hold times do not apply to neat oils/NAPL/product samples.

7. Equipment and Supplies

7.1 Gas chromatograph - Programmable, heating range from 40C to 350C; splitless-type inlet system, (Hewlett Packard 6890N Series II or similar); mass selective detector (Hewlett

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164**Revision 15
Published Date:3/4/2019 10:48:12 AM
Page 4 of 18

Packard 5973, or similar); automatic injector (Hewlett Packard 7683B or similar).

- 7.2 Chromatography Column Fused silica capillary column, 0.25mm ID x 60m length, 0.25mm film thickness RTX-5, Restek Corporation, 5% diphenyl-95% dimethyl polysiloxane, Fused silica capillary column, 0.18mm ID x 60m length, 0.18um film thickness RTX-PCB, Restek Corporation, Fused silica capillary column, 0.25mm ID x 30m length, .25mm film thickness (Zebron ZB-SemiVolatiles, Phenomenex Corporation, 5% Polysilarylene 95% Polydimethylsiloxane), or equivalent.
- **7.3 Gerstel Large Volume Injection System** Temperature programmable range from 0C to 350C; pressure programmable; capable of split or splitless injection; Injection volumes range from 1 to 50uL and Cryo cooling availability to allow for cold injections.
- 7.4 Agilent Split/Splitless injector System
- **7.5 Data Acquisition System -** Computerized system for collecting, storing, and processing detector output (Hewlett Packard Enviroquant target software) or equivalent.
- **7.6 Gases -** BIP Ultra high purity helium (99.9995%); Compressed nitrogen for N-Evap. Carbon dioxide (siphon type) for Gerstel.
- **7.7 Syringes** 10uL to 1.0mL
- 7.8 Vials- including 2ml, 4ml, 10ml, 40ml and other sizes as necessary.
- 7.9 Hamilton Gas tight Syringes varying sizes
- 7.10 Gerstel Single baffle injection port liners- packed lightly with glass wool.
- 7.11 GC Injection Port Liner: Phenomenex Direct Connect Top Hole
- **7.12 Class A Volumetric flasks:** Including 10ml, 20 ml, 50 ml, 100ml and other sizes as necessary

8. Reagents and Standards

Use reagent grade chemicals for all reagents. Deionized (DI) water is ASTM Type II laboratory reagent grade water.

- 8.1 Solvents: All solvent expirations determined as indicated by manufacturer guidelines
 - **8.1.1** Methylene Chloride, ACS approved, Pesticide grade, see SOP *Reagent, Solvent and Standard Control* (SOP 1816) for additional details regarding solvent purity. Used to extract samples and prepare instrument/analytical standards.
 - **8.1.2** Acetone, ACS approved, Pesticide grade, see SOP *Reagent, Solvent and Standard Control* (SOP 1816) for additional details regarding solvent purity. This water soluble solvent is used for surrogate and LCS/MS preparation.
 - **8.1.3** Methanol, ACS approved, Pesticide grade, *Reagent, Solvent and Standard Control SOP* (SOP 1816) for additional details regarding solvent purity.

Alpha Analytical, Inc.

Facility: Mansfield

Department: GC/MS-Semivolatiles

ID No.:2164
Revision 15
Published Date:3/4/2019 10:48:12 AM

Title: 1,4-Dioxane By GC/MS-SIM Page 5 of 18

8.2 Analytical Standards: Standards should be stored at -10C or less, away from light when not in use. They should be discarded after 1 year unless the vendor expiration date states otherwise or, if degradation is observed. Stock standards are given a 1 year expiration from the preparation date or the expiration of the primary vendor solution, whichever occurs first. Working standards are given six month expiration from the preparation date or the expiration of the primary solution whichever occurs first. All analytical standards are made up in Methylene Chloride. All prep standards are made up in Acetone.

8.3 Surrogate/Internal Standard (IS):

- 8.3.1 A 1,4-Dioxane-d8 Primary neat standard is commercially obtained from Cambridge Isotope (Cat #DLM-28-10 or equivalent). A stock surrogate/internal standard solution is prepared by weighing 0.1 g of the primary neat standard and diluting volumetrically in 10 mls methylene chloride (or equivalent preparation) to obtain a concentration of ~10,000 ug/ml.
- 8.3.2 From this stock (Section 8.3.1), the Surrogate spiking solution is made by a serial dilution (1 ml diluted up in a 50 ml volumetric flask, followed by a 0.5 ml of this solution diluted up in a 20 ml volumetric flask or equivalent preparations) in Acetone to achieve a concentration of 5 ug/mL. Of this surrogate solution, 1 mL is spiked into each water sample, and 0.4 mL is spiked into each soil/sediment sample for a surrogate/IS concentration of 500 ng/mL in samples.
- **8.3.3** The initial stock solution (Section 8.3.1) is diluted volumetrically (250 uL diluted up in a 100 ml volumetric flask or equivalent preparation) in methylene chloride to obtain an IS solution at a concentration of 25 ug/ml for the spiking of only the analytical standards (i.e. calibration curve and continuing calibration).
- **8.4 Internal Standard:** 1,4-Dichlorobenzene-d4 commercially obtained from Restek (Cat #31206 or equivalent). This primary solution is at 2000 ug/mL and contains other Semivolatile Internal standards, however this method only utilizes the 1,4-Dichlorobenzene-d4. This solution is diluted volumetrically (250 uL diluted up in a 50 ml volumetric flask or equivalent preparation) in methylene chloride to obtain a SIM-IS solution at a concentration of 25 ug/ml. All samples and standards are spiked with 20uL of internal standard before analysis. This IS is intended to be used for both quantitation (of the surrogate 1,4-Dioxane-d8) and the establishment of relative retention times.
- 8.5 Matrix Spike/Laboratory Control Spike Solutions (MS/LCS) 1,4-Dioxane primary spike solution is commercially obtained from Restek at 2000 ug/mL (Cat #31853 or equivalent). To prepare the working spike solution, syringe measure 0.25 mL of stock and bring to 100 mL in acetone (or equivalent preparation) for a 5 ug/mL concentration. From this solution, 1 mL is spiked into each MS/LCS QC water sample and 0.4 mL is spiked into each soil/sediment sample. The final concentration in a 10 mL water extract and in a 4 mL soil/sediment sample will be 500 ng/mL.
- 8.6 Stock Calibration Standards are prepared from a primary standard commercially obtained from Restek at a concentration of 2000 ug/ml (Cat #31853 or equivalent). The primary standard is used to appropriately prepare a stock standard at a concentration of 10,000 ng/ml (125 uL diluted up in a 25 ml volumetric flask or equivalent preparation). The stock solution is generally made up in 25 ml increments. The stock solution then doubles as a high level standard for the calibration curve, as well as a stock solution by which to serially dilute to prepare the other calibration standards. Except for the high level standard, the calibration curve levels are generally each made at 10 ml. Refer to the table below for example calibration curve levels. All the calibration curve levels, once aliquoted at 1 ml increments for injection,

Printouts of this document may be out of date and should be considered uncontrolled. To accomplish work, the published version of the document should be viewed online.

Alpha Analytical, Inc. ID No.:2164 Facility: Mansfield Revision 15 Department: GC/MS-Semivolatiles Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM

have 20 ul of 1,4 Dichlorobenzene IS (SIM IS - Section 8.4) added as well as 20 uL of 1,4 Dioxane-d8 IS/Surrogate (Section 8.3.3) added.

Page 6 of 18

8.7 Curve Preparation: The following 7 levels are the possible levels that can be analyzed for an ICAL.

Calibration Level	Volume of Stock Std Added
	Into 10 ml of DCM
Level 1 (10 ng/mL)	10 uL
Level 2 (50 ng/mL)	50 uL
Level 3 (100 ng/mL)	100 uL
Level 4 (500 ng/mL)	500 uL
Level 5 (1000 ng/mL) - CCV	1.0 mL
Level 6 (5000 ng/mL)	5.0 mL
Level 7 (10,000 ng/mL)	NA

8.8 The Independent Check Verification Standard is prepared from a Primary standard solution commercially obtained from Ultra (Cat NV-150-1 or equivalent) at a concentration of 100 ug/ml. The primary solution is diluted volumetrically (100 ul diluted up in a 10 ml volumetric flask or equivalent preparation) in methylene chloride to create a working ICV solution at a concentration of 1000 ng/ml. The working ICV solution is aliquoted into 1 ml increments as needed to inject with each calibration curve. Additionally, 20 ul of SIM IS (Section 8.4) and 20 ul of 1,4-Dioxane-d8 IS (Section 8.3.3) is added into the 1 ml aliquot.

9. Quality Control

The laboratory must maintain records to document the quality of data that is generated. Ongoing data quality checks are compared with established performance criteria to determine if the results of analyses meet the performance characteristics of the method.

9.1 Blank(s)

A method blank must be prepared once per every 20 samples or per extraction batch, whichever is more frequent.

Organic compounds of interest must not be detectable in the method blank at a concentration greater than the reporting limit.

Corrective Action: For contaminated blanks, all efforts must be made to identify and eliminate the source of contamination. The presence of analytes at concentrations at or above the reporting limit will warrant application of a "B" qualifier to that target compound(s) on all associated report forms, and perhaps re-extraction of all associated samples. Re-extraction of the method blank and all associated samples must be performed until the blank is in control. Surrogate recoveries must meet the QC limits for the method blank. Re-extraction must be initiated immediately so that minimum time is wasted before re-extraction can occur - if at all possible-this re-extraction should take place within holding time. Re-extraction corrective action that would exceed the sample holding time criteria should be discussed with the Organics Supervisor, Project Manager, client, and Operations Manager prior to implementation. Exceptions may be made with approval of the Organics Section Supervisor if the samples associated with an out of control method blank are non-detect for the affected compound(s) or if the concentration of the affected compound(s)in the sample is greater than 10x the blank level. In such cases, the sample results are accepted without corrective action for the high method blank result. The client must be notified, via the project narrative, of any method blank non-compliance associated with sample results

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164** Revision 15 Published Date:3/4/2019 10:48:12 AM

Page 7 of 18

9.2 Laboratory Control Sample/Laboratory Control Sample Duplicate (LCS/LCSD)

Laboratory control samples (LCS/LCSD) must be prepared once per every 20 samples or per extraction batch, whichever is more frequent, and spiked with 1,4 Dioxane spike solution (Section 8.5) and surrogate (Section 8.3.2) before extraction. The IS (Section 8.4) is spiked after extraction and before the analysis.

Acceptable Recovery limits are 40% - 140%. The relative percent difference (RPD) between the LCS/LCSD is 30%. Limits are adapted from MCP protocol and are generally monitored and documented in-house through control charts.

<u>Corrective Action</u>: Analysis must be repeated if an analytical error is suspected. If the LCS/LCSD recoveries and/or %RPD are still out of control, re-extract and re-analyze the LCS/LCSD <u>and all</u> associated samples. Samples cannot be reported until an acceptable LCS is obtained.

9.3 Initial Calibration Verification (ICV)

Refer to Section 10.2.

9.4 Continuing Calibration Verification (CCV)

Refer to Section 10.4.

9.5 Matrix Spike / Matrix Spike Duplicate (MS/MSD)

Matrix spike / matrix spike duplicate (MS/MSD) samples are performed upon project specifications. They are performed per client request. The sample is spiked with 1,4-Dioxane spike solution (Section 8.5) and surrogate (Section 8.3.2) before extraction. The IS (Section 8.4) is spiked after extraction and before the analysis. The recovery limits are 40% - 140% and the RPD limit is 30%. Limits are adapted from MCP protocol and are generally monitored and documented in-house through control charts.

<u>Corrective Action</u>: Analysis must be repeated if an analytical error is suspected. If the % recovery and/or %RPD still exceeds the control limits and the LCS/LCSD is compliant; include a project narrative with the results to client noting that there may be potential matrix effects on the accuracy or precision of the reported results as evidenced by MS/MSD recoveries and/or %RPD outside of QC limits.

9.6 Laboratory Duplicate

Duplicate analyses are performed upon client and/or workplan request. For Organic analyses, the matrix duplicate is usually in the form of the matrix spike duplicate, see Section 9.5.

Acceptable relative percent difference (RPD) of duplicates is 30%. Acceptance criterion is not applicable to sample concentrations less than 5 times the reporting limit. Calculate the RPD as follows:

$$RPD = \frac{R1 - R2}{[R1 + R2]} \times 100$$

where:

R1 = sample Replicate #1 R2 = sample Replicate #2

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164** Revision 15 Published Date:3/4/2019 10:48:12 AM

Page 8 of 18

The RPD limits should be monitored and documented in-house through control charts and updated as needed.

<u>Corrective Action</u>: Analysis must be repeated if an analytical error is suspected.. If the % RPD still exceeds the control limits; include a project narrative with the results to client noting that there may be potential matrix effects on the precision of the reported results as evidenced by the matrix duplicate % RPD exceedence.

9.7 Method-specific Quality Control Samples

9.7.1 Surrogates

Surrogate spikes (Section 8.3.2 and Section 8.3.3) must be added to QC and field samples to evaluate the extraction method performance.

The acceptable surrogate recovery limits are 15% - 110%. Limits are adapted from MCP protocol.

Corrective Action: Analysis must be repeated if an analytical error is suspected. If the % recovery still exceeds the control limits the sample must be re-extracted and re-analyzed to confirm the sample matrix. If obvious matrix interferences are noted, consultation with the Organic Supervisor or Operations Manager may be in order to confirm the need for sample re-extraction. If no re-extraction occurs, the surrogate results and reasons for the decision not to re-extract must be discussed in the project narrative to the client. Due to the isotope dilution nature of this method, the analyst must pay close attention to surrogate recoveries and areas as this recovery is then used to surrogate correct (as an Internal standard) the response of 1,4-Dioxane.

9.7.2 Internal Standards

Internal standards must be added to all sample extracts, QC samples and standards for quantitation purposes. For sample extracts, only the 1,4-Dichlorbenzene-d4 IS (SIM IS – Section 8.4) should be added since the extraction surrogate is then used as an Internal Standard upon analysis. However, all Calibration standards must be spiked with both 1,4-Dichlorobenzene-d4 as well as 1,4-Dioxane-d8 (Section 8.3.3). In the sample extracts the area counts for 1,4-Dioxane-d8 will vary based on the extraction, however the 1,4-Dichlorobenzene-d4 area counts should remain fairly constant. The 1,4-Dichlorobenzene-d4 internal standard should also remain constant with respect to the continuing calibration analyzed at the beginning of the run. Sample IS areas must be –50% to +100% of the Internal Standards in the Continuing Calibration – applies only to 1,4-Dichlorobenzene-d4. Additionally, the IS retention times should not differ more than 30 seconds from the Continuing Calibration.

<u>Corrective Action</u>: Analysis must be repeated once unless there are obvious samples matrix interferences, i.e., the sample extract was very colored and viscous, or there are obvious chromatographic interferences. If *obvious* matrix interferences are noted, consultation with the Organic Supervisor or Operations Manager may be in order to confirm the need for sample re-analysis or re-extraction.

Alpha Analytical, Inc. Facility: Mansfield

Revision 15 Department: GC/MS-Semivolatiles Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM Page 9 of 18

ID No.:2164

9.8 Method Sequence

- Tune
- CCV
- Method Blank
- LCS
- LCSD
- Samples

10. **Procedure**

10.1 Equipment Set-up

10.1.1 The instrument used for the analysis is a HP 6890N Series gas chromatograph. The HP system is equipped with a Gerstel large volume injection system, and a 7683B-type autosampler or equivalent. The mass spectrometer is an HP 5973 or 5975 with the HP Enviroquant data system. The method is modified for selective ion monitoring. The table below lists the ions monitored in one SIM window. This method must only be set up and analyzed by an experienced mass spectrometrist.

Compound	Primary Ion	Secondary Ion
1,4-Dioxane	88	58, 43
1,4-Dioxane-d8	64	96
1,4-Dichlorobenzene-d4	152	115

10.1.2 The basic GC parameters are as follows for the Gerstel Large Volume Injection System:

Oven Equib Time: 0.10 min

Oven Max: 325 Initial Temp.: 45°C Initial Time: 6.00 min

Level	Rate (°C/min)	Final Temp. (°C)	Final Time (min)
1	19.00	120	1.0
2	11.00	150	1.00
3	19.00	305	4.00

Final Time: 26.83

10.1.3 The basic injection port parameters are as follows for the Gerstel Large Volume Injection System:

"Splitless" mode Initial Temp: 46 °C Initial Time: 0.30 min Ramp Rate: 300°/second Final Temp: 300°C Final Time: 30.0

Cryo: ON

Cryo Use temp: 25 °C Cryo Timeout: 30.0 min

Cryo Fault: ON

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles Title: 1,4-Dioxane By GC/MS-SIM ID No.:**2164** Revision 15 Published Date:3/4/2019 10:48:12 AM

Page 10 of 18

Purge Flow: 25 ml/min Purge Time: 2.50 min Gas Saver: off Gas: Helium

Mode: Constant Pressure Average Velocity: 30cm/sec Initial flow: 1.3 ml/min

10.1.4 The basic GC parameters are as follows for the Agilent Split/Splitless injector System:

Oven Equib Time: 0.20 min

Oven Max: 325°C Initial Temp.: 30°C

Level	Rate (°C/min)	Final Temp. (°C)	Final Time (min)
1	0	30	2.00
2	5	50	0.00
3	50	300	2.00

Final Time: 13.00

10.1.5 The basic injection port parameters are as follows for the Agilent Split/Splitless injector System:

Mode: Splitless Temperatuire: 200 °C Flow: 1.0 mL/min Velocity: 36.074 cm/sec Septum Purge: 3mL/min

Purge Flow to Split Vent: 60 mL/min at 0.3 min

10.1.6 MS Acquisition Information:

Tune file: dftpp.u Acquisition Mode: SIM Solvent Delay: 7.70 min MS Source temp: 280°C

10.1.7 Tuning

- 10.1.7.1 Before the analytical standards are analyzed the mass spectrometer must be adjusted to meet the proper ion criteria for DFTPP. This is demonstrated by injecting into the GC/MS system 1uL of a 50ug/mL DFTPP solution. After the analysis of the DFTPP, evaluate the tune as follows:
 - Enter into the "Environmental Data Analysis" (off-line) screen.
 - · Go to "File" and select the tune data file.
 - Go into "Tuner" and select "Eval DFTPP", then select "AutoFind DFTPP to Screen," to evaluate the tune file, based on the pre-set SW-846 criteria. The

Printouts of this document may be out of date and should be considered uncontrolled. To accomplish work, the published version of the document should be viewed online.

Document Type: SOP-Technical Pre-Qualtrax Document ID: O-016

:

Alpha Analytical, Inc. Facility: Mansfield Department: GC/MS-Semivolatiles

Revision 15 Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM Page 11 of 18

> software will evaluate the tune by selecting three scans of the DFTPP peak and will display the ion intensities on the screen. That is, one scan at the apex, one scan directly preceding the apex and one scan following the apex and averages them, then takes one background subtracted scan, 20 seconds before the beginning of the DFTPP peak. If the criteria below are met, repeat. select "AutoFind to Printer", for a hardcopy of the tune evaluation for the record.

ID No.:2164

If the "AutoFind" tune evaluation does not meet the criteria below, manual evaluation of the tune can be performed by attempting either of the options below:

- Blow up the DFTPP peak on the screen and select either one single scan at the apex of the peak, or a scan immediately preceding or following the apex. Go into "Tuner" and select "Evaluate DFTPP to Screen," or "Evaluate DFTPP to Printer," as described above, OR,
- Take the average of the scans across the entire peak. Go into "Tuner" and select "Evaluate DFTPP to Screen," or "Evaluate DFTPP to Printer," as described above.
- 10.1.7.2 The following DFTPP mass intensity criteria should be used.

DFTPP KEY MASSES AND ABUNDANCE CRITERIA

Mass	m/z Abundance criteria
51	10-80 percent of mass 198.
68	Less than 2 percent of mass 69.
70	Less than 2 percent of mass 69.
127	10-80 percent of mass 198.
197	Less than 2 percent of mass 198.
198	Base peak, or >50 percent of Mass 442.
199	5-9 percent of mass 198.
275	10-60 percent of mass 198.
365	Greater than 1 percent of mass 198.
441	Present but less than 24 percent of mass 442.
442	Base Peak, or > 50 percent of mass 198.
443	15-24 percent of mass 442.

10.1.8 Tune acceptance should be verified at the beginning of every 12 hour analytical shift. The DFTPP may be combined with the calibration verification standard as long as both tuning and calibration acceptance criteria are met.

10.2 Initial Calibration

10.2.1 After the DFTPP passes criteria, a set of multi-level calibration standards listed in Section 8.7 are analyzed, from low concentration to high. A minimum of five calibration levels are analyzed. The calibration standards are stored in amber vials in the standards freezer. The labeling convention allows each standard to have a unique identifier which distinguishes it from field samples. The naming convention used throughout the laboratory identifies the standard as semivolatile, hydrocarbon, pesticide/PCB or volatile. An example of this would be SW042407E, meaning it is a semivolatile (S) working (W) standard made on April 24, 2007 and that it was the fifth standard made that day. All certificates of analysis that are shipped with standards are filed with their receipt ID written on it to insure traceability.

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Alpha Analytical, Inc.

Facility: Mansfield

Department: GC/MS-Semivolatiles

Title: 1,4-Dioxane By GC/MS-SIM

ID No.:2164

Revision 15

Published Date:3/4/2019 10:48:12 AM

Page 12 of 18

10.2.2 Once the standards have been analyzed, they are reduced by the search software of the Enviroquant data system. Once all the components are identified, a linear curve is calculated for the components. The criteria for evaluation are as follows:

- 10.2.2.1 The average RF for each compound must be greater than 0.05.
- 10.2.2.2 The %RSD of each compound must not exceed 20%. If they do this may be an indication that the chromatographic system is too reactive for analysis to begin. This indicates the instrument may need maintenance.
- Alternatively, a linear regression model may be employed, provided that the coefficient of determination (COD or r2) is ≥0.99. Otherwise, construct a nonlinear calibration of no more than a third order equation. Statistical considerations in developing a non-linear calibration model require more data than the more traditional linear approach. A quadratic (second order) model requires six standards, and a third order polynomial requires seven standards. In setting model parameters, do not force the line through the origin. The COD or r² must be greater than or equal to 0.99.
- 10.2.2.4 Once the calibration curve is reviewed, an Initial Calibration Checklist must be completed.
- 10.2.3 All samples and standards are spiked with Internal Standards (IS) before analysis. Refer to section 8.4 for specific internal standard spiking information. The IS is intended to be used for both quantitation and the establishment of relative retention times. Internal standard acceptance criteria can be found in Section 9.7.2.
- 10.2.4 Independent check standards (Section 8.8) from a separate source or different lot are analyzed after every initial calibration for evaluation against calibration standard solutions. The % Difference (%D) should not be greater than +30%.

10.3 Equipment Operation and Sample Processing

10.3.1 Tuning

A DFTPP standard must be analyzed and pass criteria before a continuing calibration verification standard or any samples are analyzed. The DFTPP may be combined with the calibration verification standard as long as both tuning and calibration acceptance criteria are met. A DFTPP tune standard must be analyzed before each 12 hour analytical shift. Please refer to Section 10.1.5 for tuning criteria and other information.

10.3.2 Daily Calibration

On a daily basis after the DFTPP has passed, a mid-level (usually 1000 ng/mL) continuing calibration standard which contains all of the analytes of interest is analyzed. The criteria for acceptance are:

- **10.3.2.1** All analytes must have response factors greater than 0.05.
- 10.3.2.2 The % D must be ±20% D from the initial calibration.

10.3.3 Sample Analysis

10.3.3.1 The prep lab staff will transfer the samples to the instrument laboratory. The samples are generally brought to a 10 mL final volume for liquids or 4 ml for soil samples; 1 mL is transferred and the remaining sample volume is archived. One

Alpha Analytical, Inc.

Facility: Mansfield

Department: GC/MS-Semivolatiles

Title: 1,4-Dioxane By GC/MS-SIM

ID No.:2164

Revision 15

Published Date:3/4/2019 10:48:12 AM

Page 13 of 18

aliquot of each sample is then placed in the sample extract holding refrigerator located in the instrument laboratory.

- 10.3.3.2 All of the samples at 1 mL (including the batch QC samples) are spiked with 20 uL internal standard (see section 8.4 for specifics regarding the internal standard). The samples are shaken briefly after the internal standard is added to ensure mixing. A sample will need a dilution for target analytes that are over calibration.
- 10.3.3.3 After the daily tune and CCAL have passed criteria, the analyst places the samples onto the autosampler tray. (Generally, the samples will be analyzed in order of color (lightest to darkest.) with QC samples being analyzed first. The instrument sequence is typed into the HP Chemstation Software. Next run "Simulate Sequence" (also under the "Sequence" dropdown list) without clicking the "Overwrite Files" box. Compare the order of the vials on the instrument versus the sequence to confirm all the samples and standards are in the right places. Next click "Run Sequence" also under the HP Chemstation "Sequence" dropdown list.
- 10.3.3.4 After the samples have been analyzed, the data files from the MS are quantitated versus the proper quantitation method. The QCPRN1.MAC macro creates a form with which to easily check internal standard and surrogate criteria are met. The following should be reviewed initially:
 - **10.3.3.4.1** Are all the surrogates within QC criteria? Please see Section 9.7.1 for surrogate information.
 - 10.3.3.4.2 Is the internal standard- 1,4 Dichlorobenzene-d4, within 50-200% of the daily CCAL? If not, the samples should be checked for matrix interferences that may be causing these issues. The IS peaks should also be evaluated for peak splitting or incorrect integration by the software. A sample may not need to be reanalyzed if it can be determined (with guidance from a supervisor) that the QC is exceeded due to matrix interference.
 - 10.3.3.4.3 Are all target analytes within calibration range? If not, the sample(s) should be diluted and re-analyzed. If a dilution is performed after the internal standard has already been added, it will be necessary to add additional IS in order to make up for the impact of the original IS added also being diluted. Conversely, if a sample has been over-diluted, it may need to be analyzed at less of a dilution to detect target analytes that may have been diluted out. Note: for this method, due to the isotope dilution, it is necessary to refer to area counts to determine if dilution is required. The area response for 1,4-Dioxane should be compared to the area of the highest level of the ICAL standard for that target. If the area of 1,4-Dioxane in the extract is greater than the area of 1,4-Dioxane in the Highest level of the Calibration Curve, then a dilution is necessary. The analyst CAN NOT use the concentration of 1,4-Dioxane found in the extract to determine if dilution is required because the concentration of 1,4-Dioxane in the sample is surrogate corrected in Chemstation based on the use of the extraction surrogate as an Internal Standard. Once it is determined based on the peak area that a sample requires a dilution, the analyst must pay attention to the concentration of IS in the dilution. Although the analyst will adjust the 1,4-Dichlorobenzene-d4 concentration by adding the appropriate amount of additional IS to the

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164**Revision 15
Published Date:3/4/2019 10:48:12 AM
Page 14 of 18

dilution, there is no way for the analyst to make up for the dilution of the extraction surrogate 1,4 dioxane-d8 (also used as internal standard). To account for this the analyst must change the concentration of the 1,4-Dioxane-d8 within the calibration table to a value which accounts for the dilution. (Example: If the analyst performs a 10x dilution then the concentration of 1,4-Dioxane-d8 must be changed within the ICAL to 50.) The analyst must ensure that only the diluted sample is calculated against this value and that the method is returned to the proper concentration once the analysis of the dilution is complete.

- **10.3.3.4.4** Are all analyses within 12 hour tune time? If a sample is analyzed outside tune time, it will need to be re-analyzed in another tune clock.
- **10.3.3.4.5** The sequence should also be printed out from Chemstation, initialed and dated, and placed in the logbook.
 - 10.3.3.4.5.1 If anything in the initial review of the data indicates that there should be a re-analysis, the reason for re-analysis should be noted on the sequence.
 - 10.3.3.4.5.2 Also, note the time the tune standard was analyzed, the time of the last sample analysis. If there are samples analyzed outside tune time, the time the last sample within tune time ran should be noted, as well as the last sample analyzed.
 - 10.3.3.4.5.3 If a re-extract is required, the "Request for Repreparation/Reclean" book should be filled out and a photocopy of the appropriate page should be given to the Preparation Group leader or the Organics Section Head.

10.4 Continuing Calibration

Continuing Calibration Criteria is outlined in section 10.3.2

10.5 Preventive Maintenance

If performing any maintenance on any piece of equipment it must be documented in the *Instrument Maintenance Logbook* located in the laboratory specific to each instrument.

Daily

Injection port maintenance: Maintenance should be done when the daily CCAL starts to demonstrate degradation either by %D outliers or area responses <50% as compared to the ICAL areas. Several tune clocks may be injected before maintenance is needed. The type of samples analyzed will have an effect on how soon maintenance should be performed. Injection port maintenance should be done as needed. General maintenance includes replacing the single baffle liner packed lightly with glass wool, cutting about 2-4 inches off the head of the column, and replacing the septa. Refer to maintenance log for more specific information.

The Gerstel Injection port should be handled with care. The liners are quite thin. Do not force the Gerstel weldmen into place as the threads are soft metal and will cross thread. Always make sure the weldmen goes on straight. The Gerstel injection port does not require significant tightening of either the weldmen or column nut, tighten enough to seal but there is no need to crank down on it.

Alpha Analytical, Inc. Facility: Mansfield Department: GC/MS-Semivolatiles

Revision 15 Published Date: 3/4/2019 10:48:12 AM Title: 1,4-Dioxane By GC/MS-SIM Page 15 of 18

11. Data Evaluation, Calculations and Reporting

11.1 Qualitative Analysis

11.1.1 The qualitative identification of compounds determined by this method is based on retention time and on comparison of mass spectrum, after background correction, with characteristic ions in a reference mass spectrum. The reference mass spectrum must be generated by the laboratory using the conditions of this method. The characteristic ions from the reference mass spectrum are defined as the two ions of greatest relative intensity, and are over 30% relative intensity. Compounds are identified when the following criteria are met.

ID No.:2164

- 11.1.2 The intensities of the characteristic ions of a compound must maximize in the same scan or within one scan of each other. A peak selected by the data system, based on the presence of target specific ions at a target specific retention time will be accepted as meeting these criteria.
- 11.1.3 The relative retention time of the sample component is within + 0.06 RRT units of the RRT of the standard component.
- 11.1.4 The relative intensities of the characteristic ions agree within 30% of the relative intensities of these ions in the reference spectrum (Example: For an ion with an abundance of 50% in the reference spectrum, the corresponding abundance in a sample spectrum can range between 20% and 80%.) The relative intensities are monitored daily. The relative intensities will be updated when they exceed established values from the reference spectrum.
- 11.1.5 Structural isomers that produce very similar mass spectra should be identified as individual isomers if they have sufficiently different GC retention times. Sufficient GC resolution is achieved if the height of the valley between two isomer peaks is less than 25% of the sum of the two peak heights. Otherwise, structural isomers are identified as isomeric pairs.
- 11.1.6 Identification is hampered when sample components are not resolved chromatographically and produce mass spectra containing ions contributed by more than one analyte. When gas chromatographic peaks obviously represent more than one sample component (i.e. a broadened peak with shoulder(s) or a valley between two or more maxima), appropriate selection of analyte spectra and background spectra is important. Selective ion monitoring eliminates this potential.

11.2 Quantitative Analysis

11.2.1 Response factors and % RSD to evaluate Initial Calibration acceptability.

$$RF = \frac{area_{cmp}}{area_{is}} \times \frac{conc_{is}}{conc_{cmp}}$$
 Calculate RF by:

where:

area cmp = Area of the characteristic ion for the compound being measured. area is = Area of the characteristic ion for the specific internal standard. conc is = Concentration of the specific internal standard. conc cmp = Concentration of the compound being measured.

Pre-Qualtrax Document ID: O-016 Document Type: SOP-Technical

Alpha Analytical, Inc. Facility: Mansfield

Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164** Revision 15 Published Date:3/4/2019 10:48:12 AM

Page 16 of 18

$$\%RSD = \frac{SD}{\overline{x}} \times 100$$

 $SD = \sqrt{\sum_{i=1}^{N} \frac{(x_i - \overline{x})^2}{N - 1}}$

Calculate % RSD by:

where:

% RSD = percent relative standard deviation

x = average of RF's

SD = standard deviation

xi = analytical results of each level in the final reporting units

N = number of results (levels)

11.2.2 Calculate % Difference (%D) by:

$$\%D = \frac{\overline{R}\,\overline{F}_i - RF_c}{\overline{R}\,\overline{F}_i} \times 100$$

where:

RFi - Initial Calibration average RF RFc = Continuing Calibration RF

11.2.3 Results of Water Analysis - calculation as performed in report form:

Concentration (ug/L) =
$$\frac{\text{(Conc) (Vf) (DF)}}{\text{(Vi)}} \times 1000$$

11.2.4 Results of Soil/Sediment Analysis – calculation as performed in report form:

Concentration (ug/Kg) =
$$\frac{\text{(Conc) (Vf) (DF)}}{\text{(Vi) \%S}} \times 1000$$

where:

Conc = Raw on-column concentration obtained from the quantitation report using Initial Calibration results.

Vf = Final volume of extract (mL)

Vi = Volume of sample extracted (mL), or weight of sample extracted in grams (g)

DF = Dilution factor, for manually prepared dilutions, not instrumental "dilutions".

%S = percent solids, as a decimal

Document Type: SOP-Technical

Pre-Qualtrax Document ID: O-016

Alpha Analytical, Inc.
Facility: Mansfield
Department: GC/MS-Semivolatiles
Title: 1,4-Dioxane By GC/MS-SIM

ID No.:**2164**Revision 15
Published Date:3/4/2019 10:48:12 AM
Page 17 of 18

12. Contingencies for Handling Out-of-Control Data or Unacceptable Data

Section 9.0 outlines sample batch QC acceptance criteria. If non-compliant organic compound results are to be reported, the Organic Section Head and/or the Laboratory Director, and the Operations Manager must approve the reporting of these results. The laboratory Project Manager shall be notified, and may chose to relay the non-compliance to the client, for approval, or other corrective action, such as re-sampling and re-analysis. The analyst, Data Reviewer, or Department Supervisor performing the secondary review initiates the project narrative, and the narrative must clearly document the non-compliance and provide a reason for acceptance of these results.

All results for the organic compounds of interest are reportable without qualification if extraction and analytical holding times are met, preservation requirements (including cooler temperatures) are met, all QC criteria defined in the table below are met, and matrix interference is not suspected during extraction or analysis of the samples. If any of the below QC parameters are not met, all associated samples must be evaluated for re-extraction and/or re-analysis.

QC Parameter	Acceptance Criteria
Method Blank	No analyte above the reporting limit The results are qualified
	with a "B" for any associated sample concentrations that are
	less than 10x the blank concentration for this analyte
Surrogate Recovery	15% - 110%
Laboratory Control Samples	40% - 140% and 30% RPD
Matrix Duplicate	30% RPD
Matrix Spike	40% - 140%
Matrix Spike Duplicate	30% - 140% and 30% RPD

13. Method Performance

13.1 Method Detection Limit Study (MDL) / Limit of Detection Study (LOD) / Limit of Quantitation (LOQ)

The laboratory follows the procedure to determine the MDL, LOD, and/or LOQ as outlined in Alpha SOP 1732. These studies performed by the laboratory are maintained on file for review.

13.2 Demonstration of Capability Studies

Refer to Alpha SOP 1739 for further information regarding IDC/DOC Generation.

13.2.1 Initial (IDC)

The analyst must make an initial, one-time, demonstration of the ability to generate acceptable accuracy and precision with this method, prior to the processing of any samples.

13.2.2 Continuing (DOC)

The analyst must make a continuing, annual, demonstration of the ability to generate acceptable accuracy and precision with this method.

Alpha Analytical, Inc.
Facility: Mansfield
Department: GC/MS-Semivolatiles

ID No.:**2164**Revision 15
Published Date:3/4/2019 10:48:12 AM
Page 18 of 18

14. Pollution Prevention and Waste Management

Refer to Alpha's Chemical Hygiene Plan and Waste Management and Disposal SOP for further pollution prevention and waste management information.

15. Referenced Documents

Chemical Hygiene Plan

Title: 1,4-Dioxane By GC/MS-SIM

SOP 1732 Detection Limit (DL), Limit of Detection (LOD) & Limit of Quantitation (LOQ)

SOP 1739 Demonstration of Capability (DOC) Generation

SOP 1731 Manual Integration & Compound Rejection

SOP 1797 Hazardous Waste and Sample Disposal

SOP 1816 Reagent, Solvent and Standard Control

16. Attachments

None

Attachment 3

Field Forms

Wittman GeoSciences, PLLC 3636 N. Buffalo Road, Orchard Park; NY 14127 michelewittmangeo@gmail.com 716-574-1513			Y 14127	Boring No:		
Project Name & Location WGS Project Number: Start Date GW Depth While Drilling GW Depth at Completion		ober: Drilling	En	WGS Representative: WGS Reviewed & Approved by: M. Wittman, P.G. Drilling Contractor Type of Drill Rig Sampler Type:	i <u>. </u>	
Sample Depth (ft)	Sample No.	Sample Depth (feet)	Recovery (%)	SAMPLE DESCRIPTION	OVM Reading (ppm)	
1						
3						
4						
5 6						
7						
8						
9						
11						
12 13						
14						
15						
16 18						
20						
22						
24		1) Organia vanar mat	or wood to field o			
N	otes:	ND - non detect or		screen and headspace soil samples.		
	1) Stratification lines represent approximate boundary between soil. Transitions may be gradual. Depths are approximate. 2) Groundwater (GW) depths approximate at time of sampling. Fluctuations in groundwater may occur. 3) f=fine; m=medium; c=coarse 4) and (36-50%); some (21-35%); little (11-20%); trace (1-10%) MC - Geoprobe Macrocore SS - Split Spoon SH - Shelby Tube BC - Bedrock Core					



Date:	Project No.:	3752 N. Buffalo Rd.
Client:		Orchard Park, NY 14127
Project:		P (716) 667-3130
Site:		F (716) 667-3156
Weather:		
	FIELD INVESTIGATION	REPORT
(Start typing he	ere making sure underline is on and text is	s justified. Hit tab at the end of the very last
row to extend the under	rline to the right margin).	
Signature	Title	

Well Data Sheet

Date:	•			Job #:			
Crew:							
Well Depti	า:						
Initial Phas	se Level:		:				
Initial Wate	er Level:						
*							
Volume Ca	alculation:						
DTB-DTW	*	=1-well vo	ol				
			Purge	Record			
	Time	Volume	рН	Cond.	Temp.	Turbidity	
-							
						•	
Purge Meth	nod:	Bailer/Sub	omersible	Pump			
Initial Wate	r Quality						
Final Wate	r Quality						
¥							
			SAMPL	E RECORD			
			1				
Date:	<u></u>		_ :	Volume:	·		
Time:			<u>.</u>	Analysis:			
Crew:		· · · · · · · · · · · · · · · · · · ·		Chain of C	Chain of Custody #:		
Method:	<u></u>		-	Sample Ty	pe:		
Sample ID:			_			•	
Water Qua	lity:		_	Diameter	Multiply by		
pH:		and the second	_	1"	0.041		
Conductivit			_	2"	0.163		
Temperatu	re:		- ··	3"	0.367		
Turbidity:	-	·	_	4"	0.653		
				6"	1.468		
	•			8"	2.61		
Comments							
			Signatu	re:			



NEW YORK CHAIN OF CUSTODY Service Centers Mahwah, NJ 07430: Albany, NY 12205: Tonawanda, NY 144

Mahwah, NJ 07430: 35 Whitney Rd, Suite 5 Albany, NY 12205: 14 Walker Way Tonawanda, NY 14150: 275 Cooper Ave, Suite 1

Page	D-4-D-4
of	in Lah

ALPHA Job#

	CUSTODY	Tonawanda, NY 14150: 275 Coo		05	0	···	9. 9 Hz	in L	ab						ALPHA JOB #	
Westborough, MA 01581 8 Walkup Dr.	Mansfield, MA 02048 320 Forbes Blvd	Project Information			Deliverables						Billing Information					
TEL: 508-898-9220 FAX: 508-898-9193	TEL: 508-822-9300	Project Name:						ASP-A	4	-		ASP	-B		Same as Client Info	
	FAX: 508-822-3288	Project Location:	entry (Maena)				EQuIS (1 File)			EQuIS (4 File)			ile)	PO #		
Client Information		Project #						Other								
Client:		(Use Project name as Pro	oject#)				Regu	latory l	Requir	emen	it				Disposal Site Information	
Address:		Project Manager:						NY TO	GS			NY P	art 375		Please identify below location of	
		ALPHAQuote #:						AWQ S	Standar	ds		NY C	P-51		applicable disposal facilities.	
Phone:		Turn-Around Time						NY Res	stricted	Use		Other			Disposal Facility:	
Fax:		Standard		Due Date:				NY Uni	restricte	ed Use					☐ NJ ☐ NY	
Email:		Rush (only if pre approved)		# of Days:				NYC S	ewer Di	ischarç	ge				Other:	
These samples have b	peen previously analyz	ed by Alpha					ANAL	YSIS						-	Sample Filtration	T
Other project specifi	c requirements/comn	nents:													Done	t
										ľ					Lab to do	а
	8														<i>Preservation</i> ☐Lab to do	
Please specify Metal	s or TAL.								ĺ							В
															(Please Specify below)	t
ALPHA Lab ID		imple ID Colle		lection Sample Sample		Sampler's			ļ							t
(Lab Use Only)	0.0	impie ib	Date	Time	Matrix	Initials			1			<u> </u>			Sample Specific Comments	е
1 × 1 × 1																
							,									
alle to the state of the state of																
			i							1						
ATA UNIXABILITY															·	
Preservative Code:	Container Code P = Plastic	Westboro: Certification No	o: MA935	-	Con	ntainer Type									Please print clearly, legibly	,
A = None B = HCl	A = Amber Glass	Mansfield: Certification No	Mansfield: Certification No: MA015			itainei Type									and completely. Samples of	
C = HNO ₃	V = Vial G = Glass	Preservativ			Dunganyakiya										not be logged in and	4
D = H ₂ SO ₄ E = NaOH	B = Bacteria Cup				reservative									turnaround time clock will not start until any ambiguities are		
F = MeOH		Relinquished By: Date			Time F		Received By:				Date/Time				resolved. BY EXECUTING	
															THIS COC, THE CLIENT	
K/E = Zn Ac/NaOH	D = BOD Bottle									\neg					HAS READ AND AGREES TO BE BOUND BY ALPHA	
O = Other															TERMS & CONDITIONS.	_
Form No: 01-25 HC (rev. 3	30-Sept-2013)				.*			7	,						(See reverse side.)	

Soil Vapor Intrusion - Structure Sampling Building Questionnaire Site Name : Site No.: Date: Time: Structure Address: Preparer's Name & Affiliation : ___ Residential ? ☐ Yes ☐ No Owner Occupied ? ☐ Yes ☐ No Owner Interviewed ? \square Yes \square No Commercial ? ☐ Yes ☐ No Industrial ? ☐ Yes ☐ No Mixed Uses ? ☐ Yes ☐ No Identify all non-residential use(s): ____) _____ - ____ Owner Name : Owner Phone : Secondary Owner Phone:) _____ - ____ Owner Address (if different) : _____ ____ Occupant Phone : () _____-Occupant Name : ___ Secondary Occupant Phone: (Number & Age of All Persons Residing at this Location : _____ Additional Owner/Occupant Information : _ Describe Structure (style, number floors, size) : _____ Approximate Year Built : Is the building **Insulated**? ☐ No ☐ Yes Lowest level: ☐ Slab-on-grade ☐ Crawlspace ☐ Basement Describe Lowest Level (finishing, use, time spent in space) : _____ Floor Type:

Concrete Slab

Dirt

Mixed: Floor Condition: \Box Good (few or no cracks) \qed Average (some cracks) \qed Poor (broken concrete or dirt) Sumps/Drains? ☐ Yes ☐ No Describe: Identify other floor penetrations & details : ___ Wall Construction: ☐ Concrete Block ☐ Poured Concrete ☐ Laid-Up Stone Identify any wall penetrations : _____ Identify water, moisture, or seepage: location & severity (sump, cracks, stains, etc): ____ Heating Fuel: ☐ Oil ☐ Gas ☐ Wood ☐ Electric ☐ Other : _____ Heating System : ☐ Forced Air ☐ Hot Water ☐ Other : ___ Hot Water System : ☐ Electric ☐ Boilermate ☐ Other: _____ ☐ Combustion Clothes Dryer: Where is dryer **vented** to? ☐ Electric ☐ Gas If combustion occurs, describe where air is drawn from (cold air return, basement, external air, etc.): Fans & Vents (identify where fans/vents pull air from and where they vent/exhaust to):

Structure ID: ___

Structure ID : _____

Describe factors that may affect indoor air quality (chemical use/storage, unvented heaters, smoking, workshop):								
Attached garage ?	☐ No Air freshene	ers? 🗆 Yes 🗆	No					
New carpet or furniture ? ☐ Yes	☐ No What/Where	e?						
Recent painting or staining?	□ Yes □ No	Where ? :						
Any solvent or chemical-like odors ?	☐ Yes ☐ No	Describe :						
Last time Dry Cleaned fabrics brought	in ?	What / Where ?						
Do any building occupants use solvents	s at work?	☐ No Describe	:					
Any testing for Radon ? ☐ Yes	□ No Result	s:						
Radon System/Soil Vapor Intrusion Mitigation System present ?								
	Lowest Building Leve	l Layout Sketch						

- Identify and label the locations of all sub-slab, indoor air, and outdoor air samples on the layout sketch.
- Measure the distance of all sample locations from identifiable features, and include on the layout sketch.
- Identify room use (bedroom, living room, den, kitchen, etc.) on the layout sketch.
- Identify the locations of the following features on the layout sketch, using the appropriate symbols:

B or F	Boiler or Furnace	0	Other floor or wall penetrations (label appropriately)
HW	Hot Water Heater	XXXXXX	Perimeter Drains (draw inside or outside outer walls as appropriate)
FP	Fireplaces	######	Areas of broken-up concrete
ws	Wood Stoves	● SS-1	Location & label of sub-slab vapor samples
W/D	Washer / Dryer	● IA-1	Location & label of indoor air samples
S	Sumps	● OA-1	Location & label of outdoor air samples
@	Floor Drains	● PFET-1	Location and label of any pressure field test holes.

Page	of	

Structure Sampling - Product Inventory

Homeowner Name & Address:			Date:	
Samplers & Company:			Structure ID:	
Site Number & Name:			Phone Number:	
Make & Model of PID:			of PID Calibration:	
Identify any Changes fro	m Original	Building Questionnaire :		
Product Name/Description	Quantity	Chemical Ingredients	PID Reading	Location



AIR/VAPOR SAMPLING FIELD DATA SHEET

Client:	Project No.:					
Site Name & Address:						
Person(s) Performing Sampling:						
Sample Identification:	_					
Sample Type: ☐ Indoor Air (ambient)	□Outdoor Air □Soil Vapor	□Sub-slab Vapor				
Date of Collection:	Setup Time:	Stop Time:				
Sample Depth:	_					
Sample Height:	_					
Sampling Method(s) & Device(s):						
Purge Volume:	_					
Sample Volume:						
Sampling Canister Type & Size (if applic	cable):					
Canister #	Regulator #					
Vacuum Pressure of Canister P	rior to Sampling:					
Vacuum Pressure of Canister A	fter Sampling:					
Temperature in Sampling Zone:						
Apparent Moisture Content of Sampling	Zone:					
Soil Type in Sampling Zone:						
Standard Chain of Custody Procedures	Used for Handling & Delivery of	Samples to Laboratory:				
□Yes □No. If	no, provide reason(s) why?					
Laboratory Name:						
Analysis:						
Comments:						
Sampler's Signature	Г)ate:				

APPENDIX C HEALTH AND SAFETY PLAN

HEALTH AND SAFETY PLAN

BROWNFIELD CLEANUP PROGRAM For **OCTAVUS STORAGE NY, LLC** 3678 Southwestern Boulevard, Orchard Park, NY 14127 BCP # TBD



Prepared For: **OCTAVUS STORAGE NY, LLC** 1810 West Kennedy Boulevard, Tampa, FL 33606 HEI Project No: e1728

Prepared By:

Hazard Evaluations, Inc. 3636 North Buffalo Road Orchard Park, New York 14127 Orchard Park, New York 14127 (716) 667-3130

Wittman GeoSciences, PLLC 3636 North Buffalo Road 716-574-1513

September 30, 2018



TABLE OF CONTENTS

	Page
1.0 INTRODUCTION	1
2.0 SITE DESCRIPTION AND HISTORY	1
2.1 Site Description	1
2.2 Site History	
3.0 ASSIGNED RESPONSIBILITIES	
3.1 Environmental Health & Safety Manager	
3.2 Project Manager	
3.3 Site Safety Officer	
3.4 Site Workers	
4.0 TRAINING and SAFTETY MEETINGS	3
4.1 Training	
4.2 Safety Meetings	
5.0 PERSONAL PROTECTIVE EQUIPMENT	4
6.0 HAZARD ANALYSIS	
6.1 Chemical Hazards	
6.2 Physical/General Hazards	
6.3 Biological Hazards	
7.0 SITE MONITORING	
7.1 Soil Borings and Monitoring Wells	7
7.2 Action Levels	
7.3 Particulate Monitoring	8
8.0 COMMUNITY AIR MONITORING PLAN	
9.0 SITE ACTIVITY AREAS AND ACCESS CONTROL	
10.0 DECONTAMINATION PROCEDURES	9
10.1 Prevention of Contamination	9
10.2 Personal Decontamination	10
10.3 Decontamination During Medical Emergencies	10
10.4 Decontamination of Equipment	
10.5 Disposal of the Contaminated Materials	
11.0 EMERGENCY RESPONSE	11
11.1 Response Procedures	
11.2 Communications	
11.3 Evacuation	
11.4 Fire or Explosion	
11.5 Personal Injury	
11.6 Adverse Weather Conditions	
11.7 Traffic, Heavy Equipment & Machinery	
11.8 Utilities	
11.9 Emergency Contingency Plan	
11.9.1 Contamination Emergency	
11.9.2 Spill/Air Release	13
11.9.3 Unknown Drums or USTs	
11.10 Additional Safety Practices	
12.0 RECORDS AND REPORTING	13



LIST OF TABLES

Table 1 Hazard Characteristics of Potential Contaminants of Concern

Table 2 Emergency Contact List

LIST OF FIGURES

Figure 1 Locus Plan

Figure 2 Map/Direction to Hospital

LIST OF ATTACHMENTS

Attachment A Community Air Monitoring Plan



1.0 INTRODUCTION

This Health & Safety Plan (HASP) has been developed for the Remedial Investigation/Alternatives Analysis Report (RI) to be completed by Hazard Evaluations, Inc. (HEI) and Wittman GeoSciences, PLLC (WGS) for Octavus Storage NY, LLC (Octavus) at 3678 Southwestern Boulevard, Orchard Park, Erie County, New York as shown on Figure 1, on behalf of Octavus (Applicant) as part of the Brownfield Cleanup Program (BCP). The proposed work will include completion of soil borings, test pits, installation of monitoring wells, soil and groundwater sampling, soil excavation and sampling, and report preparation. Such activities mandate the performance of tasks with a potential to expose remediation workers to various environmental contaminants previously identified on-site, primarily involving historical industrial fill, potentially including semi-volatile organic compounds (SVOCs) and metals. Limited exposure potential may be related to commercial substances used for equipment decontamination. A general listing of the work tasks to be completed is as follows:

- Soil sampling using a direct push method (Geoprobe) and hollow stem auger equipment;
- Test pit excavation with a track mounted excavator and bucket;
- Soil sample collection and analysis;
- Monitoring well installation, purging and development; and
- Groundwater sampling using disposable bailers, and analysis.

The intent of this HASP is to identify and present appropriate safety procedures to be followed by investigation/remediation workers involved with project activities throughout the performance of the RI. Such procedures are designed to reduce the risk of remediation worker exposure to the primary substances of concern.

The procedures also address several other physical hazards that may be encountered during the RI activities. Recommended safety procedures presented herein may be modified as the RI proceeds based upon conditions encountered at the site, with the mutual agreement of HEI, New York State Department of Environmental Conservation (NYSDEC), New York State Department of Health (NYSDOH), and Applicant. A copy of this HASP (including any modifications) will be maintained on-site throughout the duration of the RI field work to be used as a reference by HEI and their subcontractors. An initial safety meeting will be conducted at the site prior to the initiation of the sampling activities to inform all affected remediation workers of potential exposures and hazards.

2.0 SITE DESCRIPTION AND HISTORY

2.1 Site Description

The Site is addressed as 3678 Southwestern Boulevard located in Orchard Park, Erie County, New York. The Site consists of one parcel, totaling approximately 7.34 acres of land, and is bound to the north and east by Route 219 Expressway, to



the south by Southwestern Boulevard, and to the west by vacant land. The property is located within a commercial area of Orchard Park, utilized for light industrial, commercial and residential purposes.

2.2 Site History

The Site is improved with a two-story, approximate 4,756 square foot building. The parcel includes an asphalt driveway, concrete loading ramp, perimeter fence around portions, and wooded area in the northern portion which contains a former industrial process pond.

The subject site was initially utilized agriculturally, from at least 1926 through at least 1966. In 1967, the property was purchased by the County of Erie (Torrax) to be utilized for the disposal of municipal solid waste through incineration. The facility was identified as under construction in 1971, and operated by the Erie County Torrax Facility from 1972 to 1979. The building was used for an experimental pilot program for the incineration of refuse. In 1974, the experimental program was halted and all of the major equipment was removed. During its operational phase, the building housed a high temperature incinerator capable of using 13,000 cu ft/hr natural gas. The maximum waste to be burned was 8,400 lb/hr. The northern portion of the site was used as an ash landfill and scrubber process water treatment lagoon.

The building was then occupied between 1974 and 1977 by the Erie County Sewerage and Management Department as a vehicle repair and/or storage facility, and has been vacant since that time. The property was purchased in 1993 by Torrax Building Joint Venture and remained unoccupied, but was used for occasional storage of equipment. The building is currently vacant and underutilized.

3.0 ASSIGNED RESPONSIBILITIES

Specific safety responsibilities have been established for the performance of the RI as indicated below.

3.1 Environmental Health & Safety Manager

The Environmental Health & Safety Manager (EHSM) has the authority to commit any resources necessary to implement an effective RI safety program, thereby protecting the health of affected Site workers. The EHSM will delegate responsibilities, as necessary, to the Project Manager (PM) in order to facilitate various aspects of this HASP. The resolution of any on-site safety issues encountered during the RI will be coordinated by the EHSM.

3.2 Project Manager

The Project Manager (PM) will be responsible for the overall project, including implementation of the HASP. The PM will coordinate with the Site Safety Officer (SSO) to ensure that project goals of the project are met in a manner consistent with the HASP requirements.



3.3 Site Safety Officer

The Site Safety Officer (SSO) will be responsible for ensuring that the recommended safety procedures are followed during sampling activities. The SSO will supervise HEI employees and subcontractors throughout the duration of the field work. The SSO is knowledgeable of general construction safety practices and remediation worker protection techniques. Responsibilities will include:

- Ensuring day-to-day compliance with HASP safety procedures;
- Maintaining adequate personal protective equipment (PPE) supplies;
- Calibration and maintenance of monitoring instruments;
- Authority to stop work activities at any time if unsafe work conditions are identified;
- Implementing personnel decontamination procedures;
- Initiate emergency response procedures;
- Maintain a diary of activities with safety relevance; and
- Establishing and assuring adequate records of all:
 - Occupational injuries and illnesses;
 - Accident investigations;
 - Reports to insurance carrier or state compensation agencies;
 - Records and reports required by local, state and/or federal agencies;
 and
 - Property or equipment damage.

3.4 Site Workers

Affected site workers will include HEI employees and subcontractor employees. Site workers must comply with all aspects of the HASP and its safety procedures. Personnel entering the Site will have completed training requirements for hazardous waste site operations in accordance with Occupational Safety and Health Administration (OSHA) 29CFR 1910.120(c); 29CFR 1910.146(d); and 29CFR 1910.147(c). Site workers and SSO must have completed appropriate medical surveillance as required by OSHA 29CFR 1910.120(f).

3.5 Subcontractors

Various subcontractors will be utilized on the Site during RI activities, such as driller and excavation contractor. Subcontractors are responsible for development of their own HASP that is at least as stringent as this HASP. A copy of this HASP will be provided to the subcontractors for informational purposes. Subcontractors will be informed of potential health and safety hazards, as well as environmental monitoring data collected during field activities.

4.0 TRAINING AND SAFTETY MEETINGS

4.1 Training

Site personnel assigned to the Site will be in compliance with the training requirements of 29 CFR 1910 and 1926 as listed below. Site personnel will have met one of the following requirements prior to the start of on-site activities.



- A 40-hour minimum hazardous materials safety and health course, as stipulated in 29 CFR 1926.65 (e)(3); and
- An eight-hour minimum refresher course per year after the 40-hour minimum training has occurred (29 CFR 1926.65 (e)(8)).

On-site managers and supervisors must be in compliance with the additional supervisory training requirements of 29 CFR 1926.65 (e)(4). Emergency responders must be in compliance with the additional training requirements of 29 CFR 1926.65 (e)(7). Appropriate certificates of participating site workers' training programs will be maintained at HEI offices.

4.2 <u>Safety Meetings</u>

Site workers and subcontractors will be familiar with the Site and facility layout, have an understanding of known and potential hazards, and details within this HASP. On-site safety meetings will occur daily, or as needed to assist site workers and subcontractors in conducting activities safely. Attending personnel must sign an attendance sheet. Site workers must attend a safety meeting prior to being allowed to work on-site.

5.0 PERSONAL PROTECTIVE EQUIPMENT

An important aspect for site worker safety is correct selection of personal protective equipment (PPE). The levels of protection listed below are based on 29 CFR 1910.120. The majority of site activities will be conducted in Level D protection. This level of protection was selected based on the types and measured concentrations of the hazardous substances in the samples previously collected and their associated hazards and/or toxicity; and potential or measured exposure to substances in air, splashes of liquids or other indirect contact with material due to the task being performed.

Level D will generally consist of the following:

- Coveralls or long pants and long sleeve shirt to provide protection from dermal contact with soil;
- High visibility safety vest;
- Steel toe work boots;
- Safety glasses;
- Hard hat; and
- Chemical-resistant gloves.

Additional equipment may be donned at SSO requirement, including disposable boots, hearing protection, safety vest, or disposable outer chemical coveralls (Tyvek suits).



Level C will generally consist of the following:

- Full or half face air purifying respirator (APR) equipped with appropriate organic vapor canisters and/or other chemical cartridges;
- Chemical resistant clothing, such as Tyvek suit. Suits will be one piece with booties, hood, and elastic wristbands;
- High visibility safety vest (disposable);
- Outer chemical-resistant gloves (i.e. nitrile or neoprene) and inner latex gloves;
- Steel toe work boots: and
- Hard hat.

Level B will generally consist of the following:

- Self-contained breathing apparatus (SCBA) in a pressure demand mode, or supplied air with escape SCBA;
- Chemical resistant closing, such as Tyvek suit. Suits will be one piece with booties, hood, and elastic wristbands;
- High visibility safety vest (disposable);
- Outer chemical-resistant gloves (i.e. nitrile or neoprene) and inner latex gloves;
- Chemical resistant tape over PPE as needed (i.e. at glove/Tyvek location);
- Steel toe work boots; and
- Hard hat.

6.0 HAZARD ANALYSIS

Many hazards are associated with environmental work on a site. The hazards listed below deal specifically with those associated with the management of potentially contaminated soil, air, and groundwater including, physical hazards as well as environmental hazards.

6.1 Chemical Hazards

The primary chemical hazardous substances known or suspected at the subject site is semi-volatile organic compounds (SVOCs) and metals that are present within the historical industrial fill due to former industrial operations. Additional contaminants that may be present include volatile organic compounds (VOCs) associated with past petroleum storage as well past industrial usage. A summary of hazards associated with these chemicals is included on Table 1. The list has been developed based on planned activities and potential Site conditions. The most likely routes of chemical exposure during site work includes skin absorption and inhalation of airborne dust particles. The information was used to develop the levels of PPE to be used during the duration of RI field work on-site.

6.2 Physical/General Hazards

Based on the proposed scope of work to be completed, the following potential physical hazards have been identified:



- Slip/Trip/Fall Due to the timing of the project, some areas may have icy surfaces that will increase the possibility of accidental falls. Additionally, good housekeeping practices such as cleaning up garbage and stored materials from the work area are essential to reduce the occurrence of trips and falls.
- Vehicle and machinery in motion hazards A drill rig will be utilized for soil sample collection. To minimize potential hazards, the drilling subcontractor will be responsible for health and safety of its personnel, equipment and operations. Utilities must be called in via Dig Safely New York and/or Site owner. Cones and flags will be set up around each work area, as necessary. Workers must be aware of pinch points when setting the rig and lowering mast/pull rods. PPE must be worn to prevent eye injury. All body parts, clothing, and manual tools must be kept three to five feet from moving equipment when possible. Gloves and PPE must be worn when working with rods and cleaning equipment. Monitoring of the breathing zone will be completed as necessary to ensure vapors are below action levels. Each worker must have an awareness of muscle strain. All sampling liners must be opened in a motion away from the body and hands. The rig cannot be moved with the mast in a raised position.
- Electrical Heavy equipment (e.g., excavator, backhoe, drill rig) shall not be operated within ten feet of high voltage lines. Working near wet areas should also be taken into consideration when working with electrical equipment. Surge protectors and ground fault protectors must be used in such conditions.
- Noise Heavy machinery creates excessive and loud noise levels. Overexposure can result in hearing damage or loss. Proper hearing protection shall be worn during exposure to noise from heavy equipment.
- Underground utilities The proper utility clearance will be obtained before conducting any digging or drilling operations.
- Excavation and soil sampling through use of heavy equipment Excavations that are greater than four feet in depth require a protective system prior to entry into the excavation. The Project Manager will be responsible for determining if the excavation requires safety shoring. Personnel will not be permitted to work under suspended or raised loads, and shall always wear highly visible clothing. PPE, including steel-toed boots, safety glasses, and hard hats must be worn. Personnel should not walk directly in back of, or to the side of, heavy equipment without the operator's knowledge. Engineering controls can be implemented, such as water for dust particulate control.
- Cold Stress Site work is scheduled during the winter and early spring months and therefore cold weather may present hazards. Frostbite and hypothermia can occur quickly and the signs and symptoms of such should be known. Signs of hypothermia include slurred speech, confusion, and an overall warm sensation. Frostbite can be identified by red/frozen skin, numbness, and lack



of sensation on the skin. In each case, the victim should be moved to a warm place. With frostbite, the affected area should be placed in warm water and wrapped with a warm towel. Medical attention is necessary after initial treatment.

- Heat stress Although not anticipated due to the time of year operations will occur, heat stress is a severe hazard that can result in heat fatigue or even heat stroke. Signs and symptoms of heat stroke include red, dry, and hot skin as well as confusion, a rapid pulse, and nausea. Adequate shade and drinking liquids should be provided to personnel working in hot weather conditions. If a person is suspected to be suffering from heat fatigue or stroke, transport to a cool place and place cold compresses on the neck and armpits. Call 911 immediately.
- Weather (i.e. lightning storm) On-site personnel shall cease operation at the first sign of a thunderstorm/lightning strike. Workers should seek shelter within a permanent building and stay away from tall structures, trees, telephone poles, and drill rigs/equipment.

6.3 Biological Hazards

Biological hazards can be caused by contact with land animals, birds, insects, and plants. Irritation, illness, and in extreme cases, permanent disability or death can occur. The Site is located in an urban area within the City of Buffalo and field work will occur in winter/early spring. Rodents are considered the most likely biological hazards at this Site. Contact with rodents, more specifically rats, shall be avoided. If bitten or scratched by any type of rodent or fur-bearing animal, medical treatment should be sought immediately. Insect bites and stings are not considered a serious threat due to the time of year that work is planned to take place. Insect bites and stings can cause irritation and transmit disease. If stung by an insect, apply cold water and soap and immediately apply a cold compress to the area to limit swelling. If the victim is allergic to such bite or sting, immediate medical care may be necessary.

7.0 SITE MONITORING

Air monitoring will be performed on-site, using a Photoionization Detector (PID) equipped with a 10.6 electron volt (eV) lamp, in order to track contamination levels. By knowing these levels, safety is insured for personnel working on-site.

7.1 Soil Borings and Monitoring Wells

On-site monitoring will be completed by the SSO or Site worker assigned to oversee drilling operations, soil sampling and monitoring well installation/sampling. The PID will be utilized to monitor the breathing zone, the borehole, and subsurface samples for the presence of volatile organic compounds (VOCs). Auger spoils will also be monitored. Fluids produced from monitoring well development and sampling will also be monitored with the PID.



7.2 Action Levels

Work area ambient air monitoring for VOCs will be completed within the breathing zone periodically. Action levels will be based on the PID readings. The action level assumes that background level of organics is close to non-detect. Background VOC readings will be recorded daily. Action levels are listed below.

Sustained PID Reading	Action	Minimum Respiratory Protection
0 to 10 ppm*	None	None – Level D
10 to 25 ppm	Monitor for 15 minutes; if concentration does not decrease to under 10 ppm, upgrade PPE; consider venting area	respirator with organic vapor
>25 ppm	Monitor for 15 minutes; Consider venting area, upgrade PPE	Suspend work or supplied-air full face respirator – Level B

^{*} parts per million

7.3 Particulate Monitoring

Monitoring for particulates will be completed periodically in the Site worker breathing zone. The decision to upgrade levels of PPE will be made in conjunction with consideration for weather conditions, wind conditions and anticipated duration of field activity. Background particulate concentrations will be measured and recorded on a daily basis.

8.0 COMMUNITY AIR MONITORING PLAN

A Community Air Monitoring Plan (CAMP) requires monitoring of VOCs and particulates at downwind locations and is intended to provide a level of protection for neighboring residences and businesses. Continuous monitoring will be conducted during ground intrusive activities. The completed CAMP is included in Attachment A.

9.0 SITE ACTIVITY AREAS AND ACCESS CONTROL

Prior to the initiation of the RI, three work zones will be established to facilitate the implementation of the HASP. Prior to commencement of field work, a further definition of these zones will be established. Guidelines for establishing work areas is as follows.

 Exclusion Zone (EZ) – Primary exclusion zones will be established around each intrusive field activity, such as soil boring or excavation area. Locations will be identified by the placement of orange cones. Site workers in these areas must wear appropriate PPE. Upon leaving the EZ, if PPE becomes



contaminated, Site workers must remove and dispose of gloves and any other disposable PPE. After removing the PPE, Site workers should thoroughly wash their hands. Access to the EZ will be limited to Site workers only for both safety and data integrity purposes.

- Contamination Reduction Zone (CRZ) A CRZ will be established between the EZ and property limit, and provides an area for decontamination of Site equipment. The specific location of this pad will be field determined, but will be out of the way of Site activities and sampling activities. Portable wash stations will be set up in the CRZ and will consist of a potable water supply, hand soap and disposable towels. An Alconox solution will be available to decontaminate equipment used in the sampling locations. The SSO will monitor equipment cleaning procedures to ensure their effectiveness. Equipment will be adequately cleaned and Site workers will remove contaminated PPE prior to either entering the Support Zone or leaving the Site for the day once sampling activities have been completed. A fire extinguisher and first aid kit will be located in this area.
- Support Zone (SZ) The SZ is considered to be clean, and PPE are not required. The SZ will be an area on-site adjacent to the CRZ in which supplies or equipment are stored and maintained. PPE is donned in the SZ prior to entering the CRZ.

10.0 DECONTAMINATION PROCEDURES

Decontamination procedures for personnel and equipment will be implemented when exiting the work area. Decontamination involves physically removing contaminants and generally includes the removal of contamination, avoiding spreading contamination from the work zone, and avoiding exposure of unprotected personnel outside of the work zone to contaminants.

10.1 Prevention of Contamination

The first step in decontamination is to establish standard operating procedures that minimize contact with hazardous substances, and thereby the potential for contamination. Site workers should be aware of the importance of minimizing contact with hazardous substances and the use of appropriate practices and procedures for Site operations. HEI utilizes this approach by ensuring Site workers:

- Stress work practices that minimize contact with hazardous substances (e.g., do not walk through areas of obvious contamination, do not directly touch potentially hazardous substances, etc.);
- Protect sampling instruments from gross contamination by bagging and making openings in the bag for sample ports and sensors that contact site materials; and



• Wear disposable outer garments and use disposable equipment where appropriate.

10.2 <u>Personal Decontamination</u>

The degree of contamination exposure is a function of both a particular task and the physical environment in which it takes place. The following decontamination procedures will remain flexible, thereby allowing the decontamination crew to respond appropriately to changing conditions at the Site. It is expected that Site workers will be exposed to soil/fill potentially contaminated with SVOCs, metals, PCBs, and petroleum compounds. On-site sampling activities will be carried out in such a manner as to avoid gross contamination of Site workers, personal protective equipment, machinery and equipment.

Between sampling locations (or sometimes between samples at one sampling location), and upon the completion of the daily field activities, Site workers will proceed to the CRZ. Equipment (e.g., sampling tubes, shovels, tools, etc.) will be decontaminated in this area. Prior to leaving the Site for breaks, at the end of the work shift, or when PPE has been grossly contaminated, disposable boot covers, gloves, and suits will be removed and placed in a drum designated for the disposal of these materials. After removing PPE, each Site worker will wash with soap and fresh water prior to donning new PPE or leaving the site for the day. All wash-water and rinse-water will be collected and disposed of in accordance with appropriate regulations.

10.3 <u>Decontamination during Medical Emergencies</u>

In the event of a minor, non-life-threatening injury or medical problem, Site workers should follow the decontamination procedures as defined above and then administer first aid. If prompt, live-saving first aid is required, decontamination procedures should be omitted and immediate first aid should be administered, unless the environmental conditions are considered Immediately Dangerous to Life or Health (IDLH). In this case, the victim should be moved to a clean area and life-saving care should be instituted immediately without considering decontamination.

Outside garments can be removed (depending on the weather) if they do not cause delays, interfere with treatment or aggravate the problem. Respirators and backpacks must always be removed. Chemical-resistant clothing can be cut away. If the outer contaminated garments cannot be safely removed, the individual should be wrapped in plastic, rubber or blankets to help prevent contaminating the insides of ambulances and medical personnel. Outside garments will then be removed at the medical facility. No attempt should be made to wash or rinse the victim at the site. One exception would be if an individual is known to have been contaminated with an extremely toxic or corrosive material which could cause severe injury or loss of life.

10.4 <u>Decontamination of Equipment</u>

Decontamination efforts will be conducted in the CRZ. Gross contamination will first be removed with plastic scrapers or other appropriate tools. The equipment will be decontaminated at a temporary equipment decontamination pad in the CRZ via



hand washing or pressure washing. Alconox and water will then be used to wash the equipment with a cleaning brush. The equipment will then be rinsed with deionized water and allowed to air dry for a sufficient time prior to reuse or removal from the Site. Downhole tools and augers can be hand washed or pressure washed.

The decontamination of the direct push drilling rig will be undertaken (if necessary) when all on-site activities have been completed. Initially, scraping of the equipment will remove heavily caked materials prior to washing. Washing will then be accomplished with Alconox and water or pressure washing. Water generated during decontamination activities will be collected, stored and profiled for future off-site disposal.

10.5 <u>Disposal of the Contaminated Materials</u>

Potentially contaminated materials (i.e., gloves, clothing, sample sleeves etc.) will be bagged and segregated for proper disposal. Investigation derived waste will be managed in accordance with NYSDEC guidance. All fluids collected during groundwater sampling will be containerized and managed appropriately subsequent to field activities.

11.0 EMERGENCY RESPONSE

In the event of an emergency, the SSO will coordinate on-site emergency response activities. Appropriate authorities will be immediately notified of the nature and extent of the emergency. The emergency contact list is included on Table 2. The route and directions to the hospital are included as Figure 2.

11.1 Response Procedures

In the event of an emergency or acute exposure symptom, remediation workers will signal distress to the SSO. The SSO will be responsible for the response to emergencies and must:

- Have available a summary of the project's associated risk potential so it can be provided to authorities or response personnel in the event of an emergency;
- Maintain an Emergency Contact List (Table 2) and post it in a visible location;
- Maintain a map detailing directions to the nearest hospital (Figure 2); and
- Ensure appropriate safety equipment is available at the site.

11.2 Communications

Cell phones will be the primary means of communicating with emergency support services/facilities.



11.3 Evacuation

In the event of an emergency situation, such as fire, explosion, etc., all personnel will evacuate and assemble in a designated assembly area. The SSO will contact outside services (i.e. police, fire, etc.) as required. Under no circumstances will personnel be allowed to re-enter the area once the emergency signal has been given. The SSO must see that emergency equipment is available and emergency personnel are notified.

11.4 Fire or Explosion

In the event of a fire or explosion, the SSO will immediately evacuate the Site. The Buffalo Fire Department will then be notified immediately, and advised of the situation and the identification of any hazardous materials involved.

11.5 **Personal Injury**

Only basic emergency first aid will be applied on-site as deemed necessary. The SSO will supply available chemical specific information to appropriate medical personnel, as requested. First Aid kits supplied by HEI and its subcontractors will conform to Red Cross and other applicable good health standards, and will consist of a weatherproof container with individually sealed packages for each type of item. First Aid kits will be fully equipped before being sent to the Site.

11.6 Adverse Weather Conditions

In the event of adverse weather conditions, the SSO will determine if work can continue without sacrificing the safety of remediation workers. Some of the items to be considered prior to determining if work should continue are the potential for heat stress, inclement weather-related working conditions (heavy snow) and the operation of field instruments.

11.7 <u>Traffic, Heavy Equipment & Machinery</u>

Site workers must remain aware of the heavy equipment and machinery being used during RI activities. Site workers will be required to wear a high visibility safety vest during on-site work activities.

11.8 Utilities

Prior to the beginning of Site activities, all available drawings of the facility will be examined to determine the presence of underground or sub-slab utilities. HEI anticipates that a magnetic pipe and cable locator will be effective in the prevention of encountering underground utilities.

11.9 <u>Emergency Contingency Plan</u>

In the case of a spill emergency (e.g., tank/drum release, spill, fire, etc.), this section will describe the procedures to be followed during the event.

11.9.1 Contamination Emergency

It is unlikely that a contamination emergency will occur; however, if such an emergency does occur, the specific work area shall be shut down and



immediately secured. The area in which the contamination occurred shall not be entered until the arrival of trained personnel who are properly equipped with the appropriate PPE and monitoring instrumentation.

11.9.2 Spill/Air Release

In the event of a spill or air release of hazardous materials on-site, the specific area of the spill or release shall be shut down and immediately secured. The area in which the spill or release occurred shall not be entered until the cause can be determined and Site safety can be evaluated. The NYSDEC Spill Response unit shall be notified immediately. The spilled material shall be immediately contained.

11.9.3 Unknown Drums or USTs

In the event that unidentified containerized substances, including USTs, are discovered during soil sampling or soil excavation, work will be ceased immediately until hazards are addressed. The SSO will then visually assess the situation and identify any leaks or releases from the container. If leaking is identified, the spilled material shall be immediately contained. Upon visual assessment of releases and safety, properly trained personnel will then sample and remove/dispose of the waste/container.

11.10 Additional Safety Practices

The following are important safety precautions and practices that will be enforced during the field activities.

- Eating, drinking, smoking, chewing gum or tobacco or any activity that increases the probability of hand-to mouth transfer and ingestion of hazardous substances is prohibited during the RI activities;
- Remediation worker hands and face must be thoroughly washed before leaving the CRZ or before eating, drinking or other activity;
- Contact with potentially contaminated surfaces should be avoided whenever possible;
- The number of remediation workers and the amount of equipment should be minimized; and
- Alcoholic beverages will not be consumed during work hours by Site personnel.
 Personnel using prescription drugs may be limited in performing specific tasks (i.e. operating heavy equipment) without written authorization from a physician.

12.0 RECORDS AND REPORTING

The SSO will be responsible for establishing and maintaining adequate records of activities which take place at the Site. The records will pertain to Site workers involved in the project, regardless of their employer, as well as any agency personnel. A basic list of the information to be maintained is as follows:

- Occupational injuries or illnesses;
- Accident investigations;
- Reports to insurance carrier or State Compensation agencies;



- Records and reports required by local, state and federal agencies;
- Property or equipment damage;
- Third party injury or damage claims;
- Environmental testing logs;
- Explosive and hazardous substances inventories and records;
- Records of inspections and citations;
- Related correspondence; and
- Safety training level.



Tables

Table 1
Hazard Characteristics of Potential Contaminants of Concern

Contaminant	Potentially Impacted Media	Carcinogenicity/Symptoms of Acute Exposure	Occupational Exposure Values* ACGIH TLV OSHA PEL NIOSH IDLH
Benzene	Soil, Groundwater	Confirmed human carcinogen. Symptoms include irritation to eyes, skin, nose, respiratory system; headache; nausea; giddiness, fatigue.	PEL - 10 ppm; IDLH - 500 ppm; TLV - 0.5 ppm; STEL - 2.5 ppm
Chlorinated Organic Compounds	Soil, Groundwater	Exposure to the vapors of many chlorinated organic compounds such as vinyl chloride, tetrachloroethylene, 1,1,1-trichloroethane, trichloroethylene and 1,2-dichloroethylene and other chlorinated hydrocarbons may result in various symptoms including irritation of the eyes, nose and throat, drowsiness, dizziness, headache, blurred vision, uncoordination, mental confusion, flushed skin, tremors, nausea, vomiting, fatigue and cardiac arrhythmia. The liquid if splashed in the eyes, may cause burning irritation and damage. Repeated or prolonged skin contact with the liquid may cause dermatitis. Some of these compounds are considered to be potential human car-cinogens.	Refer to 29 CFR 1910.1017 for exposure values
Toluene	Soil, Groundwater	Insufficient data from carcinogenic studies to classify substance as a potential carcinogen. Symptoms include irritation to eyes, nose; fatigue; weakness; euphoria; headache; lacrimation.	PEL - 10 ppm; IDLH - 500 ppm; TLV - 20 ppm; STEL - 150 ppm
Ethyl Benzene	Soil, Groundwater	Confirmed animal carcinogen with unknown relevance to humans. Symptoms include irritation to eyes, skin, mucous membranes; headache; narcosis.	PEL - 5 ppm; IDLH - 800 ppm; TLV - 20 ppm; STEL - 30 ppm
o-, m-, and p-Xylenes	Soil, Groundwater	Insufficient data from carcinogenic studies to classify substance as a potential carcinogen. Symptoms include irritation to eyes, nose, throat; dizziness; excitement; drowsiness; nausea; vomiting.	PEL - 100 ppm; IDLH - 900 ppm; TLV - 100 ppm; STEL - 150 ppm
Polynuclear Aromatic Hydrocarbons (PAH's)	Soil, Groundwater	Many PAH's found in fuel oil and coal tar pitch volatiles (creosote) are confirmed human carcinogens. Symptoms include dermatitis and bronchitis.	Some PAH's have no established exposure values. Others considered coal tar pitch volatiles have an ACGIH TLV and OSHA PEL value of 0.2 mg/m ³ .
Cadmium	Soil	Suspected human carcinogen. Symptoms include pulmonary edema; difficulty breathing; cough; tightness in chest; substernal pain; headache; chills; nausea; vomiting; diarrhea; asnosmia.	PEL - 0.2 mg/m3; IDLH - 50 mg/m3; TLV - 0.01 mg/m3 (these limits are expressed for Cd dust)
Chromium	Soil	Hexavalent chromium compounds are confirmed human carcinogens. Symptoms include irritation to the respiratory system; nasal septum perforation; sensitization dermatitis (hexavalents). Irritation to the eyes; sensitization dermatitis (trivalents).	PEL - 0.5 mg/m3; IDLH - 250 mg/m3; TLV - mg/m3 (insoluable)
Lead	Soil	Confirmed animal carcinogen with unknown relevance to humans. Symptoms include weakness; tremor; irritation to eye; constipation; abdominal pain.	PEL - 0.05 mg/m3; IDLH - 100 mg/m3; TLV - 0.5 mg/m3
Mercury	Soil	Insufficient data from carcinogenic studies to classify substance as a potential carcinogen. Symptoms include irritation to eyes, skin; cough; chest pain; difficulty breathing; irritability; indecision; headache; fatigue; weakness; salivation.	PEL - 0.025 mg/m3 (acceptable ceiling concentration); IDLH - 2 mg/m3; TLV - 0.025 mg/m3 (elemental/inorganic)
Polychlorinated Biphenyl (PCBs)	Soil	Confirmed human carcinogen. Symptoms include dermal and ocular lesions, irregular menstrual cycles and a lowered immune response. Other symptoms included fatigue, headache, cough, and unusual skin sores	PEL - 1 mg/m3; IDLH - 5 mg/m3; TLV - 1 mg/m3

 $A CGIH\ TLV-American\ Conference\ of\ Governmental\ Industrial\ Hygienists\ Threshold\ Limit\ Value;\ Concentrations\ in\ ppm\ of\ mg/m3\ based\ on\ an\ 8-hour\ TWA$

OSHA PEL – Occupational Safety and Health Admiration Permissible Exposure Limits; Concentrations are shown in parts per million (ppm) or milligrams per cubic meter (mg/m3) based on an 8-hour time weighted average (TWA)

NIOSH IDLH - National Institute for Occupational Safety and Health Immediately Dangerous to Life or Health; Concentrations in ppm or mg/m3


OSHA STEL - Short Term Exposure Limit

Table 2
Emergency Contacts

Agency	Contact	Phone Number
Orchard Park Police	Emergency	911
Orchard Park Fire/First Aid	Emergency	911
Ambulance	Emergency	911
Poison Control Center	-	
Hospital	Mercy Hospital of Buffalo 565 Abbott Road Buffalo, NY 14220	(716) 826-7000
NYSDOH	TBD Empire State Plaza Corning Tower Rm. 1787 Albany, NY 12237	(518) 402-7860
NYSDEC	TBD 270 Michigan Ave. Buffalo, NY 14203	(716) 851-7220
NYSDEC	SPILL Hotline	(800) 457-7362
Hazard Evaluations	C. Mark Hanna, CHMM 3736 N. Buffalo Rd. Orchard Park, NY 14127	Office: (716) 667-3130 Cell: (716) 574-1513
Wittman GeoSciences, PLLC	Michele Wittman, PG 3736 N. Buffalo Rd. Orchard Park, NY 14127	Office: (716) 667-3130 Cell: (716) 574-1513
Carmina Wood Morris, DPC	R. Christopher Wood, PE 487 Main Street, Suite 500 Buffalo, NY 14203	(716) 842-3165
Octavus Storage NY, LLC	Joe McNeil 1810 W. Kennedy Tampa, Fl 33606	(813) 944-8836

Directions to Hospital - Mercy Hospital of Buffalo - Head west on Southwestern Boulevard (Hwy 20) toward Taylor Road; Turn right onto California Road; Turn left onto NY 179; Turn right at first cross street onto Abbott Road; Bear left to stay on Abbott, then next left onto Lorrain. Hospital is on left.

Figures



THIS DRAWING IS FOR ILLUSTRATIVE AND INFORMATIONAL PURPOSES ONLY AND WAS ADAPTED FROM USGS, BUFFALO SE, NEW YORK 2016 QUADRANGLE.

HAZARD EVALUATIONS, INC.		
Phase I/II Audits – Site Investigations – Facility Inspections		
SITE LOCATION		
3678 SOUTHWESTERN BOULEVARD		
ORCHARD PARK, NEW YORK		
OCTAVUS STORAGE NY, LLC		
TAMPA, FLORIDA		
DRAWN BY: SS	SCALE: NOT TO SCALE	PROJECT: e1728
CHECKED BY: MW	DATE: 09/18	FIGURE NO: 1

 \leftarrow

from 3678 Southwestern Blvd, Orchard Park, NY 14127 to Mercy Hospital of Buffalo, 565 Abbott Rd, Buffalo, ...

14 min (5.5 miles)







via Abbott Rd

Best route, lighter traffic than usual

3678 Southwestern Blvd

Orchard Park, NY 14127

† Head west on Hwy 20 W toward Taylor Rd

0.4 mi

Turn right onto California Rd

Turn left onto NY-179 W

0.1 mi

Turn right at the 1st cross street onto Abbott Rd
 Pass by Dairy Queen (Treat) (on the right in 3.2 ml)

4.1 ml -

Turn left to stay on Abbott Rd

0.2 ml

Turn left onto Lorraine Ave

233 A

Turn left


Destination will be on the right.

100 6

Mercy Hospital of Buffalo

565 Abbott Rd, Buffalo, NY 14220

These directions are for planning purposes only. You may find that construction projects, traffic, weather, or other events may cause conditions to differ from the map results, and you should plan your route accordingly. You must obey all signs or notices regarding your route.



HAZARD EVALUATIONS, INC.

Phase I/II Audits - Site Investigations - Facility Inspections

HOSPITAL DIRECTIONS 3678 SOUTHWESTERN BOULEVARD ORCHARD PARK, NEW YORK

DRAWN BY: MMW	SCALE: NOT TO SCALE	PROJECT: e1728
CHECKED BY: MMW	DATE: 09/18	FIGURE NO: 2

Attachment A Community Air Monitoring Plan

COMMUNITY AIR MONITORING PLAN

BROWNFIELDS CLEANUP PROGRAM For **OCTAVUS STORAGE NY, LLC** 3678 Southwestern Boulevard, Orchard Park, NY 14127 **BCP # TBD**



Prepared For: **OCTAVUS STORAGE NY, LLC** 1810 West Kennedy Boulevard, Tampa, FL 33606 HEI Project No: e1728

Prepared By:

Hazard Evaluations, Inc. 3636 North Buffalo Road Orchard Park, New York 14127 Orchard Park, New York 14127 (716) 667-3130

Wittman GeoSciences, PLLC 3636 North Buffalo Road 716-574-1513

September 30, 2018



TABLE OF CONTENTS

			<u>Page</u>
1.0	INTRODUC	TION	1
2.0	VOLATILE	ORGANIC COMPOUND AIR MONITORING	1
3.0	PARTICUL	ATE AIR MONITORING	2
4.0	DOCUMEN	TATION	3
5.0	WIND DIRE	CTION	3
LIST	OF ATTACH	HMENTS	
Attach	nment A	NYSDEC DER-10 Appendix 1A, New York State D Health, Generic Community Air Monitoring Plan	epartment of
Attach	nment B	NYSDEC DER-10 Appendix 1B, Fugitive Dust and Monitoring	Particulate



1.0 INTRODUCTION

This Community Air Monitoring Plan (CAMP) has been developed forRemedial Investigation/Alternatives Analysis Report (RI) to be completed by Hazard Evaluations, Inc. (HEI) and Wittman GeoSciences, PLLC (WGS) for Octavus Storage NY, LLC (Octavus) facility located at 3678 Southwestern Boulevard in Orchard Park, New York (Site), on behalf of Octavus (Applicant) as part of the Brownfield Cleanup Program (BCP).

The CAMP requires real-time monitoring of volatile organic compounds (VOCs) and particulates (dust) at downwind perimeter of each designated work area. The CAMP will be implemented during the excavation and removal of soils from the courtyard and vacant lot areas of the subject site. This CAMP will be completed in general accordance with NYSDEC DER-10 Appendix 1A, as included in Attachment A. Proposed locations include one upwind and one downwind location within the working area. Due to the large area and amount of sampling points at the site, HEI will determine monitoring points prior to intrusive activities. In addition, wind conditions will be observed during intrusive activities, which may influence the locations of the monitoring points.

2.0 VOLATILE ORGANIC COMPOUND AIR MONITORING

VOCs will be monitored at the downwind perimeter of the work are on a continuous basis and periodically during non-intrusive activities. VOC monitoring will be done using an organic vapor meter (OVM) equipped with a photoionization detector (PID) to provide real-time recordable air monitoring data.

VOCs will also be monitored and recorded at the downwind perimeter of the immediate work area(s). Upwind concentrations will be measured at the beginning of each day before activities begin and periodically throughout the day to establish background conditions. The downwind VOC monitoring device will also be checked periodically throughout the day to assess emissions and the need for corrective action. VOC monitoring action levels as per *DER-10 Technical Guidance for Site Investigations and Remediation* is as follows:

- If the ambient air concentration of total organic vapors at the downwind perimeter of the work area or exclusion zone exceeds 5 parts per million (ppm) above background for the 15-minute average, work activities must be temporarily halted and monitoring continued. If the total organic vapor level readily decreases (per instantaneous readings) below 5 ppm over background, work activities can resume with continued monitoring.
- If the organic vapor level at the perimeter of the work area persist at levels in excess of 5 ppm over background but less than 25 ppm, work activities must be halted, the source of vapors identified, corrective actions take to abate emissions,



and monitoring continued. After these steps, work activities can resume provided that the total organic vapor level 200 feet downwind of the exclusion zone or half the distance to the nearest potential receptor or residential/commercial structure, whichever is less; but in no case than that 20 feet, is below 5 ppm over background for the 15-minute average.

o If the organic vapor level is above 25 ppm at the perimeter of the work area, activities must be shut down.

3.0 PARTICULATE AIR MONITORING

The remediation crew will make all efforts to suppress dust and particulate matter during the handling of contaminated soil. Fugitive dust and particulate monitoring will be completed in accordance with DER-10 Appendix 1B, as included in Attachment B. The following techniques have been shown to be effective for the controlling the generation and migration of dust during construction activities:

- (a) Applying water on haul roads;
- (b) Wetting equipment and excavation faces;
- (c) Spraying water on buckets during excavation and dumping;
- (d) Hauling materials in properly tarped or watertight containers;
- (e) Restricting vehicle speeds to 10 mph;
- (f) Covering excavated areas and material after excavation activity ceases; and/or
- (g) Reducing the excavation size and/or number of excavations.

Care will be taken not to use excess water, which can result in unacceptably wet site conditions. Use of atomizing sprays will prevent overly wet conditions, conserve water and provide an effective means of suppressing fugitive dust.

Weather conditions will be evaluated during remedial work. When extreme wind conditions make dust control ineffective, as a last resort, remedial actions may need to be suspended.

Dust and particulate monitoring will be conducted near approximate upwind and downwind perimeters of the work area, when possible. If visual evidence of dust is apparent in other locations, monitoring equipment will be placed where necessary. Dust monitoring may be suspended during period of precipitation and snow cover.

Particulate air monitoring will be done with a DataRAM-4 (or similar), which will be capable of reading particles less than 10 micrometers in size (PM-10) and equipped with an audible alarm feature which will indicate exceedances. Dust monitoring devices will be recorded periodically throughout the day to assess emissions and the need for corrective actions. Particulate monitoring action levels as per *DER-10 Technical Guidance for Site Investigations and Remediation* is as follows:



- o If the downwind PM-10 particulate level is 100 micrograms per cubic meter (μg/m³) greater than background for the 15-minute period or if airborne dust is observed leaving the work area, then dust suppression techniques must be employed. Work may continue with dust suppression techniques provided that downwind PM-10 particulate levels do not exceed 150 (μg/m³) above the upwind level and provided that no visible dust is migrating from the work area.
- o If, after implementation of dust suppression techniques, downwind PM-10 particulate levels are greater than 150 (μg/m³) above the upwind level, work must be stopped and a re-evaluation of activities initiated. Work can resume provided that dust suppression measures and other controls are successful in reducing the downwind PM-10 particulate concentration to within 150 mcg/m³ of the upwind level and in preventing visible dust migration.

4.0 DOCUMENTATION

All 15-minute readings will be recorded and be available for or State (NYSDEC and NYSDOH) personnel to review. Instantaneous readings, if any, used for decision purposes should also be recorded.

5.0 WIND DIRECTION

Prevailing wind direction will be recorded at the beginning of each work day by visual observations of an on-site windsock. As wind direction may change throughout the work day, direction will be reestablished if a significant change in direction is observed. The wind direction results will be utilized to determine the placement of the monitoring equipment.



Attachment A

NYSDEC DER-10 Appendix 1A New York State Department of Health Generic Community Air Monitoring Plan

Appendix 1A New York State Department of Health Generic Community Air Monitoring Plan

Overview

A Community Air Monitoring Plan (CAMP) requires real-time monitoring for volatile organic compounds (VOCs) and particulates (i.e., dust) at the downwind perimeter of each designated work area when certain activities are in progress at contaminated sites. The CAMP is not intended for use in establishing action levels for worker respiratory protection. Rather, its intent is to provide a measure of protection for the downwind community (i.e., off-site receptors including residences and businesses and on-site workers not directly involved with the subject work activities) from potential airborne contaminant releases as a direct result of investigative and remedial work activities. The action levels specified herein require increased monitoring, corrective actions to abate emissions, and/or work shutdown. Additionally, the CAMP helps to confirm that work activities did not spread contamination off-site through the air.

The generic CAMP presented below will be sufficient to cover many, if not most, sites. Specific requirements should be reviewed for each situation in consultation with NYSDOH to ensure proper applicability. In some cases, a separate site-specific CAMP or supplement may be required. Depending upon the nature of contamination, chemical- specific monitoring with appropriately-sensitive methods may be required. Depending upon the proximity of potentially exposed individuals, more stringent monitoring or response levels than those presented below may be required. Special requirements will be necessary for work within 20 feet of potentially exposed individuals or structures and for indoor work with co-located residences or facilities. These requirements should be determined in consultation with NYSDOH.

Reliance on the CAMP should not preclude simple, common-sense measures to keep VOCs, dust, and odors at a minimum around the work areas.

Community Air Monitoring Plan

Depending upon the nature of known or potential contaminants at each site, real-time air monitoring for VOCs and/or particulate levels at the perimeter of the exclusion zone or work area will be necessary. Most sites will involve VOC and particulate monitoring; sites known to be contaminated with heavy metals alone may only require particulate monitoring. If radiological contamination is a concern, additional monitoring requirements may be necessary per consultation with appropriate DEC/NYSDOH staff.

Continuous monitoring will be required for all <u>ground intrusive</u> activities and during the demolition of contaminated or potentially contaminated structures. Ground intrusive activities include, but are not limited to, soil/waste excavation and handling, test pitting or trenching, and the installation of soil borings or monitoring wells.

Periodic monitoring for VOCs will be required during <u>non-intrusive</u> activities such as the collection of soil and sediment samples or the collection of groundwater samples from existing monitoring wells. "Periodic" monitoring during sample collection might reasonably consist of taking a reading upon arrival at a sample location, monitoring while opening a well cap or

Final DER-10 Page 204 of 226

overturning soil, monitoring during well baling/purging, and taking a reading prior to leaving a sample location. In some instances, depending upon the proximity of potentially exposed individuals, continuous monitoring may be required during sampling activities. Examples of such situations include groundwater sampling at wells on the curb of a busy urban street, in the midst of a public park, or adjacent to a school or residence.

VOC Monitoring, Response Levels, and Actions

Volatile organic compounds (VOCs) must be monitored at the downwind perimeter of the immediate work area (i.e., the exclusion zone) on a continuous basis or as otherwise specified. Upwind concentrations should be measured at the start of each workday and periodically thereafter to establish background conditions, particularly if wind direction changes. The monitoring work should be performed using equipment appropriate to measure the types of contaminants known or suspected to be present. The equipment should be calibrated at least daily for the contaminant(s) of concern or for an appropriate surrogate. The equipment should be capable of calculating 15-minute running average concentrations, which will be compared to the levels specified below.

- 1. If the ambient air concentration of total organic vapors at the downwind perimeter of the work area or exclusion zone exceeds 5 parts per million (ppm) above background for the 15-minute average, work activities must be temporarily halted and monitoring continued. If the total organic vapor level readily decreases (per instantaneous readings) below 5 ppm over background, work activities can resume with continued monitoring.
- 2. If total organic vapor levels at the downwind perimeter of the work area or exclusion zone persist at levels in excess of 5 ppm over background but less than 25 ppm, work activities must be halted, the source of vapors identified, corrective actions taken to abate emissions, and monitoring continued. After these steps, work activities can resume provided that the total organic vapor level 200 feet downwind of the exclusion zone or half the distance to the nearest potential receptor or residential/commercial structure, whichever is less but in no case less than 20 feet, is below 5 ppm over background for the 15-minute average.
- 3. If the organic vapor level is above 25 ppm at the perimeter of the work area, activities must be shutdown.
- 4. All 15-minute readings must be recorded and be available for State (DEC and NYSDOH) personnel to review. Instantaneous readings, if any, used for decision purposes should also be recorded.

Particulate Monitoring, Response Levels, and Actions

Particulate concentrations should be monitored continuously at the upwind and downwind perimeters of the exclusion zone at temporary particulate monitoring stations. The particulate monitoring should be performed using real-time monitoring equipment capable of measuring particulate matter less than 10 micrometers in size (PM-10) and capable of integrating over a period of 15 minutes (or less) for comparison to the airborne particulate action level. The equipment must be equipped with an audible alarm to indicate exceedance of the action level. In addition, fugitive dust migration should be visually assessed during all work activities.

Final DER-10 Page 205 of 226

- 1. If the downwind PM-10 particulate level is 100 micrograms per cubic meter (mcg/m³) greater than background (upwind perimeter) for the 15-minute period or if airborne dust is observed leaving the work area, then dust suppression techniques must be employed. Work may continue with dust suppression techniques provided that downwind PM-10 particulate levels do not exceed 150 mcg/m³ above the upwind level and provided that no visible dust is migrating from the work area.
- 2. If, after implementation of dust suppression techniques, downwind PM-10 particulate levels are greater than 150 mcg/m³ above the upwind level, work must be stopped and a re-evaluation of activities initiated. Work can resume provided that dust suppression measures and other controls are successful in reducing the downwind PM-10 particulate concentration to within 150 mcg/m³ of the upwind level and in preventing visible dust migration.
- 3. All readings must be recorded and be available for State (DEC and NYSDOH) and County Health personnel to review.

December 2009

Final DER-10 Page 206 of 226

Attachment B

NYSDEC DER-10 Appendix 1B Fugitive Dust and Particulate Monitoring

Appendix 1B **Fugitive Dust and Particulate Monitoring**

A program for suppressing fugitive dust and particulate matter monitoring at hazardous waste sites is a responsibility on the remedial party performing the work. These procedures must be incorporated into appropriate intrusive work plans. The following fugitive dust suppression and particulate monitoring program should be employed at sites during construction and other intrusive activities which warrant its use:

- Reasonable fugitive dust suppression techniques must be employed during all site activities which may generate fugitive dust.
- Particulate monitoring must be employed during the handling of waste or contaminated soil or when activities on site may generate fugitive dust from exposed waste or contaminated soil. Remedial activities may also include the excavation, grading, or placement of clean fill. These control measures should not be considered necessary for these activities.
- Particulate monitoring must be performed using real-time particulate monitors and shall monitor particulate matter less than ten microns (PM10) with the following minimum performance standards:
 - (a) Objects to be measured: Dust, mists or aerosols;
 - (b) Measurement Ranges: 0.001 to 400 mg/m3 (1 to 400,000 :ug/m3);
- (c) Precision (2-sigma) at constant temperature: +/- 10 :g/m3 for one second averaging; and +/- 1.5 g/m3 for sixty second averaging;
 - (d) Accuracy: +/- 5% of reading +/- precision (Referred to gravimetric calibration with SAE fine test dust (mmd= 2 to 3 :m, g= 2.5, as aerosolized);
 - (e) Resolution: 0.1% of reading or 1g/m3, whichever is larger;
 - (f) Particle Size Range of Maximum Response: 0.1-10;
 - (g) Total Number of Data Points in Memory: 10,000;
- (h) Logged Data: Each data point with average concentration, time/date and data point number
- (i) Run Summary: overall average, maximum concentrations, time/date of maximum, total number of logged points, start time/date, total elapsed time (run duration), STEL concentration and time/date occurrence, averaging (logging) period, calibration factor, and tag number;
- Alarm Averaging Time (user selectable): real-time (1-60 seconds) or STEL (15 minutes), alarms required;
 - (k) Operating Time: 48 hours (fully charged NiCd battery); continuously with charger;
 - (l) Operating Temperature: -10 to 50° C (14 to 122° F);
- (m) Particulate levels will be monitored upwind and immediately downwind at the working site and integrated over a period not to exceed 15 minutes.
- In order to ensure the validity of the fugitive dust measurements performed, there must be 4. appropriate Quality Assurance/Quality Control (QA/QC). It is the responsibility of the remedial party to adequately supplement QA/QC Plans to include the following critical features: periodic instrument calibration, operator training, daily instrument performance (span) checks, and a record keeping plan.
 - The action level will be established at 150 ug/m3 (15 minutes average). While conservative, 5.

this short-term interval will provide a real-time assessment of on-site air quality to assure both health and safety. If particulate levels are detected in excess of 150 ug/m3, the upwind background level must be confirmed immediately. If the working site particulate measurement is greater than 100 ug/m3 above the background level, additional dust suppression techniques must be implemented to reduce the generation of fugitive dust and corrective action taken to protect site personnel and reduce the potential for contaminant migration. Corrective measures may include increasing the level of personal protection for on-site personnel and implementing additional dust suppression techniques (see paragraph 7). Should the action level of 150 ug/m3 continue to be exceeded work must stop and DER must be notified as provided in the site design or remedial work plan. The notification shall include a description of the control measures implemented to prevent further exceedances.

- 6. It must be recognized that the generation of dust from waste or contaminated soil that migrates off-site, has the potential for transporting contaminants off-site. There may be situations when dust is being generated and leaving the site and the monitoring equipment does not measure PM10 at or above the action level. Since this situation has the potential to allow for the migration of contaminants off-site, it is unacceptable. While it is not practical to quantify total suspended particulates on a real-time basis, it is appropriate to rely on visual observation. If dust is observed leaving the working site, additional dust suppression techniques must be employed. Activities that have a high dusting potentialsuch as solidification and treatment involving materials like kiln dust and lime--will require the need for special measures to be considered.
- The following techniques have been shown to be effective for the controlling of the generation and migration of dust during construction activities:
 - (a) Applying water on haul roads:
 - (b) Wetting equipment and excavation faces;
 - (c) Spraying water on buckets during excavation and dumping;
 - (d) Hauling materials in properly tarped or watertight containers;
 - (e) Restricting vehicle speeds to 10 mph;
 - (f) Covering excavated areas and material after excavation activity ceases; and
 - (g) Reducing the excavation size and/or number of excavations.

Experience has shown that the chance of exceeding the 150ug/m3 action level is remote when the above-mentioned techniques are used. When techniques involving water application are used, care must be taken not to use excess water, which can result in unacceptably wet conditions. Using atomizing sprays will prevent overly wet conditions, conserve water, and provide an effective means of suppressing the fugitive dust.

The evaluation of weather conditions is necessary for proper fugitive dust control. When extreme wind conditions make dust control ineffective, as a last resort remedial actions may need to be suspended. There may be situations that require fugitive dust suppression and particulate monitoring requirements with action levels more stringent than those provided above. Under some circumstances, the contaminant concentration and/or toxicity may require additional monitoring to protect site personnel and the public. Additional integrated sampling and chemical analysis of the dust may also be in order. This must be evaluated when a health and safety plan is developed and when appropriate suppression and monitoring requirements are established for protection of health and the environment.

Final DER-10 Page 208 of 226 May 2010

Appendix 1C DEC Permits Subject to Exemption

In accordance with section 1.10, exemptions from the following permit programs may be granted to the person responsible for conducting the remedial programs undertaken pursuant to section 1.2:

Air - Title 5 permits

Air - State permits

Air - Registrations

Ballast Discharge

Chemical Control

Coastal Erosion Hazard Areas

Construction of Hazardous Waste Management Facilities

Construction of Solid Waste Management Facilities

Dams

Excavation and Fill in Navigatable Waters (Article 15)

Flood Hazard Area Development

Freshwater Wetland

Hazardous Waste

Long Island Wells

Mined Land Reclamation

Navigation Law - Docks

Navigation Law - Floating Objects

Navigation Law - Marinas

Non-Industrial Waste Transport

Operation of Solid Waste Management Facilities

Operation of Hazardous Waste Management Facilities

State Pollution Discharge Elimination Systems (SPDES)

Stream Disturbance

Tidal Wetlands

Water Quality Certification

Water Supply

Wild, Scenic and Recreational Rivers